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29 April 2022

OKPD2 32.50.13.190

Operating manual for the medical device OTSG Xcavator



We, **Ovesco Endoscopy AG**, at the registered address of Friedrich-Miescher-Strasse 9, 72076 Tübingen, Germany, as the legal manufacturer, hereby declare that the grasping forceps "OTSG Xcavator" with accessories meets the international standard and Russian requirements.

The product was produced according to Directive 93/42/EEC from 14.06.1993.

The attached 15 pages represent the Operating manual of the OTSG Xcavator.

Issued on behalf of Ovesco Endoscopy AG.

29 April 2022

Dipl.-Ing. (FH) Christoph Damm
Head of Quality Management
Regulatory Affairs
Friedrich-Miescher-Straße 9
72076 Tübingen
www.ovesco.com

Place of business:
Tübingen, Germany

Company Register:
Stuttgart HRB 727 461

VAT identification number:
DE 217051960

Executive Board
Prof. Dr. med. Marc O. Schurr
Prof. Dr. med. Thomas Gottwald

Chairman Supervisory Board:
Frank Ramsperger

Bank information:
Baden-Württembergische Bank
SWIFT-BIC: SOLA DEST600
Account 452 030 4
IBAN
DE62 6005 0101 0004 5203 04

Kreissparkasse Tübingen
SWIFT-BIC: SOLADES1TUB
Konto 109949
IBAN
DE43 6415 0020 0000 1099 49

OKPD2 32.50.13.190

OPERATING MANUAL

of the Medical Device

Set for foreign body retrieve OTSG Xcavator

Manufacturer: Ovesco Endoscopy AG

Meets the international standard and Russian requirements of GOST 50444-92

Manufactured according to Directive 93/42/EEC from 14.06.1993.

April 29, 2022

Ovesco Endoscopy AG

Friedrich-Miescher-Strasse 9
72076 Tübingen
Germany

Table of contents

1. Medical device design
2. Technical requirements
3. Area of application
4. Operation information
5. Operating procedure
6. Maintenance information
7. Repair information
8. Labeling
9. Safety information
10. Disposal information
11. Expiration date and period of storage
12. Transportation and storage
13. Warranties
14. List of international certificates and a declaration of conformity
15. List of the national and Russian standards
16. Engineering Design Process Steps

1. MEDICAL DEVICE DESIGN

Set for foreign body retrieve OTSG Xcavator (OTSG Xcavator – further) (REF No. 200.15) is an endoscopic attachment cap with grasping function for flexible endoscopy. The OTSG Xcavator is designed for grasping large volumes of tissue or foreign bodies. The OTSG Xcavator allows for the use of additional endoscopic instruments through the working channel of the endoscope.

The cap [A] is mounted on the distal part of the endoscope.

A flexible tube [B] extends from the forceps handle [C] through the biopsy channel of the endoscope. An endoscope sleeve [D] is included to avoid trapping of surrounding tissue between the endoscope and the shaft.

The OTSG Xcavator consists of:


Fig 1. Supply kit

The OTSG XCavator is a sterile single use medical device. OTSG XCavator does not contain medicinal products for human use and pharmaceutical substances or materials of animal or human origin.

Table 1. Specification of the device

Components according to figure 1		Indication area	Endoscope diameter [mm]	Outer diameter of cap [mm]*	Endoscope working length compatibility [mm]
A	Application cap	Gastro-intestinal tract	9.5 – 10.5	14.3 – 15.5	1650
B	Shaft				
C	Handle				
D	Endoscope sleeve				

*The maximum outer diameter of the mounted instrument depends on the actual outer diameter of the endoscope. The OTSG Xcavator cap increases the overall outer diameter by approx. 4 – 5 mm to a total of approx. 14.3 – 15.5 mm.

Supply kit:

Foreign body retrieval and biopsy forceps in the OTSG Xcavator endoscopic set, consisting of:

1. Foreign body retrieval and biopsy forceps in the OTSG Xcavator endoscopic set, consisting of:

- Foreign body retrieval and biopsy forceps OTSG Xcavator - 1 piece;

- Distal cap - 1 pc.

2. Protective tube - 1 piece;

3. Instruction for users.


2. TECHNICAL REQUIREMENTS

Requirements:

1. Total weight of the device is not more than 65 g. Weight of OTSG endoscopic foreign body retrieval and biopsy forceps is not more than 45 g., the weight of the protective tube is not more than 10 g., the weight of the distal cap is not more than 10 g.
2. The device is a sterile. An ethylene oxide sterilized medical device.
3. The device is a single use.
4. The device is sterile, non-toxic.
5. Materials:

Table 2. Materials

№	Components	System	Material / Trademark	Standard
1	Foreign body retrieve OTSG	OTSG Xcavator	Stainless steel 1.4301 ABS plastic Fluoroplastic PTFE	EN ISO 7153 AISI 303 USP class VI
2	Protective tube	OTSG Xcavator	Low Density Polyethylene (LDPE)	USP class VI FDA 177.2600

6. Tool dimensions:

Fig 2. OTSG Xcavator dimensions

- a** Total length: $18 \pm 0,1$ mm
- b** Cap depth: $6,5 \pm 0,1$ mm
- c** Jaws length: $11,5 \pm 0,1$ mm.
- d** Opening diameter: $1,1 \pm 0,05$ mm. The distance between the edges of the branches in the closed state is $1,1 \pm 0,05$ mm.
- e** Max. outer diameter: The maximum outer diameter of the mounted instrument depends on the actual outer diameter of the endoscope. The OTSG Xcavator cap increases the overall outer diameter by approx. 4–5 mm to a total of approx. 14.3–15.5 mm. Diameter of the distal part of the endoscope: 9.5 - 10.5 mm.
- f** The diameter of the endoscope channel not less than 2,8 mm. The figure shows a distal cap mounted on the distal part of the endoscope and an enlarged working part of the forceps removed from the endoscope channel with the already mounted

cap.

g Opening angle: $155 \pm 1^\circ$

h Max. diameter with fully open jaws: $28,4 \pm 0,2$ mm.

The maximum volume that can be grasped may be restricted by the conditions at the application site and patient anatomy.

7. When a force of 100 N is applied to the handle, the branches of the forceps ensure the biting of biopsy material (fragments of tissue of the upper gastrointestinal tract) with a force no more than 10 N.
8. The forceps must move in the instrumental channel of the endoscope (maximum external diameter of 10.5 mm) when applying an effort to guide the product through the instrumental channel of no more than 20 N.
9. The bendability is restricted by less than or equal to 25% of maximum angulation.
10. The force applied between the branches must be no more than 8 N.
11. Overall dimensions of the package that is subject to final sterilization (in mm) for OTSG Xcavator: $220 \pm 10 \times 360 \pm 10$ mm.
12. The material of the sterile package: PET (Polyethylenterephthalat) / PE (Polyethylen) peel, high strength paper (Tyvek) (DuPont, Luxembourg).
13. The seam strength of the sterile package is defined as ≥ 1.6 N.
14. The surface weight of the sterilization materials is 70 g/m^2 ($-5/+3$) for PET / PE and 75 g/m^2 for Tyvek.
15. Connecting dimensions of the protective tube: the tube is mounted on top of the endoscope with a diameter of the insertion part of 9.8 mm. Fastening is carried out due to the shape of the distal and proximal parts of the tube.
16. The branches of the forceps have cutting edges with a width of 1.1 ± 0.05 mm.
17. The diameter of the distal cap is 16.5 ± 1 mm, the length of the distal cap is 18 ± 1 mm.

3. AREA OF APPLICATION

The product is used exclusively in stationary conditions of hospital departments and clinics in endoscopic or surgical departments. The endoscopist using the product must be experienced in the technology of its application.

4. OPERATION INFORMATION

Do not use if packaging has been damaged.

Do not use after the expiration time.

If the packaging has been damaged or the use-by date has expired, the manufacturer does not guarantee the sterility of the device.

Conditions of using medical device:

- Temperature 10-40 °C
- Humidity 40 – 100 %

During the operation, all the tools must be resistant to the effects of climatic factors for the execution of UHL 4.2 GOST 15150, and their working parts, as products of the climatic version U6, to the effects of temperatures in the range from 32 to 42 ° C according to GOST R 50444.

Before using the device, the instructions for use must be read carefully. The user must have fully understood how to assemble, operate and use the device. Always use the forceps in strict accordance with these instructions for use.

Depends on the potential risk the device is classified as 2 a (according to the order №4 from 06.06.2012 of the Russian Ministry of Health).

OKPD2 – 32.50.13.110

The type of the device according to Russian medical device classification (according to the order №4 dated 06.06.2012 of the Russian Ministry of Health) – 179080.

INDICATION FOR USE

Set for foreign body retrieve OTSG Xcavator for flexible endoscopy is designed to grasp large volumes of tissue and foreign bodies in the gastrointestinal tract.

CONTRADICTIONS

Contraindications to the use of the system are associated with the patient's contraindications to the use of flexible endoscopy treatment.

The product can only be used with flexible endoscopes.

POSSIBLE COMPLICATIONS

- Damage or perforation of the tissue of the gastrointestinal tract
- Bleeding
- Inflammation or infection of the tissue
- Dislocation of the stent when passing through a stent

5. OPERATING PROCEDURE

5.1. Before use

Before using the device, inspect the sterile package carefully. Products with broken packing are forbidden to be used.

The physician must be sure that he is fully familiar with the instructions for use of the forceps before using the device.

Before using the device, it is always necessary to check compatibility of the distal end of the endoscope with this medical device. The forceps are used with endoscopes with a diameter of the distal part of 9.5 – 10.5 mm.

Before inserting the instrument into the working area, open and close the branches of the forceps to check their operability.

When the device is applied to patients having a transgastric access (stent) was placed before, it is always necessary to ensure first, that the diameter of the mounted system allows passing through the stent.

Due to the size of the cap, the visibility of the endoscopic image, carrying out and working with the instrument may be difficult, thereby increasing the risk of involuntary perforation.

5.2. Mounting the OTSG Xcavator cap

- The tube with the additional endoscopic channel is already mounted on the flexible part of the OTSG Xcavator system. Slide the tube with the endoscope sleeve over the endoscope tip and onto the endoscope shaft. First place the OTSG cap only lightly on the endoscope tip. Align the cap so that the working channel and endoscopic vision are not obscured.
- Thereafter, push the cap with slight lateral tilting movements until it reaches stoppers on the cap, whereupon it will be reliably fixed. Figure 3-A shows correct location of the cap; Figure 3-B shows the endoscopic perspective.


Fig 3. Mounting the device cap

When mounting the cap, it is necessary to make sure that the cap is mounted on the endoscope correctly. Visualization and functions of the endoscope must be fully accessible to the endoscopist and attending specialist.

5.3. Mounting the protective tube

- Remove the retainer and pull the protective tube from the distal part to the handle of the endoscope.
- The endoscope must be clean and dry. Fix the endoscope at a distance of 2-3 mm of the entire length of the endoscope. At the same time, do not fix the cap. Make sure that the tube of the forceps runs straight and is fixed with the endoscope.
- Pull the protective tube at the proximal end back using the tube. Make sure that the tube is stretched to the handle of the endoscope.
- Remove the remaining protective clips.

5.3. Operating the device

- Before carrying out the capture, use a guide (not included, maximum diameter 1 mm), for a more confident passage of the tool.
- A large volume of tissue or a foreign body can block the closure of the tool branches. If this happens, use an irrigator, aspirator or grab to release the branches.
- When using the instrument for manipulation in the working channel of the endoscope, it is necessary to ensure that the working part of the instrument is completely immersed in the distal cap. Otherwise, the tool may be clamped and not close completely.

When inserting, operating, and removing the forceps, the branches of the OTSG Xcavator must always be closed, and constant endoscopic visualization is recommended.

1. Inserting the device

- Guide the endoscope with mounted OTSG Xcavator over the guide wire to the site of interest.
- Target the tissue or foreign body to be grasped.

2. Grasping and removal of tissue or foreign bodies

When opening the grasper branches, make sure that there is enough space so that they do not damage the surrounding tissue. Carefully open the branches under endoscopic control.

- Push the slider of the handle forward (Fig. 2 A) to open the branches.
- Position the target tissue or foreign body between the branches.
- Pull back the slider of the handle (Fig. 2 B) to close the grasper branches. The target tissue/foreign body is clamped between the branches.
- Carefully remove the endoscope from the body.


Fig 6. Operating of the handle to open and close the branches

5.4. Dismounting of the device

- Remove all adhesive tapes.
- Carefully remove the cap together with the endoscope sleeve from the endoscope.
- Dispose all components according the hygiene regulations applicable to the user.

6. MAINTENANCE INFORMATION

The product warranty period is 3 years from the day of manufacture.

It is guaranteed that the product has no defects connected with production and materials quality. The product would be replaced if there was no physical/mechanical intervention to the product, inconsistent with the requirements of the operation manual. Make sure that any malperformance appears not due to neglect or noncompliance with instructions. The storage period is 3 years. In case of malfunction under compliance with the requirements of the operation manual, the manufacturer shall replace the product.

The device cannot be re-sterilized. Repeated use of the product may harm the patient!

For all questions, please contact official distributor in Russian Federation: Endo Stars LLC, legal address: 197101, Russia, Saint-Petersburg, Bolshaya Monetnaya, 27, lit. A, 12H ; tel./fax: +7-800-555-56-57, e-mail endo@endo-stars.com.

7. REPAIR INFORMATION

Single use. Unrepairable.

8. LABELING

Marked:

- The description of the device;
- Dimensions;
- Diameter ;

The individual packaging should contain the following symbols:

- CE mark and identification number of the notified body
- See instruction for use;
- LOT number;
- REF number;
- Manufacturer;
- Date of manufacture;

- Diameter of the compatible endoscope channel;
- Warranty period of shelf life (storage);
- EO sterilized;
- Do not re-use;
- Do not re-sterilise;
- Do not use if packaging damaged;
- Latex-free;
- Items per package.

Transport labeling should be marked on the pouch, carton boxes or on the package. Labeling is stuck or connected in other reliable ways that guarantee the safety of the goods and labeling.

9. SAFETY INFORMATION

Do not use if packaging has been damaged.

Do not use after the expiration time.

The product was produced according to Directive 93/42/EEC from 14.06.1993.

Single use: yes.

10. DISPOSAL INFORMATION

According to the Sanitary and Epidemiological Requirements for Disposal of Medical Waste: Sanitary and Epidemiological Regulations and Norms 2.1.3684-21. Class B

11. EXPIRATION DATE AND PERIOD OF STORAGE

The period of preservation of sterility is 3 years from the date of sterilization.

The guaranteed shelf life (storage) is 3 years from the date of sterilization.

12. TRANSPORTATION AND STORAGE

The product is packed in the appropriate packaging together with the operational documentation. Transportation by all modes of transport in covered vehicles and storage of the product is allowed under the following conditions:

- Temperature -50°C - +50°C
- relative humidity 45-80%,
- atmospheric pressure atmospheric pressure ranges from 86.0 kPa to 106.7 kPa (650-800 mmHg)

After transportation, please pay your special attention to possible damages.

The tools in transport packaging should be sustainable in the process of transportation under GOST R 50444.

The guaranteed shelf life (storage) is 3 years from the date of sterilization.

13. WARRANTIES

a. The manufacturer grants conformity of the product to this technical documentation subject to the condition of correct using, transporting, and storing.

b. All question should be sent to an official distributor in Russian Federation:

Endo Stars LLC, 197101, Saint-Petersburg, Bolshaya Monetnaya, 27, Lit. A, 12H,
Phone/fax 8-800-555-56-57, email: endo@endo-stars.com

14. LIST OF INTERNATIONAL CERTIFICATES AND A DECLARATION OF CONFORMITY


- Certificate ISO 13485: 2016 № 50970-14-01 from 20.05.2019
- EC certificate № 50970-16-05 from 20.05.2019
- Certificate of incorporation from 20.01.2016

15. LIST OF THE NATIONAL AND RUSSIAN STANDARDS

The list of applicable national standards of the Russian Federation and international standards to ensure safety, efficiency and quality, as well as standards with control methods: Directive 93/42/EEC, EN ISO 13485:2016.

Manufactured in accordance with Directive 93/42/EEC of 08.11.2015.

16. ENGINEERING AND MANUFACTURING DESIGN PROCESS STEPS

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