7D

ADDITION TO USER MANUAL FOR A MEDICAL DEVICE

“Electromyograph handheld Dantec CLAVIS Offline power and control”


Sean David Langan

Vice President, Global Operations

Natus Manufacturing Limited

SIGNED IN MY PRESENCE BY SEAN DAVID LANGAN WHO IS
PERSONALLY KNOWN TO ME THIS 10TH OCTOBER 2002
I BEING A NOTARY PUBLIC.

**OLIVIA LYNCH
CROW STREET
GORT, COUNTY GALWAY
NOTARY PUBLIC
COMMISSIONED FOR LIFE**

1. Purpose of the medical device

1.1 Name of the medical device

Electromyograph handheld Dantec CLAVIS self-powered and controlled

1.1 Purpose of the medical device

The DantecCLAVIS is a medical device designed to be used as a stimulator for nerve localization as well as a navigator for targeted needle insertion in intramuscular injections.

1.2 Indications

1. To control the administration of botulinum toxin (protein) and other drugs into the muscles.
2. Nerve block control.
3. Definition of muscle contracture.
4. Diagnosis of muscles involved in pathological processes (cervical, focal dystonia, spasticity).
5. Mapping of motor zones of muscles.
6. Intraoperative stimulation of nerves for the prevention of their surgical damage

1.3 Contraindications

Patients with epileptic status in the anamnesis. For patients with blood clotting disorders, prescribe procedures using the device with caution. Pregnancy and lactation of the breast are also a contraindication. There are also general physiotherapeutic contraindications.

1.4 Possible side effects when using a medical device

1. Pain and redness together injections
2. Painful muscle spasm as a result of electrical stimulation

1.5 Indication of the possibility and features of the use of a medical product for humans by medical products implanted in the human body

forbidden to perform procedures using the device for patients with implanted medical devices.

1.6 Operating conditions

Table - Operating conditions

Working environment conditions	Environment temperature	10°C to 40°C
	Relative humidity	30% to 75% (non-condensing)
	atmospheric pressure	From 700gPa to 1060gPa

1.7 Personnel qualification

This medical device is intended for use by qualified medical personnel with basic knowledge of electroneuromyography.

1.8 Target group

This medical device is used for patients in accordance with the indications for use

1.9 Conditions for the use of a medical device in medical organizations The handheld electromyograph Dantec CLAVIS with autonomous power supply and control is intended for use by qualified medical personnel who have knowledge of the basics of electroneuromyography.

2. Information about the manufacturer of the medical device

2.1 Name Natus Manufacturing Limited,
Ireland (Natus Manufacturing Limited, Ireland)

2.2 location
IDA Business Park Gort Co. Galway, Ireland

2.3 Address of the medical device manufacturing site
Natus Manufacturing Limited, Ireland
(Natus Manufacturing Limited, Ireland) IDA Business Park
Gort Co. Galway, Ireland

3. Classification of a medical device

Table - Classification of a medical device

nomenclature medical devices	classification	Class 2b
OKPD2 code		26.60.12.121
Classification electric shock	by type of protection	Medical device with internal power source
According to electric current	the degree of protection	bf
By type of software safety according to IEC62304		class B
Product class by medical non-electrical waste in accordance with SanPiN2.1.7.2790-10		ClassA
According to the degree of protection against the ingress of foreign solid objects and the ingress of water (based on the standard IEC60529)		IPX0
Operating mode		continuous

4. Description of the principles on which the operation of the medical device is based and their features

Electromyograph handheld Dantec CLAVIS with autonomous power supply and control works on the basis of a method for studying bioelectric potentials that arise in the skeletal muscles of humans and animals during excitation of muscle fibers. In EMG mode, this device is used for needle electrode examinations. During the procedure, the Dantec CLAVIS will emit a series of beeps that differ in intensity and frequency, making it easier to locate the required muscle or nerve. Flashing indicators

on the stimulation level panel allow you to control the level of current transmitted to the patient.

5. General photo of MI


Figure - General photo of MI (handheld electromyograph Dantec CLAVIS with autonomous power supply and control)



5.1 Variants of execution of a medical device

This MI is made in a water version.

5.2 Basic compose the list of components and accessories of MI

Maximum permissible deviation no more than $\pm 10\%$

Handheld electromyograph control and accessories		Dantec CLAVIS self-powered	
Name Components		Description Components	photography components
on e	Electromyograph handheld Dantec CLAVIS with self-powered and controlled	DantecCLAVIS- a medical device intended to be used as a stimulator for localizing nerves, as well as a navigator for targeted needle insertion in intramuscular injections.	
2	Management user paper on User manual electronic media on	User's manual contains the material the study of which is necessary to ensure the correct operation of the electromyograph handheld Dantec CLAVIS self -powered and controlled	 Руководство пользователя 

3	Needle electrode monopolar injection (for RU No. FSZ 2011/10502)	Disposable needle electrodes with a hypodermic needle are used for simultaneous drug administration and recording	
four	Oval adhesive electrode (according to RU No. FSZ 2011/10502)	Applies registration to a and monitoring of EMG	
5	Electrogrounding lamellar (poru No. FSZ 2011/10502)	Applies registration to a and monitoring of EMG	
6	Keyportable	Designed for transporting an electromyograph handheld Dantec CLAVIS with self-powered control.	

7	Battery9V	Designed to power the electromyograph DantecCLAVIS	
---	-----------	--	---

5.3 Possibility and methods of integration with other medical devices.

Natus Manufacturing Limited recommends the use of the following accessories with the Dantec CLAVIS Self-Powered Handheld Electromyograph:

Name	DescriptionDetails
Needle electrode (EMG input, cathode) 	Electrode needle monopolar subcutaneous Bo-ject single use with output, size: 25 x 0.30 mm, Natus Manufacturing Limited , Ireland (1pack - 12 pcs.).
	Electrode needle monopolar subcutaneous Bo-ject single use with output, size: 25 x 0.41 mm, Natus Manufacturing Limited , Ireland, (1pack - 12 pcs.).
	Electrode needle monopolar subcutaneous Bo-ject single use with output, size: 37 x 0.41 mm, Natus Manufacturing Limited , Ireland, (1pack - 12 pcs.).
	Electrode needle monopolar subcutaneous Bo-ject single use with output, size: 50 x 0.41 mm, Natus Manufacturing Limited , Ireland, (1pack - 12 pcs.).
	Electrode needle monopolar subcutaneous Bo-ject single use with output, size: 75 x 0.64 mm, Natus Manufacturing Limited , Ireland, (1pack - 12 pcs.).
	Electrode needle monopolar subcutaneous Bo-ject single use with output, size: 37 x 0.46 mm, Natus Manufacturing Limited , Ireland (1pack.- 12 pcs.).
	Electric needle monopolar injection (size: 25x0.30 mm)(Bo - ject), NatusNeurologyIncorporated , USA, (1pack-10pcs).
	Electric needle monopolar injection (size: 25x0.41 mm)(Bo - ject), NatusNeurologyIncorporated , USA, (1pack-10pcs).
	Electric needle monopolar injection (size: 37x0.41 mm)(Bo - ject), NatusNeurologyIncorporated , USA, (1pack-10pcs).
	Electric needle monopolar injection (size: 50x0.51 mm)(Bo - ject), NatusNeurologyIncorporated , USA, (1pack-10pcs).
	Electric needle monopolar injection (size: 75x0.71 mm)(Bo - ject), NatusNeurologyIncorporated , USA, (1pack-10pcs).
	Electric needle monopolar injection (size: 37x0.46 mm)(Bo - ject), NatusNeurologyIncorporated , USA, (1pack-10pcs).

Skin electrodes (anode, reference signal, grounding for the patient) 	Disposable skin electrode with gel, 50cm cable 1.5 mm TPC connector.(12 pcs.).
	Disposable skin electrode with gel, 10 cm cable TPC connector 0.7mm (12 pcs.).
	Disposable skin electrode with gel, 8cm cable TPC connector 0.7mm (12 pcs.).
	80cm unshielded cable with 0.7mm and 1.5mm TPC connectors (8pcs, colored).(Use with skin electrodes 13S021 and 13S024).
	Film electrodes adhesive electroconductive gel surface basis (100 pcs.).
	Unshielded cable with clamp connector, green, 80cm, female connector type TPC (pcs.). (Use with electrodes 13S023).
	Unshielded cable with clamp connector, black, 80cm, female connector type TPC (pcs.). (Use with electrodes 13S023).
	Unshielded cable with clamp connector, red, 80cm, female connector type TPC (pcs.). (Use with electrodes 13S023).
Ground electrodes (Patient ground) 	Grounding electrode, 10 cm Velcro, 1.5 m cable with connector TRS type 1.5 mm.
	Earth electrode, 18.5 cm Velcro, 1.5 m cable with connector TRS type 1.5 mm.
	Ground electrode, 75 cm Velcro, 1.5 m cable with connector TRS type 1.5 mm.
	Ground electrode, 45 cm Velcro, 1.5 m cable with connector TRS type 1.5 mm.
	Ground plate electrode, 1.5 m cable with connector TRS type 1.5 mm.
	Ground plate electrode, 1.2m cable with connector TRS type 1.5 mm.

5.4 Materials from which the medical device or its main parts are made, biocompatibility


Product name or parts and accessories (tested)	Materials used in the manufacture	Type of contact with the body (according to GOST ISO 10993-1-2011)
Electromyograph handheld Dantec CLAVIS with autonomous power supply and control	<i>Frame:</i> ABS plastic 1106-M-30 <i>Dye:</i> White brand SWS13-POM Weiß batch 50,	Short-term contact (less than 24 hours) with intact skin and mediated contact with the body.

6. Information about markings and symbols of the system

6.1 Labeling information

Warning signs on the body of the product		
1		Device type BF
2		Carefully!
3		This device does not contain parts that require user service.
4		Disposal of Electrical and Electronic Equipment: Compliance Information.
5		Serial number.
6		Catalog number.
7		Battery Type.
8		Caution: United States Federal Law restricts this device to sale by or on the order of a physician.
9		This device complies with EU directive 93/42/EEC on medical devices.
10		Manufacturer
11		date of manufacture
12		Batch code
13		Use before
14		Refer to instructions for use

6.2 Photomarking

Dantec CLAVIS handheld electromyograph with autonomous power supply and control

Registration certificate: RZN № XXXX/XXXX date

Country of origin: Ireland

Date of manufacture and serial number are indicated on the main labeling of the product

Purpose; normative documentation, the requirements of which the medical device complies with; weight and size characteristics and service life of the product are specified in the user documentation

Authorized representative of the manufacturer:

Infomed-Neuro LLC, 125445, Moscow, Leningradskoe shosse, 69, bldg. 1, room 1, floor 1, room. fifty.

9B; alkaline battery

IEC-6LR61, ANSI-1604A.

Figure 2 - Marking of the Dantec CLAVIS handheld electromyograph with autonomous power supply and control


Figure - Marking of electrodes according to RU No. FSZ 2011/10502

7. Information about the packaging.

The box is made of corrugated cardboard 5 mm thick, painted white on the outside. Inside the cardboard box is placed a plastic bag of air-bubble wrap (bubble wrap) to protect the case from impact.

Package dimensions and weight of the set

Dimensions, L x W x H: 359 x 264 x 86 mm

Weight: 0.967 kg

7.1 Photos packing


Figure - Packing of handheld electromyograph Dantec CLAVIS with autonomous power supply and control

User Guide Supplement - Electromyograph handheld Dantec CLAVIS Offline power and control

7.2 Location of symbols on the package


Figure - Arrangement of symbols on the front side of the package


Figure - Arrangement of symbols on the package

8. Terms of transportation and storage

Storage and transportation conditions

ambient temperature	-10 C ⁰ to +50 C ⁰
Relative Humidity	from 10% to 100%.
Atmosphere pressure	from 700 hPa to 1060 hPa

Operating conditions

ambient temperature	+10 C ⁰ to +40 C ⁰
Relative Humidity	from 30% to 75%.
Atmosphere pressure	from 700 hPa to 1060 hPa

9. Main parameters and characteristics of the medical device:

9.1 Specifications for each component of the medical device

Maximum permissible deviation no more than +/-10 %

Weight and dimensions of MD:

User Guide Supplement - Electromyograph handheld Dantec CLAVIS Offline power and control

Handheld electromyograph Dantec CLAVIS with autonomous power supply and control	
Width, mm	twenty
Height, mm	80
Length, mm	140
Weight, kg	0.14
portable case	
Width, mm	80
Height, mm	227
Length, mm	279
Weight, kg	0.527
Battery9V	
Width, mm	16
Height, mm	25.5
Length, mm	48.5
Weight, kg	0.045

Main technical characteristicsMI

Electromyograph handheld Dantec CLAVIS self-powered and controlled	
Module EMG	
Gain	45000
EMG amplifier bandwidth, Hz	600 - 5000
Sampling frequency, Hz	30kHz
RMS noise	0.7 mV
active input	black1,5mm type TRS
Reference inputs	red1,5mm type TPC
Grounding for the patient	green1,5mm type TPC
Speaker bandwidth, Hz	340 - 10000
Stimulant	
Output current, mA	1.0—15.0
Maximum excitation voltage, V	100
Electrode resistance, Ohm	200–7000
Stimulation frequency, Hz	1 or 2
Current pulse amplitude, mA	1-15
Current amplitude control step, mA	one
Pulse duration, ms	0.1 or 0.2
output signals	rectangular shape
active output	black1,5mm type TRS
Reference Outputs	red1,5mm type TPC
portable case	
Maximum load	10kg.
Battery9V	
Power supply type	alkaline battery
Type of	IEC-6LR61,ANSI-1604A.
Voltage, V	9
Battery capacity	100..125mAh
Input current, A	2
Energy consumption	2Wmax
User manual on electronic media	
media type	CD-ROM
Volume occupied by the carrier	11.2 MB

User Guide Supplement - Electromyograph handheld Dantec CLAVIS Offline power and control

Version	A2
Electric needle monopolar subcutaneous disposable with withdrawal, size 75x 0.71 mm	
Needle length in the range	74-76 mm
Outer needle diameter in the range	0.64 -0.72mm
Output length in the range	78-90 cm
Connecting the output with a needle	soldering inside one-piece plastic holder
Connector for connection to electromyograph	TRS 1.5 mm
Possible types of syringes used	Luer-Lock
Sterility	sterile
Electric needle monopolar injection	
Specifications see according to RU No. FSZ2011/10502	
Electroadhesiveoval	
Specifications see according to RU No. FSZ2011/10502	
Electrogroundingplate	
Specifications see according to RU No. FSZ2011/10502	
Needle electrode monopolar subcutaneous disposable with a lead, size 25 x 0.41 mm.	
Needle length in the range	24-26 mm
Outer needle diameter in the range	0.40 -0.42 mm
Output length in the range	78-90 cm
Connecting the output with a needle	soldering inside one-piece plastic holder
Connector for connection to electromyograph	TRS 1.5 mm
Possible types of syringes used	Luer-Lock
Sterility	sterile
Needle electrode monopolar subcutaneous disposable with a lead, size 25 x 0.30 mm.	
Needle length in the range	24-26 mm
Outer needle diameter in the range	0.29 -0.31mm
Output length in the range	78-90 cm
Connecting the output with a needle	soldering inside one-piece plastic holder
Connector for connection to electromyograph	TRS 1.5 mm
Possible types of syringes used	Luer-Lock
Sterility	sterile

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10. Service life, shelf life

Service life - 7 years, shelf life - 10 years.

User Guide Supplement - Electromyograph handheld Dantec CLAVIS Offline power and control

11. environmental protection

11.1 Requirements for the protection of the environment when using a medical device

This product does not harm the environment with proper operation, transportation and storage. In order to maintain compliance with the declared period of harmless use, work on the equipment must be carried out in strict accordance with the instructions and environmental conditions given in the operating manual, as well as with strict observance of periodic maintenance schedules.

The period of environmentally friendly operation of consumables or individual parts may be less than the value set for the product. Periodic replacement of these consumables or parts in order to ensure environmentally friendly operation of the product during the claimed safety period must be carried out in accordance with the rules of maintenance.

12. Methods and conditions of sterilization, terms of preservation of sterility

The Dantec CLAVIS is not sterilizable.

13. List of national and international regulatory documents/standards to which the medical device complies

Standard	Description
IEC/EN60601-1	Medical electrical products - Part 1 - General requirements security, taking into account the main functional characteristics
EN60601-1-6	MEDICAL ELECTRICAL PRODUCTS. Part 1 - 6. GENERAL SAFETY REQUIREMENTS. Serviceability _
IEC60601-1-4	Medical electrical devices. Part 1. General safety requirements. 4. Safety requirements for programmable medical electronic systems
IEC/EN60601-1-2	MEDICAL ELECTRICAL PRODUCTS. Part 1-2 General safety requirements, taking into account the main functional characteristics. Parallel standard. Electromagnetic compatibility. Requirements and tests
IEC62366:2008	MEDICAL PRODUCTS.Design of medical products taking into account the serviceability
EN980	Medical devices. Graphic symbols used in labeling of medical devices
EN1041	Information supplied by the manufacturer for medical appliances
EN60601-2-40	MEDICAL ELECTRICAL EQUIPMENT . Part 2-40. Particular safety requirements for electromyographs and equipment
EN62304	Medical products.Software.Processes life cycle
ENISO13485	Medical devices.Quality management systems.Regulation requirements
ENISO14971	Medical devices. Application of risk management medical devices

14. Warranty

Natus Manufacturing Limited provides a 1-year warranty for a self-powered handheld electromyograph. Any warranties and obligations of Natus Manufacturing Limited will be void in the following cases:

- non-compliance with the requirements of this manual and other accompanying documentation when using the device.

15. Complaint

Authorized representative of the manufacturer of the medical device:

Organizational and legal form and complete name of the legal entity	Limited responsibility of "Infomed-Neuro"
Address (location) of the legal faces	125445, Moscow, Leningrad highway, d. 69, bldg. 1 room 1, floor 1, room 50.
Phone number	+7495 645 47 00
An identification number taxpayer	7743785545

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