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
This document is the Instruction for Use of Genoss Co., Ltd. The purpose of notarization is to submit the Representation Letter to the Ministry of Healthcare in Russia. Genoss Co., Ltd has the responsibility of any problem which may be caused by this letter of statement.

Date : 17 Jan 2021

Company Name : Genoss Co., Ltd

Sung Min Chung

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 1

Instruction for use

GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40

Manufacturer:

Genoss Co., Ltd.

1/F, Gyeonggi R&DB Center / 226, 2F, GSBC, 105 Gwanggyo-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, Korea

Tel./fax: +82 31 888 5100; +82 31 888 5105

e-mail: mhcho@genoss.com

Authorized representative of the Manufacturer in the Russian Federation:

RAYMED TRADING GROUP PVT Llc

141001, Moscow region, Mytishchi district, Mytishchi, 1st Vokzalnaya street, bld. 8

tel: +7 (495) 956-83-16 / fax: +7 (499) 706-83-16

e-mail: raymedtgpt@mail.ru

Approved by

Sung Min Chung

Date:

02 Dec 2021

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 2

1. IDENTIFICATION OF THE MEDICAL DEVICE

1.1. Medical device name

GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40

1.2. Medical device name for registration on the territory of the Russian Federation

GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40

1.3. Models

- GENOSS Inflator B30: GBI-B30, GBI-B30S, GBI-B30V, GBI-B30VS;
- GENOSS Inflator B40: GBI-B40, GBI-B40S, GBI-B40V, GBI-B40VS;
- GENOSS Inflator B30 set: GI30V-BA1, GI30V-BA2, GI30V-BB1, GI30V-BB2, GI30-BA1, GI30-BA2, GI30-BB1, GI30-BB2, GI30V-SA1, GI30V-SA2, GI30V-SB1, GI30V-SB2, GI30V-S, GI30-SA1, GI30-SA2, GI30-SB1, GI30-SB2;
- GENOSS Inflator B40 set: GI40V-BA1, GI40V-BA2, GI40V-BB1, GI40V-BB2, GI40-BA1, GI40-BA2, GI40-BB1, GI40-BB2, GI40V-SA1, GI40V-SA2, GI40V-SB1, GI40V-SB2, GI40-SA1, GI40-SA2, GI40-SB1, GI40-SB2.

2. DEVICE DESCRIPTION

2.1. Overview

The GENOSS Inflator B30 and GENOSS Inflator B40 is an inflation/deflation device with an analog pressure gauge that registers pressures from vacuum to 30 and 40 atmospheres (vacuum to 440.88 and 587.84 psi). The GENOSS Inflator B30 and GENOSS Inflator B40 are consisted of a thread plunger and locking button which allow for the generation and sustaining of pressure. A three way stopcock should be included for use during preparation of the device and in conjunction with interventional devices.

Figure 1.
Inflator (GBI-B30)

Figure 2.
Inflator (GBI-B40)

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 3

Table 1.

Inflator (GBI-B30)	Max. Pressure	30 atm (440.88 psi)
	Barrel Capacity	20 ml (cc)
Inflator (GBI-B40)	Max. Pressure	40 atm (587.84 psi)
	Barrel Capacity	30 ml (cc)
3 way stopcock	Max. Pressure	40 atm (587.84 psi)

GENOSS Inflator, execution option GENOSS Inflator B30, set composition: Y-Connector, introducer needle and torque device, execution option GENOSS Inflator B40, set composition: Y-Connector, introducer needle and torque device

The Y-Connector with hemostasis valve is designed to be used on a guiding catheter or dilatation catheter to control backbleeding and to provide a port for introduction of fluids into the interventional system. Y-Connector has two models. GYC-A & GYC-B model of the Y-Connector has two valves that operate independently: the fixed valve and the hemostatic valve. Closing of the fixed valve allows for pressure injections up to 440 psi and also secures the diagnostic / interventional device in position within the vasculature. The hemostatic valve is a diaphragm seal that forms around diagnostic / interventional devices as they move into and out of the vasculature. This valve provides minimal fluid loss without restricting device movement. The hemostatic valve is open when the cap is pressed down, and closed when the cap is released. An open hemostatic valve allows air to be purged from the device and allows the insertion / withdrawal of diagnostic / interventional devices.

The Introducer Needle is designed to facilitate placement of a guide wire tip through the Y-connector and into the wire lumen of an interventional catheter. The introducer needle inserts the valve of Y-Connector, making a tunnel for guidewire to get through and protecting guidewire tip. The introducer needle has two models: 19G and 20G needle.

The Torque Device is designed to hold a small diameter guide wire and provide a handle for manipulating the wire.

2.2. Description of appearance and structure

GENOSS Inflator B30


Figure 3.

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	I
			Page	Page 4

Table 2.

No	Components name	Description
1	Case Upper	Main upper case
2	Case Lower	Main lower case
3	Window	Cover part on case upper
4	Button	Locking part to control the movement of handle
5	Pin	Part fix case upper and case lower
6	Barrel	Marked part to indicate the volume of fluid in the syringe
7	Piston	Cylindrical part to insert into syringe barrel
8	O-Ring	O-ring part to block the leakage of contrast media
9	High Pressure Braided Tube	Moving tube of contrast media
10	Rotating Luer	Connecting part to join the interventional devices
11	Pressure Gauge	Part to monitor the applied pressure
12	Coil Spring	Part keep the binding capability of internal thread and male thread
13	Handle	Control part to pressurize or depressurize by rotating
14	Cap	Cover part to protect the handle

GENOSS Inflator B40


Figure 4.

Table 3.

No	Components name	Description
1	Case Upper	Main upper case
2	Case Lower	Main lower case
3	Button	Locking part to control the movement of handle
4	Barrel	Marked part to indicate the volume of fluid in the syringe
5	Handle	Control part to pressurize or depressurize by rotating
6	Cap	Cover part to protect the handle
7	Window	Cover part on case upper
8	Piston	Cylindrical part to insert into syringe barrel
9	O-Ring	O-ring part to block the leakage of contrast media
10	Pressure Gauge	Part to monitor the applied pressure
11	High Pressure Braided Tube	Moving tube of contrast media
12	Rotating Luer	Connecting part to join the interventional devices
13	Pin	Part fix case upper and case lower

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 5

3-Way Stopcock


Figure 5.

Table 4.

No	Components name	Description
1	Rotating Luer (Female Luer)	Connect with Balloon catheter and Stent hub
2	Male Luer	Two part which is connecting with Rotating connector of Balloon Catheter
3	Off Lever	Decide the direction of fluid

Y-Connector, GYC-A model


Figure 6.

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 6

Table 5.

No	Name	Description
1	Main Body	Inserting channel of introducer needle
2	Rotating Luer	Connecting part to assemble with the hub of the guiding catheter
3	Outer Screw	Control part of the fixed valve to open or close by rotating the cap
4	Inner Screw	This component provides the fixation of hemostatic valve and pressure delivery on the fixed valve
5	Opener(A)	Control part of the hemostatic valve, is open when the cap is pressed down, and closed when the cap is released
6	Washer	Protecting part to keep non-separation from main body and outer screw
7	Hemostatic Valve(A)	This valve provides minimal fluid loss without restricting device movement
8	Fixed Valve	This valve allows for pressure injections up to 440 psi and also secures the diagnostic / interventional device
9	Luer Seal	Prevention ring of leakage

Y-Connector, GYC-B model


Figure 7.

Table 6.

No	Name	Description
1	Main Body	Inserting channel of introducer needle
2	Rotating Luer	Connecting part to assemble with the hub of the guiding catheter
3	Screw	Control part of the fixed valve to open or close by rotating the cap
4	Inner hub	This component provides the fixation of hemostatic valve and pressure delivery on the fixed valve
5	Opener(B)	Control part of the hemostatic valve, is open when the cap is pressed down, and closed when the cap is released
6	Coil spring	This spring provides returning of the opener when the cap is released

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 7

7	Hemostatic Valve(B)	This valve provides minimal fluid loss without restricting device movement
8	Fixed Valve	This valve allows for pressure injections up to 440 psi and also secures the diagnostic / interventional device
9	Luer Seal	Prevention ring of leakage

Introducer needle and torque device


Figure 8.

Introducer needle

A

Size	Usable length (A)	Inner diameter (B)	Outer diameter (C)
19G	100 mm	Ø 0.93 mm	Ø 1.10 mm
20G	100 mm	Ø 0.60 mm	Ø 0.90 mm

Torque device (unit: mm)


Table 7.

No	Name	Description
1	Introducer needle	The introducer needle inserts the valve of Y-connector, making a tunnel for guidewire to get through and protecting guidewire tip. The introducer needle has two models: 19G and 20G needle.
2	Torque device	The torque device holds guidewire and makes torque guidewire move easier.

INSTRUCTION FOR USE

**GENOSS Inflators in sets, execution
options GENOSS B30, GENOSS B40**

for Russia

Date 2021-12

Rev. No. 1

Page Page 8

2.3. Model classification

Table 8.

No	Model	Inflator								Y-Connector			
		Inflator B30 (Pressure gauge, max 30 atm)				Inflator B40 (Pressure gauge, max 40 atm)				GYC-A		GYC-B	
		Material of pressure gauge								GYC-A19	GYC-A20	GYC-B19	GYC-B20
		Brass		Stainless steel		Brass		Stainless steel					
		Stopcock											
Yes	No	Yes	No	Yes	No	Yes	No						
1	GBI-B30		V										
2	GBI-B30S				V								
3	GBI-B30V	V											
4	GBI-B30VS			V									
5	GBI-B40						V						
6	GBI-B40S								V				
7	GBI-B40V					V							
8	GBI-B40VS							V					
9	GI30V-BA1	V								V			
10	GI30V-BA2										V		
11	GI30V-BB1											V	
12	GI30V-BB2											V	
13	GI30-BA1		V							V			
14	GI30-BA2										V		
15	GI30-BB1											V	
16	GI30-BB2											V	
17	GI30V-SA1			V						V			
18	GI30V-SA2										V		
19	GI30V-SB1											V	
20	GI30V-SB2											V	
21	GI30-SA1				V					V			
22	GI30-SA2										V		
23	GI30-SB1											V	
24	GI30-SB2											V	
25	GI40V-BA1									V			
26	GI40V-BA2										V		
27	GI40V-BB1					V						V	
28	GI40V-BB2											V	
29	GI40-BA1									V			
30	GI40-BA2										V		
31	GI40-BB1											V	
32	GI40-BB2											V	
33	GI40V-SA1									V			
34	GI40V-SA2										V		
35	GI40V-SB1											V	
36	GI40V-SB2											V	
37	GI40-SA1									V	V		
38	GI40-SA2											V	
39	GI40-SB1												V
40	GI40-SB2												V

GYC-A and GYC-B components has set composition introducer needle and torque device.

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 9

2.4. Packaging specification

Table 9. GENOSS Inflator B30

Part Name	Material	Size
Tyvek	HDPE	305 x 113 mm
Blister	PET	305 x 113 x 67 mm
Box	Paper	361 x 240 x 320 mm

Table 10. GENOSS Inflator B40

Part Name	Material	Size
Tyvek	HDPE	345 x 113 mm
Blister	PET	345 x 113 x 67 mm
Box	Paper	361 x 240 x 360 mm

2.5. Indications

The GENOSS Inflators are used to inflate and deflate balloon catheters or other interventional devices as well as to monitor pressure within the balloon during the procedure.

2.6. Contraindications

None known.

2.7. Warnings

- Use of these devices should be restricted to those specialists and/or physicians trained to perform vascular procedures that utilize interventional devices such as balloon dilation catheter.
- Always follow the manufacturer's directions accompanying the balloon catheter for instructions for use, maximum inflation pressures, precautions, and warnings for use.
- Improper tightening of the connection between the Inflator and Y-connector and accompanying equipment may result in introduction of air into the vascular system.
- Do not inject any fluid if air bubbles are visible within fluid path of the inflator and Y-connector.
- The pressure displayed on the inflator must not exceed maximum recommended balloon catheter inflation pressure.
- Prolonged exposure to contrast medium may cause components of the device to malfunction and thereby create a safety hazard.
- Use only radiopaque inflation medium recommended by the dilation catheter manufacturer.
- To prevent damage to a balloon, it must be completely deflated during advancement / withdrawal through the Y-Connector.
- Failure to depress the opener prior to inserting a guide wire introducer could cause valve damage resulting in leakage and / or particulate embolization.
- Failure to open the hemostatic valve and the fixed valve prior to inserting / withdrawing a diagnostic / interventional device could cause damage to the diagnostic / interventional device.
- Power injection at pressures greater than 30atm (441psi) could result in leakage or detachment of components. Care should be taken to avoid overtightening the fixed valve, as this may compromise the integrity of the diagnostic / interventional device lumen.

2.8. Precautions

- The components of the package have been delivered STERILE.
- GENOSS Inflators and all accessory products are for single use only; do not reuse.

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 10

- Thoroughly inspect the components prior to use to verify that no damage has occurred during handling.
- Do not use if the package and/or its contents have been damaged, or contamination or damage is suspected. If damage is found, call your GENOSS representative.
- After use, dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Do not resterilize, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing.
- Note the "Use by" date specified on the package.
- Refer to the appropriate device labeling for indications associated with the use of interventional and/or diagnostic devices.
- Do not expose the GENOSS Inflators and accessories to organic solvents, e.g., alcohol.

2.9. Adverse effects

Not applicable.

2.10. Symbols

	Catalogue number		Batch code
	Do not reuse		Date of manufacture
	Caution		Use by
	Sterilized using ethylene oxide		Manufacturer
	CE Mark		Authorized representative
	Do not resterilize		Keep away from sunlight
	Fragile, handle with care		Do not use if package is damaged
	Consult instructions for use		Temperature limitation
	Keep dry		Quantity per package

3. ACCEPTANCE FOR ENVIRONMENTAL REQUIREMENT

GENOSS Inflators and accessories do not contain any toxic or harmful ingredients like phthalate, latex etc.

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 11

GENOSS Inflators and accessories do not adversely impact the environment.

4. STERILIZATION AND SHELF-LIFE

Sterilization of the GENOSS Inflators and accessories have been conducted in accordance with the process validation procedure and EO gas validation, and documented in the EO gas sterilization validation report in accordance with EN ISO 11135.

Revalidation shall be conducted in accordance with the revalidation period (1-year) described in validation report.

Table 11.

Sterilization conditions are as follows:

Method of sterilization	EO gas sterilization, 30% Ethylene oxide + 70% CO ₂
Amount of EO gas injection	700~800mg/L
Temperature	40 ± 5°C
Humidity	50±25%RH
Pressure	0.7±0.05kgf/cm ²
Sterilization time	4H
Flushing	7 times
Aeration temperature and time	10~40 °C, 48 hours
Sterility assurance level	10 ⁻⁶

Shelf-life of the GENOSS Inflators and accessories is 3 years from the date of sterilization.

5. DIRECTION FOR USE

Inflator

- 1) Connect the GENOSS Inflator B30 or GENOSS Inflator B40 to a 10 and 20cc syringe filled with contrast medium or mixed contrast medium with normal saline.
- 2) To fill the GENOSS Inflator B30 or GENOSS Inflator B40 with contrast medium:
 - a) Open the stopcock and press the locking button to its UNLOCK position.
 - b) Aspirate contrast medium by slowly withdrawing the plunger.
 - c) Evacuate air in the GENOSS Inflator B30 or GENOSS Inflator B40 by slowly injecting the necessary volume from the GENOSS Inflator B30 or GENOSS Inflator B40 back into the syringe.
 - d) Repeat steps 2) b) through 2) c) until all the air is removed from the GENOSS Inflator B30 or GENOSS Inflator B40 and about 6~20cc of contrast remains.
- 3) To create a closed system, release the locking button.
- 4) Connect the GENOSS Inflator B30 or GENOSS Inflator B40 to the inflation port of the interventional device such as balloon catheter hub.
- 5) When balloon inflation is desired, rotate the plunger handle clockwise to increase pressure or counterclockwise to decrease pressure.
- 6) Pressure can be quickly released by pressing the locking button to the UNLOCK position and then retracting the plunger handle. To maintain this position, release the locking button to the LOCK position.
- 7) To quickly release from the vacuum position, press the locking button to the UNLOCK position.

Devices for joint use

- 1) Prior to the interventional procedure, all equipment to be used, including the diagnostic / interventional devices should be examined carefully for defects. Do not use any defective equipment.
- 2) Attach connective tubing to a manifold and the sideport of the Y-Connector set.
- 3) Fill and flush the assembly with normal saline. Open the fixed valve and hemostatic valve, place

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 12

one finger over the luer fitting, and fill the assembly.

4) Prepare diagnostic / interventional device according to the manufacturer's instructions.

5) Connect the Y-Connector set to the guiding catheter.

6) Aspirate the system to remove any trapped air before introducing devices.

7) Insert the Introducer Needle and guide wire through an open hemostatic valve and an open fixed valve.

8) Remove the Introducer Needle.

- Insert the rapid exchange (rail) dilatation catheter through an open hemostatic valve and an open fixed valve.

- Insert the over-the-wire dilatation catheter and guide wire through an open hemostatic valve and an open fixed valve.

9) Advance the diagnostic / interventional device through an open fixed valve.

10) Once the diagnostic / interventional device is in place, if desired, close the Fixed valve to allow for pressure injections up to 30atm (441psi) and / or secure the diagnostic / interventional device in position.

11) Remove the diagnostic / interventional device through an open fixed valve. If resistance is felt, open the hemostatic valve as well.

6. TRANSPORTATION, STORAGE TERMS

Handle with care. Store in original packaging. Avoid exposure to direct sunlight or heats. Store in a cool, dark and dry place.

Devices packaged by manufacturer may be transported by all modes of transport in covered vehicles in accordance with the rules of transportation of goods for each type of transport.

The temperature at transportation should be between 1°C to +30°C with a range of $\pm 2^\circ\text{C}$ and relative humidity from 10% to 75% with a range of $\pm 5\%$.

GENOSS Inflators and accessories should be stored at controlled room temperature, nominal 1°C to 30°C and relative humidity from 10% to 75%, in an enclosed space in the original manufacturer's packaging for the labeled shelf life.

7. DISPOSAL

Products, which had contact with blood and/or other biological liquids of patients should be disposed and eliminated in accordance with the requirements for handling medical waste, established by SanPiN (Sanitary rules and regulations) 2.1.3684-21 for class B and other normative legal acts of the Russian Federation.

Unused products, including the ones which expired (which had no contact with the biological liquids of patients) should be disposed in accordance with the procedure of disposal of solid communal waste established by the legislation of the Russian Federation for class A by SanPiN 2.1.3684-21.

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	I
			Page	Page 13

8. STANDARDS

Table 12.

EN ISO 13485:2016	Medical devices - Quality management system - Requirements for regulatory purposes
EN ISO 14971:2012	Medical devices - Application of risk management to medical devices
EN ISO 10993-1:2009	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 10993-4:2017	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood
EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-7:2008	Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals
EN ISO 10993-10:2013	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
EN ISO 10993-11:2009	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
EN ISO 10993-12:2012	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
ISO 594-1:1986	Conical fittings with a 6% (LUER) taper for syringes, needles and certain other medical equipment - Part 1: General requirements
ISO 594-2:1998	Conical fittings with a 6% (LUER) taper for syringes, needles and certain other medical equipment - Part 2: Lock fittings
EN ISO 11607-1:2009	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 11607-2:2006	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
EN ISO 11135:2014	Sterilization of health-care products. Ethylene oxide. Requirements for the development, validation and routine control of a sterilization process for medical devices
EN ISO 11138-2:2017	Sterilization of health care products - Biological indicators
EN ISO 11140-1:2009	Sterilization of health care products - Chemical indicators
EN ISO 11737-1:2018	Sterilization of medical devices-Microbiological methods – Part 1: Determination of a population of microorganisms on products
EN ISO 11737-2:2010	Sterilization of medical devices-Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
EN ISO 10993-7:2008	Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals
EN ISO 15223-1:2016	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied
EN 1041:2008	Information supplied by the manufacturer with medical devices
EN 62366-1:2015	Medical devices - Application of usability engineering to medical devices

	INSTRUCTION FOR USE		for Russia	
	GENOSS Inflators in sets, execution options GENOSS B30, GENOSS B40		Date	2021-12
			Rev. No.	1
			Page	Page 14

9. LIMITED WARRANTY

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. Genoss Co., Ltd. warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2020 Under Law No.138.

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