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6. DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS OF MD AND PRINCIPLE OF OPERATION

General description of functional elements of MD

The Dual Applicator for open surgery for application surgical sealant VeraSeal (Вирасиал®) is manufactured by Ethicon and connects to the pre-filled syringes in a simple click-on design. It is a single-use delivery system and allows the VeraSeal Solutions for Sealant to be sprayed or dripped without the need for a gas assisted spraying system.

The Dual Applicator for open surgery application of surgical sealant VeraSeal (Вирасиал®) is composed of an adapter with a click-on design, as well as a malleable tip that allows the surgeon control of the position of the tip. The malleable tip consists of a wire encased in a tri-lumen tube, which is overmolded on either side. One end is ultrasonically welded to yield a manifold which can be connected via threaded Luer locks to the adapter while the other end is the drip tip. The device is provided with an Airless Spray Tip threaded onto the drip tip, the removal of which allows the surgeon to select which function they desire. The fluid path design allows the two biological components to travel from the pre-filled syringes, through the connectors in the adapter into the manifold and down the length of the tri-lumen tubing and overmold separately before mixing in the Airless Spray Tip.

The VeraSeal syringe holder has two syringe caps, which are removed when the Biologic is ready to use. Once the caps are removed, the Dual Applicator for open surgery for application surgical sealant VeraSeal

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(Вирасиал®) and Adapter are connected to the VeraSeal syringe holder and it then ready to use as illustrated in.


Figure 1 – Illustration of Attachment of the VeraSeal Dual Applicator (left) to the VeraSeal Solutions for Sealant Syringe Holder (right)

When the plunger is depressed, the VeraSeal Solutions for Sealant is driven from the glass syringes and flows through the Dual Applicator for open surgery for application of surgical sealant VeraSeal (Вирасиал®). The Thrombin (Human) and Fibrinogen (Human) travel through two Luer Connectors (11), an Adapter Shroud (9), a Manifold (8), an Overmolded Tip consisting of a Trilumen Malleable Tip (3) with an overmolded connection on one end and a Threaded Connector (4) on the other, and a removable Airless Spray Tip (6) (Figure 2).


Figure 2 - VeraSeal Dual Applicator Components

- | | |
|--|------------------------|
| 1. O-Ring | 7. Overmolded Tip |
| 2. Luer Locks | 8. Manifold |
| 3. Tri-Lumen Tubing (Malleable Wire encased) | 9. Adapter Shroud |
| 4. Threaded Connector | 10. Compression Spring |
| 5. O-Ring | 11. Luer Connector |
| 6. Airless Spray Tip | 12. Connector Retainer |

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The Luer Connector maintains a seal to the glass syringes via the Compression Springs loading the Luer Connector in the proximal direction. O-Rings maintain the seal between the Luer Connectors and Adapter Shroud. The Manifold and Adapter Shroud maintain a seal via the Luer Locks securing the Manifold to the Adapter Shroud. The Manifold and Overmolded Tip maintain a seal via an ultrasonic welding process that bonds them together. The Malleable Tip maintains a seal between the Tri-Lumen Tubing and the overmolding on either end via over-molded connections. An O-Ring on the Threaded Connector (Figure 3) ensures a seal with the Airless Spray Tip. In case of clogging, two additional Airless Spray Tips are provided with each device.


Figure 3 - Airless Spray Tip Components

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Connecting the Dual Applicator for open surgery for application of surgical sealant VeraSeal (Вирасиал®) to the VeraSeal Syringe Holder with pre-filled syringes allows the user to spray VeraSeal Solutions for Sealant. When the two components of the VeraSeal Solutions for Sealant exit through the distal tip, they are first mixed within the pre-attached Spray Tip and then exit in a conical spray pattern. If dripping of the fibrin sealant is needed for the indicated tissue site, the pre-attached Spray Tip is removed and the two components will only mix when Fluid Path of the VeraSeal Solutions for Sealant through the Dual Applicator for open surgery for application surgical sealant VeraSeal (Вирасиал®).


Figure 4 - Mixing of components of seal in VeraSeal Dual Applicator.

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*The table of technical characteristics contains the actual parameters of the MD. The parameters specified in the MD name are part of the name and presented for the convenience of the user's choice of a particular product.

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Table 2 – The technical characteristics of Dual Applicators for open surgery for application surgical sealant VeraSeal (Вирасиал®)

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Table 5 – The technical characteristics of the Dual Applicators for open surgery for application surgical sealant VeraSeal (Вирасиал®) with the airless spray tip

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8. LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Not applicable. This medical device does not contain materials of animal and / or human origin.

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9. INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD

The devices do not contain medicinal and/or pharmaceutical substances. Information on the medicinal substance used in combination with medical device:

Name: VeraSeal®

International non-proprietary name: human fibrinogen + human thrombin

Dosage form: solutions set for the preparation of a surgical sealant

Manufacturer: "Instituto Grifols, SA", Spain

10. Information on the procedure for installation, installation, adjustment, calibration, and other actions required to put medical products into operation

Circulating Nurse

Open the blister and present the VERASEAL™ Dual Applicator and two additional Airless Spray Tips to the Scrub Nurse for transfer to the sterile field. The outside of the blister should not come in contact with the sterile field.

Scrub Nurse

1. Hold the VERASEAL™ solutions for sealant syringe holder with syringe Luer caps pointed upward. Unscrew and discard the syringe Luer cap of both fibrinogen and thrombin syringes (Fig. 9).

Figure 9

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2. Hold the syringe holder with the Luers pointed upward. To remove air bubbles from syringes, strike gently the side of the syringe holder one or two times while keeping the syringe holder in an upright position and lightly depress the plunger to eject air (Fig. 10).


Figure 10

3. Attach the Dual Applicator (Fig. 11).

4. NOTE: Do not depress plunger during attachment or prior to intended use, because the two biological components will pre-mix in the Airless Spray Tip, forming a fibrin clot that prevents dispensing (Fig. 12).

5. Tighten Luer locks and ensure the Dual Applicator is firmly attached. The device is now ready to use.


Figure 11

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
Figure 12

Application by Spraying

1. Grasp and bend the Dual Applicator to the desired position. Tip will retain its shape.
 2. Position the Airless Spray Tip at least 2 cm away from the target tissue. Apply firm even pressure to the plunger to spray the fibrin sealant. Increase distance accordingly to achieve desired coverage of the target area.
 3. If expression is stopped for any reason, change the Airless Spray Tip. To change the Airless Spray Tip, remove the device from the patient and unscrew the used Airless Spray Tip (Fig. 13). Place the used Airless Spray Tip away from the spare Airless Spray Tips. Wipe the end of the applicator using dry or moist sterile surgical gauze. Then, connect a new Airless Spray Tip provided in the package and ensure it is firmly connected before use.
- NOTE: Red indicator will not be visible if Airless Spray Tip is properly connected (Fig. 14).

NOTE: Do not continue pushing the plunger in an attempt to clear the fibrin clot within the Airless Spray Tip; otherwise, the applicator may become unusable.

NOTE: Do not trim the Dual Applicator

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
Figure 13


Figure 14

Do not apply VERASEAL™ solutions for sealant intravascularly.

Application by Dripping

1. Remove the Airless Spray Tip portion of the spray and drip tip by unscrewing the Airless Spray Tip (Fig. 13).
2. Grasp and bend the drip tip to the desired position. Tip will retain its shape.
3. During dripping, keep the end of the drip tip as close to the tissue surface as possible without touching the tissue during application.
4. Apply individual drops to the surface area to be treated. To prevent uncontrolled clotting, allow the drops to separate from each other and from the end of the drip tip.

NOTE: Do not reconnect a used drip tip after it has been removed from the adapter; otherwise, a clot may form inside the drip tip and the applicator may become unusable.

Exchange of the Spray and Drip Tip (as needed)

1. Loosen the Luer locks to remove the spray and drip tip from the adapter (Fig. 15).
2. Attach the desired VERASEAL™ Applicator to the adapter by tightening the Luer locks (Figs. 16 and 17).

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Name	Material of MD Trade / Common Name
Dual Applicators for open surgery for application surgical sealant VeraSeal (Вирасиал®)	
Airless Spray Tip (Tip Assembly)	PC, coloring materials
Overmolded Threaded Connector	ABS, (Natural) Grey coloring materials
Adapter Shroud, Luer Connector, Overmolding and Manifold	ABS, coloring materials (SnoWhite)
O-Ring	SILICONE RUBBER, COLOR: RED
Tri-lumen Tubing	Thermoplastic polyurethane

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15. CONDITIONS OF TRANSPORTATION AND STORAGE

Transportation conditions

Temperature	-20°C to +50°C
Humidity	to 75%
Barometric pressure	No applicable

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	"STERILE" - Part 1: Requirements for terminally sterilized medical devices
EN ISO 11137-1: 2015 /A2: 2019	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices- Amendment 2: Revision to 4.3.4 and 11.2 (ISO 11137-1:2006/Amd 2:2018)
EN ISO 11137-2:2015	Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose
EN ISO 11737-1:2018 AMD 2021	Sterilization of health care products - Microbiological methods -Part 1: Determination of a population of microorganisms on products
EN ISO 11737-2:2020	Sterilization of health care products - Microbiological methods -Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
EN ISO 11607-1:2020	Packaging for Terminally Sterilized Medical Devices - Part 1: Requirements for Materials, Sterile Barrier Systems and packaging systems (ISO 11607-1:2019)
EN ISO 11607-2:2020	Packaging for Terminally Sterilized Medical Devices - Part 2: Validation Requirements for Forming, Sealing and Assembly Processes (ISO 11607-1:2019)
EN ISO 14644-1:2015	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness
EN ISO 14644-2:2015	Cleanrooms and associated controlled environments - Part 2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1
EN ISO 14644-3: 2019	Cleanrooms and associated controlled environments - Part 3:

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18. DATA ABOUT LABELING OF MD AND ITS PACKAGING

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Dual disposable sterile applicators for applying surgical sealant
VeraSeal (Вирасиал®):

Laparoscopic Dual Applicators for application surgical sealant
VeraSeal (Вирасиал®) 35 cm, rigid;

Laparoscopic Dual Applicators for application surgical sealant
VeraSeal (Вирасиал®) 45 cm, flexible

Subscribed and sworn to before me
this 26th September 2023

Keith R. Petersen

KEITH R. PETERSEN
Notary Public - State of New Jersey
My Commission Expires Jul 23, 2026

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1. MEDICAL DEVICE (MD) NAME AND CLASSIFICATION

MD name, variations of MD

Dual disposable sterile applicators for applying surgical sealant VeraSeal (Вирасиал®), variants:

1. Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 35 cm, rigid:
 - 1.1. Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 35 cm, rigid with the airless spray tip – 3 pcs.
 - 1.2. Airless spray tips – 6 pcs.
 - 1.3. Instructions for the use of a medical device - 1 pcs.
2. Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 45 cm, flexible:
 - 2.1. Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 45 cm, flexible with the airless spray tip – 3 pcs.
 - 2.2. Airless spray tips – 6 pcs.
 - 2.3. Instructions for the use of a medical device - 1 pcs.

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Further in the text the following names of the medical device can be used: medical device (MD), VeraSeal, applicator, Laparoscopic Dual Applicators.

The medical devices names can be mentioned in the text in abbreviated form for information perception convenience.

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2. NAME AND LEGAL ADDRESS OF MANUFACTURER AND AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION

DEVELOPER: Ethicon, Inc. 1000 Route 202 Raritan, New Jersey 08869 USA

MANUFACTURER: Ethicon, Inc. 1000 Route 202 Raritan, New Jersey 08869 USA

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DL-JNJRUMedical-Regulatory@ITS.JNJ.com

PLACE OF PRODUCTION:

1. Viant AS&O Holdings, LLC a.t.a (formerly) Lake Region Medical 45 Lexington Drive Laconia, NH 03246, USA
2. Tessa Plastics (Shanghai) Co., Ltd., Building 9#, No. 509, Ren Qing Road, East Area, Zhang Jiang High-Tech Industrial Park 201201 Pudong, Shanghai – China

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA:

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3. INTENDED USE MEDICAL DEVICE WITH AN INDICATION OF A POTENTIAL USER OF MD

Intended to be used for the simultaneous topical application of the two biological components of the VeraSeal Solutions for Sealant onto the surface.

Potential users of MD

VeraSeal Applicators are intended for use only by healthcare professionals who are trained in surgical procedures and techniques requiring the use of surgical sealants.

The healthcare professionals are advised to consult the medical literature regarding techniques, complications, and adverse reactions prior to undertaking endoscopic procedures.

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4. INDICATIONS AND CONTRAINDICATIONS FOR MD USE AND POSSIBLE ADVERSE EFFECTS

Indications

The disposable sterile Dual Applicators for applying surgical sealant VeraSeal (Вирасиал®) are intended to be used for the simultaneous topical application of the two biological components of VeraSeal Solutions for Sealant onto the surface.

Contraindication

None known.

Adverse reactions and events

See warnings and precautions

5. DIRECTIONS FOR USE

Use/operation conditions of MD

Devices are recommended for use under operation conditions indicated in Table 1.

Table 1 – Operation conditions

Operation conditions	
Humidity	not standardized
Temperature	32-42 °C
Barometric pressure	not standardized

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Removing the pre-attached Airless Spray Tip from the Spray and Drip Tip Assembly allows the user to drip VeraSeal Solutions for Sealant onto a surface. Luer Locks can be loosened to remove the pre-attached Spray and Drip Tip Assembly in order to attach another Dual Applicator tip for application of surgical sealant VeraSeal (Вирасиал®), as licensed for use with

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VeraSeal. For example, the Laparoscopic Dual Applicators for application of surgical sealant VeraSeal (Вирасиал®) 35 cm, Rigid or Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 45 cm, Flexible tip.

The Adapter that is pre-attached to the Dual Applicator for open surgery for application surgical sealant VeraSeal (Вирасиал®) is used with the VeraSeal Laparoscopic Dual Applicators also.

To attach to the Adapter:

- Remove the Spray and Drip Tip from the Adapter of the VeraSeal Dual Applicator;
- Attach the VeraSeal Laparoscopic Dual Applicator to the Adapter.

As shown in Figure 1.


Figure 1 – Removal of the Spray and Drip Tip from the Adapter of the VeraSeal Dual Applicator (left) and Attachment of the VeraSeal Laparoscopic Dual Applicator to the Adapter (right) /

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Remove the two syringe Luer caps from the VeraSeal Solutions for Sealant Syringe Holder and then connect the Laparoscopic Dual Applicators for application of surgical sealant VeraSeal (Вирасиал®) and Adapter to the VeraSeal Syringe Holder.

When the plunger is depressed, the VeraSeal Solutions for Sealant is driven from the pre-filled syringes and flows through the Laparoscopic Dual Applicators for application of surgical sealant VeraSeal (Вирасиал®).

The Thrombin (Human) and Fibrinogen (Human) travel through two Luer Connectors, an Adapter Shroud, a H-Connector, Internal Tubes, a Threaded Connector, and a removable Airless Spray Tip. Figure 2 and Figure 3 provide the Laparoscopic Dual Applicators for application of surgical sealant VeraSeal (Вирасиал®) component parts.


Figure 2 – VeraSeal Laparoscopic Dual Applicator 35 cm Rigid Component Parts

1. Threaded Connector And O-Ring
2. Stainless Steel Tube
3. Shroud
4. Internal Tubes
5. H-Connector

6. Shroud
7. Tip Insert
8. Tip Housing
9. Orifice Cup
10. Spray Tip

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
Figure 3 – VeraSeal Laparoscopic Dual Applicator 45 cm Flexible Component Parts

1. Threaded Connector And O-Ring
2. Flexible Over-Tube
3. Shroud
4. Internal Tubes
5. H-Connector
6. Grommet

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Laparoscopic Dual Applicators for application of surgical sealant VeraSeal (Вирасиал®) (35 cm Rigid or 45 cm Flexible) is comprised of a dual lumen design such that the biologics do not mix until they exit the threaded connector.


Figure 4 – Mixing of components of seal in VeraSeal Laparoscopic Dual Applicator

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When dripping, the biologics exit the dual lumen threaded connector. When spraying, the biologics enter the Spray Tip comprised of a static mixing element, a pattern of swirl chambers, and a small orifice. The biologics mix and flow through the Spray Tip and then exits through the small orifice in a conical spray pattern. If spraying is stopped and the user wishes to continue spraying, the user must replace the Spray Tip with one of the two spare Spray Tips provided in the package.

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Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 35 cm, rigid / 45 cm, flexible

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Table 2 – The technical characteristics of the Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®)

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VeraSeal Laparoscopic Dual Applicator (35 cm Rigid and 45 cm Flexible)

Circulating Nurse

Remove the VeraSeal Dual Applicator and VeraSeal solutions for sealant syringe holder with pre-filled syringes from the package and present them to the Scrub Nurse for delivery into the sterile field.

Remove the VeraSeal Laparoscopic Dual Applicator from the package and present the VeraSeal Laparoscopic Dual Applicator with two Airless Spray Tips to the Scrub Nurse for delivery into the sterile field.

Scrub Nurse

Prepare the VeraSeal solutions for sealant syringe holder with pre-filled syringes per the VeraSeal solutions for sealant Package Leaflet.

Loosen the Luer locks to remove the spray and drip tip from the adapter. (Fig. 6)

Attach the VeraSeal Laparoscopic Dual Applicator to the adapter by tightening the Luer locks. (Fig. 7)

Attach the VeraSeal Laparoscopic Dual Applicator and adapter onto the syringe holder with pre-filled syringes. (Fig. 8)

Do not depress the plunger to prime the device, because the two biological components will pre-mix in the Airless Spray Tip, forming a fibrin clot that prevents dispensing. (Fig. 9)


Figure 6

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
Figure 9

Spray Application

The device is now ready to spray. Before entering the patient, ensure that the Airless Spray Tip is firmly connected.

NOTE: Red indicator will not be visible if properly connected. (Fig. 10)

To change the Airless Spray Tip, remove the VeraSeal Laparoscopic Dual Applicator from the patient and unscrew the used Airless Spray Tip. (Fig. 10)

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Place the used Airless Spray Tip away from the spare Airless Spray Tips. Wipe the distal end of the Applicator using dry or moist sterile surgical gauze. Then connect a new Airless Spray Tip provided in the package and ensure it is firmly connected before use.

NOTE: Red indicator will not be visible if properly connected. (Fig. 10)


Figure 10

Preparation for Drip Application

Remove the Airless Spray Tip from the Applicator. The device is now ready to drip.

Do not apply VeraSeal solutions for sealant intravascularly.

For spraying, position the Airless Spray Tip at least 2 cm from the target tissue. Do not depress the plunger to remove air bubbles or prime the device. Apply firm, even pressure to the plunger to spray the fibrin sealant. Increase distance accordingly to achieve desired coverage of the target area.

NOTE: Apply continuous pressure during spraying. If expression is stopped for any reason, change the Airless Spray Tip. Do not continue pushing the plunger in an attempt to clear the fibrin clot within the Airless Spray Tip; otherwise, the Applicator device may become unusable.

For dripping, keep the end of the Applicator as close to the tissue surface as possible without touching the tissue during application. Apply individual drops to the surface area to be treated. Allow the drops to separate from each other and from the end of the Applicator.

The VeraSeal Laparoscopic Dual Applicator is not intended to lift or manipulate tissue. Do not trim the VeraSeal Laparoscopic Dual Applicator.

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11. REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD

The product is single-use, does not require maintenance and is not subject to repair. In case of damage, the product must be disposed of in accordance with the requirements of the legislation in the country of application.

12. INFORMATION ABOUT STERILITY

The VeraSeal Applicator devices are provided sterile and are provided Gamma sterilized.

VeraSeal Applicators are single-used cannot be re-sterilized.

13. INFORMATION ABOUT DESIGNING, ENGINEERING, AND VALIDATION OF SOFTWARE USED IN MD

Not applicable.

14. LIST AND DESCRIPTION OF MD MATERIALS

Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 35 cm, Rigid,
Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 45 cm, Flexible,

Shroud	Polycarbonate, coloring materials White
H-connector	Copolyester Tritan, coloring materials Black
Rigid-flex-tubing (rigid only)	Thermoplastic elastomer
Flexible-flex-tubing (flexible only)	Thermoplastic elastomer
Rigid-tube (rigid only)	Stainless Steel
Flexible -over-tube (flexible only)	Thermoplastic Polyurethane, coloring materials Gray

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15.CONDITIONS OF TRANSPORTATION AND STORAGE

Transportation conditions	
Temperature	-20°C to +50°C
Humidity	to 75%
Barometric pressure	No applicable
Storage conditions	
Temperature	-20°C to +50°C
Humidity	to 75%
Barometric pressure	No applicable

16. LIST OF RISKS

Risk analysis of instruments is fulfilled in accordance with requirements of ISO 13485 Medical Devices – Quality Management Systems – Requirements for Regulatory Purposes, EN ISO 14971 Medical Devices – Risk Management – Part 1: Application of Risk Analysis to Medical Devices, EC directive concerning medical devices 93/42/EEC.

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17. INFORMATION ON THE TESTS CARRIED OUT, TEST REPORTS, ANALYSIS OF THE OBTAINED DATA

The medical device has passed the verification and validation tests and meet the requirements of the following standards:

Applied Standard	Description
EN ISO 13485:2016	Medical device - Quality management system – Requirements for regulatory purposes
EN ISO 14971:2019	Medical devices - Application of risk management to medical devices
EN 1041:2008+A1:2013	Information supplied by the manufacturer of medical devices
EN ISO 15223-1:2021	Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied -- Part 1: General requirements
IEC 62366-1:2015 A1 2020	Medical Devices- Part 1: Application of usability engineering to medical devices
EN ISO 10993-1: 2020	Biological evaluation of medical devices – Part 1: Evaluation and testing within the risk management process
EN 556-1: 2001/AC:2006	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
EN ISO 11137-1: 2015 /A2: 2019	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices- Amendment 2: Revision to 4.3.4 and 11.2 (ISO 11137-1:2006/Amd 2:2018)
EN ISO 11137-2:2015	Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose

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

EN ISO 11737-1:2018 AMD 2021	Sterilization of health care products - Microbiological methods -Part 1: Determination of a population of microorganisms on products
EN ISO 11737-2:2020	Sterilization of health care products - Microbiological methods -Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
EN ISO 11607-1:2020	Packaging for Terminally Sterilized Medical Devices - Part 1:Requirements for Materials, Sterile Barrier Systems and packaging systems (ISO 11607-1:2019)
EN ISO 11607-2:2020	Packaging for Terminally Sterilized Medical Devices - Part 2:Validation Requirements for Forming, Sealing and Assembly Processes (ISO 11607-1:2019)
EN ISO 14644-1:2015	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness
EN ISO 14644-2:2015	Cleanrooms and associated controlled environments - Part 2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1
EN ISO 14644-3:2019	Cleanrooms and associated controlled environments - Part 3: Test methods
EN ISO 14698-1:2003	Cleanrooms and associated controlled environments - Biocontamination control - Part 1: General principles and methods

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Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information

The VeraSeal Laparoscopic Dual Applicator is a non-implantable device and therefore MRI compatibility is not applicable. No known ED substances are present at >0.1%. One or more components of this device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight (for Laparoscopic Dual Applicators for application surgical sealant VeraSeal (Вирасиал®) 35 cm, Rigid):

 #	Material Present	Residual Risk
1	Cobalt: CAS No. 7440-48-4 EC No. 231-158-0	Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless-steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

The footnote adjacent to the above symbol on the packaging refers to the relevant Material Present.

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20. MD DISPOSAL RULES AND CONDITIONS

Discard open packaging and unused devices.

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