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Tel.: 93 586 62 00*
Fax: 93 588 10 96
www.bbraun.es

TO THE
FEDERAL SERVICE ON
SURVEILLANCE IN HEALTHCARE
(ROSZDRAVNADZOR)
Department of Registration of
Medical Devices

Covering letter Operating Documentation

With this letter we

B. Braun Surgical, S.A.
Carretera de Terrassa
121, 08191 Rubí (Barcelona)
Spain

as the responsible manufacturer of medical device:

Surgical suture sterile synthetic non-absorbable monofilament Optilene®
(hereinafter referred to as "medical device")

herewith declare that we provide the following Instructions for use (IFU) in English language (as separate document) in order to update the documents of the Registration Dossier № РД-12529/45243 issued 04.08.2016 without changing the properties and characteristics, functional purpose and principle of operation of the medical device and for submission of the above document to the Federal Service for Surveillance in Healthcare (Roszdravnadzor).

To complete this declaration, this document is approved and signed in Rubí, 23rd of November 2021

For and on behalf of

B | BRAUN

B. Braun Surgical SA
Carretera de Terrassa, 121
08191 Rubí (Barcelona)
Tel.: 93 586 62 00*

Silvia Orús
Global Regulatory Affairs Manager
CoE OR Supply
B. Braun Surgical, S.A.

Martina Laporte
Quality & Regulatory Affairs Director
CoE, OR Supply
B. Braun Surgical, S.A.

**BLANK
PAGE**

Instructions for use of the medical device

Surgical suture sterile synthetic non-absorbable monofilament Optilene®

Versions*:

1. Surgical synthetic non-absorbable monofilament thread Optilene® with atraumatic needle
2. Surgical synthetic non-absorbable monofilament thread Optilene® with atraumatic needle (2 pc.) and soft PTFE pledget
3. Surgical synthetic non-absorbable monofilament thread Optilene® with atraumatic needle (2 pc.)

* - please see the full description of the product in the Registration Certificate and on the product label

Description of the device

Optilene® is a sterile, non-absorbable, monofilament surgical suture material produced from polypropylene and polyethylene (hereinafter referred to as "Optilene® sutures"). Optilene® sutures are dyed blue with copper phthalocyanine (Phthalocyaninato(2-)copper) for better visibility. Optilene® meets all requirements of the European Pharm. and the United States Pharm. for sterile, non-absorbable sutures. Optilene® sutures are available from size USP 2 (5 metric) to size USP 10/0 (0,2 metric), cut to different lengths and attached to various types and sizes of stainless steel needles.

Special presentations containing Optilene® include:

- Optilene® with Premipatch® PTFE pledgets made from 100% polytetrafluoroethylene (PTFE)
- Optilene® size USP 5/0 (1 metric) to USP 8/0 (0,4 metric) sutures for cardiovascular indications with almost square bodied CV2 Pass (square with rounded corners) needles with a thicker layer of silicone coating.

Indications for use

Optilene® is indicated for soft tissue approximation and/or ligation in general surgery, when surgical practice requires the use of a non-absorbable suture material. Optilene® sutures are also for use in cardiovascular surgery and neurosurgical procedures, where the low thrombogenic properties of polypropylene (thread material) have a positive influence on wound healing. Other indications for Optilene® include skin closure, microsurgery and ophthalmology.

The Premipatch® PTFE pledgets are intended for the use as pads between the suture and the tissue surface in order to distribute the load of the suture over a bigger tissue area.

CV Pass, CV2 Pass needles are indicated for cardiovascular surgery.

Contraindications

The implantation of this suture is contraindicated in patients with known sensitivities or allergies towards its components (polypropylene, polyethylene, copper phthalocyanine).

Usage of cutting needles is contraindicated for suturing parenchymal organs or other friable tissues where a cutting needle is not indicated.

Mode of application

Sutures should be selected depending on the tissue to be sutured, the specific duration of wound support needed, the size of the wound, the patient's condition and the specific suturing technique. Open the non-sterile outer packaging, remove the sterile inner packaging, take the needle with the needle holder and remove the thread by moving to the right.

Upon the employment of Optilene®, the suture elicits a mild acute inflammatory reaction in the tissue followed by an encapsulation of the suture by fibrous connective tissue. Optilene® sutures are not absorbed by the human body and are not known to lose a significant amount of its tensile strength after implantation.

Warnings

Optilene® and all its special presentations must not be resterilized. Open unused or damaged packs shall be discarded.

Do not reuse: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

If the needle detaches and falls into the patient and cannot be recovered, exist risk of embolism, foreign body reaction, encapsulation and potential need of additional test. Therefore, the risk is considered critical for the patient.

Precautions

- Optilene® shall be applied by trained staff and surgeons, who must be familiar with the surgical technique and the medical device characteristics in hospital. Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury.
- Contaminated or infected wounds should be treated following the appropriate surgical practice.
- Optilene® should be used applying the standard surgical suturing and knotting (flat and square ties) techniques, taking into account the surgeon's experience with the surgical procedure. Care should be taken that the knots are positioned properly and adequate knot security is given. Additional throws may be appropriate if required by the surgical circumstances.
- When working with Optilene® suture material, great care should be taken to avoid any crushing or crimping damage of the monofilament by instruments such as forceps or needle holders.

- Care should be taken to avoid damaging the needle when using the suture material. Always grasp the needle in a section 1/3 to 1/2 of the distance from the fibre attachment end to the needle point, never at the end where the fibre is attached or the needle point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage.
- Reshaping needles should be avoided and may result in a loss of their strength and resistance towards bending and breaking.
- Discard used needles in "Sharps" containers.

Side Effects

As for every other suture material, prolonged contact with salt solutions such as urine and bile can lead to lithiasis. The following side effects may be associated with the use of this product: transitory local irritation, transient inflammatory foreign body reaction, enhanced bacterial infectivity, wound dehiscence, granulation, stitch sinus, pain, palpable and/or irritating knot, haemorrhage, impaired cosmetic result, burst abdomen, incisional hernia.

Sterilization

The suture material Optilene® is sterilized in secondary (consumer, group) packaging with ethylene oxide. Sterile products are guaranteed until the expiration date of products in unopened and undamaged packaging, subject to storage conditions.

Symbols to be used on the label

Symbol	Transcript
	Do not reuse
	Use by
	Batch number
	Article number
	CAUTION (See instructions for use)
	Method of sterilisation: Ethylene oxide
Do not restirielize	Do not restirielize
 0123	CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.
	Mark of conformity with the declaration of conformity (if necessary)

Storage conditions

To maintain the sterility and expiry date the suture material Optilene® should be stored and transported at the temperature of -30 °C to + 40 °C, relative humidity not exceeding 80%..

Operating conditions

During operation, the product should withstand the effects of temperature +10 - +42 ° C and humidity up to 100%.

Packaging, scope of supply

Medical device in primary individual packaging (sterile, type BLT, RCP, .RCP) in quantity 12 / 24 / 36 pcs. and instructions for use (1 pc.) are placed in a secondary (consumer) packaging. The secondary packaging is covered from the outside with a hermetically sealed protective polypropylene film.

Symbols of surgical suture material synthetic non-absorbable sterile Optilene®

Thread diameter:

2 (5); 1 (4); 0 (3,5); 2/0 (3); 3/0 (2); 4/0 (1,5); 5/0 (1); 6/0 (0,7); 7/0 (0,5); 8/0 (0,4); 9/0 (0,3); 10/0 (0,2) - USP (EP metric (m)).

Thread length: 5CM; 12CM; 13CM; 15CM; 23CM; 30CM; 45CM; 60CM; 75CM; 90CM; 100CM; 120CM; 150CM; 200CM; 240CM, cm.

Needle angle of curvature: H - 1/2 circle; D - 3/8 circle; F - 5/8 circle; G - straight; P - progressive curved.

Body type of needle: R – taper L – lancet; S - cutting, including reverse cutting.

Point type of needle: C - short- cutting; m – microneedle; mV(MV) - micro-vascular; MP – precision point; the symbol for taper and lancet point type of the needle is absent. In this case, the point type of the body depends on the type of needle: taper and lancet respectively.

Needle length (mm)*: 4; 5; 6; 7; 8; 9; 10; 11; 12; 13; 15; 16; 17;18; 19; 20; 21; 22; 24; 26; 30; 32; 35; 36; 37; 39; 40; 43; 45; 48; 51; 60; 65; 90. * Approximate rounding to integers, the updated value of the needle length is indicated by a numeric symbol.

Type and size of pledget: PS3 - soft pledget, 3×3 mm; PS6 - soft pledget, 6×3 mm; PF6 - firm pledget, 6×3 mm.

Additional information: 2X - 2 needles; S - heavy needle; SS - extra heavy needle; F - thin bodied needle; FF - extra thin bodied needle; B - anti-reflective surface; CV - cardio vascular; LOP/LOOP - loop (both ends of thread attached to a needle); 100MIC; 140MIC; 150MI; 150MIC; 200; 200MI; 200MIC - microneedles diameter, µm

Number of pieces per pack: (First letter in the article): G - 12 pcs in packaging; B - 24 pcs in packaging; C -36 pcs in packaging.

Packaging: RCP - RacePack with substrate from polyethylene; .RCP - RacePack with substrate from cardboard; BLT – Blister.

and Optilene is intended for a broad population, adult or paediatric patient regardless of gender and who do not meet contraindications.

The medical device does not contain a medicinal product for medical use, material of animal and (or) human origin.

The shelf life of the suture material Optilene® is 5 years in an undamaged package.

The suture material Optilene® should be transported at the temperature of -30 °C to + 40 °C, relative humidity not exceeding 80% by any closed mode of transport in accordance with the requirements and regulations of the transportation of goods, the action-sponding to the appropriate mode of transport.

Open unused medical device or unused medical device un damaged packaging can not be used, they should be disposed of out in accordance with the requirements of SanPiN 2.1.3684-21 as class A waste and/or other applicable legal acts of the Russian Federation.

Disposal of the medical device after use is carried out in accordance with the requirements of SanPiN 2.1.3684-21 as class B waste and/or other applicable legal acts of the Russian Federation.

The company B. Braun Surgical, SA guarantees that the device meets the requirements of this extract from the technical documentation, subject to conditions, transporting and storage and is not responsible for the use of the product for other purposes.

The manufacturer's warranty applies to the product until the expiration date indicated on the package, and applies to unopened, undamaged, and properly stored products. Shelf life under storage conditions-5 years.

Claims for liability resulting from improper handling or use of medical devices will be rejected.

Manufacturer

B.Braun Surgical, S.A., Carretera de Terrassa, 121, 08191 Rubi (Barcelona), Spain

Information about the authorized representative of the manufacturer in Russia: LLC "B. Braun Medical", 191040, St. Petersburg, ul. Pushkinskaya, 10, Tel.: 8 (812) 320-40-04

Date of Information: 04/2021

BLANK
PAGE

The Instructions for use for the medical device
« Surgical suture sterile synthetic non-absorbable monofilament Optilene® »
is certified by the manufacturer

B. Braun Surgical, S.A.

Silvia Orús

Global Regulatory Affairs Manager
CoE OR Supply
B. Braun Surgical, S.A.

Martina Laporte

Quality & Regulatory Affairs Director
CoE, OR Supply
B. Braun Surgical, S.A.

I, **CARLOS MATEO MARTINEZ BARTOLOME**, Civil Law Notary in the City of Barcelona, Spain, certify that the signatures written on the first page and on the reverse side are the true signatures of **Mrs. SILVIA SALUD ORUS PUIGVERT**, holder of a valid Spanish D.N.I number **33.906.328-G** and **Mrs. MARTINA LAPORTE BREYSSE**, holder of a valid Spanish D.N.I number **46.239.118-X** for its coincidence with the ones they usually use and known by me. This document, written in English that I the Notary understand, consists of FIVE sheets of plain paper sealed and initialed by me, with the security seals number the present and the previous numbering in sequence. I also note that my statement does not apply to the content of the document and that this does not produce public effect to the document. Issued this, Barcelona, 29th November 2021. -----

It stays written down and filed in my register with the number: **3204**

B. Braun Surgical, S.A.
Carretera de Terrassa, 121
Почтовый адрес: Apartado 6
08191 Rubí (Barcelona)
Испания

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