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Intended Use Environment: Hospital, outpatient surgical center.

Clinical Benefit: Safe transection and resection during colorectal procedures to support effective wound closure and wound healing.

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
5. Use/operation conditions of MD


Operation conditions	
Humidity	20 % to 60 %
Temperature	20 °C to 24 °C

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11. Cartridge Anvil
12. Staple Retainer

7. DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS AND PRINCIPLE OF OPERATION OF MD

7.1. General description of MD

The CONTOUR Curved Cutter Stapler is a multifire, single patient use device with a curved head that cuts and staples. The device delivers four staggered rows of titanium staples, with a knife between the second and third row of staples, and creates a 40 mm curved transection. The device is designed with a feature which prevents closing if a used reload or no reload is in the instrument. Another feature is provided to prevent firing unless the closure trigger is latched in the closed position. A retaining pin holds tissue in place and can be positioned either manually or by squeezing the closure trigger.

The instrument may be reloaded five times, for a maximum of six firings per instrument during a single procedure. Each reload cartridge module includes a knife blade with two staggered rows of staples on each side, an anvil, a cutting washer, a retaining pin, and a staple retainer. Reload cartridges are available in two sizes.


Instrument	Reload	Description	Number of Staples	Reload Color	Tissue Thickness
CS40B	CR40B	Standard	46	Blue	1.5 mm
CS40G	CR40G	Thick	46	Green	2.0 mm

NOTE: The blue (standard) and green (thick) reloads can be used interchangeably with the instruments.

This device is intended to be used for transection and resection during colorectal procedures by healthcare professionals.

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

2



3 Grasp the instrument handle in the palm of the hand with the fingers around the closure trigger. Position tissue to be cut and stapled between the anvil and the cartridge.

NOTE: Any tissue covering the hole in the anvil will be pierced by the retaining pin. Reposition the tissue if necessary.

4 Either manually advance the retaining pin using the actuator button on top of the handle to capture the tissue within the jaw opening (Illustration 3A), or automatically advance the retaining pin by squeezing the closure trigger (Illustration 3B).

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

Note that there is an audible click halfway through the closure stroke indicating that the device is in the intermediate position. In the intermediate position, the pin is fully seated in the anvil capturing the tissue, and the jaws are partially open. Reposition tissue within the instrument if desired. The closure trigger can be released at this point, and the instrument will maintain its jaw opening to allow for final positioning of the tissue to be cut and stapled. The retaining pin is fully engaged.

Caution: Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.

NOTE: If using the automatic advancement of the retaining pin, ensure that the tissue is correctly placed to avoid accidental piercing of tissue.

NOTE: Do not grasp the firing trigger before the instrument is to be fired.

NOTE: Recommended operation is to manually advance the Retaining Pin.

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
5 Squeeze the closure trigger and handle together until the closure trigger is latched. Listen for an audible click. When the closure trigger is latched, the firing trigger moves to the ready-to-fire position (Illustration 4).

The cartridge has clamped onto the tissue which is ready to be cut and stapled. Ensure that the closing stroke is completed. All fingers should be completely removed from the closing trigger to ensure that the closing system holds.

NOTE: If the tissue needs to be repositioned within the instrument before stapling, open the jaws by squeezing the handle and closure trigger together slightly, push the release button, and slowly release grasp of the closure trigger. The closure trigger will return to the fully open position and the jaws will release the tissue but the retaining pin will remain extended. Tissue can now be repositioned.

NOTE: The CONTOUR Curved Cutter Stapler is designed with a feature which allows partial closing if a used reload, no reload, or an incorrectly inserted reload is in the instrument. In these cases, the device will close approximately one-third of the stroke and stop, but will return to the full open position when the closure

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
trigger is released. The used or misloaded cartridge module should then be replaced with a new reload or, if no cartridge was present, a cartridge should be loaded in the instrument. After following the above steps, any device which does not close either partially or completely should not be used.

NOTE: After closing, pausing for 15 seconds prior to firing may result in better tissue compression and better staple form.

NOTE: If closure trigger is difficult to lock, reposition the instrument and take a smaller amount of tissue. (See Contraindications for appropriate reload selection).

6 With the handle and closure trigger grasped in the palm of the hand, place all fingers around the firing trigger. Before firing, check to ensure the retaining pin is seated in the anvil.

Caution: Always ensure complete closing prior to firing.

Caution: Before firing, verify that the intended tissue is in the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.

Caution: When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

Caution: Care should be taken not to fire this instrument when the resulting staple line will be in tension.

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
7 Completely fire the instrument by squeezing the firing trigger until it touches the closure trigger (plastic to plastic), signifying that the instrument has fired the staples and knife blade, and the tissue has been transected. (Illustration 5)


The firing trigger automatically returns to the intermediate position as soon as the hand releases it, leaving the closure trigger in the fully closed position. Make sure, that the release button is not pressed during firing as proper formation of staples and cut line may be compromised.

Caution: Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result.

Caution: The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.

8 After firing the instrument, open the jaws of the instrument by squeezing the closure trigger and pushing the release button with

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

the thumb. For controlled release, while holding the release button down, open your hand to release the closure trigger. (Illustration 6)

The cartridge drivers are exposed as a spent cartridge indicator. If the instrument becomes locked onto tissue and will not open by pressing the release button, the device can be opened by depressing the release button and pulling the closure trigger simultaneously.

Caution: Examine the staple line for hemostasis and proper staple formation. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.


Caution: Following a division of the rectum, it is recommended that a leak test is performed prior to the creation of the anastomosis.

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
- 9 Retract the retaining pin manually. Ensure that the staple line and the cut line are complete. Remove the instrument. (Illustration 7)

Reloading the CONTOUR Curved Cutter Stapler

- 1 Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.
- 2 Depress the release button to ensure the instrument is in the open position. Ensure that the retaining pin has been retracted.
- 3 Remove the used reload from the instrument. Grasp the top of the reload and lift upward, unsnapping the reload from the jaws. (Illustration 8)

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

Properly discard the used reload. Thoroughly clean the instrument by vertically submerging completely the anvil and cartridge channel in a sterile solution and swish vigorously. After swishing, visually inspect the anvil and cartridge channel and remove all residual staples and/or foreign matter from instrument prior to reloading.

Caution: Failure to properly clean the instrument prior to reloading could damage the instrument; compromise cut line and staple formation of subsequent firings.

4 Examine the new reload for the presence of a staple retainer. If the retainer is not in place, discard the reload. Do not remove the staple retainer until the reload cartridge is inserted into the device.


NOTE: The staple retainer cannot be removed until the cartridge is correctly loaded into the device.

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
The loaded device cannot be used until the staple retainer is removed.

5 Insert the new reload into the metal housing and snap into position. (Illustration 9)


9

The tracks on each side of the reload should be used as guides to align the reload within the jaws of the instrument.


When the reload is properly aligned, push the reload into the instrument until it is fully seated. (Illustration 10)

10

Remove the staple retainer. (Illustration 11)

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
11

Check that the reload is held firmly within the jaws.

NOTE: If the cartridge is removed from the device, whether the cartridge is spent or not, and if the staple retainer has already been removed from the cartridge, the cartridge cannot be reloaded into the device.

6 After reloading and after removing the staple retainer, observe the surface of the new reload. If

drivers protrude out of the reload, replace with another reload.


7 The stapler is now reloaded and ready for use.

Muscle Stapling Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument. (Refer to **Warnings and Precautions**.)

1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.

2 Remove the staple retainer from the instrument by pushing on the pad labeled “Push to Remove”. (Illustration 2).

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


Hold the cartridge when removing the retainer. Check that the cartridge is fully loaded into the instrument (press down) after the retainer has been removed. Discard the staple retainer.

3 Grasp the instrument handle in the palm of the hand with the fingers around the closure trigger. Position tissue to be cut and stapled between the anvil and the cartridge.

NOTE: Any tissue covering the hole in the anvil will be pierced by the retaining pin. Reposition the tissue if necessary.

4 Recommended operation, in skeletal muscle, is to manually advance the Push Rod using the Retaining Pin Actuator. Squeeze the closure trigger (Illustration 3B).

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
Note that there is an audible click halfway through the closure stroke indicating that the device is in the intermediate position. In the intermediate position, the pin is fully seated in the anvil capturing the tissue, and the jaws are partially open. Reposition tissue within the instrument if desired. The closure trigger can be released at this point, and the instrument will maintain its jaw opening to allow for final positioning of the tissue to be cut and stapled. The retaining pin is fully engaged.

Caution: Make sure tissue to be stapled is properly positioned in the jaws before closing the jaws or stapling. Bunching, tension, stretching or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.

NOTE: In thick skeletal muscle tissue, two approaches are recommended for device closure:

- A single slow, controlled and continuous squeeze of the Closure Trigger
- Progressively squeezing and relaxing (pumping) the Closure Trigger

NOTE: Do not grasp the firing trigger before the instrument is to be fired.

5 Squeeze the closure trigger and handle together until the closure trigger is latched. Listen for an audible click. When the closure trigger is latched, the firing trigger moves to the ready-to-fire position (Illustration 4).

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
The cartridge is clamped onto the tissue, which is ready to be stapled and cut. Ensure that the closing stroke is completed. All fingers should be completely removed from the closing trigger to ensure that the closing system holds.

NOTE: If the closure trigger is difficult to lock, reposition the instrument and take a smaller amount of tissue.

NOTE: Do not use excessive force to close the instrument. Forcing the instrument to fully close, with excess tissue, may result in the closure trigger breaking.

NOTE: If the tissue needs to be repositioned within the instrument before stapling, open the jaws by queezing the handle and closure trigger together slightly, push the release button, and slowly release grasp of the closure trigger. The closure trigger will return to the fully open position and the jaws will release the tissue but the retaining pin will remain extended. Tissue can now be repositioned.


NOTE: The muscle stapling CONTOUR Curved Cutter Stapler is designed with a feature that allows partial closing if a used cartridge reload, no cartridge reload, or an incorrectly inserted

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cartridge reload is in the instrument. In these cases, the device will close approximately one-third of the stroke and stop, but will return to the full open position when the closure trigger is released. The used or misloaded cartridge module should then be replaced with a new cartridge reload or, if no cartridge was present, a cartridge should be loaded in the instrument. After following the above steps, any device that does not close either partially or completely should not be used.

NOTE: After closing, pausing for 15 seconds prior to firing may result in better tissue compression and better staple form.

6 With the handle and closure trigger grasped in the palm of the hand, place all fingers around the firing trigger. Before firing, check to ensure the retaining pin is seated in the anvil.

Caution: Always ensure complete jaw closure prior to firing.


Caution: Before firing, verify that the intended tissue is in the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.

Caution: When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

7 Completely fire the instrument by squeezing the firing trigger until it touches the closure trigger (plastic to plastic), signifying that the instrument has fired the staples and knife blade, and the tissue has been transected. (Illustration 5)

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


The firing trigger automatically returns to the intermediate position as soon as the hand releases it, leaving the closure trigger in the fully closed position. Make sure that the release button is not pressed during firing as proper formation of staples and cut line may be compromised.

Caution: Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result.

Caution: The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.


8 After firing the instrument, open the jaws of the instrument by squeezing the closure trigger and pushing the release button with the thumb. For controlled release, while holding the release button down, open your hand to release the closure trigger. (Illustration 6)

The cartridge drivers are exposed as a spent cartridge indicator. If the instrument becomes locked onto tissue and will not open by pressing the release button, the device can be opened by depressing the release button and pulling the closure trigger simultaneously.

Caution: Examine the staple line for hemostasis and proper staple formation. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

9. Retract the retaining pin manually. Ensure that the staple line and the cut line are complete. Remove the instrument. (Illustration 7)

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7.3. Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the CONTOUR® Curved Cutter Staplers implantable staple. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) per staple extends approximately 3-mm relative to the size and shape of the implantable staple.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular
Signal Void Size (mm ²)	17	5	42	24

7.4. The possibility and ways of integrating with other medical devices

Can be used in conjunction with MD «Cassettes with staples for stapling and cutting devices» (FSZ 2011/08985 dated 23.12.2019) (codes: CR40B, CR40G)

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10. LIST AND DESCRIPTION OF MD MATERIALS

Name	Material Marks	Type of contact
CONTOUR Curved Cutter Stapler		
Firing Bar	Stainless Steel	Limited Tissue Contact (<24 hours)
Knife Retraction Arm		
Firing Trigger		
Lockout Lever		
Anvil Plate		
Lockout Spring	Stainless Steel	
Trigger Spring		
Detent Toggle	Polycarbonate	
Release Button		
Left Shroud, Right Shroud	Polycarbonate	
Rod Guide	Polycarbonate	
Closure	Nylon	
Head Rivet	Stainless Steel	
Hook Rivet	Stainless Steel	
Retaining Pin		
Knife Guide Pin		
Push Rod		
Middle Rivet	Stainless Steel	
Washer	High Density Polyethylene	
Washer Mounting Pin	Polycarbonate	
Knife	Stainless Steel	
Knife Holder	Liquid Crystal Polymer	
Head	Aluminum alloy	
Firing Trigger Shroud	Polycarbonate	
Retaining Pin Actuator		
Closure Trigger	Polyurethane	

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While there are a number of patient and surgical factors that affect wound healing (e.g., tissue type, comorbidities, etc.), general wound healing is essentially complete by three months post treatment (Singh et al., 2017). As tissue continues to heal during this time period, the staple therapeutic function for approximating tissues post-transection / resection becomes increasingly managed by the tissue remodeling. As a result, three months is considered the therapeutic lifetime of the subject devices. Safety and performance of the devices will be demonstrated through three months post treatment. This length of time includes the wound healing period (i.e., the therapeutic lifetime) and potential later consequences associated with the healing process (e.g., late fistula or stricture) ((Reinke and Sorg, 2012; Velnar et al., 2009)).

Post Therapeutic Lifetime: Long-Term Implant

Once the staple has completed its therapeutic lifetime during the tissue healing period, the device remains in the patient for the lifetime of the patient or may be expelled after healing during normal cellular renewal (particularly in the gastrointestinal tract) or, if the physician deems it necessary, the implant may be removed either intra-operatively or post-operatively.

Finally, permanent surgical titanium alloy staples are well-established technologies with extensive history of safety over the lifetime of the patient.

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13. HOW SUPPLIED

The Curved Cutter Stapler CONTOUR are supplied sterile for single patient use. Discard after use.

Completeness of one primary packaging of the Curved Cutter Stapler CONTOUR with Reload :

- Curved Cutter Stapler CONTOUR - 1 pc.
- Reload with staples - 1 pc.

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thickness to exceed the indicated range for the selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.

- Do not fire device without a loaded cartridge.
- Be sure to remove the staple retainer prior to inserting tissue and closing the device.
- If the staple retainer is not in place in a cartridge module, discard the cartridge module and replace with a new one. Do not remove the staple retainer from the cartridge until the reload cartridge has been inserted into the instrument.
- Before closing the instrument on tissue, observe the surface of the new reload where the staples will emerge. If colored drivers protrude out of the reload, replace with another reload.
- Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.
- Ensure that unintentional tissue is not cut.
- Always ensure complete closing prior to firing.
- Before firing, verify that the intended tissue is in the closed jaws of the instrument. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.
- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc., are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.
- Make sure that the release button is not pressed during firing as proper formation of staples and cut line may be compromised.

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- Care should be taken not to fire this instrument when the resulting staple line will be in tension.

- Before firing, check to ensure the retaining pin is seated in the anvil. If the pin is not properly positioned, staples may not form properly, which may result in leakage or disruption of the staple line.

- Do not fire the instrument unless the closure trigger is properly latched against the handle.

- The firing trigger must be pulled back completely against the closure trigger to properly fire the instrument. Do not partially fire the instrument.

- Attempting to force the trigger to complete the closing stroke with too much tissue or thickened tissue may result in poor staple line integrity with possible leakage or disruption. In addition, instrument damage or failure may result.

- The firing stroke must be completed. Do not partially fire the instrument. Incomplete firing can result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device.

- After firing and release and before removing the instrument, check that the cut is complete.

- Retract the retaining pin manually prior to removing the device.

- Examine the staple lines for hemostasis and complete closure. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

- Following a division of the rectum, it is recommended that a leak test is performed prior to the creation of the anastomosis.

- The CONTOUR Curved Cutter Stapler may be reloaded during a single procedure. Do not reload the instrument more than five times for a total of six firings per instrument.

- Failure to properly clean the instrument prior to reloading could damage the instrument: compromise cut line and staple formation of subsequent firings.

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• Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.

- Dispose of all opened instruments whether used or unused.
- This device is packaged and sterilized for single use only. Multiple patient use may compromise the device integrity or create a risk of contamination that, in turn, may result in patient injury or illness.

MR Conditional

Non-clinical testing has demonstrated the implantable staple made of titanium (Ti3Al2.5V) alloy in these devices is MR Conditional. A patient with the implanted staples can be scanned safely immediately after placement of these staples, under the following conditions:

- Static magnetic field of 3.0 Tesla or less
- Spatial gradient field of 720 Gauss/cm.
- Maximum whole body averaged specific absorption rate (SAR) of 2.7 W/kg for 15 minutes of scanning.

MRI Related Heating

In non-clinical testing, the implantable staple made of titanium (Ti3Al2.5V) alloy in the CONTOUR® Curved Cutter Staplers produced a temperature rise of less than 2°C at a maximum whole body averaged specific absorption rate (SAR) of 2.7 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a (3-Tesla/128-MHz, Excite, Software HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR scanner.

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20. MD DISPOSAL RULES AND CONDITIONS

Dispose of the medical device in accordance with national legal and biohazard medical device requirements. Disposal of products must be carried out in accordance with SanPiN 2.1.3684-21 and must be disposed of properly.

Medical waste class in accordance with Russian SanPiN:
 For unused device – class A
 For used device – class B

21. LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Not applicable. This medical device does not contain materials of animal and / or human origin.

22. LIST OF APPLICABLE STANDARDS	
Standard	Title
<i>Management Systems</i>	
EN ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 14971	Medical devices - Application of risk management to medical devices
<i>Labeling</i>	
EN ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements
EN 1041	Information provided with medical device
<i>Design</i>	
EN ISO 16061	Instrumentation for use in association with non-active surgical implants

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