
TERUMOBCT

Instruction for use

Blood bag set Reveos LR

February 2019

Date: 20 February 2019

Signature: 

Position: Manager, Regulatory Affairs

Name: Karen Lowitz

State of Colorado

County of Jefferson

Subscribed and affirmed before me

this 20th day of February 2019

by Karen Lowitz

Notary Public

JENNIFER KESTER
NOTARY PUBLIC
STATE OF COLORADO
NOTARY ID # 20044014824
My Commission Expires 05-01-2020


Intended Use and Description

- The Reveos LR Set* is intended to collect a unit of whole blood and to process the whole blood unit on the Reveos Automated Blood Processing System, producing blood components.
- The Reveos LR Set has been evaluated for use only with the Reveos system.
- The blood and fluid pathways of the blood bag set are steam-sterilized and are non-pyrogenic.

CPD Anticoagulant Description Each 100 mL of CPD (citrate phosphate dextrose) anticoagulant contains: • Acidum citricum (anhydr.) 0.299 g • Natrii citras (dihydr.) 2.63 g • Mononatrii phosphas (monohydr.) 0.222 g • Dextrosum (monohydr.) 2.55 g • Aqua ad iniectionabilia q.s. 100 mL	SAG-M Additive Solution Description Each 100 mL of SAG-M (saline, adenine, glucose, mannitol) additive solution contains: • Natrii chloridum 0.877 g • Adeninum 0.0169 g • Dextrosum (anhydr.) 0.818 g • Mannitolum 0.525 g • Aqua ad iniectionabilia q.s. 100 mL
---	---

*** Please see the Annex for the Instruction for Use**

Set Diagram

- 1 Whole blood bag
- 2 Interim platelet unit (IPU) bag
- 3 Plasma bag
- 4 Red blood cell (RBC) bag
- 5 Residual leukocyte bag

Reveos LR

6 Sample/diversion bag

7 Needle

8 Needle injury protector (NIP) 9 Sample tube holder/luer adapter

10 Cross-connector

11 RBC leukoreduction filter

A, C: Y connector

B, G: Blue clamp

D, F, H, J: CLIKTIP frangible connector E: White clamp

Reveos LR Storage Conditions

- Long-term storage temperature: 1 °C to 30 °C
- Permitted temperature excursions:
 - -20 °C to 1 °C for up to 2 weeks
 - 30 °C to 50 °C for up to 1 week
- Storage humidity: 35% to 75% relative humidity (RH), non-condensing

Use blood bag sets within 28 days after you open the outer aluminum foil pouch. To store unused blood bag sets, return them to the outer aluminum foil pouch and re-close the pouch with tape or a clip. Once you open the transparent packing wrap, you must use the blood bag set within 7 days, not exceeding 28 days from when you opened the outer aluminum foil pouch.

Warnings

1. Do Not Reuse/Not for Reuse: Terumo BCT Ltd. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
 - Fluid leaks
 - Parts that are warped or deformed
 - Plastics that are brittle and discolored
 - Filters that have reduced filtration capabilities
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

2. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to

DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.

3. Inspect the packaging and the blood bag set prior to use. Do not use the set if any of the following conditions are present:
 - There are tears or holes in the outer aluminum foil packaging or in the individual transparent packaging wrap.
 - The tubing has severe kinks.
 - The blood bag set is incorrectly assembled.
 - The blood bag set is defective or damaged, or there are any leaks from the fluid-filled components of the set.
 - Any clamps are closed.
 - The needle cap is not in place.
 - The solutions are cloudy or discolored or contain particulates.

Using the blood bag set under these conditions may result in product contamination or poor performance during collection and/or processing.

Note: It is normal to have some condensation in the outer aluminum foil pouch and individual transparent packaging wrap due to sterilization.

4. The blood bag set is no longer sterile if any of the following conditions occur:
 - You disconnect the sample/diversion bag before you seal the sample/diversion line.
 - You remove blood samples before you seal the sample/diversion line.
 - The integrity of the set is compromised for any reason.

Manage the blood bag set according to your institution's standard operating procedure (SOP).

Cautions

1. Do not process whole blood less than 2 hours after collection. Processing blood too soon after collection may result in incomplete RBC leukoreduction and/or a reduced platelet yield.
2. Residual leukocytes are a by-product of the processing procedure and contain mostly white blood cells (WBC), with some plasma, platelets, and red blood cells (RBC). Residual leukocytes are not intended for transfusion.
3. Use aseptic technique during blood collection to ensure donor safety and product quality.
4. Do not vent the blood bag set.
5. Precautions for the oxygen absorber sachet:
 - Each aluminum foil pack contains a sachet that absorbs oxygen and generates heat on contact with air.
 - Handle the oxygen absorber sachet with caution.
 - Do not open the oxygen absorber sachet.
 - Do not dispose of the oxygen absorber sachet together with waste containing volatile or flammable material.

Procedures for Use

This document provides blood collection and post-processing instructions specific to the Reveos LR set. For warnings, cautions, and instructions on processing whole blood with the Reveos system, see the Reveos system operator's manual.

Blood Collection

1. If it is part of your institution's standard operating procedure (SOP), make a loose knot in the collection tubing between white clamp **E** and the whole blood bag.
2. Separate the sample/diversion bag from the other bags.
3. Load the bags onto a scale and/or blood mixing device according to your institution's SOP. Ensure that the whole blood bag and the sample bag are lower than the donor's arm.
4. Close blue clamp **B**.
5. Prepare the venipuncture site and perform the venipuncture according to your institution's SOP.
6. Open blue clamp **B**.
7. Position the sample/diversion bag with blue clamp **B** at the top and the sample tube holder/luer adapter at the bottom (see Figure 1). Allow the desired volume of whole blood to flow into the sample/diversion bag. The nominal volume of the sample bag is 60 mL. When the volume reaches the marking, approximately 40 mL has been collected.
8. Close blue clamp **B**.
9. Break CLIKTIP **D** to allow the whole blood to flow into the whole blood bag. Bend the CLIKTIP in both directions to ensure that you break it completely.
10. Seal the sample/diversion tubing as near as possible to Y connector **A**. Seal the tubing according to your institution's SOP.
11. Invert the sample/diversion bag and transfer donor blood samples from the sample/diversion bag, using evacuated blood collection tubes. Transfer samples as soon as possible after venipuncture to avoid possible clot formation in the sample/diversion bag.
 - a. Hold the sample/diversion bag and sample tube holder/luer adapter in one hand with blue clamp **B** positioned at the bottom (see Figure 2).
 - b. Using the other hand, insert a blood collection tube firmly into the sample tube holder/luer adapter. After filling the blood collection tube, remove the tube from the sample tube holder. Repeat the process to take additional samples.
12. Mix the whole blood and the anticoagulant during collection according to your institution's SOP.
13. Collect the target volume of blood $\pm 10\%$, as specified on the whole blood bag label.
14. Close white clamp **E** and perform one of the following steps:
 - a. If you made a knot in the collection tubing in step 1, tighten the knot firmly.
 - b. If you did not make a knot in step 1, seal the collection tubing near CLIKTIP **D**.
15. Slide the NIP as close to the needle hub as possible.
16. With one hand, grip the needle hub and the front of the NIP. Pull gently to withdraw the needle from the donor's arm (see Figure 3).
17. Immediately after you withdraw the needle from the donor's arm, gently pull on the tubing behind the NIP to cover the needle with the NIP. Ensure that the needle locks securely in place (see Figure 4).
18. Seal the collection tubing near CLIKTIP **D** if you did not do so in step 14. Disconnect the sealed tubing that includes the needle and the sampling assembly. Dispose of the needle, the NIP, and the sample tube holder safely according to your institution's SOP and/or local regulations.

19. Immediately after collection is complete, invert the whole blood bag several times to thoroughly mix the whole blood and the anticoagulant.
20. If required, strip the blood in the collection tubing into the whole blood bag according to your institution's SOP.
21. Seal and disconnect the collection tubing between 0.5 cm and 3.0 cm away from the whole blood bag. Dispose of the collection tubing safely according to your institution's SOP and/or local regulations.
22. Place the whole blood unit into a temperature-controlled environment according to your institution's SOP.
23. Pack and transport the whole blood unit to the processing laboratory according to your institution's SOP.

Figure 1: Sample/diversion bag


Figure 2: Collecting a whole blood sample


Figure 3: Withdrawing the needle


Figure 4: Covering the needle with the NIP

Post-Processing and Red Blood Cell Leukoreduction

1. After processing, remove the plasma bag, the residual leukocyte bag, and the IPU bag from the organizer and handle them according to your institution's SOP.

Note: You must rest and agitate the IPU prior to platelet pooling. For conditions for resting and agitating the IPU, see the platelet pooling set instructions for use.

2. Prepare a leukoreduction rack with a head height of at least 110 cm.
3. Remove the RBC bag from the organizer.
4. Hang the RBC bag from the leukoreduction rack so that the filter hangs vertically.
5. Open blue clamp **G**.
6. Break CLIKTIP **H**.
7. Break CLIKTIP **J**.
8. Drain all of the additive solution through the filter and into the whole blood bag.
9. Close blue clamp **G**.
10. Gently mix the additive solution with the RBC unit in the whole blood bag.
11. Hang the whole blood bag on the leukoreduction rack.
12. Ensure that the tubing is free of kinks or other obstructions.
13. Open blue clamp and allow the mixture of red blood cells and additive solution to flow through the leukoreduction filter.

Do not squeeze the whole blood bag to increase the draining rate. Filtration is complete when the inlet side of the leukoreduction filter collapses and the red blood cells drain from the inlet side of the filter.

14. Seal the RBC tubing below the leukoreduction filter, and disconnect the RBC bag.
Ensure that the filter remains in a vertical position while you seal the tubing, as this prevents the loss of red blood cells.
15. Dispose of the whole blood bag and filter assembly safely according to your institution's SOP and/or local regulations.
16. Seal the line for sampling segments according to your institution's SOP.
17. Store the RBC product according to your institution's SOP.

Symbols and Certification

	Number of products, if in the box indicates the number.		Indicates the date after which the medical device is not to be used.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.		Indicates the manufacturer's catalog number so that the medical device can be identified.
	Indicates a non-pyrogenic fluid pathway.		Indicates the presence of a sterile fluid path. The method of sterilization is steam.
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.		Indicates a medical device that is not to be resterilized.
	Indicates the medical device manufacturer		Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile).
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.		Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).
	Indicates a medical device that needs protection from light sources.		Indicates a medical device that needs to be protected from moisture.
	Indicates the range of temperatures to which the device can be safely exposed.		Indicates that the contents of the transport package are fragile; therefore, it shall be handled with care.
	Indicates the correct upright position of the transport package.		Indicates that the transport package is not allowed to be vertically stacked beyond the specified number.

	Indicates that the user should consult the instructions for use		Indicates that the product should not be used if the product, its sterile barrier system, or its packaging is damaged or shows any signs of deterioration.
	Indicates that the product should not be used if the package is damaged		Indicates that the product should not be vented. Do not remove the end caps from the set.
	Indicates a donor needle.		Indicates the length of time within which the product must be used after the packaging is opened.
	Indicates the nominal volume of the blood container		Indicates a whole blood bag.
	Indicates the target whole blood collection volume.		Indicates a red blood cell bag.
	Indicates a plasma bag.		Indicates a processing container.
	Indicates anticoagulant solution.		Indicates additive solution.
	Indicates a sampling site.		Indicates a leukoreduction filter.
	Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended.		
CPD	CPD-63 ml 100 ml-CPD Acidum citricum (anhydr.) 0.299 g Natrii citras (dihydr.) 2.63 g Mononatrii phosphas (Monohydr.) 0.222 g Dextrosum (monogydr.) 2.55 g Aqua ad ingiectabilia q.s. 100 ml	SAG-M	SAG-M-100 ml 100 ml – SAG-M Natrii chloridum 0.877 g Adeninum 0.0169 g Dextrosum (anhydr.) 0.818 g Mannitolium 0.525 g Aqua ad ingiectabilia q.s. 100 ml

Return of Used Product

If for any reason this product must be returned to Terumo BCT Ltd., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department. IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT. Please contact your local representative for information regarding returned goods and product complaints.

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* Further in the text Reveos LR.

FURNISHING OPTIONS OF THE MEDICAL DEVICE

There's only 1 furnishing option of the Reveos LR

INFORMATION ON MANUFACTURERS OF THE MEDICAL DEVICE

Legal Manufacturer

Terumo BCT Ltd.
Old Belfast Road
Millbrook, Larne,
BT40 2SH
United Kingdom

Developer

Terumo BCT, Inc.
10811 W. Collins Ave.
Lakewood, Colorado 80215
USA
Tel: +1 303 231 4357
Fax: +1 303 542 5215

Manufacturing Facility

Reveos LR

Terumo BCT Vietnam Co., Ltd.
Long Duc Industrial Park
Long Duc Commune
Long Thanh District
Dong Nai Province
Vietnam

MEDICAL DEVICE CLASSIFICATION

The Reveos LR is a CE-marked class IIb medical device complying with the Medical Device Directive (93/42/EEC) as amended.

GMDN Code: 46347 - Blood Donor Set – Quad Pack

Class of electrical safety, laser safety and other types of danger

NOT APPLICABLE

Type of contact with a human body

There is an invasive contact with the human body by means of a venipuncture (via the needle of the blood bag set) to collect whole blood. This is a transient contact (less than 60 minutes) to channel blood into the blood bag.

INTENDED USE

The **Reveos LR** is intended to collect a unit of whole blood and process the unit on the Reveos device, producing blood components.

Scope

Taking, processing of blood into components, component storage and transfusion after processing.

FUNCTIONAL CHARACTERISTICS

Indications

The Reveos LR Set is intended to collect a unit of whole blood and to process the whole blood unit on the Reveos Automated Blood Processing System, producing blood components

Contraindications

There are no known contraindications.

Conditions of use

This product can be used by qualified medical personnel for the collection of whole blood (i.e. phlebotomist, nurse, doctor). After whole blood collection, trained blood center personnel will process the using the Reveos Automated Blood Processing System.

Whole blood collection using the Reveos LR can occur at different locations as per the customer's operating procedures. The collected whole blood set should be stored and transported under temperature controlled conditions as per the customer's standard operating procedures. Processing of the collected whole blood contained in the Reveos LR occurs in a laboratory environment/setting.

Reveos LR

The blood bag set is no longer sterile if any of the following conditions occur:

- You disconnect the sample/diversion bag before you seal the sample/diversion line.
- You remove blood samples before you seal the sample/diversion line.
- The integrity of the set is compromised for any reason.

Manage the blood bag set according to your institution's standard operating procedure (SOP).

PRINCIPLE OF MEDICAL DEVICE OPERATION

The **Reveos Automated Blood Processing System** is an automated blood component separation system that uses centrifugal force to separate one to four units of whole blood simultaneously into two (red blood cells and plasma) or three (red blood cells, plasma and platelets), components from each unit of whole blood.

General description of how the Reveos system works:

First, whole blood is collected from a donor into the whole blood bag set of a Reveos LR.

Next, the operator loads the whole blood bag set from the Reveos LR in the centrifuge bucket of the Reveos device. The other product component bags from the Reveos LR are placed, with the labels facing towards the operator, in the blood bag component holder of the Reveos device. Once in place, this allows for routing of the cross connector and the tubings for whole blood, platelets, plasma, and residual leucocytes through the defined channels and valves in the centrifuge loading assembly.

When less than four blood bags are processed at one time, the operator loads counterweight bags until all four centrifuge buckets are filled.

The operator can select either a 2-component procedure or a 3-component procedure using the Reveos System Manager or the touch screen display on the Reveos device. The operator will input all the information required by the blood bank's standard operating procedures before the separation process can begin.

A "2 component" procedure will result in the separation of whole blood into red blood cells and plasma and a residual product. A "3 component" procedure will result in the separation of whole blood into red blood cells, plasma and platelets.

Once the procedure has been set up correctly, the operator can start the processing run. The Reveos device uses standard centrifugation methods to separate the whole blood into distinct components of differing density (red cells, plasma, leucocytes, and platelets).

There are 3 main steps:

Step 1 CENTRIFUGATION

When the procedure starts, the rotor begins to spin. The whole blood will be separated into different component layers based on the density by centrifugal force.

Step 2 EXPRESSION

After the whole blood has been centrifuged and separated into components, the hydraulic bladder fills with hydraulic fluid to express the components into their respective product bags.

Sensors in the centrifuge bucket and in the centrifuge loading assembly at the top of the rotor identify the interfaces between the various components. Valves will open and close at appropriate times during expression to route the blood components to their respective product bags depending on the protocol.

This step is repeated for each of the other loaded whole blood units.

Step 3 SEALING

When expression is complete, the valves will seal the lines to the product bags.

When the processing is complete, the outer lid of the Reveos device opens automatically.

The operator removes the product bags from the Reveos device.

Specifications for the integrated solutions

The following solutions are NOT considered medicinal substances in the current intended use for collection and preservation of whole blood and blood components. The amount of solutions are calculated strictly for preservation of blood or blood components.

➤ SAG-M Solution

This solution is a red cell additive solution for preservation of red blood cells during storage.

The red blood cell bag holds 100mL of Saline Adenine Glucose – Mannitol (SAG-M) solution.

Master Formula - Batch Contents to 1000 liters of SAG-M Solution:

Description	Mass (kg)
Adenine PhEur	0.169
Mannitol PhEur	5.25
Sodium Chloride PhEur	8.77
Dextrose Anhydrous PhEur	8.18
Water for Injection PhEur	985

➤ CPD Solution

CPD solution is an anticoagulant solution intended to prevent coagulation of whole blood and blood components in the disposable set. This solution is NOT intended for direct infusion to a donor or patient.

The whole blood collection bag holds 63mL of Citrate Phosphate Dextrose (CPD) solution.

Master Formula - Batch Contents to 1000 liters of CPD Solution:

Description	Mass (kg)
Sodium Dihydrogen Phosphate 2H ₂ O PhEur	2.51
Citric Acid Monohydrate PhEur	3.27
Dextrose Monohydrate PhEur	25.5
Sodium Citrate PhEur	26.3
Water for Injection PhEur	972

Options of executions

Device has a single execution option

Information of marking and packaging

Products supplied in packing marked with the following information:

- device name and reference number;
- LOT of the product;
- logo and address of the manufacturer;

Reveos LR

- Information about sterile;
- Number of products;
- expiry date;
- and other information and symbols.

Project of the labeling in Russian language

Project of the labeling in Russian language for the foil packaging

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Project of the labeling in Russian language for the transportation packaging

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Packaging

Each shipping carton (or case) contains 18 Reveos LR sets. Each foil pouch contains 2 individually packaged products. The Reveos LR foil pack contains two oxygen absorber packs. The shipping carton also contains the instructions for use.

Materials of which the medical device or its main components are made

Part of Device	Part name	Material /Brand of Material	Manufacturer
Игла (Needle)	Cannula	Stainless Steel (SS 304)	Terumo (Japan)
	Needle hub	PVC (AB-6)	
	Needle cap	Blue Polypropylene (PP MG02T including dye), PVC (PVC A101)	
Колпачок иглы защитный (needle injury protector)	---	Polypropylene (TitanPro 6331)	Terumo Penpol, Ltd (India)
Мешок для цельной крови (Whole blood bag); Мешок дозы промежуточной тромбомассы (IPU) (Interim platelet unit bag); Мешок для плазмы (Plasma bag); Мешок для эритроцитов (Red blood cell bag); Мешок для остаточных лейкоцитов (Residual leukocyte bag); Мешок для проб/отвода крови (sample/diversion bag);	Bag Vinyl	PVC (S-502TW9)	Terumo BCT,(Vietnam)
	Tubing Connectors	PVC (A-101)	
	Outlet Ports	PVC (A-101TW)	
Держатель пробирки для пробы/адаптер Люэра (Sample tube holder/luer adapter)	Holder and cap	Polypropylene (MG02T)	Terumo BCT Ltd, (Vietnam)
	Sample Needle:	Stainless steel (SUS 304)	Terumo Corporation (Kofu Factory)
	Sheath:	Isoprene rubber (A1G3210M)	Terumo Corporation (manufactured by Kokoku Intec Co., Ltd.) (Japan)
	Hub	Polycarbonate (S3000R)	Terumo Corporation (Kofu Factory) (Japan)
	Hub-Needle Adhesive	Urethane acrylate (UV light) Loctite 3311	Henkel, Ltd., (Japan)
	Tube Sleeve	PVC (T-640TW)	Terumo BCT Ltd, (Vietnam)
Фильтр лейкоцитарный для эритроцитов	---	PVC (Renolit 3226), PVC (AM3 W17) and	Fresenius HemoCare Italia Srl (Italy)

Part of Device	Part name	Material /Brand of Material	Manufacturer
(RBC leukoreduction filter) Bioflex (Bioflex RCC filter)		polybutylterephthalate (P22D0050)	
Соединитель хрупкий (CLIKTIP frangible connector)	Connector	PVC (IJ-500TW)	Vina Showa (Vietnam)
	Insert	Polycarbonate (Makrolon 2858)	
Зажимы трубок (Tubing Clamps)	Blue clamp	Delrin Poly Acetal Blue	Terumo Penpol, Ltd (India)
	White Clamp	Delrin Poly Acetal White	Terumo Penpol, Ltd. (India)
Трубки (Tubing) Numbered tubing	---	PVC (T-640TW) PVC (T-640TW) Printing by laser marking	Vina Showa (Vietnam)
Соединитель Y-образный (Y Connector)	---	PVC (BBM-1)* *Sample Line and donor Line "different design"	Terumo Penpol, Ltd. (India)
Перекрестный соединитель (Cross Connector)	---	Alpha Gary PVC 2212 C28/H-118 clear	Harmac (Ireland) or SMC (United States)
Организатор (Organizer)	---	Polypropylene (Pro-Fax 6523)	SMC (United States)
Tape for closing box	Adhesive Tape	Scotch® Box Sealing BOPP (Tape 375)	Galup CO, Ltd (Vietnam)
Box	brown cardboard, double wall 5 ply corrugated carton	P-230/G-175/S-150/G-175/P-230, Flute Type BC	Yuen Foong Yu Paper Enterprise Co, Ltd. (Vietnam)
Oxygen absorber	Iron Powder: Diatomaceous Earth: Moisture: Activated Carbon: Packaging Material:	60.00 % 15.00% 15.00% 10.00% PET/PAP Oil proof/PE	Jiangsu ouk packaging technology Co, Ltd (Japan)
Aluminum foil	Layer 1: Polyester	PET12/AL7/PE9	Southeast Packaging Industry Co, Ltd (Thailand)
	Layer 2: Aluminum foil	12/CPP55	
	Layer 3: Polyester	PET12/AL7/PE9	
	Layer 4: white polyethylene	RIBLENE FF30 white RAL 9016, LG CHEM	
Transparent packaging wrap	Biaxially orientated nylon	15µm BOPA (biaxially oriented polyamide film)/ 30µm CPP (cast polypropylene)	Vishakha Polyfab Pvt. Ltd (India)

Part of Device	Part name	Material /Brand of Material	Manufacturer
Bag labels	---	Thermal transfer ink ribbon for plastic (DYNIC DXR 101A 110xi 4)	Seljegat Singapore Pte Ltd. (India)

Rate Technical characteristics

General Specifications	
Bag Capacity	- Whole Blood bag: 450 mL $\pm$ 10% plus 63 mL CPD - RBC bag: Up to 600 mL, including 100 mL SAG-M storage solution - Plasma bag: Up to 600 mL if not frozen; Up to 300 ml if frozen; - Platelet bag: Up to 400 mL - Sample/diversion pouch: Up to 60 mL with marking indicating approximately 40 mL - Residual leukocytes bag: Up to 60 mL
Needle characteristics	• Outer diameter: 1.65 mm • Inner diameter: 1.39$\pm$0.03 mm • Length of the needle: – Length of the needle metal part: 33 mm, hub (connector / connector) - 30 mm and in total with the metal part-63 mm – assembled needle cap (the thickness of the cap)-7.5 mm, all 3 parts together -70.5 mm • Weight of the needle: 3.5 g • Needle shall withstand a pulling strength of 5.4 Kg • Needle shall not deform on steam sterilization at 120°C for 40 min • Wall thickness: 0.229 $\pm$ 0.0005 mm • Hardness up to 90 HRC
Needle cap	• Length 50 mm • Width 11.4 mm • Height 6.6 mm
Disposable Primary Packaging Dimension	▪ Disposable Tubing Sets with transparent packaging wrap and foil: 440 mm $\times$ 26 0 mm $\times$ 50 mm (L $\times$ W $\times$ H) ▪ Oxygen Absorber Pack (Each outer aluminum foil pouch contains a sachet that absorbs oxygen and generates heat on contact with air) Length: 50 $\pm$ 5 mm; Width: 40 $\pm$ 5 mm
Packaging dimensions (to be sterilized)	▪ Length 420 mm ▪ Width 210 mm ▪ Height 25 mm ▪ Weight Up to 0.372 kg
Whole blood bag dimensions (mm)	▪ Length: 180 ▪ Width: 125 ▪ Weight: 25.6 g
RBC bag dimensions (mm):	▪ Length: 180 ▪ Width: 125 ▪ Weight: 27.8 g

Reveos LR

Plasma bag (mm):	- Length: 180 - Width: 125 - Weight: 23.4 g
Platelet bag (mm):	- Length: 160 - Width: 125 - Weight: 18.6 g
Sample/diversion pouch (mm):	- Length: 85 (97 including welded vinyl) - Width: 64 (80 including welded vinyl) - Weight: 7.4 g
Residual leukocytes bag (mm):	- Length: 85 (97 including welded vinyl) - Width: 64 (80 including welded vinyl) - Weight: 7.4 g
Tensile properties (extension strength) of connections between bags and tubings	more than 80 N
Organizer	- Length 10.9 cm - Width 4.7 cm - Height 2.5 cm - Weight Up to 0.0067 kg
Tube	- Outer diameter of the tubes 4.4 ± 0.1 mm - Internal diameter 3.0 ± 0.1 mm - Thickness of the tube 1.4 mm (before sterilization)
Mass	- Assembled 8.20 kg - Products 372 grams - Products in the original packaging 376 grams - Tolerance for the presented values ± 10g
The hole diameter of the fitting on the bags	7,5 mm

The sachet of Oxygen Absorber Pack contains:

- Iron powder;
- Activated carbon;
- Water;
- Salt

Note: no contact with patient and medical staff.


Tubing specifications of Reveos LR

Item	Wall thickness, ± 0.1 mm	Outer diameter, ± 0.1 mm	Inner diameter, ± 0.1 mm
Access Needle Tubing	1.4 mm	4,4mm	3.0mm
Sample tube holder/luer adapter tubing (3)	2.6 mm	4,4mm	1.8mm
Tubing to Sample Bag	1.4mm	4,4mm	3.0mm
Tubing to Blood Collection Tube holder	1.4mm	4,4mm	3.0mm
Collection Line Tubing	1.4mm	4,4mm	3.0mm
Processing tubing	1.4mm	4,4mm	3.0mm
Tubing to Plasma Bag	1.4mm	4,4mm	3.0mm
Tubing to Residual Leukocytes Bag	1.4mm	4,4mm	3.0mm
Tubing to Platelet Bag	1.4mm	4,4mm	3.0mm

Item	Wall thickness, ± 0.1 mm	Outer diameter, ± 0.1 mm	Inner diameter, ± 0.1 mm
Tubing to RBC Bag	1.4mm	4,4mm	3.0mm
Tube to bag for whole blood	1.4mm	4,4mm	3.0mm
Tube of a line of a bag for Sample/diversion of a blood	1.4mm	4,4mm	3.0mm

*Note: The tubing dimensions listed on this document are nominal values and are intended for use in determining compatibility with various types of laboratory equipment. The nominal values are specified target dimensions; however, because they are based on non-sterile tubing measurements, and due to variations in the manufacturing process, the actual dimensions may be slightly different.

Tubing Pre- and Post Sterile dimensions

Position	Pre-Sterile Length		Post Sterile Length		
	Length	Tolerance	Tubing	Tolerance	Between Components
Access Needle Line					
1	260mm	± 10 mm	247 mm	± 5 %	235 mm
Primary Tubing					
2	220 mm	± 10 mm	209 mm	± 5 %	200 mm
5	1050 mm	± 20 mm	1000 mm	± 5 %	1016 mm
6	34 mm	± 2 mm	32 mm	± 5 %	19 mm
8	197 mm	± 5 mm	187 mm	± 5 %	181 mm
9	132 mm	± 3 mm	125 mm	± 5 %	117 mm

Reveos LR

10	132 mm	± 3 mm	125 mm	± 5 %	121 mm
11	132 mm	± 3 mm	125 mm	± 5 %	121 mm
12	756 mm	± 20 mm	750 mm	± 5 %	748 mm
13	294 mm	± 15 mm	279 mm	± 5 %	275 mm
Sampling Assembly Tubing					
3	30 mm	± 2 mm	28 mm	± 5 %	13 mm
4	65 mm	± 5 mm	61 mm	± 5 %	44 mm
Plasma Line					
7	405 mm	± 10 mm	390 mm	± 5 %	387 mm

Name	Length, mm	Width, mm	Thickness, mm	Weight, g
Needle injury protector	94	10.5	10.5	4.4
Sample tube holder/luer adapter	64.8	5.8 (at bottom)	5.8 (at bottom)	4.4
Cross-connector	20.7	20.7	11.7	1.1
RBC leukoreduction filter.	118.1	78.5	12	18.1
Y connector	At Needle: 33.2	24	6.5	1.3
	At sampling device: 27.7	13.1	6.7	1.2
Blue clamp	24.9	17.9	9.5	1.1
CLIKTIP frangible connector	57	8	8	1.8
White clamp	24.9	17.9	9.5	1.1

REQUIREMENTS TO ENVIRONMENT PROTECTION

Medical device does not influence on the environment.

Indication, if required, of special precaution measures to be taken when destroying unused medical devices (Utilization)

The Reveos LR Set should be disposed of in accordance with local regulations:

Class A of medical waste in accordance with SanPin 2.1.2790-10 for unused medical device

Class B of medical waste in accordance with SanPin 2.1.2790-10 for used medical device.

LIST OF HAZARDS RELATED OT THE USE OF THE MEDICAL DEVICE AND DESCRIPTION OF MEASURES/METHODS TAKEN/USED TO ENSURE ALLOWABILITY OF RESIDUAL RISKS

Terumo BCT performed a detailed Risk Analysis for the overall Reveos system including hardware, software and disposable sets to identify all potential failure modes and the associated hazards. This analysis has been performed in 2 phases. The purpose of the first phase was to

identify failure modes, causes and effects of the failure modes and to identify system and subsystem requirements to provide the required protection. The purpose of the second phase was to evaluate the effect of the detailed system and subsystem designs in alleviating the hazards. The residual hazards were estimated to assure that the system provides the desired level of protection

The purpose of the protection mechanisms is to either reduce the severity and/or likelihood (probability of occurrence or increase the detection capability of the system. The impact of the system design is indicated in the residual severity, the residual likelihood and residual detection. The process of estimating these residual severity, occurrence and detection scores were the same as was used to determine the initial scores.

Hazards for the Reveos System were identified; the risks were mitigated to the extent possible and mitigations were verified and/or validated. Residual risks have been identified and were found to be acceptable. Evaluation of the overall benefit versus the risks is found to be acceptable. Ongoing periodic review of production information and reports from the field will continue to be conducted and any impact on risk will be evaluated.

The Reveos LR Set is considered safe to use.

LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS.STANDARDS, THE MEDICAL DEVICE COMPLIES WITH

List of applicable standards used:

- EN ISO 62366:2015, Medical devices – Application of usability engineering to medical devices
- EN ISO 13485:2012, Medical Devices – Quality Management Systems – Requirements for Regulatory Purposes
- EN ISO 14971:2012, Application of Risk Management to Medical Devices
- EN ISO 3826-1:2013, Plastics collapsible containers for human blood and blood components Part 1: Conventional containers
- EN ISO 3826-3:2007, Plastics collapsible containers for human blood and blood components – Blood bag systems with integrated features
- EN ISO 1135-4:2015, section 5.4, Transfusion equipment for medical use – Part 4: Transfusion sets for single use
- EN ISO 15747:2011, Plastic containers for intravenous injections
- EN 15986:2011, Symbols for use in the Labeling of Medical Devices – Requirements for Labeling of Medical Devices containing Phthalates
- EN ISO 15223-1: 2016, Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied
- EN ISO 11607-1:2015, Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
- ISO 780:2015, Packaging – Distribution packaging – Graphical symbols for handling and storage of packages
- EN ISO 3826-2:2008, Plastics collapsible containers for human blood and blood components Part 2: Graphical symbols for use on labels and instruction leaflets
- EN ISO 1041: 2008 + A1 2013 – Information Supplied by the manufacturer with medical devices
- ISBT 128: 2016, Standard Technical Specification
- ISTA 2A: 2011, Procedure 2A – Packaged products weighing 150 lb (68 Kg) or less
- ISO 14644-1:2015 – Cleanrooms and Associated Controlled Environments – Part1: Classification of air cleanliness

- ISO 14644-2:2015 – Cleanrooms and Associated Controlled Environments – Part2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1
- EN ISO 556-1:2001 + AC:2006 – Sterilization of medical devices – requirements for medical devices to be designated “sterile” – Part 1: Requirements for Terminally sterilized medical devices.
- EN ISO 17665-1:2006 – Sterilization of Health Care Products – Moist heat – Part 1: Requirements for the development, validation and outline control of a Sterilization process for Medical Devices
- EN ISO 10993-1:2009, Overall Biological Evaluation of Medical Devices as part of Risk Management

SHELF LIFE

The product expires 2 years from date of manufacture. Refer to product labelling for expiry date.
The product keep sterility expires 2 years from date of manufacture

OPERATION CONDITIONS

The collection of whole blood using the Reveos LR typically occurs at room temperature. The processing of whole blood occurs in a laboratory setting using the Reveos Automated Blood Processing System under temperature controlled conditions.

Reveos Automated Blood Processing System Specification Type	Value
Humidity	10 to 80% RH non-condensing
Operating temperature	Recommended 18 °C to 27 °C (64.4 °F to 80.6 °F)
Altitude	Up to 2000m

The patient's body temperature is 32 °C – 42 °C, hence the operating conditions are:

Parameter	Value
Humidity	Up to 100%
Operating temperature	From 32 °C to 42 °C
Pressure	700 – 1060 gPa

TRANSPORT AND STORAGE CONDITIONS

Storage conditions

- Long-term storage temperature: 1 °C to 30 °C
- Storage humidity: 35% to 75% relative humidity (RH), non-condensing
- Protect from freezing. Avoid excessive heat.

Use blood bag sets within 28 days after you open the outer aluminum foil pouch. To store unused blood bag sets, return them to the outer aluminum foil pouch and re-close the pouch with tape or a clip. Once you open the transparent packing wrap, you must use the blood bag set within 7 days, not exceeding 28 days from when you opened the outer aluminum foil pouch.

Transportation conditions

Parameter	Value
Humidity	Permitted humidity excursions:

	- 75% to 85% ($\pm 5\%$) up to 72 hours - 15% to 35% for up to 6 months
Operating temperature	Permitted temperature excursions: - 20 °C to 1 °C for up to 2 weeks 30 °C to 50 °C for up to 1 week

Avoid freezing and high temperatures

WARRANTY LIABILITIES

Warranty shelf life: 2 years

1 When purchasing a Terumo BCT product, the general terms and conditions are applicable. This includes the following warranty statement:

2

LIMITED WARRANTY. Warranty information applicable to the equipment or disposables specified on this invoice is set forth in the price list and any contract for such equipment or disposables and shall be provided by Terumo BCT upon request. For a period of ninety (90) days from the date of the invoice, Terumo BCT will repair or replace Terumo BCT spare parts or optional feature components that Terumo BCT finds defective in material or workmanship. Terumo BCT is not responsible for, and shall have no obligation with respect to any failure caused by normal wear and tear, misuse, unauthorized alterations, accident, neglect, use of non-standard accessory attachments, and/or improper customer maintenance. Terumo BCT does not verify the safety or efficacy of, and makes no warranties, express or implied, with respect to any non-Terumo BCT components or non-Terumo BCT spare parts included at Purchaser's request or any components or products used other than as expressly intended by their manufacturer. Failure to perform preventive maintenance procedures will invalidate this warranty. Terumo BCT is not responsible for any events related to non-Terumo BCT parts installed in the equipment. THIS LIMITED WARRANTY CONTAINS THE PURCHASER'S EXCLUSIVE WARRANTIES AND REMEDIES.

CLAIM

3 When purchasing a Terumo BCT product, the general terms and conditions are applicable. This includes the following liability statement:

4 Limited liability: To the extent permitted under applicable law, Purchaser's exclusive remedy with respect to any claim, whether in contract, in tort, under statute, under warranty or otherwise shall be limited solely to the refund of the purchase price or replacement of all Products shown to be other than as warranted. Terumo BCT's liability for losses which do not affect the Product (e.g. indirect, incidental, consequential losses, damages caused by delay) shall be excluded. To the extent the foregoing limitation is not permissible under applicable law, the applicable legal limitation of liability under statute or otherwise if there shall be one, shall apply.

Representative office in Russia:

Terumo Russia LLC
123112, Moscow, Testovskaya str., 10,
Floor/premise/office 13/1/5
TEL +7-495-988-4740
FAX +7-495-988-4739

E-mail: Anna.Sabirova@terumobct.com

PRECAUTION MEASURES APPLIED WHEN USING THE MEDICAL DEVICE

Cautions

- Intended for single use only.
- Use aseptic techniques.
- The Disposable Kits described here contain phthalates (DEHP). Patient groups including pregnant or nursing women and children are considered to be the most at risk of harmful effects due to DEHP exposure. Based on five decades of successful use in medical devices without a reported adverse event related to DEHP, the risk of not doing a needed procedure is far greater than the risk associated with exposure to DEHP.

POSSIBLE SIDE EFFECTS WHEN USING THE MEDICAL DEVICE

Known adverse reactions related to blood collection

1. Local reactions related to needle insertion:
 - a. Vessel injuries: haematomas, arterial puncture, thrombophlebitis.
 - b. Nerve injuries: injury of nerve, injury of nerve by haematoma..
 - c. Other complications: tendon injury, allergic reaction (local), infection (local),
2. General reactions:
 - a. Vasovagal reaction: immediate type, delayed type.
3. Rare, significant complications :
 - a. Related to vessel injury: brachial artery pseudoaneurysm, arteriovenous fistula, compartment syndrome, axillary vein thrombosis.
 - b. Accidents: accidents or injuries related to vasovagal syncope, other kinds of accidents.
 - c. Cardiovascular reactions: angina pectoris, myocardial infarction, cerebral ischaemia.

Source: Guide to the preparation, use and quality assurance of blood components (Recommendation No. R (95) 15, 17th Edition) published by Directorate for the Quality of Medicines & HealthCare of the Council of Europe (EDQM)

INDICATION OF POSSIBILITY AND PECULIARITIES OF USING THE MEDICAL DEVICE BY PEOPLE WITH IMPLANTED MEDICAL DEVICES, PREGNANT WOMEN, BREASTFEEDING WOMEN, CHILDREN, ADULTS WITH CHRONIC DISEASES

The Disposable Kits described here contain phthalates (DEHP). Patient groups including pregnant or nursing women and children are considered to be the most at risk of harmful effects due to DEHP exposure. Based on five decades of successful use in medical devices without a reported adverse event related to DEHP, the risk of not doing a needed procedure is far greater than the risk associated with exposure to DEHP.

INFORMATION ON CAPABILITY OF THE MEDICAL DEVICE TO AFFECT THE ABILITY TO OPERATE VEHICLES, MECHANISMS

NOT APPLICABLE

Reveos LR

INDICATION, IF REQUIRED, OF SPECIAL PRECAUTION MEASURES TO BE TAKEN WHEN DESTROYING UNUSED MEDICAL DEVICES

The Reveos LR should be disposed of in accordance with local regulations.

According to hazard category of 2.1.7.2790-10 SanPiN:

Unused Reveos LR is Class A - epidemiologically safe waste, approximate in composition to solid household waste.

Used Reveos LR is Class B - epidemiologically hazardous waste.

REQUIREMENTS FOR MAINTENANCE AND REPAIR OF THE MEDICAL DEVICE

The Reveos LR is a single-use medical device and it does not need technical maintenance or repair.

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