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

Instructions for Use
DiamondTOUCH
Dated February 17, 2023
Consisting of 10 pages


ACKNOWLEDGEMENT

State of Utah)
 §
County of Salt Lake)

On this 17th day of February 2023, before me, Jennifer Fordham, a notary public, personally appeared Tom Haucter, Senior Manager, Regulatory Affairs proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand and official seal.


Notary Public

My Commission Expires: June 27, 2025

Senior Manager, Regulatory Affairs
Merit Medical Systems, Inc.
USA

(sign and seal)

17 February 2023

DiamondTOUCH angiography kit
INSTRUCTIONS FOR USE

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Merit's Angiography kit with Inflation device DiamondTOUCH is assembled with extension tubing attached to the syringe tip. A medium pressure 3-way stopcock may be packaged with the device. Merit's Angiography kit with Inflation device are available packaged with Merit Angioplasty Packs (MAP), which include a selection of hemostasis valves, extension line, guide wire insertion tools and/or torque devices:

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1. Angiography kit DiamondTOUCH, consisting of:
 - 1.1. Inflation Device DiamondTOUCH - 1 pc.
 - 1.2. Three-way stopcock - 1- pc. or 2pcs.
 - 1.3. Instruction for use.
2. Angiography kit DiamondTOUCH, consisting of:
 - 2.1. Inflation Device DiamondTOUCH - 1 pc.
 - 2.2. Three-way stopcock - 1 pc. or 2pcs.
 - 2.3. Hemostatic valve, execution options:
AccessPLUS/ Access-g/ MBA/ MBAPlus/ PhD/ Honor/ FLO50 - 1 pc.
 - 2.4. Guidewire Introducer - 1 pc.
 - 2.5. Torque Device - no more than 1 pcs.
 - 2.6. Extension line (if necessary) - 1 pc.
 - 2.7. Instruction for use.

INDICATION FOR USE

DiamondTOUCH angiography kit is designed for use during percutaneous transluminal angioplasty (PTA) and other intravascular therapeutic procedures that utilize a balloon angioplasty catheters or other interventional devices.

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
DESCRIPTION:

The DiamondTOUCH™ Inflation Device and Fluid Dispensing Syringe by Merit Medical is a 30mL disposable device with an integral pressure transducer, microcomputer, back-lit LCD, threaded plunger assembly with lock/release bar, a flexible high pressure extension tube, and a three-way stopcock. The DiamondTOUCH™ is designed to generate positive and negative pressure, and monitor positive pressures over a range of zero to +35ATM/BAR (zero to 514 PSI). The DiamondTOUCH Inflation Syringe is used to inflate and deflate balloon angioplasty catheters or other interventional devices and to measure the pressure and time of inflation within the balloon during the procedure. It is also used to dispense fluids into the body and monitor the pressure of that fluid.

PRESSURE RATING (ATM/BAR)	VOLUME (ML)	DISPLAY TYPE	TUBING LENGTH (IN/CM)
			
35 atm/bar	30 mL	Digital	13" (33.02 cm)
35 atm/bar	30 mL	Digital	36" (91.44 cm)

SPECIFICATION

Operating Range	-0.4 ATM/BAR to 35 ATM/BAR (-6 PSI to 514 PSI)
Accuracy	±3.0% of full scale at 10°C to 40°C
Fluid Dispensed	0.57ml ± 0.10ml fluid dispensed for each 360° clockwise turn of the syringe plunger handle.
Liquid Temp	10° C to 40° C
Humidity	20% to 90%
Operation Range	Non-condensing humidity
Battery Life	Fully active device, up to 10 hours

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The PhD Hemostasis Valve is intended to maintain hemostasis during the introduction/ withdrawal and use of diagnostic and interventional devices. The PhD Hemostasis Valve is indicated for maintaining a seal around diagnostic/interventional devices with an outside diameter of ≤0.096" (2.4 mm) during interventional procedures.


Honor Hemostasis Valve is a y-connector hemostasis valve with a rotating luer lock, a side port, and a push-pull hemostasis valve opener on the proximal end. The hemostasis valve is fully opened to its maximum diameter by pushing the valve opener distally. Pulling the valve opener proximally closes the hemostasis valve around the inserted interventional/diagnostic devices. With the valve in the closed position, devices can be advanced or pulled back as needed without adjustment of the hemostasis valve. The hemostasis valve is intended to maintain hemostasis during the introduction, use and withdrawal of diagnostic and interventional devices that have an outer diameter of 7 French or smaller.

The FLO50™ Hemostasis Valve is a small Y-body valve with a clear female luer side arm and male luer rotator. The Merit FLO50 Hemostasis Valve is intended to maintain hemostasis with interventional devices with outside diameters up to 9Fr providing a fluid-tight seal around interventional devices with little or no blood loss during clinical procedures.

Hemostasis Valves Characteristic: Inner Lumen					
AccessPLUS/ Access-g: 3mm/0.120 inch	MBA: 3mm/0.118 inch	MBA Plus: 3mm/0.118 inch	Honor: 2.4mm/0.094 inch	PhD: 2.4mm/0.096 inch	FLO50: 3mm/0.118 inch
					

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The main materials for the components with indirect contact: Inflation device, MAP	
Products Component	Material
Inflation Device	Rotator body & collar & housing & Retainer Cap & Barrel– Polycarbonate; Extension tubing – polyurethane; O ring – Silicone; Syringe Lubricant - Silicone Fluid; Handle body – ABS; Trigger - Nylon.
Hemostasis Valves	Clear Polycarbonate; Silicone; PTFE Teflon; Acetal homopolymer; PVC Tubing; Polyurethane; Nylon.
Extension line	PVC
Guide Wire Insertion Tool (plastic)	Polyethylene
Guide Wire Insertion Tool (metal)	Polycarbonate; Stainless Steel

INFLATION DEVICE PREPARATION:

1. Ensure the system is open to atmospheric pressure by opening the stopcock. Press the blue button behind the LCD display near the tubing to power the device on. The LCD will display "ZEro" for two seconds and then the device will be ready to use. At this point the syringe will begin its incremental time keeping.
2. To prepare syringe, turn the device with digital display facing down and aspirate up to 30 mL of contrast solution or other fluid into the syringe by squeezing the lock/release bar and pulling back on the handle. 3. Push handle against table or other solid surface to remove air in syringe.

CAUTION: Inspect the syringe tubing and stopcock (if used) to ensure that there is no air in the system.

The syringe will be set in the ATM/BAR mode when initially turned on. To change the pressure display to read in PSI, press and hold the blue button until "ATM/BAR" flashes four times. The display is now in "PSI" mode. To change back to ATM/BAR, press and hold the blue button once again.

NOTE: When in PSI mode, the tick marks on the left of the display that represent pressure will be limited to 350 PSI (23.8 ATM). If the DiamondTOUCH is pressurized past 350 PSI, the grouping of tick marks on the left will flash. The numerical digits in the center of the display will continue to show actual pressure throughout the device's pressure range (-6 to 514 PSI)

After an inflation or pressure monitored injection has been made, a graph bar or tick mark will remain to mark the highest point of pressure. Pressing the blue button once quickly will display last inflation information and a "↩" indicator on the display. After the next inflation has been started, the last inflation tick mark will disappear. **CAUTION:** If the LCD displays anything besides the pressure and time windows shown above, the syringe is defective. Please return the syringe to Merit Medical for credit.

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CAUTION: Inspect the syringe tubing and stopcock (if used) to ensure that there is no air in the system.

ATTACHING THE INFLATION DEVICE TO THE BALLOON

NOTE: Refer to the manufacturer's directions accompanying the balloon dilatation catheter or other interventional device for specific information on use, maximum inflation pressure, precautions, and warnings for that device. 1. Prepare and test the balloon catheter according to the catheter manufacturer's directions for use. 2. Create a fluid-fluid connection between the balloon and the syringe extension tube, connect the luer connectors securely. 3. Squeeze the lock/release bar and pull back on the plunger handle to apply a vacuum to the balloon.

BALLOON INFLATION AND DEFLATION

1. To inflate the balloon, squeeze the lock/release bar allowing the plunger to return to a resting position (0 ATM/BAR or PSI). Release grip on the lock/release bar, locking the plunger into position. To increase pressure, rotate handle clockwise until the desired pressure is achieved. Pressures above the maximum range will be indicated with flashing numbers. The tick mark will remain at the highest point of the last pressure reading. As the pressure decreases from the maximum pressure, the tick mark will begin to flash.

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CAUTION: To protect the threads of the lock-release bar, the pressure must be reduced to 25 ATM or lower before the quick release mechanism is used to deflate the angioplasty balloon. 2. To deflate balloon, rotate handle counterclockwise to release pressure to 25 ATM or lower. Squeeze the lock/release bar and pull back to generate a negative pressure. Release grip to lock the plunger in a negative pressure position. Pressures below the minimum range of the syringe will be indicated by flashing bars and a "NEg" in the pressure area. **DISPENSING FLUIDS AND MONITORING**

PRESSURES USING THE DIAMONDTouch SYRINGE

1. To slowly dispense fluids to the body rotate the plunger clockwise until the desired fluid is injected. To rapidly dispense fluids to the body squeeze the lock/release bar while pushing the plunger forward. The injection pressure will be displayed on the LCD

and the timer automatically starts once the device generates a positive pressure. Pressures above the maximum range of the syringe will be indicated with flashing numbers on the display. 2. A negative pressure may be generated by squeezing the lock/release bar and pulling back the plunger. Release grip on the lock/release bar to lock the plunger in the negative pressure position. 3. The timer will reset once the pressure returns to zero or less. Data associated with the last injection will be displayed when the blue button is depressed and held and the pressure is zero or less.

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RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RADIO FREQUENCY (RF) COMMUNICATIONS

EQUIPMENT AND THE DIAMONDTouch™ INFLATION DEVICE

The DiamondTOUCH Inflation Device is intended for use in an electromagnetic environment in which RF radiated disturbances are controlled. The user of the DiamondTOUCH Inflation Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the DiamondTOUCH Inflation Device as recommended below, according to the maximum output power of the communications equipment.

	Separation distance according to frequency of transmitter (in meters) m		
Rated maximum output power of transmitter (in watts) W	150 kHz to 80 MHz $d=[1.2]P$	80 MHz to 800 MHz $d=[1.2]P$	800 MHz to 2.5 GHz $d=[2.3]P$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	12	12	23

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NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should ensure that it is used in such an environment.

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1. Connect the sideport of the hemostasis valve to the Merit manifold assembly. Flush and fill the assembly with saline. To fill the valve section, open the valve, place one finger over the luer fitting and continue to fill the assembly.

2. Connect the hemostasis valve to the guiding catheter. Aspirate the valve to remove any trapped air and flush thoroughly with saline. Purge blood by opening the valve while continuing to flush the assembly. Close the valve when the blood has been purged. Inspect carefully for air bubbles and reflush if necessary.

3. Open the hemostasis valve and insert the guidewire or catheter. (If desired, insert the guidewire introducer and advance it through the valve.) Advance the guidewire or catheter an appropriate distance into the vasculature. Close the valve around the shaft of the guidewire, catheter or introducer. This forms a fluid-tight seal, yet does not inhibit movement of the guidewire, catheter or introducer.

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4. Refer to the dilatation catheter labeling for intended use, contraindications and potential complications associated with the use of dilatation catheters in PTCA and/or PTA procedure.

5. Withdraw the catheter to the point that the distal tip remains 10 to 20 cm inside the vasculature. Open the valve and withdraw the catheter completely. During insertion and removal through the valve, the balloon must be completely deflated to safely and easily traverse the valve.

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6. Once the diagnostic/interventional device is in place, if desired, close the compression seal to allow for pressure injections of up to 400 psi and/or secure the diagnostic/interventional device in position.

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5. Introduce the guiding catheter following the procedure recommended by the guiding catheter manufacturer. 6. Insert the guide wire, or guide wire and diagnostic/interventional catheter together, into the device. When inserting the guide wire by itself, use a guide wire insertion tool to protect the distal tip of the guide wire.
7. Perform the rest of the procedure following the recommendations of the respective device manufacturers.

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1. Connect the sideport of the hemostasis valve to a flush system. With the valve OPEN, flush and fill the assembly with saline. To fill the valve section, open the valve, place one finger over the luer fitting and continue to fill the assembly.
2. Connect the male luer rotator of the hemostasis valve to the catheter. Ensure the connection is airless. Aspirate the valve to remove any trapped air and flush thoroughly with saline. Purge blood by opening the valve while continuing to flush the assembly. Close the valve when the blood has been purged. Inspect carefully for air bubbles and reflush if necessary.
3. Open the hemostasis valve and insert the guidewire or catheter. (If desired, insert the guidewire introducer and advance it through the valve). Advance the guidewire or catheter an appropriate distance into the vasculature. Close the valve around the shaft of the guidewire, catheter or introducer. This forms a fluid-tight seal, yet does not inhibit movement of the guidewire, catheter or introducer.

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4. Read and follow the manufacturers IFU for devices to be used with the FLO50.
5. Withdraw the device to the point that the distal tip remains 10 to 20cm inside the vasculature. Open the valve and withdraw the device completely.
6. Disconnect the hemostasis valve from the catheter.

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- Balloon catheter with appropriate size • Contrast agent • Sterile saline solution • Luer lock syringe • Coaxial pair (sheath/dilator) with appropriate size • Guidewire with appropriate size • Appropriate guiding needle • Appropriate guiding and/or diagnostic catheter • Hemostasis valve (if not included)

TRANSPORTATION, STORAGE AND OPERATION CONDITIONS	
Transportation temperature	-25 °C to +55 °C
Storage temperature	+5 °C to +27 °C Store in a cool dry place in original packaging
Usage	15 °C to 35 °C, relative humidity ≤80%, atm. Pressure 84.0 – 106.7 kPa (630-800 mm Hg).

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WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING

WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

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
Manufacturer: Merit Medical Systems, Inc.

1600 West Merit Parkway, South Jordan, Utah 84095

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
Instructions for Use
basixCOMPAK
Dated February 17, 2023
Consisting of 10 pages

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

(sign and seal)

17 February 2023

basixCOMPAK angiography kit
INSTRUCTIONS FOR USE

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Merit's Angiography kit with Inflation device basixCOMPAK is assembled with extension tubing attached to the syringe tip. A medium pressure 3-way stopcock may be packaged with the device. Merit's Angiography kit with Inflation device are available packaged with Merit Angioplasty Packs (MAP), which include a selection of hemostasis valves, extension line, guide wire insertion tools and/or torque devices:

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1. Angiography kit basixCOMPAK, consisting of:

- 1.1. Inflation Device basixCOMPAK - 1 pc.
- 1.2. Three-way stopcock - 1- pc. or 2pcs.
- 1.3. Instruction for use.

2. Angiography kit basixCOMPAK, consisting of:

- 2.1. Inflation Device basixCOMPAK - 1 pc.
- 2.2. Three-way stopcock - 1 pc. or 2pcs.
- 2.3. Hemostatic valve, execution options:
AccessPLUS/ Access-9/ MBA/ MBAPlus/ PhD/ Honor/ FLO50 - 1 pc.
- 2.4. Guidewire Introducer - 1 pc.
- 2.5. Torque Device - no more than 1 pcs.
- 2.6. Extension line (if necessary) - 1 pc.
- 2.7. Instruction for use.

INDICATION FOR USE

basixCOMPAK angiography kit is designed for use during percutaneous transluminal angioplasty (PTA) and other intravascular therapeutic procedures that utilize a balloon angioplasty catheters or other interventional devices.

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
WARNING

- Do not inject any fluid if air bubbles are visible within the Hemostasis Valve.
- The Honor Hemostatic valve is not intended for use with pressure injections > 8 ATM/ BAR. Power injection at pressure greater than 8 ATM/BAR could result in leakage.
- Devices inserted across the hemostasis valve are not secure.
 - The device should only be used by experienced physicians trained in endovascular procedures.
 - To prevent damage to a balloon, it must be completely deflated during advancement/ withdrawal through the Hemostasis Valve.

- Any aspiration through the side port of the device should be made slowly and carefully to avoid the creation of a void at the proximal part of the connector. Failure could cause leakage.
- Failure to depress the cap prior to inserting a guide wire introducer could cause seal damage resulting in leakage and/or particulate embolization.
- Failure to open the secondary seal and compression seal prior to inserting or withdrawing a diagnostic /interventional device could cause damage to the device.
- Care should be taken to avoid over-tightening of the compression seal, as this may compromise the integrity of the diagnostic/interventional device lumen.
- Power injection at pressures greater than 400psi could result in leakage or detachment of PhD Hemostasis Valve components.
- The FLO50 hemostasis valve must be completely closed during aspiration or injection. The narrowest portion of the hemostasis valve has an inner diameter of 0.118 inches (3.00mm, approximately 9Fr).

DESCRIPTION:

The Merit basixCOMPAK Inflation Device is a 20mL disposable device capable of producing a maximum pressure of 30 ATM/bars, fitted with a threaded plunger assembly with lock/ release bar, a flexible high pressure extension tube, and a three-way medium pressure stopcock. The basixCOMPAK is designed to generate positive and negative pressure, and monitor positive pressures over a range of zero to +30ATM/BAR (zero to +441 PSI).

PRESSURE RATING (ATM/BAR)	VOLUME (ML)	DISPLAY TYPE	TUBING LENGTH (IN/CM)
			
30 atm/bar	20 mL	Analog	13" (33.02 cm)

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Adjustable Merit AccessPLUS™/ Merit Access-9™ Hemostasis Valve. The AccessPLUS/Access-9 Hemostasis Valve is recommended for maintaining a fluid-tight seal around percutaneous transluminal angioplasty catheters and guidewires.

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The PhD Hemostasis Valve™ has a 0.095" (2.44 mm) inside diameter. This device has two seals that operate independently: the compression seal and the secondary seal. The compression seal can be opened or closed by rotating the cap. Closing of the compression seal allows for pressure injections up to 400 psi and secures the diagnostic/interventional device in position within the vasculature.

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The PhD Hemostasis Valve is intended to maintain hemostasis during the introduction/ withdrawal and use of diagnostic and interventional devices. The PhD Hemostasis Valve is indicated for maintaining a seal around diagnostic/interventional devices with an outside diameter of $\leq 0.096"$ (2.44 mm) during interventional procedures.

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

The guide wire Insertion Tool is used to facilitate the passage of wires with distally shaped tips through the rotating hemostatic valve and dilatation catheters. It is manufactured in plastic and metal and accepts guide wires with diameters up to 0.020 inches, up to 0.025" for the Merit Blunt Introducer.

Characteristic	
Torque Device Range	 0.009" – 0.018" (0,22-0,45 mm)

Plastic Guide Wire Insertion Tool Inner Lumen Size	 0.028" (0,7mm)
Metal Guide Wire Insertion Tool Inner Lumen Size	 0.022", 0.025" (0,56 , 0,64 mm)

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
INSTRUCTION FOR USE: INFLATION DEVICE

This inflation device is used to inflate and deflate an angioplasty balloon or other interventional device, and to measure the pressure within the balloon.

NOTE: This device has not been cleared for dispensing fluids into the body.

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PRIMING WITH PRIMELOK™

The PrimeLok allows purging of air and excess fluid without squeezing the trigger. The basixCOMPAK is packaged with the PrimeLok and with the lock/release mechanism in the disengaged position.

1. To engage PrimeLok, squeeze trigger and slide PrimeLok into slot in plunger.

2. To prep syringe, simply aspirate up to 20mL of contrast solution into the inflation syringe by pulling back on the handle.

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4. To disengage PrimeLok, squeeze handle and slide PrimeLok out of slot. This will allow the trigger to lock in position and the device is now ready for use.

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CAUTION: If applied pressure does not indicate on gauge, discontinue use immediately and replace it with a new unit.

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NOTE: Refer to the manufacturer's directions accompanying the balloon dilatation catheter or other interventional device for specific information on use, maximum inflation pressure, precautions, and warnings for that device.

1. Prepare and test the balloon catheter according to the catheter manufacturer's directions for use.

2. Create a fluid-fluid connection between the balloon and the inflation syringe extension tube and connect the luer connectors securely.

3. Squeeze the trigger and pull back on the handle to apply a vacuum to the balloon.

BALLOON INFLATION AND DEFLATION:



1. To inflate the balloon, squeeze the trigger allowing the plunger to return to the resting position (0 ATM/BAR). Release grip on the trigger, locking the plunger into position.

To increase pressure, rotate the handle clockwise until the desired inflation pressure is reached. The lock mechanism maintains the pressure.

NOTE: Loss of pressure may indicate a leak in the system.

CAUTION: To protect the threads of the lock release handle, the device must be displaying 25 ATM or lower before quick release mechanism is used to deflate the angioplasty balloon.

2. To deflate the balloon, squeeze the trigger and pull back to generate a negative pressure. Release grip on the trigger to lock the plunger in the negative pressure position.

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PRECAUTIONS: Do not inject any fluid if air bubbles are visible within the valve. First aspirate the valve to remove the air, then flush the valve as described above. The hemostasis valve must be completely closed during aspiration or injection. The narrowest portion of the hemostasis valve has an inner diameter of 0.120 inches (0.305 cm) (approximately 9 Fr).

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The narrowest rigid portion of the hemostasis valve has an inner diameter of 0.118 inches (0.305 cm) (approximately 9 Fr).

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1. Attach the manifold to the side port of the device.
2. Connect the distal end of the device to the proximal end of the guiding catheter.
3. Flush with normal saline to remove air bubbles. To flush the valve segment, open the hemostasis valve by pushing the valve opener and continue to fill the assembly. To close the valve, pull the valve opener back to its original position.
4. Attach the pressure/infusion device to the manifold, if applicable. Check that all connections are secure to avoid air aspiration during the procedure.
5. Introduce the guiding catheter following the procedure recommended by the guiding catheter manufacturer. 6. Insert the guide wire, or guide wire and diagnostic/interventional catheter together, into the device. When inserting the guide wire by itself, use a guide wire insertion tool to protect the distal tip of the guide wire.
7. Perform the rest of the procedure following the recommendations of the respective device manufacturers.

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• Balloon catheter with appropriate size • Contrast agent • Sterile saline solution • Luer lock syringe • Coaxial pair (sheath/dilator) with appropriate size • Guidewire with appropriate size, • Appropriate guiding needle • Appropriate guiding and/or diagnostic catheter • Hemostasis valve (if not included)

TRANSPORTATION, STORAGE AND OPERATION CONDITIONS	
Transportation temperature	-25 °C to +55 °C
Storage temperature	+5 °C to +27 °C
Usage	15 °C to 35 °C, relative humidity ≤80%, atm. Pressure 84.0 – 106.7 kPa (630-800 mm Hg).

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bld 7, 1 floor, room 6, Russian Federation

+7 495 221 89 02

Applicable symbols	
Symbol	Description
	Sterilization: Ethylene-Oxide Gas
	Single use
	Do not use if package is damaged
	Non-pyrogenic
	EC mark logo - Notified body identification

Applicable symbols	
Symbol	Description
	Manufacturer
	Authorized Representative
	Caution! Consult accompanying documents

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
Instructions for Use
basixTOUCH
Dated February 17, 2023
Consisting of 11 pages

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
Merit's Angiography kit with Inflation device basixTOUCH is assembled with extension tubing attached to the syringe tip. A medium pressure 3-way stopcock may be packaged with the device. Merit's Angiography kit with Inflation device are available packaged with Merit Angioplasty Packs (MAP), which include a selection of hemostasis valves, extension line, guide wire insertion tools and/or torque devices:

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1. Angiography kit basixTOUCH, consisting of:

1.1. Inflation Device basixTOUCH - 1 pc.

1.2. Three-way stopcock - 1- pc. or 2pcs.

1.3. Instruction for use.

2. Angiography kit basixTOUCH, consisting of:

2.1. Inflation Device basixTOUCH - 1 pc.

2.2. Three-way stopcock - 1 pc. or 2pcs.

2.3. Hemostatic valve, execution options:

AccessPLUS/ Access-9/ MBA/ MBAPlus/ PhD/ Honor/ FLO50 - 1 pc.

2.4. Guidewire Introducer - 1 pc.

2.5. Torque Device - no more than 1 pcs.

2.6. Extension line (if necessary) - 1 pc.

2.7. Instruction for use.

INDICATION FOR USE

basixTOUCH angiography kit is designed for use during percutaneous transluminal angioplasty (PTA) and other intravascular therapeutic procedures that utilize a balloon angioplasty catheters or other interventional devices.

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WARNING

- Do not inject any fluid if air bubbles are visible within the Hemostasis Valve.

The Honor Hemostatic valve is not intended for use with pressure injections > 8 ATM/ BAR. Power injection at pressure greater than 8 ATM/BAR could result in leakage.


- Devices inserted across the hemostasis valve are not secure.
- The device should only be used by experienced physicians trained in endovascular procedures.
- To prevent damage to a balloon, it must be completely deflated during advancement/ withdrawal through the Hemostasis Valve.

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DESCRIPTION:

The basixTOUCH™ Inflation Syringe by Merit Medical is a 30mL disposable device with a threaded plunger assembly, a flexible high pressure extension tube. The basixTOUCH™ is designed to generate positive and negative pressure, and monitor positive pressures over a range of zero to +35ATM/BAR (zero to 514 PSI). This inflation device is used to inflate and deflate an angioplasty balloon or other interventional device, and to measure the pressure within the balloon.

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PRESSURE RATING (ATM/BAR)	VOLUME (ML)	DISPLAY TYPE	TUBING LENGTH (IN/CM)
			
35 atm/bar	30 mL	Analog	13" (33.02 cm)
40 atm/bar	30 mL	Analog	13" (33.02 cm)

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MBA Hemostasis Valve with Integral Introducer. The M3A Hemostasis Valve is designed for use during percutaneous transluminal angioplasty (PTA) and other intravascular therapeutic procedures that utilize a guiding catheter.

The PhD Hemostasis Valve™ has a 0.096" (2.44 mm) inside diameter. This device has two seals that operate independently: the compression seal and the secondary seal. The compression seal can be opened or closed by rotating the cap. Closing of the compression seal allows for pressure injections up to 400 psi and secures the diagnostic/interventional device in position within the vasculature.

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
The PhD Hemostasis Valve is intended to maintain hemostasis during the introduction/ withdrawal and use of diagnostic and interventional devices. The PhD Hemostasis Valve is indicated for maintaining a seal around diagnostic/interventional devices with an outside diameter of ≤ 0.096 " (2.44 mm) during interventional procedures.

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
Plastic Guide Wire Insertion Tool Inner Lumen Size	 0.028" (0,7mm)
Metal Guide Wire Insertion Tool Inner Lumen Size	 0.022", 0.025" (0,56 , 0,64 mm)
The main materials for the components with indirect contact: Inflation device, MAP	
Products Component	Material
Inflation Device	Rotator body & collar & housing & Retainer Cap & Barrel- Polycarbonate; Extension tubing – polyurethane; O ring – Silicone; Syringe Lubricant - Silicone Fluid; Handle body – ABS; Trigger - Nylon.
Hemostasis Valves	Clear Polycarbonate; Silicone; PTFE Teflon; Acetal homopolymer; PVC Tubing; Polyurethane; Nylon.
Extension line	PVC
Guide Wire Insertion Tool (plastic)	Polyethylene
Guide Wire Insertion Tool (metal)	Polycarbonate; Stainless Steel

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Before use, inspect the device and packaging to verify that no damage has occurred as a result of shipping. If the pressure gauge needle is not resting within the "0" box, do not use.

DEVICE PREPARATION:


1. To prepare syringe, turn the device with gauge facing down and aspirate up to 30 mL of contrast solution or other fluid into the syringe by squeezing the trigger and pulling back on the handle.
2. Push handle against table or other solid surface to remove air in syringe.

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NOTE: Refer to the manufacturer's directions accompanying the balloon dilatation catheter or other interventional device for specific information on use, maximum inflation pressure, precautions, and warnings for that device. 1. Prepare and test the balloon catheter according to the catheter manufacturer's directions for use.

2. Create a fluid-fluid connection between the balloon and the inflation syringe extension tube, connect the luer connectors securely.
3. Squeeze the trigger and pull back on the plunger handle to apply a vacuum to the balloon.

BALLOON INFLATION AND DEFLATION:

1. To inflate the balloon, squeeze the trigger to allow the plunger to return to a resting position (0 ATM/BAR or PSI). Release grip on the trigger, which will lock the plunger into position. To increase pressure, rotate handle clockwise until the desired pressure is achieved.

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CAUTION: If applied pressure does not indicate on gauge display, discontinue use immediately and replace it with a new unit 2. To deflate balloon, squeeze the trigger and pull back to generate a negative pressure. Release grip to lock the plunger in a negative pressure position.

CAUTION: To protect the threads of the lock release handle, the pressure must be reduced to 25 ATM or lower before the quick release mechanism is used to deflate the angioplasty balloon.

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4. Refer to the dilatation catheter labeling for intended use, contraindications and potential complications associated with the use of dilatation catheters in PTCA and/or PTA procedure.
5. Withdraw the catheter to the point that the distal tip remains 10 to 20 cm inside the vasculature. Open the valve and withdraw the catheter completely. During insertion and removal through the valve, the balloon must be completely deflated to safely and easily traverse the valve.
6. Disconnect the hemostasis valve from the guiding catheter.

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2. Connect the rotating male luer adapter of the hemostasis valve to the guiding catheter. Ensure the connection is airless.
3. To open the valve for catheter or guide wire insertion, with the hemostasis valve in the fully open position, push in and slightly turn the integral introducer. The valve is now ready to accept any devices.

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2. Connect the rotating male Luer adapter of the PhD Hemostasis Valve to the guiding catheter. Ensure the connection is airless.
 3. To open the PhD Hemostasis Valve for catheter or guide wire insertion, push on the back end cap to open the secondary seal. The PhD Hemostasis Valve is now ready to accept any diagnostic/interventional device.
 4. Insert the diagnostic/interventional device through an open secondary seal and an open compression seal.
- CAUTION:** PhD Hemostasis Valve will be fully open and allow for blood loss when secondary seal is open.
5. Advance the diagnostic/interventional device through an open compression seal. A fluid-tight seal, providing little to no blood loss, will form around the diagnostic or interventional device while not restricting movement of the device.
 6. Once the diagnostic/interventional device is in place, if desired, close the compression seal to allow for pressure injections of up to 400 psi and/or secure the diagnostic/ interventional device in position.
 7. Remove the diagnostic/interventional device through an open compression seal. If resistance is felt, open the secondary seal as well.

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1. Attach the manifold to the side port of the device.
2. Connect the distal end of the device to the proximal end of the guiding catheter.
3. Flush with normal saline to remove air bubbles. To flush the valve segment, open the hemostasis valve by pushing the valve opener and continue to fill the assembly. To close the valve, pull the valve opener back to its original position.
4. Attach the pressure/infusion device to the manifold, if applicable. Check that all connections are secure to avoid air aspiration during the procedure.
5. Introduce the guiding catheter following the procedure recommended by the guiding catheter manufacturer.
6. Insert the guide wire, or guide wire and diagnostic/interventional catheter together, into the device. When inserting the guide wire by itself, use a guide wire insertion tool to protect the distal tip of the guide wire.
7. Perform the rest of the procedure following the recommendations of the respective device manufacturers.

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TRANSPORTATION, STORAGE AND OPERATION CONDITIONS	
Transportation temperature	-25 °C to +55 °C
Storage temperature	+5 °C to +27 °C

Usage	15 °C to 35 °C, relative humidity ≤80%, atm. Pressure 84.0 – 106.7 kPa (630-800 mm Hg).
-------	--

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
bld 7, 1 floor, room 6, Russian Federation

+7 495 221 89 02

Applicable symbols	
Symbol	Description
	Sterilization: Ethylene-Oxide Gas
	Single use
	Do not use if package is damaged
	Non-pyrogenic
	EC mark logo - Notified body identification
	Manufacturer
	Authorized Representative
	Caution! Consult accompanying documents

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
Instructions for Use
Big60
Dated February 17, 2023
Consisting of 4 pages

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(sign and seal)

17 February 2023

BIG60 angiography kit
INSTRUCTIONS FOR USE

Only - Caution: Federal (USA) law restricts this device to sale by or on the order of a physician trained in angiography and percutaneous transluminal coronary angioplasty (PTCA), and/or percutaneous transluminal angioplasty (PTA).

Merit's Angiography kit with Inflation device BIG60 is assembled with extension tubing attached to the syringe tip. A medium pressure 3-way stopcock is packaged with the device.

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- 1. Angiography kit BIG60, consisting of:
 - 1.1. Inflation Device BIG60 - 1 pc.
 - 1.2. Three-way stopcock - 1 pc.
 - 1.3. Instruction for use.

INDICATION FOR USE

BIG60 angiography kit is designed for use during percutaneous transluminal angioplasty (PTA) and other intravascular therapeutic procedures that utilize a balloon angioplasty catheters or other interventional devices.

CONTRAINDICATIONS: None known.

ADVERSE EFFECTS: None known.


PRECAUTIONS

- Inspect the BIG60 Inflation Syringe and packaging for damage prior to use. Do not use product if opened or damaged. Confirm the product is consistent with the package label. Contact Customer Service to report and replace damaged product.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the product and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the product may lead to injury, illness, or death of the patient.
- Graduations on syringe barrel are for reference only and are not intended for exact measurement.

DESCRIPTION:

The BIG60 Inflation Syringe is a 60ml disposable inflation device capable of producing a maximum pressure of 12 ATM/BAR, fitted with a threaded plunger assembly with a lock/release bar, a flexible high pressure tube, and a three -way medium pressure stopcock. The BIG60 Inflation Syringe is designed to generate and monitor inflation pressures over a range of -.68 to 12ATM/ BAR (-10 to 176 psi) with accuracy of 3% ATM/BAR. This inflation device is used to inflate and deflate an angioplasty balloon or other interventional device, and to measure the pressure within the balloon.

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PRESSURE RATING (ATM/BAR)	VOLUME (ML)	DISPLAY TYPE	TUBING LENGTH (IN/CM)
			
12 atm /bar	60 mL	Analog	13" (33.02 cm)

INSTRUCTIONS FOR USE

1. Prior to use, free the plunger tip by twisting the syringe plunger/handle 360 degrees clockwise.
2. To prep syringe, simply aspirate up to 60ml of contrast solution into the syringe by pulling back on the syringe handle. **caution:** Inspect the syringe tubing and stopcock (if used) to ensure there is no air in the system.
3. Remove any excess air by orienting the syringe upwards, squeezing the trigger located on the syringe handle and pushing the handle forward. Push forward until the syringe plunger tip is oriented with the black arrow printed on the syringe barrel (approximately 35ml).

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NOTE: Refer to the manufacturer's directions accompanying the balloon dilation catheter, or other interventional device, for specific information on use, maximum inflation pressure, precautions, and warnings for that device.

1. Prepare and test the balloon catheter according to the catheter manufacturer's instructions for use.
2. Create a fluid-fluid connection between the balloon and the inflation syringe inspection extension tube and connect the luer connectors securely.
3. Squeeze the trigger and pull back on the handle to apply vacuum to the balloon.

BALLOON INFLATION AND DEFLATION

1. To inflate the balloon, squeeze the trigger allowing the plunger to return to the resting position (0 ATM/BAR). Release grip on the trigger, locking the plunger into position.
2. To increase pressure, rotate the handle clockwise until the desired inflation pressure is reached. The locking mechanism maintains the pressure.

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3. To deflate the balloon, (rotate the handle counter-clockwise to relieve pressures above 6 ATM) squeeze the trigger and pull the handle back to generate the desired negative pressure. Release the grip on the trigger to lock the plunger in the negative pressure position.

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- Balloon catheter with appropriate size
- Contrast agent
- Sterile saline solution
- Luer lock syringe
- Coaxial pair (sheath/dilator) with appropriate size
- Guidewire with appropriate size
- Appropriate guiding needle
- Appropriate guiding and/or diagnostic catheter
- Hemostasis valve

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WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was

not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

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

Instructions for Use
Blue Diamond
Dated February 17, 2023
Consisting of 10 pages


ACKNOWLEDGEMENT

State of Utah)
 §
County of Salt Lake)

On this 17th day of February 2023, before me, Jennifer Fordham, a notary public, personally appeared Tom Haueter, Senior Manager, Regulatory Affairs proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand and official seal.

Notary Public

My Commission Expires: June 27, 2025

**Senior Manager, Regulatory Affairs
Merit Medical Systems, Inc.
USA**

(sign and seal)

17 February 2023

Blue Diamond angiography kit

INSTRUCTIONS FOR USE

RX Only - Caution: Federal (USA) law restricts this device to sale by or on the order of a physician trained in angiography and percutaneous transluminal coronary angioplasty (PTCA), and/or percutaneous transluminal angioplasty (PTA).

Merit's Angiography kit with Inflation device Blue Diamond is assembled with extension tubing attached to the syringe tip. A medium pressure 3-way stopcock may be packaged with the device. Merit's Angiography kit with Inflation device are available packaged with Merit Angioplasty Packs (MAP), which include a selection of hemostasis valves, extension line, guide wire insertion tools and/or torque devices:

CONTENT

Angiography kit in the following versions:

1. Angiography kit Blue Diamond, consisting of:

- 1.1. Inflation Device Blue Diamond - 1 pc.
- 1.2. Three-way stopcock - 1- pc. or 2pcs.
- 1.3. Instruction for use.

2. Angiography kit Blue Diamond, consisting of:

- 2.1. Inflation Device Blue Diamond - 1 pc.
- 2.2. Three-way stopcock - 1 pc. or 2pcs.
- 2.3. Hemostatic valve, execution options:
AccessPLUS/ Access-9/ MBA/ MBAPlus/ PhD/ Honor/ FLO50 - 1 pc.
- 2.4. Guidewire Introducer - 1 pc.
- 2.5. Torque Device - no more than 1 pcs.
- 2.6. Extension line (if necessary) - 1 pc.
- 2.7. Instruction for use.

INDICATION FOR USE

Blue Diamond angiography kit is designed for use during percutaneous transluminal angioplasty (PTA) and other intravascular therapeutic procedures that utilize a balloon angioplasty catheters or other interventional devices.

CONTRAINDICATIONS: None known.

ADVERSE EFFECTS: None known.

PRECAUTIONS

- Read manufacturer's instructions for the use of catheters, guide wires, and introducers prior to use.
- Inspect the device prior to use for any damage. Do not use if packaging is opened, damaged, or broken.
- Follow local guidelines and practices regulating the disposal of the device after use.
- Contents are sterile (via ethylene oxide) and non-pyrogenic.

REUSE PRECAUTION STATEMENT


For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

WARNING

- Do not inject any fluid if air bubbles are visible within the Hemostasis Valve.
- The Honor Hemostatic valve is not intended for use with pressure injections > 8 ATM/ BAR. Power injection at pressure greater than 8 ATM/BAR could result in leakage.
- Devices inserted across the hemostasis valve are not secure.
 - The device should only be used by experienced physicians trained in endovascular procedures.
 - To prevent damage to a balloon, it must be completely deflated during advancement/ withdrawal through the Hemostasis Valve.
 - Any aspiration through the side port of the device should be made slowly and carefully to avoid the creation of a void at the proximal part of the connector. Failure could cause leakage.
 - Failure to depress the cap prior to inserting a guide wire introducer could cause seal damage resulting in leakage and/or particulate embolization.
 - Failure to open the secondary seal and compression seal prior to inserting or withdrawing a diagnostic /interventional device could cause damage to the device.
 - Care should be taken to avoid over-tightening of the compression seal, as this may compromise the integrity of the diagnostic/interventional device lumen.
 - Power injection at pressures greater than 400psi could result in leakage or detachment of PhD Hemostasis Valve components.
 - The FLO50 hemostasis valve must be completely closed during aspiration or injection. The narrowest portion of the hemostasis valve has an inner diameter of 0.118 inches (3.00mm, approximately 9Fr).

DESCRIPTION:

The Blue Diamond™ Inflation Device and Fluid Dispensing Syringe by Merit Medical is a 20ml disposable device with an integral pressure transducer, microcomputer, back-lit LCD, threaded plunger assembly with lock/release bar, a flexible high pressure extension tube, and a three-way medium pressure stopcock. The Blue Diamond is designed to generate and monitor pressures over a range of -0.4 to +30.0 ATM/BAR (-6 to 44.1 PSI). The Blue Diamond syringe dispenses 0.45ml of fluid $\pm$ 0.07ml for each 360° turn of the syringe plunger handle. This inflation device is used to inflate and deflate an angioplasty balloon or other interventional devices, and to measure the pressure within the balloon.

PRESSURE RATING (ATM/BAR)	VOLUME (ML)	DISPLAY TYPE	TUBING LENGTH (IN/CM)
30 atm/bar	20 mL	Digital	13" (33.02 cm)

SPECIFICATIONS

Operating Range:	-0.4 ATM to 30 ATM (1 ATM = 1BAR = 14.7 PSI)
Accuracy:	$\pm$ 3.0 % of full scale typical
Fluid Dispensed:	0.45ml $\pm$ 0.07ml fluid dispensed for each 360° clockwise turn of the syringe plunger handle.
Liquid Temp:	10° C to 40° C
Humidity	20% to 90%
Operation Range:	Non-condensing humidity
Battery Life:	Fully active device, up to 10 hours

The Merit Angioplasty Pack (MAP) consists of a hemostasis valve, Guide Wire Insertion Tool and/or Pin Vise torque device. The hemostasis valves, Guide Wire Insertion Tools and torque devices are all non-invasive devices.

Merit Angioplasty Pack (MAP): The MAP allows the injection of contrast and pressure monitoring through a guiding catheter, regardless of balloon catheter position. The hemostatic valve is designed to maintain hemostasis during insertion, use and withdrawal of diagnostic and interventional devices. They can be used in conjunction with inflation devices and other medical devices (introducer, guide wire, catheter, guide catheter) during angioplasty.

The hemostasis valve tightens based on the device type. The polycarbonate used to manufacture the adapter is transparent to aid in visualizing entrapped air. The entire device has no sharp edges to prevent tearing of surgical gloves.

Adjustable Merit AccessPLUS™/ Merit Access-g™ Hemostasis Valve. The AccessPLUS/Access-g Hemostasis Valve is recommended for maintaining a fluid-tight seal around percutaneous transluminal angioplasty catheters and guidewires.

MBA Hemostasis Valve with Integral Introducer. The MBA Hemostasis Valve is designed for use during percutaneous transluminal angioplasty (PTA) and other intravascular therapeutic procedures that utilize a guiding catheter.

The PhD Hemostasis Valve™ has a 0.096" (2.4 mm) inside diameter. This device has two seals that operate independently: the compression seal and the secondary seal. The compression seal can be opened or closed by rotating the cap. Closing of the compression seal allows for pressure injections up to 400 psi and secures the diagnostic/interventional device in position within the vasculature.

The secondary seal is a diaphragm seal that forms around the diagnostic/interventional devices as they move into and out of the vasculature. This seal provides minimal fluid loss without restricting device movement. The secondary seal is open when the cap is pressed down, and closed when the cap is released. An open secondary seal allows air to be purged from the device and allows insertion/withdrawal of diagnostic/interventional devices.

The PhD Hemostasis Valve is intended to maintain hemostasis during the introduction/ withdrawal and use of diagnostic and interventional devices. The PhD Hemostasis Valve is indicated for maintaining a seal around diagnostic/interventional devices with an outside diameter of $\leq$ 0.096" (2.4 mm) during interventional procedures.

Honor Hemostasis Valve is a y-connector hemostasis valve with a rotating luer lock, a side port, and a push-pull hemostasis valve opener on the proximal end. The hemostasis valve is fully opened to its maximum diameter by pushing the valve opener distally. Pulling the valve opener proximally closes the hemostasis valve around the inserted interventional/diagnostic devices. With the valve in the closed position, devices can be advanced or pulled back as needed without adjustment of the hemostasis valve.

The hemostasis valve is intended to maintain hemostasis during the introduction, use and withdrawal of diagnostic and interventional devices that have an outer diameter of 7 French or smaller.


The FLO50™ Hemostasis Valve is a small Y-body valve with a clear female luer side arm and male luer rotator

The Merit FLO50 Hemostasis Valve is intended to maintain hemostasis with interventional devices with outside diameters up to 9Fr providing a fluid-tight seal around interventional devices with little or no blood loss during clinical procedures.

Hemostasis Valves Characteristic: Inner Lumen					
AccessPLUS/ Access-9: 3mm/0.120 inch	MBA: 3mm/0.118 inch	MBA Plus: 3mm/0.118 inch	Honor: 2.4mm/0.094 inch	PhD: 2.4mm/0.096 inch	FLO50: 3mm/0.118 inch
					

Merit's Torque device provides an increased mechanical advantage in order to aid in twisting and manipulating the guide wire. The physician is able to grasp the guide wire securely with this device, allowing for greater precision and control for the physician. This torque device is brightly colored in order to improve visualization under low light conditions and is available in either green or salmon colors.

The guide wire Insertion Tool is used to facilitate the passage of wires with distally shaped tips through the rotating hemostatic valve and dilatation catheters.

Characteristic	
Torque Device Range	 0.009" – 0.018" (0.23mm – 0.46mm)
Plastic Guide Wire Insertion Tool Inner Lumen Size	 0.028" (0.71mm)
Metal Guide Wire Insertion Tool Inner Lumen Size	 0.022" (0.56mm)

The main materials for the components with indirect contact: Inflation device, MAP	
Products Component	Material
Inflation Device	Rotator body & collar & housing & Retainer Cap & Barrel – Polycarbonate; Extension tubing – polyurethane; O ring – Silicone; Syringe Lubricant - Silicone Fluid; Handle body – ABS; Trigger - Nylon.
Hemostasis Valves	Clear Polycarbonate; Silicone; PTFE Teflon; Acetal homopolymer; PVC Tubing; Polyurethane; Nylon.
Extension line	PVC
Guide Wire Insertion Tool (plastic)	Polyethylene
Guide Wire Insertion Tool (metal)	Polycarbonate; Stainless Steel

INSTRUCTIONS FOR USE: INFLATION DEVICE

Before use, inspect the device and packaging to verify that no damage has occurred as a result of shipping.

Ensure the line is open to atmospheric pressure by opening the stopcock. Press the blue button behind the LCD display near the tubing to power the device on. The LCD will display "2Ero" for two seconds and then the device will be ready to use. At this point the syringe will begin its incremental time keeping. The syringe will be set in the ATM/BAR mode when initially turned on. To change the pressure

display to read in PSI, press and hold the blue button until "ATM/ BAR" flashes four times. The user is now in "PSI" mode. To change back to ATM/BAR, press and hold the blue button once again.

NOTE: When in PSI mode, the tick marks on the left of the display that represent pressure will be limited to 300 PSI (20.4 ATM). If the Blue Diamond is pressurized past 300 PSI, the grouping of tick marks on the left will flash. The numerical digits in the center of the display will continue to show actual pressure throughout the device's pressure range (-6 to 441 PSI) After an inflation or pressure monitored injection has been made, a graph bar or tick mark will remain to mark the highest point of pressure. Pressing the blue button once quickly will display last inflation information and a "↶" indicator on the display. After the next inflation has been started, the last inflation tick mark will disappear.

CAUTION: If "ER" is displayed in the pressure area and a number appears in the time window, the syringe is defective. Please return the syringe to Merit Medical for credit.

NOTE: To conserve power the backlight will automatically turn off after ten minutes of inactivity. However, the microprocessor will continue to monitor the pressure. Pressing the blue button or inflating the balloon will reactivate the backlight. The device will power itself off after 90 consecutive minutes at zero pressure.


DEVICE PREPARATION:

To prep syringe simply aspirate up to 20 ml of contrast solution or other fluid into the syringe by squeezing the trigger and pulling back on the handle.

PRIMING WITH PRIMELOK™

The PrimeLok allows purging of air and excess fluid without squeezing the trigger. The Blue Diamond is packaged with the PrimeLok in the disengaged position.


1. To engage PrimeLok, squeeze trigger and slide PrimeLok into slot.

2. To prep syringe, simply aspirate up to 20ml of contrast solution or fluid to be dispensed into the inflation syringe by pulling back on the plunger handle.


CAUTION: Inspect the syringe tubing and stopcock (if used) to insure that there is no air in the system.

3. Push handle against table to remove air in syringe.

4. To disengage PrimeLok, squeeze trigger and slide PrimeLok out of slot. This will allow the plunger to lock in position and the device is now ready for use.

NOTE: The PrimeLok must be disengaged before pressure can be maintained by the lock/release mechanism.

ATTACHING THE INFLATION DEVICE TO THE BALLOON:

NOTE: Refer to the manufacturer's directions accompanying the balloon dilatation catheter or other interventional device for specific information on use, maximum inflation pressure, precautions, and warnings for that device. 1. Prepare and test the balloon catheter according to the catheter manufacturer's directions for use. 2. Create a fluid-fluid connection between the balloon and the syringe extension tube, connect the luer connectors securely. 3. Squeeze the trigger and pull back on the plunger handle to apply a vacuum to the balloon.

BALLOON INFLATION AND DEFLATION:

1. To inflate the balloon, squeeze the trigger allowing the plunger to return to a resting position (0 ATM/BAR or PSI). Release grip on the trigger, locking the plunger into position. To increase pressure, rotate handle clockwise until the desired pressure is achieved. Pressures above the maximum range will be indicated with flashing numbers. The tick mark will remain at the highest point of the last pressure reading. As the pressure decreases from the maximum pressure, the tick mark will begin to flash.

NOTE: Significant loss of pressure may indicate a leak in the system.

CAUTION: To protect the threads of the lock release handle, the pressure must be reduced to 25 ATM or lower before the quick release mechanism is used to deflate the angioplasty balloon. 2. To deflate balloon, rotate handle counterclockwise to release pressure to 25 ATM or lower. Squeeze the trigger and pull back to generate a negative pressure. Release grip to lock the plunger in a negative pressure position. Pressures below the minimum range of the syringe will be indicated by flashing bars and a "NEg" in the pressure area.

FLUID DISPENSING INDICATIONS AND USAGE:

This device is intended for use by healthcare professionals to dispense fluids to the body from the Blue Diamond syringe and monitor the pressure of that fluid.

DISPENSING FLUIDS AND MONITORING PRESSURES USING THE BLUE DIAMOND SYRINGE

1. To slowly dispense fluids to the body rotate the plunger clockwise until the desired fluid is injected. To rapidly dispense fluids to the body squeeze the trigger while pushing the plunger forward. The injection pressure will be displayed on the LCD and the timer automatically starts once the device generates a positive pressure. Pressures above the maximum range of the syringe will be indicated with flashing numbers on the display. 2. A negative pressure may be generated by squeezing the trigger and pulling back the plunger. Release grip on the trigger to lock the plunger in the negative pressure position. 3. The timer will reset once the pressure returns to zero or less. Data associated with the last injection will be displayed when the blue button is depressed and held and the pressure is zero or less.

CAUTION: This syringe is capable of generating high fluid pressures in a closed fluid system. The volume change of fluid dispensed may not be accurate due to compliance of the plastic components as pressure changes.

CAUTION: If applied pressure does not indicate on gauge/ digital display, discontinue use immediately and replace it with a new unit.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RADIO FREQUENCY (RF) COMMUNICATIONS EQUIPMENT AND THE BLUE DIAMOND™ INFLATION DEVICE

The Blue Diamond Inflation Device is intended for use in an electromagnetic environment in which RF radiated disturbances are controlled. The user of the Blue Diamond Inflation Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Blue Diamond Inflation Device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (in watts) W	Separation distance according to frequency of transmitter (in meters) m		
	150 kHz to 80 MHz $d = \left[1.2 \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[1.2 \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[2.3 \right] \sqrt{P}$
0.01	0.12	0.12	0.2
0.1	0.37	0.37	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 – At 80 MHz and 800 MHz, the separation distance for the higher frequency applies.

NOTE 2 – These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The EQUIPMENT is intended for use in the electromagnetic environment specified below.

The customer or the user of the EQUIPMENT should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-20	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Increasing relative humidity will reduce the potential for ESD related difficulties

INSTRUCTIONS FOR USE: AccessPLUS/ Merit Access-9 HEMOSTASIS VALVE

Inspect the device prior to use to verify that no damage has occurred during shipping.

1. Connect the sideport of the hemostasis valve to the Merit manifold assembly. Flush and fill the assembly with saline. To fill the valve section, open the valve, place one finger over the luer fitting and continue to fill the assembly.
2. Connect the hemostasis valve to the guiding catheter. Aspirate the valve to remove any trapped air and flush thoroughly with saline. Purge blood by opening the valve while continuing to flush the assembly. Close the valve when the blood has been purged. Inspect carefully for air bubbles and reflush if necessary.
3. Open the hemostasis valve and insert the guidewire or catheter. (If desired, insert the guidewire introducer and advance it through the valve.) Advance the guidewire or catheter an appropriate distance into the vasculature. Close the valve around the shaft of the guidewire, catheter or introducer. This forms a fluid-tight seal, yet does not inhibit movement of the guidewire, catheter or introducer.

WARNING: It is important that the valve be closed tight enough to prevent blood leakage, yet not so tight as to restrict function.

4. Refer to the dilatation catheter labeling for intended use, contraindications and potential complications associated with the use of dilatation catheters in PTCA and/or PTA procedure.
5. Withdraw the catheter to the point that the distal tip remains 10 to 20 cm inside the vasculature. Open the valve and withdraw the catheter completely. During insertion and removal through the valve, the balloon must be completely deflated to safely and easily traverse the valve.

6. Disconnect the hemostasis valve from the guiding catheter.

PRECAUTIONS: Do not inject any fluid if air bubbles are visible within the valve. First aspirate the valve to remove the air, then flush the valve as described above. The hemostasis valve must be completely closed during aspiration or injection. The narrowest portion of the hemostasis valve has an inner diameter of 0.120 inches (0.3 cm) (approximately 9 Fr).

INSTRUCTIONS FOR USE: MBA HEMOSTASIS VALVE

Inspect the device prior to use to verify that no damage has occurred during shipping.

1. Connect the sideport of the fully open hemostasis valve to the manifold assembly. Flush and fill the assembly with saline. To fill the valve section, ensure the thumb wheel is in the fully open position, open the slit valve by pushing in and slightly turning the integral introducer through the slit valve to fully open it. Place one finger over the male luer fitting and continue to fill the assembly until completely purged of air. Pull back to disengage the integral introducer and engage the secondary seal.

2. Connect the rotating male luer adapter of the hemostasis valve to the guiding catheter. Ensure the connection is airless.

3. To open the valve for catheter or guide wire insertion, with the hemostasis valve in the fully open position, push in and slightly turn the integral introducer. The valve is now ready to accept any devices.

Caution: Valve will be fully open and allow for blood loss when integral introducer is engaged.

4. Insert the catheter or wire into the integral introducer and advance it to an appropriate distance into the vasculature. Disengage the integral introducer.

5. A fluid tight seal, providing little or no blood loss, will form around the interventional device yet does not restrict movement of the guide wire, catheter or introducer.

6. To secure the interventional device in place, turn the adjustable thumb wheel to tighten.

Caution: Do not inject any fluid if air bubbles are visible within the valve.

The narrowest rigid portion of the hemostasis valve has an inner diameter of 0.118 inches (0.3 cm) (approximately 9 Fr).

INSTRUCTIONS FOR USE: PhD HEMOSTASIS VALVE

1. Connect the sideport of the PhD Hemostasis Valve to the manifold assembly. With the compression valve open, flush and fill the assembly with saline. To fill the valve section, ensure the secondary seal is open by pushing down on the cap. Place one finger over the male Luer fitting and continue to fill the assembly until completely purged of air. Gently release the back end cap to engage the secondary seal.

2. Connect the rotating male Luer adapter of the PhD Hemostasis Valve to the guiding catheter. Ensure the connection is airless.

3. To open the PhD Hemostasis Valve for catheter or guide wire insertion, push on the back end cap to open the secondary seal. The PhD Hemostasis Valve is now ready to accept any diagnostic/interventional device.

4. Insert the diagnostic/interventional device through an open secondary seal and an open compression seal.

CAUTION: PhD Hemostasis Valve will be fully open and allow for blood loss when secondary seal is open.

5. Advance the diagnostic/interventional device through an open compression seal. A fluid-tight seal, providing little to no blood loss, will form around the diagnostic or interventional device while not restricting movement of the device.

6. Once the diagnostic/interventional device is in place, if desired, close the compression seal to allow for pressure injections of up to 400 psi and/or secure the diagnostic/ interventional device in position.

7. Remove the diagnostic/interventional device through an open compression seal. If resistance is felt, open the secondary seal as well.

INSTRUCTIONS FOR USE: HONOR HEMOSTASIS VALVE

1. Attach the manifold to the side port of the device.

2. Connect the distal end of the device to the proximal end of the guiding catheter.

3. Flush with normal saline to remove air bubbles. To flush the valve segment, open the hemostasis valve by pushing the valve opener and continue to fill the assembly. To close the valve, pull the valve opener back to its original position.

4. Attach the pressure/infusion device to the manifold, if applicable. Check that all connections are secure to avoid air aspiration during the procedure.

5. Introduce the guiding catheter following the procedure recommended by the guiding catheter manufacturer. 6. Insert the guide wire, or guide wire and diagnostic/interventional catheter together, into the device. When inserting the guide wire by itself, use a guide wire insertion tool to protect the distal tip of the guide wire.

7. Perform the rest of the procedure following the recommendations of the respective device manufacturers.

INSTRUCTIONS FOR USE: FLO50 HEMOSTASIS VALVE

1. Connect the sideport of the hemostasis valve to a flush system. With the valve OPEN, flush and fill the assembly with saline. To fill the valve section, open the valve, place one finger over the luer fitting and continue to fill the assembly.

2. Connect the male luer rotator of the hemostasis valve to the catheter. Ensure the connection is airless. Aspirate the valve to remove any trapped air and flush thoroughly with saline. Purge blood by opening the valve while continuing to flush the assembly. Close the valve when the blood has been purged. Inspect carefully for air bubbles and reflush if necessary.

3. Open the hemostasis valve and insert the guidewire or catheter. (If desired, insert the guidewire introducer and advance it through the valve). Advance the guidewire or catheter an appropriate distance into the vasculature. Close the valve around the shaft of the guidewire, catheter or introducer. This forms a fluid-tight seal, yet does not inhibit movement of the guidewire, catheter or introducer.

WARNING: It is important that the valve be closed tight enough to prevent blood leakage, yet not so tight as to restrict function.

4. Read and follow the manufacturers IFU for devices to be used with the FLO50.

5. Withdraw the device to the point that the distal tip remains 10 to 20cm inside the vasculature. Open the valve and withdraw the device completely.
6. Disconnect the hemostasis valve from the catheter.

Components required for PTA procedure but supplied separately:

• Balloon catheter with appropriate size • Contrast agent • Sterile saline solution • Luer lock syringe • Coaxial pair (sheath/dilator) with appropriate size • Guidewire with appropriate size • Appropriate guiding needle • Appropriate guiding and/or diagnostic catheter • Hemostasis valve (if not included)

TRANSPORTATION, STORAGE AND OPERATION CONDITIONS	
Transportation temperature	-25 °C to +60 °C
Storage temperature	+5 °C to +27 °C Store in a cool dry place in original packaging
Usage	15 °C to 35 °C, relative humidity ≤80%, atm. Pressure 84.0 – 106.7 kPa (630-800 mm Hg).

Use Conditions: in medical and preventive treatment facility.

SHELF LIFE

Shelf life – 3 years.

STERILIZATION

Ethylene Oxide

DISPOSAL

Dispose the devices according to local environmental requirements.

Environmental Monitoring

Dispose the devices according to local environmental requirements.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.3684-21) Class B.

WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

Russian Authorized Representative:


Merit Technologies LLC

119017, Moscow,

1st Kazachiy pereulok,

bld 7, 1 floor, room 6, Russian Federation

+7 495 221 89 02

Applicable symbols	
Symbol	Description
	Sterilization: Ethylene-Oxide Gas
	Single use
	Do not use if package is damaged
	Non-pyrogenic
	EC mark logo - Notified body identification
	Manufacturer
	Authorized Representative
	Caution! Consult accompanying documents
	Interference may occur in the vicinity of equipment marked with this symbol
	Contains Batteries - Do Not Remove

Manufacturer: Merit Medical Systems, Inc.
 1600 West Merit Parkway, South Jordan, Utah 84095
 U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

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