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Instructions For Use

for relieve dry eye symptoms, eye discomfort, and to protect the cornea as a tear substitute

Read all of this instructions for use carefully before you start using this product.

Single use device

Keep this instructions for use; you may need to read it again.

If you have any further questions, ask your doctor or pharmacist.

If symptoms persist and side effects appear, see your doctor.

Name of Medical Device

Ocutears® Hydro+ Eye drops, solution in consist of:

1. Bottle 10 ml – 1 pc;
2. Instruction for use – 1 pc.

What are eye drops in the form of Ocutears® Hydro+ Eye drops, solution?

The medical device is a lubricating ophthalmic solution based on highly concentrated cross-linked sodium hyaluronate (0.2%), designed to treat dry eye symptoms, eye discomfort and corneal protection.

The product is sterile and does not contain preservatives.

Cross-linked sodium hyaluronate imparts high viscoelasticity to the product, providing a high degree of hydration and retention to the ocular surface for long-term relief of dry eye symptoms.

Information about the manufacturer of the medical device

MANUFACTURER:

OFFHEALTH, S.P.A.

Via Giovanni Paisiello, 10-50144 FIRENZE (FL)- Italy

SITE OF MANUFACTURING:

1. **COC FARMACEUTICI S.R.L.**

IT-40019 Sant'Agata Bolognese (BO) – Via Modena 15, Italy-

2. **COC FARMACEUTICI S.R.L.**

IT-41016 Novi Di Modena, Fraz. Rovereto S/S (MO) – Via Chiesa Sud 156 CDEFGHI, Italy

Usage conditions

The product is used at home or in a medical facility.

Potential user

The intended patient population is patients suffering from dry eye symptoms and/or to protect the cornea as a tear substitute (artificial tear).

Purpose of the medical device

Ocutears® Hydro + is recommended for relieve dry eye symptoms, eye discomfort, and to protect the cornea as a tear substitute (artificial tear), providing long-lasting hydration.

It can also be used when wearing contact lenses to moisturize the eyes.

Indications:

Ocutears® Hydro + indications:

- for the symptomatic treatment of dry eye syndrome, eye discomfort, as well as for the protection of the cornea as a substitute for tear fluid, providing long-term hydration;
- Can be instilled to lubricate the eyes even when using contact lenses.

Ocutears® Hydro + is beneficial for and after:

- prolonged exposure to the sun;
 - stay in air-conditioned rooms;
 - prolonged use of a phone, tablet, laptop;
 - contact with cleaning agents.
- In addition, drops can relieve discomfort after eye surgery.

Contraindications

Hypersensitivity to device components or other substances closely related to it from a chemical point of view.

Possible Adverse Reactions

As with other eye drops, there may be mild temporary discomfort associated with instillation of the device.

ADVERSE REACTIONS, REPORT THE DOCTOR OR PHARMACIST.

Warnings

- Do not use the device if you are allergic or hypersensitive to any of the ingredients.

- Avoid swallowing.
- Do not allow the pump-cap to touch the surface of the eye.

Recommendations for the use of Ocutears® Hydro+ Eye drops, solution

Directions for use: 1-2 drops in each eye, one or more times a day, or as recommended by your doctor. Wash your hands thoroughly before opening the bottle.

- 1) Before opening for the first time, pull on the tamper evident strip and remove it completely.
- 2) Put 1 drop into the eye by lightly pressing on the bottle.
- 3) Close the bottle well after use.

ATTENTION

When used simultaneously with other ophthalmic drugs, the time interval between the instillation of each subsequent eye drops should be at least 5 minutes. It is recommended to use Ocutears® Hydro+ Eye drops, solution last.

Key Accessories, Medical Devices or Non-Medical Devices, to Be Used in Combination With Medical Device Submitted For Registration

Not applicable for this medical device.

Main technical and functional characteristics

Tests	Specifications	Conducted during validation / release
Exterior	Transparent, slightly viscous, colorless, odorless solution	V./R.
pH	6.90–7.50	V./R.
Osmolality	270–330 mOsm/kg	V./R.
Viscosity (at +20 °C, measured on manufactured product)	6.00–14.00 cSt	V./R.
Sterility	Sterile	V./R.

Operating conditions

The product is used in the internal environment of the human body.

Storage and transportation conditions

- Store vials in consumer packaging at 2-25 °C at relative humidity range up to 80% non-condensing (at 25 °C).
- Do not freeze.
- Protect from heat, sunlight and moisture.
- Keep out of the reach of children.
- Products are intended for transportation under storage conditions.
- Do not use the product if the packaging is damaged.
- Do not use the product after the expiration date indicated on the package.
- The bottle should be used within 90 days after opening.

Package contents

Eye drops are packaged in a white linear low density polyethylene bottle, sealed with a white OSD cap-pump with a blue cap. During normal operation, the pump-cap cannot be removed from the bottle, only the cap is removed after removing the tamper evident strip. Drops in a bottle are packed in consumer packaging - a cardboard box.

Explanation of symbols used in labelling the medical product

The symbol	decoding of the symbol
	Manufacturer
	Store 2°C-25°C
	Keep away from sunlight
	Sterile via aseptic methods
	Batch number
	Expiry date
	Read the instruction for use

	Do not use if the product sterile barrier system or its packaging is compromised
	EC symbol
	Date of manufacture
	Catalog number
	Keep dry

Eye drops in the form Ocutears® Hydro+ Eye drops, solution contains: Cross-linked sodium hyaluronate 0.2%, Disodium tetraborate decahydrate, Sodium chloride, Boric acid, Water for injection

Pharmaceutical substances

Not applicable. The product does not contain any pharmaceutical substances

Information about Pharmaceutical substances

The medical device does not contain active pharmaceutical ingredients. The medical device contains substances that can be used in the manufacture of pharmaceutical products or act as medicines, however, these substances in the given concentrations do not affect the human body in addition to the effects of the medical device and do not have pharmacological, immunological, genetic or metabolic effects on the human body.

Information about materials of animal/human origin

Not applicable. The product does not contain materials of animal/human origin.

Product Expiration / Shelf life

2 years. (24 months).

The product expiration after the first opening is 90 days

Description of sterilization method

STERILIZATION OF THE FIRST PACKAGE

Primary packaging material is sterilized by gamma radiation in a qualified external sterilization workshop. A sterilizing dose for bottle and cap-pump is 12.5 kGy.

COC periodically conducts dose audits and dose mapping if they meet the requirements of EN 552 and ISO 11137 and to verify the effectiveness of the sterilization dose applied.

STERILIZATION OF THE SOLUTION

The solution is double sterilized filtration under aseptic conditions through a suitable membrane with a pore size of 0.2 µm.

To ensure that the sterilizing filtration of the medical device is carried out correctly, safely and effectively, COC systematically checks the integrity of the filters used before and after filtration.

The aseptic process of manufacturing and filling the product into containers is validated by COC using the culture medium filling method. This is done by simulating the preparation and filling steps typically performed at COC using a specific culture broth and checking for sterility after filling the strips. This activity is carried out in accordance with the special procedure "Aseptic Filling Validation Standards". Validation protocols and reports are stored in COC's Q&A.

Disposal

The product that has not come into contact with biological fluids is subject to disposal and destruction in accordance with the requirements for handling medical waste as epidemiologically safe waste (Class A) in accordance with SanPiN 2.1.3684-21.

Information about the initial release or the last revision of the operational documentation

Operational documentation is issued for the first time.

Warranty

The manufacturer guarantees the quality and safety of the product during the shelf life, subject to the instructions for use, transportation and storage.

Reclamation

For additional information related to the use of the medical product, please contact an authorized representative of the manufacturer:

Limited Liability Company "SANTEN" (LLC "SANTEN")

105064, Moscow, Nizhny Susalny lane, 5, building p. 19, floor 2, space XII, room 16

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