20

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COOK MEDICAL EUROPE LTD.

O'Halloran Road,

National Technology Park

Limerick, Ireland.

Retrievers-baskets Reusable

Company Registration No. 490000

INSTRUCTION FOR USE VAT Registration No. 9776237X

Retrievers-baskets Reusable used to remove polyps, foreign bodies, stents and prostheses of the bile ducts and gastrointestinal tract.

Retrievers-baskets (MPB-2X4-8; MPB-3X6-8)

INTENDED USE

This device is used for the endoscopic removal of polyps after polypectomy.

NOTES

This device is supplied non-sterile and reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in conjunction with removal of polyps.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to medication, respiratory depression or arrest, cardiac arrhythmia or arrest.

Others associated with basket removal of a polyp include, but are not limited to: localized inflammation.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

INSTRUCTIONS FOR USE

1. With basket fully retracted into sheath, insert device into endoscope accessory channel. Advance device through channel in short increments until basket sheath exits endoscope.
2. Endoscopically visualize polyp to be removed and advance basket to appropriate position. Place sheath of basket beyond target, then extend basket to entrap polyp. **Note:** It may be necessary to move basket back and forth along target in order to capture polyp.
3. Withdraw basket into sheath while maintaining entrapment. Support polyp against endoscope tip, and withdraw endoscope from patient. Maintain endoscopic visualization to ensure continued entrapment of polyp. While withdrawing device from endoscope, wipe secretions from outer sheath.
4. Prepare extracted specimen per institutional guidelines.
5. Upon completion of procedure, clean and prepare device for future use following cleaning instructions. During cleaning, inspect integrity and function of device to determine advisability of reuse. **If kinks, bends or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.**

CLEANING INSTRUCTIONS

1. Immediately after use, attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaning solution to Luer lock fitting on device, then irrigate.
2. Extend basket, then submerge solution-filled device in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
3. Beginning at handle and working toward tip, wipe device with a textured cloth soaked in enzymatic solution.
4. Irrigate inner sheath by attaching a syringe containing a minimum of 20 cc of enzymatic solution to Luer lock fitting on handle of device. Repeat this step until all visible secretions are removed.
5. Gently brush basket with enzymatic solution and a soft nylon bristle brush to remove all foreign matter.
6. Remove device from enzymatic solution, then thoroughly rinse inner sheath using a syringe filled with clean water. Rinse handle, outer sheath and extended basket thoroughly with clean, running water.

2. Block.

7. Ultrasonically clean device, with basket extended, for minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.
8. Remove device from ultrasonic cleaner, then thoroughly rinse inner sheath with a syringe containing a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse external portions of device and basket with clean, running water.
9. Gently wipe external portions of device dry with a lint-free cloth. Air dry lumen and basket by hanging device, with basket extended, or use medical grade (oil-free) forced air through Luer lock fitting on handle.
10. Package device for sterilization following AAMI recommended practices.
11. Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4°C (130°F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4°C (130°F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space.

AUTOCLAVE CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	132.2°C (270°F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

INSTRUCTION FOR USE

Retrievers-baskets (MSB-1.5X3.5; MSB-2.5X5; MSB-2X4; MSB-2X4-6; MSB-3X6; MSB-3X6-6; MSB-1.5X3.5-F; MSB-2.5X5-F; MSB-2X4-F; MSB-3X6-F)

INTENDED USE

This device is used for the endoscopic removal of biliary stones and foreign bodies.

NOTES

Select Memory Soft Wire Baskets are compatible with the Soehendra Lithotripter.

This device is reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to the primary endoscopic procedure to be performed in conjunction with the removal of the foreign body.

Perforation.

If this device is to be used for the removal of biliary stones, contraindications include an ampullary opening inadequate to allow for the unimpeded passage of the stone and basket.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy and ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction, respiratory depression or arrest, cardiac arrhythmia or arrest.

Others associated with basket extraction include, but are not limited to: impaction of the object, aspiration of foreign body, localized inflammation, pressure necrosis

PRECAUTIONS

Refer to package label for minimum channel size required for this device and to determine basket compatibility.

If this device is to be used for removal of biliary stones, the following caution applies: Assessment of stone size and ampullary orifice must be made to determine necessity of sphincterotomy. If sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

WARNINGS

When retrieving a foreign body from upper GI tract, extreme caution must be exercised to avoid occlusion of airway.

Surgical intervention may be required if stone impaction and/or basket fragmentation occurs.

If a stone cannot be removed endoscopically with this basket, the Soehendra Lithotripter may be used with select Memory Soft Wire Baskets (see package label) to mechanically crush stone and aid in removal. Due to mechanical pressure generated with Soehendra Lithotripter, basket fragmentation and/or stone impaction in common bile duct may occur and require surgical intervention. Risk of basket fragmentation or stone impaction must be weighed against potential benefit of using lithotripter.

INSTRUCTIONS FOR USE

1. With basket fully retracted into sheath, insert device into endoscope accessory channel. Advance device through channel in short increments until basket sheath exits endoscope. Note: When using a basket that is removable from sheath, basket may be removed to facilitate use of a wire guide. When sheath is in proper position, remove wire guide, then reintroduce basket through sheath.
2. Endoscopically or fluoroscopically visualize extraction target and advance basket to appropriate position. Fluoroscopic visualization may be enhanced by injecting contrast via side arm fitting. Place sheath of basket beyond target before extending basket to entrap extraction target. Note: Tighten connector cap on handle before injection of contrast. (See fig. 1)
3. Confirm desired position of basket sheath relative to target. With one hand holding side arm fitting, gently push forward on pin vise handle with other hand to advance basket out of sheath. Caution: Pulling on sheath while advancing or retracting basket may damage device, rendering it inoperable.
4. Place basket around extraction target or foreign body and slowly retract handle until entrapment is achieved. Note: It may be necessary to move basket back and forth along target in order to capture object. Note: Connector cap on handle may be loosened in order to more freely manipulate basket. (See fig. 1)
5. Withdraw basket into channel until foreign body is supported against endoscope tip, then withdraw endoscope from patient. Maintain endoscopic visualization to ensure continued entrapment of object. While withdrawing device from endoscope, wipe secretions from outer sheath. Caution: If difficulty is encountered during removal of basket from duct, moderate force may be applied by pulling on handle to manually fracture stone. If this proves unsuccessful, a compatible basket may be used with Soehendra Lithotripter to mechanically crush stone. If passage is still restricted, surgical intervention may be necessary.
6. Prepare or dispose of extracted specimen per institutional guidelines.
7. Upon completion of procedure, clean and prepare device for future use following cleaning instructions below. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, bends, or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.

CLEANING INSTRUCTIONS

Baskets Not Removable From Sheath

1. Immediately after use, attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting on device, then irrigate.
2. Extend basket, then submerge solution-filled device in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
3. Beginning at handle and working toward tip, wipe device with a textured cloth soaked in enzymatic solution.
4. Irrigate inner sheath by attaching a syringe containing a minimum of 20 cc of enzymatic solution to Luer lock fitting on handle of device. Repeat this step until all visible secretions are removed.
5. Gently brush basket with enzymatic solution and a soft nylon bristle brush to remove all foreign matter.
6. Remove device (or components, if applicable) from enzymatic solution, then thoroughly rinse inner sheath using a syringe filled with clean water. Rinse handle, outer sheath and extended (or separate) basket thoroughly with clean, running water.

7. Ultrasonically clean device (or components, if applicable) with basket extended, for a minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.
8. Remove device (or components, if applicable) from ultrasonic cleaner, thoroughly rinse inner sheath with a syringe containing a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse external portions of device and basket with clean, running water.
9. Gently wipe external portions of device dry with a lint-free cloth. Hang device with basket extended to air dry lumen and basket, or use medical grade (oil-free) forced air through Luer lock fitting on handle.
10. Package device for sterilization following AAMI recommended practices.
11. Sterilize device using set points for desired sterilization method shown in tables.

Baskets Removable From Sheath

1. Immediately after use, disconnect handle from catheter. While withdrawing basket, gently wipe secretions from wires. Attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting on catheter, then irrigate.
2. Submerge solution-filled catheter and basket in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.

REFER TO STEPS 3-9 OF "BASKETS NOT REMOVABLE FROM SHEATH" ABOVE, THEN CONTINUE WITH STEP 3 BELOW:

3. Gently re-assemble basket and sheath.
4. Package device for sterilization following AAMI recommended practices.
5. Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4°C (130°F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4°C (130°F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

1 To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

2 Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.


Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilization, the product should be taken out of service immediately.

AUTOClave CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	132.2°C (270°F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

Connector Cap

INSTRUCTION FOR USE

Retrievers-baskets (MWB-1.5X3.5; MWB-2.5X5; MWB-2X4; MWB-3X6).

INTENDED USE

This device is used for the endoscopic removal of biliary stones and foreign bodies.

NOTES

This device is reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, **do not use**. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. **Please** notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in conjunction with removal of foreign body.

Perforation.

If device is to be used for biliary stone extraction, contraindications include an ampullary opening inadequate to allow for unimpeded passage of stone and basket.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy and ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction, respiratory depression or arrest, cardiac arrhythmia or arrest.

Others associated with basket extraction include, but are not limited to: impaction of object, aspiration of foreign body, localized inflammation, pressure necrosis.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

If device is to be used for removal of biliary stones, assessment of stone size and ampullary orifice must be made to determine necessity of sphincterotomy. If sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

WARNINGS

This device is not compatible with the SoehendraR Lithotripter or any other mechanical lithotripter.

When retrieving a foreign body from the upper GI tract, extreme caution must be exercised to avoid occlusion of airway.

Surgical intervention may be required if impaction occurs.

INSTRUCTIONS FOR USE

1. **With basket fully retracted into sheath, insert device into endoscope accessory channel. Advance device through channel in short increments until basket sheath exits endoscope.**
2. **Endoscopically or fluoroscopically visualize extraction target and advance basket to appropriate position.** Fluoroscopic visualization may be enhanced by injecting contrast via side arm fitting. Place sheath of basket beyond target before extending basket to entrap extraction target.
3. **After confirming desired position of basket sheath relative to target, advance basket out of sheath by pushing forward on handle. Place basket around extraction target and slowly retract handle until entrapment is achieved.**

Note: It may be necessary to move basket back and forth along target in order to capture object.

4. Withdraw basket into channel while maintaining entrapment. Support foreign body or biliary stone against endoscope tip, and withdraw endoscope from patient. Maintain endoscopic visualization to ensure continued entrapment of object. While withdrawing device from endoscope, wipe secretions from outer sheath. **Caution:** If difficulty is encountered when removing basket from duct, moderate force may be applied by pulling on handle to manually fracture stone. If passage is still restricted, surgical intervention may be necessary.

5. Prepare or dispose of extracted specimen per institutional guidelines.

6. Upon completion of procedure, clean and prepare device for future use following cleaning instructions below. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, bends, or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.

CLEANING INSTRUCTIONS

1. Immediately after use, attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting on device, then irrigate.

2. Extend basket, then submerge solution-filled device in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.

3. Beginning at handle and working toward tip, wipe device with a textured cloth soaked in enzymatic solution.

4. Irrigate inner sheath by attaching a syringe containing a minimum of 20 cc of enzymatic solution to Luer lock fitting on handle of device. Repeat this step until all visible secretions are removed.

5. Gently brush basket with enzymatic solution and a soft nylon bristle brush to remove all foreign matter.

6. Remove device from enzymatic solution, then thoroughly rinse inner sheath using a syringe filled with clean water. Rinse handle, outer sheath and extended basket thoroughly with clean, running water.

7. Ultrasonically clean device, with basket extended, for a minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.

8. Remove device from ultrasonic cleaner, then thoroughly rinse inner sheath with a syringe containing a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse external portions of device and basket with clean, running water.

9. Gently wipe external portions of device dry with a lint-free cloth. Air dry lumen and basket by hanging device, with basket extended, for three hours, or by using medical grade (oil free) forced air through Luer lock fitting on handle.

10. Package device for sterilization following AAMI recommended practices.

11. Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER

Prevacuum Pressure

Relative Humidity

Humidity Dwell

Chamber Temperature

Exposure Pressure

Exposure Dwell

ETO Concentration

Post Vacuums

Post Vacuum Pressure

Air Wash

Aeration Temperature

Aeration Dwell

Air Changes During Aeration

SET POINT

2.0 psia (25.9" Hg Vac)

60% RH

20 minutes

54.4°C (130°F)

To Be Determined

2 hours

600 mg/L

3

2.0 psia (25.9" Hg Vac)

2.5 minutes

54.4°C (130°F)

12 hours

90 per hour

1 To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

2 Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilizing, it should be taken out of service immediately.

AUTOClave CYCLE PARAMETERS

PARAMETER

Cycle Type

Prevacuum Pressure

Chamber Exposure Temperature

Chamber Exposure Pressure

Exposure Dwell

SET POINT

Prevacuum

2.5 psia (24.8" Hg Vac)

132.2°C (270°F)

41.8 psia (27.1 psig)

5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Wilson Cook Medical, Inc. manufacture all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

Wilson Cook Medical, Inc. and Cook Incorporated facility is *also certified to ISO 13485 (2003) Quality Management System*.

SHELF LIFE

The "use by" date as denoted by symbols on the labelling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

The Shelf-life of all Retrievers – baskets reusable, which fall under the scope of this Technical Data Summary – is three (3) years.

STERILIZATION DETAILS

Retrievers-baskets (MSB-1.5X3.5; MSB-2.5X5; MSB-2X4; MSB-2X4-6; MSB-3X6; MSB-3X6-6; MSB-1.5X3.5-F; MSB-2.5X5-F; MSB-2X4-F; MSB-3X6-F; MWB-1.5X3.5; MWB-2.5X5; MWB-2X4; MWB-3X6) are supplied sterile in peel-open packages. To guarantee sterility, an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} applies. Reusable.

Retrievers-baskets (MPB-2X4-8; MPB-3X6-8) are non-sterile and reusable.

STORAGE RULES

These devices should be stored in a dry location, away from temperature extremes.

TRANSPORTATION RULES

Transportation of Retrievers – baskets reusable shall account for the type of transportation employed, the destination, the vendor used for transportation.

SUPPLY DETAILS

Retrievers – baskets reusable (MPB-2X4-8; MPB-3X6-8; MSB-1.5X3.5; MSB-2.5X5; MSB-2X4; MSB-2X4-6; MSB-3X6; MSB-3X6-6; MSB-1.5X3.5-F; MSB-2.5X5-F; MSB-2X4-F; MSB-3X6-F; MWB-1.5X3.5; MWB-2.5X5; MWB-2X4; MWB-3X6) are supplied in individual packages. Individual packages are put in the corrugated board boxes. Boxes are being taped with the control tape to cover the sides or along the perimeter.

Terms of use

Retrievers – baskets reusable used in medical institutions.

LABELING



Retrievers – baskets reusable are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

- company-manufacturer and manufacturing address;
- device name (trade mark);
- serial number;
- marking about CE certification;
- marking about sterilization;
- address of the EU representative;
- manufacturer's trade mark
- manufacturing date;
- sizes;
- graphic picture;
- lot number;
- marking "use with handle" and "do not use with handle"
- marking for "Attention. Consult accompanied documentation"
- shelf life.

MSB-1.5X3.5
MEMORY SOFT WIRE BASKET

G22136

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste. There is no environmentally harmful effect as the result of the device disposal in accordance with the relevant national guidelines.

Regulatory Affairs Manager

COOK MEDICAL EUROPE LTD.

O'Halloran Road,
National Technology Park
Limerick, Ireland.

Company Registration No. 4021

VAT Registration No. 977621

COOK MEDICAL EUROPE LTD
O'Halloran Road,
National Technology Park,
Limerick, Ireland.

Company Registration No. 483272
VAT Registration No. 977623-01

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