80

CALIFORNIA ALL-PURPOSE ACKNOWLEDGMENT

CIVIL CODE § 1189

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.


State of California)
County of Riverside)

On 05 OCT 2022 before me, Steve Oberlander, Notary Public
Date Here Insert Name and Title of the Officer
personally appeared Nishitaben S. Patel
Name(s) of Signer(s)

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature Steve Oberlander
Signature of Notary Public

Place Notary Seal Above

OPTIONAL

Though this section is optional, completing this information can deter alteration of the document or fraudulent reattachment of this form to an unintended document.

Description of Attached Document

Title or Type of Document: IFU - Solitaire X (max)
Document Date: _____ Number of Pages: _____
Signer(s) Other Than Named Above: _____

Capacity(ies) Claimed by Signer(s)

Signer's Name: _____
☐ Corporate Officer — Title(s): _____
☐ Partner — ☐ Limited ☐ General
☐ Individual ☐ Attorney in Fact
☐ Trustee ☐ Guardian or Conservator
☐ Other: _____
Signer Is Representing: _____

Signer's Name: _____
☐ Corporate Officer — Title(s): _____
☐ Partner — ☐ Limited ☐ General
☐ Individual ☐ Attorney in Fact
☐ Trustee ☐ Guardian or Conservator
☐ Other: _____
Signer Is Representing: _____

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Solitaire™ X Revascularization Device

CAUTION


- The Solitaire™ X Revascularization Device should only be used by physicians trained in interventional neuroradiology and treatment of ischemic stroke.

DESCRIPTION

The Solitaire™ X Revascularization Device is designed to restore blood flow in patients experiencing ischemic stroke due to large intracranial vessel occlusion. The device is designed for use in the neurovasculature such as the internal carotid artery, M1 and M2 segments of the middle cerebral artery, basilar, and the vertebral arteries.

Figure 1:

Solitaire™ X Revascularization Device

- | | |
|-----------------------|---|
| 1. Proximal marker | A. Stent diameter |
| 2. Distal markers | B. Usable length |
| 3. Introducer sheath | C. Stent length |
| 4. Fluoroscope Marker | D. Length from distal tip to fluoroscope marker |
| 5. Push wire | E. Marker to marker length |
| | F. Push Wire Length |

DEVICE COMPATIBILITY

The following parts are required to use the Solitaire™ X Revascularization Device:

- All Solitaire™ X Revascularization Devices should be introduced only through a microcatheter with an inside diameter of 0.017 inches (0.43 mm)-0.027 inches (0.69 mm). Refer to Table 1 for recommended microcatheter sizing guidelines.

Compatibility testing for all Solitaire™ X devices has been performed with Phenom™ 21/27, Rebar™ 18, and Marksman™ microcatheters.

Refer to the instructions supplied with all interventional devices and materials to be used in conjunction with the Solitaire™ X Revascularization Device for their intended uses, contraindications and potential complications.

Other accessories for performing a procedure are NOT supplied; to be selected based on the physician's experience and preferences:

- Minimum 5F (1.67 mm) guide catheter
(NOTE: Users should select a guide catheter with appropriate support to deliver interventional devices. For a guide catheter a minimum inside diameter 0.061 inches (1.55 mm) is recommended. For a balloon guide catheter ensure to use 8F (2.67 mm)- 9F (3.00 mm) with minimum inside diameter 0.075 inches (1.91 mm))
- Microcatheter (refer to Table 1)
- Guidewire
- Aspiration Device (60 cc syringe or aspiration pump/system)
- Saline solution/heparin-saline continuous flush set
- Rotating Hemostasis Valve (RHV)
- Infusion stand
- Femoral arterial lock

INTENDED PURPOSE / INDICATIONS FOR USE

The Solitaire™ X Revascularization Device is designed for use in the flow restoration of patients with ischemic stroke due to large intracranial vessel occlusion. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

CONTRAINDICATIONS

Use of the Solitaire™ X Revascularization Device is contraindicated under these circumstances.

- Patients with known hypersensitivity to nickel-titanium.
- Patients with stenosis and/or pre-existing stent proximal to the thrombus site that may preclude safe recovery of the Solitaire™ X Revascularization Device.
- Patients with angiographic evidence of carotid dissection.

PREPARATIONS FOR USE

- Administer anti-coagulation and anti-platelet medications per standard institutional guidelines.
- Aided by angiographic radiography, determine the location and size of the area to be revascularized.
- Select a Solitaire™ X Revascularization Device per Table 1.
- To achieve optimal performance of the Solitaire™ X Revascularization Device and to reduce the risk of thromboembolic complications, maintain continuous flushing action between a) the femoral arterial sheath and the guide catheter, b) the microcatheter and the guide catheter and c) the microcatheter and the push wire and the Solitaire™ X Revascularization Device. Check all connections to make sure that during the continuous flush no air enters the guide catheter or the microcatheter.
- Position a suitable guide catheter as close to thrombus site as possible employing a standard method. The guide catheter should be appropriately sized to retrieve clot if so desired in subsequent steps. Connect a RHV to the fitting of the guide catheter, and then connect a tube to the continuous flush.
- With the aid of Table 1, select a microcatheter suitable for advancing the Solitaire™ X Revascularization Device.
- Connect a second RHV to the fitting of the microcatheter and then connect a tube to the continuous flush.
- Set the flush rate per standard institutional guidelines.
- With the aid of a suitable guide wire, advance the microcatheter until the end of the microcatheter is positioned distal to the thrombus so that the usable length portion of the Solitaire™ X Revascularization Device will extend past each side of the thrombus in the vessel when fully deployed. Tighten the RHV around the microcatheter. Remove guide wire.

DIRECTIONS FOR USE

Delivering the Solitaire™ X Revascularization Device

- Insert the distal end of the introducer sheath partially into the RHV connected to the microcatheter. Tighten the RHV and verify that fluid exits the proximal end of the introducer sheath.
- Loosen the RHV and advance the introducer sheath until it is firmly seated in the hub of the microcatheter. Tighten the RHV around the introducer sheath to prevent back flow of blood, but not so tight as to damage the Solitaire™ X Revascularization Device during its introduction into the microcatheter. Confirm that there are no air bubbles trapped anywhere in the system.
- Transfer the Solitaire™ X Revascularization Device into the microcatheter by advancing the push wire in a smooth, continuous manner. Once the flexible portion of the push wire has entered the microcatheter shaft, loosen the RHV and remove the introducer sheath over the proximal end of the push wire. Once completed, tighten the RHV around the push wire. Leaving the introducer sheath in place will interrupt normal infusion of flushing solution and allow back flow of blood into the microcatheter.
- Visually verify that the flushing solution is infusing normally. Once confirmed, loosen the RHV to advance the push wire.
- Once the Solitaire™ X Revascularization Device has been advanced to the fluoroscope marker band begin fluoroscopic imaging. With the aid of fluoroscopic monitoring, carefully advance the Solitaire™ X Revascularization Device until its distal markers line up at the end of the microcatheter. The Solitaire™ X Revascularization Device should be positioned such that the usable length portion of the device will extend past each side of the thrombus in the vessel when the device is fully deployed.

Table 1.
Solitaire™ X Revascularization Device:
Product Specifications and Recommended Sizing Guidelines

Model	Recommended Vessel Diameter ¹	Recommended Microcatheter Inner Diameter	Push Wire Length	Stent Diameter	Usable Length ²	Stent Length	Length from Distal tip to fluoroscope marker	Distal Radio-paque Markers	Proximal Radio-paque Markers	Radiopaque Stent Markers Spacing
SFR4-3-20-10	1.0 mm-3.0 mm	0.017" (0.43 mm)-0.027" (0.69 mm)	200 cm	3.0 mm	20.0 mm	30.6 mm	<150 cm	3	1	10 mm
SFR4-3-40-10	1.0 mm-3.0 mm	0.017" (0.43 mm)-0.027" (0.69 mm)	200 cm	3.0 mm	40.0 mm	51.6mm	<150 cm	3	1	10 mm

¹ Select a Solitaire™ X Revascularization Device based on the sizing recommendations in Table 1 and based on the smallest vessel diameter at thrombus site.

² Select a Solitaire™ X Revascularization Device usable length that is at least as long as the length of the thrombus.

WARNINGS


If excessive resistance is encountered during the delivery of the Solitaire™ X Revascularization Device, discontinue the delivery and identify the cause of the resistance. Advancement of the Solitaire™ X Revascularization Device against resistance may result in device damage and/or patient injury.

Deploying the Solitaire™ X Revascularization Device

- Loosen the RHV around the microcatheter. To deploy the Solitaire™ X Revascularization Device, fix the push wire to maintain the position of the device while carefully withdrawing the microcatheter in the proximal direction.

Figure 2:

Solitaire™ X Revascularization Device Deployment

- Retract the microcatheter until it is just proximal to the proximal marker of the Solitaire™ X Revascularization Device. Tighten the RHV to prevent any movement of the push wire. The usable length of the deployed Solitaire™ device should extend past each side of the thrombus.
- Tighten the RHV around the microcatheter. Angiographically assess the revascularization status of the treated vessel.

Solitaire™ X Revascularization Device Re-Sheathing


If resheathing of the Solitaire™ X Revascularization Device is necessary (e.g. repositioning), follow these steps:

WARNINGS

- Advancing the microcatheter while the device is engaged in clot may lead to embolization of debris.
 - Do not advance the microcatheter against any resistance.
 - Do not reposition more than two times.
- Loosen the RHV around the microcatheter and around the push wire. With the aid of fluoroscopic monitoring, hold the push wire firmly in its position to prevent the Solitaire™ X Revascularization Device from moving.
 - Carefully re-sheath the Solitaire™ X Revascularization Device by advancing the microcatheter over the Solitaire™ X Revascularization Device until the distal markers of the Solitaire™ X Revascularization Device line up at the end of the microcatheter as illustrated in Figure 3 below.
- If significant resistance is encountered during the re-sheathing process, stop immediately and proceed to the section below entitled "Revascularization Device Recovery".

Figure 3:

Solitaire™ X Revascularization Device Re-sheathing

Revascularization Device Recovery

- If using balloon guide catheter, inflate guide catheter balloon to occlude vessel as specified in Balloon Guide Catheter labeling.
- Loosen the RHV around the microcatheter enough for movement while still maintaining a seal. If desired remove the microcatheter completely.
- To retrieve thrombus, withdraw the Solitaire™ X Revascularization Device to the guide catheter tip while applying aspiration to the guide catheter with the aspiration device. If applicable withdraw the microcatheter and the Solitaire™ X Revascularization Device as a unit. Never advance the deployed Solitaire™ X Revascularization Device distally.

- Apply aspiration to the guide catheter using the aspiration device and recover the Solitaire™ X Revascularization Device inside guide catheter. Continue aspirating guide catheter until the Solitaire™ X Revascularization Device is nearly withdrawn from the guide catheter.

NOTE: If withdrawal into the guide catheter is difficult, deflate balloon (if using balloon guide catheter) and then simultaneously withdraw guide catheter, microcatheter and Solitaire™ X Revascularization Device as a unit through the sheath while maintaining aspiration. Remove sheath if necessary.

WARNINGS

- If excessive resistance is encountered during recovery of the Solitaire™ X Revascularization Device, discontinue the recovery and identify the cause of the resistance.
 - For vessel safety, do not perform more than three recovery attempts in the same vessel using the Solitaire™ X Revascularization Device.
- Open the guide catheter RHV to allow the Solitaire™ X Revascularization Device to exit without resistance. Use care to avoid interaction with the site of the intervention and to prevent air from entering the system. If applicable allow the microcatheter and the Solitaire™ X Revascularization Device to exit without resistance.
 - Aspirate the guide catheter to ensure the guide catheter is clear of any thrombus material.
 - Deflate guide catheter balloon if using balloon guide catheter.
 - If additional flow restoration attempts are desired with:
 - A new Solitaire™ X Revascularization Device, then:**
 - Repeat the steps described above starting with the "Preparation for Use" section.
 - The same Solitaire™ X Revascularization Device, then:**
 - Clean the device with saline solution.
 - Note: Do not use solvents or autoclave.**
 - Carefully inspect the device for damage.
 - If no damage, load the introducer sheath over the proximal end of the push wire. For 3 mm devices, one end of the introducer sheath is tapered and should be (1) aligned with the distal end of the Solitaire™ X Revascularization Device and (2) interface with the microcatheter hub. If there is any damage, do not use the device and use a new Solitaire™ X Revascularization Device for subsequent flow restoration attempts following the steps described above starting with the "Preparation for Use" section. Use of a damaged device could result in additional device damage or patient injury.

WARNINGS

- For device safety, do not use each Solitaire™ X Revascularization Device for more than three flow restoration recoveries.

POTENTIAL COMPLICATIONS

Potential complications of the device and the endovascular procedure include or are synonymous with, but may not be limited to:

- Access site complications such as hematoma, pain, nerve damage
- Adverse reactions to antiplatelet/anticoagulation agents, contrast media, or anesthesia such as organ failure, hypertension, hypotension, hemorrhage
- Cardiac complications such as arrhythmia, myocardial infarction, angina
- Complications of radiation exposure such as alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia
- Device malfunctions such as friction, foreign body in patient, device separation, fracture
- Hematologic complications such as coagulopathy, hemolysis, intracranial hemorrhage
- Neurological deficit/dysfunctions such as cerebral infarction, transient ischemic attack, ischemia
- Systemic complications such as infection, fever, inflammation, edema, hypersensitivity, foreign body reaction, toxicity, shock
- Therapeutic Response Decreased
- Vascular complications such as dissection, perforation, rupture, embolism (including distal embolization to previously uninvolved territory), thromboembolism, thrombosis, fistula, necrosis, stenosis, vasoconstriction, occlusion
- Vision complications such as impaired vision and vision loss
- Death

*Consult instructions for use for other therapy devices and medications for additional potential complication information. If a serious incident related to the device occurs, contact your Medtronic representative and the competent authority in your respective country/region.

WARNINGS

- The appropriate anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.
- Administer IV t-PA as soon as possible for all patients who are indicated to receive the drug. Do not cause delays in this therapy.
- Do not torque the Solitaire™ X Revascularization Device.
- For vessel safety, do not perform more than three recovery attempts in the same vessel using Solitaire™ X Revascularization Devices.
- For device safety, do not use each Solitaire™ X Revascularization Device for more than three flow restoration recoveries.
- For each new Solitaire™ X Revascularization Device, use a new microcatheter.
- Solitaire™ X Revascularization Device does not allow for electrolytic detachment.

WARNINGS

To prevent device separation:

- Do not oversize device.
- Do not recover (i.e. pull back) the device when encountering excessive resistance. Instead, resheath the device with the microcatheter and then remove the entire system under aspiration. If resistance is encountered during resheathing, discontinue and remove the entire system under aspiration.
- Do not treat patients with known stenosis proximal to the thrombus site.

This device is intended for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

Visually inspect all sterile barrier systems, that are labeled as sterile, immediately prior to use. Do not use the device if breaches in sterile barrier system integrity are evident.

Nickel (CAS 7440-02-0) is present in the Solitaire™ X Revascularization Device, and may cause sensitization or an allergic reaction.

For additional Materials of Concerns information such as REACH, CA Prop 65 or other product stewardship programs, go to www.medtronic.com/productstewardship

PRECAUTIONS

- Operators should take all necessary precautions to limit X-ray radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors whenever possible.
- Carefully inspect the sterile package and the Solitaire™ X Revascularization Device prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components.
- The Solitaire™ X Revascularization Device is not to be used after the expiration date imprinted on the product label.
- Refer to the appropriate intravenous tissue plasminogen activator (IV t-PA) manufacturer labeling for indications, contraindications, warnings, precautions, and instructions for use.
- Initiate mechanical thrombectomy treatment as soon as possible.

HOW SUPPLIED

This device is supplied STERILE using ethylene oxide. This device is Non-pyrogenic.

STORAGE AND DISPOSAL

- This device should be stored in a dry place, away from sunlight.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

SAFETY AND CLINICAL PERFORMANCE

- Expected clinical benefit - The Solitaire™ X Revascularization Device is designed for use in the flow restoration of patients with ischemic stroke due to large intracranial vessel occlusion to attain the following clinical benefits:
 - removal of the thrombus from the occluded intracranial vessel segments
 - restoration of blood flow in the affected vessels
- Performance characteristics of the device - The Solitaire™ X Revascularization Device is designed to restore blood flow in the neurovasculature such as the internal carotid artery, M1 and M2 segments of the middle cerebral artery, basilar, and the vertebral arteries. The Solitaire™ X Revascularization Device removes the thrombus from the occluded intracranial vessel segments and restores blood flow in the affected vessels. Solitaire™ X Revascularization Device is intended to treat ischemic stroke due to large intracranial vessel occlusion.
- A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, the manufacturer has no control over the conditions under which this product is used. The manufacturer therefore disclaims all warranties, both expressed and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The manufacturer shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind the manufacturer to any representation or warranty with respect to the product. The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

This page is intentionally left blank.

SYMBOL GLOSSARY

	Sterilized using ethylene oxide		Authorized representative in the European Community
	Single Sterile Barrier System		Catalogue number
	Do not re-use		Manufacturer
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician		Use-by date
	Do not re-sterilize		Batch code
 www.medtronic.com/manuals	Consult instructions for use at this website		Contents of Package
	Do not use if package is damaged		Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union Acts.
	Caution		Medical device
	Non-pyrogenic		Importer
	Keep away from sunlight		Inner diameter (compatible micro catheter)
	Keep dry		

Medtronic, Medtronic with logo and Medtronic logo are U.S. and internationally registered trademarks of Medtronic.

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
Figure 1. Labelled Schematic of Solitaire™ X Revascularization Device

- | | |
|----------------------|--|
| 1. Proximal marker | A. Stent diameter |
| 2. Distal markers | B. Usable length |
| 3. Introducer sheath | C. Stent length |
| 4. Push wire | D. Length from distal tip to fluorosafe marker |
| 5. Fluorosafe marker | E. Marker to marker length |
| | F. Push Wire Length |

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Indications for use

- The Solitaire X Revascularization Device is designed for use in the flow restoration of patients with ischemic stroke due to large intracranial vessel occlusion. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

Contraindications

- Use of the Solitaire X Revascularization Device is contraindicated under these circumstances.
- Patients with known hypersensitivity to nickel-titanium.
- Patients with stenosis and/or pre-existing stent proximal to

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calibration and other actions necessary for putting a medical device into operation

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no air enters the guide catheter or the microcatheter.

5. Position a suitable guide catheter as close to thrombus site as possible employing a standard method. The guide catheter should be appropriately sized to retrieve clot if so desired in subsequent steps. Connect a RHV to the fitting of the guide catheter, and then connect a tube to the continuous flush.
6. With the aid of Table 1, select a microcatheter suitable for advancing the Solitaire™ X Revascularization Device.
7. Connect a second RHV to the fitting of the microcatheter and then connect a tube to the continuous flush.
8. Set the flush rate per standard institutional guidelines.
9. With the aid of a suitable guide wire, advance the microcatheter until the end of the microcatheter is positioned distal to the thrombus so that the usable length portion of the Solitaire™ X Revascularization Device will extend past each side of the thrombus in the vessel when fully deployed. Tighten the RHV around the microcatheter. Remove guide wire.

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1. Insert the distal end of the introducer sheath partially into the RHV connected to the microcatheter. Tighten the RHV and verify that fluid exits the proximal end of the introducer sheath.
2. Loosen the RHV and advance the introducer sheath until it is firmly seated in the hub of the

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microcatheter. Tighten the RHV around the introducer sheath to prevent back flow of blood, but not so tight as to damage the Solitaire™ X Revascularization Device during its introduction into the microcatheter. Confirm that there are no air bubbles trapped anywhere in the system.

3. Transfer the Solitaire™ X Revascularization Device into the microcatheter by advancing the push wire in a smooth, continuous manner. Once the flexible portion of the push wire has entered the microcatheter shaft, loosen the RHV and remove the introducer sheath over the proximal end of the push wire. Once completed, tighten the RHV around the push wire. Leaving the introducer sheath in place will interrupt normal infusion of flushing solution and allow back flow of blood into the microcatheter.
4. Visually verify that the flushing solution is infusing normally. Once confirmed, loosen the RHV to advance the push wire.
5. Once the Solitaire™ X Revascularization Device has been advanced to the fluorosafe marker band begin fluoroscopic imaging. With the aid of fluoroscopic monitoring, carefully advance the Solitaire™ X Revascularization Device until its distal markers line up at the end of the microcatheter. The Solitaire™ X Revascularization Device should be positioned such that the usable length portion of the device will extend past each side of the thrombus in the vessel when the device is fully deployed.


WARNINGS

- If excessive resistance is encountered during the delivery of the Solitaire™ X Revascularization Device, discontinue the delivery and identify the cause of the resistance. Advancement of the Solitaire™ X Revascularization Device against resistance may result in device damage and/or patient injury.

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
6. Loosen the RHV around the microcatheter. To deploy the Solitaire™ X Revascularization Device, fix the push wire to maintain the position of the device while carefully withdrawing the microcatheter in the proximal direction.

Figure 1:

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7. Retract the microcatheter until it is just proximal to the proximal marker of the Solitaire™ X Revascularization Device. Tighten the RHV to prevent any movement of the push wire. The usable length of the deployed Solitaire™ device should extend past each side of the thrombus.
8. Tighten the RHV around the microcatheter. Angiographically assess the revascularization status of the treated vessel.

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WARNINGS

• Advancing the microcatheter while the device is engaged in clot may lead to embolization of debris.

• Do not advance the microcatheter against any resistance.


• Do not reposition more than two times.

9. Loosen the RHV around the microcatheter and around the push wire. With the aid of fluoroscopic monitoring, hold the push wire firmly in its position to prevent the Solitaire™ X Revascularization Device from moving.

10. Carefully re-sheath the Solitaire™ X Revascularization Device by advancing the microcatheter over the Solitaire™ X Revascularization Device until the distal markers of the Solitaire™ X Revascularization Device line up at the end of the microcatheter as illustrated in Figure 2 below.


If significant resistance is encountered during the re-sheathing process, stop immediately and proceed to the section below entitled "Revascularization Device Recovery".

Figure 2:
Solitaire™ X Revascularization Device Re-sheathing

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11. If using balloon guide catheter, inflate guide catheter balloon to occlude vessel as specified in Balloon Guide Catheter labeling.
 12. Loosen the RHV around the microcatheter enough for movement while still maintaining a seal. If desired remove the microcatheter completely.
 13. To retrieve thrombus, withdraw the Solitaire™ X Revascularization Device to the guide catheter tip while applying aspiration to the guide catheter with the aspiration device. If applicable withdraw the microcatheter and the Solitaire™ X Revascularization Device as a unit. Never advance the deployed Solitaire™ X Revascularization Device distally.
 14. Apply aspiration to the guide catheter using the aspiration device and recover the Solitaire™ X Revascularization Device inside guide catheter. Continue aspirating guide catheter until the Solitaire™ X Revascularization Device is nearly withdrawn from the guide catheter.
- NOTE: If withdrawal into the guide catheter is difficult, deflate balloon (if using balloon guide catheter) and then simultaneously withdraw guide catheter, microcatheter and Solitaire™ X Revascularization Device as a unit through the sheath while maintaining aspiration. Remove sheath if necessary.

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WARNINGS

- If excessive resistance is encountered during recovery of the Solitaire™ X Revascularization Device, discontinue the recovery and identify the cause of the resistance.
- For vessel safety, do not perform more than three recovery attempts in the same vessel using the Solitaire™ X Revascularization Device.
- 15. Open the guide catheter RHV to allow the Solitaire™ X Revascularization Device to exit without resistance. Use care to avoid interaction with the site of the intervention and to prevent air from entering the system. If applicable allow the microcatheter and the Solitaire™ X Revascularization Device to exit without resistance.
- 16. Aspirate the guide catheter to ensure the guide catheter is clear of any thrombus material.
- 17. Deflate guide catheter balloon if using balloon guide catheter.
- 18. If additional flow restoration attempts are desired with:
 - a. A new Solitaire™ X Revascularization Device, then:
 - i. Repeat the steps described above starting with the "Preparation for Use" section.
 - b. The same Solitaire™ X Revascularization Device, then:
 - i. Clean the device with saline solution.

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Symbols

Symbol	Description
STERILE EO	Sterilized using ethylene oxide
	Do not re-use
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Do not re-sterilize
	Single sterill barrier system
	Inner diameter (compatible micro catheter)
	Consult instructions for use
	Caution
	Non-pyrogenic
	Keep away from sunlight
	Keep dry
REF	Catalogue number

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© 2021 Medtronic.