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«I certify»

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Managing Director
Botiss Biomaterials GmbH
Dr. Drazen Tadic

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Addresses of the places of production of medical products:

Cells+Tissuebank Austria gemeinnützige GmbH, Magnesitstr. 1, AT- 3500 Krems
Tel.: +43 (0) 2732 76954 0
Fax: +43 (0) 2732 76954 40
www.ctba.at

Composition

Cancellous, cortical, or cortico-cancellous portions of human bone tissue

Product description

Bone transplant of human origin, lyophilized, in cancellous, cortical or corticocancellous form as blocks, granules, struts, rings and cylinders.

Intended use

maxgraft[®] is intended to replace missing or inadequate bone tissue, for filling and stabilizing bone defects. maxgraft[®] should only be used by doctors who are familiar with surgical intervention techniques in dental, oral, and maxillofacial surgery.

Indications

maxgraft[®] granules:

- Regeneration of periodontal osseous defects
- Regeneration after cysts and root tip resection
- Regeneration of extraction sockets (socket preservation)
- Regeneration of missing bone tissue around dental implants
- Regeneration of gaps around block grafts
- Horizontal augmentation of the alveolar ridge
- Sinus augmentation
- Three-dimensional (horizontal and/or vertical) augmentation of the alveolar ridge

maxgraft[®] blocks:

- Horizontal augmentation of the alveolar ridge
- Sinus augmentation
- Three-dimensional (horizontal and/or vertical) augmentation of the alveolar ridge

maxgraft[®] cortico:

- Horizontal augmentation of the alveolar ridge
- Three-dimensional (horizontal and/or vertical) augmentation of the alveolar ridge

maxgraft[®] bonering:

- Horizontal augmentation of the alveolar ridge
- Sinus augmentation
- Three-dimensional (horizontal and/or vertical) augmentation of the alveolar ridge

74337
18.10.2021

For Russia

maxgraft® bonebuilder:

- Horizontal augmentation of the alveolar ridge
- Sinus augmentation
- Three-dimensional (horizontal and/or vertical) augmentation of the alveolar ridge

Contraindications

Necrosis in the oral cavity is a contraindication for maxgraft®. Following the complete removal of infected bone tissue, the use of maxgraft® should be considered carefully due to a limited healing rate.

Side effects

The evaluation of side effects is based on the following frequencies: very common ($\geq 1/10$); common ($\geq 1/100$, $< 1/10$); occasional ($\geq 1/1,000$, $< 1/100$); rare ($\geq 1/10,000$, $< 1/1,000$); very rare ($< 1/10,000$); unknown (frequency cannot be estimated from the available data).

As with all bone implants, the use of surgery makes it difficult to assess whether side effects, such as infections, complications during healing, or loosening of the implant are caused by the product or by the surgery.

A retrospective analysis demonstrated the safety of human lyophilized allografts. No microbiological or viral infections were caused by the allograft. Side effects known to be due to the surgery include haematomas and pain. Donors of bone material that is processed into maxgraft® are subject to strict selection criteria with respect to donor selection and serological testing. A validated inactivation procedure is used for the processing. In spite of these extensive safety measures, it is not possible to exclude the presence of unknown or undetectable pathogens.

Use during pregnancy and lactation

There are no adequate data on use in pregnancy and lactation. The risks that are generally linked to an operation including drug administration, blood loss etc., must be especially taken into careful consideration for pregnant or lactating patients. Each case should be considered carefully.

Precautions for use and warnings

The application in host bed with poor blood circulation is subject to closer examination and consideration by the attending physician due to poor healing rates.

Although no case of immunological rejection was observed in the long-standing use of human bone allografts, the possibility of such a rejection cannot be ruled out completely. After the initial testing of each donor, the batch processing takes place in accordance with all technical safety specifications, whereby various donors can be combined.

Interaction with other products

In concomitant medication with the substance group diclofenac, an impairment of bone growth is possible. No interactions with other medicinal products or medical devices are currently known.

Dosage forms and contents

maxgraft® products have double packaging. The primary packaging (granules: vial closed with a lid or blister with sealing film, others: sealed peel-off packaging) keeps the product sterile and is suitable for use in the surgery. The secondary packaging is a sealed peel-off package.

Form	Composition	Dimensions (mm)	Content
Granules	cortico-cancellous	< 2	0.5ml – 1.0ml – 2.0ml – 4.0ml
	cancellous	< 2	0.5ml – 1.0ml – 2.0ml – 4.0ml
Block	uncortical	20 x 10 x 10	One item each
		10 x 10 x 10	

74337
18.10.2021

For Russia

	cancellous	20 x 10 x 5	
		20 x 10 x 10	
		10 x 10 x 10	
		20 x 10 x 5	
		15 x 15 x 5	
Strut	cortical	30 x 15 x 1.5	One item each or 3x1 piece
		25 x 10 x 1.5	
		25 x 10 x 1	
Bonebuilder	cancellous	34 x 13 x 13 mm (max.)	One item
		23 x 13 x 13 mm (standard)	One item

Form	Composition	Outside diameter, height (mm)	Content
Bonering 3.3	cancellous	6, 10	One item each
		7, 10	
Bonering 4.1	cancellous	7, 10	

Method and time of application

maxgraft® is only intended for a single use. Aseptic conditions are essential when opening the maxgraft® packaging and for subsequent handling. If the implant is contaminated during the operation, it must be discarded.

Before application in patients, bigger material forms in particular should be rehydrated in a suitable physiological medium for at least 10 minutes (e.g. physiological saline). maxgraft® bonebuilder should be implanted soon after rehydrating it. It must be juxtaposed in the area of use and fixed with osteosynthesis screws.

The area to be filled must always be completely filled with bone substitute material. The final wound closure must be complete and saliva proof. Wound closure should be free of tension. The allogenic bone tissue must be stabilized in order to avoid any micromovement. A summary of clinical studies with human allogenic bone tissue showed its good incorporation, a lack of rejection and good remodelling.

The material is dense to X-rays, and therefore radiological control presents no problem.

Ambient temperature 25±10°C

Relative humidity not controlled

Atmospheric pressure not controlled

Implantable product: internal environment of the body

Temperature from 32 °C to 42 °C

A relative humidity of 100 %

Deleting and replacing

Not applicable as maxgraft® will be replaced by bodies own bone tissue over time and does not need to be removed.

Shelf life

maxgraft® has a shelf life of 5 years from the radiation date in the closed packaging. The expiry date is stated on the folding box as well as on the secondary and primary packaging. maxgraft® must not be used after the expiration date printed on the packaging.

74337
18.10.2021

For Russia

Limited or nullified use, even before the expiry date

If the primary or the secondary packaging is damaged, the sterility of the product is no longer guaranteed. In this case, maxgraft® may no longer be used. The contents of unused but opened or damaged packages may not be re-sterilised and are to be discarded.

Storage and transport

maxgraft® needs to be stored between +5°C and +30°C.

Store in the original packaging to protect the contents from humidity. The transport of maxgraft® is to be carried out protected from humidity at ambient temperature.

Sterilisation procedures

maxgraft® is supplied sterile in a closed vial, a sealed blister or a sealed peel-off packaging (primary packaging) and in a sealed peel-off packaging (secondary packaging). The content of both the primary packaging and of the secondary packaging is sterile and contained in a non-sterile outer cardboard box. Gamma rays are used for sterilisation and have no negative effect on stability, healing success or incorporation of the transplant.

Safe disposal

Remaining residues are to be discarded. Disposals are to be carried out according to local regulations.

Products with expired shelf life are classified as medical waste of class according to SanPiN 2.1.3684-21 and must be disposed of accordingly. After contact with the patient's body, the components of the product should be considered medical waste of class B according to SanPiN 2.1.3684-21 and disposed of accordingly.

Date of revision of the text

04/2021

Symbols used

Manufacturer

Use-by date

Batch number

Sterilized by irradiation

Do not re-sterilize

Do not use if package is damaged

Keep dry

Upper and lower limit of temperature

Do not re-use

Caution, consult the instructions for use

Registered trademark

Explanation of manufacturer symbols:

<https://www.ctba.at/symbole>

74337
18.10.20

For Russia

Any serious adverse reactions observed during or after clinical application, which may be linked to maxgraft[®], should be immediately reported to:

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Email: vigilanz@ctba.at

Warranty

botiss biomaterials GmbH guarantees the Buyer that "maxgraft[®]" (the product) purchased by the Buyer does not have material and quality defects at the time of delivery to the Buyer. The company "botiss biomaterials GmbH" does not give any guarantee (direct, indirect or statutory) to products that underwent a non-standard mechanical loading, storage in improper conditions, inappropriate usage/misuse, incorrect or negligent application, improper testing, usage in conjunction with other products or components other than those for which "maxgraft[®]" (the product) has been developed, or underwent a usage in a way or in medical procedures which are not intended for this product. Guarantee claims must be made without delay. Any guarantee claim is limited to the replacement of the defective product(s) and payment of shipping charges to the Buyer of the replaced product(s). In particular, "botiss biomaterials GmbH" does not compensate any additional components, tools, treatment costs and other costs and expenses happened due to or in connection with the replacement of the defective product to the Buyer or the patient.

Manufacturer's authorized representative:

Botiss Plus Ltd.,

Russia, St Peterburg, 194291, Lunacharskogo pr. 47, Lit. A

Tel. +7812558-15-15

maxgraft[®] Bone implant for dental surgery complies with the following standards as far as applicable

MEDDEV 2.7.1	Clinical evaluation: A guide for manufacturers and notified bodies
European directive 2006/17/EC	Commission Directive implementing Directive 2004/23/EC of the European Parliament and of the Council as regards certain technical requirements for the donation, procurement and testing of human tissues and cells
European directive 2006/86/EC	Commission Directive implementing Directive 2004/23/EC of the European Parliament and of the Council as regards traceability requirements, notification of serious adverse reactions and events and certain technical requirements for the coding, processing, preservation, storage and distribution of human tissues and cells
EN ISO 13485	Medical devices — Quality management systems — Requirements for regulatory purposes
EN ISO 14971	Medical devices. Application of risk management to medical devices
DIN EN 1041	Information supplied by the manufacturer of medical devices
EN ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 14476	Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of virucidal activity in the medical area - Test method and requirements
EN ISO 11137-2	Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose
ASTM F1929	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration

74337
18.10.2021

For Russia

ASTM F2096	Standard Test Method for Detecting Gross Leaks in Packaging by Internal Pressurization (Bubble Test)
ISO 10993	Medical devices. Biological evaluation of medical devices.

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