
True Copy Statement

As a responsible official of New World Medical, Inc., I hereby certify that the «Instructions Application at Medical Device» is a true and correct to best of my knowledge. The Instructions Application at Medical Device will be submitted to the Ministry of Health of the Russian Federation for registration:

Ahmed Glaucoma Valve, models FP7, FP8 composed of:

- Valve Glaucomatous Ahmed - 1 pc.
- Instructions for use 1 - pcs.
- Label for the patient's history - 1 pc.
- Implant registration card - 1 pc.
- Patient's personal card 1 - pcs.

This copy contains no alterations from the original. A photocopy of the document listed above is acceptable with attachment of this original statement.

Victor Arellano
Sr. Global Regulatory Affairs Manager
New World Medical, Inc.

08/02/2022
Date

10763 Edison Court • Rancho Cucamonga, California 91730
Telephone (909) 466-4304 • Fax (909) 466-4305 • Website www.newworldmedical.com

New World Medical, Inc.


INSTRUCTIONS APPLICATION AT MEDICAL DEVICE

The Ahmed Glaucoma Valve, models FP7, FP8 composed of:

- Valve Glaucomatous Ahmed - 1 pc.
- Instructions for use 1 - pcs.
- Label for the patient's history - 1 pc.
- Implant registration card - 1 pc.
- Patient's personal card 1 - pcs.

“APPROVE”

«New World Medical, Inc.», USA

 Victor Arellano

(signature)


07.09.2015

«day» month (in numbers)

L.S.

Date of last revision:

05.07.2022

NAME OF MEDICAL DEVICE

The Ahmed Glaucoma Valve, models FP7, FP8 composed of:

- Valve Glaucomatous Ahmed - 1 pc.
- Instructions for use 1 - pcs.
- Label for the patient's history - 1 pc.
- Implant registration card - 1 pc.
- Patient's personal card 1 - pcs.

(further in the text The Ahmed Valve).

PURPOSE OF THE MEDICAL DEVICE, INSTALLED BY THE MANUFACTURER

The Ahmed Valve is intended for treatment of refractory glaucoma when previous surgical treatment has been unsuccessful, or, based on the experience, it is known that it has not yielded satisfactory results with the use of antimetabolic agents.

Such refractory glaucoma may include neovascular glaucoma, initial and advanced open-angle glaucoma, refractory to drug therapy, advanced congenital and infantile glaucoma. Also refractory (repeatedly operated) glaucomas caused by aphakia and uveitis, angle-closure glaucoma, with glaucoma induced by keratoprosthetics, penetrating keratoplasty, trauma, with Sturges-Weber syndrome and other types of secondary glaucomas. Traumatic, absolutely painful, glaucoma in the single eye, congenital glaucoma, iridocorneal endothelial syndrome, aniridia, open-angle glaucoma, with angiomas, with combined hypotensive surgery with phacoemulsification, organic angle blockade with arphakia, IOP from 30 and higher, absolute glaucomas for purpose of organ saving.

The Ahmed valves are used as a primary intervention in cases of neovascular, congenital, uveal and open-angle glaucoma, in the presence of perilimbal scarring or anterior cyanosis.

It is used as a secondary intervention in the absence of intraocular pressure compensation with maximum drug therapy after laser intervention.

APPLICATION CONDITIONS:

Application conditions: $20^{\circ}\text{C} \pm 5^{\circ}\text{C}$.

Relative humidity: $\leq 75\%$.

Atmosphere pressure: 84,0-106,7 kPa (630-800 mm. Hg).

CONTRAINDICATIONS

Bacterial conjunctivitis

Bacterial corneal ulcers.

Endophthalmitis.

Orbital cellulitis.

Bacteremia or septicemia.

Active scleritis.

COMPLICATIONS AND ADVERSE REACTIONS

Edema of the cornea of the eye.

Contact with the cornea.

Contact of tubing with the iris.

Fusion.

Open transplant of the sclera.

Choroidal detachment.

Inflammation of the iris.

Hyphema.

Tubing occlusion.

Retraction of the tubing.

Upon proper implantation, it is possible to avoid complications: hypotension, small eye chamber, choroidal effusion, suprachoroidal hemorrhage, retinal detachment, cataracts, loops in the conjunctival cavity, phthisis of the eyeball, bullous keratopathy, uveitis and end ophthalmitis.

DESCRIPTION OF OPERATING PRINCIPLES AND THEIR FEATURES

The basic principle of the Ahmed valve is the ability to remove aqueous humor from the anterior chamber without the development of hypotension. Another criterion is the duration of the hypotensive effect with minimal side effects on the surrounding tissues.

The Ahmed valve represents itself a silicone membrane folded together and elongated at the top. The membrane is based on a polypropylene plate (the body), (Fig. 1), which having additional protrusions and grooves, which pull the membrane (Fig. 2). When the IOP becomes higher than 8-10 mm. Rt., the silicone membrane shall be opened for the outflow of liquid.


Figure 1. Tubing inserted into a membrane that is clamped by upper and lower plateaus of the body.


Figure 2. Cross section of the Ahmed Valve body with membrane.
The stretched type membrane.

The resistance of the Ahmed valve is self-regulated by level of intraocular pressure. Upon high initial pressure in the anterior chamber is high, the Ahmed valve automatically opens in full. As soon as the pressure decreases, the membrane stretches, closing the lumen, reducing the outflow of fluid (Fig. 3). Since the Ahmed valve begins to work only at a pressure higher than 8-10 mm Hg., it prevents the development of hypotension.


Figure 3. Comparison of fluid flow.

A special feature of the Ahmed Valve is its implantation in one step, which instantly provides a decrease in intraocular pressure after surgery. Among all tubular implants, the Ahmed Valve has the lowest postoperative level of hypotension and fewer cases of complications.

The peculiarity of the internal resistance to the outflow of fluid in the implant system ensures the principle of operation of the Venturi mechanism of the Ahmed Valve (at the valve, the inlet area is larger than the outlet area) (Fig. 4).


Figure 4. The Venturi Mechanism of the Ahmed Valve.

The Ahmed valve is a product of non-obstructive type. Since the transverse size of the Ahmed Valves is larger than the circumference of the narrow tubing, each particle passing through this tubing also passes through the valve mechanism.

TECHNICAL SPECIFICATIONS

Table 1. Characteristics of the Glaucomatous Ahmed Valve, Models FP7, FP8

Product	Component parts	Specification	
Model FP-7	Tubing	I.D. O.D. Length	0,305 mm $\pm$ 0,076 mm 0,635 mm $\pm$ 0,076 mm 25,4 mm $\pm$ 1,524 mm
	Body	Thickness Length Width	1,524 mm $\pm$ 0,254 mm 4,851 mm $\pm$ 0,127 mm 8,636 mm $\pm$ 0,254 mm
	Plate	Thickness Length Width	2,108 mm $\pm$ 0,076 mm 16,002 mm $\pm$ 0,127 mm 13,005 mm $\pm$ 0,127 mm
	Valve Parameters Glaucomatous Ahmed	Weight Outflow speed Pressure control Plate's Surface Area	0,11 g $\pm$ 0,006 g 12 μ l/min $\pm$ 3 μ l/min 8 -12 mm. Hg. 186 mm ² $\pm$ 9,3 mm ²
Model FP-8	Tubing	I.D. O.D. Length	0,305 mm $\pm$ 0,076 mm 0,635 mm $\pm$ 0,076 mm 25,4 mm $\pm$ 1,524 mm
	Body	Thickness Length Width	1,524 mm $\pm$ 0,254 mm 4,851 mm $\pm$ 0,127 mm 8,636 mm $\pm$ 0,254 mm
	Plate	Thickness Length Width	2,108 mm $\pm$ 0,076 mm 10,846 mm $\pm$ 0,127 mm 10,922 mm $\pm$ 0,254 mm
	Valve Parameters Glaucomatous Ahmed	Weight Outflow speed Pressure control Plate's Surface Area	0,09 g $\pm$ 0,006 g 12 μ l/min $\pm$ 3 μ l/min 8 -12 mm. Hg. 96 mm ² $\pm$ 4,8 mm ²

OPERATING TECHNIQUE

OPERATION PROGRESS

To prepare the Ahmed Valve for operation, it should be flushed with a balanced saline solution to make sure it is completely bled and the Valve is functioning (fig. 5). For this, a syringe with a 27G cannula is inserted into the valve tubing, which is held in place with the tweezers. Press sharply on the syringe plunger with a force sufficient to expel the liquid from the valve, till the moment when 2-3 ml of solution shall not be visible behind the valve device or at the top of its plate. Then, valve is then ready for implantation.

Figure 5. Flushing the Ahmed Valve with Balanced Salt Solution.

The preferred site for implantation is the upper outer quadrant. Do not implant an Ahmed valve in the upper inner quadrant, because this is associated with a high risk of Pseudo Brown syndrome and damage to the tendon of the superior oblique muscle.

Make an incision in the conjunctiva at the limbus or in the conjunctival fornix.

Form a conjunctival pocket between the conjunctiva and tenon capsule between the two adjacent rectus muscles (Fig. 6). For this purpose, it is suitable to use Stevens scissors for this purpose.


Figure 6. Formation of the conjunctival flap with the plate to the fornix.

The Ahmed valve should not be grasped by the body of the valve with tweezers (Fig. 7), as damage to the body will result to damage of membrane that controls the outflow of fluid.

It is need, to grasp the Ahmed valve should be grasped by the plate near the small opening and then to insert into the formed subconjunctival space.


Figure 7. Insertion of the Ahmed valve into the subconjunctival space.
Do not grasp by the plate as shown here.

Before suturing the valve's plate to the sclera, it is needed to make sure that the valve is positioned between the rectus muscles and the tubing is directed in radial to the limb. The front edge of the Ahmed Valve should be situated at 8-10 mm from the limb (Fig. 8).


Figure 8. The Ahmed valve is placed as possible from the edge of the cornea as possible (8-10 mm from the limb).

Next, sew the Ahmed valve should be sutured to the sclera (Fig. 9) with non-absorbable sutures (for example, nylon suture 8/0, 9/0).


Figure 9. Suturing of the Valve with non-absorbable sutures.

Cutting off the tubing.

Place the tubing over the cornea. Calculate the length of the drain tubing so that the length of it should be sufficient to be inserted into the anterior chamber.

The length of the tubing in the anterior chamber should not be more than 2-3 mm.

Then the drainage tubing should be cut with sharp scissors at an angle of 30° so that when inserted into the anterior chamber, the beveled end shall be directed towards the corneal endothelium.

Creation of a hole for the tubing.

Option 1. Creation of a conventional scleral flap as in trabeculectomy.

Using a diamond knife, a scleral flap of 4 mm by 6 mm in size, 2/3 thick, is cut with the plate to the limb in order to cover the tubing after it has been inserted into the anterior chamber (Fig. 10).


Figure 10. Scleral flap formation and piercing with a 23G needle.

Lift the scleral flap with forceps to continue the operation.

Use a sharp 23G needle to puncture the anterior chamber in the corneoscleral area (Fig. 11). The needle travel should be in parallel to the iris.


Figure 11. Puncture with a 23G needle using your own sclera

Viscoelastic can be applied to the tubing cut for easier insertion of the tubing into the anterior chamber.

Using forceps, the tubing is passed into the anterior chamber through a puncture made with a 23G needle (Figure 12). The tubing should not touch the cornea, iris or lens. The length of the tubing in the anterior chamber should be 2-3 mm.


Figure 12. Inserting the tubing into the anterior chamber using forceps through the puncture.

The tubing may appear a bit longer. When the anterior chamber recovers, the tubing may move slightly backward. Make sure the beveled edge is facing up. It is important. In such position, the tubing is not blocked by the iris when excess fluid outflows from the anterior chamber.

Option 2. Using a donor sclera or perekarda.

Using a sharp 23G needle, start puncturing at distance of 2-3 mm from the limb. The needle travel

should be in parallel to the iris (Fig. 13). Viscoelastic can be applied to the tubing cut for easier insertion of the tubing into the anterior chamber. Firmly, push the tubing into the anterior chamber. Make sure so that the beveled edge is facing up.


Figure 13. Cutting and tunneling of tubing with a 23G needle while covering the tubing with donor material.

The donor flap must be not less than 4x6 mm in order to close the tubing.

The donor sclera or pericardium is sutured to the sclera with 9/0 polyglactin or nylon suture on a flat needle (Fig. 14).

Make sure so that the tubing is well closed with your own tissue or donor material. It permits to avoid infectious complications.


Figure 14. Using the donor sclera to cover the tubing.

To reconstruct the anterior chamber, inject saline solution or hyaluronic viscoelastic through paracentesis outside the area of surgical intervention. At this time, the tubing should be situated at a distance of 2-3 mm in the anterior chamber.

After the end of the surgery, administer antibiotics and corticosteroids at 180° from the site of intervention subconjunctivally.

If tubing blockage occurs at any time after surgery, use a long 27G needle, inserting it through the paracentesis from the opposite side of the anterior chamber into the valve tubing (Figure 15).


Figure 15. Flushing the valve tubing through paracentesis.

After inserting the needle, push with force 0.5 ml of saline solution through the valve tubing.

MODE OF APPLICATION

Preparation of Valve Ahmed. Before application / use The Ahmed valve must be carefully inspected and prepared. Preparation is carried out by injecting a balanced saline solution through the drain tubing and valve. Insert a 27G Balanced Saline Syringe with a 27G cannula into the open end of the drain tubing. Slowly, inject about 1.0 cm³ of the solution into the implant. The Ahmed valve is now ready for use. *General surgical intervention.* The plate of Ahmed Valve is placed above the equator of the eyeball between the rectus muscles, in a previously prepared pocket formed by the dissection of the conjunctiva and tenon's capsule and by separating the tenon's capsule from the sclera. The plate of Ahmed Valve is sutured to the sclera by fenestration of its anterior edge, approximately 8-10 mm from the limbus. The drainage tubing is cut so that it is inserted into the anterior chamber of the eye by 2-3 mm. The edge of the tubing should be cut at a 30-degree angle towards the cornea. A two-thirds thick scleral flap is created and the plate cut to the limb. A puncture is made under the flap and into the anterior chamber with a sharp 23G gauge needle. The drainage tubing is inserted through the channel formed by the needle and into the anterior chamber in parallel to the iris plane. The scleral flap is sutured.

Warnings, cautions.

Before using the Ahmed valve, the implanting surgeon must be proficient in the technique of filtration surgery for glaucoma using drainage devices and postoperative care.

MANUFACTURING MATERIALS

Tubing, membrane, plate - silicone elastomer.

The body is made of polypropylene.

The adhesive is organosilicon compounds (silicone).

STERILIZATION METHODS

The sterilization method is radiation.

STORAGE CONDITIONS

- Store the medical device in its original packaging at temperatures from +10°C up to +35°C.
- Air relative humidity: no more than 75%.
- Atmosphere pressure: 84,0-106,7 kPa (630-800 mm Hg).



TRANSPORTATION CONDITIONS

Transportation of a medical device of Valve Ahmed must comply with the standards for transportation of medical devices and must be carried out at a temperature from -20 ° C to + 35 ° C, humidity from 60 to 80%.

The medical device is transported by all types of transport in covered vehicles in accordance with the rules on carriage of cargos, in force for each type of transport.

MARKING

Symbols used:

	Manufacturer		Manufacturing date		Не стерилизовать повторно
	Batch code		Use-by date Marking		Do not re-use
	Serial number		Catalog number		Sterilized using irradiation
	CE labeling approved		Do not use if the package is damaged		Authorized Representative in the European Community
	Refer to instructions for use		Safe for MRI		Caution: Federal (USA) law restricts this device to sale of this product or on the order of a physician
	Temperature Range		Range humidity		Fragile, handle with care
	Keep away from moisture		Do not expose to direct sunlight		Carefully! Refer to instructions for use

Manufacturer codes

(01)... .. barcode number
 (11)... .. manufacturing date
 (17)... .. valid until
 (21)... .. serial number
 (10)... .. batch code

Sample of Russian-language labeling (stickers)

The Ahmed Glaucoma Valve, model FP7 composed of:

- Valve Glaucomatous Ahmed - 1 pc.
- Instructions for use 1 - pcs.
- Label for the patient's history - 1 pc.
- Implant registration card - 1 pc.
- Patient's personal card 1 - pcs.

New World Medical, Inc. Made in the USA,
 country of origin USA.

Good (sterile) until 04/22/2023

RU No. _____ dated _____ year.

DS/CC No. _____ dated _____ year.

Storage conditions:

temperature: from + 10 ° C to + 35 ° C, humidity: no more than 75%

Authorized representative in territory of the RF:

LLC "Vision Technology"

room number 6, premise I. 15, Astradamskaya str., Moscow, 127206

Postal address: office 4, 46, 3rd Tverskaya-Yamskaya str., Moscow, 125047,

rustem@ahmedglaucoma.ru phone – 7 (495) 600-44-14

SHELF LIFE

Shelf life 2 years. The expiration date is indicated on the package.
Do not use after the expiration date and if a package defect is found.

PACKAGING

Blister - three-layer package.

-**Blister** packaging.

The valve is placed in an inner and outer blister pack of 1 pc.

-**Consumer** packaging.

The blister package is packed into a consumer package (box) - 1 unitary unit.

Consumer packaging also includes:

- **Instructions** for use - 1 pc.
- **label** for the patient's history - 1 pc.
- **Implant** registration card - 1 pc.
- **patient's** personal card - 1 pc.

-**Transportation** packaging.

Consumer packages are placed in shipping cartons for transportation of the product in quantities according to the order (from 1 piece or more). The instructions for use are enclosed in an amount corresponding to the number of consumer packages. Empty spaces are filled with wrapping paper to prevent movement and damage to consumer packages.

Completeness

- **Valve Glaucomatous Ahmed** - 1 pc.
- **Instructions** for use - 1 pc.
- **Patient** history label - 1 pc.
- **Implant** registration card - 1 pc.
- **Patient's** personal card - 1 pc.

REQUIREMENTS FOR ENVIRONMENTAL PROTECTION WHEN USING A MEDICAL DEVICE

This medical device during use, transportation and storage does not have a negative impact on humans and the environment.

REQUIREMENTS FOR ENVIRONMENTAL PROTECTION WHEN USING A MEDICAL DEVICE

This medical device during use, transportation and storage does not have a negative impact on humans and environment.

PROCEDURE FOR DISPOSAL OR DESTROYING OF A MEDICAL DEVICE

Dispose of the product after use in accordance with applicable national practice.

Disposal of the product must be carried out in accordance with the documents regulating the handling of medical waste.

An unused medical device with an expired expiration date, with damaged packaging, as well as packaging, are disposed of in accordance with SanPin 2.1.3684-21, class A medical waste hazard.

After use, the medical device is disposed of according to the hazard class of medical waste - B in accordance with SanPin 2.1.3684-21.

MANUFACTURER'S WARRANTIES

New World Medical, Inc. guarantees the quality of the medical device Glaucomatous Ahmed Valve, models FP7, FP8 throughout the shelf life from the date of release, subject to the proper conditions of release, transportation, storage and use, read the instructions before use.

Manufacturer:

New World Medical, Inc., 10763 Edison Court, Rancho Cucamonga, CA, 91730, USA
+1(909)4664304, +1(909)4664305
info@newworldmedical.com

Quality claims are accepted by the Manufacturer's Authorized Representative / Regional Distributor in the Russian Federation:

Limited Liability Company "Vision Technology"
(Vision Technology LLC),

Legal address:

room number 6, premise I, 15, Astradamskaya str., Moscow, 127206,

Actual:

office 4, 46, 3rd Tverskaya-Yamskaya str., Moscow, 125047,
8 (495) 600-44-14, 8 (909) 910-10-70

Patient's personal card layout (front side)

AHMED™ GLAUCOMA VALVE	
PATIENT IDENTIFICATION CARD	
NAME: _____	
DATE OF IMPLANTATION: _____	
EYE: ☐	RIGHT ☐ LEFT
IN CASE OF EMERGENCY CONTACT DOCTOR ON REVERSE SIDE	

Patient's personal card mock-up (back side)

DOCTOR: _____	
TELEPHONE: _____	
Part #50-0013 REV A 06/15	
MR MR NAME	SN G691471
REF FP7	USE-BY 2023-04-22
LOT C1921	

Implant registration card mockup

NO POSTAGE
NECESSARY
IF MAILED
IN THE
UNITED STATES

BUSINESS REPLY MAIL
FIRST-CLASS MAIL PERMIT NO. 9188 RANCHO CUCAMONGA, CA

New World Medical Inc
 19763 EDISON COURT
 RANCHO CUCAMONGA, CA 91730-8754

REMOVE TAPE, SEAL & MAIL - DO NOT STAPLE

REMOVE TAPE, SEAL & MAIL - DO NOT STAPLE

IMPLANT SURGEON:
 Please complete and mail this form immediately upon use of implant.

DOCTOR: _____

HOSPITAL/CLINIC NAME: _____

CITY, STATE: _____

DR. TELEPHONE #: _____

DATE OF SURGERY: _____

DISPOSITION ☐ IMPLANTED
 OF IMPLANT ☐ RETURNED
☐ DESTROYED

SN	G691471	 USE-BY  2023-04-22
REF	FP7	
LOT	C1921	

Part 2 50-0014 Rev C 07/16

Patient history label mockup

SN	G691471	USE-BY	2023-04-22	SN	G691471	USE-BY	2023-04-22
REF	FP7			REF	FP7		
LOT	C1921			LOT	C1921		
SN	G691471	USE-BY	2023-04-22	SN	G691471	USE-BY	2023-04-22
REF	FP7			REF	FP7		
LOT	C1921			LOT	C1921		
SN	G691471	USE-BY	2023-04-22	SN	G691471	USE-BY	2023-04-22
REF	FP7			REF	FP7		
LOT	C1921			LOT	C1921		

CALIFORNIA ALL-PURPOSE ACKNOWLEDGMENT

CIVIL CODE § 1189

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California)

County of San Bernardino)On August 2, 2022 before me, Arlene Renee Maldonado, Notary Public

Date

Here Insert Name and Title of the Officer

Personally appeared Victor Arellano

Name(s) of Signer(s)


who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature

Signature of Notary Public

Place Notary Seal Above

OPTIONAL

Though this section is optional, completing this information can deter alteration of the document or fraudulent reattachment of this form to an unintended document.

Description of Attached Document

Title or Type of Document Instructions Application at Medical Device Document Date August 2, 2022Number of Pages n/a Signer(s) Other Than Named Above n/a

Capacity(ies) Claimed by Signer(s)

Signer's Name n/a

- ☐ Corporate Officer—Title(s) n/a
☐ Partner ☐ Limited ☐ General
☐ Individual ☐ Attorney in Fact
☐ Trustee ☐ Guardian or Conservator
☐ Other n/a

Signer Is Representing n/aSigner's Name n/a

- ☐ Corporate Officer—Title(s) n/a
☐ Partner ☐ Limited ☐ General
☐ Individual ☐ Attorney in Fact
☐ Trustee ☐ Guardian or Conservator
☐ Other n/a

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