ALVO®

APPROVE

Agata Olszewska

Proxy

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
25.2 International standards.....	26
Appendix 1. Record card for repairs and inspections.....	26
Appendix A. Declaration of Electromagnetic Compatibility.....	27

Index	4
1. Name of the medical device	4
2. Manufacturer of the medical device	4
3. Purpose of the medical device	4
Potential consumers:.....	4
Users of the medical device:.....	4
4. Indications and contraindications for the use of a medical device	4
4.1 Indications:.....	4
4.2 Contraindications.....	4
4.3 Side effects when following the manual.....	4
5. Description of the main functional elements and accessories of the medical device	5
5.1 Functional characteristics of the ALVO Conveyor:.....	5
5.2 ALVO Conveyor control panels.....	9
5.3 Built-in Euro Power EV2212 rechargeable battery.....	10
5.4 Power supply cable.....	11
5.5 Potential equalization cable.....	11
6. Sterility	12
7. Safety requirements during maintenance	12
7.1. Safety requirements.....	12
7.2 Safety measures.....	13
8. Sequence of unpacking and first turn on of the ALVO Conveyor	13
8.1. Unpacking.....	13
8.2. First turn on.....	14
8.3 ALVO Conveyor move order.....	14
9. Charging batteries	14
10. Antistatic properties of the ALVO Conveyor	16
11. Assessment of the serviceability of the ALVO Conveyor	16
12. Crashes and defects	16
13. Care and disinfection	17
14. Maintenance, inspection, repair	17
15. Monitoring the technical condition and inspection the technical condition	17
16. Elimination of possible malfunctions	19
17. Marking of the ALVO Conveyor	20
17.1. Factory identification plate of the ALVO Conveyor and its decoding.....	20
17.2. ALVO Conveyor identification plate for Russian Federation.....	20
17.3 Marking of the transport packaging of the ALVO Conveyor.....	21
17.4 Symbols used.....	21
18. Packing of the ALVO Conveyor	22
19. Conditions of transportation and storage	22
19.1 Transport conditions.....	22
19.2 Storage conditions.....	22
20. Maintenance conditions	22
21. Expiry date	23
22. Warranty and rules for its provision	23
22.1 Warranty conditions.....	23
22.2 Exclusions from warranty.....	23
22.3 Post-warranty repair.....	23
23. Disposal	24
24. Authorized representative of the manufacturer	24
25. List of national and international regulatory documents / standards that the medical device complies with	24
25.1 National standards (compliance with RF standards).....	24

1. Name of the medical device

Patient transfer system: Model ALVO Transpa 4-12 (hereinafter referred to as the ALVO

Conveyor, product).

Accessories:

1. Built-in rechargeable battery EURO POWER EV22-12-2 pcs.
2. Power cable 3 m -1 pc.
3. 5 m potential equalization cable-1 pc.
4. User manual 1 piece.

2. Manufacturer of the medical device

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Ul. Poludneva 21A,
64-030 Smigiel, Poland

3. Purpose of the medical device

The ALVO Conveyor is designed for transporting and safely moving (shifting) low-mobility and immobilized patients on operating tables, hospital beds and other trolleys.

Potential consumers:

The ALVO Conveyor is used in the emergency rooms of clinical hospitals, operating rooms, departments of surgery and intensive care, in the departments of traumatology, radiology, cardiology, maternity hospitals and other departments of hospitals.

Users of the medical device:

Qualified medical personnel.

4. Indications and contraindications for the use of a medical device

4.1 Indications

The ALVO Conveyor is designed for transporting and transferring low-mobility and immobilized patients to operating tables, hospital beds and other trolleys. There are no age restrictions for moving.

4.2 Contraindications

There are no contraindications to use.

4.3 Side effects when following the manual

No side effects.

5. Description of the main functional elements and accessories of the medical device.

The ALVO Conveyor is equipped with an automatic conveyor for moving a sedentary or immobilized patient with a weight of up to 250 kg.

Large wheels provide high reliability when transporting inside a medical facility.

The process of shifting a patient is as simplified as possible, it takes place with the help of an electromechanical belt conveyor.

The movement (shifting) of the patient with the help of the device is provided by an electromechanical method. The ALVO Conveyor is based on the principle of lateral movement of the patient between two horizontal surfaces without the use of lifts. Control of the process of moving, lowering and lifting the table is carried out using the control buttons located on the main handrail of the product.

The elements of the ALVO Conveyor are made of steel, plastic and aluminum alloys.

5.1 Functional characteristics of the ALVO Conveyor:

1. No risk of injuring the patient;
2. The process of shifting a patient is as simplified as possible, it takes place with the help of an electromechanical belt transporter.
3. The movable belt of the electromechanical belt transporter is resistant to disinfection and washing;
4. There is an electromechanical adjustment of the height of the table (the upper part of the conveyor);
5. High load capacity up to 250 kg, which allows you to shift even very large patients without the physical effort of the staff;
6. Two specialists from the medical staff are enough to transport and transfer the patient: caregiver for supporting the patient during the transport and second one for operating the device;
7. The system can be operated for a long time by a rechargeable battery, which is important for transporting the patient between departments, and wards;
8. The conveyor is powered by electricity from the built-in battery, in offline mode;
9. There is an indication of the battery charge level;
10. Antistatic wheels provide high traffic inside the medical facility;
11. Wheels are protected from winding of hoses and cables;
12. The conveyor is equipped with a drive wheel.
13. The wheels of the ALVO Conveyor are equipped with a mechanical brake;
14. Easy to manage processes: the function of the re-laying system is carried out using the control buttons located on the main handrail of the product.


Illustration 1. ALVO Conveyor

Table 1. Description of ALVO Conveyor parts

Nº	Part Name	Description
1	Footing	Metal structure equipped with wheels and a mechanical wheel lock system
2	The upper part	Electromechanical device for shifting patients
3	Protective limiter	Metal structure adjusted with gas spring that prevents patient fall during shifting
4	Drive rack I	Special device equipped with the drive that moves in the vertical plane the upper part (2)/ Electromechanical height adjustment of the upper part (2)
5	Upper corrugated casing	Protective element. Corrugated plastic construction. Closes the rack drives (4 and 9)
6	Frame	The main element of the trolley, containing the drive of the device for shifting patients in a horizontal plane and electronic components in the casing
7	Motorized wheel	Drive wheel with DC motor preassembled which is an electromechanical system. The movement of the conveyor is provided by an electromechanical system enclosed in a casing
8	Main grab handle	Metal structure with built-in control panels I and II. Designed for maneuvering the trolley when shifting patients to the bed/operating table. Equipped with a motorized wheel control panel (7) and a device for shifting patients
9	Drive rack II	A special device equipped with the drive that allows to move the upper part (2) in the vertical plane. Electromechanical height adjustment of the upper part (2)
10	Side grab handle	Metal construction to move the ALVO Conveyor more efficiently. Designed for maneuvering the trolley when moving patients between different points in the room
11	Bottom corrugated casing	Protective element. Corrugated plastic construction. Closes the rack drives
12	Rolling belt	Canvas with a width of ... in a looped-back form. Designed for horizontal movement (shifting) of the patient
13	Mechanical wheel lock	Pedal with a special device that has a braking mechanical effect on the wheel. Designed to block the movement of the wheel


Illustration 2. Photo of the ALVO Conveyor

Table 2. Technical, weight and size and electrical characteristics

Total length	2150 ± 50 mm
Overall height	1320 ± 50 mm
Width of the upper part	650 ± 50 mm
The length of belt table	1835 ± 50 mm
Table length	1950 ± 50 mm
Overall width	1100 ± 50 mm
Height adjustment range	750-1000 mm ± 20 mm
Wheel diameter	150 mm
Range of movement (movement) of the rolling belt	650 mm ± 50 mm
ALVO Conveyor weight	250 kg
Maximum permissible load	250 kg
Maximum speed when moving with a motorized wheel	4 km / h
Number of batteries	2 pcs.
Supply voltage	24 V
Battery parameters (maintenance-free)	24 V, 15 A / h
Battery charging time	< 12 h
Operating cycle	2/18 min.
Battery life without chargings	About 60 operations of moving the patients.
Charger Type	Built-in
Charger power supply voltage	230 V, ~50 Hz
Power consumption	350 VA
Power cable length	At least 3 m.
Service life (Warranty period)	10 years (2 years)
The service life of the battery	5 years
The degree of water protection IEC 60529:2013 (GOST 14254-2015)	IP X4
Protection class against electric shock according to GOST R MEK 60601-1-2010	I
Working part according to GOST R MEK 60601-1-2010	Type B
The force of movement	With the brakes released, no more than 50 N. With the brakes on, no more than 295 N.


Illustration 3. Overall dimensions of the ALVO Conveyor


Illustration 4. Placement of the motorized wheel on the Conveyor.

Illustration 5. The motorized wheel photo.

Table 3. Technical and electrical characteristics of the motorized wheel

№	Characteristics	Value
1	Wheel size	125 mm
2	Tread width	35 mm
3	Maximum lifting force / downforce	400H
4	Rated voltage	24 V
5	Rated moment	4 Hm
6	Rated power	73.1 V
7	Rated speed of rotation	174.50 min ⁻¹
8	Rated current	7.0 A
9	Max speed	4 km/H
10	Protection class	IP X4
11	Direction of movement	forward / backward
11	Weight	6 kg

5.2 ALVO Conveyor control panels.

The control panels I and II are located on the main handrail designed to move the ALVO Conveyor. Panels I (Illus. 7) and panels II (Illus. 8,9) contain six function buttons that allow you to adjust the height of the upper part of the trolley, control the device for shifting patients and the motorized wheel.

To select the desired function, press and hold the button marked with the corresponding symbol. The function is not disabled while the button is still pressed.


Illustration 6. Placing the control panel on the main handrail


Illustration 7. Designation of the control buttons on panel I.


Illustration 8. Designation of the control buttons on panel II

Illustration 9. Photo of the control buttons on panel II

5.3 Built-in EURO POWER EV22-12 rechargeable battery

The ALVO Conveyor is equipped with two built-in rechargeable batteries and a EURO POWER EV22-12.

Illustration 10. EURO POWER EV22-12 battery pack.

Table 4. EURO POWER EV22-12 Battery Specifications

No	Name of the characteristic	Indicator
1	Output voltage (V)	12
2	Capacity (A / h)	15
3	Height (mm)	166
4	Width (mm)	76
5	Depth (mm)	181
6	Weight (kg)	6.5
7	Service life	5 years
9	Battery type	Lead-acid, maintenance-free
10	Charge current	$I_{max} = 0.3 C_{nom}$ (the recommended charge current is $0.1 C_{nom}$)
11	Charge voltage	Buffer mode: $U = 2.25-2.30 V / cell$ Cyclic mode: $U = 2.40-2.50 V / cell$
12	Discharge current:	$I_{max} = 15 C_{nom}$ (up to 40Ah); $I_{max} = 10 C_{nom}$ (over 40 A / h)
13	Operating temperature	Charge: $0^{\circ}C$ to $+40^{\circ}C$ Discharge: $-20^{\circ}C$ to $+60^{\circ}C$

5.4. Power supply cable

For enabling charging from AC 230 V; 50/60 Hz; 3.5 A.

The power cord:

- Length: 3.0 m;
- Weight: 500 g

Illustration 11. Power supply cable.

5.5. Potential equalization cable

If it is necessary to equalize the potentials, the ALVO Conveyor is connected via a cable to the elements of electrical installations that have zero potential, the location of the corresponding connector on the conveyor is shown in Illus. 12.

Specifications

- length 5 m;
- cross-section $5 mm^2$;
- yellow-green color.

Manufacturer: HAGER POLO company, Poland.

Illustration 12. Potential equalization cable.

6. Sterility

The ALVO Conveyor is a non-sterile product. Sterilization before use is not required.

7. Safety requirements during maintenance

7.1. Safety requirements

1. Carefully read the user manual, familiarize yourself with the working methods and warnings;
2. Before performing functional movements, block the free movement of the ALVO Conveyor by pushing the mechanical wheel lock in the direction marked with a locked padlock icon on a red background.
3. When adjusting the height of the upper part of the ALVO Conveyor, secure the ALVO Conveyor elements against collisions with other objects.
4. When using the ALVO Conveyor near high-frequency medical devices or defibrillators, the instructions in the manuals for those devices must be followed.
5. When moving the ALVO Conveyor, protect it from collision with other objects.
6. The ALVO Conveyor must be transported on a level floor free of any obstacles.
7. Do not move the ALVO Conveyor across electric cables.
8. The parameters of the power supply used to connect the ALVO Conveyor must correspond to the rating data given on the nameplate.
9. Do not use power cords that appear damaged.
10. Do not connect the ALVO Conveyor to a power supply in hazardous conditions, for example, in explosive environments.
11. Do not store the ALVO Conveyor in a state with discharged batteries.
12. Do not disinfect the ALVO Conveyor in the disinfection chamber and place it in a high pressure chamber.
13. Do not use bleaches containing active chlorine and oxygen to disinfect the ALVO Conveyor.

14. Do not use any substances that destroy the plastic structure to wash the plastic parts of the ALVO Conveyor.

15. Do not use alcohol-based disinfectants to clean and disinfect the ALVO Conveyor.

7.2 Safety measures

1. Use the ALVO Conveyor for its intended purpose.
2. Use only compatible accessories recommended by the manufacturer;
3. Before switching on, let the ALVO Conveyor stand for three hours at room temperature;
4. Do not use the ALVO Conveyor if the power cord is damaged, inoperative, or malfunctioning. Contact the service center;
5. The power cord should be away from heat sources;
6. To reduce the risk of electric shock, the conveyor should only be connected to a grounded outlet;
7. Service personnel are obliged to work in medical gloves;

WARNING:

1. Improper use of the grounding plug may result in electric shock;
2. Do not use any adapters with the charger of the cart, if the plug does not fit into the outlet, do not try to change it. Contact a qualified electrician to install the correct outlet;
3. Maintenance of the ALVO Conveyor must be carried out by ALVO representatives or the manufacturer's representatives;
4. Warning: danger of electric shock. Do not attempt to carry out maintenance by yourself. Opening the product may result in electric shock or other damage and will void the warranty;
5. To avoid damage to the ALVO Conveyor and additional equipment, the ALVO Conveyor must be protected from mechanical collisions with other objects.

Be careful when moving the cart! Remember that you are responsible for the correct handling of the equipment!

8. Sequence of unpacking and first turn on of the ALVO Conveyor

8.1. Unpacking

The ALVO Conveyor together with the packaging container in which it is delivered weighs 400 kg, the weight of the conveyor without the packaging container is 250 kg.

Do not unpack the ALVO Conveyor outdoors.

To prepare the ALVO Conveyor, proceed as follows:

- Make sure that the ALVO Conveyor in its transport packaging container has been stored for a sufficiently long time in the room (at least three hours) where it is intended to be used.
- Open the packaging and remove the materials protecting the cart.

- Remove and set aside all the supplied accessories.
- Remove the ALVO Conveyor mechanically from the pallet for transport. Do not lift the cart by holding the upper parts.
- Carefully read the operating instructions for the cart.
- Observe the instructions in point 8.2 when starting the cart for the first time.

8.2. First turn on

1. After unpacking, place the ALVO Conveyor on a level surface.
2. Using the power cable, connect the ALVO Conveyor to a power source to charge the batteries.
3. When the batteries are charged (the corresponding LED should turn green), disconnect the conveyor from the power supply and move the ALVO Conveyor to the working position.
4. Move the mechanical wheel lock in the direction marked with an open padlock icon on a green background.
5. Perform the checkout procedures described in Section 15 "Technical condition Monitoring and Inspection" in the following sequence in the presence of a responsible technician:
 - Visual inspection;
 - Functionality check;
 - Do not use the cart if the inspection reveals the deviations specified in Section 15 "Technical condition Monitoring and Inspection";
 - Report any irregularities in the ALVO Conveyor to the manufacturer or his representative;
 - Using a faulty ALVO Conveyor can lead to damage for which the manufacturer cannot be held liable, as a result of erroneous actions of the user.

8.3 ALVO Conveyor move order

On the base of the ALVO Conveyor, there is a set of wheels that facilitate its movement in any direction. The base is also equipped with a mechanical wheel lock to block movement.

To move the ALVO Conveyor, remove the lock by moving the mechanical wheel lock in the direction marked with an open padlock icon on a green background and by holding the handrail with your hands, move the ALVO Conveyor, depending on the required direction of movement, using buttons 1 or 2 (illustration 8, 9) intended for control by drive wheel.

When moving the ALVO Conveyor, be careful not to forget to control the surrounding area. To lock the trolley again, move the mechanical wheel lock in the direction marked with a locked padlock icon on a red background.

9. Charging batteries

The ALVO Conveyor is equipped with a battery charger. The charger is activated by connecting the special socket in the ALVO Conveyor body to a power source using the power cable and switching the power switch on the ALVO Conveyor body from position 0 to position 1.


Illustration 13. Charger function control


Illustration 14. Photo of the power socket, potential equalization slot and switch on and LED

To monitor the battery charging status, an LED indicator is provided on the trolley body.

The LED indicator gives the following signals:

Green light - batteries are charged.

Orange light - the batteries are not fully charged and can be recharged.

Red light - batteries are discharged, they need to be charged immediately.

Charging can also be done when the LED indicator is orange.

When the ALVO Conveyor is connected to a power source, the LED indicator turns green. The batteries should be charged for at least 8 hours. If the batteries are fully charged in less than 8 hours, charging stops automatically. In case of interruption of charging by the user, the battery status is monitored by the LED indicator.

The batteries are rated to operate without recharging for approximately 60 patient movements. Intensive use of electromechanical drive elements contributes to its reduction.

ATTENTION!

It is forbidden to store the ALVO Conveyor in a state with discharged batteries.

When installing new batteries, the entire set must be replaced. Recharging batteries too often will shorten battery life.

10. Antistatic properties of the ALVO Conveyor

The ALVO Conveyor is equipped with special structural elements that ensure safe neutralization of electrostatic discharges, namely antistatic wheels and an antistatic patient transfer device.

If potential equalization is required, the ALVO Conveyor is connected with a cable to potential-free electrical elements, the location of the corresponding connector on the conveyor you can find on the illustration 14.

11. Assessment of the serviceability of the ALVO Conveyor

Before starting work with the ALVO Conveyor, it is necessary to check its condition every time.

Procedure for checking the serviceability of the cart:

1. Place the ALVO Conveyor on a level surface and block its movement. Try to move the ALVO Conveyor with force in any direction. The ALVO Conveyor must remain stationary.
2. Make sure there is no air gap in the mechanical connections that could interfere with the proper functioning of the conveyor.
3. Check the operation of the electromechanical drive by activating the appropriate functions from the control panel.
4. Check the battery charge level with the alarm LED indicators.

The ALVO Conveyor is usable without any abnormalities or faults. If you find them, refer to the troubleshooting section or to the manufacturer, the manufacturer's official representative.

Do not use the cart if it is not fully functional, that is, if it does not function as described in the instruction manual included in the delivery.

The use of a defective ALVO Conveyor can lead to damage for which the manufacturer cannot be held liable as a result of the user's erroneous actions.

12. Crashes and defects

It is the operator's responsibility to promptly report any product malfunctions or faults found in the product to the person responsible for maintaining the equipment. After assessing the nature of the fault and its causes, the person responsible for maintaining the equipment should contact an authorized technician or

the manufacturer for instructions on how to proceed. A product whose operation may be unsafe (with mechanical or electrical failure) should not be used until repair is carried out.

13. Care and disinfection

1. Only use detergents that do not contain active oxygen and chlorine to disinfect the ALVO Conveyor.
 2. After disinfection, the product must be rinsed with distilled water to remove traces of the disinfectant solution and wiped dry with a soft cloth.
 3. Do not disinfect the ALVO Conveyor in the disinfection chamber.
 4. Before disinfection, disconnect the power cable from the ALVO Conveyor.
 5. Do not wash the ALVO Conveyor with a water jet.
 6. Do not use any substances that destroy the structure of the plastic to wash the plastic parts of the
 7. Do not use alcohol-based disinfectants to clean and disinfect the rolling belt.
- Failure to comply with the above requirements will void the product warranty.

14. Maintenance, inspection, repair

All repairs must be carried out by: the manufacturer, authorized technicians or the manufacturer's direct representative.

Users should not attempt to modify or repair the product without special training and authorization. If a written consent is given for the repair to be carried out by the buyer's technicians, the manufacturer must provide the buyer with all the information necessary to carry out the repair.

Use only original spare parts provided by the manufacturer for long-term and trouble-free operation of the cart. Some parts of the product are hazardous to the environment.

Disposal of used parts must be carried out in accordance with current environmental protection regulations.

All operations for repair, inspection and maintenance should be reflected in the repair and maintenance record attached to the operation manual (see Appendix 1).

15. Monitoring the technical condition and inspection the technical condition

To maintain the product in good working order, the user is obliged to periodically inspect it.

Inspection must be performed by an authorized technician or manufacturer's representative.

Further maintenance of the ALVO Conveyor is only permitted if the inspection has been completely passed.

Below in tables 5-7 there is a list of operations that must be performed with an interval of 6 months:

Table 5. Visual inspection.

No.	Action	Description
1.	Checking device markings	All markings are clean, clear and free from possible signs of damage.
2.	Checking the belt surface	Checking the integrity of the surface and the absence of mechanical damage.
3.	Checking moving parts	Pay special attention to places subject to the greatest pressures during maintenance.
4.	Check all removable items	Check whether the product is complete.
5.	Checking surfaces most exposed to disinfectants	In case of visible corrosion, inform staff about the use of prohibited substances containing oxygen or active chlorine.
6.	Checking the potential equalization cable	Cable and clench insulation without mechanical damage.
7.	Checking the condition of gears and toothed belts	Belts and gears are without visible mechanical damage and excessive wear.
8.	Checking drive bolts and bearings	Drive bolts and bearings are without visible mechanical damage. The functioning is quiet, smooth and stable; lubricate as needed.
9.	Visual assessment of the device	Individual parts are without damage that could affect the safety and usability of the device.

Table 6. Functionality check

No.	Action	Description
1.	Checking mobility and brakes	The wheels run smoothly, making a full revolution. The brake levers must be able to lock and unlock without significant resistance.
2.	Checking the electric wheel drive	Activate the corresponding function from the keypad on the handle – quiet, smooth, steady movement in both directions.
3.	Up / Down Movement Check	Activate the corresponding function from the keypad on the handle – smooth movement. The table top is straight when moving. Clear button operation, keyboard without damage
4.	Checking belt movement	Activate the corresponding function from the keypad on the handle - quiet, smooth movement. No lateral movement. The belt is evenly tensioned over the entire range of motion, correct if necessary. Clear operation of buttons, keyboard without damage. LED signaling check.
5.	Checking the safety stop	Checking the operation of gas springs. Checking the holding of a given position.

Table 7. Checking electrical systems.

No.	Action	Description
1.	Checking the status of the protective installation	Test: Protection system test: {limit: <0.200 Ohm} Insulation resistance test: {limit: > 2.000 MOhm} Ground current: {limit: <0.500 mA} SFC ground current: {limit: <1.000 mA} Cover leakage current: {limit: <0.100 mA} SFC cover leakage current: {limit: <0.500 mA} Patient leakage current AC: {limit: <0.100 mA} Patient leakage current AC SFC: {limit: <0.500 mA} Patient leakage current DC: {limit: <10.00 µA} Patient leakage current DC SFC: {limit: <050.0 µA} [Device: Secutest SHI + MH A08 D00 / SR.Nr. YI 522132 0001 IEC 60601.1]
2.	Checking wires	Cable and clench/ plug insulation without mechanical damage. Separate electronics and system components without visible damage / overheating.
3.	Checking the charging system and batteries	Checking the charging voltage and the correct operation of the charger. Battery voltage check. [Device: BRYMEN BM859CF Multimeter]
4.	Checking all electrical modules	Separate electronics and system components without visible damage / overheating.
5.	Checking the condition of electric motors	Separate electronics and system components without visible damage / overheating.
6.	Checking the tape rotation potentiometer	Correct endpoints, adjustments if necessary.
7.	Checking the timer system	Activation button operation, check the shutdown time.

16. Elimination of possible malfunctions

1. The ALVO Conveyor does not perform movements controlled from the control panels:

- Check the battery charge level;
- Check the power cable and control panel cable.

2. The ALVO Conveyor does not move or is in an unstable position:

- Check if the base has bumped into any obstacle;
- Check if any object has fallen under the wheels.

3. If you have any questions, contact the manufacturer for advice.

17. Marking of the ALVO Conveyor

17.1. Factory identification plate of the ALVO Conveyor and its decoding

Explanation of designations:

1. Name, logo and address of the manufacturer;
2. Date of manufacture;
3. Working cycle of the lifting rack;
4. Degree of protection against penetration of foreign objects and water;
5. Name of the medical device, model;
6. CE mark of conformity;
7. Warning - refer to the instruction manual;
8. Input electrical voltage;
9. Electrical safety type B;
10. Disposal instructions in the form of a symbol;
11. Made in Poland.


Illustration 15. Identification plate

17.2. ALVO Conveyor identification plate for Russian Federation

1. Name, logo and location of the manufacturing organization;
2. Name and location of the authorized representative organization;
3. Name of the medical device, model;
4. Serial number;
5. Symbol: "Batch number";
6. Electric voltage, current frequency, power;
7. Safety class;
8. Symbol: "Date of manufacture (month, year)";
9. Registration certificate number;
10. Disposal instructions in the form of a symbol;
11. Symbol "refer to the operating instructions";
12. Degree of protection against penetration of foreign objects and water;
13. Symbol: "Electrical safety type B";
14. Symbol: Warning! Electric shock possible;
15. CE mark of conformity;

Illustration 16. Identification plate layout for deliveries to Russia

16. Inscription: Made in Poland.

17.3 Marking of the transport packaging of the ALVO Conveyor

The transport packaging must contain the following information:

1. Manufacturer: name, logo and location of the manufacturer's organization;
2. Name and location of the authorized representative organization (Supplier);
3. Name of the medical device, model;
4. Serial number;
5. Registration certificate date and number;
6. Storage conditions;
7. Transport conditions;
8. Symbol "Batch number";
9. Symbol: "Keep away from moisture";
10. Symbol: "Top";
11. Overall dimensions of the package indicating the unit of measurement;
12. Net weight;
13. Gross weight;
14. Inscription: Made in Poland.

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17.4 Symbols used

Table 8. Symbols.

Symbols designation			
Serial number	Date of manufacture	Warning, refer to the User manual	Working part type B
Refer to instruction manual	Batch code	Use until (expiration date)	Waste electrical and electronic equipment
CE mark of conformity	Manufacturer	Keep away from	Protect from sunlight

		moisture	
			
Alternating current	Caution. Potential electric shock	Up	Potential equalization connector
Electric shock protection class			Water protection degree
I			IPX4

18. Packing of the ALVO Conveyor

Assembled ALVO Conveyors are packed in a wooden box, accessories and removable parts are packed in individual packaging (polyethylene) inside the main box and protect the products from damage during transportation and storage.

19. Conditions of transportation and storage

19.1 Transport conditions

Transportation is carried out in covered vehicles in accordance with the requirements and rules for the carriage of goods in force for a certain type of transport at a temperature not lower than minus 15 °C and not higher than plus 45 °C and humidity 0-80% without condensation.

When transporting, the ALVO Conveyor must not be stacked in more than one row.

After transportation in freezing temperatures, the ALVO Conveyor in the transport packaging must be kept at room temperature for at least 24 hours. Relative humidity: 10 - 90% without condensation. Pressure: 700 - 1060 hPa.

19.2 Storage conditions

The ALVO Conveyor is stored in a clean, dry place in the manufacturer's packaging at a temperature from - 15 °C to + 45 °C, relative humidity 0 - 80% without condensation, atmospheric pressure 700 hPa 1060 hPa.

It is forbidden to store the ALVO Conveyor near heating devices and in places with high temperatures.

During storage, the ALVO Conveyor must not be stacked in more than one row.

The ALVO Conveyor should be kept out of the reach of children.

The criterion of the limiting state is the economic inexpediency of restoration.

20. Maintenance conditions

Temperature: + 5 °C + 40 °C.

Humidity: 10 - 90%.

Pressure: 700 - 1060 hPa.

21. Expiry date

The service life of the ALVO Conveyor is 10 years from the date of production.

22. Warranty and rules for its provision

22.1 Warranty conditions

The warranty period for the ALVO Conveyor is 2 years from the date of purchase without payment for repairs, if the damage occurred while observing the rules of storage, operation, maintenance and transportation.

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guarantees the conformity of the ALVO® MEDICAL conveyor, model ALVO TRANSPA 4-12, with the declared properties and characteristics during the warranty period, subject to the rules of storage, operation, maintenance and transportation.

The guaranteed shelf life of the ALVO® MEDICAL Mobile conveyor, model ALVO TRANSPA 4-12 in the manufacturer's packaging - 2 years from the date of production.

In the event of a warranty case, please contact the manufacturer directly:

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or an authorized representative: LLC "MNPB Biomir XXI", Russia, 109129, Moscow, 8th Tekstilschikov st., 11, off. 101, E-mail: info@biomir21.ru, Tel. +7 495 228 13 59.

Free warranty service is carried out for 2 years from the date of purchase.

22.2 Exclusions from warranty

The manufacturer is not responsible for:

- 1) Damage due to careless handling or user error.
- 2) Damage caused by the product being updated or repaired by others other than the after-sales service.
- 3) Damage caused by the use of parts or consumables other than the equipment specified by us.
- 4) Damage due to power problems.
- 5) Damage caused by using the device for other purposes.
- 6) Damage caused by natural disaster, for example fire, earthquake, flood, etc.


22.3 Post-warranty repair

After the end of the warranty period, repairs are paid by the buyer. Exceptional issues that are not mentioned in the warranty card are considered on the basis of the current legislation of the Russian Federation.

23. Disposal

When taking the product out of service, the user must disinfect it (in accordance with current environmental protection regulations, a non-disinfected product is classified as hazardous waste). There are three ways to dispose of the product:

1. Disposal by the manufacturer.
2. Disposal by a company licensed to collect equipment for destruction or disposal, observing all necessary measures to protect human life and health and the environment.
3. Disposal by the user, if he has the appropriate technology and resources.

Do not throw this product into the trash can.

24. Authorized representative of the manufacturer

LLC "MNPK Biomir XXI", Russia, 109129, Moscow, 8th Tekstilshchikov st., 11, off. 101, E-mail: info@biomir21.ru, Tel. +7 495 228 13 59.

25. List of national and international regulatory documents / standards that the medical device complies with

25.1 National standards (compliance with RF standards)

- 1) GOST R 50444-2020 "Medical devices, apparatus and equipment. General technical conditions";
- 2) GOST R IEC 60601-1-2010 "Medical electrical equipment. Part 1. General safety requirements considering basic functional characteristics";
- 3) GOST R IEC 60601-1-6-2014 "Medical electrical equipment. Part 1-6. General safety requirements taking into account the main functional characteristics";
- 4) GOST R IEC 62366-2013 "Design of medical devices taking into account serviceability";
- 5) GOST ISO 14971-2011 "Medical devices - Application of risks to medical devices";
- 6) GOST R ISO 15223-1-2020 "Medical devices. Symbols used in marking on medical devices, labels and accompanying documentation."

25.2 International standards

- 1) Council Directive 93/42 / EEC (amended 2007/47 / EC) of 14 June 1993 concerning medical devices;
- 2) EN ISO 13485: 2012 Medical devices - Quality control systems - Regulatory requirements;
- 3) EN ISO 14971: 2012 Medical devices - Application of risk management to medical devices;
- 4) EN ISO 1041: 2008 Information provided by the manufacturer of medical devices;
- 5) ISO 15223-1: 2006 Medical devices - Symbols used for labeling medical devices, labeling and information provided - Part 1: General requirements;

- 6) GOST R IEC 60601-1-2010 "Medical electrical equipment. Part 1. General safety requirements considering the main functional characteristics";
- 7) IEC 60601-1-2: 2007 Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Efficiency Rules - Subsidiary Standard - Electromagnetic Compatibility - Requirements and Testing;
- 8) EN 55011: 2009 / A1: 2010 Electronic equipment for industrial, scientific and medical purposes. High frequency interference characteristics. Limits and methods of measurement;
- 9) GOST 30804.3.2-2013 (IEC 61000-3-2: 2009) "Electromagnetic compatibility of technical equipment. Emission of harmonic current components of technical equipment with a current consumption of not more than 16 A (in one phase). Standards and test methods"
- 10) GOST 30804.3.3-2013 (IEC 61000-3-3: 2008) "Electromagnetic compatibility of technical equipment. Limitations of voltage changes, voltage fluctuations and flicker in low-voltage general-purpose power supply systems. Technical equipment with a current consumption of no more than 16 A (in one phase) connected to the electrical network if certain connection conditions are not met. Standards and test methods";
- 11) IEC 60601-2-2: 2009 Medical Electrical Equipment - Part 2-2: Specific Requirements for General Safety and Performance Practices for High Frequency Electrosurgical Equipment and High Frequency Electrosurgical Accessories;
- 12) ISO 10993-1: 2009 Evaluation of the biological effects of medical devices - Part 1: Evaluation and testing in the risk management process;
- 13) ISO 10993-5: 2009 Evaluation of the biological effects of medical devices - Part 5: In vitro cytotoxicity tests;
- 14) ISO 10993-10: 2002 Evaluation of the biological effect of medical devices - Part 10: Tests for irritability and hypersensitivity of delayed-type skin;
- 15) GOST 30804.4.2-2013 (IEC 61000-4-2: 2008) "Electromagnetic compatibility of technical equipment. Resistance to electrostatic discharge. Requirements and test methods";
- 16) GOST 30804.4.3-2013 (IEC 61000-4-3: 2006) "Electromagnetic compatibility. Immunity to radio-frequency electromagnetic field. Requirements and test methods";
- 17) IEC 61000-4-5: 2005 Electromagnetic compatibility (EMC) - Part 4-5: Testing and measurement methods - Surge immunity test;
- 18) IEC 61000-4-6: 2008 Electromagnetic compatibility (EMC) - Part 4-6: Testing and measurement methods - Immunity to conducted disturbances due to high frequency fields;
- 19) GOST 31204-2013 * "Electromagnetic compatibility. Part 4-8. Test and measurement methods. Resistance to magnetic field of industrial frequency";
- 20) GOST 30804.4.11-2013 (IEC 61000-4-11: 2004) "Electromagnetic compatibility of technical equipment. Immunity to dips, short-term interruptions and changes in power supply voltage. Requirements and test methods";
- 21) EN 60601-1-6: 2010 Medical Electrical Equipment - Part 16 -: General Requirements for Basic Rules and Effectiveness - Secondary Standard: Serviceability;
- 22) EN 62366: 2008 Medical devices - Application of serviceability design to medical devices.

Appendix 1. Record card for repairs and inspections.

Conveyor type Serial number Purchase date

Inspection number	Date of inspection or repair	Inspection type (annual, semi-annual)	Person who performed the inspection or repair	Signature of the person who performed the inspection or repair	Comments on the inspection or repair performed
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					

Appendix A. Declaration of Electromagnetic Compatibility

Medical device Mobile conveyor ALVO® MEDICAL, model: ALVO TRANSPA 4-12. In an electrical device. Every electrical device generates electromagnetic radiation. In addition, all electrical devices are affected by it.

The use of electrical devices requires appropriate precautions related to electromagnetic compatibility.

The use of electrical devices requires appropriate precautions related to electromagnetic compatibility.

The tables in Appendix A "Characteristics of the electromagnetic environment" describe the electromagnetic environment in which the ALVO TRANSPA 4-12 product is designed for.

Follow the instructions and the recommendations given in the tables.

The use of accessories, service equipment, cables and spare parts that are not offered and recommended by the manufacturer may increase the level of electromagnetic radiation and / or reduce the immunity of the product to electromagnetic interference.

Table A 1. Instructions and manufacturer's declaration - electromagnetic emissions.

The ALVO TRANSPA 4-12 medical device is intended for use in the electromagnetic environment specified below. The user of the ALVO TRANSPA 4-12 medical device must ensure that it is used in such an environment.

Radiation type *	Classification	* Explanations and instructions
RF emissions, CISPR 11	Group 1	RF energy is generated by the ALVO TRANSPA 4-12 medical device for internal use only. Therefore, this RF energy remains very low and will not interfere with nearby electronic equipment.
RF emissions, CISPR 11	Class B	The ALVO TRANSPA 4-12 medical device is suitable for use in residential premises and premises connected directly to the low-voltage network through which electricity is supplied to residential buildings.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flickering emissions, IEC 61000-3-3	Compatible	

Table A 2. Instructions and manufacturer's declaration – immunity

The ALVO TRANSPA 4-12 medical device is intended for use in the electromagnetic environment specified below. The user of the ALVO TRANSPA 4-12 medical device must ensure that it is used in such an environment.

Immunity	Test level IEC 60601-1-2	Immunity level	Explanations and instructions
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV, contact ± 8 kV, in air	± 6 kV, contact ± 8 kV, in air	The floors in the place of use of the product must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Surge current IEC 61000-4-5	± 1 kV, "wire-wire" in the power line ± 2 kV, "wire-land" in the power line	± 1 kV, "wire-wire" in the power line ± 2 kV, "wire-land" in the power line	The power supply parameters and the existing level of interference should be comparable to those typical for commercial or hospital settings.

Rapid electrical transient, IEC 61000-4-4	± 2 kV, for power lines ± 1 kV, for signal lines	± 2 kV, for power lines ± 1 kV, for signal lines	The power supply parameters and the existing level of interference should be comparable to those typical for commercial or hospital settings.
Voltage dips, off-peaks and changes in power lines, IEC 61000-4-11	$<5\% U_T$ (U_T dip $> 95\%$) for 0.5 period $40\% U_T$ (U_T dip)		The power supply parameters and the existing level of interference should be comparable to those typical for commercial or hospital settings.
	60% for 5 periods 70% U_T (U_T dip 30%) for 25 periods $<5\% U_T$ (U_T dip $> 95\%$) for 5 periods		

NOTE: U_T is the rated AC voltage.

Table A 3. Instructions and manufacturer's declaration – immunity

The ALVO TRANSPA 4-12 medical device is intended for use in the electromagnetic environment specified below. The user of the ALVO TRANSPA 4-12 medical device must ensure that it is used in such an environment.			
Immunity	Test level IEC 60601-1-2	Immunity level	Explanations and instructions
Conducted emissions in radio frequency electromagnetic fields, IEC 61000-4-6	3 V rms 150 kHz - 80 MHz	3 V rms	Portable RF communications equipment should not be used near any part of the cart, including cables, if the distance is less than the recommended separation distance calculated using the appropriate formula based on the frequency of the transmitter. Recommended removal distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{PP}$ 80 MHz – 800 MHz $d = 2,3 \sqrt{PP}$ 800 MHz – 2,5 GHz where: P is the maximum output power of the transmitter in watts (W) according to the manufacturer's specifications. d is recommended removal distance in meters (m). ^a The intensity of electromagnetic fields generated by radio transmitters measured at the installation site of the product should not exceed the product's immunity level for each frequency range. ^b Interference may occur with a medical device in the vicinity of equipment marked with the following symbol. It is necessary to control the correct functioning of the medical device near the equipment marked with the following symbol. 
RF electromagnetic fields, IEC 61000-4-3	3 V / m 80 MHz - 2.5 GHz	3 V / m	

Table A 4.

Recommended separation distance between the product and portable radio-communication devices

Maximum rated output power of the transmitter, W	150 kHz - 80 MHz $d = 1,2 \sqrt{PP}$ Distance, m	150 kHz - 800 MHz $d = 1,2 \sqrt{PP}$ Distance, m	800 MHz - 2,5 GHz $d = 2,3 \sqrt{PP}$ Distance, m
0,01	0,1	0,1	0,2
0,1	0,4	0,4	0,7
1	1,2	1,2	2,3
10	4	4	7
100	12	12	23

For transmitters whose maximum output power is not listed above, the separation distance is calculated using the presented formulas. P is the transmitter power in watts (W) according to the manufacturer's data.

ATTENTION!

Under certain conditions, the above guidelines may not apply. The propagation of electromagnetic radiation is influenced by the reflective and absorption properties of the objects through which it passes.

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