
“APPROVED”
By the Director of
Vacutest Kima S.r.l.
Monica Vianello

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October 01, 2020

INSTRUCTIONS FOR USE OF THE MEDICAL DEVICE

VACUTEST Transfer Straw for Urine Collection, single use

2020

NAME OF THE MEDICAL DEVICE

VACUTEST transfer straw for urine collection, single use, in design variants (hereafter referred to as the VACUTEST transfer straw for urine collection, transfer straw(s))

DESIGN VARIANTS

1. VACUTEST transfer straw for urine collection, single use, short straw - 100 pcs/pack.
2. VACUTEST transfer straw for urine collection, single use, long straw - 50 pcs/pack.

DESCRIPTION

VACUTEST transfer straw for urine collection is a disposable non-sterile device made of semitransparent polypropylene. The transfer straw is designed for a fast transfer of a sample from any container to a vacuum tube. It can be used with urine collection containers and any other tanks. The transfer straw can be used with one or more vacuum tubes for the same patient. The straw of the device is placed in a container, and due to the fact that there is a needle closed with a rubber valve inside the device, when the tube cap is pierced, urine is quickly, safely, without direct contact with the medical staff and environment, transferred into the vacuum tube.

INTENDED USE

VACUTEST transfer straw is intended for transferring the patient's urine into a vacuum tube for the subsequent laboratory tests, as an auxiliary for in vitro diagnostics.

VACUTEST transfer straw can be used as an auxiliary in conjunction with the corresponding vacuum tubes for common (clinical) urine analysis or for microbiological urine analysis.

SCOPE OF APPLICATION

The transfer straw is intended for use in medical treatment and prevention institutions, and laboratories as intended.

INDICATIONS FOR USE

Use as intended for all patient groups to transfer a urine sample from a container or any other tank to a vacuum tube for further studies.

LIST OF POSSIBLE STUDIES FOR WHICH THE DEVICE IS USED

The list (including but not limited) is provided in Table 1.

Table 1

No.	Design variant	Ref. No.	Studies
1.	VACUTEST transfer straw for urine collection, single use, short straw - 100 pcs/pack.	21810	For clinical and microbiological studies, in conjunction with suitable vacuum tubes
2.	VACUTEST transfer straw for urine collection, single use, long straw - 50 pcs/pack.	21840	For clinical and microbiological studies, in conjunction with suitable vacuum tubes

LIMITATIONS OF USE (CONTRAINDICATIONS)

- Use only as intended;
- Do not use for transferring other biological fluids;
- Do not use if there are signs of use and damage to any part of the transfer straw.

POTENTIAL USERS

Medical staff.

SAFETY MEASURES

- Do not use after expiry date.
- Do not use if there are signs of use and damage to any part of the transfer straw.
- Do not reuse. Single use only!
- Do not use if there is a sharp obtrusive edge, protrusion or roughness on the surface of the device that could cut, pierce or tear the gloves or the user's finger.
- Do not ignore cases of contact with biological samples in order to avoid the transmission of infectious agents.

CONDITIONS OF USE

The device should be used at temperatures from +5°C to +25°C and relative humidity not higher than 75%.

COMPATIBILITY

The transfer straws are compatible with containers without integrated urine collection device and any other tanks. Requirements to the urine collecting container: height from the edge to the bottom of the container is no more than 8 cm (for ref. No. 21810) and no more than 16 cm (for ref. No. 21840).

The transfer straws are compatible with vacuum tubes supplied by the main manufacturers. The list of compatible tubes is provided in table 2.

Table 2

Seq. No.	Ref. No.	Description	Manufacturer	Data on registration
VACUTEST tube for urine collection, sterile, single use			Vacutest Kima S.r.l., Italy	It is registered in parallel with the device. Information about the RC will be added at the end of the registration procedure.
1.	14900	Sterile VACUTEST tube for urine collection, single use, 13x75 mm, 4 ml, round bottom, no additive		
2.	14911	Sterile VACUTEST tube for urine collection, single use, 13x100 mm, 6 ml, round bottom, no additive		
3.	148201	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 6 ml, round bottom, no additive		
4.	14930	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 9 ml, round bottom, no additive		
5.	14936	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 9 ml, round bottom, no additive		
6.	14850	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 9.5 ml, conical bottom, no additive		
7.	14851	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 9.5 ml, conical bottom, no additive		

Seq. No.	Ref. No.	Description	Manufacturer	Data on registration
8.	14750	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 11 ml, round bottom, no additive		
9.	14762	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 11 ml, round bottom, no additive		
10.	149610	Sterile VACUTEST tube for urine collection, single use, 13x75 mm, 4 ml, round bottom, with boric acid		
11.	149600	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 9 ml, round bottom, with boric acid		
12.	149685P	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 9 ml, round bottom, with boric acid		
13.	148600	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 9 ml, conical bottom, with boric acid		
14.	147700	Sterile VACUTEST tube for urine collection, single use, 16x100 mm, 11 ml, round bottom, with boric acid		
Devices for urine collection, transportation and analysis, single-use BD Vacutainer in kits and single packs			Becton Dickinson and Company, USA	RC No. FSZ 2007/00888 dd. 28.01.2016
1	368500	Vacutainer urine tube, no additive, 4 ml, 13x75 mm		
2	368501	Vacutainer urine tube, no additive, 6 ml, 13x100 mm		
3	364992	Vacutainer urine tube, with additive, 8 ml, 16x100 mm		
4	365017	Vacutainer urine tube, with additive, 8 ml, 16x100 mm		
5	365000	Vacutainer urine tube, no additive, 9.5 ml, 16x100 mm		
6	364938	Vacutainer urine tube, no additive, 10 ml, 16x100 mm		
7	364915	Vacutainer urine tube, no additive, 11 ml, 16x100 mm		
8	364958	Vacutainer urine tube, with additive, 4 ml, 75 mm		
9	364955	Vacutainer urine tube, with additive, 10 ml, 16x100 mm		
Vacurette system for vacuum urine collection, with accessories			Greiner Bio-One GmbH, Austria	RC No. FSZ 2011/09315 dd. 23.01.2015
1	456007	Vacurette urine tube, no additive, 6.5 ml, 13x100 mm		
2	455028	Vacurette urine tube, no additive, 9.5 ml, 16x100 mm		
3	455007	Vacurette urine tube, no additive, 10 ml, 16x100 mm		
4	456009	Vacurette urine tube, with additive, 6.5 ml, 13x100 mm		
5	454239	Vacurette urine tube, with additive, 4 ml, 13x75 mm		
6	455003	Vacurette urine tube, with additive, 10 ml, 16x100 mm		
7	455029	Vacurette urine tube, with additive, 9.5 ml, 16x100 mm		
8	455049	Vacurette urine tube, with additive, 9.5 ml, 16x100 mm		
9	455048	Vacurette urine tube, with additive, 10 ml, 16x100 mm		

Seq. No.	Ref. No.	Description	Manufacturer	Data on registration
		mm		
Devices for laboratory studies made of polymer materials			Ningbo Greetmed Medical Instruments Co., LTD., Chine	RC No. FSZ 2012/11857 dd. 28.03.2012
1	12007050	Urine tube, no additive, 9 ml, 16x100 mm		
2	12007051	Urine tube, with additive, 9 ml, 16x100 mm		
3	12007052	Urine tube, no additive, 9 ml, 16x100 mm		
Vacurette system for vacuum urine collection, with accessories			Greiner Bio-One GmbH, Austria	RC No. FSZ 2011/09315 dd. 23.01.2015
1	456007	Vacurette urine tube, no additive, 6.5 ml, 13x100 mm		
2	455028	Vacurette urine tube, no additive, 9 ml, 16x100 mm		
3	455007	Vacurette urine tube, no additive, 10 ml, 16x100 mm		
4	455049	Vacurette urine tube, with additive, 9 ml, 16x100 mm		
5	455048	Vacurette urine tube, with additive, 10 ml, 16x100 mm		
Medical polymeric single-use devices for collection and analysis of biological materials in laboratories V-Monovette			Sarstedt AG&Co., Germany	RC No. FSZ 2009/04699 dd. 18.10.2017
1	11.2252.0 01	V-Monovette Vacuum Urine Collection Tubes, no additive, 4 ml		
2	11.2352.0 01	V-Monovette Vacuum Urine Collection Tubes, no additive, 6 ml		
3	11.2452.0 01	V-Monovette Vacuum Urine Collection Tubes, no additive, 10 ml		
4	11.2552.0 01	V-Monovette Vacuum Urine Collection Tubes, no additive, 9.5 ml		
5	11.2253.0 01	V-Monovette Vacuum Urine Collection Tubes, with additive, 4 ml		
6	11.2453.0 01	V-Monovette Vacuum Urine Collection Tubes, with additive, 10 ml		
Devices for laboratory studies made of polymer materials			Ningbo Greetmed Medical Instruments Co., LTD., China	RC No. FSZ 2012/11857 dd. 28.03.2012
1	12007050	Urine tube, no additive, 9 ml, 16x100 mm		
2	12007051	Urine tube, with additive, 9 ml, 16x100 mm		
3	12007052	Urine tube, no additive, 9 ml, 16x100 mm		

INDICATIONS FOR USE

- Check the integrity and absence of damage to the group package by visual inspection, if the package is damaged, the integrity of the contents should be checked before use.
- Check the expiration date of the device.
- Open the cardboard box, unwrap the plastic bag.
- Take one transfer straw at a time directly from the box as needed.

CAUTION! The transfer straws are disposable devices.
They can be used with one or more vacuum tubes for the same patient.

USE OF THE DEVICE



Recommendations for use

Use of the transfer straw

- Thoroughly clean your hands before use.
- Use personal protective equipment (gloves, mask).
- Prepare all necessary items (transfer straw, tube, container with urine).

Transfer of the urine to the tube using the transfer straw

1. Before transfer of urine to the tube, gently homogenize the urine sample (mix the urine in the closed container or other tank by circular motions).
2. Open the container's/tank's cap. Place the transfer straw into the container/tank.
3. Insert the vacuum tube into the urine transfer straw.
4. Gently press the tube to allow the cap piercing. Wait for the complete filling of the tube (stop of the flow).
5. Pull the tube out of the device. For tube with additive: Gently invert the sample tube 6-8 times to dissolve and mix the additive. (If necessary, insert the next tube into the transfer straw).
6. Close the container's/tank's cap.

GENERAL VIEW OF THE DEVICES, DESCRIPTION

General view of the device and principles of operation

The transfer straw is a holder for a vacuum tube, there is a needle inside the holder, on the outer side the holder's luer adapter is equipped with a short or long immersion tube, which is placed in a container or tank with urine. Due to the fact that there is a needle closed with a rubber valve inside the device, when the tube cap is pierced, urine is quickly, safety, without direct contact with the medical staff and environment, transferred into the tube under vacuum.

See Figure 1.


Figure 1.

General view of the devices is given in table 3.

Table 3

	Design variant	Ref. No.	Weight of the transfer straw, g	View
1.	VACUTEST transfer straw for urine collection, single use, short straw	21810	2.6	
2.	VACUTEST transfer straw for urine collection, single use, long straw	21840	2.8	

REQUIREMENTS TO MANUFACTURE PROCESS

- On the surface of the transfer straw, there should be no sink marks, burrs, traces from the mold connector with a height of more than 0,5 mm and foreign inclusions protruding above the surface.
- Expected expiry date must be 3 years from the manufacturing date.
- The devices should be resistant to mechanical and climatic effects during operation, transportation and storage.
- Materials, safety of which is confirmed in the prescribed manner, must be used during manufacture of the devices.
- The devices should be resistant to the effects of biological fluids and secretions of body tissues with which they come in contact during operation.

Material

Materials used in the manufacture are presented in table 4.

Table 4

Intended use	Material
Luer Adapter	Polypropylene (yellow) Note-dyes: Yellow -YELLOW PE 3/R067-AT
Cannula (18G):	Siliconized with single tip & beige rubber back-flow protection Cannula material: Stainless steel (AISI 304) Rubber protection material: 100% isoprene rubber Adhesive: Epoxy bonding Lubricant: silicone fluid
Holder	Polypropylene, uncolored
Straw collection	Polypropylene, uncolored

TECHNICAL SPECIFICATIONS

Weight and main dimensions of the devices are provided in table 5.

Table 5

Parameter	Data			
	21810		21840	
Dimensions, mm	Length, mm	Diameter, mm	Length, mm	Diameter, mm
Total length of transfer straw	138.2±5.0	-	212.2±5.0	-
Holder	52.20 ± 0.5 (total)	20.00 ± 0.5 (cylindrical part) 30.00± 0.5 (base)	52.20± 0.5 (total)	20.00 ± 0.5 (cylindrical part) 30.00± 0.5 (base)
Luer Adapter	14.1± 0.5	6.00± 0.15	14.1± 0.5	6.00± 0.15
Cannula (18G)	18.00± 0.5	18G (1,2 mm)	18.00± 0.5	18G (1.2 mm)
Straw collection	86.00± 5.0	2.40± 0.15	160.00± 5.0	2.40± 0.15

SCOPE OF SUPPLY

	Design variant	Ref. No.	Packaging	
1.	VACUTEST transfer straw for urine collection, single use, short straw	21810	100 pcs	Instructions for Use, 1 pc per shipping case
2.	VACUTEST transfer straw for urine collection, single use, long straw	21840	50 pcs	Instructions for Use, 1 pc per shipping case

STERILIZATION

Non-sterile.

SHELF LIFE

Shelf life is 3 years from the manufacturing date.

STORAGE AND TRANSPORTATION CONDITIONS

- Transportation and storage can be carried out by any means of transport in closed vehicles in accordance with the rules applicable to each mean of transport and approved in the prescribed manner. Transportation conditions: from -5 to +50 °C.
- Storage and preservation of transfer straws for a long time is carried out in heated spaces at temperatures from +5 to +25 °C in a dry place.

DISPOSAL

- The used devices shall be disposed of in accordance with current legislation.
- The used devices are classified as epidemiologically hazardous waste (class B) and must be collected and disposed of in accordance with current legislation (for the Russian Federation in accordance with requirements of SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public

premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" dd. 29.01.2021).

- The devices unused for their intended purpose due to the expiration of the shelf life or other reasons are classified as epidemiological safe waste (class A), similar in composition to solid household waste, and disposed of as household waste.

WARRANTIES OF THE MANUFACTURER

- The manufacturer warrants that the transfer straws conform to the requirements if the use observes conditions of transportation and storage, as well as operating rules.
- Warranty shelf life is 3 years from the manufacturing date.

TECHNICAL SERVICE AND REPAIR

The transfer straws are disposable devices. They are not a subject to technical service and repair.

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AUTHORIZED MANUFACTURER REPRESENTATIVE IN THE RUSSIAN FEDERATION

Limited Liability Company "BioSystems" ("BioSystems" LLC)

Legal Address: Russia, 197022, Saint Petersburg, inner-city municipality Aptekarsky Ostrov, 23 Professora Popova Str., Let. E, premises 3H/5H, room 207

Post Address: Russia, 197022, Saint Petersburg, 23 Professora Popova Str., Let. E

Tel.: (812) 320-49-49

E-mail: main@bioline.ru

LIST OF STANDARDS APPLIED BY THE MANUFACTURER OF THE MEDICAL DEVICE

UNI EN ISO 9001:2015 Quality Management System. Requirements;

UNI EN ISO 13485:2016 Medical devices - Quality management systems - Requirements for regulatory purposes;

CE: quality guarantee system through the CE Declaration of Conformity issue after preparation of technical files available for competent authorities, and according to CE 98/79/CE Directive (Italian legislative decree 08/09/2000 No.332)

UNI EN 928 In Vitro Diagnostic Systems - Directions for application of EN 29001, EN 46001 and EN 29002 E EN 46002

UNI EN 980 Symbols used for labeling medical devices.

UNI EN 1041 Information supplied by the manufacturer of medical devices

UNI EN 14971 – Application of risk management to medical devices

PACKAGING

Consumer package

The transfer straws in plastic bag are placed in pieces in a cardboard box (100 pcs/pack. (ref. No. 21810) and 50 pcs/pack. (ref. No. 21840)). There are warning symbols on the package.

Shipping case

Cardboard box containing 10 consumer packaging boxes.

Dimensions and weight of the package are given in table 6.

Table 6

Package type / ref. No.	Dimensions, L x W x H, mm		Weight, g	
	21810	21840	21810	21840
Consumer package	210*155*155	210*155*155	96 +/- 7%	96 +/- 7%
Shipping case	330*330*443	330*330*443	874 +/- 7%	874 +/- 7%

MARKING OF THE DEVICE

The marking should be legible, clear and resistant to external exposure.

Marking on a consumer package

Manufacturer marking:

The marking include: reference number, lot and expiry date, quantity, barcode, manufacturer.

For marking in the Russian Federation, consumer packaging has an additional marking in the form of a sticker:

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
Shipping case marking:

Manufacturer marking:

For marking shipping case in the Russian Federation, additional marking (located on the shipping case) is used:

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

IMPORTANT SYMBOLS ON MARKING

Symbol	Meaning
	Manufacturer
	Best before
EXP.	Best before
	Lot number
	Reference number
	IVD medical device
	Mandatory conformity mark for products placed on the market in the European Economic Area
	No reuse

Symbol	Meaning
	Manufacturing date
	Keep dry
	Avoid sunlight
	Top
	Bottom
	Side
	Take care! See Instructions for Use
	Temperature range
	Fragile Handle with care
	Country of origin: Italy
	Keep out of the reach of children
	Take care! Sharp needle inside

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