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

«I certify»

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Managing Director
botiss biomaterials GmbH
Dr. Drazen Tadic

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Instructions for use
cerabone® plus natural bovine bone substitute material with hyaluronate for filling of bone
defects

cerabone® plus**INSTRUCTIONS FOR USE**

Manufacturer:
 botiss biomaterials GmbH
 Hauptstraße 28
 15806 Zossen, Germany
 Tel.: +49 33769 / 88 41 985
 Fax: +49 33769 / 88 41 986
www.botiss.com

Rev: B01CBp-04RU/2023-10

Material

Hydroxyapatite ceramic (cerabone® granules) with the excipient sodium hyaluronate.

Composition

1 cm³ ceramic is equivalent to an average of 0.6 - 1.1 g hydroxyapatite (pentacalcium hydroxide [tris] phosphate), depending on the porosity of the ceramic.

The granules are provided in the form of an interconnecting macroporous and microporous hydroxyapatite ceramic, which is manufactured from bovine spongiosa in a high-temperature process (T > 1200 °C) lasting several hours. The macroporosity of the ceramic is in the range of 65 – 80 Vol.% and the pore width is in the range of approximately 100 - 1500 µm.

Other components

According to the specification and ASTM standard (F 1185-03), the granules contain permissible proportions of, e.g., other calcium phosphates and trace elements. Sodium hyaluronate of non-animal origin is used as an excipient. It is obtained biotechnologically by fermentation.

Dosage forms and contents

REF	Granule size	Volume
1810	0,5 - 1,0 mm	1 x 0,5 ml
1811	0,5 - 1,0 mm	1 x 1,0 ml
1820	1,0 - 2,0 mm	1 x 0,5 ml
1821	1,0 - 2,0 mm	1 x 1,0 ml

The sterile bone substitute material is delivered in a double blister.

Completeness:

The cardboard box contains a material of the appropriate size and volume – 1 pc., instructions for use - 1 pc. and patient labels – 4 pcs.

Material class

Bone substitute material made of calcium phosphate ceramic

Responsible manufacturer

botiss biomaterials GmbH, Hauptstraße 28, 15806 Zossen, Germany

Intended use

cerabone[®] plus is intended for the permanent filling or reconstruction of aseptic bone defects and for volume augmentation in autologous spongiosa transplantation in dental, oral and maxillofacial surgery.

Indications

- Alveolar ridge augmentation/reconstruction
- Filling of bone defects (including after root resection, apicoectomy or cystectomy)
- Filling of extraction sockets to support alveolar ridge preservation
- Sinus floor elevation
- Filling of periodontal bone defects
- Filling of extraction sockets as part of immediate implantations
- Filling of peri-implant bone defects

Mode of action

cerabone[®] plus consists of a hydroxyapatite ceramic (cerabone[®] granules) and the excipient sodium hyaluronate. The contained sodium hyaluronate is a very hygroscopic material. Due to this property of the material, a viscous solution forms upon contact with aqueous solutions, which, following hydration and careful mixing, leads to the granules being bound into a connected mass with malleable texture. This improves handling and facilitates the application of the material to the bone defect. After application, the sodium hyaluronate is rapidly degraded in the human body (within less than 2 (two) weeks). The granules have a structure of interconnecting macro and micro pores. Due to the porous structure, it serves as a scaffold for blood vessels and regenerative cells after implantation, thereby supporting bone regeneration. The bone substitute material is integrated into the newly forming bone within about 6-9 months. The granulate is a volume-stable material that is not or only very slowly resorbed by the body. It is therefore very well suited for the permanent filling or reconstruction of aseptic bone defects.

Dosage instructions, method and duration of administration

cerabone[®] plus must be hydrated before use (see table below). **Approx. 0.5 ml of liquid (corresponds to about 10 - 12 drops) must be added to 1 ml of bone substitute material.** Hydration can be performed with sterile saline solution or patient blood.

Hydration Protocol

REF	Granule size	Volume	Hydration with
1810	0.5 - 1.0 mm	0.5 ml	approx. 0.25 ml liquid
1811	0.5 - 1.0 mm	1.0 ml	approx. 0.5 ml liquid
1820	1.0 - 2.0 mm	0.5 ml	approx. 0.25 ml liquid
1821	1.0 - 2.0 mm	1.0 ml	approx. 0.5 ml liquid

After adding the liquid, it should be mixed thoroughly with the material until it forms a connected mass with a malleable texture that can be applied. Insufficient or excessive liquid can result in a failure to form the described texture or a loss thereof. If the product appears dry following hydration and no connected mass has formed, additional saline solution can be added dropwise and mixed with the product until the desired texture is obtained. When working with autologous bone, harvested bone chips can be added to cerabone[®] plus before hydration is performed. In this case, the amount of liquid to be added may vary. The material does not harden.

Amount and frequency of administration

The amount of bone substitute material depends on the size of the bony defect. cerabone[®] plus is intended for single use and is implanted during the surgical procedure. The hydroxyapatite ceramic is intended for permanent implantation in the patient.

The maximum dose is 17 ml.

Method and time of administration

The prerequisite for successful use is the stable implantation of cerabone® plus into a vital, well-perfused bone bed with moderate biomechanical stress. cerabone® plus is particularly suitable for filling cancellous bone defects and should always be implanted in direct, preferably all-round contact with the natural bone. Excessive liquid at the augmentation site, i.e. in case of profuse bleeding, may reduce the handling of the bone substitute material during application and its adaptation to the defect. It is recommended to remove excess liquid from the defect site prior to the application.

In order to ensure stable implantation, the bone defect should be filled up with cerabone® plus as completely as possible. The implant bed is prepared by freshening the adjacent bone. Each surgical procedure is determined by the location and the type and extent of the bone defect. The goal of stable intraosseous implantation is a direct bony contact between the ceramic and the natural bone without any intermediate connective tissue layer. The interconnecting porous structure of the cerabone® plus granules ensures their function as an osteoconductive matrix and/or scaffold for the ingrowth of new bone, so that a solid bond can develop between the granules and the adjacent bone.

To support bony regeneration, cerabone® plus and the augmentation area should be covered with a barrier membrane for guided bone regeneration (GBR) and guided tissue regeneration (GTR). Any macro- or microparticles that may have split off into the soft tissue at implantation are to be removed.

After augmentation with cerabone® plus, a healing time of about 6-9 months should be allowed before inserting the implant.

Contraindications

Absolute contraindications

- Acute and chronic flord infections in the implant area

Relative contraindications

- Severe endocrine bone diseases
- Severe metabolic bone disorders
- Pre-existing therapy with glucocorticoids, mineralocorticoids and drugs that interfere with the calcium metabolism (e.g. calcitonin)
- Severe or poorly controlled diabetes mellitus
- Immunosuppressive therapy
- Malignant tumour diseases (due to more difficult X-ray diagnostics of tumour recurrences in the implant area)

Side effects and potential complications

A successful result may not always be achieved. A surgical error, certain medical circumstances or failure of the bone substitute material may require a secondary operation to remove or replace the introduced bone substitute material.

Potential complications which may occur with any surgery include transient discomfort or pain, swelling at the surgical site, flap sloughing, bleeding, wound dehiscence, local inflammation, infection, bone loss, or dental implant loss. Potential complications of an augmentation procedure may include material displacement, perforation of mucous and soft tissues, fast material resorption, insufficient material integration or insufficient bone formation.

The hydroxyapatite ceramic in cerabone® plus is manufactured in a high temperature sintering process lasting several hours at temperatures of > 1200°C. This high temperature processing results in all organic components, including potential infectious agents such as prions¹, being eliminated (deproteinization), so that the risk of occurrence of allergic side effects or rejection reactions is

reduced to a minimum.

Allergic reactions or symptoms of intolerance to the sodium hyaluronate in cerabone[®] plus are extremely rare but cannot be ruled out.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authorities.

Notes

Hydroxyapatite ceramics have proven to be very well-tolerated (biocompatible) bone substitute materials, both in vitro (cell tests) and in vivo (experimental, clinical). Experimental mutagenicity/genotoxicity studies did not indicate the presence of any mutagenic/genotoxic properties.

In the scientific literature, hydroxyapatite ceramics are generally considered to be non-resorbable. They are radiopaque. The ceramic does not pose any safety risk in a nuclear magnetic resonance imaging environment.

It gives no signal in magnetic resonance imaging and can, if required, be imaged only by means of special methods (STIR sequence, contrast medium administration).

cerabone[®] plus consists of a hydroxyapatite ceramic (cerabone[®] granules) and hyaluronate, which is added as an excipient to improve the handling of the bone substitute. This characteristic of cerabone[®] plus, which is caused by the water-binding activity of hyaluronate, has been verified by usability tests and clinical trying. Based on the clinical data of cerabone[®] plus, there is no evidence for a metabolic, immunological or pharmacological effect of the contained hyaluronate. Nevertheless, published preclinical data on the device has demonstrated a potential beneficial effect on regenerative cells types and angiogenesis, which is attributed to the contained hyaluronate, however, with the exact underlying mechanisms being unknown.

Use during pregnancy and lactation

No studies are available on the use of cerabone[®] plus during pregnancy or while breastfeeding, or on its effect on human fertility. Before using cerabone[®] plus, the treating physician must therefore evaluate the possible risks to the child against the benefits for the mother on a case by case basis.

Use in children and the elderly

There have been no investigations into the use of cerabone[®] plus in children. No special precautions for use have to be considered in elderly people.

Precautions for use and warnings

It is recommended to check the correct implantation with appropriate medical imaging procedures.

Interactions with other medicinal products

Apart from those listed under Contraindications, no interactions with medicinal products or other medical devices are known to date.

Duration of use

The hydroxyapatite ceramic is intended for permanent (life-long) implantation. The excipient is broken down enzymatically after the implantation.

Operational conditions

cerabone[®] plus should only be used by doctors who are familiar with surgical intervention techniques in dental, oral and maxillofacial surgery.

Ambient temperature 22±10°C

Relative humidity not controlled

Atmospheric pressure not controlled
 Implantable product: internal environment of the body
 Temperature from 32 °C to 42 °C
 A relative humidity of 100 %

Removal and replacement

Replacement is not required, the auxiliary substance is degraded after implantation.
 If removal is necessary (for example, in case of infection or other complications of the healing process), complete removal of the implant should be ensured and thorough surgical treatment of adjacent tissues should be performed to obtain a well vascularized surface.

Shelf life

The expiry date is printed on the folding box as well as on the secondary and primary packaging. cerabone[®] plus must no longer be used after the stated expiry date.

Restrictions or loss of usability even prior to expiration date

If the primary packaging (inner blister) or the secondary packaging (outer blister) is damaged, the sterility of the bone substitute material is no longer guaranteed. In such case, cerabone[®] plus must no longer be used.

The contents of unused, but opened or damaged packaging must not be re-sterilized and are to be discarded.

Storage

cerabone[®] plus should be stored in a dry place, protected from light exposure, at a temperature between +5°C and +25°C. Avoid any contact with heat sources and/or storage in direct sunlight, as this can cause damage to the packaging and thus a loss of sterility.

Transportation

Packed by the manufacturer devices are transported by all modes of transport in closed vehicles in accordance with the rules of carriage of goods, operating in each mode of transport.

- Lower temperature limit: - 20°C, for a maximum of 5 days
- Upper temperature limit: + 40°C, for a maximum of 72 h.

Sterilization procedure

cerabone[®] plus is provided sterile in a double blister. The contents of both the inner blister and those of the outer blister are sterile and are contained in a non-sterile outer cardboard box. Gamma rays are used for sterilization.

Safe disposal

Any remaining product residues are to be disposed of with the practice or hospital waste. Products with expired shelf life are classified as medical waste of class according to SanPiN 2.1.3684-21 and must be disposed of accordingly. After contact with the patient's body, the components of the product should be considered medical waste of class B according to SanPiN 2.1.3684-21 and disposed of accordingly.

Date of revision of the text

10/2023

Symbols used

STERILE R	Sterilized using irradiation
LOT	Batch code

	Catalogue number
	Manufacturer
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Consult instructions for use
	Use-by date
	Temperature limit
	Keep away from sunlight
	Keep dry
	Conformity mark
	Registered trademark

Any serious adverse reactions observed during or after clinical application, which may be linked to cerabone[®] plus, should be immediately reported to:

botiss biomaterials GmbH

Hauptstraße 28, 15806 Zossen, Germany

Tel.: +49 33769 88 41 985

Fax: +49 33769 88 41 986

Email: Siba@botiss.com

Warranty

botiss biomaterials GmbH guarantees the Buyer that "cerabone[®] plus" (the product) purchased by the Buyer does not have material and quality defects at the time of delivery to the Buyer. The company "botiss biomaterials GmbH" does not give any guarantee (direct, indirect or statutory) to products that underwent a non-standard mechanical loading, storage in improper conditions, inappropriate usage/misuse, incorrect or negligent application, improper testing, usage in conjunction with other products or components other than those for which "cerabone[®] plus" (the product) has been developed, or underwent a usage in a way or in medical procedures which are not intended for this product. Guarantee claims must be made without delay. Any guarantee claim is limited to the replacement of the defective product(s) and payment of shipping charges to the Buyer of the replaced product(s). In particular, "botiss biomaterials GmbH" does not compensate any additional components, tools, treatment costs and other costs and expenses happened due to or in connection with the replacement of the defective product to the Buyer or the patient.

Manufacturer's authorized representative:

Botiss Plus Ltd.,

Russia, St Peterburg, 194291, Lunacharskogo pr. 47, Lit. A

Tel. +7812558-15-15

cerabone[®] plus natural bovine bone substitute material with hyaluronate for filling of bone defects complies with the following standards as far as applicable

DIN EN 556-1	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
DIN EN 1041	Information supplied by the manufacturer of medical devices
DIN EN 62366-1 (VDE 0750-241-1)	Medical devices - Part 1: Application of usability engineering to medical devices
DIN EN ISO 7405	Dentistry - Evaluation of biocompatibility of medical devices used in dentistry
DIN EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
DIN EN ISO 10993-2	Biological evaluation of medical devices - Part 2: Animal welfare requirements
DIN EN ISO 10993-3	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity
DIN EN ISO 10993-5	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
DIN EN ISO 10993-6	Biological evaluation of medical devices - Part 6: Tests for local effects after implantation
DIN EN ISO 10993-9	Biological evaluation of medical devices - Part 9: Framework for identification and quantification of potential degradation products
DIN EN ISO 10993-10	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
DIN EN ISO 10993-11	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
DIN EN ISO 10993-12	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
DIN EN ISO 10993-17	Biological evaluation of medical devices - Part 17: Establishment of allowable limits for leachable substances
DIN EN ISO 10993-18	Biological evaluation of medical devices - Part 18: Chemical characterization of medical device materials within a risk management process
DIN EN ISO 11137-1	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
DIN EN ISO 11137-2	Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose
DIN EN ISO 11137-3	Sterilization of health care products - Radiation - Part 3: Guidance on dosimetric aspects of development, validation and routine control
DIN EN ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
DIN EN ISO 11607-2	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
DIN EN ISO 11737-1	Sterilization of health care products - Microbiological methods - Part 1: Determination of a population of microorganisms on products
DIN EN ISO 11737-2	Sterilization of health care products - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
DIN EN ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
DIN EN ISO 14630	Non-active surgical implants - General requirements

DIN EN ISO 14644-1	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness by particle concentration
DIN EN ISO 14644-2	Cleanrooms and associated controlled environments - Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration
DIN EN ISO 14644-3	Cleanrooms and associated controlled environments - Part 3: Test methods
DIN EN ISO 14644-5	Cleanrooms and associated controlled environments - Part 5: Operations
DIN EN ISO 14644-8	Cleanrooms and associated controlled environments - Part 8: Classification of air cleanliness by chemical concentration
DIN EN ISO 14698-1	Cleanrooms and associated controlled environments - Biocontamination control - Part 1: General principles and methods
DIN EN ISO 14698-2	Cleanrooms and associated controlled environments - Biocontamination control - Part 2: Evaluation and interpretation of biocontamination data
DIN EN ISO 14971	Medical devices - Application of risk management to medical devices
DIN EN ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
ISO 22442-1	Medical devices utilizing animal tissues and their derivatives - Part 1: Application of risk management
ISO 22442-2	Medical devices utilizing animal tissues and their derivatives - Part 2: Controls on sourcing, collection and handling
DIN EN ISO 22442-3	Medical devices utilizing animal tissues and their derivatives - Part 3: Validation of the elimination and/or inactivation of viruses and transmissible spongiform encephalopathy (TSE) agents
ISO/TS 10993-19	Biological evaluation of medical devices - Part 19: Physico-chemical, morphological and topographical characterization of materials
ISO/TS 10993-20	Biological evaluation of medical devices - Part 20: Principles and methods for immunotoxicology testing of medical devices
ISO 13779-3	Implants for surgery - Hydroxyapatite - Part 3: Chemical analysis and characterization of crystallinity ratio and phase purity
ISO 20417	Medical devices - Information to be supplied by the manufacturer
ASTM D4169	Standard Practice for Performance Testing of Shipping Containers and Systems
ASTM D7386	Standard Practice for Performance Testing of Packages for Single Parcel Delivery Systems
ASTM F1929	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM F1980	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
ASTM F 2347	Standard Guide for Characterization and Testing of Hyaluronan as Starting Materials Intended for Use in Biomedical and Tissue Engineered Medical Product Application
ASTM F88/F88M	Standard Test Method for Seal Strength of Flexible Barrier Materials

Specification device:

Specification	Acceptance criteria
---------------	---------------------

Appearance and content of packaging	Ok, complete
Appearance product	White to yellowish or reddish granules in a white / almost white powder or granules
Grain size hydroxy apatite ceramic	0.5-1.0 mm 1.0-2.0 mm
Quantity hydroxy apatite ceramic	0.5 ± 0.2 ml 1.0 ± 0.1 ml
Quantity sodium hyaluronate	118.8 ± 18.8 mg 237.5 ± 37.5 mg
Molecular weight of sodium hyaluronate	1.90–2.10 MDa
Concentration of sodium hyaluronate	95.0-105.0 %
Intrinsic viscosity of sodium hyaluronate	2.430-2.600 m ³ /kg
Consistency, texture (after hydration)	Connected mass with malleable texture
Volumetric porosity	45%-85%
Pore width	100-1500 µm
Compressive strength	> 1.5 MPa
pH value	<11
pH value (in TRIS buffer)	Around 7.4
Sterility	Sterile
Endotoxin content	<20 I.U. / product
Bioburden (prior sterilisation)	≤100 CFU / product

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