AGBY

Notarius Publicus

MARKS DOMSTOLE

INSTRUCTIONS FOR USE

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Phone +45 48 21 64 00

Uncontrolled Copy

Document Title **IFU: Auxiliary components for dental prosthetics**
Document ID **DEV-02261**
Version **1.0**
Release Date **01-Dec-2021**
Lifecycle State **Approved**
Internal Legacy Number

Content Approval

The Electronic signature of Content Approval means that this document has been reviewed and approved for compliance with the end user's and owner's requirements.

Signed by	ELOSMEDTECH\henrik.andersen
Decision date	30-nov-2021 21:54

QA Authorisation

The Electronic Signature of QA Authorisation means that this document has been reviewed and approved for compliance with business quality standards and SOPs.

Signed by	ELOSMEDTECH\azfar.rehman
Decision date	01-dec-2021 07:51

Table of contents

1. Medical device name.....	3
2. Information about the developer, manufacturer and manufacturing places of a medical device.....	3
3. The use of the medical device intended by the manufacturer.....	4
Intended use.....	4
Indications for Use.....	4
Contraindications.....	4
Scope.....	5
Medical device application method:.....	5
4. The main parameters and characteristics of the medical device.....	7
Compatible medical devices.....	10
Scope of supply of a medical device.....	10
5. Safety requirements.....	11
6. Sterilization.....	12
7. List of international regulatory documents.....	13
Scope.....	13
8. Operation conditions.....	13
9. Data on labeling and packaging.....	13
Data on labelling.....	13
Symbols on the packaging of a medical device:.....	14
Design of medical device labeling in Russian.....	15
Data on packaging.....	15
10. Storage and transportation conditions.....	18
11. Risk Analysis.....	18
12. Warranty liabilities.....	18
13. Expiration date.....	18
14. The procedure for the disposal or destruction of a medical device.....	19
15. Claims.....	19

1. Medical device name.

Auxiliary components for dental prosthetics:

1. Elos Accurate Scan Bodies:

Elos Accurate® Scan Body Nobel Replace NP
Elos Accurate® Scan Body Nobel Replace RP
Elos Accurate® Scan Body Nobel Replace WP
Elos Accurate® Scan Body Nobel CC NP
Elos Accurate® Scan Body Nobel CC RP
Elos Accurate® Scan Body Nobel CC WP
Elos Accurate® Scan Body Nobel CC NP (Single Abutment)
Elos Accurate® Scan Body Nobel CC RP (Single Abutment)
Elos Accurate® Scan Body Nobel CC WP (Single Abutment)
Elos Accurate® Scan Body Nobel CC 3.0 (Single Abutment)
Elos Accurate® Scan Body Nobel Multi-unit NP/RP
Elos Accurate® Scan Body Nobel Multi-unit WP
Elos Accurate® Scan Body Brånemark NP
Elos Accurate® Scan Body Brånemark RP
Elos Accurate® Scan Body Brånemark WP

2. Accurate Analog for Printed Models

Elos Accurate® Analog for Printed Models NobelReplace® NP
Elos Accurate® Analog for Printed Models NobelReplace® RP
Elos Accurate® Analog for Printed Models NobelReplace® WP
Elos Accurate® Analog for Printed Models Nobel Biocare® CC NP
Elos Accurate® Analog for Printed Models Nobel Biocare® CC RP
Elos Accurate® Analog for Printed Models Nobel Biocare® CC WP
Elos Accurate® Analog for Printed Models Nobel Biocare® Multi-unit NP/RP
Elos Accurate® Analog for Printed Models Nobel Biocare® Multi-unit WP
Elos Accurate® Analog for Printed Models Brånemark NP
Elos Accurate® Analog for Printed Models Brånemark RP
Elos Accurate® Analog for Printed Models Brånemark WP

Accessories:

Elos Accurate® Scan Body Driver Screwdriver
Elos Accurate® Scan Body Driver/Long Screwdriver
Elos Accurate® Analog Insertion Pin
Elos Accurate® Analog Insertion Screw
Elos Accurate® Analog Pliers

2. Information about the developer, manufacturer and manufacturing places of a medical device.

Developer:

Elos Medtech Pinol A/S, Engvej 33, 3330 Gørlose, Denmark.

Manufacturer:

Elos Medtech Pinol A/S, Engvej 33, 3330 Gørlose, Denmark.

Manufacturing site:

Elos Medtech Pinol A/S, Engvej 33, 3330 Gørlose, Denmark.

The manufacturer's authorized representative in the Russian Federation:

Nobel Biocare Russia Limited Liability Company, Moscow, Stanislavsky St., 21, bldg. 2, 109004.

3. The use of the medical device intended by the manufacturer.

Intended use

The Scan Body is intended to be used together with a dental scanner to locate the exact position of an implant/analog or abutment in relation to prosthetic restorations. The scan body is intended to be mounted directly on the implant/analog or abutment using the Scan Body Driver. The Scan Body is intended to be used directly in a patient's mouth when using an Intra oral scanner, or on a model of the patients' teeth setup when using a desktop scanner.

The Elos Accurate Analog for Printed Models is a dental implant replica intended to be mounted in a dental laboratory working model in order to duplicate the location and restorative platform orientation of the final dental implant.

Indications for Use

The Elos Accurate Scan Body (Q) is used in conjunction with a dental scanning device to determine the exact position of an implant, its analog or abutment during prosthetic restoration procedures. The Elos Accurate Scan Body is inserted directly onto the implant, implant analog or abutment using the Scan Body Driver screwdriver. The Elos Accurate Scan Body is intended for use directly in the patient's oral cavity in conjunction with an intraoral scanner, as well as on the patient's teeth set-up model when using a desktop scanner.

The Elos Accurate Analog for Printed Models is a dental implant replica intended to be mounted in a dental laboratory working model in order to duplicate the location and restorative platform orientation of the final dental implant.

- The Analog is intended for dental laboratory use only.
- The Analog is compatible with dental laboratory working models manufactured by additive 3D printing only.

Elos Accurate Analog Pliers & Elos Accurate Analog Insertion Pin

The Elos Accurate Analog Pliers (Pliers) and the Elos Accurate Analog Insertion Pin (Pin) shall be used to fix the Implant Analog in the working model in a dental laboratory environment.

Elos Accurate Analog Screw

In case of a loose fit due to inaccurate printing the optional Elos Accurate Analog Screw (Screw) can be mounted in the Analog by using the screwdriver interface at the end of the Pliers. This will help secure firm positioning of the Implant Analog in the working model in a dental laboratory environment.

Elos Accurate Driver Screwdriver

The special screwdriver for the Elos Accurate Scan Body Driver used to secure the fixed screw is made of biocompatible stainless steel.

The Pliers, Pin and Screw are compatible to all the variants of the Implant Analogs.

Contraindications

None.

Warnings and precautions

- The Elos Accurate Scan Body shall only be used with original branded implants, their analogs or abutments shown on the product label.

- When installing the device, make sure the Elos Accurate Scan Body fitting surface is clean and free of sediment, liquids, bone or soft tissue.
- The Elos Accurate Scan Body is small and must be handled carefully to prevent swallowing by the patient or aspiration.
- Tightening the screw shall be done manually using a screwdriver for the Elos Accurate Scan Body Driver, using a slight torque (no more than 5 Ncm). Too strong tightening torque may damage the device and cause inaccurate scan results.
- DO NOT use damaged or modified Elos Accurate Scan Bodies. Using damaged or modified Elos Accurate Scan Bodies may result in inaccurate scan results.
- Do not attempt to disassemble the Elos Accurate Scan Body.
- Do not attempt to rotate the Elos Accurate Scan Body while attached to the implant.
- The Elos Accurate Scan Body shall only be used with the compatible digital library file installed in the design software.
- The Elos Accurate Scan Body shall only be used by dental professionals.
- Allergic reactions to materials and the components they contain can occur in very rare cases.
- Inspect Scan Body/Scan Body Driver under normal lightning for the removal of visible soil and for damage or wear: Any Scan Body/Scan Body Driver Showing visible damage or wear must be disposed of according to section 14.
- The Analog shall only be used with a compatible digital library file installed in the dental design software.
- The Analog is for single-use only.
- Reuse of the Analog can result in loss of functionality.
- The Analog, Pin, Pliers or Screw must never be changed or modified.
- The Analog, Pin, Pliers and Screw must only be used and handled by dental professionals.
- The torque applied to a Screw for tightening a component (e.g. abutment) placed into the Analog must not exceed the recommended torque defined by the original component manufacturer.
- The torque applied to Screw must not exceed 5 Ncm.

Scope

Dentistry, used in the dental laboratory to reproduce the location and orientation of the final dental implant restorative platform; the device is not intended for intraoral use.

Potential risk class of the medical device use is 1.

Medical device application method:

Scanned Analog Procedure

Scanning


1. Clean and sterilize the Elos Accurate Scan Body according to the procedure in Section 7. Sterilization is only necessary for intraoral scanning.
2. Remove the Cover Screw, Healing Abutment, or similar device from the implant or abutment on which the Elos Accurate Scan Body is to be placed.
3. Place the Elos Accurate Scan Body on the implant / abutment in the patient's oral cavity or on its analog if a model is being scanned.
4. Tighten the attached screw manually using the Scan Body Driver Screwdriver applying slight torque.
5. Perform the scanning according to the scanner manufacturer's instructions to obtain a 3D digital model of the patient's oral cavity.
6. After scanning, the Elos Accurate Scan Body must be removed using the Scan Body Driver Screwdriver.

Design


1. Open the resulting 3D digital model in the design software provided with the scanner.
2. Select the corresponding digital library based on the prosthetic solution used in this case.
3. In the library, select the implant platform that matches the scanned implant, analog, or abutment.
4. Design the prosthetic restoration according to the instructions of the design software. Please note that for each of the prosthetic solutions, there may be restrictions on the design possibilities according to the maximum value of angulation and minimum value of wall thickness.

Implant Analog Procedure

1. Place the Insertion Pin in the dental laboratory working model from underneath.


2. Screw the Analog onto the Insertion Pin by hand (max. 5Ncm)


3. Pull down and rotate clockwise the Insertion Pin until the Analog is rotational fixated in the dental laboratory working model (one position only).

4. Slide the groove of the Pliers onto the Insertion Pin until the Pliers rest on the Insertion Pin and in the dental laboratory working model. Rotate the Pliers 90 degrees according to the Insertion Pin.

5. Apply compression force to the handle of the Pliers until the Analog is pulled into position.

6. The Screw can be mounted in the Analog by using the screwdriver interface at the end of the Pliers.

4. The main parameters and characteristics of the medical device

Elos Accurate® Scan Body Nobel Replace NP:

Total length 15.4 mm
Effective length: 10 mm
Diameter: 3.5 mm
Platform type: NP

Elos Accurate® Scan Body Nobel Replace RP:

Total length 17.4 mm
Effective length: 10 mm
Diameter: 4.3 mm
Platform type: RP

Elos Accurate® Scan Body Nobel Replace WP:

Total length 17.6 mm
Effective length: 10 mm

Diameter: 5 mm
Platform type: WP

Elos Accurate® Scan Body Nobel CC NP:

Total length 14.21 mm
Effective length: 10 mm
Diameter: 3.5 mm
Platform type: NP

Elos Accurate® Scan Body Nobel CC RP:

Total length 15.2 mm
Effective length: 10 mm
Diameter: 3.9 mm
Platform type: RP

Elos Accurate® Scan Body Nobel CC WP:

Total length 15 mm
Effective length: 10 mm
Diameter: 5 mm
Platform type: WP

Elos Accurate® Scan Body Nobel CC NP (Single Abutment):

Total length 12.7 mm
Effective length: 10 mm
Diameter: 3.5 mm
Platform type: NP

Elos Accurate® Scan Body Nobel CC RP (Single Abutment):

Total length 15.9 mm
Effective length: 10 mm
Diameter: 3.9 mm
Platform type: RP

Elos Accurate® Scan Body Nobel CC WP (Single Abutment):

Total length 16 mm
Effective length: 10 mm
Diameter: 5 mm
Platform type: WP

Elos Accurate® Scan Body Nobel CC 3.0 (Single Abutment):

Total length 15.7 mm
Effective length: 10 mm
Diameter: 3 mm
Platform type: 3.0

Elos Accurate® Scan Body Nobel Multi-unit NP/RP:

Total length 10 mm
Effective length: 10 mm
Diameter: 4.8 mm
Platform type: NP / RP

Elos Accurate® Scan Body Nobel Multi-unit WP:

Total length 10 mm
Effective length: 10 mm
Diameter: 6 mm
Platform type: WP

Elos Accurate® Scan Body Brånemark NP:

Total length 13.7 mm
Effective length: 10 mm
Diameter: 3.5 mm
Platform type: NP

Elos Accurate® Scan Body Brånemark RP:

Total length 12.7 mm
Effective length: 10 mm
Diameter: 4.1 mm
Platform type: RP

Elos Accurate® Scan Body Brånemark WP:

Total length 13.6 mm
Effective length: 10 mm
Diameter: 5 mm
Platform type: WP

2. Accurate Analog for Printed Models

Elos Accurate® Analog for Printed Models NobelReplace® NP:

Total length 6.8 mm
Platform type: NP

Elos Accurate® Analog for Printed Models NobelReplace® RP:

Total length 8 mm
Platform type: RP

Elos Accurate® Analog for Printed Models NobelReplace® WP:

Total length 8.3 mm
Platform type: WP

Elos Accurate® Analog for Printed Models Nobel Biocare® CC NP:

Total length 6.5 mm
Diameter: 2.9 mm
Platform type: NP

Elos Accurate® Analog for Printed Models Nobel Biocare® CC RP:

Total length 6.5 mm
Diameter: 3.35 mm
Platform type: RP

Elos Accurate® Analog for Printed Models Nobel Biocare® CC WP:

Total length 6.5 mm
Diameter: 4.4 mm
Platform type: WP

Elos Accurate® Analog for Printed Models Nobel Biocare® Multi-unit NP/RP:

Diameter: 1.73 mm

Platform type: NP / RP

Elos Accurate® Analog for Printed Models Nobel Biocare® Multi-unit WP:

Total length 8.7 mm

Diameter: 2.8 mm

Platform type: WP

Elos Accurate® Analog for Printed Models Brånemark NP:

Total length 7.3 mm

Diameter: 1.7 mm

Platform type: NP

Elos Accurate® Analog for Printed Models Brånemark RP:

Total length 7.2 mm

Diameter: 2.1 mm

Platform type: RP

Elos Accurate® Analog for Printed Models Brånemark WP:

Total length 7.0 mm

Diameter: 2.6 mm

Platform type: WP

Accessories:

Elos Accurate® Scan Body Driver Screwdriver:

Length: 12.5 mm

Elos Accurate® Scan Body Driver/Long Screwdriver

Length 80 mm

Elos Accurate® Analog Insertion Pin:

Total length 36.5 mm

Diameter m: 1.55 mm

Diameter b: 5 mm

Elos Accurate® Analog Insertion Screw:

Diameter 2.4 mm

Elos Accurate® Analog Pliers:

Length 145 mm

7 mm width

Compatible medical devices

Medical device "Auxiliary components for dental prosthetics" are compatible with implants of various platforms and sizes.

Scope of supply of a medical device

Elos Accurate® Scan Body Nobel Replace NP – 1 pc.

Elos Accurate® Scan Body Nobel Replace RP– 1 pc.

Elos Accurate® Scan Body Nobel Replace WP– 1 pc.

Elos Accurate® Scan Body Nobel CC NP– 1 pc.

Elos Accurate® Scan Body Nobel CC RP– 1 pc.
Elos Accurate® Scan Body Nobel CC WP– 1 pc.
Elos Accurate® Scan Body Nobel CC NP (Single Abutment) – 1 pc.
Elos Accurate® Scan Body Nobel CC RP (Single Abutment) – 1 pc.
Elos Accurate® Scan Body Nobel CC WP (Single Abutment) – 1 pc.
Elos Accurate® Scan Body Nobel CC 3.0 (Single Abutment) – 1 pc.
Elos Accurate® Scan Body Nobel Multi-unit NP/RP– 1 pc.
Elos Accurate® Scan Body Nobel Multi-unit WP– 1 pc.
Elos Accurate® Scan Body Brånemark NP– 1 pc.
Elos Accurate® Scan Body Brånemark RP– 1 pc.
Elos Accurate® Scan Body Brånemark WP– 1 pc.

2. Accurate Analog for Printed Models

Elos Accurate® Analog for Printed Models NobelReplace® NP– 1 pc.
Elos Accurate® Analog for Printed Models NobelReplace® RP– 1 pc.
Elos Accurate® Analog for Printed Models NobelReplace® WP– 1 pc.
Elos Accurate® Analog for Printed Models Nobel Biocare® CC NP– 1 pc.
Elos Accurate® Analog for Printed Models Nobel Biocare® CC RP– 1 pc.
Elos Accurate® Analog for Printed Models Nobel Biocare® CC WP– 1 pc.
Elos Accurate® Analog for Printed Models Nobel Biocare® Multi-unit NP/RP– 1 pc.
Elos Accurate® Analog for Printed Models Nobel Biocare® Multi-unit WP– 1 pc.
Elos Accurate® Analog for Printed Models Brånemark NP– 1 pc.
Elos Accurate® Analog for Printed Models Brånemark RP– 1 pc.
Elos Accurate® Analog for Printed Models Brånemark WP– 1 pc.

Accessories:

Elos Accurate® Scan Body Driver Screwdriver - 1 pc.
Elos Accurate® Scan Body Driver/Long Screwdriver - 1 pc.
Elos Accurate® Analog Insertion Pin– 1 pc.
Elos Accurate® Analog Insertion Screw– 1 pc.
Elos Accurate® Analog Pliers– 1 pc.

Materials used in the instruments

1. The Elos Accurate Scan Body - Titanium alloy and polyetheretherketone - short-term contact with skin and oral mucosa
2. Elos Accurate Analog for Printed Models - Titanium Grade 5– no contact, laboratory device.
3. Elos Accurate Driver Screwdriver - Stainless Steel - no contact, laboratory device.
4. Elos Accurate Analog Insertion Pin - Stainless Steel - no contact, laboratory device.
5. Elos Accurate Analog Insertion Screw - Stainless Steel - no contact, laboratory device.
6. Elos Accurate Analog Pliers - Titanium Grade 5– no contact, laboratory device.

The medical device does not contain and is not a source of medicines and pharmaceutical substances.

The medical device does not contain and is not manufactured with the addition of any derivatives of blood, tissues and / or materials of human and / or animal origin.

5. Safety requirements

Auxiliary components for dental prosthetics can only be used by dentists and clinicians who are familiar with dental surgery in accordance with the indications for its use and the general rules of dental / surgical procedure, as well as in compliance with occupational health and safety regulations. Before starting a procedure, please make sure all the required components, instruments and

auxiliary materials are available in a sufficient quantity and fully operational. All system components for intraoral use should be properly isolated to prevent accidental aspiration or swallowing. Therefore we recommend consulting an experienced specialist regarding the use of these devices. If there is any doubt regarding an indication or a method of system use, the device shall not be used until all subtleties have been clarified.

6. Cleaning and Sterilization Information

Cleaning

- Prepare a solution of a mild alkaline detergent 0.5%, temperature 15-20 °C.
- Completely submerge the instruments in the cleaning solution for five (5) minutes. Whilst submerged, brush the instrument with attention to the critical and hard-to reach areas with a soft-bristled brush for a minimum of 30 seconds and flush hollow instruments using a syringe with 0.5 mL of cleaning solution.
- Rinse one (1) minute under cold running tap water (15-20°C)
- Place the devices into the washing racks without overloading them. Then place the racks on the wash supports of the washer-disinfector. Avoid any contact between the devices, which may cause damage during washing
- Using a valid washer disinfector (according to ISO 15883 requirements) and an alkaline enzymatic detergent use the minimal cycle parameters set points mentioned in the following table.

Table 1 Procedure for cleaning

Cycle	Minimum Time	Minimum Temperature	Type of Detergent/Water
Pre-Cleaning	2 minutes	Cold (< 45°C)	Tap water
Wash	10 minutes	Heated (50-60°C)	Alkaline enzymatic detergent (0.2-1%)
Rinse	2 minutes	Cold (< 45°C)	Critical water ¹
Thermal Rinse	5 minutes	Heated (90°C)	Critical water ¹
Dry	25 minutes	Heated (90°C)	Not applicable

- At the end of the program remove the instruments.
- If not completely dry, dry instruments with single use, non-linting wipe

Sterilization

- Use a steam sterilizer confirming to EN 13060, DIN 285, and /or ST79
- Use the sterilization parameters provided in the table below, which have been validated

Table 2 Dynamic-air-removal Steam Sterilization Cycle, Exposure Time, Temperature & Drying Time

Procedure	Dynamic-air-removal-sterilization Cycle	
Exposure time	4 min.	3 min.
Temperature	132°C (270°F)	134°C (273°F)

¹ As per AAMI TIR34, water extensively treated usually by a multistep treatment process that may include a carbon bed, softening, DI, and RO or distillation, to ensure that the microorganisms and the inorganic and organic material are removed from the water.

Drying time	20 min.
--------------------	---------

Warning. Since the Elos Accurate Scan Body is a highly accurate instrument, its use immediately after sterilization in an autoclave is not permissible, as thermal expansion can lead to incorrect scan results. Cool the Elos Accurate Scan Body to 18-37 °C prior to use.

See Appendix # 1 to the Technical Documentation for details.

7. List of international regulatory documents.

Elos Medtech Pinol A/S implements and maintains a Quality Management System according to ISO 13485:2016 and MDD 93/42/EEC, Annex II, establishing limits for usability of Elos Medtech Pinol A/S products.

Document number	Issue date	Title
Directive 93/42/EEC	1993	Council Directive 93/42/EEC of 14 June 1993 on medical devices
MEDDEV 2.7.1 Rev. 4	2016	Clinical Assessment - A guide for manufacturers and notified bodies
ISO 13485:2016	2016	Medical devices. Quality management systems. Requirements for regulatory purposes
ISO 15223-1:2016	2016	Medical devices. Symbols to be used with medical device labels, labeling and information to be supplied — Part 1. General requirements
EN 1041:2008 + A1:2013	2013	Information supplied by the manufacturer of medical devices
IEC 62366	2015	Medical devices. Application of usability engineering to medical devices (IEC 62366:2007)
ISO 14971	2019	Medical devices — Application of risk management to medical devices
ISO 10993-1	2018	Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing within a risk management process
ISO 10993-5	2009	Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity
ISO 10993-12	2012	Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference materials
ISO 10993-18	2020	Medical devices. Biological evaluation of medical devices. Part 18. Chemical characterization of materials
ISO 11737-1:2018	2018	Sterilization of medical devices. Microbiological methods. Part 1. Determination of a population of microorganisms on products
ISO 11737-2:2019	2019	Sterilization of medical devices. Microbiological methods. Part 2: Tests of sterility performed in the validation of a sterilization process
EN/ISO 11607-1:2019	2019	Packaging for terminally sterilized medical devices. Part 1. Requirements for materials, sterile barrier systems and packaging systems
EN/ISO 11607-2:2019	2019	Packaging for terminally sterilized medical devices. Part 2: Validation requirements for forming, sealing and assembly processes
EN/ISO 14937	2009	Sterilization of health care products. General requirements for the characterization of a sterilizing agent and for the development, validation and routine control of a sterilization process for medical devices

EN 1639	2009	Dentistry. Medical devices for dentistry. Instruments
ISO 14644-1	2015	Cleanrooms and associated controlled environments – Part 1. Classification of air cleanliness by particle concentration
ISO 14698-1	2003	Cleanrooms and associated controlled environments – Biocontamination control. Part 1. General principles and methods
ASTM F899	2019	Standard Specification for Forged Stainless Steel for Surgical Instruments
ASTM D4169	2016	Standard Practice for Conducting Shipping Tests of Containers and Systems

Scope

Design, manufacture, and sales of general non-active, non-implantable medical devices, non-active implants, dental devices and accessories and components to In Vitro Diagnostic Medical Devices (IVD).

8. Operation conditions.


During operation, Elos Accurate Scan Bodies can withstand temperatures ranging from +32 °C to +42 °C and exposure to biological fluids and effusions of body tissues they come in contact with during operation.

The instruments can be used in public or private dental clinics and only by personnel who have been trained and have sufficient experience in this area.

9. Data on labeling and packaging.

Data on labelling

The following original labeling must be applied to the packaging of a medical device: medical device name, catalog number, lot number, expiration date, warning symbols, CE marking, QR code, Manufacturer information (name, address and website of the company).

Elos Accurate Analog for Printed Model

Elos Accurate®
ELOS

Analog for Printed Models

Nobel Biocare® Multi-unit™ NP/RP

REF

PMA-MUA45-1

LOT

123456

[EN] Implant analog [BG] Аналог на имплант [CS] Analog implantátu [DA] Implantatreplik [ET] Implanntaadi koopia [FI] Implanttiansalogi [FR] Analogie d'implant [DE] Implantatgegenstück [EL] Ανάλογο εμφυτεύματος [ES] Análogo de implante [HR] Analog implantata [HU] Analóg implantátum [IS] Plantafléstaða [IT] Analogo di impianto [LV] Implanā analogs [LT] Implantų analogas [NL] Implanntaat-analoog [NO] Implantatanalog [PL] Analog implant [PT] Análogo do implante [RO] Copie implant [SE] Implantat analog [SK] Analóg implantátu [SL] Analog vsadka [SR] Analog implanta [TR] Implant analog

www.elosmedtech.com/IFU

Elos Medtech Pinol A/S, Engvej 33
3330 Gørlose, Made in Denmark

2019-01-01

MD

UDI

(01)05712821005847
(11)190101
(10)123456

C29721_01

Symbols on the packaging of a medical device:

REF	Part number
LOT	Lot
	Refer to the instructions for use
	Date of manufacture
	Non-sterile
	Do not reuse!
	Attention! This device can only be sold to a clinician or supplied upon the order of a clinician!
	Do not use if package is damaged
	CE Marking

	Manufacturer
	Medical device
	Unique numeric or alphanumeric code
	QR code
	Avoid direct sunlight
	Keep dry
Made in	-

Design of medical device labeling in Russian

The Russian-language labeling must contain the following information:

- name according to the Registration Certificate;
- the number of the Registration Certificate will be indicated on the labeling after the state registration of the medical device and the receipt of the corresponding Registration Certificate;
- The date of the state registration of the medical device will be indicated on the labeling after the state registration of the medical device and the receipt of the corresponding Registration Certificate;
- data of the authorized representative on the territory of the Russian Federation;

Part No.: _____ Lot: _____

Elos Accurate® Scan Body Nobel CC NP Single Abutment

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Registration Certificate: No. RZN ____/____/____ dated ____/____/____

Auxiliary components for dental prosthetics

Authorized representative, seller: Nobel Biocare Russia LLC 109004 21 Stanislavskogo Street, Building 2, 109004 Moscow Phone +7 (495) 974-77-55.

Do not use if package is damaged Refer to the instructions for use

Non- Sterile Date of manufacture, see on the package

Do not reuse

FOR PROFESSIONAL DENTISTRY USE ONLY

Data on packaging

The Elos Accurate Scan Bodies are supplied non-sterile in a blister packaging. The Elos Accurate Scan Body is designed for dental clinic sterilization prior to intraoral use.

Primary packaging

The primary packaging is being in the direct contact with the product (device) and consists of PET-G (Polyethylene Terephthalate Glycol) blister tray with the paper/thin cardboard) onlay/label performing the function of the barrier.

The blister tray and the paper onlay are binded together with the method of the heat sealing.

Materials of the primary packaging

Blister tray	PET-G grade 6763, manufacturer PRENT THERMOFORMING, Rørvangsvej 37, DK-4300, Holbæk
Paper/thin cardboard onlay	Medical grade paper 330 g/m ² Medical paper grade R114 29GP BL, manufacturer Oliver™ Healthcare Packaging, Spurkt 1, 5804 AR Venray, Niðurlond

Primary packaging

Analog Insertion Pin, Analog Insertion Screw, Screwdriver are also supplied non-sterile in a blister packaging.

Materials of the primary packaging

Blister tray	PET-G grade 6763, manufacturer PRENT THERMOFORMING, Rørvangsvej 37, DK-4300, Holbæk
Paper/thin cardboard onlay	Medical grade paper 330 g/m ² Medical paper grade R114 29GP BL, manufacturer Oliver™ Healthcare Packaging, Spurkt 1, 5804 AR Venray, Niðurlond

Primary packaging

Analog for Printed models, Pliers and a Long Screwdriver are supplied non-sterile and packaged in a welded bag according to specification

Materials of the primary packaging

Welded bag (66x162mm)	LD/LLDPE – Linear low-density polyethylene Manufacturer, A/S Christian Berner Maglebjergvej 5A DK-2800 Kgs. Lyngby
-----------------------	---

Transportation packaging

During transportation, blisters are placed in group shipping boxes made of white, durable corrugated cardboard with a capacity of up to 50 blisters. This delivery of boxes avoids damage to boxes and devices during transportation and storage.

The density of the cardboard shipping box is 230 g/m².

The labels are printed on the outer box and on each of the blisters.

Design of the transportation packaging labelling

The typical transportation packaging labelling is indicated on the image below

Ship-to Country: RU
 Parent Reference: CON14082
 Delivery note number: 891413536
 Weight: 0.62 KG
 Handling unit number: 105438920

Information concerning the transportation packaging labelling in Russian is indicated below:

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Definitions of the margins of the transportation packaging labelling:

Ship to Country is stated individually


Patient Reference – number of the main shipment which might include several orders. The number is generated in the system SiS (ShipItSmarter).


Delievery note number is the number of the specific order. The number is generated in the system SiS (ShipItSmarter). The Delivery note is applied to each shipment. The Delivery note contains information regarding: the manufacturer, cargo shipper, cargo receiver, the list and names of the medical devices.

Weight of the shipment is indicated individually.

Handling unit number is a physical unit that includes the packaging (cargo packaging/wrapping material) and goods for shipping that are included in the packaging. The handling unit always represents the combination of products and packaging. The number is generated in SAP.

Symbols on the transportation packaging

10. Storage and transportation conditions.

Devices shall be stored and transported in their original packaging at room temperature in a dry place, protected from direct sunlight, at a room temperature of +18 to + 25 °C. Failure to comply with the storage or transportation conditions may affect the properties of the component and result in damage to it.

Storage, transportation and handling of devices must be carried out in their original packaging. Transportation is carried out by all means of transport in accordance with local legislation.

Transportation conditions: from +10 °C to +40 °C (50 °F to 104 °F), and the recommended relative humidity range is 0% to 80%.

11. Risk Analysis.

Elos implements risk management in compliance with EN ISO 14971.

The Risk Management Report is an overview of all types of risk management activities in the following aspects:

- the risk management plan has been properly implemented;
- a brief overview of the design failure mode and effect analysis (DFMEA) and process failure mode and effect analysis (PFMEA);
- the overall residual risk is acceptable;
- appropriate methods are in place to collect relevant production and post-production information;

The analysis of the risk management process was carried out, and the implementation of the measures provided for in the risk management plan was ensured. Risks associated with design, application (operation) and manufacturing process have been identified, analyzed and assessed. In addition, appropriate mitigation measures have been identified and implemented, their effectiveness has been verified or validated, and residual risks have been assessed.

Appropriate techniques have been developed and applied to obtain relevant post-production information.

12. Warranty liabilities.

If the packaging of the instrument is not opened, has no mechanical damage or other defects, the products can be accepted for return or exchange.

Return or exchange is subject to the date the instrument was purchased by the customer. Within 30 days from the date of purchase, a refund can be made, in terms of a later date of purchase - the instrument can be exchanged.

If the packaging of the device was opened, but a defect was detected, the instrument can be replaced.

Regional distributor in Russia:

Nobel Biocare Russia Limited Liability Company

Address: 21 Stanislavsky Street, bldg. 2, Moscow 109004

Tel.: +7 495 974 77 55, fax: +7 495 974 77 66

13. Expiration date

Expiry date N/P

14. The procedure for the disposal or destruction of a medical device.

Medical devices that are potentially contaminated or unusable for further use must be disposed of as medical waste in accordance with local and federal legislation requirements.

The separate collection, recycling and disposal of packaging material must be carried out in accordance with applicable packaging and packaging waste legislation.

After use, medical devices are classified as medical waste and are class B medical waste - epidemiologically hazardous waste. Measures for the disinfection and disposal of used devices must be carried out in accordance with the requirements of the sanitary rules SanPiN 2.1.3684-21 Sanitary and Epidemiological Requirements for Medical Waste Management.

Safely discard potentially contaminated or no longer usable medical devices as healthcare (clinical) waste in accordance with local healthcare guidelines, country and government legislation or policy.

Separation, recycling or disposal of packaging material shall follow local country and government legislation on packaging and packaging waste, where applicable.

15. Claims

Contact information for contacting on issues related to the operation, maintenance, quality and compliance of medical devices with regulatory requirements on the territory of the Russian Federation:

Nobel Biocare Russia Limited Liability Company

Legal address: Moscow, st. Stanislavsky, 21, bldg. 2, 109004;

Mailing address: Moscow, st. Stanislavsky, 21, bldg. 3, 109004;

Contact phone: +7 (495) 974-77-55;

Email address: olga.kalach@nobelbiocare.com

Website: www.nobelbiocare.ru

This is to certify that **Mr Søren Olesen** who today in my presence at the Notarial Office approved and signed this document. No conspicuous corrections or addenda were found in the document.

The notary has not reviewed or approved the contents of the document and only confirms the identity, the accuracy of the dating and the authenticity of the signature.

Mr Søren Olesen has proved his identity by presenting Danish Passport No 212724997.

District Court of Lyngby, Denmark, on December 1, 2021

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