Bensheim, August 03rd 2023

Dentsply Sirona
Implants
Dentsply Implants
Manufacturing GmbH
Rodenbacher Chaussee 4
63457 Hanau, Germany

Dentsply Sirona

For Instructions for Use
and symbols glossary refer to
ifu.dentsplysirona.com

Instructions for Use

PrimeTaper EV

Dr. Aline Frare Mocruha
Regulatory Affairs Manager

Description of the System

The DS Implants system assortment comprises dental implants, cover screws, transmucosal abutments, other prosthetic components, as well as surgical and prosthetic instruments.

Intended Users

Users must be licensed dentists/physicians, be trained in implant dentistry, and should be familiar with the DS Implants system.

Clinical Benefit and Patient Target Group

The DS Implants assortment is designed for replacing one or more missing teeth in patients to restore chewing function. The patient target group are persons with adult jaws to whom none of the contraindications apply that are related to dental implantations.

Compatibility

The implant is compatible with the restorative assortment of DS system products and the Astra Tech Implant System with the EV connection. The components may only be used together with original products and instruments.

Purpose

The DS Implants assortment is designed for replacing one or more missing teeth in patients to restore chewing function.

Description of the Products

PrimeTaper EV

The PrimeTaper EV Implant is a screw-shaped dental implant with a defined surface achieved through grit blasting, followed by a process involving treatment in diluted hydrofluoric acid. There is a conical connection between the implants and the abutments for a tight and stable connection.

The PrimeTaper EV Implant has an inner geometry with a conical and indexing configuration.

Different diameters and lengths of the implants are available providing versatility according to the thickness and shape of the bone.

Implants with a diameter of 3.0 mm, are not indicated for the replacement of a central incisor, a canine, a premolar, or a molar in the maxilla as well as for replacement of a canine, a premolar, or a molar in the mandible.

PrimeTaper Drills

The PrimeTaper Drills are used during implant surgery according to the drilling protocol. The drills have laser markings showing the drilling depth, are marked with the respective diameters and have an RA-connection to fit into a contra angle.

PrimeTaper Taps

The PrimeTaper Taps are used for widening the osteotomy in very dense bone.

PrimeTaper Radiographic Templates

PrimeTaper Radiographic Templates are transparent plastic sheets with printed contours representing PrimeTaper EV implants in different magnifications. The magnifications range is from 1.0 to 1.5. The templates are used as planning aids.

Intended Use

PrimeTaper EV Implants

The implants are intended for both one- and two-stage surgical procedures in the following situations and with the following clinical protocols:

- Replacing missing teeth in single or multiple unit applications in the mandible or maxilla.
- Immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- Especially indicated for use in soft bone applications where implants with other implant surface treatments may be less effective.
- Immediate and early loading for all indications, except in single tooth situations on implant shorter than 8 mm or in soft bone (type 4)¹ where implant stability may be difficult to obtain and immediate loading may not be appropriate.
- PrimeTaper EV implants Ø 3.0 are only indicated for the replacement of lateral incisors in the maxilla and the replacement of incisors in the mandible.

PrimeTaper Drills and PrimeTaper Taps

The PrimeTaper Drills and PrimeTaper Taps are used for the preparation of an appropriate osteotomy for the implant and provide preparation variations to meet different bone conditions at the individual implant site.

PrimeTaper Radiographic Templates

The PrimeTaper Radiographic Templates are planning tools used prior to implant surgery. The templates are intended to be used together with the X-ray of the patient's jaw and helps the clinician in planning the optimal position, direction, and size of the implant.

PrimeTaper Direction Indicator

The PrimeTaper Direction Indicator is intended for measuring the distance between osteotomies, controlling the direction of the osteotomy and facilitates parallelization of osteotomies.

Indications

Partial and total edentulism

Contraindications

- Allergic or hypersensitivity reactions to materials used (see section "Material")
- Recent myocardial infarction (< 3 months)
- Recent cerebrovascular accident
- Recent cardiac-valvular prosthesis placement
- Hemorrhagic diathesis
- Immunosuppression
- Active treatment of malignancy
- Severe liver dysfunction
- Florid infection
- Drug abuse
- Psychiatric illness
- Intravenous bisphosphonate use
- Implants with a diameter of 3.0 mm, are not indicated for the replacement of a central incisor, a canine, a premolar, or a molar in the maxilla as well as for replacement of a canine, a premolar, or a molar in the mandible.

Warnings

Oral surgery and rehabilitation entail general risks for complications:

- Components that accidentally are dropped in the patient's mouth may be swallowed or aspirated which can result in suffocation or physical injury. Care must be taken to have control over small devices.
- Allergy or hypersensitivity to chemical ingredients of material used. Patient health history must be carefully assessed.
- Failure to plan and prepare implant site position and depth correctly can result in harm to nerves, blood vessels, teeth, jawbone, or maxillary sinus. It is strongly recommended that clinicians always undergo special training before starting up new treatment methods.
- Highly angled abutments on small diameter implants are not recommended for use in the posterior region.

Precautions

- Prior to each procedure it must be ensured that all necessary components, instruments and materials are available in the required quantities and are fully operational.
- Position the patient so that the danger of aspiration of components is minimized.
- All components that are used intraorally must be secured to prevent aspiration or swallowing.
- Improper patient selection/analysis or failure to comply with the step-by-step instructions can lead to failed osseointegration or loss of osseointegration.
- During preparation with drills the recommended speed must be used.
- Replace reusable drills and taps after a maximum of 10 reprocessing cycles/patients or if they are damaged or dull.
- Fracture of mandible may occur in a very atrophic mandible, during surgery or during routine oral function after surgery.
- Fracture or damage of implants, healing components, prosthetic components or prosthesis requiring removal and replacement of an implant may occur due to unfavorable loading conditions. Mastication forces and restoration design must therefore be evaluated carefully.
- When treatment planning for implants shorter than 8 mm in length consider the widest possible implant, a two-stage surgical approach and splinting of implants. Closely monitor the patient for peri-implant bone loss or change in the implants' response to percussion. If the implant shows greater than 50 % bone loss or mobility, consider possible removal of the implant.
- Routine implant treatment is not recommended until the end of the jaw-bone growth and for pregnant women.

CE/CE 0123

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.
Other languages are available upon request.

Not all products are available in all countries.

Dentsply Sirona

Status of information:
IFU 7002-EN_RU Rev. 002 F 2023-07

Instructions for Use

- Inadequate cooling, incorrect implant site preparation, incorrect instrument settings or excessive installation torque can result in failed implant.
- All products intended for single use must not be re-used. Failure to observe this precaution may prevent osseointegration of an implant or cause local or systemic infection.
- Bone quality and quantity, local infections, disturbed initial healing or premature loading are all factors, which may influence implant survival. Successful healing and osseointegration may be negatively affected by the following factors:

- Inflammations, tumors, abscess or cysts in the jaw area
- Insufficient height and/or width of bone
- Insufficient soft tissue coverage
- History of therapeutic radiation in the area
- Bone and wound healing disorders (e.g. uncontrolled diabetes)
- Autoimmune diseases
- Uncontrolled para-functional habits (e.g. bruxism)
- Inadequate oral hygiene
- Alcohol abuse
- Smoking

Adverse Reactions

Adverse reactions during or after implant treatment can be:

- Allergic or hypersensitivity reactions to materials used (see section "Material")
- Aspirated or swallowed components
- Permanent harm to nerve or tooth
- Harm to jawbone or maxillary sinus
- Necrosis due to inadequate cooling or excessive torque
- Fractured mandible
- Uncontrolled post-operative bleeding (e.g. damage to arteria lingualis)
- Postoperative discomfort such as:
 - Edema
 - Hematoma
 - Temporary local swelling
 - Temporary anesthesia or paresthesia and masticatory impediment
 - Pain
- Implant loss
- Peri-implantitis
- Local or systemic infection
- Bone loss as a result of peri-implantitis or mechanical overload
- Fracture or damage of implants, healing components, prosthetic components or prosthesis
- Removal and replacement of an implant
- Esthetic or functional complications

MRI Safety Information

All DS Implants products which remain in the patient's body have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the DS Implants in the MR environment is unknown. Scanning a patient who has such a product may result in patient injury.

Step-by-Step Instructions

Preoperative planning

Each patient must be carefully evaluated from a general medical and from a dental/oral point of view. Before any treatment procedure begins, all pathology in the jaws and remaining teeth must be resolved or be in the process of healing.

Clinical, radiographic as well as model analysis, are all essential prerequisites for the successful implant treatment. Pre-operative planning should be based on the expected restorative treatment outcome.

PrimeTaper Radiographic Templates can be used as a planning aid in conjunction with an X-ray. The X-ray is placed on a light box and the appropriate template is placed over the X-ray. The clinician can then select the most suitable implant position, direction and size.

NOTE:

Double lines under the images of implants mean the maximum additional drilling depth due to Drill Tap (1 mm).

Preparation of the implant site

Prepare the implant site according to chosen surgical procedure and implant size.

Removal of the sterile drills and taps from the blister

- Open the outer blister packaging.
- Pour out the inner blister onto a sterile area.
- Secure the drill or tap by squeezing the blister.
- Expose the shaft by bending the top of the blister.
- Engage the drill or tap with the contra angle.
- Follow the recommended drilling protocol below.

Recommended drilling protocol

NOTE:

- Drill at a maximum speed of 1500 rpm under profuse irrigation with physiological saline solution.
- Use the taps at a maximum speed of 25 rpm.
- The appropriate implant site preparation requires sharp drills and taps.
- The 2 mm mark on the drill indicates the height of the micro-threaded portion of the implant.

Recommended Protocol	Ø 1.9	Ø 2.35	Ø 2.85	Ø 3.55	Ø 4.15	Ø 4.75	Ø 5.35	Very dense	Very soft
Implant	1	2	3	4	5	6	7	Tap	
Ø 3.0	X S	1	2	3				+	1 3
Ø 3.6	S	1		3	4			+	1 2 4
Ø 4.2	S	1		3	4	5		+	1 3 5
Ø 4.8	L	1		3		5	6	+	1 3 4 6
Ø 5.4	L	1		3		5	6	7	1 3 6 7

- Figures 1 – 7 refer to drill numbers
- "]" refers to cortical preparation only (mandatory)
- + Very dense bone (type 1)¹ – screw tap available, used after cortical preparation
- Very soft bone (type 4)¹ – In case of very soft trabecular bone, as final full depth drill, use one size smaller drill than in the recommended protocol, and finalize with cortical preparation

Considerations for special clinical situations

Finetuning			Ø 2.65	Ø 3.25	Ø 3.85	Ø 4.45	Ø 5.05			
Recommended Protocol		Ø 1.9	Ø 2.35	Ø 2.85	Ø 3.55	Ø 4.15	Ø 4.75	Ø 5.35	Very dense	Very soft
Implant		1	2	3	4	5	6	7	Tap	
Ø 3.0	X S	1	2	3	4	5	6	7	+	1 3
Ø 3.6	S	1		3	4				+	1 2 4
Ø 4.2	M	1		3	4	5			+	1 3 5
Ø 4.8	L	1		3		5	6		+	1 3 4 6
Ø 5.4	L	1		3		5	6	7	+	1 3 5 7

5 additional drill diameters (Ø 2.65, Ø 3.25, Ø 3.85, Ø 4.45, Ø 5.05) are available for finetuning the osteotomy diameters compared to the recommended protocol, if required:

- a slight under-preparation
 - a slight widening of the osteotomy
- Finalize the osteotomy with cortical preparation]

Extraction sockets

Use drilling sequence for very soft bone (final full depth drill one size smaller than in the recommended protocol), or alternatively, one half size smaller than the recommended protocol

Narrow crest

A tap may be used in case of contact with lingual/buccal cortical bone during osteotomy preparation.

Short implants

For 6.5 mm implants consider using, as final full depth drill, a drill that is one size or half a size smaller than in the recommended protocol

In dense trabecular bone

Use a tap or a final full depth drill one half size larger than in the recommended protocol.

CE/CE 0123

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.
Other languages are available upon request.

Not all products are available in all countries.

Dentsply
Sirona

Status of information:
IFU 7002-EN_RU Rev. 002 F 2023-07

Instructions for Use

PrimeTaper Direction Indicator

Place the PrimeTaper Direction Indicator in the osteotomy after preparation, with drill number 1 and 3 respectively, to check direction and to facilitate parallelism between osteotomies.

Depth Gauge

Check the osteotomy depth using the Depth Gauge.

Implant Installation

1. Open the packaging including the blister.
2. Drop the sterile inner container onto a sterile area.
3. Remove the cap from the container using a twisting motion to expose the top of the implant.
4. Attach the appropriate Implant Driver EV to a contra-angle or a surgical driver handle.
5. When picking up the implant, make sure the implant driver is fully seated into the implant. Press downwards to activate carrying function before picking up the implant.
6. Install the implant with a contra-angle at low speed (25 rpm) under profuse irrigation with physiological saline solution.
7. Allow the implant to work its way into the osteotomy until the top of the micro threaded part is flush with the bone level.

NOTE: Avoid applying unnecessary pressure.

8. For the final seating, use the Torque Wrench EV with the Torque Wrench EV Surgical Driver Handle and the appropriate Implant Driver EV.

NOTE: Do not exceed 45 Ncm when installing the implant. If not completely seated before reaching 45 Ncm, reverse/remove the implant and widen the osteotomy appropriately. It is recommended to have a titanium forceps available in case the implant driver does not provide sufficient carrying function during the removal procedure.

Post-operative care

Instruct the patient regarding the need for regular oral hygiene. For the first 7 – 10 days post-op, oral hygiene can be supplemented by a suitable anti-microbial rinse. Mechanical loading of the implant site should be avoided in the post-operative period. A refined provisional restoration can be worn for esthetic purposes.

Follow-up

After treatment is completed, regular follow-up examinations are important to maintain health.

Delivery – Storage – Disposal

NOTE: Inspect sterile packed products visually for breaches of packaging integrity prior to use. Do not use the products if any breaches of packaging integrity are detected.

Delivery

The PrimeTaper EV Implants are delivered sterile and are intended for single use. The PrimeTaper EV Drills and Taps are delivered sterile. The products are sterilized by e-beam irradiation. The Depth Gauge is delivered non-sterile.

The Depth Gauge, the PrimeTaper Direction Indicator and the PrimeTaper Radiographic Templates are delivered non-sterile.

Packaging

PrimeTaper EV Implant

The PrimeTaper EV Implants package comprises of a cardboard outer box containing a plastic blister which is sealed by a foil. The foil of the blister serves as a sterile barrier and the content of the blister is sterile (single sterile-barrier system). Inside the blister there is an inner package containing the PrimeTaper EV Implant.

Do not use the PrimeTaper EV Implants after expiry date or if the blister is damaged. Do not use PrimeTaper EV Implants if the blister was unintentionally opened before use.

PrimeTaper Drills and Taps

The PrimeTaper Drills and Taps package comprises of an inner blister that serves as a holder for the drill or tap and an outer sealing blister. The foil of the outer blister serves as a sterile barrier and the content of the blister is sterile (single sterile-barrier system).

The drills and taps can be used beyond the expiry date or in case the blister was unintentionally opened if the cleaning and sterilization instructions are followed. However, do not use the drills and taps if the blister is damaged.

Storage

Store the products at room temperature in the original packaging under conditions normal for dental practices (from +18 to +25 °C and relative humidity

from 10 % to 90 %). PrimeTaper Drills should be stored at temperature from +10 °C to +45 °C and relative humidity from 10 % to 90 %.

Disposal

For disposal of the packaging and components comply with the currently applicable national waste disposal regulations in your country. Instruments and components that were used in the patient's mouth may be potentially infectious and must be disposed of safely in containers approved for this purpose. Dispose of sharps immediately after use in designated puncture-proof containers according to ISO 23907-1 to prevent infection and injury.

Note: Dispose of reusable drills and taps after a maximum of 10 reprocessing cycles/patients or if they are damaged or dull.

Material

Product	Material
PrimeTaper EV Implants	Titanium Grade 4 acc. to ASTM-F67
PrimeTaper Drills	Stainless Steel acc. to ASTM-F899 with TiN Coating
PrimeTaper Taps	Stainless Steel acc. to ASTM-F899
PrimeTaper Radiographic Templates	Polycarbonate Foil
PrimeTaper Direction Indicator	Titanium Alloy (Ti6Al4V-ELI) acc. to ASTM-F136

Cleaning and Sterilization

PrimeTaper EV Implant

Not applicable.

The product is sterile.

The product is for single use only.

PrimeTaper Drills and Taps

All drills and the taps can be used for a maximum of 10 cases.

PrimeTaper Radiographic Templates

The templates can be cleaned with common cleaning liquids.

Automated Procedure

IMPORTANT:

- Prior to clinical use, all components and instruments delivered non-sterile must be cleaned, disinfected and sterilized according to ISO 17664. Instruments that can be disassembled must be disassembled for cleaning and may be assembled again for sterilization only in a dry condition.
- Do not use ultrasonic cleaning for the Implant Driver EV.
- Drills must be kept separated from each other when using ultrasonic cleaning, e.g. keep the drills hanging in a holder during the ultrasonic cleaning.

Manual pre-cleaning

Coarse soiling must be removed immediately (within a maximum of two hours) after use with a soft nylon brush until visibly clean.

Ultrasonic cleaning

The action times and concentrations specified by the manufacturer of the cleaning agent must be observed.

Mechanical processing

Manual pre-cleaning is necessary to ensure efficient mechanical cleaning. Suitable washer disinfectors as well as cleaning and neutralization agents are to be used for mechanical cleaning. The instructions from the manufacturer of the washer disinfectant must be followed. Cleaning and neutralization agents are to be dosed and used in accordance with the manufacturer's instructions. A cleaning program with thermal disinfection ($A_0 \geq 3000$) is recommended. Either demineralized water or water which satisfies this level of purity is recommended for cleaning.

The mechanical processing method was validated with the following materials and devices (acc. to ISO 15883 and AAMI TIR 30):

- Washer disinfectant, type G 7836 CD, Miele & Cie. KG, Gütersloh
- Cleaning program VARIO TD
Pre-rinsing with cold water 4 min
Cleaning 10 min from 45–55 °C
Neutralization 6 min
Intermediate rinsing with cold water 3 min
Thermal disinfection 5 min at 90 °C ($A_0 \geq 3000$)
Drying 30 min at 105 °C
- Cleaning agent neodisher® MediClean forte (0.5 %)/neutralization agent neodisher® Z (0.1 %), Dr. Weigert, Hamburg.

Maintenance

Before sterilization, all products should be inspected with the naked eye for visible soil, impairments and/or corrosion. Particular attention should be paid to design features such as channels, blind holes, threads, undercuts and mating surfaces.


Check all markings on products for visibility and readability by inspection with the naked eye. Check matching components for proper assembly and functionality. Check products with movable parts for correct operation. Check

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.
Other languages are available upon request.

Not all products are available in all countries.

Status of information:
IFU 7002-EN_RU Rev. 002 F 2023-07

Instructions for Use

pre-mounted components for completeness and integrity.
Dispose of defective products.

Sterilization

Sterilization is to be performed corresponding to the following scheme:

- Preparation for sterilization:** Place components in a sterilization pouch which is legally marketed for use with the recommended sterilization parameters. Packaging must comply with the requirements according to EN ISO 11607, ANSI/AAMI ST79 and AAMI TIR12. Every sterilization package must have a sterilization indicator and sterilization date.
- Parameters (acc. to ISO 17665, EN 13060, EN 285 and AAMI TIR12)**

Method	Cycle	Temperature	Exposure Time*	Dry Time
Steam (EU**)	Dynamic air removal (prevacuum)	134 °C	3 min	20 min
Steam	Dynamic air removal (prevacuum)	132 °C	4 min	20 min
Steam	Dynamic air removal (prevacuum)	135 °C	3 min	20 min
Steam	Gravity displacement	121 °C	30 min	20 min

* Minimum exposure times; the operating times are longer and may vary depending on the device.
** EU customers are required to use these steam sterilization parameters.

- Storage:** Store the sterilized components dry and dust-free at room temperature.

Transportation Conditions

Products must be transported and carried with care in accordance with the special precaution measures applying to each type of a medical product, at temperatures ranging from -50 to +50 °C and a relevant humidity from 10 % to 95 %.

Operating Conditions

Dental practices and laboratories.

Can only be used by dentists / physicians who have been trained accordingly.

For medical products which remain in the mouth (for example, implants, abutments).

Products may be operated at temperatures from +10 to +45 °C.

Application Area

Stomatology. Dental implantology.

Shelf Life

Shelf life is set at five years in the original packaging.

Shelf life is indicated on the labeling.

Implant System Warranty

For full implant systems warranty information, refer to the detailed terms and conditions.

Product Availability

ATTENTION: Not all products described in this IFU may be regulatory cleared/released/licensed in all markets. See current product catalog and brochures for product assortment and representative product photos or contact the local Dentsply Sirona sales office for product availability.

Documentation

To ensure traceability, we recommend that the removable label on the blister is attached to the patient record.

Trademarks

All trademarks and company names are the property of their respective owners. To improve readability for our customers, DS Implants does not use ® or ™ in

body copy. However, Dentsply Sirona Implants does not waive any right to the trademark and nothing herein shall be interpreted to the contrary.

Authorized Representative on the Territory of the Russian Federation

Dentsply Sirona LLC

115035, Moscow, Ovchinnikovskaya embankment, 18/1, bld.2, premises 3N
Tel. +7 (495) 725 10 87

Valid for registration (Russia)

Dental implants PrimeTaper EV

Dental drills, sterile, multiple use, for dental implant systems, with accessories

Explanation of Symbols

Symbol

Title

Manufacturer and Date of manufacture

Use-by date

Batch Code

Catalog Number

Sterilized using irradiation

Do not re-sterilize

Do not re-use/ Single use/ Use only once

Do not use if package is damaged

Consult Instructions for Use

European conformity (CE) mark

Prescription only.
U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician.

Medical Device

Single sterile barrier system

CE/CE 0123

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hensu, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.
Other languages are available upon request.

Not all products are available in all countries.

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