
STATE OF MINNESOTA

COUNTY OF HENNEPIN

I, Susan Elizabeth Talmadge, Notary for the State of Minnesota, hereby certify that the attached Appendix to the Instructions for Use for the Prevail™ Paclitaxel-coated PTCA Balloon Catheter is a true copy of the original document filed at Medtronic, Inc.

Dated: 9/28/2022

(stamp)

Susan Talmadge
Signature of Notary Public

Title: Notary Public
My commission expires: January 31, 2027

Medtronic

Medtronic, Inc.

710 Medtronic Parkway

Minneapolis, MN 55432

USA

www.medtronic.com

tel 001-763-514-4000

fax 001-763-514-4879

To Whom It May Concern

As an authorized representative of Medtronic CRDN business unit, I hereby certify that the attached Appendix to the Instructions for Use for the Prevail™ Paclitaxel-coated PTCA Balloon Catheter is a true and accurate copy of the original document.

Yours sincerely,

Catherine Kelly

Digitally signed by Catherine Kelly
DN: cn=Catherine Kelly, o, ou,
email=catherine.kelly@medtronic.com, c=US
Date: 2022.09.28 13:44:31 -07'00'


Signed on behalf of Medtronic Inc.

Catherine Kelly

Senior Regulatory Affairs Manager

Date: September 28, 2022

Susan Talmadge 28-SEP-2022

Appendix to the Instruction for use

Prevail™ Paclitaxel-coated PTCA Balloon Catheter

DESCRIPTION OF MEDICAL DEVICE

Intended for percutaneous transluminal coronary angioplasty (PTCA) in the coronary arteries with a vessel diameter from 2.0 mm to 4.0 mm to treat de novo lesions, in-stent restenosis (ISR), and small vessel disease (SVD).

PACKAGE CONTENTS EACH PACKAGE CONTAINS THE FOLLOWING ITEMS:

Catheter - 1 pc.,
Flushing cannula - 1 pc.,
Looper device - 1 pc.,
Compliance chart - 1 pc.,
Instructions for use - 1 pc.

INFORMATION ABOUT STICKER

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Warning and informing symbols and inscriptions (if applicable)

Sample of the sticker:

Name of medical device	
Manufacturer:	Code UPN
Country of origin:	Code CFN
Manufacturer's authorized representative::	
Registration Certificate № _____ от _____	
Date of manuf./Lot/UBD: See on pack	See IFU See packaging symbols

STERILIZATION

Available upon request

Appendix to the Instruction for use
Prevail™ Paclitaxel-coated PTCA Balloon Catheter

1. This device is single-use only.
2. Device is sterile.
3. Sterilized by EO.
4. Sterility Assurance Level (SAL) - 10^{-6}
5. Term preservation of sterility for 2 years.
6. Primary packaging is sterilized in accordance with ISO 11607.

MAIN TECHNICAL DESCRIPTION

Product Characteristic	Specification with tolerance (if applicable)
Active pharmaceutical ingredient	Paclitaxel, 3.5 µg/mm ² target dose
Excipient type	Urea
Balloon Coating	FreePac™ (Paclitaxel with Urea as the excipient)
Length of balloon coating	The drug coating is applied as a single layer along the entire working length of the balloon.
Paclitaxel Drug Dose	3.5 µg/mm ² based on expanded balloon surface
Balloon diameters	(2.00, 2.25, 2.50, 2.75, 3.00, 3.50, 4.0) mm ± 10 %
Balloon Working Lengths	(10, 15, 20, 25, 30 mm) ± 1 mm NOTE: The 30 mm length is not offered on the 2.25 mm and 2.75 mm diameter balloon
Balloon Wrap Configurations	3 Folds: 2.0, 2.25, 2.50, 2.75, 3.00, 3.5 mm
	5 Folds: 4.0 mm
Catheter effective length	142 cm + 4 cm / – 5 cm
Guidewire compatibility	0.36 mm (0.014 in)
Guide catheter compatibility	5.0 Fr (0.056 in) (1.42 mm)
Shaft Diameter	Distal Shaft Outer Diameter:
	2.5 Fr (0.83 mm) for 2.0 to 3.5 mm diameters
	2.7 Fr (0.89 mm) for 4.0 mm diameter
	Transition Shaft Outer Diameter
	2.5 Fr (0.83 mm) for 2.0 to 4.0 mm (all sizes)
	Proximal Shaft Outer Diameter
	2.1 Fr (0.69 mm) for 2.0 to 4.0 mm (all sizes)
Shaft Design	Rapid Exchange (RX)

Available upon request

Product Characteristic	Specification with tolerance (if applicable)		
Nominal Balloon Pressure	8 ATM (810.6 kPa)		
Rated Burst Pressure	14 ATM (1418.6 kPa)		
Introducer Sheath Compatibility	5Fr (1.42 mm)		
Luer	Female locking luer specifications		
	Length, cm	less than 2,0	
	Width, cm	less than 1,5	
	Connecting dimensions: Inside Diameter, mm	4,315 (+0 / - 0,005)	
	Connecting dimensions: Outside Diameter, mm	7,83 (+0 / -0,1)	
	The device must meet ISO 594-1 and ISO 594-2 (ISO 80369-7) requirements		
Radiopaque Markers location	Two radiopaque marker bands swaged onto the inner shaft inside the balloon define the working length, enabling accurate balloon placement under fluoroscopy.		
Maximum Tip ID/ Inner Diameter (ID) distal section	0.0165 in $\pm$ 0.0007 in / 0.42 mm $\pm$ 0.02 mm		
Maximum Tip OD/ Outer Diameter (OD) distal section	0.0195 in $\pm$ 0.0008 in / 0.495 mm $\pm$ 0.020 mm		
Tip Length	2.5 mm $\pm$ 0.5mm		
Balloon Deflation & Inflation Time	Inflation Time to RBP shall be $\leq$ 15 seconds. Deflation Time from RBP shall be $\leq$ 30 seconds.		
Weight	3 g $\pm$ 10 %		

CONDITIONS OF TRANSPORTATION AND STORAGE

The recommendation is to “store in a cool, dry, and dark place.”

Based on the temperatures the device was subjected to during testing, storage between 15°C and 30°C is recommended for shelf life of the product. Temperature excursions down to -35°C for 2 days and up to 55°C for 75 days are supported by testing.

OPERATING CONDITIONS

The products can be used under the standard conditions of a hospital operating room environment. There are no additional assembly, installation or maintenance requirements.

The product is intended for use by trained medical personnel in the operating room: at room temperature (18°C) not higher than 25°C, relative humidity of 75%, a pressure of from 86.0 to 106.0 kPa.

Available upon request

Appendix to the Instruction for use
Prevail™ Paclitaxel-coated PTCA Balloon Catheter

SHELF LIFE 2 years.

WARRANTY SHELF LIFE 2 years.

PROCEDURE AND CONDITIONS FOR THE DISPOSAL OR DESTRUCTION OF A MEDICAL DEVICE

After use, this product may be a potential biohazard (hazard class B in accordance with SanPiN 2.1.3684-21). Handle and dispose of all such device in accordance with accepted medical practice and applicable local, regional and national laws and regulations.

THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS, WHICH MATCHES MEDICAL DEVICE

EN ISO 14971, EN ISO 13485, EN ISO 14644-1, EN ISO 14644-2, EN ISO 10555-1, EN ISO 10555-4, EN ISO 80369-7, EN ISO 12417-1, EN ISO 62366-1, EN ISO 11607-1, EN ISO 11607-2, ASTM D4169, ASTM D5276, ASTM D642, ASTM D999, ASTM D4728, ASTM F88/F88M, ASTM F1980, ASTM F2096, ASTM F1929, ASTM D4332, ISO 2233, EN ISO 15223-1, EN 1041, EN 556-1, EN ISO 11135, EN ISO 11737-1, EN ISO 11737-2, EN ISO 11138-1, EN ISO 11138-2, EN ISO 11138-7, EN ISO 10993-7, ANSI/AAMI ST72, EN ISO 10993-1, EN ISO 10993-3, EN ISO 10993-4, EN ISO 10993-5, EN ISO 10993-10, EN ISO 10993-11, EN ISO 10993-12, EN ISO 10993-17, EN ISO 10993-18, EN ISO 10993-23, EN ISO 14155

CUSTOMER SERVICE

Contact for all questions, including complaints of medical devices the Authorized Representative in Russia:

Medtronic, LLC,

123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41

Tel.: (495) 580-73-77, Fax: (495) 580-73-78,

www.medtronic.ru

Signed on behalf of Medtronic Inc.

Catherine Kelly

Senior Regulatory Affairs Manager

Signature: **Catherine Kelly**
Date: September 28, 2022

Digitally signed by Catherine Kelly
DN: cn=Catherine Kelly, o.ou,
email=catherine.kelly@medtronic.com, c=US
Date: 2022.09.28 14:03:31 -07'00'

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