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www.brightwaygroup.ru

General Manager
Bloomage Biotechnology Corporation Limited
Date: June 6, 2024
Signature: Zhao Yan

Seal of the Company: Bloomage Biotechnology Corporation Limited (Sealed)

Instructions for Use

Product Name

Hyaluronic acid solution "Miawell" (20mg/2.0mL) for intra-articular injection, filled in a syringe

International Nonproprietary Name (INN):

Sodium Hyaluronate for Injection

Intended Use:

To alleviate the pain caused by osteoarthritis of the joints, relieve pain syndrome, and improve joint mobility.

Description:

"Miawell" is a sterile, pyrogen-free, colorless, and transparent viscous solution of highly purified sodium hyaluronate, produced by biofermentation. "Miawell" contains 10mg of sodium hyaluronate per milliliter, dissolved in a physiological saline-phosphate buffer solution (pH 6.8-7.5).

Sodium hyaluronate is the salt form of hyaluronic acid, a glycosaminoglycan consisting of D-glucuronic acid and N-acetyl-D-glucosamine, widely present in the extracellular matrix of animals and humans.

Indications:

Osteoarthritis of the joints, degenerative and nutritional changes, traumatic and post-traumatic joint changes, including knee joints, hip joints, and other large synovial joints;

Relieve pain syndrome and improve joint mobility;

Used as an adjunct treatment in orthopedic surgery.

Contraindications:

Contraindicated in patients allergic (sensitive) to sodium hyaluronate. Do not inject the "Miawell" product into the knee joint in areas with infections or skin diseases.

Conditions of Use:

The use of this medical device is restricted to trained individuals within medical institutions.

Specifications:

Content: 10mg/mL sodium hyaluronate (1.0%).

Source of sodium hyaluronate: Biofermentation.

pH: 6.8-7.5.

Ingredients: Sodium hyaluronate, disodium hydrogen phosphate, sodium dihydrogen phosphate, sodium oxide, and water for injection.

Warnings:

Do not use disinfectants containing quaternary ammonium salts simultaneously when preparing the skin, as precipitation may occur in the presence of sodium hyaluronate. Do not administer this product intravenously as it may cause systemic adverse reactions. Do not use "Miawell" during intraocular surgery.

Precautions:

General Precautions:

The safety and efficacy of using this product concomitantly with other intra-articular injections have not been investigated.

The efficacy of a single treatment cycle comprising less than 5 injections has not been determined.

Therefore, all requirements of aseptic injection technique must be strictly followed.

In the presence of joint effusion, it should be removed before administering "Miawell."

Sterile contents. Pre-filled syringes are for single use only. Once the container is opened, its contents should be used immediately. Any unused "Miawell" product must be discarded.

Do not use "Miawell" if the blister pack is opened or damaged.

Storage: Keep the product in a cool, dry place between 2°C and 25°C, away from direct sunlight and extreme temperatures.

Shelf Life: 36 months.

Patient Instructions:

Following the injection of "Miawell," temporary pain and/or swelling in the injected joint may occur.

Like any invasive procedure involving a joint, patients are advised to avoid any physical activities and prolonged heavy lifting activities exceeding 1 hour (e.g., jogging and tennis) within 48 hours after the injection.

The safety and efficacy of repeated treatment cycles with "Miawell" have not been determined.

Instructions for Special Group:

The safety and efficacy of "Miawell" in pregnant women have not been determined.

The safety and efficacy of "Miawell" in lactating women have not been determined.

The safety and efficacy of "Miawell" in children have not been determined.

Side Effects:

The most common side effects include joint pain (arthritis), joint diseases (arthropathy)/arthritis, headache, pain at the injection site, and joint swelling. Allergies and hypersensitive reactions have rarely been reported. Adverse reactions at the injection site are usually mild to moderate and can be managed symptomatically.

Treating physicians should also be aware of the potential allergy risks associated with the use of any biological material.

Possible Adverse Reactions:

Topical pain, swelling, heat, and redness may occur at the injection site, and in rare cases, these symptoms are typically mild and temporary.

More severe inflammatory reactions, occasionally associated with sodium pyrophosphate crystals, have been reported following intra-articular hyaluronic acid injections.

Just like any intra-articular treatment, there is a rare possibility of developing suppurative arthritis if general injection precautions are not followed or if the injection site is not sterile.

Product Form:

"Miawell" is provided in the form of a single-use glass syringe with a capacity of 2.25mL, enclosed in a blister pack. The content of the syringe, namely sodium hyaluronate solution, is sterile. Pre-filled syringes are for single use only. Once the container is opened, its contents should be used immediately. Do not use the product if

the blister pack or the syringe is opened or damaged. For intra-articular injection into the knee joint, a needle of 22-23G and appropriate length should be employed. Please note that the needle is not included in the product packaging.

Instructions for Use:

1. "Miawell" should be injected into the joint once a week (with a one-week interval between each injection), for a total of 5 injections. It is advisable to perform a subcutaneous injection of lidocaine or a similar local anesthetic prior to the administration of "Miawell."
2. Only qualified physicians or healthcare personnel within appropriately equipped medical institutions are allowed to administer "Miawell" into the affected knee joint.
3. Do not use disinfectants containing quaternary ammonium salts simultaneously when preparing the skin, as precipitation may occur in the presence of sodium hyaluronate.
4. Do not use "Miawell" if the blister pack is opened or damaged. The product should not be used beyond the expiration date, which is indicated on the packaging.
5. All requirements for aseptic drug injection must be strictly observed.
6. If necessary (in the presence of joint effusion), the joint effusion should be removed before administering "Miawell."

The spiral cap and needle cover of the syringe should be twisted open under aseptic conditions. "Miawell" should be injected into the knee joint using a 22-23G needle. Please note that the needle is not included in the product packaging. A complete dose of 2.0mL should be injected into each knee. If directed to inject the product into both knees, a separate syringe should be used for each knee.

7. The pre-filled syringe is for single use only. It must be used immediately once removed from the packaging. Any remaining "Miawell" must be discarded and should not be used again.

Basic Parameters

S/N	Analytical Method	Acceptance Criteria (Specifications)
1	Description	Colorless, transparent, and viscoelastic (elastoviscous) gel, free from visible foreign matters
	Authenticity	a) White flake precipitate b) Phosphate test: positive
2	Quantitative determination	90.0%-120.0% of labeled amount
3	pH	6.8-7.5
4	Molecular weight	800-T200kDa (dry Dalton)
5	Minimum extractable	≥90% of labeled volume

	amount	
6	Sterility	Sterile
7	Osmotic pressure	250-350mOSM/L
8	Visible particles	No visible particles
9	Intrinsic viscosity (apparent)	4.5-19.5dl/g
10	Bacterial endotoxin	<0.5EU/ml
11	Optical density (ultraviolet absorption)	A _{257nm} 0.2 ≤ A _{280nm} 0.2 ≤

Storage Conditions:


Keep the product in a cool, dry place between 2°C and 25°C, away from direct sunlight, high temperatures and freezing.

Shelf Life:

When stored as recommended, the product has a shelf life of 36 months.

Symbols Used in Labeling

	No Reuse
	Valid Until
	Steam or Dry Heat Sterilization
	Batch No.
	Manufacture Date
	Attention! Please carefully read the instructions for use
	Manufacturer
	Do not use if packaging is damaged
	Please read the instructions for use

	No Resterilization
	Avoid Direct Sunlight
	Temperature Limit

Manufacturer: Bloomage Biotechnology Corporation Limited

Address: No. 678 Tianchen Street, Jinan Innovation Zone, Shandong Province, 250101,
P.R. China

Tel.: 86-531-82685988

Fax: 86-531-82685988

All rights reserved

Disposal and Destruction Procedures

After use, expiration, or when the packaging is damaged, the product shall be disposed of and destroyed in accordance with the local laws of the consuming jurisdiction.

Manufacturer's Warranty

The manufacturer guarantees the product will remain functional until the expiration date, provided that the original packaging is kept intact and undamaged.

Authorized representative in the Russian Federation:

Bright Way LLC


Office 16, Building 1, No. 83 Koptevskaya Street, Moscow, 125239, Russia

www.brightwaygroup.ru

Bloomage Biotechnology Corporation Limited (Sealed)

本附加证明书仅证明文书上的签名、签署人姓名和身份，不证明文书内容。
This Apostille only certifies the authenticity of the signature, the capacity of the person who has signed the public document and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued.
Visit <https://consular.mfa.gov.cn/VE-RHS/> and scan the QR code to verify the issuance of this Apostille.

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