Niall Sheehy
Notary Public

TO ALL TO WHOM THESE PRESENTS SHALL COME, I, Niall Sheehy, Notary Public, duly authorised, admitted, and sworn, practicing at 6/7 Glenthworth St., Limerick, Republic of Ireland, DO HEREBY CERTIFY that I am familiar with certificates as issued from time to time by the **Cook Medical Europe Ltd** and that the copy certificate hereunto annexed namely, *Instruction for use for Dilation Catheters, Sterile*, which I have signed my name and affixed my Seal of Office is a true and faithful photo statically reproduced copy of the original and that the original appears to me to be genuine.

AND I HAVE CERTIFIED as herein at the request of Cook Ireland Limited having it's place of business at O'Halloran Road, National Technology Park, in the city of Limerick, Ireland, and for the benefit of the proper authorities in **Russia** and for all to whom it shall come.

IN FAITH AND TESTIMONY whereof I have hereunto set my hand and affixed my Seal of Office.

Niall Sheehy

Niall Sheehy
Notary Public

29-11-2013


Date

NIALL SHEEHY
Notary Public
6/7 Glenthworth St., Limerick.

increatic D
NTEND
This

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Pays/País:		IRELAND	
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por		Niall Sheehy	
3. acting in the capacity of agissant en qualité de quien actúa en calidad de		Notary Public	
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de		-----	
Certified Attesté / Certificado			
5. at à / en	Cork	6. the le / el día	02/12/2013
7. by par / por	Department of Foreign Affairs and Trade		
8. No sous no bajo el número	0633522013		
9. Seal / stamp: 	10. Signature: Signature:  Firma:		
This Apostille certifies the authenticity of the signature and the capacity of the person who has signed the public document and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.austlii.edu.au/au/other/dfat/ie			

COOK MEDICAL EUROPE LTD.

O'Haboran Road,
National Technology Park
Limerick, Ireland.

Company Registration No. 492272

VAT Registration No. 9776234K

pancreatic Dilation Catheters GGDC

INTENDED USE

This device is used to dilate pancreatic strictures or stenotic minor papilla.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTION

Refer to package label for minimum channel size required for this device.

INSTRUCTIONS FOR USE

1. Introduce device into accessory channel of endoscope over pre-positioned wire guide. Advance in small increments until device is endoscopically visualized exiting scope.
2. Proceed with cannulation.
3. Under fluoroscopic monitoring, advance catheter to appropriate position to dilate stricture.
4. Proceed with dilation.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

Biliary Dilation Catheter FS-BDC

INTENDED USE

This device is used to dilate biliary strictures.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

device restricted to a trained healthcare professional.

INDICATIONS

Indications include those specific to ERCP.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

INSTRUCTIONS FOR USE (See fig. 1)

I. If using Zip port and a pre-positioned short wire.

1. Advance tip of catheter over pre-positioned wire guide ensuring wire guide exits catheter at zip port, then unlock short wire from Wire Guide Locking Device.
2. Introduce catheter into endoscope accessory channel and relock wire guide. Then advance catheter until it is endoscopically visualized exiting endoscope.
3. Proceed with cannulation.

4. Under fluoroscopic monitoring, advance catheter to appropriate position and dilate stricture or papillary orifice. **Note:** Dilation is achieved once distal radiopaque band has passed through stricture or papillary orifice.

Note: Previously placed wire guide may be left in position, in order to facilitate introduction of other wire guided devices.

If wire guide is to remain in place while device is withdrawn, utilize following steps:


5. Upon completion of dilation, remove the catheter from accessory channel until resistance is achieved. Unlock wire guide from Wire Guide Locking Device and perform wire guide exchange. Relock wire guide and remove catheter from wire guide.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide. (See fig. 1)

Note: For best results, wire guide should be kept wet.

1. Remove stylet.
2. Advance tip of catheter over pre-positioned wire guide until endoscopically visualized exiting endoscope. REFER TO STEPS 3-4 IN "SECTION I" THEN RESUME WITH STEP 3 BELOW.
3. Remove the catheter using standard long wire exchange technique.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

Biliary Dilation Catheters CGDC; SBDC; T7

SoehendraR/CottonR/Siegel-Cohen/van Andel

INTENDED USE

This device is used to dilate the papillary orifice, biliary strictures or pancreatic strictures.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTION

Refer to package label for minimum channel size required for this device.

INSTRUCTIONS FOR USE

1. Introduce device into accessory channel of endoscope over a pre-positioned wire guide. Advance in small increments until it is endoscopically visualized exiting endoscope and proceed with cannulation.

2. Under fluoroscopic monitoring, advance device to appropriate position to dilate stricture or papillary orifice and proceed with dilation.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

Biliary Dilation Balloon FS-BDB

INTENDED USE

This device is used to dilate strictures of the biliary tree. Supplied sterile and intended for single use.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

ions include those specific to ERCP and any procedures to be performed in conjunction with
ation.

ications to dilation include, but are not limited to: uncooperative patient, asymptomatic strictures,
y to advance the dilation balloon through the strictured area, coagulopathy, known or suspected
ration, severe inflammation or scarring near the dilation site.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis,
aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication,
hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Do not inflate the balloon prior to introduction into the scope, as this may cause damage to the scope.

The duodenoscope must remain as straight as possible with the elevator in the open position when advancing
or withdrawing the dilation balloon.

The entire dilation balloon must be extended beyond the tip of the duodenoscope and be fluoroscopically
visualized and positioned before inflation.

The dilation balloon is used in conjunction with an inflation device or a manometer and fluidfilled syringe.
DO NOT use air or gaseous substances to inflate the balloon, as this will result in reduced balloon
effectiveness.

During withdrawal of the dilation balloon from the duodenoscope, negative pressure is mandatory to
maintain balloon deflation.

A compromised balloon may prohibit removal from the endoscope accessory channel. Removal of the
endoscope along with the compromised balloon may be required.

WARNINGS

Do not exceed the recommended balloon inflation pressure as listed on the catheter tag of the dilation
balloon.

Do not advance the dilation balloon if resistance is encountered. Assess the cause of resistance to determine
if dilation should be re-attempted.

SYSTEM PREPARATION

1. Attach inflation device to dilation balloon.
2. Create and maintain a vacuum with inflation device. Dilation balloon is now ready for placement through
accessory channel of duodenoscope.
3. Remove protective sheath from balloon. Apply a water-soluble lubricant to balloon to allow easier
passage through accessory channel.

INSTRUCTIONS FOR USE (SEE FIG.1)

I. If using Intra Ductal Exchange (IDE) port and short wire:

1. Unlock short wire from Wire Guide Locking Device and advance tip of balloon catheter over pre-
positioned wire guide ensuring that wire guide exits catheter at IDE port.
2. Introduce deflated balloon into endoscope accessory channel and relock wire guide. Then advance device
in short increments until endoscopically visualized exiting scope.
3. Using fluoroscopic monitoring, continue advancement into duct until desired position is achieved. **Note:**
Position may be verified using fluoroscopic markers located at either end of balloon.
4. Inflate balloon to specified balloon pressure. Recommended pressures are listed on catheter tag of device.
5. Monitor balloon fluoroscopically during procedure.
6. To deflate balloon, decrease balloon pressure below 20 psi. Then create and lock in negative pressure.
Note: Indicator on pressure gauge should be in vacuum portion of gauge during deflation.
7. Maintain negative pressure and observe balloon for complete deflation.
8. Observe balloon size. If balloon is not completely deflated, create negative pressure again as indicated in
steps 6 & 7. Repeat these steps as necessary until balloon is completely deflated.

balloon must be thoroughly deflated, with all fluid removed, before withdrawal. A compromised balloon may prevent removal from accessory channel, requiring removal of endoscope along with balloon. Previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices.

If wire guide is to remain in place while device is withdrawn, utilize following steps:

● Prior to removing device, utilize reference marks on catheter to ensure that IDE port is within ductal system.

10. Fluoroscopically visualize radiopaque band at IDE port. Retract wire guide until radiopaque distal tip of wire guide passes band; a disengagement from wire guide lumen will occur.

11. Advance disengaged wire guide to maintain ductal access.

12. Lock wire guide into Wire Guide Locking Device.

13. To remove deflated balloon from accessory channel, straighten endoscope tip and open elevator. Maintain negative pressure while withdrawing balloon using a continuous twisting motion.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide:

Note: For best results, wire guide should be kept wet.

1. Remove stylet.


2. Advance tip of balloon catheter over pre-positioned long wire guide, ensuring that wire guide exits catheter at PWP.

3. Advance deflated balloon in short increments through accessory channel using standard long wire exchange technique until visualized exiting endoscope.

REFER TO STEPS 3-8 IN "SECTION I" THEN RESUME WITH STEP 4 BELOW:

4. To remove deflated balloon from accessory channel, straighten endoscope tip and open elevator. Maintain negative pressure while withdrawing balloon using a continuous twisting motion.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

Biliary Balloon Dilators QBD / Biliary Inflation Device

INTENDED USE

This device is used to dilate strictures of the biliary tree. The Biliary Inflation Device is used for inflation, deflation and pressure monitoring of the Biliary Balloon Dilator.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

device restricted to a trained healthcare professional.
Plasty is associated with increased incidence of pancreatitis.

INDICATIONS

Indications are specific to ERCP and any procedures to be performed in conjunction with balloon dilation, including, but not limited to: asymptomatic strictures; inability to advance the balloon dilator through the strictured area; known or suspected perforation; severe inflammation or scarring near the dilation site.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Do not inflate the balloon prior to introduction into the scope, as this may cause damage to the scope.

The duodenoscope must remain as straight as possible with the elevator in the open position when advancing or withdrawing the balloon dilator.

The entire balloon must be extended beyond the tip of the duodenoscope and be fluoroscopically visualized and positioned before inflation.

The balloon dilator is used in conjunction with the inflation device or a manometer and fluidfilled syringe.

DO NOT use air or gaseous substances to inflate the balloon, as this will result in reduced balloon effectiveness.

During withdrawal of the balloon dilator from the duodenoscope, negative pressure is mandatory to maintain balloon deflation.

A compromised balloon may prohibit removal from the endoscope accessory channel. Removal of the endoscope along with the compromised balloon may be required.

WARNINGS

Do not exceed the recommended balloon inflation pressure as listed on the catheter tag of the balloon dilator.

Do not advance balloon dilator if resistance is encountered. Assess the cause of resistance to determine if dilation should be re-attempted.

SYSTEM PREPARATION

1. Ensure the plunger handle of the inflation device is rotated completely clockwise and is pushed all the way in. To move the plunger, squeeze and hold the plunger lever located to the side of the syringe cylinder. Release the lever to lock the plunger in position.

2. Rotate the stopcock and the connecting tube fittings until tight. The stopcock is provided to facilitate purging of air.

3. Verify the sideport of the stopcock is closed by turning the “off” arm of the stopcock to the sideport.

4. While holding the inflation device, submerge the connecting tube in a minimum of 30 ml/cc of sterile water, saline or dilute contrast. Fill the device with a minimum of 15 ml/cc by compressing the plunger lever and pulling the plunger handle. Release the plunger lever to lock the plunger in place. (See fig. 1)

5. To completely purge the inflation device of air, hold the inflation device upright and compress the plunger lever while simultaneously grasping and pushing the plunger handle until all air is expelled. (See fig. 2) While maintaining the position of the plunger handle, release the plunger lever to lock the plunger in place. Release the plunger handle. To add fluid to the syringe, repeat steps 4 & 5.

6. Attach the male Luer lock fitting of the stopcock to the hub of the balloon, and secure by rotating the fittings. (See fig. 3)

7. While compressing the plunger lever, pull the plunger handle until vacuum is indicated on pressure gauge. (See fig. 4) Maintain a vacuum with the handle, then release the plunger lever. The balloon dilator is now ready for placement through the accessory channel of the duodenoscope.

INSTRUCTIONS FOR USE

Remove the protective sheath from the balloon. Apply a water soluble lubricant to the balloon to allow passage through the accessory channel.

Introduce the deflated balloon into the accessory channel of the duodenoscope, then advance in short increments until the dilator is completely visualized endoscopically.

Note: The balloon may be introduced over a pre-positioned .035 inch or smaller wire guide.

Prior to advancing the device over a pre-positioned wire guide, flush the wire guide lumen with sterile water or saline to facilitate advancement.

3. Following cannulation, contrast may be injected to fluoroscopically confirm the position of the device. If a wire guide is in place, remove the wire guide and inject contrast through the wire guide port.

4. Using fluoroscopic monitoring, continue advancement into the duct until the desired position is achieved.

Note: The position may be verified using the metal bands located at either end of the balloon as fluoroscopic markers. If the wire guide was removed for contrast injection, flush the contrast from the catheter with sterile water or saline before reintroducing the wire guide, observing the necessary cautions.

Balloon Inflation

Note: The balloon dilator must be inflated with **LIQUID ONLY**.

5. For rapid inflation, compress the plunger lever, then push the plunger handle forward until significant resistance is met. (*See fig. 5*) Release the plunger lever while maintaining plunger handle position. (*See fig. 6*)

6. Rotate the plunger handle clockwise to increase to the specified balloon pressure for the balloon dilator in use. Recommended pressures are listed on the catheter tag of the balloon dilator. Pressure may be decreased by rotating the plunger handle counterclockwise.

7. Monitor the balloon fluoroscopically during the procedure.

Balloon Deflation

8. To decrease balloon pressure and deflate, rotate the plunger handle counterclockwise until the pressure on the inflation device decreases below 20 psi. Compress the plunger lever, then pull plunger handle until negative pressure, as seen on the pressure gauge, is created. (*See fig. 7*) Release the plunger lever to lock in negative pressure. Release the plunger handle. **Note:** The indicator on the pressure gauge should be in the vacuum portion of the gauge during deflation.

9. Maintain negative pressure and observe the balloon for complete deflation. Negative pressure may be increased by continued counterclockwise rotation of the plunger handle.


10. Observe balloon size. If the balloon is not completely deflated, open the stopcock by rotating the "off" arm to close the balloon catheter port. Purge the inflation device of air. To continue balloon deflation, turn the stopcock to close the open sideport. Create negative pressure again as indicated in steps 8 & 9. Repeat these steps as necessary until the balloon is completely deflated.



Balloon Withdrawal

Caution: The balloon must be thoroughly deflated, with all fluid removed, before withdrawal.

11. To remove the deflated balloon from the accessory channel, straighten the duodenoscope tip and open the elevator. Maintain negative pressure while withdrawing the balloon using a continuous twisting motion. The previously placed wire guide may be left in position, in order to facilitate introduction of other wire guided devices.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

Esophageal Dilation Balloon HBD

INTENDED USE

This device is used to endoscopically dilate strictures of the esophagus.

PRODUCT DESCRIPTION

This balloon can be inflated to three distinct diameters, as indicated on package and catheter tag.

This device is used in conjunction with an inflation device.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Do not use this device for any purpose other than stated intended use. Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications to dilation include, but are not limited to: asymptomatic rings, webs, or strictures; known or suspected perforation; severe inflammation or scarring near the dilation site.

ADDITIONAL COMPLICATIONS

complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, hypotension, respiratory depression or arrest, cardiac arrest, thrombosis or embolism.

RECAUTIONS

Refer to package label for minimum channel size required for this device.

Do not pre-inflate the balloon.

WARNINGS

During dilation, do not inflate balloon beyond the maximum indicated inflation pressure.

The balloon must be inflated with water only.

Do not advance the balloon dilator if resistance is encountered. Assess the cause of resistance to determine if dilation should be reattempted.

INSTRUCTIONS FOR USE

1. Attach the balloon to a 60 ml (cc) inflation device with gauge to monitor balloon pressure.
2. To facilitate passage through the endoscope, apply negative pressure to the catheter.
3. Remove the protective sheath from the balloon.
4. Apply a lubricating agent to the balloon to facilitate passage through the endoscope accessory channel.
5. Maintain balloon deflation with negative pressure and introduce into the accessory channel of the endoscope, advancing in short increments until the dilator is completely visualized endoscopically. Monitor endoscopically until the balloon is in the desired position within the stricture.
6. Inflate the balloon to the pressure corresponding to the smallest balloon diameter and maintain until desired dilation is achieved. To achieve increasingly larger balloon diameters, increase pressure as indicated on catheter tag. **Note:** Do not exceed the maximum indicated inflation pressure.
7. Apply a vacuum and remove all fluid from the balloon while observing the balloon endoscopically.
8. Remove the deflated balloon from the accessory channel. **Caution:** A compromised balloon may prohibit removal from the endoscope accessory channel. Removal of the endoscope along with the compromised balloon may be required.
9. Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

Esophageal • Pyloric • Colonic

Wire Guided Balloon Dilator **HBD-W**

INTENDED USE

This device is used to dilate strictures of the gastrointestinal tract, including strictures of the esophagus, pylorus, duodenum and colon.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Balloon can be inflated to three distinct diameters, as indicated on package and catheter tag.

This balloon is used in conjunction with an inflation device or manometer and may be filled with sterile water, sterile saline or up to a 1:1 contrast medium mixture consisting of contrast and sterile saline.

CONTRAINDICATIONS

itions to dilation include, but are not limited to: uncooperative patient; asymptomatic rings, strictures; known or suspected perforation; acute corrosive injury of less than one week duration; advance device through strictured area.

ENTIAL COMPLICATIONS

ential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Do not pre-inflate the balloon.

WARNINGS


During dilation, do not inflate balloon beyond the maximum indicated inflation pressure, as this could result in overextension or bursting of the balloon.

Do not advance the balloon dilator if resistance is encountered. Assess cause of resistance to determine if dilation should be reattempted.

INSTRUCTIONS FOR USE

1. Attach the balloon to a 60 ml (cc) inflation device with gauge to monitor balloon pressure.
2. To facilitate passage through the endoscope, apply negative pressure to the device.
3. Remove the protective sheath from the balloon. **Note:** Wire guide can be removed from balloon catheter prior to use and replaced with a standard .035" wire guide to aid in bridging difficult strictures or device can be back loaded over a pre-positioned wire guide. **Note:** Prior to advancing device over a wire guide, flush wire guide lumen.
4. Apply a lubricating agent to the balloon to facilitate passage through the endoscope accessory channel.
5. Maintain balloon deflation with negative pressure and introduce into endoscope accessory channel, advancing in short increments until balloon is completely visualized endoscopically. Positioning indicator (*fig. 1*) on catheter indicates balloon clearance from endoscope accessory channel. **Note:** If using wire guide, move slide to "on" position and advance wire guide to desired location. Approximate wire guide advancement can be determined using markers located at 5-cm increments on wire guide. **Note:** Fluoroscopy is recommended if advancing wire guide without direct visualization.
6. Monitoring endoscopically, advance device until balloon is in desired position within stricture. **Note:** When balloon and wire guide are in desired location, wire guide can be secured to maintain its position by moving slide to "off" position.
7. Inflate balloon to pressure corresponding to smallest balloon diameter and maintain until desired dilation is achieved. To achieve increasingly larger balloon diameters, increase pressure as indicated on catheter tag.
8. To deflate balloon, apply negative pressure and remove all fluid from balloon while observing balloon endoscopically.
9. Remove deflated balloon from accessory channel. **Caution:** A compromised balloon may prohibit removal from endoscope accessory channel. Removal of endoscope along with compromised balloon may be required.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Please notify Cook for return authorization.

Cook Medical Europe Ltd
O'Halloran Road
National Technology Park
Limerick, Ireland
Tel: +353 61 334440 FAX: +353 61 334441
WWW.COOKMEDICAL.COM

SHELF LIFE

The "use by" date as denoted by symbols on the labelling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Dilation catheters have a shelf-life of 3 years

STERILIZATION DETAILS

Dilation Catheters are supplied sterile

STORAGE RULES

The devices are to be stored in a dry location, at room temperature.

DISPOSAL RULES

The devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

Regulatory Affairs Manager

COOK MEDICAL EUROPE
O'Halloran Road,
National Technology Park
Limerick, Ireland.
Company Registration No. 492271
VAT Registration No. 9776234K

COOK MEDICAL LTD.
O'Halloran Road,
National Technology Park
Limerick, Ireland.
Company Registration No.
VAT Registration No. 97742477

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