CERTIFICATION

State of Maryland, Montgomery County, Sct.

In the Office of the Clerk of the Circuit Court for Montgomery County

I, Barbara H. Meiklejohn, Clerk of the Circuit Court for Montgomery County, Maryland, a court of record, hereby certify that **DINNA KENCANASARI** was a commissioned/ appointed and qualified Notary Public commencing on the 27th day of April, 2018.

In Testimony Whereof, I have hereunto set my hand and affixed the seal of the Circuit Court for Montgomery County this 22nd day of January, 2020.

Barbara H. Meiklejohn

Clerk of the Circuit Court for Montgomery County

CERTIFICATION

Montgomery county:

I, Budi Isyono, hereby declare that the attached document is satisfactory to the best of my knowledge and belief.

Signatory

Budi Isyono

SWORN TO AND SUBSCRIBED before me, a notary public for the State of Maryland, this 21st day of January, 2020.

Dinna Kencanasari

Notary public, State of Maryland

My Commission Expires May 30, 2022

RVS.

INSTRUCTION FOR USE FOR MEDICAL DEVICE


Only for use by qualified specialists in dentistry

for medical device

Universal Stain and Overglaze Stain Kit for individualization of ceramic restorations

Manufactured by Dentsply Prosthetics – Also trading as: Dentsply Trubyte, Dentsply Austenal, Dentsply Ceramco, Dentsply International Inc. and Dentsply International, USA


Approved by:

Khyati Dave
Regulatory Affairs Manager


COMMONWEALTH OF PENNSYLVANIA
COUNTY OF YORK

On this 13th day of January 2020, before me a notary public, the undersigned officer, personally appeared Khyati Dave, Regulatory Affairs Manager, known to me to be the person whose name is subscribed to the within instrument, and acknowledged that she executed the same for the purposes therein contained. In witness whereof, I hereunto set my hand and official seals

Anne E. Shaffer, Notary Public

Issue date: January 2020

1	Name of the medical device	3
2	Name and address of the Manufacturer, Developer and Official Representative	3
3	Medical device classification.....	3
4	Performance characteristics of the medical device.....	3
5	Contraindications.....	3
6	Precautions	3
7	Possible adverse effects	4
8	Special warnings/Adverse reactions:.....	4
9	Product composition:.....	4
10	Description, appearance of the kit.....	5
11	General physical characteristic/properties of the device:	8
12	Composition and materials	10
13	Secondary packaging.....	11
14	Instruction for use.....	12
15	Shelf life:	13
16	Handling, transportation and storage.....	13
17	Labelling with symbol description..	15
18	Pyrogenicity.....	16
19	Information regarding sterility and frequency of use	16
20	Environmental impact	16
21	List of applicable national/international standards.	16
22	Disposal	16
23	Warranties.....	16
24	Authorized representative in the Russian Federation	16

1 Name of the medical device

Universal Stain and Overglaze Stain Kit for individualization of ceramic restorations (refer to section 9).

2 Name and address of the Manufacturer, Developer and Official Representative

Official Manufacturer and Developer: Dentsply Prosthetics – Also trading as: Dentsply Trubyte, Dentsply Austenal, Dentsply Ceramco, Dentsply International Inc. and Dentsply International, 570 West College Avenue, York, Pennsylvania, 17401, USA

Manufacturing site: DENTSPLY International Inc. Prosthetics Division – Also trading as: Dentsply Prosthetics, Dentsply Trubyte, Dentsply Austenal, Dentsply Ceramco and Dentsply Sirona, 570 West College Avenue, York, Pennsylvania, 17401, USA

Official Representative: Dentsply Sirona LLC, Andropov prospect 18 build. 6, rooms 1, 39, 40, 115432, Tel. +7 (495) 725 10 87

3 Medical device classification

Potential risk class: 2a

- GMDN code – 16187 Dental Ceramic
- Tissues of animal origin are not used during manufacture of the MD
- Human blood substances and derivatives are not the part of the medical device
- The medical device is not sterile

4 Performance characteristics of the medical device

Indication: the device is intended for individualization of ceramic restorations, incl. staining and glazing, as well as correction of contact points and occlusal surface.

Scope: dentistry

The device is intended only for use by qualified specialists in dentistry.

Patient group: for patients requiring dental restorations according to the indication.

Indication: Universal Stain and Overglaze Stain Kit for individualization of ceramic restorations is intended for staining and glazing the dental restorations made of the following materials:

- Celtra Press zirconia reinforced lithium silicate (ZLS)
- Celtra DUO zirconia reinforced lithium silicate (ZLS)
- Cercon ht, Cercon xt zirconia
- Ceramco PFZ porcelain for zirconia
- Cercon Ceram Kiss porcelain
- Ceramco3 porcelain
- Ceramco iC porcelain
- Duceram Kiss porcelain

5 Contraindications

Only the indications listed above are suitable.

6 Precautions

- Wear suitable protective eyewear, clothing and gloves

- Do not inhale abrasive dusts
- Do not ingest any of the materials
- Frameworks that show evidence of cracking or surface voids must be discarded and not processed further

7 Possible adverse effects

There were no complaints regarding the adverse effects resulting from application of the material. Follow the instructions for use and the corresponding requirements specified in material data sheet when using the material.

8 Special warnings/Adverse reactions:

Consult the Universal Stains and Overglaze Safety Data Sheets (SDS) for safe handling and usage. Use only according to the indications listed above.

Undesired side effects of these medical products are extremely rare with proper processing and use. In case of skin sensitization or rash, discontinue use and seek medical attention. Immunoreactions, such as allergies and/or local irritations (taste or oral mucosa), cannot completely be excluded on principle.

For patients with hypersensitivity to any of the ingredients, this medical device must not be used or only used under strict supervision of a treating physician/dentist. Similarly, known cross-reactions or interactions of this medical product with other materials existing in the mouth must be considered by the physician/dentist.

9 Product composition:

Universal Stain and Overglaze Stain Kit for individualization of ceramic restorations, including:

1. Universal Stain in 5 g jars; shades: Body S0, Body S1, Body S2, Body S3, Body S4, Incisal i1, Incisal i2, White, Crème, Sunset, Copper, Khaki, Olive, Mahogany, Violet, Raspberry, Purple, Grey, Chestnut, Blue, Pink – 21 pcs;
2. Universal Overglaze in 5 g jar – 1 pcs;
3. Universal Overglaze High Flu in 5 g jar – 1 pcs;
4. Universal Stain and Glaze Liquid in 15 ml vial – 1 pc.
5. Brush – 1 pc.
6. Shade guide – 1 pc.
7. Instruction for use – 1 pc.
8. Illustrated technical manual – 2 pcs.

10 Description, appearance of the kit

Kit contents:

1. Universal Stain in 5 g jars; shades: Body S0, Body S1, Body S2, Body S3, Body S4, Incisal i1, Incisal i2, White, Crème, Sunset, Copper, Khaki, Olive, Mahogany, Violet, Raspberry, Purple, Grey, Chestnut, Blue, Pink – 21 pcs;
2. Universal Overglaze in 5 g jar – 1 pcs;
3. Universal Overglaze High Flu in 5 g jar – 1 pcs;
4. Universal Stain and Glaze Liquid in 15 ml vial – 1 pc.
5. Brush – 1 pc.
6. Shade guide – 1 pc.
7. Instruction for use – 1 pc.
8. Illustrated technical manual – 2 pcs (see picture below)

Instruction for use – 1 pc.

10.1 Description of the kit components and their purpose

Component	Purpose	Patient exposure
Universal Stain/Universal Stain Body/Universal Stain Incisal 	The material is intended for staining of the dental prostheses, inlays, onlays, veneers and bridges made from the following materials: lithium silicate, zirconium dioxide, ceramics.	Permanent contact with hard tooth tissues, enamel and dentine
Universal Overglaze 	Material is intended for glazing of the dental prostheses made from lithium silicate, zirconium dioxide, ceramics.	Permanent contact with hard tooth tissues, enamel and dentine
Universal Overglaze High Flu 	Material is intended for glazing of the dental prostheses made from lithium silicate, zirconium dioxide, ceramics.	Permanent contact with hard tooth tissues, enamel and dentine
Universal Stain and Glaze Liquid 	The liquid used to obtain creamy consistency of the glaze and stain	Permanent contact with hard tooth tissues, enamel and dentine

Brush 	Brush for application of stain and/or glaze.	Possible contact with intact skin
Shade guide 	Intended to facilitate the selection of the required shade for staining of the ceramic restorations.	Possible contact with intact skin

11 General physical characteristics/properties of the device:

11.1 Stains and glaze

- Appearance: glaze – Opaque paste; stains – opaque paste in different colours (depending on the shade).
- Flexural strength: ≥ 100 MPa
- Chemical solubility: $12 \mu\text{g}/\text{cm}^2 (<0.05\%)$ for stains; $6 \mu\text{g}/\text{cm}^2 (<0.05\%)$ for glaze.
- Radioactivity: <22 м.д. ppm for stains and <40 ppm for glaze
- Coefficient of Thermal Expansion (CTE): $7,5 \pm 0,5 \times 10^{-6} \text{K}^{-1}$ (25-500°C),
- Glass transition point: $535 \pm 20^\circ\text{C}$;
- Softening point: $576 \pm 20^\circ\text{C}$
- Linear shrinkage during firing: Less than 16 %
- Porosity: The number of pores with diameter $> 30 \mu\text{m}$ should not exceed 16, and the number of pores with diameter in the range from 40 to $150 \mu\text{m}$ should not exceed 6 on the 1 mm^2 surface area. There should be no pores with diameter $> 150 \mu\text{m}$.
- Staining resistance, number of stains on samples upon exposure in staining solution: 0.
- Producibility: The material should be brought to desired consistency using special liquid, while there should be no lump formation. Excess of moisture should be well removed during condensation. Sample must hold well

11.2 Universal Stain and Glaze Liquid

Physical indicator	Values
Appearance	Transparent colourless liquid
Odour	Odourless
Freezing point	-77°C
Initial boiling point	$207,5^\circ\text{C}$
Solubility	Completely soluble in water
Viscosity, +/-10%	104 mPa.s
Vial dimensions	Height: $70 \text{ mm} \pm 10\%$ Diameter: $25 \text{ mm} \pm 10\%$
Liquid volume	Not less than 15 ml

11.3 Weight and dimensions of the kit components

11.3.1 Stains and glaze

Name	Weight, g	Weight in primary packaging, g.
Universal Stain/Universal Stain Body/Universal Stain Incisal	Not less than 5	Not less than 27
Universal Overglaze, Universal Overglaze High Flu	Not less than 5	Not less than 27

11.3.2 Liquid

Name	Volume, ml, not less than	Weight, g, not less than	Weight in primary packaging, g, not less than
Stain and Glaze Liquid	15	15	24

11.3.3 Other materials

Name	Weight, g	Dimensions, mm
Shade guide	21±10%	10.7 x 8.1 x 3±10%
Brush	brush: 2.94 ± 0.3 g Protective cap for brush: 0.16 g ± 10%	166.5± 10% Handle: 132 mm ± 10%; Tip: 29 mm ± 10%; Brush (bristles): 5.5 mm ± 10%, Protective cap for brush: 21 mm x 220 mm ± 10%

12 Composition and materials

12.1 General composition of the material (stains Universal Stain/Universal Stain Body/Universal Stain Incisal and Universal Overglaze, Universal Overglaze High Flu)

Stains and glaze are in the form of paste. Paste consists of powder (70%) and gel (30%).

12.1.1 Powder composition is presented below.

Component	Formula
Silicon oxide	SiO ₂
Aluminium oxide	Al ₂ O ₃
Sodium oxide	Na ₂ O
Potassium oxide	K ₂ O
Lithium oxide	Li ₂ O
Calcium oxide	CaO
Barium oxide	BaO
Boron trioxide	B ₂ O ₃
Cerium oxide	CeO ₂
Terbium oxide	Tb ₂ O ₃
Antimony oxide	Sb ₂ O ₃
Tin oxide	SnO ₂
Titanium oxide	TiO ₂
Stains	-

12.1.2 Composition of the gel used for paste preparation is presented below.

Ammonium nitrate, butylene glycol, propylene glycol, silicone dioxide, colloid, tylose H.

12.2 Universal Stain and Glaze Liquid, 15ml

Universal Stain and Glaze Liquid, 15ml	Contents
	Butylene glycol

12.3 Brush materials

Bristles	Finest Kolinsky Hair
Tip	CuZn MS
Handle	Wood (birch)
Protection cap for bristles	LDPE


12.4 Shade guide materials

Polycarbonate.

12.5 Drugs

Not applicable. The product does not contain any drug substances or drugs.

13 Secondary packaging

The device is supplied in a carton box with IFU included.

14 Instruction for use

14.1 Framework preparation:

Follow the current corresponding Directions for Use for preparing the substructure material.

Instructions for use:

- Celtra Press zirconia reinforced lithium silicate – item 50539994
- Celtra DUO zirconia reinforced lithium silicate – item 50539985
- Cercon ht, Cercon xt – item 50537617
- Ceramco PFZ* porcelain for zirconia – item 090862
- Cercon Ceram Kiss, porcelain – item 28078
- Ceramco3, porcelain – item 301700
- Ceramco iC, porcelain s – item 745-215
- Duceram Kiss, porcelain – item 22147

14.2 Application of stains and glaze:

General steps:

1. Take required amount of Universal Stain using plastic or zirconium spatula. Put on the mixing pallet. If a thinner consistency is required – dilute the material using Stain and Glaze Liquid
2. Apply thin layer of stain or glaze on the whole surface of the ingot.
3. Perform annealing depending on the ingot material.

(detailed information is given in the Instruction for Use).

Note: If the stain or glaze inside the jar has separated, mix thoroughly with a glass, plastic or zirconia spatula. Withdraw desired amount of Universal Stain or Glaze from the jar and place it on the mixing palette. If a thinner consistency is desired, dilute the material with the Dentsply Sirona Stain and Glaze Liquid (REF 601315/15mL, 601350/50mL). (See Figure 1)

Figure 1.

Apply a thin layer of the stain or glaze material all over the crown surface in the usual manner using a brush. Make sure not to apply the glaze material too thick or too thin. Too thin of a layer will result in a more matte finish. Too thick of a layer can result in puddling and/or pitting of the material. (See Figure 2)

Figure 2.

If a more intensive shade effect is desired, it is achieved by several staining applications and firings. However, the application of too many stain layers can result in an unnatural shade effect.

Stains and glaze firing:

Please use the suggested furnace cycle/program shown below:

	Drying (min)	Preliminary heating (min)	Low temperature	Heating rate	High temperature	Vacuum start temperature	Vacuum stop temperature	Hold Time Vac	Hold Time Air (min)	Cool time (min)
Celtra Press (full contour)	2	2	400C	55C/min	760C	0	0	0	2:00	5
Celtra Press (decrease or secondary annealing of full contour)	2	2	400C	55C/min	750C	0	0	0	2:00	5
Celtra Duo	2	2	500C	55C/min	820C	0	0	0	1:30	3
Cercon ht	3	3	450C	55C/min	820C	0	0	0	1:30	6
Ceramco PFZ*	5	5	450C	60C/min	850C	0	0	0	0:30	6
Cercon Ceram Kiss	2	3	450C	55C/min	800C	0	0	0	1:00	6
Ceramco3	3	3	650C	70C/min	935C	0	0	0	0:30	0
Ceramco iC	3	3	500C	35C/min	750C	500C	750C	1	1:30	0
Duceram Kiss	3	1	575C	55C/min	890C	0	0	0	1:00	0

* Drying time depends on the type of furnace

* Additional Stain/glaze and characterization firing cycles can be conducted with the same firing parameters.

Do not go below or above the temperature range 740 °C to 950 °C.

15 Shelf life:

Stability: close the jars and bottles tightly right after use. Use the material by expiration date.

The difference between the shelf life indicated on the packaging and opening date: shelf life corresponds to the actual period of product validity.

Shelf life of the kits:

- Universal Stain and Overglaze Stain kit – 4 years

Shelf life of the kit components:

Name	Shelf life (years)
Universal Overglaze, Overglaze High Flu	4
Universal Stain	4
Stain and Glaze Liquid	5

16 Handling, transportation and storage

16.1 Transport and storage conditions

Avoid direct exposure to sunlight and heat, protect from moisture, keep the lid tightly closed when not in use.

Storage and transportation temperature from 10 to 30 °C

Humidity: less than 70%

After each use, immediately close the jars and bottles with the original cap.

16.2 Conditions of use:

Avoid direct exposure to sunlight and heat, protect from moisture, keep the lid tightly closed when not in use.

Temperature: from 10 to 30 °C

Humidity: less than 70%.

Materials are also used inside the oral cavity as intended, in the presence of biological fluids at nominal temperature values: from 32 to 42° C

17 Labelling with symbol description.

17.1 Symbols used on primary and secondary packaging:

	Medical device manufacturer
	Official representative in EU
	Catalog no.
	Lot No.
	Expiration date
	Keep dry
	Date of manufacture
	See Instruction for Use
	Caution! Refer to SDS
	CE mark – indicates compliance with EU standards
	Only for professional use.

18 Pyrogenicity

Not applicable. The materials of the kit are non-sterile.

19 Information regarding sterility and frequency of use

Not applicable. The kit is non-sterile and is intended for single use. Secondary sterilization and disinfection are not applicable.

20 Environmental impact

The product doesn't have any negative environmental effects during life cycle.

21 List of applicable national/international standards.

- ISO 6872:2015 – Dentistry. Ceramic materials
- ISO 14971 – MEDICAL DEVICES – Application of risk management to medical devices
- ISO 9693-2:2016 – Dentistry. Compatibility tests. Part 2. Ceramic-ceramic systems

22 Disposal

The product should be disposed of according to local and medical facility regulations.

Disposal according to sanitary regulations and standards (SanPiN) of the Russian Federation 2.1.7.2790-10 "Public health requirements for handling of medical waste". Hazard class of the medical waste – B.

The product should be used only by qualified staff.

23 Warranties

Manufacturer guarantees compliance of the medical device characteristics to specifications given in operating instructions provided the product is transported, stored and used according to Instruction for Use supplied by manufacturer.

Manufacturer takes no responsibility and does not pay compensation for possible damage or accidents caused by:

- Usage of components not included in the kit and which may interfere with normal operation.
- Non-compliance with instructions for use.

The user is responsible for testing the material for its suitability and use for any purpose that is not explicitly specified in the instructions!

Warranties are valid within the shelf life of the product.

24 Authorized representative in the Russian Federation

Dentsply Sirona LLC, Andropov Avenue 18, build. 6, rooms 1 39 40, Moscow, 115432
Tel. +7 (495) 725 10 87

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