
Instructions for Use

**Dentsply
Sirona**
implants
Dentsply Implants
Manufacturing GmbH
Rodenbacher Chaussee 4
63457 Hanau, Germany

Bensheim, November 21st 2022

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General information

Xive® threaded dental implants are intended for surgery implantation in the alveolar and/or basal bones of the mandible or maxilla to provide support and retention of the denture (e. g. bridge, single tooth, cover prosthesis) to restore the masticatory function of the patient.

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Indications

The Xive® S plus Implants D 3.4 – D 5.5 are for single-stage or two-stage surgical procedures and cemented, removable or screw retained restorations. The Xive® S plus Implants D 3.4 – D 5.5 may be used for immediate placement and function on single tooth and/or multiple tooth applications when adequate primary stability is achievable, with appropriate occlusal loading, in order to restore chewing function. Multiple tooth applications may be splinted.

TempBase

Placement head for implant insertion; serves as a basis for temporary restoration and index registration in conjunction with the TempBase Cap

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- Clotting disorders e.g. anticoagulant therapy, inherited or acquired clotting disorder
- Bone and wound healing disorders e.g. uncontrolled diabetes, metabolic or systemic disorders which affect bone or wound healing; general restrictions on elective surgery apply to tobacco and alcohol abuse
- Immunosuppressive therapy e.g. chemotherapy and radiation therapy
- Intraoral infection or inflammation e.g. periodontitis
- Uncontrolled parafunctional habits e.g. bruxism
- Inadequate oral hygiene
- Patient not prepared to undergo total oral rehabilitation
- Untreatable occlusal/articulation disorders and insufficient interarch space
- Insufficient height and/or width of bone / insufficient soft tissue coverage

Warnings

Please read the instructions for use prior to using Xive® implants. Xive® implants may only be used for their intended purpose according to the general rules for dental/surgical treatment, occupational safety and accident prevention.

Xive®/Friadent® components must only be used for medical/dental procedures with the Xive® Implant system. If the indication and use are unclear, treatment should not be carried out until all of these points have been clarified. The following instructions are not adequate to allow inexperienced providers to provide professional implant dentistry. We recommend that an experienced user provide proper use instructions. Xive® implants must only be used by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.

Dentsply Implants Manufacturing GmbH is not liable for damage resulting from use outside the intended use of the product.

Oral surgery may cause general complications. These include:

- Components in the patient's mouth can be swallowed or aspirated.
- Nerve injury
- Implant dislocation in the maxillary sinus
- Perforation of the maxillary sinus, mandibular base, lingual plate and the inferior alveolar canal
- Fracture of the jaw bone
- Necrosis due to inadequate cooling or excessive torque
- Excessive bone loss requiring intervention
- Edema
- Loosening after ingrowth/implant loss due to insufficient osseointegration

- In rare cases, implant breakage can occur under extreme unfavorable loading conditions (prosthetic overload, severe bone resorption)
- Periimplantitis
- Infections
- Post-operative bleeding
- Iatrogenic trauma
- Esthetic or functional complications
- Suture dehiscence

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- Jammed or over-tightened placement head (TempBase)
- Loss of implant due to an over-tightened cover screw

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- Prior to each procedure it must be ensured that all necessary components, instruments and materials are available in the required quantities and are fully operational.
- Always wear protective clothing for your own safety.
- Position the patient so that the danger of aspiration of components is minimized.
- All components that are used intraorally must be secured to prevent aspiration or swallowing.
- Observe the maximum speed.
- Comply with the recommended speed during preparation with drills; in particular, avoid a reduction in speed.
- Cutting instruments should be replaced after about 20 uses or if they are damaged or dull.
- Ensure adequate cooling.
- Extensive manipulation of the implant (e.g. change of prosthetic components) should be carried out at the earliest after the osseointegration phase, unless absolutely necessary.
- No product intended for single use may be reused. Noncompliance generally leads to the risk of infection and complications due to a loss of accuracy of fit and the precision of the components. This may lead to fractures through to implant loss.

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MRI Safety Information

All Xive® products, which remain in the patient's body, have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the Xive® products in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

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There are 2 types of instruments:

- Type 1 with shaft in accordance with ISO 1797-1

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Not all products are available in all countries

Status of information
IFU 2601-en, RU Rev. 032 A, 2022-11

**Dentsply
Sirona**
Implants

Instructions for Use

- Type 2 with shaft in accordance with ISO 1797-1 and hexagon clamping system (HXSS)

Instruments with the HXSS used with the appropriate contra-angle offer improved transmission.

NOTE: Use the implant driver for D 3.8-D 5.5 implants when seating the D 3.8-D 5.5 screw implant (without TempBase).

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The laser ring markings and grooves correspond to the implant lengths. The laser band markings are for orientation purposes. In addition, the implant site depth can be checked with a depth gauge.

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Implant diameter	Drill sequence (Drill without designation = Twist Drill Anterior or Posterior) The round drill can be used to smoothen the ridge or to mark the cortical bone.	Rec. drill speed
D 3.4 mm	D2.0 → D3.0 → D3.4 → Twist Drill Crestal D3.4	800 rpm
D 3.8 mm	D2.0 → D3.0 → D3.4 → D3.8 → Twist Drill Crestal D3.8	800 rpm
D 4.5 mm	D2.0 → D3.0 → D3.4 → D3.8 → D4.5 → Twist Drill Crestal D4.5	800 rpm
D 5.5 mm	D2.0 → D3.0 → D3.4 → D3.8 → D4.5 → D5.5 → Twist Drill Crestal D5.5	800 rpm

ATTENTION: The Twist Drill Crestal must generally be used (see the following information). The only exception is missing cortical bone.

Final preparation is accomplished with a Twist Drill Crestal at a recommended speed of 800 rpm. The depth is dependent on the bone quality/density:

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The depth of the Twist Drill Crestal should not exceed 4 mm for implants 8 mm in length. The Friadent® Drill extension can be used in case of limited space. It must always be used with the long internal irrigation tube.

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In case of a very dense cortical bone (type D I), the Xive® Tap must be used at a speed of 15 rpm after implant site preparation with the Twist Drill Crestal. If the thread region of the tap is countersunk, then the site preparation depth is 6 mm. To use the tap without an integrated contra-angle, the implant driver for screw implants D 3.4 and TempBase must be used (with O-ring to improve TempBase seat. Replace the O-ring if the clamping effect decreases).

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Use the Xive® implant driver with ISO shaft to drive in Xive® implants with the corresponding ratchet insert and Xive® ratchet or contra-angle. Set the arrow on the switch button of the Xive® ratchet according to the direction of insertion (clockwise).

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1. Open the box and remove the non-sterile outer blister.
2. Open the outer blister and hand the sterile inner blister to the dentist or drop it onto the surgical tray.
3. Hold the inner blister horizontal while opening it to prevent the cover screw from falling out.
4. Remove the implant holder (shuttle) from the inner blister.
5. Place the implant driver for contra-angle or ratchet into the internal hex of the TempBase. Check that it is firmly in place and if necessary secure the implant against rotation by pressing the opposing wings of the shuttle together gently.
6. Gently bend open the shuttle and remove the implant. (Avoid touching the osseous threaded section).
7. Place the implant into the prepared site.
8. Screw the implant into the prepared site at 15 rpm until the top of the implant is 0.4 mm supra-crestal. (The transition between implant and TempBase can be determined by the Friadent® color-coding.)

ATTENTION: If the torque reaches 50 Ncm during seating then the TempBase must be removed and the implant screwed in using the implant driver for screw implants D 3.8-D 5.5. After implant insertion, the TempBase is screwed into the implant for index registration or to fabricate a temporary restoration.

The following options are available at the time of implant placement

Submerged healing phase

The TempBase is removed by loosening the retaining screw with the 0.9 mm hex driver. The cover screw is screwed into the implant and the soft tissue is sutured over it. Make sure no soft tissue is trapped between the implant and the cover screw. When using the platform switch concept (PS), use of a cover screw with the next smallest diameter is recommended.

You will find more information on this in the PS Step-By-Step brochure.

Index registration

The position of the implant is transferred to the master model with the aid of the TempBase. The TempBase Cap is press-fitted onto the TempBase until it snaps into position. The Cap is then luted to the surgical stent with light-cured resin.

Remove the stent with the TempBase Cap and remove the TempBase from the implant. Cover the implant with the cover screw and suture the soft tissue back into place. A provisional restoration can be fabricated and delivered at the second stage. Model fabrication is accomplished in the laboratory by mounting the TempBase on an implant analog.

Fabricating a temporary restoration with the TempBase

Functional and non-functional immediate loading

Functional immediate loading is only possible if all Xive® implants have sufficient primary stability. The functional immediate loading may only be performed on single tooth replacement, on primarily splinted or on implants secondarily connected by stable rigid components.

The TempBase remains in the implant after implant insertion. Retighten the retaining screw with the hex driver without torque control. Shorten the TempBase Cap and fabricate a provisional with a template. Fiber material is fastened to the TempBase Cap groove for complex restorations.

Uncovering the implant after the healing phase

After the healing phase (3-4 months), the implant is uncovered and the cover screw removed. The healing phase should be adequately extended if augmentation procedures were performed. The inside chamber of the implant is cleaned and a gingiva forming component of matching diameter and soft tissue thickness is screwed into position. Ensure that the gingiva forming component is seated correctly. The soft tissue is sutured tightly around the gingiva forming component. Prosthetic restoration can begin once the soft tissues have completely healed.

Post-operative care

Instruct the patient regarding the need for regular oral hygiene. For the first 7-10 days post-op, oral hygiene can be supplemented by a suitable anti-microbial rinse. Mechanical loading of the implant site should be avoided in the post-operative period. A relined provisional restoration can be worn for esthetic purposes.

Delivery – Storage – Disposal

Xive® S plus implants are manufactured from commercially pure grade 2 titanium under GMP standards. The implants are available in various lengths and diameters. For ease of identification, the various implant diameters are colour-coded on the respective packaging (D3.4 mm = gray, D3.8 mm = yellow, D4.5 mm = blue, D5.5 mm = red). For type and size see product label.

Xive® S plus implants are gamma sterilized and are intended for single use only. **Xive® S plus implants must not be re-sterilized.** Do not use the Xive® S plus implants after expiration of the sterilization date or if the seal flap is damaged or missing or if the outer blister is damaged. Store Xive® S plus implants at room temperature from +15°C to +30°C in the original packaging under conditions normal for dental practices. Do not expose to direct sunlight. For disposal of the packaging and components comply with the currently applicable national waste disposal regulations in Russian Federation: Implants and cover screws for implants after use - as waste class B (epidemiologically hazardous waste).

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Not all products are available in all countries

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Mechanical processing

Manual pre-cleaning is necessary to ensure efficient mechanical cleaning. Suitable washer disinfectors as well as cleaning and neutralization agents are to be used for mechanical cleaning. The instructions from the manufacturer of the washer disinfectant must be followed. Cleaning and neutralization agents are to be dosed and used in accordance with the manufacturer's instructions. A cleaning program with thermal disinfection ($A_0 \geq 3000$) is recommended. Either demineralized water or water which satisfies this level of purity is recommended for cleaning.

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- Preparation for sterilization:** Place components in a sterilization pouch which is legally marketed (for the US market: FDA-cleared) for use with the recommended sterilization parameters. Every sterilization package must have a sterilization indicator and sterilization date.
- Sterilization:**

Method	Cycle	Temperature	Exposure Time*	Dry Time
Steam	Dynamic air removal (prevacuum)	132 °C (270 °F)	4 Min.	20 Min.
		135 °C (275 °F)	3 Min.	
Steam	Gravity displacement	121 °C (250 °F)	30 Min.	20 Min.

* Minimum exposure times, the operating times are longer and may vary depending on the device.

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Stomatology. Dental implantology. The products should be operated at a temperature of +32 to +42 °.

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Copyright and Trademark

All rights reserved. No part of these instructions for use covered by the copyright hereon may be reproduced or used in any form or by any means graphic, electronic or mechanical, including photocopying, recording, taping, or information storage and retrieval systems without the permission of Dentsply Implants Manufacturing GmbH.

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Explanation of Symbols

Symbol	Explanation
	Manufacturer
	Date of manufacture
	Use-by date
	Batch Code
	Catalog Number
	Sterilized using irradiation
	Do not resterilize
	Non-Sterile
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Do not reuse
	Electronic Instruction for Use
	European conformity (CE) mark
	Prescription only: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician.
	Double sterile barrier system
	Medical device

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Not all products are available in all countries

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Indications

Prosthetic concept

Single-tooth restoration anterior region (tight gaps) 12 + 22 and 32 - 42
Splinted single-tooth restoration region 32 - 42

Prosthetic treatment

Non-functional immediate loading
Two-stage procedure

Time of implant placement

Delayed immediate implant placement
Late implant placement

Healing method

Submerged
Transgingival with gingival forming components

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- Clotting disorders e.g. anticoagulant therapy, inherited or acquired clotting disorder
- Bone and wound healing disorders e.g. uncontrolled diabetes, metabolic or systemic disorders which affect bone or wound healing; general restrictions on elective surgery apply to tobacco and alcohol abuse
- Immunosuppressive therapy e.g. chemotherapy and radiation therapy
- Intraoral infection or inflammation e.g. periodontitis
- Uncontrolled parafunctional habits e.g. bruxism
- Inadequate oral hygiene
- Patient not prepared to undergo total oral rehabilitation
- Untreatable occlusal/articulation disorders and insufficient interarch space
- Insufficient height and/or width of bone
- Insufficient soft tissue coverage
- Proven hypersensitivity to one of the metals in the alloy

Warnings

Please read the instructions for use prior to using Xive® implants. Xive® implants may only be used for their intended purpose according to the general rules for dental/surgical treatment, occupational safety and accident prevention. If the indication and use are unclear, treatment should not be carried out until all of these points have been clarified. The following instructions are not adequate to allow inexperienced clinicians to administer professional implant dentistry. We recommend that an experienced user provide proper use instructions. Xive® implants must only be used by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning. Dentsply Implants Manufacturing GmbH is not liable for damage resulting from use outside the intended use of the product.

ATTENTION: The size D 3.0 refers to the prosthetic diameter of the implant, the thread diameter is **3.25 mm!**

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- Components in the patient's mouth can be swallowed or aspirated.
- Nerve injury
- Implant dislocation in the maxillary sinus
- Perforation of the maxillary sinus, mandibular base, lingual plate and the inferior alveolar canal
- Fracture of the jaw bone
- Necrosis due to inadequate cooling or excessive torque
- Excessive bone loss requiring intervention
- Edema
- Loosening after ingrowth/implant loss due to insufficient Osseointegration
- In rare cases, implant breakage can occur under extreme unfavorable loading conditions (prosthetic overload, severe bone resorption)
- Periimplantitis
- Infections
- Post-operative bleeding
- Iatrogenic trauma
- Esthetic or functional complications
- Suture dehiscence

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- Prior to each procedure it must be ensured that all necessary components, instruments and materials are available in the required quantities and are fully operational.
- Observe the maximum speed.
- Comply with the recommended speed during preparation with drills; in particular, avoid a reduction in speed.
- Cutting instruments should be replaced after about 20 uses or if they are damaged or dull.
- Ensure adequate cooling.
- Always wear protective clothing for your own safety.
- Position the patient such that the danger of aspiration of components is minimized.
- All components that are used intraorally must be secured to prevent aspiration or swallowing.
- No product intended for single use may be reused. Noncompliance generally leads to the risk of infection and complications due to a loss of accuracy of fit and the precision of the components. This may lead to fractures through to implant loss.

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During preoperative planning, it must be established that there is an adequate height and width of bone to allow placement of the implant. The width of the oral and vestibular bone lamella should be at least 1.5 mm. Appropriate radiographic preoperative analysis is necessary to determine the position of important anatomical structures. Augmented areas must show complete and mechanically stable regeneration prior to receptor site preparation. Using Friadent® Select components improves surgical planning.

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NOTE: Use the implant driver for D 3.0 implants when seating the D 3.0 screw implant (without TempBase).

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The Xive® depth stops ensure that the correct depth is reached. Pull the depth stop, with the colour-coding towards the top, over the drill tip. Additional instruments are not required for depth stop removal.

Implant diameter	Drill sequence (Drill without designation = Twist Drill Anterior or Posterior) The round drill can be used to smoothen the ridge or to mark the cortical bone.	Rec. drill speed
D 3.0 mm	D 2.0 → D 3.0 → Twist Drill Crestal D 3.0	800 rpm

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The Friadent® Drill extension can be used in case of limited space. It must always be used with the long internal irrigation tube. Ensure there is sufficient internal and external irrigation when drilling.

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1. Open the box and remove the non-sterile outer blister.
2. Open the outer blister and hand the sterile inner blister to the dentist or drop it into the surgical tray.
3. Hold the inner blister horizontal while opening it to prevent the cover screw from falling out.
4. Remove the implant holder (shuttle) from the inner blister.
5. Place the implant driver for contra-angle or ratchet into the internal hex of the TempBase. Check that it is firmly in place and if necessary secure the implant against rotation by pressing the opposing wings of the shuttle together gently.
6. Gently bend open the shuttle and remove the implant. (Avoid touching the osseous threaded section).
7. Place the implant into the prepared site.

8. Screw the implant into the prepared site at 15 rpm until the top of the implant is 0.4 mm supra-crestal. (The transition between implant and TempBase can be determined by the colour-coded TempBase).

ATTENTION: If the torque reaches 50 Ncm during seating, then the TempBase **must** be removed and the implant screwed in using the implant driver for screw implants D 3.0. After implant insertion, the TempBase is screwed into the implant for index registration or to fabricate a temporary restoration.

The following options are available at the time of implant placement:

Submerged healing phase

The TempBase is removed by loosening the retaining screw with the 0.9 mm hex driver. The cover screw is screwed into the implant and the soft tissue is sutured over it.

Index registration

The position of the implant is transferred to the master model with the aid of the TempBase. The TempBase Cap is press-fitted onto the TempBase until it snaps into position. The Cap is then luted to the surgical stent with light-cured resin. Remove the stent with the TempBase Cap and remove the TempBase from the implant. Cover the implant with the cover screw and suture the soft tissue back into place. A provisional restoration can be fabricated and delivered at the second stage. Model fabrication is accomplished in the laboratory by mounting the TempBase on an implant analog.

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The TempBase remains in the implant after implant insertion. Retighten the retaining screw with the hex driver without torque control. Shorten the TempBase Cap and fabricate a provisional with a template. If the single-tooth restoration needs to be splinted, fibre material is fastened to the TempBase Cap groove for reinforcement.

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Post-operative care

Instruct the patient regarding the need for regular oral hygiene. For the first 7-10 days post-op, oral hygiene can be supplemented by a suitable anti-microbial rinse. Mechanical loading of the implant site should be avoided in the post-operative period. A relined provisional restoration can be worn for aesthetic purposes.

Delivery - Storage - Sterilization

Xive® S plus implants are manufactured from commercially pure grade 2 titanium under GMP standards. The implants are available in various lengths and diameters. For ease of identification, the various implant diameters are colour-coded on the respective packaging (D3.0 mm = bronze-colored). For type and size see product label.

Xive® S plus implants are **gamma sterilized** and are intended for **single use only**. **Xive® S plus implants must not be re-sterilized**. Do not use the Xive® S plus implants after expiration of the sterilization date or if the seal flap is damaged or missing or if the outer blister is damaged. Store Xive® S plus implants at room temperature +15 to +30 °C in the original packaging under conditions normal for dental practices. Do not expose to direct sunlight.

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Mechanical processing

Manual pre-cleaning is necessary to ensure efficient mechanical cleaning. Suitable washer disinfectors as well as cleaning and neutralization agents are to be used for mechanical cleaning. The instructions from the manufacturer of the washer disinfectant must be followed. Cleaning and neutralization agents are to

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Dentsply
Sirona
implants

Instructions for Use

be dosed and used in accordance with the manufacturer's instructions. A cleaning program with thermal disinfection (A₀ ≥ 3000) is recommended. Either demineralized water or water which satisfies this level of purity is recommended for cleaning.

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- Washer disinfector, type G 7836 CD, Miele & Cie. KG, Gütersloh
- Cleaning program VARIO-TD
Pre-rinsing with cold water 4 min
Cleaning 10 min from 45-55 °C
Neutralization 6 min
Intermediate rinsing with cold water 3 min
Thermal disinfection 5 min at 90 °C (A₀ ≥ 3000)
Drying 30 min at 105 °C
- Cleaning agent Neodisher MediClean (0.5 %) / neutralization agent Neodisher Z (0.1 %), Dr. Weigert, Hamburg.

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1. **Preparation for sterilization:** Place components in a sterilization pouch which is legally marketed for use with the recommended sterilization parameters. Every sterilization package must have a sterilization indicator and sterilization date.
2. **Sterilization:**

Method	Cycle	Temperature	Exposure Time*	Dry Time
Steam	Dynamic air removal (prevacuum)	132 °C	4 Min.	20 Min.
		135 °C	3 Min.	
Steam	Gravity displacement	121 °C	30 Min.	20 Min.

* Minimum exposure times, the operating times are longer and may vary depending on the device

3. **Storage:** Store the sterilized components dry and dust-free at room temperature.

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Copyright and Trademark

All rights reserved. No part of these instructions for use covered by the copyright hereon may be reproduced or used in any form or by any means graphic, electronic or mechanical, including photocopying, recording, taping, or

information storage and retrieval systems without the permission of Dentsply Implants Manufacturing GmbH.

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
Instructions for Use

Bensheim, November 21st 2022

Dr. Aline Frare Mocruha
Regulatory Affairs Manager

**Dentsply
Sirona**
Implants
Dentsply Implants
Manufacturing GmbH
Rodenbacher Chaussee 4
63457 Hanau, Germany

**Dentsply
Sirona**
Implants

For Instructions for Use
and symbols glossary refer to
ifu.dentsplysirona.com

Xive® Implant System Xive® TG plus Implants

Description of the System

The Xive® implant system is designed for replacing one or more missing teeth in patients to restore chewing function. It comprises dental implants, cover screws, transmucosal abutments, other prosthetic components, as well as surgical and prosthetic instruments.

General information

Xive® threaded dental implants are intended for surgery implantation in the alveolar and/or basal bones of the mandible or maxilla to provide support and retention of the denture (e.g. bridge, single tooth, cover prosthesis) to restore the masticatory function of the patient.

Ordering information and Dimensions

See the current product catalog for updated and detailed information about the products.

Material

See the current product catalog for material information of the products.

Indications

The Xive® TG implant system is for single-stage surgical procedures and cemented, removable or screw retained restorations. The Xive® TG implant system may be used for immediate placement and function on single tooth and/or multiple tooth applications when adequate primary stability is achievable, with appropriate occlusal loading, in order to restore chewing function. Multiple tooth applications may be splinted.

Contraindications

Contraindications customarily associated with elective oral surgery should be observed when selecting patients. These include:

- Clotting disorders e.g. anticoagulant therapy, inherited or acquired clotting disorder
- Bone and wound healing disorders e.g. uncontrolled diabetes, metabolic or systemic disorders which affect bone or wound healing, general restrictions on elective surgery apply to tobacco and alcohol abuse
- Immunosuppressive therapy e.g. chemotherapy and radiation therapy
- Intraoral infection or inflammation e.g. periodontitis
- Uncontrolled parafunctional habits e.g. bruxism
- Inadequate oral hygiene
- Patient not prepared to undergo total oral rehabilitation
- Untreatable occlusal/articulation disorders and insufficient interarch space
- Insufficient height and/or width of bone
- Insufficient soft tissue coverage
- Proven hypersensitivity to one of the metals in the alloy

Warnings

Please read the instructions for use prior to using Xive® TG implants. Xive® TG implants may only be used for their intended purpose according to the general rules for dental/surgical treatment, occupational safety and accident prevention. If the indication and use are unclear, treatment should not be carried out until all of these points have been clarified. The following instructions are not adequate to allow inexperienced clinicians to administer professional implant dentistry. We recommend that an experienced user provide proper use instructions. Xive® TG implants must only be used by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning. Dentsply Implants Manufacturing GmbH is not liable for damage resulting from use outside the intended use of the product.

The following complications have been observed after using endosteal implants:

- Components in the patient's mouth can be swallowed or aspirated.
- Nerve injury
- Implant dislocation in the maxillary sinus
- Perforation of the maxillary sinus, mandibular base, lingual plate and the inferior alveolar canal
- Fracture of the jaw bone
- Necrosis due to inadequate cooling or excessive torque
- Excessive bone loss requiring intervention
- Edema
- Loosening after ingrowth/implant loss due to insufficient osseointegration
- In rare cases, implant breakage can occur under extremely unfavorable loading conditions (prosthetic overload, severe bone resorption)
- Perimplantitis

- Infections
- Post-operative bleeding
- Iatrogenic trauma
- Esthetic or functional complications
- Suture dehiscence

Possible complications with the product or application include:

- Jammed or over-tightened placement head
- Loss of implant due to an over-tightened cover screw

Precautions

The following precautions must be observed prior to and during treatment:

- Prior to each procedure it must be ensured that all necessary components, instruments and materials are available in the required quantities and are fully operational.
- Always wear protective clothing for your own safety.
- Position the patient such that the danger of aspiration of components is minimised.
- All components that are used intra-orally must be secured to prevent aspiration or swallowing.
- Products designed for single use must not be re-used. Ignoring this requirement may cause infection!
- Observe the maximum speed.
- Comply with the recommended speed during preparation with drills; in particular, avoid a reduction in speed.
- Cutting instruments should be replaced after about 20 uses or if they are damaged or dull.
- Ensure adequate cooling.
- Extensive manipulation of the implant (e.g. change of prosthetic components) should be carried out at the earliest after the osseointegration phase, unless absolutely necessary.

Adverse Reactions

- Dysesthesia or paresthesia due to nerve injury
- Necrosis due to inadequate cooling or excessive torque
- Excessive bone loss requiring intervention
- Hyperplasia of the oral mucosa
- Esthetic or functional complications
- Edema
- Hematoma
- Temporary local swelling
- Temporary anesthesia or paresthesia/impaired mastication
- Pain

MRI Safety Information

All Xive® devices, which remain in the patient's body, have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the Xive® device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Application

Preoperative prosthetic planning

During preoperative planning, it must be established that there is adequate height and width of bone to allow placement of the implant. The width of the oral and vestibular bone lamella should be at least 1.5 mm. Appropriate radiographic preoperative analysis is necessary to determine the position of important anatomical structures. Augmented areas must show complete and mechanically stable regeneration prior to receptor site preparation. Using Friadent® Select components improves surgical planning.

Preparation of the implant site

Adequate bone quantity must be available.

There are 2 types of instruments:

- Type 1 with shaft in accordance with ISO 1797-1
- Type 2 with shaft in accordance with ISO 1797-1 and hexagon clamping system (HXSS)

Instruments with HXSS used with the appropriate contra-angle offer improved transmission.

CE/CE 0123

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice

Not all products are available in all countries.

Status of information:
IFU 2611-en_RU Rev 023 A 2022-11

**Dentsply
Sirona**
Implants

Instructions for Use

Checking the depth: there are two types of depth markings on the drills.

- Type 1 has laser ring markings only.
- Type 2 has grooves and laser band markings.

ATTENTION: The laser ring markings and grooves correspond to the implant lengths. The laser band markings are for orientation purposes. In addition, the implant site depth can be checked with a depth gauge.

The Xive® depth stops ensure that the correct depth is reached. Pull the depth stop, with the color-coding towards the top, over the drill tip. Additional instruments are not required for depth stop removal.

Implant diameter	Drill sequence (Drill without designation = Twist Drill Anterior or Posterior) The round drill can be used to smoothen the ridge or to mark the cortical bone.	Rec. drill speed
D 3.4 mm	D2.0 → D3.0 → D3.4 → Twist Drill Crestal D3.4	800 rpm
D 3.8 mm	D2.0 → D3.0 → D3.4 → D3.8 → Twist Drill Crestal D3.8	800 rpm
D 4.5 mm	D2.0 → D3.0 → D3.4 → D3.8 → D4.5 → Twist Drill Crestal D4.5	800 rpm

ATTENTION: The Twist Drill Crestal must always be used. The only exception is missing cortical bone (see the following information).

Final preparation is accomplished with a crestal twist drill at a recommended speed of 800 rpm. The depth is dependent on the bone quality/density:

- 6 mm preparation in type D I bone.
- 2 mm preparation in type D IV bone.

The depth of the Twist Drill Crestal should not exceed 4 mm for implants 8 mm in length. The Friadent® Drill extension can be used in case of limited space. It always has to be used with the long internal irrigation tube.

Ensure there is sufficient internal and external irrigation when drilling.

In case of very dense cortical bone (type D I), the Xive® Tap must be used at a speed of 15 rpm after implant site preparation with the Twist Drill Crestal. If the thread region of the tap is countersunk, then the site preparation depth is 6 mm. To use the tap without an integrated contra-angle, the implant driver for screw implants D 3.4 and TempBase must be used (with O-ring to improve TempBase seat. Replace the O-ring if the clamping effect decreases).

Implant placement

Use the Friadent® implant driver without ISO shaft with the Friadent® ratchet.

Use the Friadent® implant driver with ISO shaft without HXSS with the contra-angle.

Use the Xive® implant driver with ISO shaft to drive in Xive® implants with the corresponding ratchet insert and Xive® ratchet without contra-angle.

Use the Xive® extension instrument without internal cooling to extend instruments as required. **Do not** use as a drill extension! Align the flattened shaft of the instrument with the dot markings on the extension instrument/ratchet insert and snap into position.

1. Open the box and remove the non-sterile outer blister.
2. Open the outer blister and hand the sterile inner blister to the dentist or drop it onto the surgical tray.
3. Remove the implant holder (shuttle) from the inner blister.
4. Place the Xive® implant driver for implants and TempBase D 3.4 into the internal hex of the placement head checking that the instrument is firmly seated in place.
5. Gently bend open the shuttle and remove the implant. Avoid touching the osseous threaded section.
6. Place the implant into the prepared site and screw it into place. Ensure that the entire threaded section is in bone.
7. Remove the implant driver from the placement head.

NOTE: The placement head must not be used in prosthetic treatment.

8. Loosen the retaining screw with the hex driver 0.9 mm and remove the placement head from the implant.

ATTENTION: If during seating the torque exceeds 50 Ncm the placement head must be removed and the implant has to be screwed in directly with the implant driver for Xive® TG implants.

Under normal soft tissue conditions, the implant is screwed in place so that the apical section of the gingival collar is 0.5 mm below the crest of the bone. Thick gingival tissue should be reduced. After implant insertion, suture the soft tissue tightly around the gingival collar.

Impressions and Closure

Impressions and closure for bar-supported immediate dentures

Transfer Technique: Remove the cover screw and mount the impression coping on the implant using the 0.9 mm hex driver. Block out the hex opening in the coping with wax. Make the impression with a custom tray or use a relieved immediate denture. After the impression is made, remove the impression coping and mount it on a Xive® TG implant analog.

PickUp Technique: Mount the PickUp impression coping on the implant with a long PickUp guide pin using the 1.22 mm hex driver. A custom tray or a perforated denture with an opening for the guide pin must be used. The impression coping remains in the impression. A Xive® TG implant analog is mounted on the impression coping and the guide pin is tightened. Until the prosthesis is delivered, seal the top of the implant with the cover screw using the 0.9 mm hex driver. Make sure no soft tissue is trapped between the implant and the cover screw.

Impressions and closure during transgingival healing

During the 3-4 month healing phase, the top of the implant is protected with the cover screw. Relieve the removable provisional restoration around the implant and relapse as necessary.

Transfer Technique: Remove the cover screw and mount the impression coping to the implant using the 0.9 mm hex driver. Make the impression with a custom tray. After the impression is made, remove the impression coping and mount it on a Xive® TG implant analog.

PickUp Technique: Mount the PickUp impression coping on the implant with a long PickUp guide pin using the 1.22 mm hex driver. A custom tray or a perforated denture with an opening for the guide pin must be used. The impression coping remains in the impression. A Xive® TG implant analog is mounted on the impression coping and the guide pin is tightened. Until the prosthesis is delivered, seal the top of the implant with the cover screw using the 0.9 mm hex driver. Make sure no soft tissue is trapped between the implant and the cover screw.

Prosthetic restoration fabrication

Detailed information is available in the Xive® manual.

Delivery of the functional prosthetic restoration

Functional immediate loading is only possible if primary stability of all Xive® TG implants is ensured. Functional immediate loading may only be performed on single tooth replacement, on primarily splinted or on implants secondarily connected by stable rigid components. For bar-supported dentures the implants should be splinted with a passively fitting bar within 48 h up to a maximum of 5 days. The screws should be tightened to 14 Ncm. If the bar does not seat passively on the implants, it should be sectioned and luted together until a passive fit is obtained. Check for loose screws regularly during the osseointegration phase. The screws should be tightened to a torque of 24 Ncm after integration. Loading the implants with a removable prosthesis should be avoided if an immediate restoration is not possible.

Delivery of the prosthetic restoration after the healing phase

The prosthesis is mounted onto the implants and the abutment screws are tightened to a torque of 24 Ncm. Passive fit and occlusion must be checked and verified.

Post-operative care

Instruct the patient regarding the need for regular oral hygiene. For the first 7-10 days post-op, oral hygiene can be supplemented by a suitable anti-microbial rinse. Mechanical loading of the implant site should be avoided during the post-operative period. A relined provisional restoration can be worn for aesthetic purposes.

Delivery – Storage – Sterilisation

Xive® TG plus implants are manufactured from commercially pure grade 4 titanium under GMP standards. The implants are available in various lengths and diameters. For ease of identification, the various implant diameters are colour-coded on the respective packaging (D3.4 mm = gray, D3.8 mm = yellow, D4.5 mm = blue). For type and size see product label.

Xive® TG plus implants are gamma sterilised and are intended for **single use only**. **Xive® TG plus implants must not be re-sterilised.**

Do not use the Xive® TG plus implants after expiration of the sterilisation date or if the seal of the external blister is damaged or missing, or if the blister packaging is damaged. Store Xive® TG plus implants at room temperature +15 to +30 °C in the original packaging under conditions normal for dental practices. Do not expose to direct sunlight.

IMPORTANT: Prior to clinical use, all components and instruments delivered non-sterile must be cleaned, disinfected and sterilized according to a validated method. Instruments that can be disassembled must be disassembled for cleaning and may be assembled again for sterilization only in a dry condition.

Manual pre-cleaning

Coarse soiling must be removed immediately after use (soft nylon brush until visibly clean).

Ultrasonic cleaning (optional for very heavy soiling)

The action times and concentrations specified by the manufacturer of the cleaning agent must be observed.

Mechanical processing

Manual pre-cleaning is necessary to ensure efficient mechanical cleaning.

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice

Not all products are available in all countries.

Status of information
IFU 2611-en_RU Rev 023.A 2022-11

Instructions for Use

Suitable washer disinfectors as well as cleaning and neutralization agents are to be used for mechanical cleaning. The instructions from the manufacturer of the washer disinfectant must be followed. Cleaning and neutralization agents are to be dosed and used in accordance with the manufacturer's instructions. A cleaning program with thermal disinfection ($A_0 \geq 3000$) is recommended. Either demineralized water or water which satisfies this level of purity is recommended for cleaning.

The mechanical processing method was validated with the following materials and devices:

- Washer disinfectant, type G 7836 CD, Miele & Cie. KG, Gütersloh
- Cleaning program VARIO-TD
 - Pre-rinsing with cold water 4 min
 - Cleaning 10 min from 45-55 °C
 - Neutralization 6 min
 - Intermediate rinsing with cold water 3 min
 - Thermal disinfection 5 min at 90 °C ($A_0 \geq 3000$)
 - Drying 30 min at 105 °C
- Cleaning agent Neodisher MediClean (0.5 %) / neutralization agent Neodisher Z (0.1 %), Dr. Weigert, Hamburg.

For detailed information on manual and mechanical processing, please refer to the Dentsply Implants Instructions for Sterilization – Instrument Care.

Sterilization

Sterilization is to be performed corresponding to the following scheme:

- Preparation for sterilization:** Place components in a sterilization pouch which is legally marketed for use with the recommended sterilization parameters. Every sterilization package must have a sterilization indicator and sterilization date.

- Sterilization:**

Method	Cycle	Temperature	Exposure Time*	Dry Time
Steam	Dynamic air removal (prevacuum)	132 °C	4 Min.	20 Min.
		135 °C	3 Min.	
Steam	Gravity displacement	121 °C	30 Min.	20 Min.

*Minimum exposure times, the operating times are longer and may vary depending on the device.

- Storage:** Store the sterilized components dry and dust-free at room temperature.

Transportation Condition

Products should be transported and carried carefully in a manner that corresponds to the special transport precautions for each medical device's nature at a temperature range from -30 to +60 °C, with a humidity up to 80 %.

Terms of Use

Dental clinics and laboratories. May only be used by dentists/physicians who have had the appropriate education and training.

Application Area and Operating Conditions

Stomatology. Dental implantology.

Disposal of Products

For disposal of the packaging and components comply with the currently applicable national waste disposal regulations in Russian Federation: Implants and cover screws for implants after use - as waste class B (epidemiologically hazardous waste).

Templates, packaging of products, as well as products with an expired shelf life and not in use - according to waste class A (epidemiologically safe waste, similar with MSW).

Influence on the environment

Xive® implant system products with accessories, during storage, operation and disposal do not adversely affect the human body and the environment.

Shelf Life

Unsterile: 10 years in its original packaging
Sterile: 5 years in its original sterile packaging.
Expiration date indicated on the labeling.

Implant System Warranty

For full implant systems warranty information, refer to the detailed terms and conditions in the Implant Systems Warranty leaflet.

Product Availability

ATTENTION: Not all products described in this IFU may be regulatory cleared/released/licensed in all markets. See current product catalog and brochures for product assortment and representative product photos or contact the local Dentsply Sirona sales office for product availability.

Copyright and Trademark

All rights reserved. No part of these instructions for use covered by the copyright hereon may be reproduced or used in any form or by any means

graphic, electronic or mechanical, including photocopying, recording, taping, or information storage and retrieval systems without the permission of Dentsply Implants Manufacturing GmbH.

Xive® and Friadent® are registered trademarks of Dentsply Implants Manufacturing GmbH.

Authorized Representative on the Territory of the Russian Federation

"SIMKO BIOTECH" LLC
Olkhovskaya str. 45, bld. 1, office 4, Moscow, 105066, Russia
Tel. +7 (495) 737-80-03

Valid for registration (Russia)

Xive® dental implants with accessories

Explanation of Symbols

Symbol	Explanation
	Manufacturer
	Date of manufacture
	Use-by date
	Batch Code
	Catalog Number
	Sterilized using irradiation
	Do not resterilize
	Non-Sterile
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Do not reuse
	Consult Instructions for Use
	European conformity (CE) mark
	Prescription only: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician.
	Double sterile barrier system
	Medical device

CE/CE 0123

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.

Not all products are available in all countries.

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