Lynne N. McCormack

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MD classification

Pump and shaver for arthroscopic interventions medical device.

Health country	Authority,	Class	
FDA, US		II	
European Commission, EU	Commis-	IIa	
Ministry of Health, Russia	Health,	IIa	

INTENDED USE MEDICAL DEVICE

The FMS VUE II Fluid Management and Tissue Debridement System is intended to provide controlled fluid distention and suction, controlled cutting, burring, shaving and abrading of bone and tissue for use in gery of the shoulder, knee, ankle, elbow, wrist and hip joints

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INDICATIONS AND CONTRAINDICATIONS FOR MD USE AND POSSIBLE ADVERSE EFFECTS

Indications

The FMS VUE II Fluid Management and Tissue Debridement System is intended to provide controlled fluid distention and suction, controlled cutting, burring, shaving and abrading of bone and tissue for use in arthroscopic gery of the shoulder, knee, ankle, elbow, wrist and hip joints.

Contraindications

The FMS VUE II System is contraindicated in any non procedure and is for use with saline irrigation only. The System is not appropriate for patients for whom an arthroscopic procedure is contraindicated for any reason.

- DO NOT use the FMS VUE II in instances where capsular integrity is suspect.
- DO NOT use the FMS VUE II with a gas distension medium. Use only with sterile irrigation solution for distension of the operative site.
- DO NOT use the FMS VUE II to administer drugs, IV substitutes

Adverse effects

The following list of adverse effects may be associated with the use of this device:

- Electrical Shock;
- Infection;
- Soft Tissue Irritation;
- Burns (device overheating);
- Surgical Delay;

Adverse effects associated with the

sues, and or infection.

Local or systemic tissue reaction including allergic reaction.

Excess fluid extravasation may lead to edema near the surgical site.

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INFORMATION ON POTENTIAL USERS OF MD AND DIRECTIONS FOR USE

Directions for use

Operator should wear medical gloves when using the medical device

Unpacking

To unpack the pump:

1 Place hands in the open areas between the two pieces of packaging foam on both sides of the console until you can securely grasp the unders pump.

2. Carefully remove the pump from the box and place it on a level and sturdy surface.

3. Verify that all items on the packing list match the products you have received.

4. Report damage immediately to a DePuy Synthes Sales Representative or call DePuy Synthes Customer Service in the USA at +1 all packing materials to use in the event you must return product(s) to DePuy Synthes. Outside the United States, contact your local affiliate.

5. Recycle all FMS VUE II packaging materials and antistatic bag.

Pump Placement

Position the FMS VUE II pump outside the sterile field, and where required depending on operator requirements. Ensure that there is enough space for Circulating Nurses to access the pump controls. Position the pump so that access to power main is not obstructed. There are two options for setting up pump elevation:

- Place the pump (fill line of Fill Chamber) at the same level as the patient joint.
- Place the pump at a convenient cart elevation and adjust the joint elevation from the front panel at startup.

Pump Assembly and System Check

To assemble the FMS VUE II pump:

1. Connect the power cord to the power source.
2. Connect the other end of the power cord to the rear panel of the pump.
3. Connect the accessories to the pump (Figure 1).

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4. Turn the pump ON by pressing “|” on the power switch on the front of the pump (Figure 1). The software version appears on the display followed by the initial settings.


Fig 1. Connect medical devices to the pump
Connecting and Disconnecting Pump Tubing
To connect the Inflow Tubing:
Use aseptic technique.

Fig 2. Remove Tube from Package

1. Remove the Inflow tubing from the package and close three (3) white clamps on the Inflow tubing (Fig 2.).

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Fig 3. Position the Green Autolock on the Pump

2. Lift the door and position the green autolock on the tube in the green autolock slot on the pump (Fig 3.).
3. Place and center the tube around the roller.


Fig 4. Place and Center the Inflow Tube on the Pump

4. Pull gently on the tube to engage the tube fitting lock while maintaining tension on the tube (Fig. 4).

CAUTION: Make sure the transparent tubing connection is securely positioned in the auto-lock to avoid damage to the tubing. Make sure the tubing is centered in the pump roller groove.

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
Fig 5. Position Fill Chamber

CAUTION: Make sure the Quick Connect port is open before inserting the Quick Connect tube. To open, press down on the lever at the top of the port (4).

5. Connect the pressure sensor Quick Connect tube (1) to the fill chamber port on the FMS VUE II pump (Fig 5.). The connector clicks when it is locked.
6. Position the Fill Chamber in the Fill Chamber bracket (2).
7. Close the pump door.

WARNING: Do not invert fill chamber at any time.

Position the chamber in the bracket below the pressure sensor connection. Make sure the white clamp on the pressure chamber is closed.

CAUTION: When placing the fill chamber in the bracket, make sure that the fill chamber symbol (3) is facing toward the user and tubing is not twisted. Improper positioning could cause the tube to twist and obstruct the flow of fluid.

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Fig 6. Insert Irrigation Spikes into the Irrigation Bags

8. Remove the caps from the irrigation spikes and insert the spikes into the irrigation bags (Fig 6.).
9. Open the two (2) clamps on the Inflow tubes.
10. Press the RUN/STOP pad to turn on the pump

Fig 7. Press FILL Chamber Pad to Prime Pump

11. To prime the pump, press and hold the FILL Chamber Pad until the chamber is 1/3 full. Refer to the mark on the fill chamber shown in Fig 7.

WARNING: The pressure chamber must remain vertical during the procedure. Make sure that the pump is positioned at the same level as the patient.

To connect the Intermediary tubing:
Use aseptic technique.

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Fig 8. Remove Blue Caps from Intermediary and Inflow Tubing and Connect Intermediary Tube

1. **STERILE OPERATOR:** Make sure that the luer ened and pass the blue end to the non
2. **NON-STERILE OPERATOR:** Remove the blue caps from the intermediary tube and Inflow tubing and immediately connect the Inflow Tubing luer-lock to the Intermediary Tube luer

Fig 9. Open Clamp under Pressure Chamber

3. Prime the Intermediary Tube chamber (Fig 9.)

To install the Outflow Tube:

Use aseptic technique.

Fig 10. Position Outflow Tubing on Pump Roller

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1. Remove the Outflow Tubing from the package.
2. **STERILE OPERATOR** – Pass the waste end of the suction non-sterile operator (circulating nurse).
3. Lift the right pump door and position the orange autolock in the orange autolock slot on the pump (Fig 10).

Fig 11. Close Pinch Valve

4. Place and center the tube around the roller.
5. Pull gently on the tube to engage the tube fitting opposite the auto while maintaining tension on the tube. **CAUTION:** Make sure the transparent tubing connection is securely positioned in the auto avoid damage to the tubing. Make sure tubing is centered in the p roller groove.
6. Close the right pump door.
7. Rotate the Pinch Valve to the vertical position to open.
8. Insert both tubes into the Pinch Valve slots.
9. Rotate the Pinch Valve to horizontal position to close (Fig 11.).

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Fig 12. Connect Waste

10. Insert the waste end of the outflow tubing into waste collection (Fig 12.).

Connecting the Handpiece – STERILE ZONE

Sterile operator: Use this procedure to connect the handpiece in the sterile zone.

Use aseptic technique.

Insert the tube labelled “SHAVER” onto the h

Connecting the Cannula – STERILE ZONE

Sterile operator: Use this procedure to connect the Cannula in the sterile zone.

Use aseptic technique.

Insert the tube labelled “CANNULA” onto the Cannula suction port.

Tubing Connections for SOLO Mode

To operate the pump in SOLO Mode (irrigation and handpiece control only):

Before operation:

1. Make sure the pump tubing is connected.
2. Connect the Inflow Tubing.
3. Connect the Intermediary Tubing.

End of Operation

To end operation:

Fig 13. Close the Clamp under the Pressure Chamber

1. Close the clamp under the pressure chamber (Figure 13). When another surgery is scheduled for the same operating day, leave the Inflow Tubing.

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Fig 14. Disconnect the Luer-lock Below One

2. Disconnect the blue luer-lock of the Inflow tube below the One Valve (Fig. 14).

Fig 15. Press on Blue FLOW+ button (FLOW+)

3. Press the blue FLOW+ button on the 5 patient joint (optional) (Fig. 15).
4. Rotate the Pinch Valve to the vertical position to release the tube.
5. Remove the orange suction tubing and discard according to hospital regulations.

Set Up for Next Operation

Fig 16. Remove Intermediary Tube

1. STERILE OPERATOR: Check that the luer and pass the blue end to a non
2. NON-STERILE OPERATOR: Disconnect the used One and immediately connect the new Intermediary Tube with the One Way Valve.

End of Operating Day

At the end of the operating day, remove and discard all tubing according to approved hospital procedures. Turn off the pump.

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Potential users of MD

This equipment may be operated only by arthroscopic surgeons. The surgeon using this device must be trained in arthroscopic surgical procedures, be aware of the risks associated with those procedures, and have current knowledge of technological advances in surgical products and techniques.

Use/operation conditions of MD

Operation conditions	
Humidity	
Temperature	
Pressure	

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Description of the main MI nodes

The FMS VUE II pump controls irrigation and suction independently which can help minimize the

risk of tissue extravasation by eliminating the need to increase fluid pressure in the joint. The pump may be operated in either DUO or SOLO mode:

- DUO Mode (irrigation, handpiece control, and suction). When operating in DUO mode, the pump provides irrigation, handpiece control and suction.
- SOLO Mode (irrigation and handpiece control). When operating in SOLO mode, the pump provides irrigation and handpiece control. SOLO mode is activated by pressing the SOLO pad on the control panel


Fig.17. MD connection diagram

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Element	
1	Transparent safety doors Doors must be closed for pump operation.
2	Irrigation roller pump Provides fluid inflow.
3	Pressure adjustment pads (mmHg) Controls the pressure level in increments of 5 mmHg.
4	Pressure display (mmHg) Displays both the preset pressure when the base pressure is being set, and the dynamic pressure following the release of the pressure adjustment pads.

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5	FLOW+ Suction setting	Sets the flow rate through the cannula.
6	Shaver speed display (RPM)	Indicates direction of rotation.
7	Shaver Speed Selection pads (RPM)	Increases and decreases handpiece speed.
8	Suction roller pump	Provides fluid outflow.
9	Suction Auto-locking slot	Locks the tubing in place around the suction roller pump.
10	RFID antenna	Detects correct pump tubing.
11	Pinch valve	Moves automatically when the handpiece rotates and selects the appropriate suction tube.
12	Shaver Suction setting	Increases and decreases handpiece suction.
13	SOLO pad	Press to
14	Power Switch	Provides power to the pump.
15	Foot Pedal connection	Connection port for 5
16	RUN/STOP pad	Starts and stops the pump.
17	Remote hand control connection	Connection port for the remote hand control.
18	Handpiece interface connection	Connection port for handpiece or FMS CONNECT Interface Cable (separately, representative for availability)
19	Fill Chamber pad	Press to fill the pressure chamber.

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Fuse ratings: 5 Amps, 250 Volts, **Breaking Capacity:** Н

Conforms to standard IEC 60601-1 Medical Device Electrical Safety and IEC 60601-1-2 (Electromagnetic Compatibility)

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The FMS VUE II pump is a Class 1 device according to IEC 606061-1. (No electrical contact with patient). BF isolation applies to both pumps and shaver.

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Joint Pressure Range: 20–150mmHg (increments of 5).

When pressure is stabilized, levels will be within 20% of the set point.

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Cannula flow rate: 100–600ml/min.

Pre-selected Precision ($\pm 20\%$) of maximum flow rate.

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Handpiece Utilization Cycle: 1 minute every 6 minutes

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Aspiration (Outflow)

An integrated pinch valve on the aspiration side of the pump has a two way valve that allows two different outflow modes (suction modes): shaver suction or cannula suction. This feature is available on the predicate device DUO+. In the DUO mode, the pinch valve is always cycling between the two outflow modes to maintain intra pressure. When the shaver is activated the pinch valve opens allowing for the removal of fluid and debris from the joint through the shaver. When the shaver is not activated the pinch valve allows for removal of fluid and debris from the joint via the outflow cannula.

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Name	Material of MD
Console FMS VUE II	AlMg3 (powder coated) 6060 Al (powder coated) 303 Stainless Steel 6082 Al (powder coated) Polyphenylsulfone Radel R 5100 GY 1137 Acetal

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


Definitions of symbols used on the label

	Consult instructions for use
	Refer to instruction manual/booklet
	Catalogue number
	Serial number
	Quantity
	Date of manufacture

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


NONSTERILE	Nonsterile
	Type B. Conforms to standard IEC 60601
	Type BF. Conforms to standard IEC 60601
	Per Directive 2002/96/EC, product must be collected separately. Do not dispose of as unsorted municipal waste. Contact your DePuy Mitek Sales Representative for disposal information.
IP40	Ingress Protection Rating
Made in	Made in
Rx Only	Federal Law (United States of America) restricts this Device to sale by, or on the order of a physician.
	Manufacturer

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	EC Representative
	Intertek Semko evaluated the FMS VUE System and its accessories per standards IEC-60601 2005 and IEC-60601 pump and accessories were found to be compliant.
	European Conformity
	Temperature
	Humidity
	Atmospheric pressure
	Caution
IP58	Ingress Protection Rating

Device service life and storage

There is no shelf life, as the FMS VUE II

The Lifetime of these products is determined by their service warranty which is 1 year.

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The shipping box contains:

- FMS VUE II
- Power Cord
- Operator's Manual

Registered medical devices have their own product codes and are shipped separately.

LIST OF STANDARDS APPLICABLE TO A

Device meet requirements EN ISO 13485 Medical device system – Requirements for regulatory purposes. A full list of applicable standards is available upon request.

LIST OF MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

Not applicable. This medical device does not contain materials of animal and / or human origin.

12. INFORMATION ON THE TESTS CARRIED OUT, TEST REPORTS, ANALYSIS OF THE OBTAINED DATA

FMS VUE II

Guidance and Manufacturer's Declaration

— Emissions

The FMS VUE II is intended for use in the electromagnetic environment specified below. The customer or user of the FMS VUE II should ensure that it is used in such an environment.

Emissions Test	Compliance	
RF Emissions CISPR 11	Group 1	

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RF Emissions CISPR 11	Class A	
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Complies	
The FMS VUE II is intended for use in the electromagnetic environment specified below. The customer or should ensure that it is used in such an environment.		

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	30% Dip for: - 25 Cy- cles at 50Hz - 30 Cycles at 60Hz 100% Dip for: - 250 Cycles at 50Hz - 300 Cycles at 60Hz	30% Dip for: - cles at 50Hz - 30 Cycles at 60Hz 100% Dip for: Cycles at 50Hz Cycles at 60Hzs		
Power Fre- quency 50/60Hz Magnetic Field IEC 61000-4-8	30 A/m	30 A/m		
The FMS VUE II is intended for use in the electromagnetic environment specified below. The customer or user of the FMS VUE II should ensure that it is used in such an environment.				
Immunity Test	IEC 60601 Test Level	Compliance Level		
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6Vrms in ISM bands between	(V1) = 3 Vrms (V2) = 6 Vrms		

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mended Separations Distances for FMS-

The FMS-VUE II is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the FMS help prevent

electromagnetic interference by maintaining a minimum distance between and mobile RF Communications Equipment and the FMS

below, according to the maximum output power of the communications equipment.

Max Output Power (Watts)	Separation (m) 150kHz to 80MHz $D = (3,5/V1)(\sqrt{P})$	Separation (m) for ISM Bands 150kHz to 80MHz $D = (3,5/V2)(\sqrt{P})$	Separation (m) 80 to 800MHz $D = (3,5/E1)(\sqrt{P})$	
0.01	0.11667	0.05833	0.11667	
0.1	0.36894	0.18447	0.36894	
1	1.1667	0.58335	1.1667	
10	3.6894	1.8447	3.6894	
100	11.667	5.8335	11.667	

RFID Transmitter/Receiver Specifications

Frequency	
RF Output Power (EIRP)	Field Strength = 60.73 dBμV/m @ 10meters
Modulation	ASK (Amplitude Shift Keying)
Тип антенны	Spiral Loop PCB Antenna (Gain = 0 dBi)

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Electrical Safety

- Electrical safety testing should be performed by a biomedical engineer or other qualified person.

- Risk of electrical shock: Do not remove the pump housing. Refer servicing to authorized Repair Service Center.

- Do not connect the device to a power source that is not properly earthed (grounded). If the power source is not grounded, use the equipotential connector.
- Disconnect the device from the main power source when cleaning, servicing or inspecting.

- Inspect all equipment and cables periodically for wear. If damage is noted, replace or return to Repair Service Center.

- Avoid fluid contact with the system and its electrical connectors.

- To avoid risk of fire, replace fuses with same type and rating

- Do not restrict access to the power switch. Ensure that the pump is positioned so that access to power main is not obstructed. To remove all mains power from the unit, disconnect the power cord from the appliance inlet.

- The FMS VUE II pump power supply isolation between the power supply and the power main. Ensure that the placement of the equipment does not interfere with access to the power supply cord.

Operating Warnings

- This equipment may be operated only by using this device must be trained in arthroscopic surgical procedures, be aware

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of the risks associated with those procedures, and have current knowledge of technological advances in surgical products and techniques.

- While using the tubing with one-way valve, if additional joint pressure is created by the surgeon during joint manipulation, the pressure displayed at the pump will not reflect the pressure change, which could possibly lead to extravasation.

- Do not use the pump adjacent to or stacked with other electrical equipment. If adjacent or stacked use is necessary, observe both the equipment and pump to verify normal operation in the configuration in which it will be used.

- Portable and mobile RF communications

Electrical Equipment.

- Do not use in the presence of a flammable anesthetic mixture with air, or in the presence of a flammable anesthetic mixture with oxygen or nitrous oxide.

- Factory settings should only be used as guide

for settings pertaining to the surgical procedure.

- Failure to follow all applicable instructions may result in serious surgical consequences to the patient.

- Use the FMS CONNECT Interface Cable to connect handpieces o

those manufactured by DePuy Synthes (

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supported handpieces, refer to FMS CONNECT Cable for FMS VUE II Instructions for Use. (P/N: IFU-111105).

- Any modifications to the system must be authorized by DePuy Synthes. Modifying the device in a manner other than specified by DePuy Synthes may result in injury. Contact Customer Service at 1

requests. Outside the United States, contact your local affiliate.

- The FMS VUE II pump is equipped with a terminal for the connection of a Potential Equalization Connector as required by Electrical Safety Standard IEC 60601-1 (поставляется отдельно, обратитесь

чтобы узнать о наличии). To install the connector, connect one end to the terminal on the back of the pump. Attach the other end to the bus bar of the electrical installation. The Potential Equalization Connector is provided to prevent electric shock hazards by maintaining a constant potential (voltage) to all exposed metal parts so a person cannot simultaneously make contact with two different potentials. Use of the potential equalization connector should be determined by the facility where it is in distribution system.

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- This device has been tested and found in compliance with limits of medical devices according to IEC 60601-1-2 standard. These limits are designed to provide reasonable protection in a type interference when devices are located in close proximity. If the pump is causing harmful interference, turn the main power off, then on again to determine the cause of the interference, relocate the system or separate the pump. If you cannot resolve the problem contact DePuy Customer Service at 1-800-382-4682. Outside the United States, contact your local affiliate.
- Accessories labeled "REUSABLE" must be processed only according to the recommended procedures. Refer to the applicable Instructions for Use.
- Do not wrap handpiece, foot pedal or pump power cables around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire and/or injury to patient.
- The use of accessories, transducers and cables other than those specified, with the exception of transducers and cable sold by the Manufacturer of the ME Equipment or ME System, may result in increased emissions or decreased immunity of the ME Equipment or ME System.
- No parts of the FMS VUE II system can be serviced or maintained while the pump is in use.
- The product shall be used in accordance with instructions and standard aseptic techniques, otherwise cross-contamination may occur.

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INFORMATION ABOUT MEDICINES AND (OR) SUBSTANCES INCLUDED IN MD

Not applicable. Devices do not interact with medicinal and / or pharmaceutical substances.

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INFORMATION ABOUT DESIGNING, ENGINEERING, AND VALIDATION OF SOFTWARE USED

Software development plan for this device has been developed in accordance with the overall quality assurance effort as defined in the FMS VUE II Software Development Plan (103094123).
The software version will be released as 103271866 rev. A, in pr according to CO-103346495.
Software safety class - B.

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This section provides functional tests that comply with the IEC 60601 rective. The following tests are provided:

- Functional test of active elements
- Functional test of pinch valve

This test is intended for use by readers who are familiar with the material in this Manual or biomedical engineers/technicians. All output, safety, product conformity, acceptance testing, commissioning acceptance testing, must be performed by personnel trained and authorized by DePuy Mitek. There are no user servicable parts inside the pump.

WARNING: Do not modify this equipment without authorization of the manufacturer.

WARNING: Electrical shock hazard. Do not open pump housing.

Materials required for test: Service Test Tubing Ref.284511

- **Functional test of active elements**

Verify Motors and Electro-valve

1. Turn the pump ON by pressing “|” on the power switch on the front of the pump.
2. Install the QC test tubing Ref.284511 as described in Chapter 5.
3. Press RUN/STOP to activate the pump and verify that the following components operate correctly:
 - Irrigation Motor: With the left door down and the right door raised, verify that the motor rotates.
 - Suction Motor: With the right door down and the left the motor rotates.
 - Electrovalve: With both doors raised press and hold the FILL Chamber pad. The electrovalve should engage for 5 seconds.

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• Functional Test of Pinch Valve

Rotate the pinch valve vertically and verify that both si

REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD Factory Settings

Factory settings are based on experience and should only be used as guidelines. The surgeon is responsible for settings depending on the surgical procedure performed. **WARNING!** Risk of electrical shock: do not remove the pump housing. Refer servicing to your DePuy Synthes Sales Representative, or contact DePuy Synthes Customer Service in the USA at +1

United States, contact your local affiliate. Do not connect the device to a power source that is not properly earthed (grounded). Dis main power source when cleaning, servicing, or inspecting. Inspect equipment and cables periodically for wear. Avoid fluid contact with the FMS VUE II, and its electrical connectors. To avoid risk of fire, replace fuses with sam rating

Cleaning the FMS VUE II Pump

At the end of the operating day, switch the FMS VUE II off, unplug the power cord and discard all tubing. Wipe off the instrument's sides, front panel and doors with a soft cloth that has been slightly moiste

Alcohol) or equivalent At the end of the operating day, switch the FMS VUE II off, unplug the power cord and discard all tubing. Wipe off the instrument's sides, front panel and doors with a soft cloth that has been slightly moi

IPA (Isopropyl Alcohol) or equivalent

DO NOT IMMERSE. DO NOT AUTOCLAVE.

Cleaning

Refer to the Instructions for Use for cleaning instructions.

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The following procedures may be performed by a biomedical engineer or biomedical technician.

Replacing Fuse

The fuse door is located on the back panel below the power cord receptacle. Use a flat head screw driver to open the door. Replace the fuse with 5A Time Delay Fuse rated for 250VAC.

Physical Inspection

Check all equipment to assess its functionality and technical safety. These inspections may be performed on a less frequent basis depending on the duration and frequency of use.

To inspect:

1. Disconnect the console from the power supply.
2. Inspect all exterior components, surfaces, and corners for damage. If the unit was dropped, contact your DePuy Synthes Sales Representative or call DePuy Synthes Customer Service in the USA at +1 replacement. Outside the United States, contact your local
3. Inspect the foot pedal and handpiece for obvious signs of mechanical damage or wear.
4. Inspect the power cable, verify that the cable is intact.
5. Look for cuts on the cables and missing or bent pins on the connectors.
6. Remove the fuse holder (located below the power cord receptacle) and verify that both fuses are intact and match the current, time delay and voltage ratings per the back panel label.

Weekly inspection

- Ensure that the transparent safety door and hinges are working
- Inspect power cord and all cables for wear or cuts.
- Check all connections.
- Ensure that the foot pedal cable and handpiece interface cable connector pins are not missing or bent.
- Plug the power cord to the power connection.
- Switch the I/O switch on "I".
- Press RUN/STOP pad.
- Press on all pads. You should hear an audible tone for each pad pressed. When the pump is primed, the screen displays the actual pressure.

Annual inspection

Check all equipment annually to assess its functionality a

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Repairs

For repairs, contact Customer Service Center in the USA at
For repairs outside the United States, contact your local affiliate. Relevant service instructions are provided in DePuy Mitek work instructions. DePuy provide circuit diagrams, component part instructions, or other information that will assist service personnel authorized by the manufacturer in repairing the pump.

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POSSIBLE FAULTS AND TROUBLESHOOTING.

Problems	Potential Causes	
The pump stops. Red ALARM LED illuminates.	The HIGH PRES-SURE safety alert has activated.	

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CONDITIONS OF TRANSPORTATION AND STORAGE

Designation	Range
Temperature	-10°C
Relative humidity	10%
Barometric pressure	500-

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MD DISPOSAL RULES AND CONDITIONS

The FMS VUE contains electronic printed circuit assemblies. At the end of the useful life of the equipment disposed of in accordance with any applicable national or obsolete electronic equipment.

Environmentally Friendly Use Period (EFUP) China. The during which the toxic contained in electronic information products do not leak or mut potentially causing serious pollution to the to human life and property.

This symbol indicates that the product contains less than the maximum concentration value of all six hazardous substances.

Dispose of any system accessories relating to potentially contaminated items.

Any opened but unused device, as well as an expired device should be utilized in accordance with accepted medical practice and applicable local and national requirements.

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WARRANTY

DePuy Mitek products are guaranteed to be free from defects in materials and workmanship for the warranty period of a particular product, beginning from date of invoice.

All spare parts are used only by a qualified engineer when repairing in a service center.

Replacement of spare parts and repairs outside the service center prohibited.

Contact your DePuy Mitek representative for specific product warranty information.

DePuy Mitek offers an expedited service replacement program for its products in the United States to minimize downtime in your operating room. For authorization number or for additional information on this program, call DePuy Mitek Customer Service at +1 (800) 382

your DePuy Mitek representative. Outside the United States, affiliate.

NAME AND LEGAL ADDRESS OF MANUFACTURER AND AUTHORIZED REPRESENTATIVE

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325 Paramount Dr.,

Raynham, Massachusetts 02767, USA

**AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN
RUSSIA:**

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Refer to IFU FMS VUE II Fluid Management and Tissue Debridement System

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



Name	Material of MD
Remote control with connection cable	Silicone Medical Grade Polyester Medical Grade Stainless steel 301 PEEK

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



Warnings

- Prior to use, check the system components for damage. Check the cable integrity carefully. If there are signs of damage, do not use.
- Failure to follow all applicable instructions may result in serious surgical consequences to the patient.
- Do not use flammable agents when cleaning and disinfecting.
- Do not wrap Remote Hand Control cable around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire and/or injury to patient or surgical personnel.
- This device has been tested and found in compliance with limits of medical devices according IEC 60601-1-2 standard. These limits are designed to provide reasonable protection in a typical medical installation, against harmful interference when devices are located in close proximity. If the pump is causing harmful interference with other devices, turn the main power off, then on again to determine if it is causing the interference, relocate the system or separate it from other devices. If you cannot resolve the problem contact DePuy Customer Service at 1-800-382-4682. In the European Union, contact your local affiliate.
- Disconnect the device from the pump when

PRECAUTIONS

- Read the instructions, cautions, and warnings provided with all FMS VUE System accessories before use. Your sales representative can advise which accessories are used with the FMS VUE System.
- Electrical safety testing should be performed by a biomedical engineer or other qualified person.

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- Wash at 60°C with pH neutral cleaning agent for 3 min.
- Rinse with hot water for 20 seconds.
- Dry at 95°C for 5 min.

9. After automated cleaning, proceed to inspection instructions below.

Manual Cleaning

Clean before and after each use:

1. Place the protective cap securely onto the Remote Hand Control cable connector.
2. Carry out a first soaking bath in a pH for 15 minutes (or an equivalent decontaminant solution) and refer to decontaminant instructions for immersion time and bath temperature
3. Carry out 5 (five) to-and-fro motions with a swab paying particular attention to the buttons, the screws on the back side and the seam around the sides.
4. Check for remaining debris and remove as needed.
5. Carry out a second soaking bath in a pH CIDEXPLUS for a minimum of 15 minutes while using an ultrasound bath frequency of 50–60Hz (or an equivalent decontaminant solution) and refer to cleaning agent solution manufacturer's instructions for immersion time and bath temperature.
6. Rinse with deionized water for 30 seconds at ambient temperature.
7. Dry with medical air at less than 100°C or 212°F for 30 seconds.

Inspection after cleaning

1. Inspect before sterilization or storage to ensure the complete removal of soil from surfaces and holes.
2. If areas are difficult to inspect visually, check for blood by immersing or flushing in a 3% hydrogen peroxide solution. If bubbling is observed, blood is present. Rinse thoroughly after using hydrogen peroxide solution.
3. If soil is still present, re-clean.

Sterilization

- The Remote Hand Control is factory sealed and can be completely sterilized using a steam autoclave.
- Use a validated, properly maintained and calibrated steam sterilizer.
- Effective steam sterilization can be achieved for 4 minutes.

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Forward
Shaver direction

Oscillate
Shaver direction

Fig 2. Foot Pedal and Handpiece Functions

The possibility and ways of integrating

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LIST AND DESCRIPTION OF MD

Name	Material of MD
5- Way Foot Pedal with connection cable	Acetal Polyethylene Polycarbonate

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




DATA ABOUT LABELING OF MD AND ITS PACKAGING

	Consult instructions for use
	Refer to instruction manual/booklet
	Catalogue number
	Serial number
	
	Date of manufacture
	
	Type B. Conforms to standard IEC

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	type BF. Conforms to standard IEC
	Per Directive 2002/96/EC, product must be collected separately. Do not dispose of as unsorted municipal waste. Contact your DePuy Mitek Sales Representative for disposal information.
IP40	Ingress
Made in	
	Federal Law (United States of America) restricts this Device to sale by, or on the order of a physician.
	Manufacturer
	EC Representative
	Intertek Semko VUE System and its accessories per standards IEC 2005 and IEC pump

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- Prior to use, check the system components for damage. Check the cable integrity carefully. If there are signs of damage, do not use.
- Failure to follow all applicable instructions may result in serious surgical consequences to the patient.
- Factory settings should be used only as guidelines. The surgeon is responsible for settings pertaining to his
- Do not wrap Foot pedal power cables around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire and/or injury to patient or surgical personnel.
- This device has been tested and found in compliance with limits of medical devices according to IEC 60601-1-2 standard. These limits are designed to provide reasonable protection in a typical medical installation, against harmful interference when devices are located in close proximity. causing harmful interference with other devices, turn the main power off, then on again to determine if it is causing the interference, relocate the system or separate it from other devices. If you cannot resolve the problem contact DePuy Customer Service at 1-800-382-4682 your local affiliate.
- Disconnect the device from the pump when cleaning or inspecting

INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD
Not applicable. Devices do not interact with medicinal and / or pharmaceutical substances.

FUNCTIONAL TEST

Refer to IFU FMS VUE II Fluid Management and Tissue

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After each operation, clean the foot pedal and the power cord with a Glutaraldehyde disinfectant such as CIDEX® OPA solution manufacturer's instructions.

1. Prepare the disinfectant solution per manufacturer instructions.
2. Wipe the surface of the Foot Pedal with a cloth, swab or gauze pad soaked in disinfectant solution.
3. Let the disinfectant stay on all device areas for 19
4. Rinse thoroughly in running water.
5. Allow the Foot Pedal to drain and dry thoroughly

Caution:

DO NOT AUTOCLAVE

Do not twist the cable when storing. Do not unplug instrument by pulling on cable.

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REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD MAINTENANCE

Inspect all equipment and cables periodically for wear. If damage is noted, replace or return to Repair Service Center. Any repairs must be performed by a qualified Biomedical Engineer.

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Instruction for use
Tornado Micro II Shaver Handpiece No-Button
Tornado Micro II Shaver Handpiece Hand Controlled

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DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS OF MD AND PRINCIPLE OF OPERATION

Tornado Micro II Shaver Hand Controlled

The Tornado Micro II handpiece is an accessory device to the FMS VUE™ systems FMS VUE II. The handpiece consists of a locking and rotating device for holding shavers and burrs, and includes three buttons for controlling shaver or burr rotation direction and BLOODSTOP™ f
foot pedal also may be used to control shaver or burr rotation direction. The device provides suction through a suction port when connected to the FMS VUE II pump.

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Part	Function
1 Locking ring	Holds shaver or burr in place.
2 Suction lever	Controls suction through the suction port.
2a BLOODSTOP button	Activates the BLOODSTOP mode.
2b Handpiece rotation mode button	Press to to oscillate mode
2c ON/OFF button	Turns the handpiece on and off.
3 Serial number	Handpiece Identification.
4 Suction port	Connects to the FMS tubing.
5 Alignment slots	Aligns shaver or burr face up or face down.

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Tornado Micro II Shaver No-Button

The Tornado Micro II handpiece is an accessory device to the FMS VUE systems FMS VUE II. The handpiece consists of a locking and rotating device for holding shavers and burs. A foot pedal also may be shaver or burr rotation direction. The device provides suction through a suction port when connected to the FMS VUE II pump.

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










Dedicated button control mode (Only for FMS VUE II)

To...	Press...	
Activate/deactivate Oscillate function		
Activate/deactivate Forward/Reverse		
Change direction between Reverse/Forward	 1 sec	
Start BLOODSTOP mode	 2 sec	
Activate/deactivate override feature when the Foot Pedal is connected to the pump	 +  2 sec	
Activate/deactivate Dedicated button control mode to Default Button Control mode	 +  2 sec	

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Tornado Micro II Shaver Handpiece Hand Controlled

To...	Движение	
Increase/turn on suction		
Reduce/stop suction		
Remove protective cap at end of cable	Unscrew... 	
Attach Handpiece Cable Connector to the FMS connection port using Connector alignment arrow	Align and press 	
Attach FMS tubing to Suction port		

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BLOODSTOP

Feature The BLOODSTOP feature provides for a timed (2 minute) cycle of increased flow and pressure to stop bleeding and clear the field of view. Once the cycle is complete, the FMS VUE II will return to its original settings. When the BLOODSTOP feature is activated, the pump will have flashing brackets in the pump pressure display (


Fig. 3. Flashing brackets in the pump pressure display on FMS VUE II Dedicated button control mode (Only for FMS VUE II)

A configuration for your controls that forward/reverse button. The controls will return to default button control mode every time the handpiece is reconnected to the pump or when the pump is re-started.

Inserting and Removing Shavers & Burrs

To insert shavers and burrs into the handpiece, align the shaver/burr's locking pin with the handpiece alignment slot (face up or face down, depending on user preference) (Fig. 4) and push in until an audible "click" is heard. To remove the shaver/burr, use your thumb to the left, then pull out shaver/burr (

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Handpiece Controls when Foot Pedal is Connected

When the foot pedal is connected to the pump, the handpiece MODE button and the ON/OFF button are automatically deactivated, but the BLOODSTOP button remains functional. The Override feature allows both foot pedal and handpiece controls when pressing (BLOODSTOP Feature) and (ON/OFF) buttons for 2 seconds simultaneously (handpiece control is activated when Foot Pedal is connected (

Fig. 6 a. To override feature press (BLOODSTOP Feature) and (ON/OFF) buttons.


Fig. 6 b. Hand Control de-activated (left) Hand Control activated (right)

The possibility and ways of integrating with other medical devices

Refer to IFU FMS VUE II Fluid Management and Tissue Debridement System.

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Recommended Handpiece Utilization Cycle: 1 minute every 6 minutes

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






	Consult instructions for use
	Refer to instruction manual/booklet
REF	Catalogue number
SN	Serial number
QTY	Quantity
	Date of manufacture
NONSTERILE	Nonsterile
	Type B. Conforms to standard IEC 60601
	Type BF. Conforms to standard IEC 60601

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	Per Directive 2002/96/EC, product must be collected separately. Do not dispose of as unsorted municipal waste. Contact your DePuy Mitek Sales Representative for disposal information.
IP40	Ingress Protection Rating
Made in	Made in
	Federal Law (United States of America) restricts this Device to sale by, or on the order of a physician.
	Manufacturer
	EC Representative
	Intertek Semko evaluated the FMS VUE System and its standards IEC-60601 2005 and IEC-60601 pump and accessories were found to be compliant.
	European Conformity

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MD SAFETY AND EFFECTIVENESS REQUIREMENTS

This equipment may be operated only by personnel personnel must also be aware of the risks associated with those procedures and have current knowledge of technological advances in surgical products and techniques.

- Prior to use, check the system components for damage. Check t carefully. If there are signs of damage, do not use.
- Failure to follow all applicable instructions may result in serious surgical consequences to the patient.
- Do not wrap handpiece cable around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire and/or injury to patient or surgical personnel.
- This device has been tested and found in compliance with limits of medical devices according to IEC 60601-1-2: 2014 standard. These l provide reasonable protection in a typical medical installation, against harmful interference when devices are located in close proximity. If the pump is causing harmful interference with other devices, turn the main power off, then o determine if it is causing the interference, relocate the system or separate it from other devices. If you cannot resolve the problem, contact DePuy Synthes Customer Service at +1-(800)-382-4682 in the United States or contact your DePuy Synthes sales representative. Outside of the United States, contact your local affiliate.
- Disconnect the device from the pump when cleaning or inspecting.
- Do not insert anything other than approved DePuy Mitek shavers and burrs into the distal end of the handpiece.
- For patients with, or suspected with, Creutzfeldt CJD or other transmissible spongiform encephalopathy (TSE) and related infections, it is recommended to treat the patient using single Safely dispose of all devices used in accordance with local procedures and guidelines.
- Failure to follow the Cleaning and Sterilization instructions may result in an infection.

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INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD

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DESCRIPTION OF REPROCESSING AND STERILIZATION

The Tornado Micro II handpiece

Sterilization

Package trays/instruments with a barrier wrap procedures, using standardized wrapping techniques such as those described in ANSI/AAMI ST79. In the United States, use an FDA

Sterilize in a validated sterilizer according to the guidelines bel

- Steam (moist heat) sterilization is recommended in an approved, pre vacuum (force air removal) cycle. The steam sterilizer should be validated to the requirements of any local standards and guidance such as EN 285 or ANSI/AAMI ST8. Ensure that the steam sterilizer cycle is chosen that is designed to remove air from porous loads and/or devices with lumens in accordance to manufacturer's instructions.
- Effective steam sterilization can be achieved using the following cycles:

Cycle Type	Minimum Exposure Time	Minimum Sterilization Temperature	

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Pre-vacuum	4 minutes	132°C (270°F)	
Pre-vacuum (outside of US)	3 minutes	134°C (273°F)	

The total combined sterilization tray system should not exceed 22 lbs. (weight limit). Do not stack tray systems for sterilization.

Additional Information:

- Cleaning Agent Information: DePuy Orthopaedics, Inc. used the following cleaning agents during validation of these reprocessing recommendations. These cleaning agents are not listed in preference to other available cleaning agents which may perform satisfactorily: Valsure® and Enzol®.
- The cleaning and sterilization information is provided in accordance with ANSI/AAMI ST81 and EN ISO 17664.

REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD HANDPIECE CLEANING

The Tornado Micro II handpiece shall be cleaned using either the Manual or the Automated instructions and sterilized before the first use and in Clean instruments as soon as possible after use. If cleaning must be delayed, immerse instruments in a compatible detergent solution, spray with an instrument pre-soak solution, or cover instruments with a towel moistened with purified water to prevent drying and encrustation of surgical soil.

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
Cleaning agents with a used dilution pH of within 7 alkaline conditions (pH > 10) can damage components/ devices, such as aluminum materials. Do not use saline, environmental disinfection (including chlorine solutions) or surgical antiseptics (such as iodine products). Do not use a cleaning aid that can damage the surface of instruments such as steel wool, abrasive cleaners or wire brushes.

Carefully inspect instruments between uses to verify proper functioning. Return damaged instruments to a DePuy Synthes sales

These reprocessing instructions have been validated as being capable of preparing reusable DePuy Orthopaedics, Inc. and DePuy International Ltd. instruments for reuse. It is the responsibility of the reprocessor to ensure that the reprocessing is performed using appropriate equipment, materials, and personnel to achieve the desired result. This normally requires validation and routine monitoring of the process. Any deviation from these instructions should be evaluated for effectiveness and potential adverse consequences.

Automated Cleaning

1. Screw the protective cap securely onto the shaver handpiece connector


Fig. 7. Screw the protective cap on.

2. Pre-clean through rinsing in cold tap water for a minimum of 1 minute. While rinsing, actuate all moving parts at least five (5) times to clean all surfaces. Rinse the suction port with the suction lever in open position. Moveable parts:

- The Locking Ring: (Inner and Outer parts) Open and close the locking system to ensure that it is free of foreign bodies

Fig. 8. Locking ring

- The Suction Port: Pass running water through the suction port, with the suction lever switch in the fully open position.

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
Fig. 9. Suction port

- The Suction Lever: Open and close the suction lever the inner and outer parts thoroughly. (Fig

Fig. 10. Suction lever

3. Leave the suction lever in the open «MAX» position (Fig instruments per manufacturer's instructions solution. Immerse the instrument in the neutral pH cleaning (pH 7 or may not contain enzymes) solution to avoid aerosol bristle brush to remove all traces of blood and debris; pay close attention to hard-to-reach areas, textured surfaces, or crevices. While soaking, use a soft bristled, tight-fitting lumen brush to brush both the distal and device's lumen for 1 minute. Flush all lumens with the detergent solution for a minimum of 1 minute.


Fig. 11. Suction lever in open position

4. Check for remaining debris and remove as needed.
5. Ultrasonically clean for a minimum of 4 minutes in a second soaking bath using a pH-neutral detergent prepared in accordance with instructions. While submerged, use a syringe to flush all lumens and articulating areas.
6. Rinse the device components with warm tap water for a minimum of one minute and until visual evidence of debris, soil, and cleaning. Include actuation of all moving parts and rinse the suction lumen with the lever in the open position.
7. Confirm the protective cap is securely screwed onto the connector.
8. Load the device components so that the lumens can drain.

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9. Automated cleaning can be conducted in a validated washer compliance to ISO 15883-1 and -2. The following automated cycle is an example of a validated cycle:

Phase	Recirculation Time (Minimum)	Water Temp (Set point)	
Pre-Wash	2:00	Cold Tap Water	
Enzyme wash	1:00	Hot Tap water	
Wash	3:00	60°C	
Rinse	1:00	Hot Tap water	
Rinse	0:20	Ambient	
Thermal Disinfection*	1:00	>90°C	
Dry	7:00	66°C	

*A minimum A_0 of 600 is recommended in accordance to ISO 15883. Higher levels of A_0 can be achieved by increasing the temperature (e.g., A_0

of 3000 at >90°C for 5 min, in accordance with local

**DI/ PURW/ Critical water (high purity water generated by processes such as RO, deionization or distillation)

Manual Cleaning

1. SCREW the protective cap securely onto the shaver handpiece connector 12)

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2. Leave the suction lever in the open «MAX» position


Fig.13. Suction lever in open position

3. Pre-clean through rinsing in cold tap water for a minimum of 1 minute. While rinsing, actuate all moving parts at least five (5) times to clean all surfaces. Rinse the suction port with the suction lever in open position.

4. Carry out a first soaking bath in a pH-neutral detergent solution (pH 7 may or may not contain enzymes) for instruments per manufacturer's instructions (minimum of 5 minutes) in the pH neutral solution. Immerse the instrument in the cleaning solution to avoid aerosol generation.

5. Carry out brushing with a swab-type brush on all surfaces, paying particular attention to lumens, holes and spring moveable mechanisms, such as hinges, box locks, or spring features. While soaking, use a soft-bristled, tight the distal and proximal ends of device's lumen for 1 minute the detergent solution for a minimum of 1 minute.

Moveable parts

- The Locking Ring: Open and close the locking ring to ensure that it is free of foreign bodies (Inner and Outer parts)

Fig. 14. Locking ring

- The Suction Port: Pass running water through the suction port, with the suction lever switch in the fully open position

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
Fig. 15. Suction port

- The Suction Lever: Open and close the suction lever several times. Irrigate the inner and outer parts thoroughly (Fig

Fig. 16. Suction lever

6. Check for remaining debris and remove as needed.

7. Ultrasonically clean for a minimum of 10 minutes in a second soaking bath using a pH-neutral detergent prepared in accordance instructions. While submerged, use a syringe to flush all lumens and articulating areas.

8. Remove the devices and rinse under running ambient (minimum of 1 minute. While rinsing, actuate all times and rinse the suction port with the suction lever in the open position.

9. Submerge each device in PURW/ Critical water for a minimum of 30 seconds with actuation of all moving parts and flushing of the Critical water is high purity water deionization or distillation.

10. Dry using clean lint-free cloths and/or with filtered 100°C or 212°F for a minimum of 30 seconds.

CAUTION: Do not direct air into the shaver/burr seal.

Inspection after Cleaning

1. Inspect before sterilization or storage to ensure the complete removal of soil from surfaces and holes.
2. If soil is still present, re-clean.

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