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
I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG., hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: 

Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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Created November 2016.

INTENDED USE

The ARCHITECT B-R-A-H-M-S PCT Calibrators are for the calibration of the ARCHITECT iSystem when used for the quantitative determination of procalcitonin (PCT) in human serum and plasma.

Refer to the reagent package insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

6 Bottles (2.0 mL each) of ARCHITECT B-R-A-H-M-S PCT Calibrators. Calibrator A contains normal human plasma. Preservatives: ProClin 950 and sodium azide. Calibrators B–F contain different concentrations of recombinant PCT in phosphate buffer. Preservatives: ProClin 950 and sodium azide.

The calibrators are at the following concentrations:

Calibrators	Concentration	
	(ng/mL)	(µg/L)
CAL A	0.00	0.00
CAL B	0.10	0.10
CAL C	0.50	0.50
CAL D	12.10	12.10
CAL E	20.50	20.50
CAL F	100.00	100.00

STANDARDIZATION

There is currently no internationally recognized reference method or reference material for standardization. The ARCHITECT B-R-A-H-M-S PCT Calibrators are traceable to internal reference standards which underwent a one-time value assignment to align with another commercially available assay.

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use
- **Rx ONLY**

- **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human plasma used in Calibrator A is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HIV-1/HIV-2, and anti-HCV.

The following warnings and precautions apply to: **CAL A** **CAL F**

WARNING	Contains methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbott diagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

STORAGE

- Calibrators are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- Calibrators must be stored frozen at -10°C or colder when not in use.

PREPARATION FOR ANALYSIS

Thaw calibrators at room temperature (15–30°C) until completely thawed (30–60 minutes).

Prior to use, mix by gentle inversion (10 times).

After each use, tightly close the caps and return to -10°C or colder storage.

It is suggested to record each thaw date on the carton or the bottles as an aid in tracking the number of times the calibrators are thawed.

Avoid more than 3 freeze/thaw cycles.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CAL A	Calibrator (A,B,C,D,E or F)
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
PRODUCED FOR ABBOTT BY	Produced for Abbott by
PRODUCT OF USA	Product of USA
REF	List Number
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).
WARNING: SENSITIZER	Warning: May cause an allergic reaction.

B-R-A-H-M-S PCT is a trademark and property of Thermo Fisher Scientific Inc. and its subsidiaries.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Fisher Diagnostics,
A Div. of Fisher Scientific Company, LLC
A Part of Thermo Fisher Scientific, Inc.
8365 Valley Pike
Middletown VA 22645 USA

**Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com**

Created November 2016.
©2016 Abbott Laboratories

ADDENDUM

to the Instruction for Use ARCHITECT B-R-A-H-M-S Calibrators for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for device are -10°C or lower.

Medical device must be stored in freezers to ensure regulated temperature conditions with a daily temperature record.

After thaw, must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. The maximum shelf life of the medical device is 9 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010

"Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

ISO 9001: 2008

ISO 13485: 2003 (including CMDCAS)

EN ISO 13485:2012

EN ISO 14971

EN 13641

EN 13612

EN ISO 15223-1

EN 980

EN ISO 18113-1

EN ISO 18113-2

EN ISO 18113-3

EN ISO 23640

EN ISO 20776-1

EN 13975

EN ISO 17511

ISO 80000-1

POTENTIAL USER

Professional technician.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs.

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage.

TECHNICAL CHARACTERISTICS

- Calibrator A: pH-value at 20 °C: 7,4
- Calibrator B-F: pH-value at 20 °C: 6,9
- Density at 20 °C - 1,01 g/cm³

MATERIALS OF ANIMAL/HUMAN ORIGIN

Materials of animal/human origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottiagnostics.com or contact your local representative.

SAFETY PRECAUTIONS

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

CALIBRATION

- Test Calibrators A-F in duplicate. The calibrators should be priority loaded. A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Calibration Range: 0.00 - 100.00 ng/mL (0.00 - 100.00 μ g/L).
- Once an ARCHITECT B·R·A·H·M·S PCT calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used.
 - Daily quality control results are outside of statistically-based quality control limits, as described in the Quality Control Procedure section of this package insert, used to monitor and control system performance.
 - If statistically-based quality control limits are not available then the calibration should not exceed a 30-day limit for recalibration frequency.

The ARCHITECT B·R·A·H·M·S PCT assay may also need to be recalibrated after specified service procedures have been performed or maintenance to critical part or subsystems that might influence the performance of the assay.

MATERIALS REQUIRED BUT NOT PROVIDED

- ARCHITECT B·R·A·H·M·S PCT Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbottiagnostics.com.
- 6P22-10 ARCHITECT B·R·A·H·M·S PCT Controls
- 6P22-25/35 ARCHITECT B·R·A·H·M·S PCT Reagent Kit
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

INFORMATION ON THE MULTIPLICITY OF DEFROST

After thaw and use, store at -10° C or colder. Discard after 3 freeze/thaw cycles.

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