30

APPROVAL

Ethicon, LLC
Organization

Highway 183 Km 8.3, San Lorenzo, 00754, Puerto-Rico, USA
Address


Director, Regulatory Affairs, Ethicon Biosurgery
Occupation


Subscribed and sworn to before me

Joseph Chmielewski
Surname, Name


this 22nd day of September 2022

22nd September, 2022
Date DD.MM.YYYY

Signature

Notary Public

Instruction for use for the medical device

Surgicel Endoscopic Applicator for applying Surgicel Powder

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can be used on its own, or in conjunction with the rigid cannula provided.

INDICATIONS AND CONTRAINDICATIONS FOR MD USE AND POSSIBLE ADVERSE EFFECTS

Indications

The SURGICEL Endoscopic Applicator is intended for use in delivering SURGICEL Powder, absorbable haemostat to bleeding surgical sites through a 5 mm or larger trocar.

Contraindications

- Do not inject or place SURGICEL™ Powder into an open blood vessel. Do not use to treat bleeding from large defects in arteries or veins.
- SURGICEL™ Powder should not be used to control hemorrhage from large arteries or veins.
- The SURGICEL™ Endoscopic Applicator was not designed for intraluminal procedures.

Warnings

The SURGICEL™ Endoscopic Applicator is supplied sterile and as the material is not compatible with autoclaving or ethylene oxide sterilization, the SURGICEL™ Endoscopic Applicator should not be resterilized.

- To prevent clogging, do not touch the tip to wet surface. Be careful to avoid damaging tissue with the rigid tip.
- Do not attempt to trim the applicator tip. Replace the tip if it becomes clogged.
- Do not use the SURGICEL Endoscopic Applicator if package is opened or damaged.
- Apply SURGICEL Powder according to the product's labeling.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

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Precautions

- The SURGICEL Endoscopic Applicator is supplied with a flexible inner tip inside a rigid cannula. The rigid cannula cannot be used independently.
- The SURGICEL Endoscopic Applicator should only be used by persons having adequate training and familiarity with endoscopic techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any endoscopic procedure.
- To prevent inadvertent device spillage, or unintended contact with tissue, organs, or blood, maintain visualization of the SURGICEL Endoscopic Applicator tip at all times.
- Do not compress or excessively bend the flexible inner tip which could obstruct the application of the powder. It is possible that powder accumulated in the applicator could disperse beyond the target bleeding site upon compression of the bellows, which may require additional irrigation and aspiration.
- After use, dispose of the SURGICEL Endoscopic Applicator in accordance with biohazardous material protocol.

Adverse reactions and events



None stated.

Directions for use

1. Open package and transfer SURGICEL Endoscopic Applicator to sterile field. The SURGICEL Endoscopic Applicator is composed of two parts, a flexible inner tip and a rigid cannula.

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



2. Remove the cap from the SURGICEL Powder Device (Fig. 1).


3. Twist to open (Figs. 2, 3).


4. Attach the SURGICEL Endoscopic Applicator to the SURGICEL Powder Device without pumping the bellows and ensure the connection is secure (Fig. 4).

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

5. The device with the rigid tip is now ready for use (Fig. 5).

6. If a completely flexible tip is desired, remove the rigid cannula (Fig. 6-optional).

7. Introduce the SURGICEL Endoscopic Applicator into the trocar port. **CAUTION:** Insert the SURGICEL Endoscopic Applicator using caution not to inadvertently express powder. Carefully position the distal end of the SURGICEL Endoscopic Applicator at the site where the powder is to be delivered. To prevent clogging, do not touch the tip to wet surface. Do not trim the applicator tip (see WARNINGS). Be careful to avoid damaging tissue with the rigid tip.

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8. There is potential for clogging with use of the SURGICEL Endoscopic Applicator. In the event of clogging, the tip can be wiped off using dry, sterile surgical gauze to remove clog. If clog cannot be removed using dry gauze, dispose of the clogged tip and use a new SURGICEL Endoscopic Applicator.

9. While holding the SURGICEL Endoscopic Applicator in place, pump to apply powder to treatment site. NOTE: To view remaining powder in the reservoir or flexible inner tip while using the rigid cannula, partially retract the SURGICEL Powder device and flexible inner tip from the rigid cannula.

10. Carefully remove the SURGICEL Endoscopic Applicator from the trocar port when sufficient powder has been delivered to the treatment site. Apply direct pressure over the treated site as needed without disturbing the clot.

NOTE: Distance for application in laparoscopic settings may vary between procedure types and individual patient anatomy. In preclinical animal studies, safety and efficacy were demonstrated at application distances ranging from 1cm to 10cm from the bleeding site.

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INFORMATION ON POTENTIAL USERS OF MD AND DIRECTIONS FOR USE

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Directions for use

Devices are intended for use in healthcare facilities.

Devices are recommended for use under operation conditions indicated below.

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Potential users of MD

SURGICEL Endoscopic Applicator is intended for use during surgical procedures, by health care professionals.

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Use/operation conditions of MD

The device is used in the internal environment of the body.

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Once the “open” applicator tip is removed from the SURGICEL Powder device, the Endoscopic Applicator tip can be attached to the SURGICEL Powder device.

Figure 9 shows the SURGICEL Powder device with the “Open” Applicator Tip fully removed, and the flexible inner tip of the Endoscopic Applicator now attached to the device. Figure 4 shows the complete Endoscopic Ap-plicator device attached to the SURGICEL Powder device.

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The possibility and ways of integrating with other medical devices
SURGICEL Endoscopic Applicator is used with the SURGICEL Powder Device for endoscopic or laparoscopic powder application.

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Table 2 - Technical characteristics of rigid cannula

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CONDITIONS OF TRANSPORTATION AND STORAGE

No special storage conditions required. Do not use after expiration date.

Transportation conditions	
Humidity	No special conditions
Temperature	
Pressure	
Storage conditions	
Humidity	No special conditions
Temperature	
Pressure	

MD DISPOSAL RULES AND CONDITIONS

After use, dispose of the SURGICEL Endoscopic Applicator in accordance with biohazardous material protocol.

WARRANTY

Manufacturer guarantees to the end user absence of defects and conformity to product specifications for the entire shelf life of the product under the conditions of adherence to storage, transportation and package integrity requirements. There are no warranties that extend beyond the description on the face of instruction for use. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use.

NAME AND LEGAL ADDRESS OF MANUFACTURER AND AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION

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PLACE OF PRODUCTION:

Alltrista Plastics LLC, 20 Setar Way, Reedsville, Pennsylvania, 17084, USA

**AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN
RUSSIA:**

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