Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879

To Whom It May Concern

It is hereby certified that the document attached is the 'Instructions for Use' (0110119 Rev. 0) for the Admiral Xtreme PTA Balloon Catheter.

Yours sincerely,

Signed on behalf of Medtronic Inc.
Sara Bakker
Regulatory Affairs Manager

15 DEC. 2021

Date

Medtronic

Admiral™ Xtreme

PTA Balloon Catheter	4
PTA балонен катетър	5
Balónkový katétr pro PTA	7
PTA-balloonkateter	8
PTA-Balloonkatheter	10
Καθετήρας μπαλονιού PTA (διαθερμικής διαυλικής αγγειοπλαστικής)	12
Catéter de balón para ATP	14
PTA balloonkateeter	15
PTA-pallokatetri	17
Cathéter à ballonnet pour ATP	18
Balonski kateter za PTA	20
PTA-balloonkatéter	21
Catetere a palloncino per PTA	23
PTA 풍선 카테터	25
PTA balioninis kateteris	26
PTA balonkatetrs	28
Балон катетер за перкутана транслюминална ангиопластика (ПТА)	29
PTA-balloonkatheter	31
PTA-balloonkateter	33
Cewnik balonowy do PTA	34
Cateter balão para PTA	36
Cateter de balão para PTA	38
Cateter cu balon pentru ATP	39
Баллонный катетер для ЧТА	41
Balónkový katéter na PTA	43
Balonski kateter PTA	44
PTA balon-kateter	46
PTA-balloonkateter	48
PTA Balon Kateteri	49
Балонний катетер для ЧТА	51

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning •
Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet •
Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • 사용 지침 •
Naudojimo instrukcijos • Lietošanas pamācība • Упатства за употреба • Gebruiksaanwijzing •
Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização •
Instrucțiuni de utilizare • инструкция по применению • Pokyny na používanie • Navodila za uporabo •
Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з експлуатації

2017-03-06
0110119 Rev. 0

2017/APR/06 at 10:39 a.m. Doc number: 0110118 Rev. 0 [Global]
Printspec a5-55 - 148 x 210 mm

global:20170221, alldocModCol:20170216, alldocStyleCol:20150528, master_pageset:20150528, languageLookup:20161110, table:
20170228, characterJoin:20170404, 3of9:20110211, 3of9_svg:20110211

Symbols / Символы / Symboly / Symboler / Symbole / Σύμβολα / Símbolos / Symbold / Symbolit / Symboles / Oznake / Szimbólumok / Simboli / 기호 / Simboliai / Simboli / Символи / Symbolen / Symboler / Symbole / Símbolos / Símbolos / Simboluri / Символы / Symboly / Simboli / Simboli / Symboler / Semboller / Символы

Consult instructions for use / Разгледайте инструкции за употреба / Víz návod k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar instrucciones de uso / Lugege kasutusjuhendit / Katso käyttöohjeet / Consulter le mode d'emploi / Pogledajte upute za upotrebu / Läs a hanvändar-
ömnututót / Consultare le istruzioni per l'uso / 사용 지침 참고 / Skaitkite naudojimo instrukcijas / Skatīt lietošanas pamācību /
Прочитайте го Упатството за употреба / Zie gebruiksaanwijzing / Se i brugsanvisningen / Należy zapoznać się z instrukcją
użytkowania / Consultar as instruções de uso / Consultar as instruções de utilização / Consultati Instrucțiunile de utilizare /
См. инструкцию по эксплуатации / Prečítajte si návod na použitie / Glejte navodila za uporabo / Pogledajte uputstva za
upotrebu / Läs bruksanvisningen / Kullanim talimatlarına bakın / Див. інструкцію з експлуатації

Catalogue number / Каталоген номер / Katalogové číslo / Katalognummer / Katalognummer / Αριθμός καταλόγου / Número
de catálogo / Katalooginumber / Luettelonumero / Numéro de catalogue / Kataloški broj / Katalógusszám / Numero di
catalogo / 카탈로그 번호 / Katalogo numeris / Kataloga numurs / Каталогски број / Catalogusnummer / Artikelnummer /
Nimer katalogowy / Número de catálogo / Número de catálogo / Număr de catalog / Номер по каталору / Katalogové číslo /
Kataloška številka / Kataloški broj / Katalognummer / Katalog numarası / Номер у каталогу

Lot number / Номер на партида / Číslo šarže / Partinummer / Chargennummer / Αριθμός παρτίδας / Número de lote / Parti
number / Eränumero / Numéro de lot / Broj serije / Gyártási szám / Numero di lotto / 로트 번호 / Partijos numeris / Partijas
numurs / Број на серија / Partinummer / Lotnummer / Numer partii produkcyjne / Número de lote / Número de lote / Număr
de lot / Номер партии / Číslo šarže / Serijska številka / Broj serije / Lotnummer / Lot numarası / Номер партии

Date of manufacture / Дата на производство / Datum výroby / Fremstillingsdato / Herstellungsdatum / Ημερομηνία
κατασκευής / Fecha de fabricación / Tootmiskuupeä / Valmistuspäivämäärä / Date de fabrication / Datum proizvodnje /
Gyártási ideje / Data di fabbricazione / 제조일자 / Pagaminimo data / Ražošanas datums / Датум на производство /
Produktionsdatum / Produksjonsdato / Data produkci / Data de fabricação / Data de fabrico / Data fabricării / Дата
изготовления / Datum výroby / Datum izdelave / Datum proizvodnje / Tilverkningsdatum / Üretim tarihi / Дата изготовления

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Valmistaja /
Fabricant / Proizvođač / Gyártó / Fabricante / 제조업체 / Gamintojas / Ražotājs / Производител / Fabrikant / Producent /
Prodcent / Fabricante / Fabricante / Producător / Производител / Výrobca / Izdelovalec / Proizvođač / Tilverker / Üretici /
Виробник

Manufactured in / Произведено в / Vyrobeno v / Fremstillet i / Hergestellt in / Κατασκευάστηκε σε / Fabricado en /
Tootmistukohd / Valmistuspaikka / Lieu de fabrication / Zemlja proizvodnje / A gyártás helye / Fabricato in / 제조 국가 /
Pagaminia / Ražošanas vieta / Произведено во / Vervaardigd in / Produziert in / Miejsce produkcji / Fabricado em / Fabricado
em / Fabricat în / Произведено в / Vyrobené v / Izdelano v / Proizvedeno u / Tillverkad i / Imalat yeri / Micile de fabricație

Use-by date / Използвайте не по-късно от / Datum použitelnosti / Kan anvendes til og med / Verwendbar bis / Ημερομηνία
«Χρήση έως» / Fecha de caducidad / Köblik kuni / Viimeinen käyttöäivämäärä / Date de péremption / Rok upotrebe / Lejārti
idā / Utilizzare entro / 사용 기한 / Naudoti iki / Deriguma termiņš / Употреблять до / Uiterste gebruiksdatum / Siste
forbrukstid / Termin przydatności do użycia / Data de vencimento / Utilizar antes da data / A se utiliza până la data de / Срок
годности / Datum najneskoršej spotreby / Uporabno do / Datum „Upotrebljivo do“ / Utgångsdatum / Son kullanna tanhi /
Термін придатності

Sterilized using ethylene oxide / Стерилизовано с етиленов оксид / Sterilizováno ethylenoxidem / Steriliseret med
ethylenoxid / Sterilisiert mit Ethylenoxid / Αποστειρωμένο με αιθυλνοξείδιο / Esterilizado mediante óxido de etileno /
Steriliseeritud etüleenoksiidiga / Steriloitu etyleenoksidilla / Stérilisé à l'oxyde d'éthylène / Sterilizarono etilen oksidom / Etüleen-
oksidall sterilizálva / Sterilizzato a ossido di etilene / 산화 에틸렌을 사용하여 멸균함 / Sterilizuota etileno oksidu / Sterilizēts,
izmantot etilēnoksidu / Sterilizirano so etilen oksid / Gesteriliseerd met ethyleenoxide / Steriliseret med etylenoksid /
Steryliowany tlenkiem etylenu / Esterilizado com óxido de etileno / Esterilizado utilizando óxido de etileno / Sterilizat cu oxid
de etilenă / Sterilizovano etilenoksidom / Sterilizované pomocou etylénoksidu / Sterilizirano z etilenoksidom / Sterilizano
konččenjem etilen-okside / Steriliserad med etylenoxid / Etülen oksit ile sterilize edilmiştir / Стерилизовано за допомогою
етиленоксида

Caution / Внимание / Pozor / Forsigtig / Vorsicht / Προσοχή / Precaución / Ettevaatus / Varoitus / Attention / Oprez /
Figyelem / Attenzione / 주의 / Demesio / Uzmanību / Внимание / Let op / Forsiktig / Przestroga / Advertência / Atenção /
Atenție / Предупреждение / Varovanie / Pozor / Oprez / Var forsiktig / Dikkat / Увага!

Do not reuse / Да не се използва повторно / Nepoužívejte opakovaně / Må ikke genbruges / Nicht wiederverwenden / Мнү
статочупотребляти / No reutilizar / Mitte kasutada korduvalt / Älä käytä uudelleen / Ne pas réutiliser / Nemojte upotrebljavati
više puta / Kizárólag egyszer használatra / Non riutilizzare / 재사용 금지 / Nenaudoti pakartotinai / Nelietot atkārtoti / Само за
еднократна употреба / Niet opnieuw gebruiken / Skal ikke brukes flere ganger / Nie stosować ponownie / Não reutilizar / Não
reutilizável / De unică folosință / Не использовать повторно / Nepoužívejte opakovaně / Za enkratno uporabo / Nije za
ponovnu upotrebu / Får inte återanvändas / Yeniden kullanmayın / He zastrocovувати повторно

Do not re-sterilize / Да не се стерилизира повторно / Neprovádějte sterilizaci / Må ikke reesteriliseres / Nicht reesterilisieren /
Мнү не повторювати стерилізацію / No reesterilizar / Mitte steriliseenda uuesti / Älä steriloi uudelleen / Ne pas reesteriliser / Nemojte
iznova sterilizirati / Nem újrasztelizálható / Non riesterilizzare / 재멸균 금지 / Nestenilizuokite pakartotinai / Nesterilizēt atkārtoti /
Да не се стерилизира повторно / Niet opnieuw steriliseren / Skal ikke reesteriliseres / Nie resterylizować / Não reesterilizar /
Não reesterilizar / Nu resterilizați / Не стерилизовать повторно / Opakovane nesterilizuje / Ne sterilizirajte znova / Nije za
ponovnu sterilizaciju / Får inte omsteriliseras / Yeniden sterilize etmeyin / He стерилизувати повторно

Keep dry / Поддържайте сухо / Udržujte v suchu / Skal opbevares tørt / Trocken aufbewahren / Διατηρήστε στεγνό /
Mantener seco / Hoida kuivas / Säilytä kuivassa / À conserver dans un endroit sec / Držite na suhom / Szárazon tartandó /
Mantenere asciutto / 건조한 상태 유지 / Laikyti sausioje vietoje / Uzglabāt sausā vietā / Да се чува на сухо място / Droog
bewaren / Skal holdes tørt / Chronić przed wilgocią / Manter seco / Manter seco / Pástrati într-un spațiu uscat / Храните в
сухом месте / Uchovávejte v suchu / Hranite v suhem prostoru / Držati na suvom mestu / Förvaras tørt / Kuru yerde
saklayın / Зберігати у сухому місці

Keep away from sunlight / Пазете от слънчева светлина / Chraňte před slunečním světlem / Må ikke udsættes for direkte
sollys / Vor Sonnenlicht schützen / Διατηρείται μακριά από το ηλιακό φως / Mantener alejado de la luz del sol / Hoida
pääkesevalguse eest / Suojaa auringonvalolta / Ne pas exposer à la lumière du soleil / Držite podalje od sunčeve svetlosti /
Napfénytől védendő / Tenere lontano da fonti di calore / 직사광선을 피할 것 / Saugoti nuo saulės šviesos / Sargāt no tiešiem
saules stariem / Да се чува подальше от солнечной светлина / Uit de buurt van zonlicht houden / Må ikke utsettes for sollys /
Chronić przed światłem słonecznym / Manter afastado da luz do sol / Manter afastado da luz do sol / Nu expuneți la radiația
solară / Не допускать попадания солнечных лучей / Uchovávejte mimo slunečního žiarenia / Zaščitite pred sončno
svetlobo / Držati daleko od sunčeve svetlosti / Får inte utsättas för solljus / Güneş ışığından uzak tutun / Зберігати у місці,
захищеному від прямих сонячних променів

Do not use if package is damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, pokud je balení
poškozeno / Må ikke anvendes, hvis emballagen er beskadiget / Nicht verwenden, wenn die Verpackung beschädigt ist / Μη
χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημία / No utilizar si el envase está dañado / Mitte kasutada, kui pakend on
kahjustatud / Älä käytä, jos pakkaus on vaurioitunut / Ne pas utiliser si l'emballage est endommagé / Nemojte upotrebljavati
ako je pakiranje oštećeno / Ne használja fel, ha a csomagolás sérült / Non utilizzare se l'imballaggio non è integro / 포장에 손
상된 경우 사용 금지 / Ne naudoti, jei pakuotė pažeista / Nelietot, ja iepakojums ir bojāts / Да не се употребява доколку
пакването е оштетено / Niet gebruiken als de verpakking beschadigd is / Skal ikke brukes hvis pakningen er skadet / Nie
stosować, jeśli opakowanie jest uszkodzone / Nao utilizar se a embalagem estiver danificada / Nao utilizar se a embalagem
estiver danificada / Nu folosiți dacă ambalajul este deteriorat / Не използвайте, если упаковка повреждена / Nepoužívejte,
ak je obal poškozený / Ne uporabljajte, če je embalaža poškodovana / Ne koristiti ako je pakovanje oštećeno / Får inte
användas om förpackningen är skadad / Ambalaj hasar gördüyse kullanmayın / Не використовувати, якщо упаковка
пошкоджена.

In Conformity with the European Council Directive 93/42/EEC concerning medical devices / В съответствие с Европейската
Директива 93/42/ЕО на Съвета относно медицинските изделия / Splňuje požadavky směrnice Evropské rady 93/42/EHS
týkající se zdravotnických prostředků / I overensstemmelse med Rådets Direktiv 93/42/EØF om medicinsk udstyr / Entspricht
den Vorschriften der Richtlinie 93/42/EWG des Europäischen Rats über Medizinprodukte / Σχ. συμμόρφωση με την
Οδηγία 93/42/ΕΟΚ του Ευρωπαϊκού Συμβουλίου περί ιατροτεχνολογικών προϊόντων / De conformidat cu la Directiva del
Consejo Europeo 93/42/CEE relativa a productos sanitarios / Vastab Euroopa Liidu Nõukogu meditsiiniseadmete direktiivi
93/42/EMU nõuetele / Euroopan unionin neuvoston lääkitäilaitteita koskevan direktiivin 93/42/ETY mukainen / Conforme à la
directive 93/42/CEE du Conseil de l'Union européenne concernant les dispositifs médicaux / Usklađeno s direktivom
Evropskog vijeća 93/42/EEZ o medicinskim uređajima / Megfelel az Európai Tanács orvostechnikai eszközökre vonatkozó
93/42/EGK irányelvének / In conformità con la Direttiva del Consiglio Europeo 93/42/CEE in materia di dispositivi medici / 의료
기기에 관한 유럽이사회 지침 93/42/EEC 준수 / Atitinka Europos Tarybos direktyvą 93/42/EEB dėl medicinios prietaisų /
Atbilst Eiropas Padomes Direktīvai 93/42/EEK par medicīnas ierīcēm / Во соответствии со Директивата 93/42/ЕЕЗ на
Европейскиот совет за медицински помагала / Voldoet aan de Europese Richtlijn 93/42/EEG met betrekking tot medische
apparatuur / I samsvar med rådsdirektiv 93/42/EØF om medisinske udstyr / Zgodność z dyrektywą Rady Unii Europejskiej
93/42/EWG dotyczącą wyrobów medycznych / Em conformidade com a Diretiva do Conselho Europeu 93/42/CEE relativa a
dispositivos médicos / Em conformidade com a Diretiva do Conselho Europeu 93/42/CEE relativa a dispositivos médicos / În
conformitate cu Directiva 93/42/CEE a Consiliului European privind dispozitivele medicale / Соответствие требованиям
директивы Европейского Совета 93/42/ЕЕС относительно устройств медицинского назначения / Zariadenie je v súlade
s ustanoveniami smernice Európskej rady 93/42/EHS o zdravotníckych pomôckach / V skladu z Direktivo Sveta 93/42/EGS o
medicinskih pripomočkih / Usklađeno sa Direktivom Saveta Evrope 93/42/EEC u vezi sa medicinskim uređajima /
Överensstämmer med rådets direktiv 93/42/EEG om medicintekniska produkter / Tibbi cihazları konulu alan 93/42/EEC sayılı
Avrupa Konseyi Direktifi ile Uyumludur / Згідно з вимогами Директиви Ради Європи 93/42/ЕЕС стосовно пристроїв
медичного призначення

Over the wire / Over the wire (Катетър с метален водач) / Over The Wire (po drátě) / Over The Wire / Over-The-Wire / Eml
tou οδηγού σύρματος (Over The Wire: OTW) / Sobre guía / Üle traadi sisestatav / Väijena pitkin asetettava / Sur guide /
Putem țice / Vezetődróttal beültethető / Over The Wire (su filo guida) / 오버 더 와이어(Over The Wire) / Istumiamas viela /
Izmantojams kopā ar vadītājsīgu / Преву жица / Over-the-wire / Over the wire (over mandrengen) / Nadprowadnikowy
(OTW) / Sobre o fio / Over the wire (sobre o fio) / Peste firul de ghidare / Доставка по проводнику / Zariadenie typu Over The
Wire (s vodiacim drôtom) / Za vstavljanje z žico / Preko žice / Over The Wire (OTW, over ledaren) / Tel úzen / Система
доставки за допомогою проводника

Nominal pressure / Номинално налягане / Nominální tlak / Nominell tryk / Nenndruck / Ονομαστική πίεση / Presión nominal /
Nimirohk / Nimellinen paine / Pression nominale / Nominalni tlak / Névleges nyomás / Pressione nominale / 호칭 압력 /
Nominālais slēgis / Nominālais spiediens / Номинален притисок / Nominale druk / Nominell trykk / Ciśnienie znamionowe /
Pressão nominal / Pressão nominal / Presiune nominală / Номинальное давление / Nominálny tlak / Nominalni tlak /
Nominální pritisak / Nominell tryck / Nominal basınç / Номинальный тиск

Rated burst pressure / Номинално налягане на спукване / Maximální spolehlivý plnicí tlak / Nominell sprængningstryk /
Maximaler Nenndruck / Μέγιστη πίεση ρήξης / Presión de rotura / Nimilõhkemürõhk / Nimellinen puhkeamisaine / Pression
théorique de rupture / Navedeni tlak pucanja / Névleges repedési nyomás / Pressione nominale di scoppio / 규정 파열 압력 /
Nominālais spriegums / Nominālais pārraušanas spiediens / Номинален притисок на пуканье / Vastgesteld barstdruk /
Nominell sprengtrykk / Znamionowe ciśnienie rozrywania / Pressão de ruptura nominal / Pressão de ruptura nominal / Presiune
nominală de rupere / Расчетное давление разрыва / Menovity tlak prasknutia / Ocenjeni razpočni tlak / Procenjeni pritisak
pucanja / Nominell bristningstryck / Anna patlama basıncı / Розрахунковий тиск розриву

Maximum guidewire diameter / Максимален диаметър на водач / Maximální průměr vodičho drátu / Maksimal guidewirediameter / Maximaler Durchmesser des Führungsdrahts / Μέγιστη διέκτρος οδηγού σύρματος / Diametro máximo de la guía / Maksimaalne juhtetraadi läbimõõt / Ohjainvajerin enimmäisläpimitta / Diamètre maximum du fil-guide / Najveći promjer žice vodilice / A vezetődírt legnagyobb átmérője / Diametro massimo del filo guida / 가이드 와이어 최대 직경 / Didžiausias kreipiamosios vielos skersmuo / Maksimālais vadītājsfigas diametrs / Максимален диаметър на жица-водилка / Maximale diameter voerdraad / Sterste diameter på ledevaier / Maksymalna średnica przewodnika / Diâmetro máximo do fio-guia / Diâmetro máximo do fio-guia / Diametru maxim fir de ghidare / Максимальный диаметр проводника / Maximálny priemer vodiaceho drôtu / Najveći premer vodične žice / Maksimalni prečnik vodič-žice / Ledarens maximala diameter / Maksimum kilavuz tel çapı / Максимальный диаметр проводника

Balloon diameter / Диаметър на балона / Průměr balónku / Ballondiameter / Ballondurchmesser / Διέκτρος μπαλονιού / Diametro del balón / Balloni läbimõõt / Pallon läpimitta / Diamètre du ballonnet / Promjer balona / Lunghezza del palloncino / Diametro del palloncino / Balonelo skersmuo / Balona diametrs / Диаметар на балонот / Ballondiameter / Ballongdiameter / Średnica balonu / Diâmetro do balão / Diâmetro do balão / Diametru balon / Диаметр баллона / Premier balónika / Premier balona / Prečnik balona / Ballongens diameter / Balon çapı / Диаметр балона

Balloon length / Дължина на балона / Délka balónku / Ballonlänge / Länge des Ballons / Μήκος μπαλονιού / Longitud del balón / Balloni pikkus / Pallon pituus / Longueur du ballonnet / Duljina balona / A ballon hossza / Lunghezza del palloncino / Balonelo ilgis / Balona garums / Должина на балонот / Ballonlengte / Ballonglengde / Długość balonu / Comprimento do balão / Comprimento do balão / Lungime balon / Дължина балона / Dłzka balónika / Dolžina balona / Dužina balona / Ballongens längd / Balon uzunluğu / Довжина балона

Minimum sheath inner diameter / Минимален вътрешен диаметър на катетъра / Minimální vnitřní průměr pouzdra / Sheathens mindste indvendige diameter / Minimaler Innendurchmesser der Einführschleuse / Ελάχιστη εσωτερική διέκτρος θήκας / Diametro interno minimo de la vaina / Minimale hülsl siseläbimõõt / Holkin pienin sisäläpimitta / Diamètre interne minimum de la gaine / Najmanjši unutarnji promjer obloge / Bevezetőhüvely minimális belső átmérője / Diametro interno minimo della guaina / 와이어 최소 내경 / Mažiausias vidinis movos skersmuo / Ievadslēdzas minimālais iekšējais diametrs / Минимален вътрешен диаметър на обвивката / Minimale binnendiameter sheath / Minste indre diameter på hylse / Minimalna średnica wewnętrzna koszulki / Diâmetro interno mínimo da bainha / Diâmetro interno mínimo da bainha / Diametru intern minim al tecii / Минимальный внутренний диаметр оболочки / Minimálny vnútorný priemer puzdra / Najmanjši notranji premer tulca / Minimalni unutrašnji prečnik omotača / Hylsans minsta innerdiameter / Minimum kilif iç çapı / Минимальный внутренний диаметр оболочки

Nonpyrogenic / Непирогенно / Nепyroгenni / Ikke-pyrogen / Pyrogenfrei / Μη πυρρογόνο / Aπирόγено / Mitterpyrogeenne / Pyrogeenit / Apyroгène / Nezapaljivo / Nem piroгén / Apirogeno / 비발열성 / Nedegus / Apiroгēns / Алирогено / Niet-pyrogen / Pyrogenfri / Produkt niepiroгenny / Não piroгénico / Apiroгénico / Apyrogen / Алирогено / Apyroгēne / Apyrogeno / Nije pirogeno / Ikke-pyrogen / Pirogenik Değildir / Алирогено

Quantity / Количество / Množství / Antal / Menge / Ποσότητα / Cantidad / Kogus / Määr / Quantité / Količina / Mennyiség / Quantità / 수량 / Kiekis / Daudzums / Коничина / Aantal / Antall / Ilość / Quantidade / Quantidade / Cantitate / Количество / Množstvo / Količina / Količina / Antal / Miktar / Кількість

Do not exceed rated burst pressure / Не превишавайте номиналното налягане на спукване / Nepřekračujte maximální spolehlivý plnicí tlak / Det nominelle sprængningstryk må ikke overskrides / Maximalen Nenndruck nicht überschreiten / Μη υπερβαίνετε τη μέγιστη πίεση ρήξης / No superar la presión de rotura / Mitte ületada nimilõhkemisrõhku / Älä ylittää nimellistä puhkeamispainetta / Ne pas dépasser la pression théorique de rupture / Tlak uređaja ne smije premašiti navedeni tlak pucanja / Ne lēpje tūl a nēvleges repedési nyomást / Non superare la pressione nominale di scoppio / 규정마열압력 초과 금지 / Neviršykite nominalaus plyšimo slėgio / Nepārsnīgt nominālo pārsaušanas spiedienu / Да не се надминава номиналниот притисок на пукање / Vastgestelde barstdruk niet overschrijden / Nominelt sprengrtrykk må ikke overstiges / Nie należy przekraczać znamionowego ciśnienia rozerwania / Não exceder a pressão de ruptura nominal / Não exceda a pressão de ruptura nominal / Nu depășiți presiunea nominală de rupere / Не превышать расчетное давление разрыва / Nepražaukite menovytį tlak prasknuta / Ne prekaročite ocenjenega razpočnega tlaka / Nemojte premašiti procenjeni pritisak pucanja / Överskrid ej nominellt bristningstryck / Anna patlama basincini aşmayın / Не превышать розрахункове значення тиску розриву

Inflation pressure / Налягане при разгъване / Plnicí tlak / Inflationstryk / Inflationsdruck / Πίεση εμφύσησης / Presión de inflado / Täitmisrõhk / Täyttöpaine / Pression de gonflage / Tlak napuhavanja / Feltöltési nyomás / Pressione di gonfiaggio / 팽창 압력 / Pnplidymo slėgis / Uzplidšanas spiediens / Притисок на дуење / Vuldruk / Fyllingstrykk / Ciśnienie napełniania / Pressão de insuflação / Pressão de insuflação / Presiune de umflare / Давление заполнения / Inflačný tlak / Tlak polnjenja / Pritisak naduvavanja / Fyllningstryck / Şişirme basıncı / Тиск роздування

Usable catheter length / Използваема дължина на катетъра / Použitelná délka katétru / Anvendelig kateterlængde / Nutzlänge des Katheters / Οφέλιμο μήκος καθετήρα / Longitud útil del catéter / Kasutatav katetri pikkus / Katetrin käyttöpituus / Longueur utile du cathéter / Upotrebjliva duljina katetera / Hasznos katéterhossz / Lunghezza utile del catetere / 사용 가능한 카테터 길이 / Naudojamas kateterio ilgis / Lietojamais katetra garums / Unotrebjliva dužina na kateterot / Bruikbare lengte van de katheter / Nyttbar kateterlengde / Użyteczna długość cewnika / Comprimento útil do cateter / Comprimento útil do cateter / Lungime utilizabilă a cateterului / Рабочая длина катетера / Použitelná dĺžka katétra / Uporabna dolžina katetra / Upotrebjliva dužina katetera / Användbar kateterlängd / Kullanlabilir kateter uzunluğu / Рабочая довжина катетера

English

Admiral™ Xtreme

PTA Balloon Catheter

1. Description

The Admiral™ Xtreme is an Over the Wire (OTW) peripheral balloon catheter designed for percutaneous transluminal angioplasty in atherosclerotic obstructed vessels. The catheter has a dual-lumen shaft with a balloon welded onto the distal tip. The dual-lumen shaft is branched at the proximal end. One tube forms the entrance to the central lumen for the guidewire, while the other tube is used to inflate and deflate the dilatation balloon with a mixture of contrast medium and saline solution. The catheter construction and the special balloon material are designed so that a specific balloon diameter can be reached, depending on the balloon size and defined pressure. The length of each balloon is specified. The Admiral Xtreme is available in different balloon sizes. Nominal balloon diameter and lengths are printed on the hub.

2. Indications

The Admiral Xtreme is indicated for percutaneous transluminal angioplasty (PTA) in patients with obstructive disease of peripheral arteries.

3. Contraindications

- The Admiral Xtreme PTA balloon catheter is contraindicated for use in coronary arteries and supra-aortic vessels.
- The Admiral Xtreme PTA balloon catheter is contraindicated for cases where the guidewire cannot cross the target lesion.

4. Warnings

- This device is designed and intended for single use only. DO NOT RSTERILIZE OR REUSE. Reuse or resterilization may create a risk of contamination of the device or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. Reuse or resterilization may compromise the structural integrity of the device or lead to device failure which, in turn, may result in patient injury, illness, and death. Medtronic is not responsible for any direct, incidental, or consequential damages resulting from resterilization or reuse.
- To reduce the potential for vessel damage, inflate the balloon to a diameter approximating that of the vessel just proximal and distal to the stenosis.
- When the PTA balloon catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation.
- Do not manipulate the PTA balloon catheter in the inflated state. The position of the PTA balloon catheter can only be changed with the guidewire in place.
- If resistance occurs during manipulation, ascertain the cause via fluoroscopy, road mapping, or DSA before the balloon catheter is moved backward or forward.
- Do not under any circumstance move the guidewire during inflation of the Admiral Xtreme PTA balloon catheter.
- The balloon must be completely deflated before retrieving the catheter from the vascular system.
- Do not exceed the Rated Burst Pressure (RBP). The RBP is based on the results of in vitro testing. Use of a pressure-monitoring device is recommended to prevent over pressurization.
- Do not use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium (equal parts contrast medium and saline solution).
- Use the catheter prior to the Use-by date specified on the package.

5. Precautions

- Only interventionalists who have sufficient experience in percutaneous transluminal angioplasty should use this PTA balloon catheter. A thorough understanding of the technical principles, clinical applications, and risks associated with percutaneous transluminal angioplasty is necessary before using this product.
- Administer appropriate drug therapy (anticoagulant, vasodilator, etc.) to the patient according to standard protocols for PTA before inserting the Admiral Xtreme PTA balloon catheter.
- Exercise care during handling in order to avoid possible damage to the catheter. Do not use a catheter that has been damaged.
- To minimize the introduction of air, aspirate and flush the system and keep a tight catheter connection throughout the procedure.

- Use the device with caution for procedures involving calcified lesions or synthetic vascular graft due to the abrasive nature of these lesions.
 - Identify allergic reactions to contrast medium before treatment.
 - The general technical requirements for catheter insertion must be observed at all times. This includes flushing the components with sterile, isotonic saline solution prior to use and the usual prophylactic, systemic heparinization.
 - Catheter applications vary. Select the technique on the basis of the patient's condition and the experience of the interventionalist.
 - The maximal diameter of the guidewire must not exceed 0.035 in (0.89 mm). Certain coatings on guidewires may lead to increased friction between the guidewire and the catheter.
- Caution:** Larger sizes of the Admiral Xtreme PTA balloon catheter may exhibit slower deflation times particularly on long catheter shafts.

6. Potential complications / Adverse effects

Possible complications associated with intravascular balloon catheter implantation may include, but are not limited to:

Puncture related:

- Local hematoma
- Local hemorrhage
- Local or distal thromboembolic episodes
- Thrombosis
- Arteriovenous fistula
- Pseudoaneurysm
- Local infections

Dilatation related:

- Dissection in the dilated artery wall
- Perforation of the artery wall
- Prolonged spasms
- Acute reocclusion necessitating surgical intervention
- Restenosis of the dilated artery
- Total occlusion of the peripheral artery.

Angiography related:

- Hypotension
- Pain and tenderness
- Arrhythmias
- Sepsis/infection
- Systemic embolization
- Endocarditis
- Short-term hemodynamic deterioration
- Death
- Drug reactions
- Allergic reaction to contrast medium
- Pyrogenic reaction

7. Preparation technique

Inspect the Admiral Xtreme PTA balloon catheter before to the procedure to verify that the product is intact and functional. Do not use if the outer or the inner packaging is damaged or opened. After moistening the protective sheath with physiological saline solution, remove the sheath with caution so as not to damage the balloon. In case of resistance, twist the sheath with one hand while holding the shaft with the other.

Caution: Do not use with Lipiodol™ or Ethiodol™ contrast media, or other such contrast media that incorporate the components of these agents.

TO DISPLACE AIR:

- Flush out the lumen for the guidewire by connecting a syringe filled with saline solution. Start flushing until liquid rinses out of guidewire lumen at the distal tip of the balloon catheter.
- Avoid manipulation of the balloon catheter during removal from packaging and flushing of guidewire lumen.

8. INSERTION TECHNIQUE

- Access to the vessel must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention.
- In case the guidewire has not already been positioned across the lesion, under fluoroscopy, insert a 0.035 in (0.89 mm) guidewire across the lesion, following standard PTA techniques.
- Insert the proximal end of the guidewire into the distal tip of the catheter.
- Carefully insert the PTA balloon catheter through the introducer sheath. Under fluoroscopic guidance, carefully advance it through the peripheral vasculature following the guidewire toward the lesion.

Caution: If resistance is met while advancing the balloon catheter, determine the cause of resistance and take remedial action before further advancing the guidewire or the balloon catheter.

- Position the balloon within the lesion using the 2 radiopaque markers on the shaft as reference points for precise placement across the target lesion.
- To avoid kinking, advance the dilatation catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.

9. Balloon inflation

Caution: Do not use air or any gaseous substances to inflate the balloon.

- Connect a syringe containing contrast medium to the Luer connector located at the proximal hub of the balloon catheter. Apply negative pressure for about 15 seconds until no bubbles appear in the contrast medium solution. Return to atmospheric pressure allowing contrast medium flow into the catheter lumen. Remove the syringe, leaving a meniscus of contrast in the hub of the balloon lumen.
- Remove air bubbles from the inflation device following the manufacturer's instructions.
- Using the stopcock, attach an inflation device to the PTA balloon catheter. Do not allow air to enter the system.
- Open the stopcock on the inflation device. Inflate the dilation catheter gradually to dilate the lesion to the calculated diameter in accordance with the compliance chart. Apply a constant pressure for about 30 seconds.

Caution: Do not exceed Rated Burst Pressure (RBP)

10. Balloon deflation and removal

Deflate the balloon in accordance with standard PTA procedures. Apply negative pressure to the balloon for between 30 and 120 seconds until the balloon is empty. Ensure the balloon is fully deflated before carefully pulling the catheter out of the vessel. A contrast mix of 70% saline and 30% contrast media is recommended with balloon lengths ≥ 200 mm in order to reduce deflation times.

After use, this product may be a biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

11. How supplied

The Admiral Xtreme PTA balloon catheter is supplied sterile and intended for single use only. The Admiral Xtreme is sterilized by ethylene oxide gas. It will remain sterile as long as the packaging remains unopened and undamaged. Use the product before the labeled Use-by date.

The Admiral Xtreme PTA catheter does not contain natural rubber latex. However, during the manufacturing, assembly, or packaging process, it may have incidental contact with materials, tools, or equipment containing or coming in contact with latex.

Caution: Do not use if the inner package is open or damaged.

12. Storage

Store at controlled room temperature in a dry place. Keep away from sunlight. Do not expose to organic solvents (such as alcohol), ionizing radiation, or ultraviolet light. Rotate inventory so that catheters are used prior to the expiration date on package label.

13. WARRANTY / LIABILITY

The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in these instructions for use are expressly to be considered as an integral part of this provision. Invatec guarantees the product until the expiry date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. Invatec disclaims any warranty of merchantability or fitness for a particular purpose of the product. Invatec is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on Invatec's part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and compliance, and excludes any other possible liability of Medtronic, however originating, from the product supplied. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects Invatec's legitimate interest in limiting its liability or warranty. No person has any authority to bind Invatec to any warranty or liability regarding the product.

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STATE OF MINNESOTA


COUNTY OF RAMSEY

I, Danielle Jean Rasmussen, Notary for the State of Minnesota, hereby certify that the attached Cover letter for Appendix to the Instructions for Use, Admiral Xtreme is a true copy of the original document filed at Medtronic, Inc.

Dated: 20th day of December 2021

(stamp)

Signature of Notary Public

Title: Notary Public
My commission expires: January 31, 2023

Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879

To Whom It May Concern

It is hereby certified that the document attached is the Appendix to the Instructions for Use for the Admiral Xtreme PTA Balloon Catheter.

Yours sincerely,

Sara Bakker

Digitally signed by Sara Bakker
Date: 2021.12.17 14:02:24
-06'00'

Signed on behalf of Medtronic Inc.

Name: Sara Bakker

Title: Regulatory Affairs Manager

Date

Appendix to the Instruction for use Admiral Xtreme PTA Balloon Catheter

PRESCRIPTION OF MEDICAL DEVICE

for percutaneous transluminal angioplasty (PTA) of peripheral arteries.

INFORMATION ABOUT STICKER

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Warning and informing symbols and inscriptions (if applicable)

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STERILIZATION

1. This device is single-use.
2. Device is sterile and non-pyrogenic.
3. Sterilized by EO.
4. Sterility Assurance Level (SAL)- 10^{-6}
5. Term preservation of sterility during 3 years.

CONDITIONS OF TRANSPORTATION AND STORAGE

Storage Conditions:

Keep dry and away from sunlight. Store at controlled room temperature (from +18°C to +25°C). Do not expose to organic solvents (e.g. alcohol), ionizing radiation or ultraviolet light. Rotate inventory so that catheters are used according to labeled "use by" date.

Transportation:

Care should be taken when transporting the devices. Packed by the manufacturer devices for peripheral vascular stenting are transported by all modes of transport in closed vehicles in accordance with the rules of carriage of goods, operating in each mode of transport.

Temperature of transporting: -18°C up to 55°C, humidity from 30 to 85 %.

OPERATING CONDITIONS

Do not autoclave. Do not use if package is opened or damaged. Use at controlled room temperature.

The product is intended for use by trained medical personnel in the operating room: at room temperature (18°C) not higher than 25°C, relative humidity of 75%, a pressure of from 860 to 1060 kPa.

SHELF LIFE

The shelf life of the device is 3 years.

PROCEDURE AND CONDITIONS FOR THE DISPOSAL OR DESTRUCTION OF A MEDICAL DEVICE

After use, this product may be a potential biohazard (hazard class B). Handle and dispose of all such device in accordance with accepted medical practice and applicable local, regional and national laws and regulations.

THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS, WHICH MATCHES MEDICAL DEVICE

Directive 93/42/EEC, EN ISO 13485, ISO 10993-1, ISO 10993-4, ISO 10993-5, ISO 10993-10, ISO 10993-11, ISO 10993-12, ISO 11135, ISO 10993-7, EN 556-1, EN ISO 11737-2, ANSI/AAMI ST72, ISO 11607-1, ISO 11607-2, ASTM D 4169, ASTM F 1980, EN 1041, ISO 15223-1, MEDDEV.2.7.1 Rev 4, ISO 14971, ISO 10555-1, ISO 10555-4, EN 62366-1, EN 1707, ISO 80369-7, ISO 80369-20.

CUSTOMER SERVICE

Contact for all questions, including complaints of medical devices the Authorized Representative in Russia:

Medtronic, LLC,

123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41

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© 2017 Medtronic