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Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG., hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: 
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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Created November 2016.

INTENDED USE

The ARCHITECT B-R-A-H-M-S PCT Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT iSystem when used for the quantitative determination of prolactin (PCT) in human serum and plasma.

Refer to the reagent package insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

2 x 3 Bottles (3.0 mL each) of ARCHITECT B-R-A-H-M-S PCT Controls. Low Control, Medium Control, and High Control contain recombinant PCT prepared in phosphate buffer. Preservatives: ProClin 950 and sodium azide.

The controls are at the following ranges and approximate concentrations:

Control	Target Concentration (ng/mL or µg/L)	Control Range (ng/mL or µg/L)
CONTROL L	0.20	0.14–0.26
CONTROL M	2.00	1.38–2.62
CONTROL H	70.00	42.00–98.00

Each laboratory should establish its own concentration ranges for new control lots at each control level.

Note: The insert ranges for the controls are not lot specific and represent the total range of values which may be generated throughout the life of the product. It is recommended that each laboratory establish its own means and acceptable ranges which should fall within the package insert ranges. Sources of variation that can be expected include:

- Calibration
- Calibrator lot
- Control lot
- Reagent lot
- Instrument

PRECAUTIONS

- IVD
- For In Vitro Diagnostic Use
- Rx ONLY

The following warnings and precautions apply to: CONTROL L / CONTROL M / CONTROL H	
WARNING	Contains methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

STORAGE

- Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- Controls must be stored frozen at -10°C or colder. After thaw, store at 2-8°C for up to 30 days.

PREPARATION FOR ANALYSIS

Thaw controls at room temperature (15–30°C) until completely thawed (30–60 minutes).

Prior to each use, mix by gentle inversion (10 times).

It is suggested to record the thaw date on the carton or the bottles.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL: L	Control Low, Medium, High (L,M,H)
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
PRODUCED FOR ABBOTT BY	Produced for Abbott by
PRODUCT OF USA	Product of USA
REF	List Number
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).
WARNING: SENSITIZER	Warning: May cause an allergic reaction.

B-R-A-H-M-S PCT is a trademark and property of Thermo Fisher Scientific Inc. and its subsidiaries.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

 Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY
Fisher Diagnostics,
A Div. of Fisher Scientific Company, LLC
A Part of Thermo Fisher Scientific, Inc.
8365 Valley Pike
Middletown VA 22645 USA

**Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com**

Created November 2016.
©2016 Abbott Laboratories

to the Ins

SHIPPING
Recomm
Medical
After the
record.
Device
have b
carrier.

WARN
Manu
life of
requi

MAIN
Not

REF
Not

DE
De
"S
Tr

L
IS
IS
E
I

ADDENDUM

to the Instruction for Use ARCHITECT B-R-A-H-M-S Controls for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for device are -10°C or lower.

Medical device must be stored in freezers to ensure regulated temperature conditions with a daily temperature record.

After thaw, must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. The maximum shelf life of the medical device is 9 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPiN 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

ISO 9001: 2008

ISO 13485: 2003 (including CMDCAS)

EN ISO 13485:2012

EN ISO 14971

EN 13641

EN 13612

EN ISO 15223-1

EN 980

EN ISO 18113-1

EN ISO 18113-2

EN ISO 18113-3

EN ISO 23640

EN ISO 20776-1

EN 13975

EN ISO 17511

ISO 80000-1

POTENTIAL USER

Professional technician.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs.

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage.

TECHNICAL CHARACTERISTICS

- pH-value at 20 °C - 6,9
- Density at 20 °C - 1,01 g/cm³

MATERIALS OF ANIMAL/HUMAN ORIGIN

Materials of animal/human origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available www.abbottiagnostics.com or contact your local representative.

VIAL TO VIAL VARIATION

Vial to vial variation parameter falls into overall uncertainty specification limit for test and is:

≤ 7% for the range > 0.1 ng/mL to ≤ 45 ng/mL

≤ 10% for the range > 45 ng/mL to ≤ 100 ng/mL

SAFETY PRECAUTIONS

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

QUALITY CONTROL PROCEDURES

- The recommended control requirement for the ARCHITECT B·R·A·H·M·S PCT assay is that a single sample of each control level be tested:
- Once every 24 hours each day of use
- After performing calibration
- After instrument service procedures or maintenance that may affect assay performance have been performed. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.
- Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.
- Each laboratory should establish control ranges to monitor the acceptable performance of the assay. If a control is out of its specified range, the associated sample results are invalid and the samples must be retested. Recalibration may be indicated.
- To establish statistically-based control limits, each laboratory should establish its own concentration target and ranges for new control lots at each control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days and using the reported results to establish the expected average (target) and variability about this average (ranges) for the laboratory. Sources of variation that should be included in this study to be representative of future system performance include:
 - Multiple stored calibrations
 - Multiple reagent lots
 - Multiple calibrator lots
 - Multiple processing modules
- These results should be applied to your laboratory's quality control practices. In addition, the laboratory must ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.
- Refer to Clinical and Laboratory Standards Institute (CLSI) Document C24-A3 or other published guidelines for general quality control recommendations.

MATERIALS REQUIRED BUT NOT PROVIDED

- ARCHITECT B·R·A·H·M·S PCT Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbottiagnostics.com.
- 6P22-01 ARCHITECT B·R·A·H·M·S PCT Calibrators
- 6P22-25/35 ARCHITECT B·R·A·H·M·S PCT Reagent Kit
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

INFORMATION ON THE MULTIPLICITY OF DEFROST

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