
STATE OF MINNESOTA


COUNTY OF HENNEPIN

I, Ashley Sophia Peterson, Notary for the State of Minnesota, hereby certify that the attached Cover Letter for Appendix to the IFUs-Amplatz Goose Neck Snares and MicroSnares is a true and accurate copy of an original document in possession of Medtronic, Inc.

Dated: May 17, 2023

(stamp)

Signature of Notary Public

Title: Notary Public
My commission expires: January 31, 2027

Your endovascular company.

4600 Nathan Lane N
Plymouth, MN 55442
USA
www.medtronic.com

Date: 17 May 2023

To Whom It May Concern


It is hereby certified that the document attached contains the Appendix to the Instructions for Use for the "Amplatz Goose Neck Snare kits extraction and manipulation, versions".

Signed on behalf of ev3, Inc.

Rebecca McCormick

Associate International Regulatory Affairs Specialist

Date

Appendix to the Instructions for Use

Amplatz Goose Neck Snare kits extraction and manipulation, versions

(When the device is delivered to the Russian Federation, the information contained in this Appendix to the Instructions for use is a priority)

NAME OF MEDICAL DEVICE

Amplatz Goose Neck Snare kits extraction and manipulation, versions:

1. Amplatz Goose Neck Snare Kit, models: GN 500, GN1000, GN1001, GN 1500, GN2000, GN2500, GN2501, GN3000, GN 3500, as part of:

- Snare with loop – 1 pc.
- Snare catheter-1 pc.
- Snare introducer – 1 pc.
- Torque device – 1 pc.
- Snare catheter, models: MC4000, MC4001, MC6000, MC6001 – 1 pc. (if necessary);
- Instructions for use

2. Amplatz Goose Neck Microsnare Kit, models: SK200, SK 201, SK 400, SK 401, SK700, SK 701, as part of:

- Microsnare with loop-1 pc.
- Microcatheter-1 pc.
- Microsnare introducer-1 pc.
- Torque device-1 pc.
- Instructions for use

PRESCRIPTION OF MEDICAL DEVICE

The Amplatz Goose Neck snare is intended for use in the cardiovascular system or hollow viscus to retrieve and manipulate foreign objects.

The Amplatz Goose Neck microsnare kit is intended for use in the retrieval and manipulation of atraumatic foreign bodies located in the coronary and peripheral cardiovascular system and the extra-cranial neurovascular anatomy.

INFORMATION ABOUT STICKER

Information on Russian Sticker:

The final packaging for the Russian market is accompanied with additional sticker:

Devices supplied in packaging marked with the following possible information:

Name of medical device in Russian language is written in accordance with the registration certificate of the medical device;

Name of the manufacturer;

Country of origin;

Code UPN (unique product number);

Code CFN (customer facing number);

Date and number of Registration Certificate;

Sign of conformity in Russia;

Warning and informing symbols and inscriptions (if applicable).

Appendix to the Instruction for use
Amplatz Goose Neck Snare kits extraction and manipulation, versions

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Appendix to the Instruction for use

Amplatz Goose Neck Snare kits extraction and manipulation, versions

Introducer	Aids in the insertion of the snare into the catheter
Torque device (pin vise)	Aids in the manipulation of the snare once placed in the vasculature

Snare Introducer, mm	5, 10	15, 20, 25, 30, 35
Inner diameter, in (mm)	0.040 (1,02)	0.063 (1,60)
Torque Device		
Compatibility with wire, in (mm)	0.038 (0,96)	0.052 (1,32)

MicroSnare Introducer	
inner diameter, in (mm)	0.022 (0,56)
Torque Device	
Compatibility with wire, in (mm)	0.018 (0,46)

Torque Device	
Inner diameter, in (mm)	0.106 (2,69)
Compatibility with wire, in (mm)	0.014 (0,36) Minimum wire for functional requirement for tightening on wire

STERILIZATION

Amplatz Goose Neck Snare kits extraction and manipulation are provided sterile for single use. They are sterilized using an ethylene oxide method that results in a sterility assurance level (SAL) minimum of 10^{-6} . The devices are provided to the user labeled sterile and single-use. The devices are not intended for sterilization by the user.

CONDITIONS OF TRANSPORTATION AND STORAGE

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Conditions of transportation and storage:

Temperature	from +10°C to +35°C
Humidity	30-70 %
Pressure	70-106 kPa

OPERATING CONDITIONS

The implantable medical devices are intended for use into a human body.

Temperature	from +32°C to +42°C
Humidity	100 %
Pressure	70-106 kPa

SHELF LIFE

The shelf life of the Amplatz Goose Neck Snare Kit, Amplatz Goose Neck Microsnare Kit - 5 years from the date of sterilization.

The shelf life of the Snare catheter, models: MC4000, MC4001, MC6000, MC6001 – 3 years from the date of sterilization.

INDICATIONS FOR USE

Snare Kit

Appendix to the Instruction for use

Amplatz Goose Neck Snare kits extraction and manipulation, versions

The Amplatz Goose Neck snare is intended for use in the cardiovascular system or hollow viscus* to retrieve and manipulate foreign objects.

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MicroSnare Kit

The Amplatz Goose Neck microsnare kit is intended for use in the retrieval and manipulation of atraumatic foreign bodies located in the coronary and peripheral cardiovascular system and the extra-cranial neurovascular anatomy.

Snare catheter

The Snare catheter is intended for use in the cardiovascular system or hollow viscus* to retrieve and manipulate foreign objects.

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* The hollow viscera include the digestive system, specifically the esophagus, stomach, small intestine, and large intestine; the renal system, specifically the kidneys, ureters, bladder and urethra; the biliary system, specifically the gallbladder and bile ducts; and the heart/central circulatory system.

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Do not use substances that are incompatible with any of the product's components (see «Product components section» Instruction for Use medical device).

Microsnare Kit

This device is not intended for the removal of foreign bodies entrapped by tissue growth.

Do not use substances that are incompatible with any of the product's components (see «Product components section» Instruction for Use medical device).

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Do not use substances that are incompatible with any of the product's components (see «Product components section» Instruction for Use medical device).

PROCEDURE AND CONDITIONS FOR THE DISPOSAL OR DESTRUCTION OF A MEDICAL DEVICE

After use, this product may be a potential biohazard. Dispose the medical device in accordance with the requirements of the national legislation on environment and health and/or hospital biohazard requirements.

Devices should be utilized according to SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of

Appendix to the Instruction for use

Amplatz Goose Neck Snare kits extraction and manipulation, versions

industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures)".

Classical medical waste:

Medical waste class for unused device: A.

Medical waste class for used device: B.

THE LIST OF INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS, WHICH MATCHES MEDICAL DEVICE

EN ISO 13485, EN ISO 14971, EN ISO 11607-1, EN ISO 11607-2, EN 1041, ISO 15223-1, EN ISO 10993-7, ISO 11135, ISO 11138-1, EN ISO 11138-2, ISO 11737-1, ISO 14644-1, ISO 14644-2, ANSI/AAMI ST72, EN ISO 10993-1, EN ISO 10993-4, EN ISO 10555-1, EN ISO 20594-1.

CUSTOMER SERVICE

Contact for all questions, including complaints, repair and service of medical devices the Authorized Representative in Russia:

Medtronic, LLC,

123112, Russia, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III, room 41.

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


To Whom It May Concern

It is hereby certified that the document attached is the "Instructions for Use" (510600-001, 510601-001 for the Amplatz Goose Neck extraction and manipulation kits).

Yours sincerely,

Signed on behalf of Ev3, Inc.
Ryan Zenanko
Associate Regulatory Affairs Specialist

30 NOV 2021
Date

Mary Kathryn Stagni
11/30/21

Amplatz Goose Neck™

Snare Kit / Snaresæt / Schlingenset / Kit βρόχου / Kit de lazo / Kit de lasso / Kit del dispositivo di recupero a cappio / Teugelset / Snaresett / Zestaw pętli / Conjunto de ansa / Набор ловушки / Snarsats / Kement Seti

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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2. Indications for use

The Amplatz Goose Neck snare is intended for use in the cardiovascular system or hollow viscous to retrieve and manipulate foreign objects.

Manipulation procedures include indwelling venous catheter repositioning, indwelling venous catheter fibrin sheath stripping, and central venous access venipuncture procedure assistance.

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4. Potential adverse events

Potential adverse events associated with the use of the snare include, but are not limited to, those listed in section 4.1 to section 4.4.

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- Inspect the device by extending and retracting the snare loop through the distal end of the catheter 2 or 3 times, while carefully examining the catheter, the radiopaque marker band, and the snare for defects or damage.

- Use a guiding catheter or guidewire to insert the snare and snare catheter together.

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9.4. Procedure

9.4.1. Snare assisted manipulation or retrieval

- Remove the indwelling balloon or delivery catheter, if present. It may be necessary to change or extend the indwelling guidewire to facilitate balloon removal and to upsize the indwelling guide catheter to accommodate the snare.
- Position the snare just proximal to the foreign body with one of the methods described in section 9.2 or section 9.3.
- Gently push the snare shaft forward within the snare catheter to move the snare forward and completely open the loop. Advance the loop forward, around the proximal end of the foreign body.
- While maintaining the position of the snare shaft, advance the snare catheter to cinch the loop of the snare around the foreign body.
Note: Attempting to close the loop by pulling the snare shaft within the snare catheter moves the loop from its position around the foreign body.
- To manipulate a foreign body, maintain tension of the snare catheter on the snare loop to retain the hold on the foreign body, and move the snare and snare catheter together to manipulate the foreign body to the desired position.
- To retrieve a foreign body, maintain tension of the snare catheter on the snare loop and move the snare and snare catheter together proximally to or into the guide catheter or sheath. Then withdraw the foreign body through or together with the guiding catheter or vascular sheath. To withdraw large foreign bodies, if necessary insert larger sheaths or guiding catheters or perform a cut down at the peripheral site.

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- Insert a snare and, from a common femoral vein approach, insert the snare and snare catheter with one of the methods described in section 9.2 or section 9.3. Position the snare in the vasculature at the desired venous access site.
- Open the snare loop by advancing the snare shaft

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Instructions for use

1. Device description

The Amplatz Goose Neck™ snare kit contains 1 Amplatz Goose Neck snare and 1 Amplatz Goose Neck snare catheter. The snare is made of a nitinol cable and a gold-plated tungsten loop. The preformed snare loop can be introduced through catheters without risk of snare deformation because of the snare's superelastic construction. The snare catheter contains a platinum-iridium, radiopaque marker band.

2. Indications for use

The Amplatz Goose Neck snare is intended for use in the cardiovascular system or hollow viscous to retrieve and manipulate foreign objects. Manipulation procedures include indwelling venous catheter repositioning, indwelling venous catheter fibrin sheath stripping, and central venous access venipuncture procedure assistance.

3. Contraindications

- This device is not intended for removing foreign bodies that have become entrapped by tissue growth.
- This device should not be used for fibrin sheath stripping in the presence of atrial or ventricular septal defects or persistent foramen ovale.
- This device is not intended for removal of implanted pacing leads.

Do not use substances that are incompatible with any of the product's components. (See Product components section.)

4. Potential adverse events

Potential adverse events associated with the use of the snare include, but are not limited to, those listed in from section 4.1 to section 4.4.

4.1. Use of a snare retrieval device in arterial vasculature

- Embolization
- Myocardial infarction (depending on placement)
- Stroke

4.2. Stripping fibrin sheath from an indwelling venous catheter

- Indwelling catheter damage on small-bore catheters
- Pulmonary embolization

4.3. Use of a snare retrieval device in venous vasculature

- Pulmonary embolization

4.4. Other potential adverse events associated with snare retrieval devices

- Device entrapment
- Infection
- Lengthened procedure
- Vessel damage
- Vessel perforation

5. Warnings

- This device is supplied sterile, for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing can increase the risks of patient infection and compromised device performance.
- Do not use excessive pull force when attempting fibrin sheath stripping of indwelling catheters less than or equal to 2.0 mm (6 Fr) in diameter to avoid the risk of damaging, stretching, or breaking them.
- To minimize risk of pulmonary embolism, use caution when removing a large fibrin sleeve with this device.
- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.

6. Product components

Each Amplatz Goose Neck snare kit contains 1 snare and 1 snare catheter. The snare is packaged with a snare introducer and a torque device. Table 1 lists the components that may contact skin or blood.

Table 1: Product components

Snare	Nitinol, gold-plated tungsten, polytetrafluoroethylene, adhesive
Snare introducer	Low-density polyethylene, high-density polyethylene, barium sulfate
Snare catheter	High-density polyethylene, platinum-iridium marker band
Torque device	Aluminum anodized to a yellow-gold color, polysulfone

7. How supplied

Warning: This device is supplied sterile, for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing can increase the risks of patient infection and compromised device performance.

8. Storage conditions

Store away from sunlight at room temperature in a dry place.

9. Directions for use

9.1. Preparation

- Select the appropriate size snare for the location of the foreign body. In general, the snare loop should approximate the size of the vessel in which the snare will be used.
- Remove the snare and snare catheter from their protective holders and inspect for any damage.

9.2. To insert the snare and snare catheter together as a unit

- Remove the introducer and torque device from the proximal end of the snare shaft.
- Load the snare into the snare catheter by inserting the proximal (nonlooped) end of the snare into the distal (nonhubbed) end of the snare catheter, until the proximal end of the snare shaft exits the hub and the loop end can be retracted into the distal end of the snare catheter.

- Inspect the device by extending and retracting the snare loop through the distal end of the catheter 2 or 3 times, while carefully examining the catheter, the radiopaque marker band, and the snare for defects or damage.
- Use a guiding catheter or guidewire to insert the snare and snare catheter together.

9.2.1. With an indwelling guide catheter

- Pull the snare just within the distal end of the snare catheter and advance both together through the indwelling guide catheter until the distal end of the snare catheter is positioned as desired.

9.2.2. With a guidewire

- Place the loop of the snare over the proximal end of the wire and cinch the snare loop over the wire by advancing the snare catheter over the snare shaft.
- Advance the snare and the snare catheter together over the wire and into the guide catheter or sheath until the distal end of the snare catheter is positioned as desired.

9.3. To insert the snare after the snare catheter is positioned within the vasculature

- Move the snare introducer (provided on the proximal end of the snare) distally until the loop of the snare is enclosed within the tubing portion of the introducer.
- Insert the distal end of the introducer into the hub of the indwelling snare catheter until resistance is felt. Resistance indicates that the tip of the introducer is properly aligned with the inner lumen of the catheter.
- While keeping the introducer as straight as possible, grasp the shaft of the snare just behind the hub of the introducer and advance the snare until it is secure within the lumen of the catheter. To remove the introducer, first remove the torque device from the snare shaft and then pull the introducer off the proximal end of the snare shaft.

9.4. Procedure

9.4.1. Snare assisted manipulation or retrieval

- Remove the indwelling balloon or delivery catheter, if present. It may be necessary to change or extend the indwelling guidewire to facilitate balloon removal and to upsize the indwelling guide catheter to accommodate the snare.
- Position the snare just proximal to the foreign body with one of the methods described in section 9.2 or section 9.3.
- Gently push the snare shaft forward within the snare catheter to move the snare forward and completely open the loop. Advance the loop forward, around the proximal end of the foreign body.
- While maintaining the position of the snare shaft, advance the snare catheter to cinch the loop of the snare around the foreign body.

- Note:** Attempting to close the loop by pulling the snare shaft within the snare catheter moves the loop from its position around the foreign body.
- To manipulate a foreign body, maintain tension of the snare catheter on the snare loop to retain the hold on the foreign body, and move the snare and snare catheter together to manipulate the foreign body to the desired position.
 - To retrieve a foreign body, maintain tension of the snare catheter on the snare loop and move the snare and snare catheter together proximally to or into the guide catheter or sheath. Then withdraw the foreign body through or together with the guiding catheter or vascular sheath. To withdraw large foreign bodies, if necessary insert larger sheaths or guiding catheters or perform a cut down at the peripheral site.

9.4.2. Snare assisted removal of fibrin sheaths from indwelling catheters

- Select a snare and, from a common femoral vein approach, insert the snare and snare catheter with one of the methods described in section 9.2 or section 9.3. Advance the snare and the snare catheter together to the inferior vena cava or right atrium.
- Advance a 0.89 mm (0.035 in) guidewire through the port (the distal or venous port if more than one lumen) of the indwelling catheter and into the right atrium or inferior vena cava.
- Encircle the guidewire with the snare loop.
- Advance the snare loop over the distal end of the indwelling catheter to a position proximal to the fibrin sheath.
- Close the snare loop around the indwelling catheter by advancing the snare catheter while maintaining the position of the snare shaft. Apply light traction continuously while gently pulling the snare down toward the distal end of the indwelling catheter, over the port or ports.
- Repeat steps 4 and 5 until the indwelling catheter is free of fibrin sheath.

9.4.3. Snare assisted venous cannulization

- At a patient venous access site, introduce the snare and snare catheter with one of the methods described in section 9.2 or section 9.3. Position the snare in the vasculature at the desired venous access site.
- Open the snare loop by advancing the snare shaft through the snare catheter to provide a target to guide a needle puncture into the desired venous access site.
- Introduce a guidewire through the needle and through the snare loop.
- Cinch the snare loop over the needle by advancing the snare catheter over the snare shaft while maintaining the position of the snare shaft.
- While maintaining tension of the snare catheter on the snare loop, move the snare shaft and the snare catheter together and pull the guidewire into the inferior vena cava.

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

STATE OF MINNESOTA

COUNTY OF RAMSEY

I, Danielle Jean Rasmussen, Notary for the State of Minnesota, hereby certify that the attached Cover letter for Appendix for IFUs Signed, Snares is a true copy of the original document filed at Medtronic, Inc.

Dated: 22nd day of December 2021

Signature of Notary Public


Title: Notary Public
My commission expires: January 31, 2023

(stamp)

To Whom It May Concern

It is hereby certified that the document attached is the Appendix to the Instructions for use 510600-001, 510601-001) for the Amplatz Goose Neck extraction and manipulation kits.

Yours sincerely,

Digitally signed by Ryan
Zenanko
Date: 2021.12.20 09:59:13
+06'00'

Signed on behalf of ev3, Inc.
Ryan Zenanko
Associate Regulatory Affairs Specialist

20Dec2021
Date

Appendix to the Instructions for use Amplatz Goose Neck extraction and manipulation kits

PRESCRIPTION OF MEDICAL DEVICE

This is intended for manipulating and extracting foreign objects.

INFORMATION ABOUT STICKER

The final packaging for the Russian market is accompanied with additional sticker.
Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language is written in accordance with the registration certificate of the medical device;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Sign of conformity in Russia;
- Warning and informing symbols and inscriptions (if applicable).

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STERILIZATION

The Amplatz Goose Neck Snare Kit and Amplatz Goose Neck Microsnare Kit and Snare catheter are provided sterile for single use.

They are sterilized using an ethylene oxide method that results in a sterility assurance level (SAL) minimum of 10^{-6} . The devices are provided to the user labeled sterile and single-use. The devices are not intended for sterilization by the user.

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Appendix to the Instruction for use 510600-001 Rev C, 510601-001 Rev B
Amplatz Goose Neck extraction and manipulation kits

Pressure	70 - 106 kPa
----------	--------------

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Snare catheter - 3 years from the date of sterilization.

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Devices should be utilized according to SanPiN

Classical medical waste:

Medical waste class for unused device: A

Medical waste class for used device: B

THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS, WHICH MATCHES MEDICAL DEVICE

EN ISO 13485:2016, EN ISO 14971:2012, EN ISO 11607-1:2009 ISO 11607-1:2006/Amd-1:2014, EN ISO 11607-2:2019, EN 1041:2008+A1:2013, ISO 15223-1:2016 Corr 2017, EN ISO 10993-7:2008/AC:2009, ISO 11135:2014 +A1:2019 CORR 2020, ISO 11138-1:2017, EN ISO 11138-2:2017, ISO 11737-1:2018, ISO 14644-1:2015, ISO 14644-2:2015, ANSI/AAMI ST72:2011, EN ISO 10993-1:2020, EN ISO 10993-4:20017, EN ISO 10555-1:2013+A1:2017, EN ISO 20594-1:1993/AC:1996

CUSTOMER SERVICE

Contact for all questions, including complaints of medical devices the Authorized Representative in Russia:

Medtronic, LLC,

123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41

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