COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, Michelle Mirzaee, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

Phenom IFU

Description of documents

Is a true and exact copy of the correct and complete original document.

Michelle Ball

Signature of document custodian (affiant)


A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 10th day of November,
2021 by Michelle Mirzaee —, proved
to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Steve Oberlander
Signature (Seal)

Medtronic

Phenom™ Delivery Catheter

fr	Cathéter de pose Phenom™
de	Phenom™ Abgabekatheter
it	Catetere di rilascio Phenom™
es	Catéter de liberación Phenom™
sv	Phenom™ införingskateter
nl	Phenom™-plaatsingskatheter
pt-pt	Cateter introdutor Phenom™
fi	Phenom™-asetuskatetri
da	Phenom™-indføringskateter
el	Καθετήρας τοποθέτησης Phenom™
cs	Zaváděcí katétr Phenom™
hu	Phenom™ bevezetőkatéter
ru	Доставочный катетер Phenom™
pl	Cewnik wprowadzający Phenom™
tr	Phenom™ Taşıma Kateteri
no	Phenom™-innføringskateter
sk	Zavádzací katéter Phenom™
ro	Cateter de poziționare Phenom™
bg	Доставящ катетър Phenom™
et	Sisestuskateeter Phenom™
hr	Kateter za uvođenje Phenom™
lt	„Phenom™“ įterpimo kateteris
lv	Phenom™ ievadkatetrs
sl	Dovajalni kateter Phenom™

TABLE OF CONTENTS

. *henom*[™] Delivery Catheter

English.....	3
Français.....	5
Deutsch.....	7
Italiano.....	9
Español.....	11
Svenska.....	13
Nederlands.....	15
Português.....	17
Suomi.....	19
Dansk.....	21
Ελληνικά.....	23
Česky.....	25
Magyar.....	27
Русский.....	29
Polski.....	31
Türkçe.....	33
Norsk.....	35
Slovenčina.....	37
Română.....	39
Български.....	41
Eesti.....	43
Hrvatski.....	45
Lietuvių.....	47
Latviski.....	49
Slovenščina.....	51
Symbol Glossary.....	53

English Instructions for Use

en

Phenom™ Delivery Catheter

CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.
This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures.

DESCRIPTION

The Phenom™ Delivery Catheters are variable stiffness, single lumen catheters designed to access small, tortuous vasculature. They are available in a variety of lengths, stiffness and inner and outer diameters. The outer surface of the catheter is coated to enhance navigation in the vessel. The Phenom™ 17, 21, and 27 Delivery Catheters have a hydrophilic coating that spans the distal 100cm. The Phenom™ Plus Delivery Catheters have a hydrophilic coating that spans the distal 90cm. The catheter also incorporates a liner to facilitate movement of introduction devices passing through its lumen. The distal tip has radiopaque marker(s) to aid visualization and positioning under fluoroscopy. The Phenom™ Delivery Catheter is packaged with a shaping mandrel and may be accompanied with a split introducer sheath.

DEVICE FEATURE AND COMPATIBILITY

Phenom™ 17 Delivery Catheter			
	French	Inch	mm
Proximal OD	2.2	0.029	0.74
Distal OD	1.8	0.024	0.61
Proximal ID	-	0.017	0.43
Distal ID	-	0.017	0.43
Guiding Catheter or Catheter ID	-	≥ 0.035	≥ 0.89
Guidewire OD	-	≤ 0.014	≤ 0.36

Phenom™ 21 Delivery Catheter			
	French	Inch	mm
Proximal OD	2.6	0.034	0.86
Distal OD	2.3	0.030	0.76
Proximal ID	-	0.021	0.53
Distal ID	-	0.021	0.53
Guiding Catheter or Catheter ID	-	≥ 0.038	≥ 0.97
Guidewire OD	-	≤ 0.018	≤ 0.46

Phenom™ 27 Delivery Catheter			
	French	Inch	mm
Proximal OD	3.1	0.040	1.02
Distal OD	2.8	0.036	0.91
Proximal ID	-	0.027	0.69
Distal ID	-	0.027	0.69
Guiding Catheter or Catheter ID	-	≥ 0.045	≥ 1.13
Guidewire OD	-	≤ 0.025	≤ 0.64

Phenom™ Plus Delivery Catheter			
	French	Inch	mm
Proximal OD	4.7	0.061	1.55
Distal OD	4.2	0.055	1.40
Proximal ID	-	0.045	1.13
Distal ID	-	0.045	1.13
Guiding Catheter or Catheter ID	-	≥ 0.070	≥ 1.78
Guidewire OD	-	≤ 0.041	≤ 1.04

INTENDED PURPOSE / INDICATIONS FOR USE

Phenom™ Delivery Catheters are intended for the introduction of interventional devices or diagnostic agents into the neuro, peripheral, and coronary vasculatures.

CONTRAINDICATIONS

There are no known contraindications.

DIRECTIONS FOR USE

1. The Phenom™ Delivery Catheters may be packaged in a protective coil dispenser fitted with a flushing luer connector or in a packaging tray. The catheters have a hydrophilic coating and require hydration prior to use.
2. For catheters packaged in protective coil dispenser, flush the dispenser with saline/heparinized saline through the luer fitting connected to the end of the coil dispenser. Repeat if friction is felt when attempting to remove the catheter.
3. For catheters packaged in a tray, hydrate the catheter in the tray with saline/heparinized saline.

CAUTIONS: Do not attempt to use catheters without flushing or hydrating first with saline/heparinized saline. Failure to do so, may compromise the coating and lubricity of the catheter.

4. After removing catheter from coil dispenser or tray, inspect to verify that it is undamaged.

WARNING

- Do not use catheters that are damaged in any way. Damaged catheters may break or rupture causing vessel damage or tip detachment during the procedure.

5. Shaping mandrel:

- Remove shaping mandrel from mandrel card and insert into distal tip of the catheter.
- Bend the catheter tip and shaping mandrel to desired shape.
- Hold shaping mandrel and catheter tip 1" from a steam source for approximately 10 seconds to set shape. Remove from heat source and allow cooling in either air or saline prior to removing the mandrel.
- Discard shaping mandrel.

WARNING

- Shaping mandrel is not for in-vivo use.

6. Prior to use, flush the lumen of the catheter with saline/heparinized saline solution by attaching a saline/heparinized saline filled syringe to the catheter hub.
7. Carefully insert an appropriate guidewire into the hub and advance into the catheter lumen.
8. Place an appropriate guiding catheter using standard percutaneous technique. Attach a rotating hemostasis valve to the guiding catheter luer connector and maintain a continuous flush.
9. A split introducer sheath is included with Phenom™ 17 models to aid in insertion into hemostatic valve. To use, slide the split introducer sheath from the proximal end until it covers the distal end of the catheter.
10. Introduce the catheter and guidewire into the guiding catheter through the hemostasis adaptor. Tighten the valve around the catheter to prevent backflow but allowing movement of the catheter through the valve.
11. Advance the catheter and guidewire to the selected target site by alternately advancing the guidewire and then tracking the catheter over the guidewire.
12. If using split introducer, slide the introducer back to the proximal end of the catheter next to the hub.
13. Withdraw the guidewire from the catheter. Connect the desired syringe with infusate or interventional device and continue the procedure according to the manufacturer's instructions.
14. Withdraw the catheter and discard after completing the procedure. Follow your institutional procedures for disposal of biohazards.

NOTE: Flow rates for given pressures and infusates are provided below as a reference and are approximate values.

Flow Rate Tables				
Phenom™ 17 Delivery Catheter				
Length (cm)	Pressure (psi)	100% 76% Ionic Contrast (cc/sec)	50/50 76% Ionic Contrast/Saline (cc/sec)	100% Saline (cc/sec)
150	100	0.03	0.2	0.4
	300	0.1	0.5	1.0

Phenom™ 21 Delivery Catheter				
Length (cm)	Pressure (psi)	100% 76% Ionic Contrast (cc/sec)	50/50 76% Ionic Contrast/Saline (cc/sec)	100% Saline (cc/sec)
150	100	0.1	0.4	1.0
	300	0.2	1.0	1.8
160	100	0.1	0.4	0.9
	300	0.2	0.9	1.7

Phenom™ 27 Delivery Catheter				
Length (cm)	Pressure (psi)	100% 76% Ionic Contrast (cc/sec)	50/50 76% Ionic Contrast/Saline (cc/sec)	100% Saline (cc/sec)
150	100	0.2	1.0	1.9
	300	0.6	2.5	3.0
160	100	0.1	0.9	1.7
	300	0.5	2.4	3.0

Phenom™ Plus Delivery Catheter				
Length (cm)	Pressure (psi)	100% 76% Ionic Contrast (cc/sec)	50/50 76% Ionic Contrast/Saline (cc/sec)	100% Saline (cc/sec)
120	100	1.5	7.4	6.8
	300	4.4	12	12

This product does not contain raw materials made with DEHP or BPA.

POTENTIAL COMPLICATIONS

Potential complications of the devices and the endovascular procedure include or are synonymous with, but may not be limited to the following:

- Pain at insertion site
- Hematoma
- Hemorrhage
- Infection
- Sensitization or an allergic reaction
- Complication of radiation exposure (e.g., alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia).
- Distal emboli (air, foreign material or thrombus)
- Vessel spasm, dissection, perforation or injury
- Ischemia
- Stroke that could lead to Neurological deficit and death

* Consult instructions for use for other therapy devices and medications for additional potential complication information. If a serious incident related to the device occurs, contact your Medtronic representative and the competent authority in your respective country/region.

WARNINGS

- THIS DEVICE IS INTENDED FOR SINGLE USE ONLY. Discard after one procedure. Reuse and/or re-sterilization can potentially result in compromised device performance and cross contamination.
- The Phenom™ Delivery Catheter should be manipulated under fluoroscopy only. Do not attempt to move the catheter without observing the resultant tip response.
- Never advance or withdraw an intraluminal device against resistance until the cause of resistance is determined by fluoroscopy. If the cause cannot be determined, withdraw the catheter. Movement of the catheter against resistance may result in damage to the catheter, or the vessel.
- The infusion pressure should not exceed the listed pressure as indicated in Flow Rate Table (above).
- Nickel (CAS 7440-02-0) is present in the Phenom™ Delivery Catheter and may cause sensitization or an allergic reaction.
- The catheter has not been tested for use with automated high-pressure contrast injection equipment, do not use this equipment with the device because it may damage the device.
- Visually inspect all sterile barrier systems, that are labeled as sterile, immediately prior to use. Do not use the device if breaches in sterile barrier system integrity are evident.
- For additional Materials of Concerns information such as REACH, CA Prop 65 or other product stewardship programs, go to www.medtronic.com/productsstewardship

PRECAUTIONS

1. Inspect the Phenom™ Delivery Catheter prior to use for any kinks or bends. Any catheter damage may decrease the desired performance characteristics.
2. Do not use if package is open or damaged. Contents of unopened, undamaged package are sterile.
3. Use prior to the "Use By" Date.
4. Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.
5. Phenom™ Delivery Catheter is sterilized with ethylene oxide gas. Do not re-sterilize.
6. Phenom™ Delivery Catheter is non-pyrogenic.
7. Operators should take all necessary precautions to limit x-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

STORAGE AND DISPOSAL

- This device should be stored in a dry place, away from sunlight.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

SAFETY AND CLINICAL PERFORMANCE

A copy of the Summary of Safety and Clinical Performance can be viewed by searching the device name on the EUDAMED website <https://ec.europa.eu/tools/eudamed>

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CALIFORNIA ALL-PURPOSE ACKNOWLEDGMENT

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State of California

County of Los Angeles

On 03/25/2022 before me, Steve Oberlander, Notary Public

Date

Here Insert Name and Title of the Officer

personally appeared


Nishitaben Patel

Name(s) of Signer(s)

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature

Steve Oberlander
Signature of Notary Public

Place Notary Seal Above

OPTIONAL

Though this section is optional, completing this information can deter alteration of the document or fraudulent reattachment of this form to an unintended document.

Description of Attached Document

Title or Type of Document: Appendix to IFU for Phenom

Document Date: _____ Number of Pages: _____

Signer(s) Other Than Named Above: _____

Capacity(ies) Claimed by Signer(s)

Signer's Name: _____

☐ Corporate Officer — Title(s): _____

☐ Partner — ☐ Limited ☐ General

☐ Individual ☐ Attorney in Fact

☐ Trustee ☐ Guardian or Conservator

☐ Other: _____

Signer Is Representing: _____

Signer's Name: _____

☐ Corporate Officer — Title(s): _____

☐ Partner — ☐ Limited ☐ General

☐ Individual ☐ Attorney in Fact

☐ Trustee ☐ Guardian or Conservator

☐ Other: _____

Signer Is Representing: _____

Medtronic

9775 Toledo Way

Irvine, CA 92618

USA

Tel: 949-837-3700

www.medtronic.com

M. S. Patel, 03/25/2022

Nishitaben Patel
Regulatory Affairs Associate

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© 2021 Medtronic.
Medtronic, Medtronic with logo and Medtronic logo are U.S. and internationally registered trademarks of Medtronic.