
MILESTONE S.r.l.
Via Fatebenefratelli, 1/5
24010 SORISOLE (BG)
P. IVA 01879330163

MILESTONE S.r.l.
Via Fatebenefratelli, 1/5
24010 SORISOLE (BG)
P. IVA 01879330163

MILESTONE
HELPING
PATIENTS

PrestoCHILL
Flash Chiller for Frozen Sections

Operator Manual
MM138

Thank you for having selected our system and welcome to the ever growing world club of users for Milestone laboratory instrumentation.

We are sure that you will be completely satisfied with this new tool entering your laboratory.

We invite you to read carefully this operator manual and to keep it in reach for convenient and fast consulting.

For any possible clarification or any request for assistance please contact either our Representative in your country:

For the Russian Federation the authorized representative:

BioLine LLC

Pinsky lane 3 letter A,

197101, Saint Petersburg,

the Russian Federation

E-mail: service@bioline.ru

Tel: (812) 320-49-49

Fax: (812) 320-49-40

Service Hotline: 8-800-333-00-49.

Or contact

Milestone s.r.l.

Via Fatebenefratelli, 1/5

24010 Sorisole (BG) Italy

Tel. +39.035.412 8264

Fax +39.035.575498

web site www.milestonemedsrl.com

e-mail marketing@milestonemedsrl.com

Please read the user manual before using the device.

SUMMARY

SUMMARY

1. INTRODUCTION

- 1.1. Symbol used
- 1.2. Biological and hygiene warnings
- 1.3. Designated Uses
- 1.4. The operating principle
- 1.5. Information about consumers and intended users
- 1.6. Indications and contraindications to the use, side effects
- 1.7. The type of tested specimens
- 1.8. Technical specifications
 - 1.8.1. Touch control terminal
- 1.9. Environmental conditions
- 1.10. Transportation and storage conditions
- 1.11. Warning information
- 1.12. Labeling explanation
- 1.13. Compatible reagents
- 1.14. Decommissioning and waste disposal
 - 1.14.1. Disposal of equipment
 - 1.14.2. Disposal of filter
 - 1.14.3. Disposal of frozen blocks
 - 1.14.4. Disposal of condensed liquid

2. INSTALLATION

- 2.1. Space Requirements
- 2.2. Complete set
- 2.3. Moulds
 - 2.3.1. Standard different sizes moulds (if required)
 - 2.3.2. Mould for freezing in cryomould (if required)
 - 2.3.3. Mould for 55mm chuck (if required)
- 2.4. Materials or equipment required for testing, but not included in the set of supply
- 2.5. Electrical installation
 - 2.5.1. Power supply
 - 2.5.2. Back side connection plate
- 2.6. Final placement
- 2.7. Unit overview
- 2.8. Plate conformation

3. OPERATE WITH PrestoCHILL

- 3.1. Set Date and Time
- 3.2. Help
- 3.3. Events log
- 3.4. Start up MANUAL MODE
- 3.5. Start up AUTOMATIC MODE
 - 3.5.1. Set automatic programs during the week
- 3.6. Set up freezing parameters
- 3.7. Defrost
 - 3.7.1. HEPA Filter
- 3.8. Freezing a sample for Frozen Section
 - 3.8.1. Freezing time
 - 3.8.2. Using multiple timers
- 3.9. Safe switching of

4. MAINTENANCE

- 4.1. Introduction
- 4.2. Before/After each run
- 4.3. Daily
- 4.4. Every 250 defrost pump hours
 - 4.4.1. How to change the HEPA filter

5. DECONTAMINATION AND CLEANING

- 5.1. Safety precautions
- 5.2. PrestoCHILL decontamination and cleaning procedure
 - 5.2.1. Detergents
 - 5.2.2. Procedure

6. Warranty

- A. Appendix – Tips and Tricks
 - A.1. Fatty tissues
 - A.2. Muscle tissue
 - A.3. Small Biopsies
- B. Appendix – Freezing time
 - B.1. Recommended freezing time for moulds
 - B.2. Freezing time changes
- C. Appendix – Use the extractor for chuck
 - C.1. How to use the extractor for chucks
- D. Appendix – How to change a mould
 - D.1. How to change a mould
- E. Appendix – Troubleshooting
 - E.1. Block detachment from the chuck
 - E.2. Power failure
- F. Appendix – How to manage an alarm
 - F.1. How to manage an alarm
- G. Appendix – Information about Milestone Cryoembedding Compound (MCC)
- H. Appendix – Marking of PrestoCHILL and accessories
 - H.1. Shield of instrument
 - H.2. Label of PrestoCHILL package
 - H.3. Labels and marking of accessories
- I. Appendix – The list of national standards applied by the manufacturer

1. INTRODUCTION

1.1. Symbol used

An instruction accompanied by this symbol provides important information and requires more attention.

An instruction accompanied by this symbol provides a cautionary statement: failure to follow the instruction may endanger the operator or cause damage to the instrument.

An instruction accompanied by this symbol provides a caution against electric shock. Be sure to follow the instruction to prevent accidents resulting from electric shock.

Biohazard.

Be careful: parts identified with this symbol are zones with risks of crushing hand and/or finger(s).

Surfaces identified with this symbol are at danger temperature. Be careful. Cold.

Be careful. Hot surface.

CE logo: this instrument complies to European Community directives.

IVD medical device according to 98/79/EC directive.

Catalog number

Lot number

Serial number

Switched on

Switched off

Alternating current

Nominal voltage, V

Maximum nominal power, W

USB slot

Best before

The content is sufficient for n-number of tests.

For single use only

WEEE European directive symbol (2002/96/EC): Electric/electronic device do not throw out in the environment.

Manufacturer.

Fragile.

1.2. Biological and hygiene warnings

Please, carefully read the following note.

This instrument is an in-vitro medical device and it is installed in a location where there is presence of blood, specimen and other biological tissues. For your safety, it is required to clean and disinfect the device before entering in contact with it.

Also you should wear gloves when you touch it.

In addition, it is required that the unit is cleaned and disinfected before sending it to MILESTONE S.R.L.

In any case, international health rules forbid shipments in standard packages of biohazard materials (risk of sanctions).

1.3. Designated Uses

Presto CHILL is an excellent choice for supporting tissue during cryotomy.

PrestoCHILL is specifically designed to produce frozen blocks of samples of biological section or frozen tissues in a fast, reliable and standardized way, before their microtoming, as part of intraoperative histological study.

Milestone PrestoCHILL device and its components are designed for use in laboratories and therefore able to withstand the most severe laboratory conditions associated with continuous use when freezing the sections of biological samples.

The device is stationary, intended **only for laboratory use** (indoor use). For in vitro diagnostic.

Milestone PrestoCHILL device is to be used with reagents listed in the "Compatible Reagents" of Chapter 1.13.

Any other use of the instrument is considered improper and may void the warranty provided by the manufacturer.

1.4. The operating principle

The operating principle of a medical product: the PrestoCHILL allows producing frozen tissue blocks with a standardized protocol using the "face down" cryoembedding technique.

Standardized freezing method: the PrestoCHILL freezing way is based on the "face down" technique. The sample without pre-treatment is transferred into the cold mould with the interesting side in contact with the mould base. Then the mould is filled with medium and covered with the relative chuck, in order to create a frozen block to be cut in the cryostat. The "face down" technique allows simultaneously the fast freezing of sample and its embedding into the block.

1.5. Information about consumers and intended users

Only for the professional users (Pathologist or medical laboratory technician (medical laboratory assistant)).

The PrestoCHILL may be operated by trained laboratory personnel only.

All laboratory personnel designated to operate this instrument must read the Operator Manual carefully and must be familiar with all technical features of the instrument before attempting to operate it.

1.6. Indications and contraindications to the use, side effects

Indications for use: producing of frozen blocks of samples of biological section or frozen tissues (for example, fat, muscle tissue, small biopsies and other) in a fast, reliable and standardized way, before their microtoming, as part of intraoperative histological study.

Contraindications: not applicable.

Side effects: not applicable.

1.7. The type of tested specimens

Tested samples: specimens of biopsies and other biological samples/tissues.

1.8. Technical specifications

TYPE OF INSTRUMENT: Tabletop immovable flash chiller for high-quality frozen sections

MODEL: PrestoCHILL

REF 100850 (230V~ 50/60Hz)

CLASSIFICATION: In-vitro medical device professional use only (European directive 98/79/EC)

SAFETY MARKS: CE – IVD

The Quality System Management of Milestone s.r.l. is ISO 9001:2015 and ISO13485:2016 certified.

TECHNICAL CHARACTERISTICS:

- Refrigerant: Helium
- The acceptable time for setting the operating mode (the time is required to reach the set temperature (-40°C) from room temperature (approximate 20°C)) varies from 40 minutes to 1 hour.
- Use temperature range: -24°C to -60°C, tolerances (± 2) °C
- Setting increment of operating temperature: 1 °C
- Range of values of the set freezing time: from 1 sec to 999 sec
- Setting increment of freezing time: 1 sec
- Freezing time for each mold and heat removal unit: see Appendix B.1.
- Device defrost time: depends on the set operating temperature. Standard defrost program at working temperature of -40 °C takes 3 hours, plus 1 hour to reach the defrost temperature;
- Number of moulds placed in the working area: 6 pcs.
- Force required for lifting the cover: max. 22 N
- Corrected sound power level of noise: not more than 63 dB
- Sound power level of alarm: not more than 63 dB
- Device operation modes: manual and automatic.
- Multiple timer function: 6 independently adjustable timers
- **DIMENSIONS:** W 300mm x D 540mm x H 450mm / 730mm with open cover (with open cover) (± 5) %
- **WEIGHT:** 25kg (± 10) %
- **POWER SUPPLY:** single phase 230V~ 50/60Hz, Power Input 250W – 854 BTU/h Max, tolerances (± 10) %
- **FUSES:** two main fuses 1.6A, curve T, size 5x20, voltage 250 V
- Maximum altitude: 2000m

Product classification in accordance with the Russian standard GOST R 50444-92:


Class of potential failure effects during the using – class V for the device.

Group of perceived mechanical actions – group 1 for the device.

Enclosure protection degree (according to GOST 14254-2015) – IP20.


1.8.1. Touch control terminal

The following picture shows the Touch Control TERMINAL of the PrestoCHILL unit. It is possible to operate the unit by simply touching the icons on the screen. The inscriptions on the touch terminal are in English.

	[Settings]
 Start chiller	[Start chiller]
 Start defrost	[Start defrost]
 AUTO CYCLE	[AUTO CYCLE]
Automatic mode	[Automatic mode]
	[Home page]

One USB PORT is located on the lateral sides of the Touch Control Terminal.

When not in use, it's required to leave the USB protection inserted into the slot, as shown in the picture above.

- Characteristics:
- TFT-LCD 4.3" widescreen color display
Resolution 480X272,
Resistive Touchscreen
Mass storage 512MByte
1CPU 700MHz
- Software:
- Windows™ Compact Embedded 6.0
- Programs:
- Customizable time and temperature
- I/O-Interface:
- 1x USB port 2.0

Software security class in accordance with the Russian standard GOST R IEC 62304-2013: Class A: No injury or damage to health is possible.

Since the installed software is only the controlling software/ firmware of the equipment. No damage or injury is possible. There is no likelihood of causing harm by destruction. Software makes no clinical decisions.

1.9. Environmental conditions

Operating conditions on the request of the manufacturer.

WORKING TEMPERATURE: 15 °C to 25 °C

RELATIVE HUMIDITY: up to 80%

Ambient pressure: (760 ± 30) mm Hg

Indoor use only.

1.10. Transportation and storage conditions

Transport and storage conditions on the request of the manufacturer.

Temperature: -20°C up to +60°C

Humidity: up to 80% (up to 50% at 40°C) non-condensing

Transportation by all types of roofed transport. Storage in a dry places.

Before turning PrestoCHILL on (when coming from a storage room), let the unit reach working environmental conditions (at least half an hour).

It is very important to ship and move the PrestoCHILL system in a vertical position to avoid that the cooling plate could be damaged and unscrewed from its support.

1.11. Warning information

Should the equipment be used in different ways than those specified by the manufacturer, the warranty provided by the manufacturer may be void.

Main power plugs are considered disconnecting device. Disconnect all plugs from the socket outlet before assembling of the system, before connection of the accessories (if present) and before cleaning.

Equipment shall be connected to a main power socket outlet which is accessible and visible by the operator. Before starting work, while the device is in a safe state, it is necessary to visually make sure that the device is evenly placed on a horizontal surface, and check the absence of external damage to the device and power cable.

Equipment is provided with replaceable fuses, however their replacement is to be performed by properly authorized service personnel. If the system is not properly working, please contact the local authorized and trained service representative of the manufacturer.

Parts of this instrument may attain temperatures as low as -42°C. It is important to avoid allowing bare skin to touch such cold surfaces – when in doubt, wear gloves.

The equipment conforms to the requirements for noise resistance and electromagnetic emission, and does not affect other equipment.

It is recommended to place the device in places where the intensity of external electromagnetic radiation and electrostatic discharge is below the maximum allowable intensity.

Do not use this equipment near sources of strong electromagnetic radiation, ionizing radiation, or electrostatic discharges, as they may disrupt its proper operation.

When connecting to the instrument of items or equipment which is not included in the delivery, the electromagnetic emission exceeding the permissible level may be occurs.

It is recommended to use the device in the specified environmental conditions (see chapter 1.9.) in order to avoid adverse effects on its operation.

HEPA filter and residues of frozen blocks of sections or tissues should be treated as a biohazardous material after use. It should be disposed in accordance with federal, municipal and local government regulations (see chapters 1.14.2 and 1.14.3).

All parts of the equipment and all accessories must be supplied only by the manufacturer.

DO NOT OPEN THE CHASSIS. Failure to follow the instruction in the documentation might lead to a reduction of device safety. If the equipment is not working properly, please contact your supplier or the manufacturer.

Do not remove the instrument, accessories, panels or cover. Only authorized and qualified **SERVICE** personnel may repair the instrument and access the instrument's internal components.

It is important that normal standards of safety and good laboratory practices are employed. Always use common sense and the best known practice when operating the instrument.

The institution owning the instrument has primary responsibility for accident-free operation together with designated personnel who operate, service or repair it.

To avoid damage to the instrument, use only listed in chapter 1.13 Compatible Reagents. Failure to follow the instructions in the documentation might lead to a reduction of device safety.

Low Temperatures are present in this equipment.

Extreme care should be taken.

DO NOT let bare skin come into contact with metal surfaces.

When using equipment with violations of the requirements established by the manufacturer, the level of equipment protection may be reduced.

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For symbols explanation see chapter 1.1

1.13. Compatible reagents

1. Embedding medium: gel consisting of polyethylene glycol and polyvinyl alcohol which when frozen has the same density as frozen tissue.
More detailed information of the embedding medium composition is shown in Appendix G.
2. Alcohol-based detergents and polyphenolic based disinfectants to clean internal work plan and PrestoCHILL accessories.

1.14. Decommissioning and waste disposal

1.14.1. Disposal of equipment

This instrument is an in-vitro medical device and it is usually installed in a laboratory where specimens and other biological tissues are present. For your safety, it is therefore required to clean and disinfect before entering in contact with it. You must also wear gloves when operating the system.

In case of return it is required that the unit is cleaned and disinfected before sending it back to Milestone. International health rules require that shipments of biohazard materials are not done in standard packages (risk of sanctions).

Do not remove the side panels of the unit during the cleaning operations.

**Before doing any operation disconnect the power supply.
Do not use direct or high-pressure water to clean the instrument.**

For safety reason, any operation without the individual protective devices must be avoided.

It is recommended to use alcohol or polyphenolic based disinfectants to clean and disinfect the instrument. For further information, please contact the manufacturer at application@milestonemedsrl.com.

Removal of equipment from use for disposal:

- 1) Switch OFF the unit. Disconnect the power supply.
- 2) Disconnect all the accessories, if available.
- 3) Clean and disinfect all parts of the equipment according to the manufacturer specification.
- 4) Dispose all electric and electronic parts of the equipment.

Complying with directive 2002/96/EC of the European Parliament and of the Council of 27 January 2003 on waste electrical and electronic equipment (WEEE), the separate collection environmental managed of equipment is mandatory.

The in-observance of Directive 2002/96/EC or of the local law which acknowledge it can have potential effects on the environment and human health.

**This symbol indicates separate collection
for electrical and electronic equipment**

In case of installation outside the European Community territory, the instrument and/or its part after disinfecting must be disposed in accordance with the existing applicable, local regulations. For the Russian Federation it is recommended to be disposed in accordance with the rules implemented by "SanPin 2.1.7.2790-10 (Sanitary and epidemiological requirements for the medical waste handling)" for the medical waste referred to class A (non-hazardous waste), excluding electronic parts. Electronic parts must be disposed in accordance with local specified regulations.

1.14.2. Disposal of filter

HEPA filter (Code 100624) must be treated as a biohazard material after its use.

The used filter must be disposed in accordance with federal, municipal and local government regulations. For the Russian Federation, it is recommended to dispose in accordance with the rules established by "SanPiN 2.1.7.2790-10 (Sanitary and epidemiological requirements for medical waste handling)" for medical waste related to Class B (epidemiologically hazardous waste).

1.14.3. Disposal of frozen blocks

Residues of frozen blocks of sections or tissues after use must be disposed in an approved incinerator in accordance with federal, municipal and local government regulations.

For the Russian Federation, it is recommended to dispose in accordance with the rules established by "SanPiN 2.1.7.2790-10 (Sanitary and epidemiological requirements for medical waste handling)" for medical waste related to Class B (epidemiologically hazardous waste).

1.14.4. Disposal of condensed liquid

Condensed liquid is water formed during the defrosting of the working plate (defrosting process is described in section 3.7).

For the Russian Federation, it is recommended to dispose in accordance with the rules established by "SanPiN 2.1.7.2790-10 (Sanitary and epidemiological requirements for medical waste handling)" for medical waste related to Class A (epidemiologically hazardous waste).

2. INSTALLATION

Milestone PrestoCHILL is an instrument which requires the utmost care when it is unpacked and installed.


Provide a stable level surface of at least 70cm in width and 70cm in depth and that can stand the PrestoCHILL weight of $25\text{kg} \pm 10\%$.

At least 2 people should remove the unit from its box and locate it on the bench.

Then place PrestoCHILL on the bench, make sure that the unit is level and the bench is made of non-flammable material and is vibration free. Ensure that the instrument has been positioned away from hot, direct sunlight and is in a location completely free from drafts.

Make sure the unit is not placed against any side surface limiting the air cooling inlet/exhaust: allow at least a distance of 20cm for each side and of 30cm behind.

2.1. Space Requirements

PrestoCHILL is delivered in a **wood box** with these **dimensions**:

Width: 420mm $\pm 5\%$	Height: 680mm $\pm 5\%$	Depth: 650mm $\pm 5\%$
------------------------	-------------------------	------------------------

Unit dimensions:

Width: 300mm $\pm 5\%$	Height: 450mm/ 730 mm (with open cover) $\pm 5\%$	Depth: 540mm $\pm 5\%$
------------------------	--	------------------------

Space dimensions required to locate the unit:

Width: 700mm	Height: 800mm	Depth: 700mm
--------------	---------------	--------------

2.2. Complete set


When unpacking, check that all parts correspond to the packing list included in the shipment. Utmost care must be taken during unpacking of PrestoCHILL, to avoid scratching or damage of the external coating.

The PrestoCHILL unit (230V version **Code 100850**) may include the following items:

- Qty.1 Main power cable (Code 50036 for 230V version)

- Qty.1 Magnetic pen to work with the touch-screen terminal Code 50596/S
- Qty.1 USB – pen with Operator manual
USB storage capacity - 8 Gb.
- Qty: 1 pcs. Operation manual on paper
- Dispensing slides:
 - Qty. 3 Code 100551, large tip (if required)
 - Qty. 3 Code 100545, narrow tip (if required)

- Qty.1 Key for the changing of the moulds (Code 100575) (if required)


- Tweezers:
 - Qty.1 Code 100593: the small one (if required)
 - Qty.1 Code 100741: the big one (if required)


- Qty.2 Heat extractor standard (Code 100533) (if required)


- Qty.4 Heat extractor slim (Code 100677) (if required)

- Qty: 1 pcs. Heat extractor rectangular (Code 100592) (if required)

- Qty: 1 pcs. Heat extractor with bigger hole (Code 100736) (if required)

- Qty.1 Chuck Extractor (Code 100737L) (if required)

- Qty: 5 Box with 100 paper discs, diameter 16 mm (**Code 100623**) (if required)
- Qty: 5 Box with 100 paper discs, diameter 22mm (**Code 100622**) (if required)
- Qty: 5 Box with 100 paper discs, diameter 30 mm (**Code 100615**) (if required)

Paper discs are for single use only.

- Qty.2 Kit HEPA filter (**Code 100624**) (if required)

- Qty.12 Milestone Cryoembedding Compound (MCC) 125ml (Code 51420) (if required)

- Qty: 6 pcs. Moulds (if required), see chapter 2.3.







2.3. Moulds









We advise against the use of accessories other than those provided by Milestone.




2.3.1. Standard different sizes moulds (if required)




1	MLD3228 2926-5 rectangular shape chuck of 26 x 31mm depth well - 5mm well length – 29 mm well width – 26 mm max sample size 26,5 x 23,5mm outer diameter – 50mm total height – 19,5 mm weight – 0,062 kg Code 100590	
---	---	--





2	<u>MLD3228 2926-7</u> rectangular shape chuck of 26 x 31mm depth well - 7mm well length – 29 mm well width – 26 mm max sample size 26 x 23mm outer diameter – 50mm total height – 19.5 mm weight – 0.059 kg Code 100694	
3	<u>MLD3030 2525-5</u> square shape chuck of 28 x 28mm depth well - 5mm well length – 25 mm well width – 25 mm max sample size 22.5 x 22.5mm outer diameter – 50mm total height – 19.5 mm weight – 0.063 kg Code 100693	
4	<u>MLD3030 2525-7</u> square shape chuck of 28 x 28mm depth well - 7mm well length – 25 mm well width – 25 mm max sample size 22 x 22mm outer diameter – 50mm total height – 19.5 mm weight – 0.060 kg Code 100698	
5	<u>MLD3030 2525-10</u> square shape chuck of 28 x 28mm depth well - 10mm well length – 25 mm well width – 25 mm max sample size 21 x 21mm outer diameter – 50mm total height – 19.5 mm weight – 0.056 kg Code 100699	
6	<u>MLD3939 3333-5</u> square shape chuck of 37 x 37mm depth well - 5mm well length – 33 mm well width – 33 mm max sample size 31,5 x 31,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.054 kg Code 100692	
7	<u>MLD3939 3333-7</u> square shape chuck of 37 x 37mm depth well -7mm well length – 33 mm well width – 33 mm max sample size 31 x 31mm outer diameter – 50mm total height – 19.5 mm	


	weight – 0.048 kg Code 100696	
8	<u>MLD3939 3333-10</u> square shape chuck of 37 x 37mm depth well - 10mm well length – 33 mm well width – 33 mm max sample size 30 x 30mm outer diameter – 50 mm total height – 19.5 mm weight – 0.040 kg Code 100697	
9	<u>MLD4129 3324-5</u> rectangular shape chuck of 40 x 28mm depth well - 5mm well length – 33 mm well width – 24 mm max sample size 30,5 x 21,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.058 kg Code 100666	
10	<u>MLD4129 3324-7</u> rectangular shape chuck of 40 x 28mm depth well - 7mm well length – 33 mm well width – 24 mm max sample size 30 x 21mm outer diameter – 50mm total height – 19.5 mm weight – 0.055 kg Code 100700	
11	<u>MLD4129 3324-10</u> rectangular shape chuck of 40 x 28mm depth well - 10mm well length – 33 mm well width – 24 mm max sample size 29 x 20mm outer diameter – 50mm total height – 19.5 mm weight – 0.050 kg Code 100701	
12	<u>MLD4234-5</u> round shape chuck of 35mm depth well - 5mm well diameter – 34 mm max sample size 31,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.054 kg Code 100702	
13	<u>MLD4234-7</u> round shape chuck of 35mm depth well - 7mm well diameter – 34 mm max sample size 31mm outer diameter – 50mm total height – 19.5 mm	

	weight – 0.050 kg Code 100703	
14	MLD4234-10 round shape chuck of 35mm depth well - 10mm well diameter – 34 mm max sample size 30mm outer diameter – 50mm total height – 19.5 mm weight – 0.044 kg Code 100704	
15	MLD3022-5 round shape chuck of 25mm depth well-5mm well diameter – 22 mm max sample size 19,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.064 kg Code 100605-A	
16	MLD3022-7 round shape chuck of 25mm depth well -7mm well diameter – 22 mm max sample size 19mm outer diameter – 50mm total height – 19.5 mm weight – 0.062 kg Code 100705	
17	MLD3426-5 round shape chuck of 30mm depth well - 5mm well diameter – 26 mm max sample size 23,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.061 kg Code 100606-A	
18	MLD3426-7 round shape chuck of 30mm depth well - 7mm well diameter – 26 mm max sample size 23mm outer diameter – 50mm total height – 19.5 mm weight – 0.059 kg Code 100607-A	
19	MLD3426-10 round shape chuck of 30mm depth well - 10mm well diameter – 26 mm max sample size 22mm outer diameter – 50mm total height – 19.5 mm weight – 0.055 kg Code 100608-A	


20	<u>MLD4436-5</u> round shape chuck of 40mm depth well - 5mm well diameter – 36 mm max sample size 33,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.052 kg Code 100609-A	
21	<u>MLD4436-7</u> round shape chuck of 40mm depth well -7mm well diameter – 36 mm max sample size 33mm outer diameter – 50mm total height – 19.5 mm weight – 0.047 kg Code 100610-A	
22	<u>MLD4436-10</u> round shape chuck of 40mm depth well -10mm well diameter – 36 mm max sample size 32mm outer diameter – 50mm total height – 19.5 mm weight – 0.040 kg Code 100611-A	
23	<u>MLD3330 2419-5</u> rectangular shape chuck of 26 x 31mm depth well - 5mm well length – 24 mm well width – 19 mm max sample size 21,5 x 16,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.065 kg Code 100713	
24	<u>MLD3330 2419-7</u> rectangular shape chuck of 26 x 31mm depth well - 7mm well length – 24 mm well width – 19 mm max sample size 21 x 16mm outer diameter – 50mm total height – 19.5 mm weight – 0.063 kg Code 100714	
25	<u>MLD3330 2419-10</u> rectangular shape chuck of 26 x 31mm depth well - 10mm well length – 24 mm well width – 19 mm max sample size 20 x 15mm outer diameter – 50mm total height – 19.5 mm weight – 0.060 kg Code 100715	

26	<u>MLD4442 4436-5</u> oval shape round chuck, diameter 40mm depth well - 5mm well length – 44 mm well width – 36 mm max sample size 33,5 x 39,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.050 kg Code 100716	
27	<u>MLD4442 4436-7</u> oval shape round chuck, diameter 40mm depth well - 7mm well length – 44 mm well width – 36 mm max sample size 33 x 39mm outer diameter – 50mm total height – 19.5 mm weight – 0.044 kg Code 100717	
28	<u>MLD4442 4436-10</u> oval shape round chuck, diameter 40mm depth well - 10mm well length – 44 mm well width – 36 mm max sample size 32 x 38mm outer diameter – 50mm total height – 19.5 mm weight – 0.036 kg Code 100718	
29	<u>MLD4410-5</u> round shape chuck of diameter or diagonal 20 - 40mm depth well - 5mm well diameter – 10 mm max sample size 7,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.064 kg Code 100707	
30	<u>MLD4410-7</u> round shape chuck of diameter or diagonal 20 - 40mm depth well - 7mm well diameter – 10 mm max sample size 7mm outer diameter – 50mm total height – 19.5 mm weight – 0.064 kg Code 100708	
31	<u>MLD4410-10</u> round shape chuck of diameter or diagonal 20 - 40mm depth well - 10mm well diameter – 10 mm max sample size 6mm outer diameter – 50mm total height – 19.5 mm weight – 0.064 kg Code 100709	

32	<u>MLD4415-5</u> round shape chuck of diameter or diagonal 20 - 40mm depth well - 5mm well diameter – 15 mm max sample size 12,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.063 kg Code 100710	
33	<u>MLD4415-7</u> round shape chuck of diameter or diagonal 20 - 40mm depth well - 7mm well diameter – 15 mm max sample size 12mm outer diameter – 50mm total height – 19.5 mm weight – 0.063 kg Code 100711	
34	<u>MLD4415-10</u> round shape chuck of diameter or diagonal 20 - 40mm depth well - 10mm well diameter – 15 mm max sample size 11mm outer diameter – 50mm total height – 19.5 mm weight – 0.062 kg Code 100712	
35	<u>MLD2618-5</u> round shape chuck of 22mm depth well -5mm well diameter – 18 mm max sample size 15,5mm outer diameter – 50mm total height – 19.5 mm weight – 0.066 kg Code 100727	
36	<u>MLD2618-7</u> round shape chuck of 22mm depth well -7mm well diameter – 18 mm max sample size 15mm outer diameter – 50mm total height – 19.5 mm weight – 0.065 kg Code 100728	
37	<u>MLD4326 1818-5</u> square shape chuck of 25 x25mm depth well - 5mm well length – 18 mm well width – 18 mm max sample size 15.5 x 15.5mm outer diameter – 50mm total height – 19.5 mm	


	weight – 0.067 kg Code 100730	
38	MLD4326 1818-7 square shape chuck of 25 x 25mm depth well – 7 mm well length – 18 mm well width – 18 mm max sample size 15 x 15mm outer diameter – 50mm total height – 19.5 mm weight – 0.065 kg Code 100731	

2.3.2. Mould for freezing in cryomould (if required)


MLD50 2520-5 rectangular shape cryomould of 22 x20mm compatible with chuck of diameter or diagonal 20- 55mm depth well - 5mm well length – 25mm well width – 20 mm outer diameter – 50mm total height – 19.5 mm weight – 0.070 kg Code 100726	
--	---

Use of moulds for freezing in cryomould


1. Prepare the specimen in cryomould at room temperature and fill the latter with MCC.

2. Transfer to mould MLD50 2520-5 and add chuck and standard heat extractor. Freeze for 90 seconds and then check that the block is completely frozen.


3. Remove the cryomould from the block and cut in cryostat or detach the block from the chuck and store it in the freezer.

Freezing with cryomould is slower and could lower the end quality.

2.3.3. Mould for 55mm chuck (if required)

Usually these moulds are used for big specimens.

MLD5951-5 round shape chuck of 55mm depth well -5mm well diameter – 51 mm max sample size 48.5mm total height – 19.5 mm weight – 0.083 kg Code 100722	
MLD5951-7 round shape chuck of 55mm depth well -7mm well diameter – 51 mm max sample size 48mm total height – 19.5 mm weight – 0.073 kg Code 100723	

The tolerance: for all sizes of moulds (± 1.0 mm)
for the weight ($\pm 10\%$)


You can use a single mould MLD5951 in the working area of the PrestoCHILL.

If using a MLD5951 mould the PrestoCHILL conformation allows the presence of 3 moulds: one MLD5951 mould in position 4 and two standard sized moulds in position 1 and 5.

With this type of moulds it is not advisable to use the slim heat extractor. Always use the standard heat extractor, see next paragraph.

Avoid several consecutive freezing cycles with this type of mould.

2.4. Materials or equipment required for testing, but not included in the set of supply

Frozen blocks of biological tissue samples or sections obtained with PrestoCHILL can be sectionalized on any cryostat or rotary microtome available to the user, using a disk or cryostat holder.

Based on the size of the disc or cryostat holder available to the user, Milestone recommends choosing the appropriate molds from the list of molds.

The cryostat holder or disc is an integral part of any cryostat, serving as an auxiliary device to the PrestoCHILL, providing the convenience of transferring the frozen block to the cryostat for subsequent sectioning, and does not affect the ability to produce frozen tissue blocks or sections.

2.5. Electrical installation

2.5.1. Power supply

Power supply line must have a ground connection (yellow/green, green or bare wire).

Do not use the blue/white wire (neutral of power line). Power supply neutral has to be referred to ground.

Carefully check that:

- A Ground Connection is made available in the lab to assure the unit works correctly. Do not connect a unit without ground connection.
- Power supply line wires (their size and distance from the switch board) used in the laboratory are compliant with the label positioned on the back side of units.
- Plug/socket is visible during the normal use and easily accessible.

Power line protection:

- 230V: Miniature circuit breaker curve C, 6A, interruption more than 6000A.

Residual current device: 30mA (Class I device).

230V version: a schuko terminated power cord is provided.


Before connecting the plug to the laboratory's socket CHECK that the power supply voltage of the unit matches the power supply line of the laboratory. Damages occur when wrong power supplies the instrument.

2.5.2. Back side connection plate

On the back side of the PrestoCHILL are positioned:

- the HEPA filter
- the electrical plug
- the label with unit serial number
- the vent
- the spare fuses box


See also the label of the instrument that provides some important information. Always communicate the S.N. of the unit to Customer Support.

2.6. Final placement

PrestoCHILL unit can now be positioned in its final place of use. Make sure that the unit is level.

Then connect:

- the main power cable into the electrical plug that is located on the back of the unit (see the picture below),

- and the plug of main power cable into the wall socket.


- Now it is possible to turn the PrestoCHILL on.

2.7. Unit overview

1. Work area cover
2. Touch control terminal
3. Main switch
4. USB port
5. Embedding medium and spare accessories storage area
6. Cool trap


1. Heat extractors standard
2. Heat extractors slim
3. Moulds for 30mm chucks
4. Moulds for 40mm chucks

We advise against the use of accessories other than those provided by Milestone.

Connector and Switch Marking

On the side of the device


1. Main switch

I – switched on

O – switched off

2. USB port 

On the back panel

3. Connector for a power cable

- Warning of the potential electric shock. Care must be taken to avoid the effects of an electric shock.

2.8. Plate conformation

In the standard configuration in the work area of the PrestoCHILL houses two different types of heat extractors: standard and slim, and six 5mm deep moulds.

The standard heat extractor can be used with any type of mould among those listed in chapter 2.3.

The slim heat extractor can only be used with 5mm deep moulds.

The PrestoCHILL moulds can be replaced according to the laboratory requirements: size of frozen specimens and type of chuck used. To change the moulds see appendix D and for the list of mould models available see chapter 2.3.

The conformation used by the PrestoCHILL that ensures optimal performances is: four 5mm deep moulds and two moulds with depth of your choice.

If the PrestoCHILL conformation features moulds deeper than 5mm (7 or 10mm), we recommend positioning the slim heat extractors near the 5mm deep moulds and the standard heat extractors near the deeper moulds.

For optimal results, Milestone recommends using 5mm deep moulds.


Do not work with a different conformation to that indicated above.


3. OPERATE WITH PrestoCHILL

3.1. Set Date and Time

Press the  "Settings" icon from the status bar.


Then select  "Date/time" and set.

3.2. Help

In the  "Help" section in the "Settings" menu, it is possible to insert useful info: serial number (only for service), technical, application and sales assistance contacts and website.

3.3. Events log

Clicking on “Events log” icon , allows access to the PrestoCHILL actions list.


3.4. Start up MANUAL MODE

- 

The time required to reach the set temperature (-40°C) from room temperature (approximate 20°C) varies from 40 minutes to 1 hour.
- 

Before starting, wipe the work plate of the PrestoCHILL with paper or similar absorbent material to dry it.


Switch on the unit. This page will appear:


Touch the button “Start chiller” to set the working temperature and the sample freezing time and to initiate the cooling power. Software records the last set temperature and time data.

Set the required temperature with   in Temperature and Set areas.

Set the freezing time of the sample with   in Time area on the touch terminal.

The spot Status near the actual temperature value is red if the plate has not reached the set temperature. In this case it is better to not freeze any samples. When the actual temperature is $\pm 4^{\circ}\text{C}$ the set temperature, the spot becomes green and the plate is ready for freezing.

Push Start  , to start the freezing.

3.5. Start up AUTOMATIC MODE


The time required to reach the set temperature (-40°C) from room temperature (approximate 20°C) varies from 40 minutes to 1 hour.

Before starting, wipe the work plate of the PrestoCHILL with paper or similar absorbent material to dry it.

Switch on the unit. This page will appear:

Activate Automatic mode pushing the icon  to light:  on the main page. If you turn the unit on in a moment between Chiller start and Stop time, in 10 seconds the plate cooling starts and the unit will be usable when the temperature reaches the set one. If you turn the unit on in another moment, cooling will start at set time. It is always possible to modify the set temperature and the sample freezing time, by entering in the chiller page.


For safety reasons, if the cover of the PrestoCHILL is open when the beginning of a chiller program is set, the chiller program does not start.

3.5.1. Set automatic programs during the week

Automatic cycle is a function, which allows the user to decide the periods during the week in which the unit is at the working temperature, and the ideal moment during each day to do the unit defrost.


Press the "Automatic cycle" icon from the Settings page. It is necessary to insert the Administrator password (issued as a separate document inserted in the Operator Manual).


Edit the parameters on the base of laboratory needs:

Set Chiller time:

Select the day to set up.
Select the start chiller hour, by touching on it: the hour box will have a grey border. Press on the “Start chill” button, and select the precise starting time: at beginning, half or end of the selected hour. The selected hour will flash.


Now select the chilling stop time, touching on it: it will become highlighted. Press the “Stop chill” button to confirm the chiller ending time and choose the precise starting time: at beginning, half or end of the selected hour. The period in which chiller is on will be underlined with a blue line.

The chiller stop time can be selected also the day after the chiller start.
The chiller time can't be plus than 24h, it is preferable to do a defrost cycle once a day.

Set Defrost time:


Select the start defrost hour, by touching on it: it will have a grey border. Press on the "Start defrost" button, and select the precise starting time: at beginning, half or end of the selected hour. The selected hour will flash. The software automatically calculates and reserves the needed hours for the defrost program. A standard defrost program requires 3 hours plus 1 hour to reach the defrost temperature. The period reserved for the defrost cycle will be underlined with a yellow line.


Defrost cycle duration depends on Milestone default settings, accessible only to Service personnel. Contact application@milestonemedsrl.com for further information.




The software does not allow Chiller and Defrost programs in the same time. For example, if you try to start a defrost program during the chiller time you will have the following message:


Or if you start a defrost program that will end during a chiller program, the message will be:


To reset all the week programs, press on the "Reset program" button, while to cancel only one program, select it and press the "Reset selected" button, answering YES to the relative question.


Save settings, by clicking on YES to the exit question:

Activate the automatic mode by pressing the


icon to render it green (active):


in the main page.

Then the main screen shows icons for use in manual mode OFF (grey color):

During an automatic cycle, from the **Home page** if the current time is between the plate cooling start and end time, the chiller screen is automatically displayed after 10 seconds:

With automatic mode on, if the “Settings” page or a page other than the “Home page” is accessed, after 10 minutes of inactivity the unit returns to the home page in order to allow the implementation of the weekly program set.

On the left status bar, it is present the PrestoCHILL manual/automatic mode icon:


Automatic mode active

Automatic mode inactive
(Manual mode active)

3.6. Set up freezing parameters


IN MANUAL MODE: To set up freezing temperature and duration, enter in the Chiller page, through the "Start chiller" icon.

IN AUTOMATIC MODE: To set up freezing temperature and duration must be in the working chiller period (pre-set in automatic cycle settings- see chapter 3.5.1) otherwise you can not access the screen below.

When you are in the screen displayed below, modify the set temperature with the right arrows (red arrow to increase and blue arrow to decrease) or by touching the set temperature box (in this case a keypad will appear).

Modify the set freezing time in with the lateral arrow or by touching the Time box (in this case a keypad will appear).

Milestone suggests freezing samples at -40°C.

The freezing time depends on mould size used, refer to chapter 3.8.1.


3.7. Defrost

The defrost procedure allows to obtain, in a short time, the complete elimination of frost that forms during normal PrestoCHILL use.

Defrost procedure consists of plate heating, while an aspiration removes aqueous vapors derived from frost evaporation.

Since samples prepared in PrestoCHILL are fresh and so potentially dangerous, defrost vapors are filtered through an HEPA12 filter (blue arrow on the picture below), which removes biological pathogens.


Then vapors condense in the cool trap, positioned in the left angle of the PrestoCHILL frontal side, yellow arrow on the picture below.

Check the cool trap daily. If needed, empty it. Refer to chapter MAINTENANCE 4.3. It is strongly suggested to do a defrost cycle once a day.

When in manual mode, press the **“Start defrost”** icon on main screen to start defrost procedure. At the end of this procedure the plate returns to room temperature.

In automatic mode, defrost start is already set and it will not need to be initiated by the operator.

Defrost settings are under Service settings, accessible only to Service personnel. Contact application@milestonemedsrl.com for further information.

3.7.1. HEPA Filter

For safety reasons, between the working area and the defrost pump, there is a HEPA filter that filters possibly dangerous pathogens in the vapors.

Filter class H13 (in accordance with Russian standard GOST R EN 1822-1-2010), the filter retains 99.97% of all particles contained in the air with diameter $\geq 0.3 \mu\text{m}$, skip $\leq 0.05\%$, maximum pressure of 4.1 bar, filtration area of $1,300 \text{ cm}^2$, and nominal air flow rate of max. 3.9 l/min .

This filter must be changed after the defrost pump has worked for 250 hours.

The software advises you that the filter life is finished with the message “Filter needs to be replaced. Please do it”, at the starting of the first freezing program of the day.

Avoiding HEPA filter change, can be dangerous for operator safety. To change the filter, refer to chapter 4.4.1.


3.8. Freezing a sample for Frozen Section

Use PrestoCHILL unit to prepare blocks for frozen section.

- 1) Prepare the mould by cleaning with paper soaked in alcohol and placing the appropriate chuck on it (this procedure will prevent ice inside the mould).


Cleaning of the mould is very important to avoid the block from sticking in the mould.


- 2) Apply a thin layer of embedding medium (MCC) to the tip of the dispensing slide. This thin layer of MCC creates a clean plane of separation between the tissue and the mould floor and will prevent the sticking of any tissue residue.

- 3) Place tissue FACE DOWN on the dispensing slide. The transparency of the dispensing slide permits orientation of the tissue into the desired position.


- 4) Look through the underside of the slide while adjusting the tissue so desired margin or other aspect is visible in position.

- 5) Let sample protrude from the end of the slide and touch the leading edge of the tissue to the mould floor. It will adhere.

- 6) Pull the dispensing slide out from under the tissue. Tissue requiring precise positioning can be flattened or manipulated and guided into desired position by carrying out this step and carefully.

- 7) Fill the mould as quick as possible so that a meniscus of medium is bulging above the brim. Too much medium is not a problem. Too little medium can cause failure of the specimen to anchor to the chuck

- 8) Place the chuck over the mould as quick as possible and press it.

- 9) Place the heat extractor over the chuck stem. The heat extractor is a nickel-plated brass block designed to fit over the chucks. The extractor is maintained at a cold temperature to provide additional instant and rapid heat removal from the chuck, embedding medium, and specimen.

- 10) Close PrestoCHILL cover.

If the cover is left open for more than 3 minutes will trigger the alarm "cover open". Just close the cover to reset the alarm.


Alternatively press the icon

that appears in the status bar to the left of the screen, to silence the alarm.

- 11) Press "Start" button in "Time" screen area to begin the timer countdown.

- 12) Remove the frozen specimen and chuck when timer alarms. It will be flat and the surface will require less trimming and wastage of tissue before the entire tissue face is available for sectioning.

3.8.1. Freezing time

The freezing time depends on the size of the mould used and on the sample type.

The mould for 30 mm chuck provides a standard freezing time of 40 seconds. In case of delicate specimens which are cut in cryostat at temperatures equal to or greater than -20°C , the freezing time can be reduced up to 20 seconds, so as to have a less hard block and easier to cut. In case of fatty tissues we recommend extending the freezing time to up to 80 seconds, so as to have a harder block that helps to cut the adipose tissue.


The mould for 40 mm chuck provides a standard freezing time of 60 seconds. In case of delicate specimens which are cut in cryostat at temperatures equal to or greater than -20°C , the freezing time can be reduced up to 30 seconds, so as to have a less hard block and softer to cut. In case of fatty tissues we recommend extending the freezing time to up to 120 seconds, so as to have a harder block that helps to cut the adipose tissue.

For freezing times of differently shaped or sized moulds, see appendix B.

3.8.2. Using multiple timers

The “multiple timers” button on the Chiller page opens a screen with 6 timers:

More than one timer can be useful in some cases:

- When it is necessary to freeze more samples simultaneously.
- When different size moulds are present in the PrestoCHILL.
- When wanting to freeze different tissues in the PrestoCHILL with specific times.

The multiple timers page cannot be accessed if the timer of the main page is running.

The 6 timers are modified and used like the individual timer (see chapter 3.6).

The name of the timer can be modified: click on the name to display the alphanumeric keyboard and enter the desired name (maximum 13 digits).

3.9. Safe switching off

In order to safely switch the device off, move the main switch on the lateral side of the device to the position **O**.

- Device is switched off.
- Disconnect from power supply.
- Disconnect accessories, if necessary.

4. MAINTENANCE

4.1. Introduction

Maintenance and current repairs of devices Milestone is available anywhere in the world in one of the international service centers, as well as in your area. To arrange service or repair of your instrument, please contact your nearest office; updated contact information can be found at www.milestonemedsrl.com. For service organizations and technical support see details below.

Customer's Support

Tel.: +39 035 572735

E-mail: customersupport@milestonemedsrl.com

For Russian Federation:

BioLine LLC

197101, Saint Petersburg,

Pinsky lane 3 letter A,

Tel: (812) 320-49-49

Fax: (812) 320-49-40

Service Hotline: 8-800-333-00-49.

E-mail: service@bioline.ru

Maintenance of PrestoCHILL is carried out at two levels:

The first is at laboratory level: a simple routine maintenance can prevent the likelihood of more serious problems. Minimal routine maintenance is required by laboratory personnel that use PrestoCHILL.

The second level of maintenance is technically more detailed but can be carried out from normal laboratory personnel.

4.2 Before/After each run

Ensure that PrestoCHILL is generally kept clean after each use. Do not use abrasives or sharp instruments that may damage the coating.

Before each run cleaning of the mould is very important to avoid the block from sticking in the mould, cleaning must be performed with paper soaked in alcohol.

4.3 Daily

It's required to check the level of water condensed in the condensation trap, daily.

As you can see below, the maximum level of liquid is labeled on condensation trap. When the level is near the limit, remove the cool trap pulling it and empty any collected fluids. Then reinsert the condensation traps. Make sure that the condensed liquid does not go over the limit indicated. Information on disposal of condensed liquids, see the chapter 1.14.4.

Inscription on the condensation trap

**EMPTY THE CONTAINER IF
LEVEL IS OVER THE LINE**

A decontamination procedure must be performed daily, please refer to chapter 5 DECONTAMINATION AND CLEANING.

4.4. Every 250 defrost pump hours

For safety reasons, between the working area and the defrost pump, there is a HEPA filter that filters possibly dangerous pathogens in the vapors.

This filter must be changed after the defrost pump has worked for 250 hours.

The software advises you that the filter life is finished with the message "Filter needs to be replaced. Please do it", at the starting of the first freezing program of the day.


Avoiding HEPA filter change, can be dangerous for operator safety.

Dispose of the HEPA filter according to current regulations.

4.4.1. How to change the HEPA filter

For safety reasons, after the defrost pump has worked for 250 hours, the HEPA filter must be changed. The filter is on the rear of the unit.


1. Extract the metal connection of the filter hose from the unit cover.


2. Take out the second filter hose from its fitting.


3. Remove the filter from its housing, leveraging with the finger.

4. Insert the new filter into the housing, being sure that the filter direction (see the blue arrow on the filter) corresponds with the black arrow on the unit.

5. Now insert the metal part of the long hose into the cover hole then connect the other filter hose to its place.

After changing the filter, reset the defrost pump hour counter.

Click on the icon from the Settings page. It is necessary to insert the Administrator password (issued as a separate document inserted in the Operator Manual). Then press YES to the below question:

Have you changed the HEPA filter?

Yes

No

5. DECONTAMINATION AND CLEANING

The procedure described in this chapter must be performed daily.

For safety reasons, any operation without the required personal protection equipment is to be avoided.

PrestoCHILL decontamination is mandatory due to the biological hazard, in view of the use of fresh specimens.

All types of moulds, heat extractors, tweezers and dispensing slides are resistant to the maximum freezing temperature, cryoembedding compound, samples of biological tissues and sections, disinfectants and detergents.

5.1. Safety precautions

The personnel involved in the decontamination procedure of PrestoCHILL and its accessories (if present) must wear the appropriate individual protective devices, in particular:

- Single use, HI-RISK protective gloves (100% hi-risk nitrile or equivalent). Individual protection device that meets the standard EN 374 and the Directive 89/686/CEE.
- Antibacterial mask face for surgery room and white coat. Individual protection device, meets the standard EN 149 Grade FFP3D and the Directive 89/686/CEE.

5.2. PrestoCHILL decontamination and cleaning procedure

5.2.1. Detergents

The following instructions must be strictly complied with in order to perform proper decontamination and cleaning procedure.

Do NOT remove the unit's side panels during cleaning operations.

For decontamination and cleaning do not use procedures other than the one explained in this chapter.

Before any operation, disconnect from the power supply mains. Do not use direct or high pressure water to clean the instrument.

It is recommended to use disinfectants containing alcohols and polyphenols, e.g. Cryofect, Puli-Jet plus or any other suitable according to their parameters.

For further information, contact our specialists at the following email address:
application@milestonemedsrl.com.

Germicidal efficacy of polyphenolic based disinfectant agents:

- minutes: Immunodeficiency Virus Type 1 (HIV-1), Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), Rotavirus, Coronavirus SARS.
- 15 minutes: Staphylococcus Aureus, Candida Albicans.
- 30 minutes: Pseudomonas Aeruginosa, Mycobacterium Bovis (TB).

5.2.2. Procedure

1) Cold Plate, moulds, heat extractors, dispensing slide and tweezers

- Disinfect plate, mould and heat extractor with polyphenol-based reagent, leaving it to work for the time necessary as explained above and, if necessary, remove any organic traces using a non-abrasive tool.
- Rinse with a cloth moistened with water
- Clean with ethanol (concentration in the range 70 – 100%)
- Dry with adsorbent paper or a clean cloth.

2) Chucks

- Rinse chucks under tap water until all freezing medium (OCT) traces go away.
- Wash with an appropriate alcohol-based detergent.
- Dry with adsorbent paper or a clean cloth.

3) Final steps

- Place the positioning tool, chucks and heat extractors on the cold plate.
- Before removing gloves, wash them with soap, water and disinfectant.
- Remove the individual protection devices.
- Wash your hands.
- Make sure the contaminated material is packed and disposed in an approved incinerator in compliance with all federal, provincial and local government regulation.

On completion of the cleaning operations, wait at least 5 minutes before connecting the power supply cable.

6. WARRANTY

Warranty applies to medical instruments manufactured by Milestone S.r.l., Italy.

Milestone S.r.l. guarantees that the contractual product delivered has been subjected to a comprehensive quality control procedure based on international standards, and that the product is faultless and complies with all technical specifications.

The warranty period of device storage is not less than 6 months.

The warranty period is 13 months from the date of shipment from the date of sale.

This warranty becomes invalid in case of mishandling, damage and negligence and further in case of usage of inappropriate spare parts and/or accessories or unauthorized modification of the unit.

The instrument service life is ten years.

A. Appendix – Tips and Tricks


A.1. Fatty tissues

Fatty tissues are difficult to section using a cryostat because they are non-aqueous materials and simply do not freeze. The problem is clear when the fat begins to smear and impedes cutting of any tissue in its path. To make fat hard enough to section, it must be cooled to a very low temperature that often is too cold for non-fatty components present in the same sample. Another problem is that fat does not adhere well to the embedding medium.


The PrestoCHILL allows to freeze fat samples for an easier cut while preserving the morphology also of non-fatty tissues possibly present in the sample. Our suggestion is to freeze fatty tissue into PrestoCHILL with a temperature of -40°C . At this temperature fat freezes very fast and became hard enough to make thin slices.

Other important tips are:

- Prolong the freezing time in order to have a harder block.
- Set cryostat temperature at $-30 / -35^{\circ}\text{C}$ because this is the best to cut fat. If some non-fatty parts of the sample shatter, it is possible to warm these parts with gloved thumb. In the picture below the sample has a fatty part (arrows) and a glandular part, which is warmed by the thumb.


- Orient the tissue so that fat (arrows in the cut slide) hits the blade last or by itself.

- Utilize a new blade because sharpness helps to section this soft material.
- Freeze the sample in order to have always a handle of embedding medium for use as a hold for brush.

If, despite these measures, it is impossible to obtain a slice with the entire sample, it is possible to cut a thicker section up to $10\mu\text{m}$. Due to the transparency of adipose tissue, it is still quite interpretable even at this thickness (see picture below).


Last problem about fat is that this material tends to adhere less to the glass and then to come off during staining. In this case, it is possible to try warming the glass microscope slide before picking up the slice or cut a thinner slice.

A.2. Muscle tissue

Muscle tissues are very rich in water and this causes the formation, during freezing, of large ice crystals within and between the muscle fibres. In addition, it is important in order to morphologically evaluate the muscle tissues to have a cross section.

The PrestoCHILL is used to freeze the muscle specimen quickly and in the correct position:


1. Prepare a drop of MCC on the tip of the spatula. Take the specimen with tweezers so as to have the section of interest free and wet it with the drop of MCC.

2. Place the specimen in the mould so that the muscle fibres are vertical, this will give their cross section to the cut.

3. Holding the specimen firmly, fill the mould with MCC. Then place the chuck and heat extractor and wait for the freezing time to elapse.

4. Extract the block and cut and stain the slide.

At the end of the freezing and cutting processes, it is possible to follow the BOP 105 procedure in order to further improve the quality of the muscle tissue sections.

A.3. Small Biopsies


The difficult part in small biopsies freezing is to put them in the same level. In fact, they tend to collapse inside the medium and consequently it is necessary to do many slices to view the entire biopsy or all the biopsies. In this case, we suggest the paper embedding:

- 1) Prepare the mould by cleaning with paper soaked in alcohol and placing the appropriate chuck on it (this procedure will prevent ice inside the mould).

- 2) Apply a thin layer of embedding medium (MCC) to the tip of the dispensing slide. Place a circle of lens paper on the MCC drop and press it in order to wet the paper, then flip the paper over and place it on the drop of medium, wetting both sides. It is better to leave an edge of the paper without OCT in order to have a not sticky handle for the tweezers during sliding.

We suggest the use of blue lens paper because the blue colour gives the right contrast to better understand if the sample is correctly positioned and oriented.

- 3) Place tissue FACE DOWN and in the needed position and orientation. More than one piece of tissue can be put on one paper square according to the size, as shown in the photo above.


- 4) Touch the leading edge (preferably without MCC) of the paper to the mould floor, block the paper to the mould floor with the forceps. It will adhere to the cold mould. Then pull the dispensing slide out from under the paper.


Paper discs are resistant to the maximum freezing temperature, cryoembedding compound, samples of biological tissue and sections.

Avoid taking the paper disk with tweezers to drag it into the mould because this makes the positioning of the disc in the center of the base mould more difficult.


During this phase it is very important to use tweezers with cleaned and thin tips, otherwise they may stick to the paper, making difficult its correct positioning into the mould.

- 5) Fill the mould so that a meniscus of medium is bulging above the brim. Too much medium is not a problem. Too little medium can cause failure of the specimen to anchor to the chuck.

- 6) Place the chuck over the mould as quickly as possible and press it. Place the heat extractor over the chuck stem. Now close PrestoCHILL cover and press on countdown Start button and wait pre-set freezing time.

- 7) Take the block out. It will be flat with the paper over the specimen. The entire sample, in case of long thin biopsy, or all samples, in case of multiple short biopsies, will be at the same level and specimen will be visible through the blue lens paper due to its semi-transparency. Moreover trimming will be easier because when the paper goes away, it means that the sample is ready to cut.

B. Appendix – Freezing time

B.1. Recommended freezing time for moulds

	Freezing times (seconds)		
	Heat extractor standard Heat extractor rectangular Heat extractor with bigger hole		Heat extractor Slim
	Non-intensive use	Intensive use	Intensive use
MLD3022-5	10	30	40
MLD3022-7	30	40	
MLD3426-5	30	45	60
MLD3426-7	60	70	
MLD3426-10	80	90	
MLD4436-5	45	50	60
MLD4436-7	75	75	
MLD4436-10	180	180	
MLD3330 2419-5	30	45	50
MLD3330 2419-7	60	70	
MLD3330 2419-10	80	90	
MLD4234-5	45	50	60
MLD4234-7	75	75	
MLD4234-10	180	180	
MLD3228 2926-5	40	50	60
MLD3228 2926-7	70	75	
MLD3030 2525-5	40	50	60
MLD3030 2525-7	70	75	
MLD3030 2525-10	130	135	
MLD3939 3333-5	45	50	60
MLD3939 3333-7	75	75	
MLD3939 3333-10	180	180	
MLD4129 3324-5	40	50	60
MLD4129 3324-7	70	75	
MLD4129 3324-10	130	135	
MLD4442 4436-5	45	50	60
MLD4442 4436-7	75	75	
MLD4442 4436-10	180	180	
MLD4410-5	10	10	10
MLD4410-7	15	15	
MLD4410-10	15	15	
MLD4415-5	15	15	15
MLD4415-7	20	20	
MLD4415-10	25	25	
MLD50 2520-5	30	45	50
MLD2618-5	10	30	40
MLD2618-7	30	40	
MLD4326 1818-5	10	30	40
MLD4326 1818-7	30	40	
MLD5951-5	90		
MLD5951-7	150		

Intensive use of the PrestoCHILL means a use that involves the possibility of executing at least 5 consecutive freezing cycles.

Non-intensive use means a use in which each freezing cycle is spaced out by at least 15 minutes of inactivity.

The freezing times in the table are indicative and can vary according to the conditions of the laboratory and of the embedding medium used.

Milestone recommends using Milestone Cryoembedding Compound (MCC).

The freezing times must be confirmed by each laboratory. For further information, please contact application@milestonemedsrl.com.

B.2. Freezing time changes

The ideal freezing time depends on the size of the mould used but also on the type of frozen tissue.

In case of delicate specimens which are cut in cryostat at temperatures equal to or greater than -20°C , the freezing time can be reduced to halve the time suggested in the above table, so as to have a less hard block and better suited to the cut.

In case of fatty tissues we recommend extending the freezing time to double the freezing time suggested in the above table, so as to have a harder block that helps to cut the adipose tissue.

C. Appendix – Use the extractor for chucks

At the end of freezing the block could stick to the mould, as a result of one of the following factors:

- Frost inside the mould before its use. In this case ice forms between the block and the mould that acts as glue between the two parts.

Cleaning the mould is very important to prevent the block from sticking to the mould.

- Use of warm embedding medium (temperature above 25°C). In this case ice may form in the embedding medium layer in contact with the mould.
- Using moulds deeper than 5mm: when preparing specimens in 7 and 10mm deep moulds it is more difficult to extract the block because the contact surface between block and mould is larger.
- Using embedding medium other than Milestone Cryoembedding Compound.

Milestone recommends using embedding medium Milestone Cryoembedding Compound (MCC - Code 51420).

For optimal results, Milestone recommends using 5mm deep moulds.

C.1. How to use the extractor for chucks

- Loosen, by rotating anti-clockwise, the jaw stop, so that the jaws widen more than the diameter of the chuck handle

- Insert the jaws under the chuck handle
- Close and secure the jaws by rotating the stop clockwise

- Pull the extractor for chucks towards the other, so as to perform a vertical movement. Do not tilt the extractor for chucks.

- Loosen, by rotating anti-clockwise, the jaw stop to free the chuck with the block.

The chuck extractor is used to extract the chuck if this has a handle: in case of flat chucks, you can use the connection to the cryostat, in case of chucks with pin, there is no handle unless the o-ring is removed in order to use the o-ring slot as connection for the chuck extractor.


D. Appendix – How to change a mould

D.1. How to change a mould

On the basis of the blocks size needed in the laboratory, it is possible to change the mould on the cold plate.

- Using the key for the changing of the mould, unscrew the wrong size mould.

- Screw the desired mould into place and tighten with the key.

E. Appendix – Troubleshooting

E.1. Block detachment from the chuck

At the end of freezing, it may happen that the MCC block with the specimen, when detaching from the chuck, remains inside the mould.

In this case, add MCC on the block and re-fit the chuck and heat extractor. Wait 30 seconds then remove the block by pulling the chuck vertically.

Tips to avoid this problem:

1. Clean the mould with paper soaked in ethanol before its use.
2. Make sure to place a suitable amount of embedding medium to reach the clamping engravings of the chuck: better to put excess embedding medium rather than not enough.
3. Extend the freezing time.
4. Make sure that the embedding medium residues are removed from the used chuck and that the engravings for clamping to the MCC are not worn. If the chuck does not have optimum clamping, it is best to keep it at room temperature (outside the working area of the PrestoCHILL) before its use. This allows the MCC to have more time, before freezing, to penetrate the clamping engravings.
5. Prefer the standard heat extractor that induces faster freezing in the clamping area between embedding medium and chuck.
6. Prefer embedding medium Milestone Cryoembedding Compound (MCC) which, being more fluid, penetrates better the clamping engravings of the chuck.
7. Do not use the embedding medium at a temperature above 25°C.
8. Avoid extracting the chuck by tilting it, pull with a vertical movement. Use the extractor for chucks provided with the unit. For its use, refer to appendix C.
9. Prefer the use of 5mm deep moulds to deeper ones.


E.2. Power failure

In the event of a power failure during a chiller program, when the power returns, PrestoCHILL turns on again and goes back to the home page.

If in manual mode (see chapter 3.4) the chiller program must be restarted by the user.

If in automatic mode (see chapter 3.5) the chiller program restarts after 10 seconds.

In the event of a power failure during a defrosting process in manual mode, when the power returns, the defrosting process does not restart and the following screen appears:


It is not possible to start with a chiller program: the software informs that it is necessary to run a defrosting program to clean the working area of the PrestoCHILL.

Press OK and start a defrosting program. Alternatively, dry and thoroughly clean the working area of the PrestoCHILL and all accessories before starting a chiller program.

Before starting, wipe the work plate of the PrestoCHILL with paper or similar absorbent material to dry it.

In the event of a power failure during a defrosting process in automatic mode, when the power returns, the defrosting process does not restart and any chiller programs are blocked. The following screen will appear:

Press OK and start a manual defrosting program. Alternatively, dry and thoroughly clean the working area of the PrestoCHILL and all accessories and wait for the PrestoCHILL to run the preset automatic cycle (paragraph 3.5.1).

Before starting, wipe the work plate of the PrestoCHILL with paper or similar absorbent material to dry it.

F. Appendix – How to manage an alarm**F.1. How to manage an alarm**

CODE	DESCRIPTION	WHAT TO DO
1	Communication fault	1. Press 'Silent' to hide the alarm. 2. Reboot the system. (wait 2 minutes before the rebooting) 3. Call your customer service if the issue can't be fixed.
2	Chiller temperature not reached	
3	Defrost temperature not reached	
6	Chiller system error	
7	Chiller temperature out of range	

G. Appendix – Information about Milestone Cryoembedding Compound (MCC)


No	Feature description	Values
1	Content	polyvinyl alcohol 16%, polyethyleneglycol 6,2%, benzalkonium chloride 0,26%, distilled water
2	Appearance	Liquid
3	Color	colorless
4	Odor	odorless
5	Density	1,046 kg/l
6	Freezing temperature	- 5 °C
7	Validity	24 months

H. Appendix – Marking of PrestoCHILL and accessories


H.1. Shield of instrument

See chapter 1.12.

H.2. Label of PrestoCHILL package

H.3. Labels and marking of accessories

* - marking of Operator Manual on USB-pen.

MILESTONE **MILESTONE s.r.l.**
 Via Fatebenefratelli, 1-5
 24010 Sorisole (BG) - Italy
 Tel. +39.035.4128284
 Fax. +39.035.575498
 e-mail: medical@milestonesrl.com
 www.milestonemedsrl.com

REF 100545
SMALL TIP SPATULA

LOT

IVD CE

MILESTONE **MILESTONE s.r.l.**
 Via Fatebenefratelli, 1-5
 24010 Sorisole (BG) - Italy
 Tel. +39.035.4128284
 Fax. +39.035.575498
 e-mail: medical@milestonesrl.com
 www.milestonemedsrl.com

REF 100551
LARGE TIP SPATULA

LOT

IVD CE

P/N: 100575

Descr.: PRESTOCHILL COMPLETE MOULD KEY

Lot N°:

Serial N°: **dated** **Q.ty: 1**

Origin: Italy

P/N: 100593

Descr.: PRESTOCHILL SMALL FORCEP

Lot N°:

Serial N°: **Q.ty: 1**

Origin: Italy

P/N: 100741

Descr.: TITANIUM FORCEPS L145

Lot N°:

Serial N°: **Q.ty: 1**

Origin: Italy

P/N: 100533**Descr.: HEAT EXTRACTOR STANDART****Lot N°:****Serial N°:** **dated:** **Q.ty: 1****Origin: Italy****P/N: 100677****Descr.: HEAT EXTRACTOR SLIM****Lot N°:****Serial N°:** **dated:** **Q.ty: 1****Origin: Italy****P/N: 100592****Descr.: HEAT EXTRACTOR RECTANGULAR****Lot N°:****Serial N°:** **dated:** **Q.ty: 1****Origin: Italy****P/N: 100736****Descr.: HEAT EXTRACTOR WITH BIGGER HOLE****Lot N°:****Serial N°:** **dated:** **Q.ty: 1****Origin: Italy****P/N: 100737L****Descr.: CHUCK EXTRACTOR****Lot N°:****Serial N°:** **dated:** **Q.ty: 1****Origin: Italy**

REF 100623
LENS PAPER DISC 16

LOT

MILESTONE s.r.l.
Via Fatebenefratelli, 1-5
24010 Sorisole (BG) - Italy
Tel. +39.035.4128284
Fax. +39.035.575488
e-mail: medical@milestoneart.com
www.milestonemedart.com


REF 100622
LENS PAPER DISC 22

LOT

MILESTONE s.r.l.
Via Fatebenefratelli, 1-5
24010 Sorisole (BG) - Italy
Tel. +39.035.4128284
Fax. +39.035.575488
e-mail: medical@milestoneart.com
www.milestonemedart.com

REF 100615
LENS PAPER DISC 30

LOT

MILESTONE s.r.l.
Via Fatebenefratelli, 1-5
24010 Sorisole (BG) - Italy
Tel. +39.035.4128284
Fax. +39.035.575488
e-mail: medical@milestoneart.com
www.milestonemedart.com

P/N: **100624**


Descr.: KIT HEPA FILTER

Lot N°:

Serial N°:

dated:

Q.ty: **1**

Label and marking for Moulds


XXXXXX* - mould catalog number, e.g. 100590

MLD0000 0000-0* - abbreviation name of mould, e.g. MLD3228 2926-5 or MLD2618-7.


For symbols explanation see chapter 1.1

I. Appendix – The list of national standards applied by the manufacturer

- 1) 98/79/EC In vitro diagnostic medical devices Directive
- 2) 2014/35/EC Low Voltage Directive
- 3) 2014/30/EC Electromagnetic compatibility EMC Directive
- 4) 2011/65/EC Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment
- 5) CEI EN 61010-1:2001 Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1. General requirements.
- 6) CEI EN 61010-2-101:2002 Safety requirements for electrical equipment for measurement, control and laboratory use. Particular requirements for in vitro diagnostic (IVD) medical equipment.
- 7) CEI EN 61326-1:2013 Electrical equipment for measurement, control and laboratory use. EMC requirements. General requirements.
- 8) CEI EN 61326-2-6:2013 Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6: particular requirements - In vitro diagnostic (IVD) medical equipment.
- 9) EN IEC 62304:2006 Medical devices. Software. Life cycle process.
- 10) IEC/ CEI EN 61010-2-010 Safety of electric control measuring and laboratory equipment. Part 2-010. Particular requirements to laboratory equipment for heating materials.
- 11) IEC/ CEI EN 55011 Industrial, scientific and medical equipment. RFI parameters. Limit values and measurement methods.

Approved 
(VISINONI Francesco - PRESIDENT)

MILESTONE S.r.l.
Via Fatebenefratelli, 1/5
24010 SORISOLE (BG)
P. IVA 01879330163

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At the request of the person concerned, we declare, without responsibility on the content, that copy of this document has been lodged at Bergamo Chamber of Commerce.

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