APOSTILLE

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 Eurocor
JOKER®035
BALLOON DILATATION CATHETER

Instructions for use

STERILE. Sterilized with ethylene oxide gas. Pyrogen-free. For single use only. Do not autoclave.

The **JOKER®035 PTA Balloon Dilatation Catheter** is an over the wire (OTW) catheter developed for the percutaneous transluminal angioplasty (PTA) of peripheral vessels.

The **JOKER®035 PTA Balloon Dilatation Catheter** has a double-lumen coaxial shaft with a balloon close to the distal tip and a Y-connector on the proximal end. The straight port of the Y-connector is intended for the guide wire insertion. The catheter is compatible with a 0.035" inch (0.89mm) guide wire - the compatibility is indicated on the product label.

The side port of the connector is intended for the inflation and deflation of the balloon. The working length of the instrument is 120cm or 150cm.

Due to the controlled compliance of the balloon material, specific balloon diameters can be achieved at specific pressure (see enclosed compliance table). The **JOKER 035 PTA Balloon Dilatation Catheter** is available in different balloon sizes. The nominal balloon diameter and the nominal balloon length are printed on the proximal shaft.

The **JOKER®035 PTA Balloon Dilatation Catheter** is recommended for the percutaneous transluminal angioplasty of peripheral vessels, e.g. iliac, femoralis, iliac-femoral, poplitea, infrainguinal, renal and carotid as well as for treating obstructive lesions of native or artificial arteriovenous dialysis fistulas. This balloon is also intended for stent post-dilatation in the peripheral vascular system.

No contraindications are known for percutaneous transluminal angioplasty (PTA). The use of the **JOKER 035** in the coronary arteries is contraindicated.

- This catheter is intended for single use only. DO NOT resterilize and / or reuse

resterilization and / or reuse may potentially affect the product's performance and increase the risk of improper resterilization and cross-contamination. This type of treatment may result in a non-sterile product and might result in infection of the patient and the risk of a rupture or bursting of the balloon.

- Do not bend or squeeze the balloon area of the **JOKER 035 PTA Balloon Dilatation Catheter** before the procedure.
- To reduce the risk of vascular damage, the diameter of the inflated balloon should correspond to the diameter of the artery distal and proximal to the stenosis.
- As soon as the catheter has been introduced into the vascular system, a control of the manipulations using very advanced X-ray equipment is required. Advance and pull back the catheter only when the balloon is completely deflated under vacuum.

If some resistance is felt during manipulation, determine its cause before continuing the examination.

- Balloon pressure must not exceed the rated burst pressure (RBP). The RBP is

Eurocor GmbH
In den Dauen 6a
53117 Bonn - Germany
eurocor.de

- At least 99.9% of the balloons will not burst (with the RBP). We strongly recommend the use of a pressure monitoring device in order to avoid excessively high pressure.
- Use the recommended balloon inflation agent. Do not use air or another gaseous agent for balloon inflation.
- Use the catheter before the end of the expiration month ("Use by") printed on the packing.
- Do not use the product if the use by date on the packing has been exceeded.
- Store the packed product in a dry place at a temperature between +5 °C and +40° C.
- Protect the packed product against direct sunlight.

Complications Possible complications in connection with the use of the **JOKER®035 PTA Balloon-Dilatation Catheters** are similar to those connected with standard PTA. Possible complications include, without being limited to:

- Death
- Haemorrhage or haematoma
- Restenosis of the dilated vessel
- Dissection, perforation, rupture, spasms or injuries of the treated vessel
- Allergic reactions to drugs and / or contrast medium
- Hypotension
- Infection
- Arteriovenous fistula
- Pain
- Arrhythmia
- Pseudo-aneurysms
- Thrombosis and / or embolisation
- Sepsis
- Endocarditis

Preparation of the balloon catheter

1. The balloon filling fluid is usually a mixture of 50:50 contrast medium and sterile saline solution or a commonly used contrast medium with 200 mg/ml iodine.
2. Selection of the dilatation catheter
The nominal balloon size should be selected to correspond to or be smaller than the internal diameter of the vessel proximal and distal to the stenosis when the balloon is inflated. If the desired balloon diameter cannot pass the stenosis, a balloon catheter of a smaller diameter can be used to predilate the stenosis.
3. Insert approx. 3 ml contrast medium solution into a 10 ml or 20 ml syringe. Connect the syringe to the Luer connector of the balloon catheter. Hold the two devices pointed downwards. Use the syringe to draw a vacuum (aspirate). Repeat, until no more air bubbles rise in the syringe. When deflated, the catheter contains small air bubbles that have to be removed before inserting the dilatation catheter. To remove the air bubbles, connect the 3-way cock to the lateral port of the connector and flush the stopcock. Connect the inflation device with a mixture of contrast medium / saline solution to the 3-way cock. Never dilate the balloon using air or a gaseous medium. Hold the inflation device and the catheter pointing downwards and aspirate for approx. 15 seconds. Close the stopcock and evacuate all the air from the inflation device. Repeat the previous step until no more air bubbles appear when aspirating.
4. Connection of the balloon catheter to the manometer / syringe system. It has been filled with a sufficient quantity of filling fluid before (see step 1).
5. Gently pull the catheter out of the protective dispenser.
6. Remove the balloon cover and the guard wire from the lumen of the guiding wire. The catheter is now ready for the procedure.

Instructions for use / Application technique

Note: Do not use Ethiodol* or lipiodol or a contrast medium having the same ingredients.

Note: Do not expose the catheter to organic solutions (e.g. alcohol).

Note: Do not use when the sterile individual packing has been opened or damaged.

* Ethiodol is a registered trademark of Guerbet S.A.

Note: Do not autoclave. Do not resterilize.

Note: Store in a cool, dark and dry place

1. Percutaneous insertion technique or arteriotomy suitable for the use of insertion instruments and guiding wires.
 2. Insert a guiding catheter or a sheath with haemostatic adapter using standard technique.
 3. If required, the lumen of the guiding wire can be rinsed using physiological saline solution.
 4. Slowly advance the balloon catheter with completely deflated balloon through the haemostatic valve. It is important to close the haemostatic valve firmly enough to avoid the leaking of blood, but not to impede the mobility of the balloon catheter. If a resistance is felt, do not further advance the balloon catheter.
 5. Gentle insertion of the guiding wire and placing the tip distal of the stenosis
- Note:** In order to avoid kinking, advance the balloon catheter slowly in small steps until the proximal end of the guiding wire leaves the catheter.
6. The **JOKER 035 PTA Balloon Dilatation Catheter** is placed in the area of stenosis under contrast medium and X-ray control.
 7. After dilatation, completely empty the balloon catheter and pull it back into the sheath and / or the guiding wire under X-ray control and check the result of dilatation by means of angiography.
 8. As the case may be, it is also possible to replace the catheter used through the guiding wire left in the vessel.
 9. Connect the lateral connection of the haemostatic adapter for the guiding catheter with the proximal pressure indicator / infusion line or the distributor in order to make registering of the proximal pressure or the infusion through the guiding wire possible
 10. Place the balloon in the lesion to be dilated under X-ray control and using the radiopaque markers and inflate to the prescribed pressure (see compliance table). Maintain low pressure at the balloon between balloon inflations.
 11. Completely deflate the balloon catheter. Remove the deflated balloon catheter and the guiding wire through the haemostatic valve from the sheath / the guiding wire. Tighten the knurled knob at the haemostatic valve.
- In case of insufficient primary result, continue with further dilatations, as the case may be, postdilate using a balloon of a larger diameter.

Compliance Chart

Page.
JOKER 035.

The tolerance for the balloon diameter is $\pm 5\%$ for the entire operating

Inflation Pressure		Balloon Diameter (mm)				
atm/bar	kPa	4.00	5.00	6.00	7.00	8.00
4	400	3.78	4.80	5.68	6.67	7.45
6*	600*	4.00*	5.00*	6.00*	7.00*	8.00*
8	800	4.15	5.15	6.16	7.19	8.31
10	1000	4.28	5.26	6.30	7.33	8.53
11	1100	4.33	5.31	6.37	7.41	8.64
12***	1200	4.40	5.37	6.44	7.48***	8.77***
13	1300	4.46	5.43	6.52	7.55	8.89
14** ***	1400	4.53	5.49***	6.59***	7.63**	9.00**
15	1500	4.60	5.58	6.67	7.70	9.12
16** ***	1600	4.67** ***	5.64**	6.79**	7.78	9.25
17	1700	4.75	5.72	6.84		
18	1800	4.83	5.80			

* Nominal pressure

** Rated burst pressure (RBP) for Balloon length 20 -60 mm

*** Rated burst pressure (RBP) for Balloon length 80 -150 mm

- This product should only be used by physicians who are experienced in angiography, percutaneous transluminal angioplasty (PTA) and, as the case may be, a vascular surgical emergency team should be on call.
- When taking the **JOKER®035 PTA Balloon-Dilatation Catheter** out of the packaging and passing the haemostatic valve, you have to take care not to bend, damage or make the catheter unsterile.
- It is possible to rinse the lumen of the guiding wire with a sterile, isotonic saline solution.
- The balloon must not be inflated prematurely. The balloon pressure recommended by the manufacturer must not be exceeded. It is recommended to use a pressure manometer for monitoring the inflation pressure.
- If resistance is felt when advancing the catheter in the vessel, utmost precaution is required since this obstacle might indicate a damage of the balloon catheter. If the resistance occurs on the way through the guiding wire, the complete insertion system should be pulled back, as the case may be.
- **Attention:** Non-ionic contrast medium differs from ionic contrast medium with respect to viscosity and precipitation and may result in long in- and deflation times.
- The smallest compatible sheath French size is printed on the product label. Do not attempt to push the PTA catheter through a sheath smaller than that indicated on the label.
- Not suitable for precise blood pressure monitoring.

Handling and precautionary measures

The inflated diameter of the balloon should correspond to the reference diameter of the peripheral artery, never use a balloon of a larger diameter.

Do not shake the balloon during preparation in order to remove bubbles.

Do not insert the guide wire into the balloon catheter when there is resistance without having determined its cause and taken remedial action. This PTA catheter is pressure tested in a way that, with 95% probability, 99.9% of the balloons or the balloon connections do not show any defects or leaks under the rated burst pressure (RBP).

It is not recommended to exceed the rated burst pressure, since this pressure may result in the rupture of the balloon or may damage the joints of the catheter and/or balloon. If this happens, the balloon must be emptied and removed. Note: Do not use when the sterile packing is opened or damaged. Do not resterilize. Inspect the **JOKER®035 Balloon PTA Balloon-Dilatation Catheter** before use in order to ensure that it has not been damaged during transit. The catheter should only be used by physicians who have experience in the implementation of angiography and are trained in percutaneous transluminal angioplasty.

Attention

Eurocor individual products and kits are compatible with one another. Before use, the user must ensure the intercompatibility of the products. This applies in particular, if Eurocor products are used in connection with products of other manufacturers.

References

The user should be familiar with the most recent publications concerning the current medical practice of balloon dilatation.

Medication scheme

Upon admittance to hospital, the values for prothrombin time (PT) and the partial thromboplastin time (PTT) as well as the platelet count must be determined for each patient for documentation of the baseline values. The following medication scheme serves as a possible guideline only and cannot be considered a strict medical requirement. Prior to balloon inflation:

- ASA 100 up to 325 mg 1 x daily
- Clopidogrel at a loading dose of 300 mg, or alternatively at 75 mg/day (start 3 days prior to implantation)

During balloon inflation:

- 2500 - 5000 I. E. Heparin i.u. After

balloon inflation:

- ASA 100 up to 325 mg 1 x daily
- Clopidogrel 75 mg/day for at least 4 weeks

Note: The doses, notes and values stated are to be examined very carefully by the user prior to the dilatation of the **JOKER®035 Balloon PTA Balloon-Dilatation Catheter**. Given the further development of the corresponding procedures during the last years, other approaches, differing from those described above, are absolutely possible.

The packed products have to be protected from direct light. When stored properly, the product may be used until the stated date (see packing). The packed products should not be exposed to temperatures under +5°C or above +40°C. Use by date: 24 months after sterilization.

The product has tight individual packing, which ensures the preservation of its operational and medical characteristics during the shelf life. Each packing contains 1 item. Packaging material: plastic package put in a box.

After use After their use, medical products and their accessories may constitute potential biological hazards. Therefore, the medical products used and their accessories must be handled and disposed of in compliance with accepted medical procedures and in accordance with the respective statutory regulations and local provisions.

Environmental control During usage, transportation and storage the product must not have an adverse impact on human and environment.

Transportation Catheters can be transported by all means of transport in compliance with haulage rules under the temperature control (from 5°C up to 45°C) which applies to the particular mean of transport. The product should be shipped in roofed transport.

Warranty / Liability The product and every component of the system have been designed, manufactured, checked and packed as carefully as possible. The warnings in the instructions for use issued by Eurocor shall expressly apply as integral part of this clause. Eurocor's warranty for the product shall be valid until its expiry date. The warranty will only be valid if the product has been used in accordance with the instructions for use. Eurocor does not assume any warranty for the product in terms of marketability or suitability for a certain purpose. Eurocor shall not be liable for direct, indirect, collateral or consequential damage caused by the product. The indemnification for damage suffered by the buyer shall in no case exceed the invoice price of the products except for the case of fraud or gross misconduct. The warranty stated in this clause shall contain and replace the legal warranties for defects and compliance and shall exclude any other possible liability of Eurocor for the delivered product, irrespective of their origin. It is not intended that these restrictions with respect to liability and warranty conflict with the compulsory provisions of applicable laws. If a clause of the warranty disclaimer is declared as invalid or conflicting with applicable law by a competent court, this shall not affect its remaining provisions, which will remain fully effective. The invalid provision will be replaced by a valid provision that reflects Eurocor's legitimate interests regarding the limitation on liability or warranty best. No person shall be authorized to commit Eurocor to a warranty or liability for the product.

Subject to technical changes.

Manufacturer

Eurocor GmbH
In den Dauen 6a, 53117 Bonn, Germany
Phone +49 (0) 228 20150-0
Fax +49 (0) 228 20150-5
info@eurocor.de
www.eurocor.de

**Representative of the
manufacturer in Russia**

LLC "FORLIFE" 194295, St. Petersburg, street of Ivan Fomin, 5, Building 2, ofis 143
E-mail: forlife2012@mail.ru
Tel (fax): (812) 327-13-90, 328-89-69

REF

Catalogue No.

LOT

Batch Code

STERILE EO

Sterilized using ethylene oxide

Do not re sterilize

Do not reuse

Consult instructions for use

Do not use if package is damaged

Keep away from sunlight

Keep away from rain

Date of manufacture

Use by date

Attention

Eurocor GmbH

In den Dauen 6a, 53117 Bonn, Germany

Phone +49 (0) 228 20150-0

Fax +49 (0) 228 20150-5

www.eurocor.de

CE

0434

Rev. No. 0712-L6 J035
Version 03
Last update July 2012

JOKER®014

BALLOON DILATATION CATHETER

Instructions for use

STERILE. Sterilized with ethylene oxide gas. Pyrogen-free. For single use only. Do not autoclave.

Designation

The **JOKER®014 PTA Balloon Dilatation Catheter** is an over the wire (OTW) catheter developed for the percutaneous transluminal angioplasty (PTA) of peripheral vessels.

Description

The **JOKER®014 PTA Balloon Dilatation Catheter** has a double-lumen

shaft with a balloon close to the distal tip and a Y-connector on the proximal end. The straight port of the Y-connector is intended for the guide wire insertion. The catheter is compatible with a 0.014" inch (0.36mm) guide wire - the compatibility is indicated on the product label.

The side port of the connector is intended for the inflation and deflation of the balloon. The working length of the instrument is 120cm or 150cm.

Due to the controlled compliance of the balloon material, specific balloon diameters can be achieved at specific pressure (see enclosed compliance-table). The **JOKER 014 PTA Balloon Dilatation Catheter** is available in different balloon sizes. The nominal balloon diameter and the nominal balloon length are printed on the proximal shaft.

Indications

The **JOKER 014 PTA Balloon Dilatation Catheter** is recommended for the percutaneous transluminal angioplasty of peripheral vessels, e.g. iliac, femoralis, iliac-femoral, poplitea, infringuinal, renal and carotid as well as for treating obstructive lesions of native or artificial arteriovenous dialysis fistulas. This balloon is also intended for stent post-dilatation in the peripheral vascular system.

Contraindications

No contraindications are known for percutaneous transluminal angioplasty (PTA). The use of the **JOKER 014** in the coronary arteries is contraindicated.

Warning

- This catheter is intended for single use only. DO NOT resterilize and / or reuse since resterilization and / or reuse may potentially affect the product's performance and increase the risk of improper resterilization and cross-contamination. This type of treatment may result in a non-sterile product and might result in infection of the patient and the risk of a rupture or bursting of the balloon.
 - Do not bend or squeeze the balloon area of the **JOKER®014 PTA Balloon Dilatation Catheter** before the procedure.
 - To reduce the risk of vascular damage, the diameter of the inflated balloon should correspond to the diameter of the artery distal and proximal to the stenosis.
 - As soon as the catheter has been introduced into the vascular system, a control of the manipulations using very advanced X-ray equipment is required. Advance and pull back the catheter only when the balloon is completely deflated under vacuum.
- If some resistance is felt during manipulation, determine its cause before continuing the examination.
- Balloon pressure must not exceed the rated burst pressure (RBP). The RBP is based on in-vitro tests. At least 99.9% of the balloons will not burst (with

 Eurocor GmbH
In den Dauen 6a
53117 Bonn - Germany
eurocor.de

...is not suitable for use with hypertensive high pressure.
...Do not use air or another
...through catheter.
...at the end of the expiration month ("use by") printed on the
...the use by date on the packing has been exceeded.
...in a dry place at a temperature between +5 °C and +40°

- Protect the packed product against direct sunlight.

Complications Possible complications in connection with the use of the **JOKER 014 PTA Balloon-Dilatation Catheters** are similar to those connected with standard PTA. Possible complications include, without being limited to:

- Death
- Haemorrhage or haematoma
- Restenosis of the dilated vessel
- Dissection, perforation, rupture, spasms or injuries of the treated vessel
- Allergic reactions to drugs and / or contrast medium
- Hypotension
- Infection
- Arteriovenous fistula
- Pain
- Arrhythmia
- Pseudo-aneurysms
- Thrombosis and / or embolization
- Sepsis
- Endocarditis

**Inflation of the
balloon catheter**

1. The balloon filling fluid is usually a mixture of 50:50 contrast medium and sterile saline solution or a commonly used contrast medium with 200 mg/ml iodine.
2. Selection of the dilatation catheter
The nominal balloon size should be selected to correspond to or be smaller than the internal diameter of the vessel proximal and distal to the stenosis when the balloon is inflated. If the desired balloon diameter cannot pass the stenosis, a balloon catheter of a smaller diameter can be used to predilate the stenosis.
3. Insert approx. 3 ml contrast medium solution into a 10 ml or 20 ml syringe. Connect the syringe to the Luer connector of the balloon catheter. Hold the two devices pointed downwards. Use the syringe to draw a vacuum (aspirate). Repeat, until no more air bubbles rise in the syringe. When deflated, the catheter contains small air bubbles that have to be removed before inserting the dilatation catheter. To remove the air bubbles, connect the 3-way cock to the lateral port of the connector and flush the stop cock. Connect the inflation device with a mixture of contrast medium / saline solution to the 3-way cock. Never dilate the balloon using air or a gaseous medium. Hold the inflation device and the catheter pointing downwards and aspirate for approx. 15 seconds. Close the stop cock and evacuate all the air from the inflation device. Repeat the previous step until no more air bubbles appear when aspirating.
4. Connection of the balloon catheter to the manometer / syringe system. It has been filled with a sufficient quantity of filling fluid before (see step 1).
5. Gently pull the catheter out of the protective dispenser.
6. Remove the balloon cover and the guard wire from the lumen of the guiding wire. The catheter is now ready for the procedure.

**Instructions for use /
Application technique**

Instructions for use. Note: Do not use Ethiodol* or Iopiodol or a contrast medium having the same ingredients.

Note: Do not expose the catheter to organic solutions (e.g. alcohol).

Note: Do not use when the sterile individual packing has been opened or damaged.

Note: Do not autoclave. Do not resterilize.

Note: Store in a cool, dark and dry place

1. Percutaneous insertion technique or arteriotomy suitable for the use of insertion instruments and guiding wires.

* Ethiodol is a registered trademark of Guerbet S.A.

2. Insert a guiding catheter or a catheter with haemostatic adapter using standard technique.
 3. If required, the lumen of the guiding wire can be rinsed using physiological saline solution.
 4. Slowly advance the balloon catheter with completely deflated balloon through the haemostatic cock. It is important to close the haemostatic valve firmly enough to avoid the leaking of blood, but not to impede the mobility of the balloon catheter. If a resistance is felt, do not further advance the balloon catheter.
 5. Gentle insertion of the guiding wire and placing the tip distal of the stenosis
Note: In order to avoid kinking, advance the balloon catheter slowly in small steps until the proximal end of the guiding wire leaves the catheter.
 6. The JOKER® 014 PTA Balloon-Dilatation Catheter is placed in the area of stenosis under contrast medium and X-ray control.
 7. After dilatation, completely empty the balloon catheter and pull it back into the sheath and / or the guiding wire under X-ray control and check the result of dilatation by means of angiography.
 8. As the case may be, it is also possible to replace the catheter used through the guiding wire left in the vessel.
 9. Connect the lateral connection of the haemostatic adapter for the guiding catheter with the proximal pressure indicator / infusion line or the distributor in order to make registering of the proximal pressure or the infusion through the guiding wire possible.
 10. Place the balloon in the lesion to be dilated under X-ray control and using the radiopaque markers and inflate to the prescribed pressure (see compliance table). Maintain low pressure at the balloon between balloon inflations.
 11. Completely deflate the balloon catheter. Remove the deflated balloon catheter and the guiding wire through the haemostatic valve from the sheath / the guiding wire. Tighten the knurled knob at the haemostatic valve.
- In case of insufficient primary result, continue with further dilatations, as the case may be, postdilate using a balloon of a larger diameter.

The tolerance for the balloon diameter is $\pm 5\%$ for the entire operating range.

Inflation Pressure		Balloon Diameter (mm)				
atm / bar	kPa	2.0	2.5	3.0	3.5	4.0
4	400	1.88	2.37	2.79	3.27	3.82
6	600	2.00	2.50	2.92	3.40	3.95
8	800	2.09	2.60	3.17	3.72	4.12
10	1000	2.16	2.70	3.29	3.87	4.19
12	1200	2.22	2.77	3.39	4.00	4.27
14	1400	2.27	2.83			
16	1600			3.54	4.18	4.42
18	1800	2.42	2.94	3.61	4.26	4.50
20	2000	2.53	3.00			

I-I nominal pressure I-I rated burst pressure

Warning / Precautions for use

- This product should only be used by physicians who are experienced in angiography, percutaneous transluminal angioplasty (PTA) and, as the case may be, a vascular surgical emergency team should be on call.
- When taking the JOKER®014 PTA Balloon-Dilatation Catheter out of the packaging and passing the haemostatic valve, you have to take care not to bend, damage or make the catheter unsterile.
- It is possible to rinse the lumen of the guiding wire with a sterile, isotonic saline solution.
- The balloon must not be inflated prematurely. The balloon pressure recommended by the manufacturer must not be exceeded. It is recommended to use a pressure manometer for monitoring the inflation pressure.
- If resistance is felt when advancing the catheter in the vessel, utmost precaution is required since this obstacle might indicate a damage of the balloon catheter.

If the resistance occurs on the way through the guiding wire, the complete insertion system should be pulled back, as the case may be.

- **Attention:** Non-ionic contrast medium differs from ionic contrast medium with respect to viscosity and precipitation and may result in long in- and deflation times.
- The smallest compatible sheath French size is printed on the product label.
Do not attempt to push the PTA catheter through a sheath smaller than that indicated on the label.
- Not suitable for precise blood pressure monitoring.

The inflated diameter of the balloon should correspond to the reference diameter of the peripheral artery, never use a balloon of a larger diameter.
Do not shake the balloon during preparation in order to remove bubbles.
Do not insert the guide wire into the balloon catheter when there is resistance without having determined its cause and taken remedial action. This PTA catheter is pressure tested in a way that, with 95% probability, 99.9% of the balloons or the balloon connections do not show any defects or leaks under the rated burst pressure (RBP).

It is not recommended to exceed the rated burst pressure, since this pressure may result in the rupture of the balloon or may damage the joints of the catheter and/or balloon. If this happens, the balloon must be emptied and removed. Note: Do not use when the sterile packing is opened or damaged. Do not resterilize. Inspect the **JOKER 014 Balloon PTA Balloon-Dilatation Catheter** before use in order to ensure that it has not been damaged during transit. The catheter should only be used by physicians who have experience in the implementation of angiography and are trained in percutaneous transluminal angioplasty.

Eurocor individual products and kits are compatible with one another. Before use, the user must ensure the intercompatibility of the products. This applies in particular, if Eurocor products are used in connection with products of other manufacturers.

The user should be familiar with the most recent publications concerning the current medical practice of balloon dilatation.

Upon admittance to hospital, the values for prothrombin time (PT) and the partial thromboplastin time (PTT) as well as the platelet count must be determined for each patient for documentation of the baseline values. The following medication scheme serves as a possible guideline only and cannot be considered a strict medical requirement.

Prior to balloon inflation:

- ASA 100 up to 325 mg 1 x daily
- Clopidogrel at a loading dose of 300 mg, or alternatively at 75 mg/day (start 3 days prior to implantation)

During balloon inflation:

- 2500 - 5000 I. E. Heparin i. u.

After balloon inflation:

- ASA 100 up to 325 mg 1 x daily
- Clopidogrel 75 mg/day for at least 4 weeks

Note: The doses, notes and values stated are to be examined very carefully by the user prior to the dilatation of the **JOKER®014 Balloon PTA Balloon-Dilatation Catheter**. Given the further development of the corresponding procedures during the last years, other approaches, differing from those described above, are absolutely possible.

the product may be used until the stated date (see packing). The packed products should not be exposed to temperatures under +5°C or above +40°C. Use by date: 24 months after sterilization.

The product has tight individual packing, which ensures the preservation of its operational and medical characteristics during the shelf life. Each packing contains 1 item. Packaging material: plastic package put in a box.

after use

After their use, medical products and their accessories may constitute potential biological hazards. Therefore, the medical products used and their accessories must be handled and disposed of in compliance with accepted medical procedures and in accordance with the respective statutory regulations and local provisions.

environmental control

During usage, transportation and storage the product must not have an adverse impact on human and environment.

transportation

Catheters can be transported by all means of transport in compliance with haulage rules under the temperature control (from 5°C up to 45°C) which applies to the particular mean of transport. The product should be shipped in roofed transport.

warranty / Liability

The product and every component of the system have been designed, manufactured, checked and packed as carefully as possible. The warnings in the instructions for use issued by Eurocor shall expressly apply as integral part of this clause. Eurocor's warranty for the product shall be valid until its expiry date. The warranty will only be valid if the product has been used in accordance with the instructions for use. Eurocor does not assume any warranty for the product in terms of marketability or suitability for a certain purpose. Eurocor shall not be liable for direct, indirect, collateral or consequential damage caused by the product. The indemnification for damage suffered by the buyer shall in no case exceed the invoice price of the products except for the case of fraud or gross misconduct. The warranty stated in this clause shall contain and replace the legal warranties for defects and compliance and shall exclude any other possible liability of Eurocor for the delivered product, irrespective of their origin. It is not intended that these restrictions with respect to liability and warranty conflict with the compulsory provisions of applicable laws. If a clause of the warranty disclaimer is declared as invalid or conflicting with applicable law by a competent court, this shall not affect its remaining provisions, which will remain fully effective. The invalid provision will be replaced by a valid provision that reflects Eurocor's legitimate interests regarding the limitation on liability or warranty best. No person shall be authorized to commit Eurocor to a warranty or liability for the product. Subject to technical changes.

Manufacturer

Eurocor GmbH

In den Dauen 6a, 53117 Bonn, Germany
Phone +49 (0) 228 20150-0
Fax +49 (0) 228 20150-5
info@eurocor.de
www.eurocor.de

**Representative of the
manufacturer in Russia**

LLC "FORLIFE" 194295, St. Petersburg, street of Ivan Fomin, 5, Building 2, ofis 143
E-mail: forlife2012@mail.ru
Tel (fax): (812) 327-13-90, 328-89-69

	Catalogue No.
	Batch Code
	Sterilized using ethylene oxide
	Do not resterilize
	Do not reuse
	Consult instructions for use
	Do not use if package is damaged
	Keep away from sunlight
	Keep away from rain
	Date of manufacture
	Use by date
	Attention

Eurocor GmbH
In den Dauen 6a, 53117 Bonn, Germany
Phone +49 (0) 228 20150-0
Fax +49 (0) 228 20150-5
www.eurocor.de

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