State of Colorado Secretary of State

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4. bears the seal / stamp of

est revêtu du sceau / timbre de
y está revestido del sello / timbre de

Deanna Jacobs, Notary Public State of Colorado

Certified

Attesté / Certificado

5. at

à / en

Denver, Colorado

6. the

le / el día

20th day of March, 2017

7. by

par / por

Wayne W. Williams, Secretary of State of the State of Colorado

8. Number

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6027374563

9. Seal / stamp:

Sceau / timbre :
Sello / timbre:

10. Signature:

Signature :
Firma:

Wayne W. Williams

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

To verify the issuance of this Apostille, see: www.sos.state.co.us/auth

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
CERTIFICATE OF AUTHENTICITY
THIS DOCUMENT IS FOR RUSSIA

In the State of Colorado, County Jefferson, I, Deanna Jacobs, a Notary Public in and for said state, do certify that on this 17th day of March, 2017, I carefully compared the original with the attached copy of the Instruction for Use, 41910, Platelet Pooling Set, the copy I now hold in my possession. It is a complete, full, true, and exact copy of the document it purports to reproduce.

Notary Public

Seal:

DEANNA JACOBS
NOTARY PUBLIC
STATE OF COLORADO
NOTARY ID 20164047609
MY COMMISSION EXPIRES DECEMBER 16, 2020

Terumo BCT, Inc.
10811 West Collins Ave.
Lakewood, Colorado 80215-4440
USA
USA Phone: 1.877.339.4228
Phone: +1.303.231.4357
Fax: +1.303.542.5215

Terumo BCT Europe N.V.
Europe, Middle East and Africa
Ikarostraat 41
1930 Zaventem
Belgium
Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

Terumo BCT (Asia Pacific) Ltd.
Room 3903-3903A, 39/F
ACE Tower, Windsor House
311 Gloucester Road
Causeway Bay, Hong Kong
Phone: +852.2283.0700
Fax: +852.2576.1311


Terumo BCT Latin America S.A.
La Pampa 1517 - 12th Floor
C1428DZE
Buenos Aires
Argentina
Phone: +54.11.5530.5200
Fax: +54.11.5530.5201

Terumo BCT Japan, Inc.
20-14, 3-chrome
Higashi Gotanda, Shinagawa-ku
Tokyo 141-0022
Japan
Phone: +81.3.6743.7890
Fax: +81.3.6743.9800

UNLOCKING THE POTENTIAL OF BLOOD | TERUMOBCT.COM

Instructions For Use

41910, Platelet Pooling Set

Warnings

1. Platelet pooling sets are intended for single use only.

Cautions

1. The blood and fluid pathways of the platelet pooling set are sterilized with ethylene oxide and are nonpyrogenic.

If any of the following occurs, do not use the platelet pooling set:

- The tubing has severe kinks.
- The platelet pooling set is incorrectly assembled.

Terumo BCT
10811 W. Collins Ave
Lakewood, CO 80215
USA

BART ROMANUS
Terumo BCT Senior Regulatory Affairs
Specialist
29/08/2016
B. Romanus

- The platelet pooling set is defective or damaged.
 - Any clamps are closed.
 - The end caps are not in place.
2. The Disposable Kits described here contain phthalates (DEHP). Patient groups including pregnant or nursing women and children are considered to be the most at risk of harmful effects due to DEHP exposure. Based on five decades of successful use in medical devices without a reported adverse event related to DEHP, the risk of not doing a needed procedure is far greater than the risk associated with exposure to DEHP.
 3. Protect platelet pooling sets from freezing and excessive heat during storage.
 4. When filtering platelets (and an additive solution if necessary), make sure the tubing is free of kinks or other obstructions. Do not squeeze the bags to increase the draining rate.
 5. Do not seal along the red-striped tubing.

Intended Use

The Platelet Pooling Set is intended for use with whole blood derived platelets, including pooling and filtration, resulting in a platelet product for transfusion. The pooling kit has only been evaluated on the Atreus system and the Reveos system. Therefore, if it will be used with other systems an evaluation will need to be performed according to the institution's SOP.

Interim Platelet Unit Handling

To create a final Pooled Platelets product that is ready for storage, pool four to six interim platelet units (IPUs), and then filter the pooled IPUs. The Platelet Pooling Set allows for simultaneous pooling and filtration of the IPUs and sterile connection of a bag of PAS (platelet additive solution) or plasma, if needed.

There are three options for pooling and filtering the IPUs, which are based on the volume of plasma in the individual interim platelet units.

- Sterile connect a bag of PAS to the platelet pooling set.
- Sterile connect a bag of plasma to the platelet pooling set.
- Use plasma-suspended IPUs, which do not require sterile connection of any additional solutions to the platelet pooling set.

Note: Adjusting the volume of plasma in the IPU requires the assistance of a Terumo BCT implementation specialist.

Follow your institution's SOP to determine whether you should use PAS, plasma, or no additional solutions for platelet pooling and filtration.

Before you begin pooling and filtering the IPUs, you must rest and agitate the individual interim platelet units. The time needed for resting and agitation of the IPUs varies based on when you process the whole blood on the machine. If you process the whole blood unit on the machine the day of collection (2 to approximately 14 hours from collection) follow the IPU handling instructions in Table 1. However, if you process the whole blood unit on the machine the day after collection (approximately 14 to 24 hours from collection) follow the IPU handling instructions in Table 2.

collection (2 to approximately 14 hours from collection) follow the IPU handling instructions in Table 1. However, if you process the whole blood unit on the machine the day after collection (approximately 14 to 24 hours from collection) follow the IPU handling instructions in Table 2.

Note: 14 hours is an approximate value. You should validate the appropriate time range in your own setting.

Table 1: IPU Handling for System Processing the Same Day as WB Collection

Action	Temperature Range	Time
Rest the IPU.	18 °C to 28 °C	Approximately one hour.
Agitate the IPU.	18 °C to 28 °C	Overnight.
Pool and filter the IPU.	18 °C to 28 °C	Within 36 hours of processing the whole blood unit on the machine.

Table 2: IPU Handling for System Processing the Day After WB Collection

Action	Temperature Range	Time
Rest the IPU.	18 °C to 28 °C	Approximately one hour.
Agitate the IPU.	18 °C to 28 °C	A minimum of two hours.
Pool and filter the IPU.	18 °C to 28 °C	Within 36 hours of processing the whole blood unit on the machine.

Note: You will need to validate specific time ranges for overnight, same day, and day after platelet handling in your own setting.

Table 3: Storage Recommendations for the Pooled Platelets Storage Bag

Storage Parameter	Recommended Ranges for Storage in Plasma	Recommended Ranges for Storage in PAS
Volume per bag	100 to 400 mL	100 to 400 mL
Platelet concentration	0.5 to 2.1 x 10 ⁶ /μL	0.5 to 2.1 x 10 ⁶ /μL
Maximum total platelets per bag [(vol. x conc.) x 10 ³]	≤ 5.1 x 10 ¹¹ platelets	≤ 4.5 x 10 ¹¹ platelets
Temperature	22 °C ± 2 °C	22 °C ± 2 °C

Resting and Agitating the Interim Platelet Units

1. Allow each interim platelet unit bag to rest at room temperature for approximately one hour.

Note: Protect the interim platelet units from cool temperatures and air flow vents during resting.

2. After resting, place the interim platelet unit bags on a flatbed rotator and agitate gently. See Table 1 and Table 2.

Note: Do not store interim platelet units in the interim platelet unit bags attached to the processing set. The interim platelet units require additional pooling and filtration before final storage.

Pooling and Filtering the Interim Platelet Units

In the following example instructions, PAS is used to pool and filter the interim platelet units; however, the additive solution you use (if any at all) may vary per your institution's SOP.

1. Prepare a leukoreduction rack.
2. Close all the clamps on the platelet pooling set.
3. Sterile connect the PAS bag to the blue-striped tube; use the minimum amount of tubing needed to complete the sterile connection procedure.
4. Sterile connect 4 to 6 interim platelet units, that have been rested and agitated, to the platelet pooling set using the pooling tubes.
5. Hang all the interim platelet unit bags from a single filtration hook.
6. Hang the PAS bag from a separate filtration hook.

Note: The fluid level of the PAS must be higher than the fluid level of the interim platelet units.

7. Open sterile connection seals on each pooling tube and solution line tube.
8. Ensure that all lines are free of kinks or other obstructions and the leukoreduction filter hangs vertically.
9. After all air bubbles flow into the interim platelet unit bags, open the colored, storage solution line clamp and gently squeeze the PAS bag.
This helps remove any additional air from all the lines and forces the PAS to flow into the interim platelet unit bags.
10. Allow the PAS to flow into the interim platelet unit bags for at least 30 seconds.
11. Ensure that all air bubbles (larger than 1 cm) clear all the lines.
12. Open the clamp on the line to the filter, and allow the mixture of interim platelet units and PAS to flow through the leukoreduction filter and into the pooled platelets storage bag.

Caution: When filtering platelets, make sure the tubing is free of kinks or other obstructions. Do not squeeze the bags to increase the draining rate.

13. After filtration is complete, close all the clamps on the platelet pooling set.

14. Permanently seal the line to the pooled platelets storage bag where the normal tubing meets the red-striped restriction tubing (approximately 7.5 cm from the storage bag).

Caution: Do not seal along the red-striped tubing.

15. Disconnect the pooled platelets storage bag that contains the Pooled Platelets product from the platelet pooling set. (The sample pouch attached to the pooled platelets storage bag should remain intact.)
16. Discard the empty interim platelet unit bags, the empty PAS bag, and the leukoreduction filter assembly according to your institution's SOP.
17. Open the clamp on the line between the pooled platelets storage bag and the sample pouch, and express residual air from the storage bag into the sample pouch, if necessary.
18. Store the Pooled Platelets product according to your institution's SOP.
19. Follow your institution's SOP for sampling the Pooled Platelets product using the sample pouch.

Note: Before dispensing a product sample from the sample pouch, permanently seal the tubing between the pouch and the Pooled Platelets storage bag.

20. Disconnect the sample pouch from the pooled platelets storage bag.
21. Discard the sample pouch according to your institution's SOP.

Symbols and Certification

Symbol	Definition
	Indicates the product quantity when the quantity is placed in the square.
	Indicates the expiration date of the product when accompanied by a specific date.
	Indicates the product lot number when accompanied by a number.
	Indicates the product catalog number when accompanied by a number.
	Indicates that the product is non-pyrogenic.
	Indicates a sterile fluid pathway using ethylene oxide (EtO).

Symbol	Definition
	Indicates that the user must read the instructions for use before application.
	Not for Reuse: Tubing sets manufactured by Terumo BCT bearing the not for reuse symbol are intended for single use only and are not intended to be reprocessed or reused. Terumo BCT cannot ensure the functionality or sterility of the set if it is reprocessed or reused.
	Indicates that the product should not be resterilized.
Rx Only	Indicates that the use of product is by prescription only.
	Device is manufactured in conformance with the European Council Directive 93/42/EEC of June 14, 1993, concerning medical devices. This mark with the Notified Body Number 0086 indicates the approval by the British Standards Institute for these devices within the European Community.
	Indicates the manufacturer of the product when accompanied by the name of the manufacturer.
	Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).
	Indicates that the product must be kept away from sunlight.
	Indicates that the product must be kept dry.
	Indicates the European Authorized Representative of the product when accompanied by the name of the European Authorized Representative.

Return Of Used Product

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department. **IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT.** Please contact your local representative for information regarding returned goods and product complaints.

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142

State of Colorado
Secretary of State

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Pays / País:		United States of America	
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por		Deanna Jacobs	
3. acting in the capacity of agissant en qualité de quien actúa en calidad de		Notary Public, State of Colorado	
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de		Deanna Jacobs, Notary Public State of Colorado	
Certified Attesté / Certificado			
5. at à / en		Denver, Colorado	6. the le / el día
			15th day of February, 2017
7. by par / por		Wayne W. Williams, Secretary of State of the State of Colorado	
8. Number sous n° bajo el número		4743371970	
9. Seal / stamp: Sceau / timbre : Sello / timbre:			10. Signature: Signature : Firma:


This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.
This Apostille does not certify the content of the document for which it was issued.
This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.
To verify the issuance of this Apostille, see: www.sos.state.co.us/auth

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
In the State of Colorado, County Jefferson, I, Deanna Jacobs, a Notary Public in and for said state, do certify that on this 13th day of February, 2017, I carefully compared the original with the attached copy of the Addition to Instructions For Use file for Platelet Pooling Set, the copy I now hold in my possession. It is a complete, full, true, and exact copy of the document it purports to reproduce.

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Position:

Name:

Bart Romanus
Senior Regulatory Affairs Specialist
BARTROMANUS

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
Design and appearance

The **Platelet Pooling Set** is an ethylene oxide sterilized single use disposable set. It contains 6 pooling arms and separate tubing for sterile connection of platelet additive solution or plasma, a platelet leukoreduction filter, a pooled platelet storage bag and a sampling pouch.

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
**Note: The tubing dimensions listed on this document are nominal values and are intended for use in determining compatibility with various types of laboratory equipment. The nominal values are specified target dimensions; however, because they are based on non-sterile tubing measurements, and due to variations in the manufacturing process, the actual dimensions may be slightly different.*

Position	Pre-Sterile Length (mm)		Post Sterile Length	
	Length	Tolerance	Tubing	Between Components
Platelet Pooling Lines				
1	135 mm	+10 / - 5 mm	132 mm	120 mm
Primary Tubing				
2	80 mm	+10 / - 5 mm	80 mm	62 mm
4	80 mm	+10 / - 5 mm	75 mm	57 mm
6	280 mm	+10 / - 5 mm	266 mm	255 mm
7	110 mm	+10 / - 5 mm	104 mm	83 mm
Blue Striped PAS Line				
3	385 mm	+10 / - 5 mm	380 mm	363 mm
Red Striped Restrictive Tubing Line				
5	229 mm	± 6 mm	230 mm	216 mm

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OPERATION CONDITIONS

Operating conditions of a MD

The pooling and leukoreduction of platelet units into a platelet unit for transfusion using the Platelet Pooling Set occurs in a laboratory setting under ambient room temperature conditions.

Requirements for installation and set-up

Not applicable

Information on technical maintenance

The Platelet Pooling Set is a single-use medical device and it does not need technical maintenance

Information on verification

Not applicable, the Platelet Pooling Set does not incorporate a measuring function

Requirements for maintenance or repair

Not applicable, the Platelet Pooling Set is a single-use medical device and it does not need maintenance or repair

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Type of storage, temperature, moisture, illuminance etc.

Protect from freezing. Avoid excessive heat.

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Representative office in Russia:

Terumo Russia LLC

BC Northern Tower, 13th Floor,
10 Testovskaya Street, Moscow 123317,
Russia

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Addition to Instruction For Use

Platelet Pooling Set

Date: July 5, 2017
Signature: Mark Herzberg
Position: Senior Specialist, Regulatory Affairs
Name: Mark Herzberg

Terumo BCT
10811 W. Collins Ave
Lakewood, CO 80215

State of Colorado
County of Jefferson
Subscribed and affirmed before me
this 5th day of July 2017
by Mark Herzberg
Notary Public [Signature]

DEANNA JACOBS
NOTARY PUBLIC
STATE OF COLORADO
NOTARY ID 20184047609
MY COMMISSION EXPIRES DECEMBER 16, 2020

Addition to Instruction for Use

INFORMATION ON MANUFACTURERS OF THE MEDICAL DEVICE

Manufacturer

Terumo BCT, Inc.
10811 W. Collins Avenue
Lakewood, Colorado 80215
USA
Phone: +1.303.231.4357
Fax: +1.303.542.5215

Manufacturing facility

Terumo BCT, Inc.
10811 W. Collins Avenue
Lakewood, Colorado 80215
USA
Phone: +1.303.231.4357
Fax: +1.303.542.5215

THE USE OF THE MEDICAL DEVICE FOR WHICH INTENDED BY THE MANUFACTURER

Intended use

The **Platelet Pooling Set** is intended for use with whole blood derived platelets, including pooling and filtration, resulting in a platelet product for transfusion.

The pooling kit has only been evaluated on the Atrius system and the Reveos system. Therefore, if it will be used with other systems an evaluation will need to be performed according to the institution's SOP.

Indications

The Platelet Pooling Set is to be used for pooling and leukoreduction of whole blood derived platelets resulting in a platelet product for transfusion.

Contraindications

There are no known contraindications

Conditions of use

This product is used by trained and qualified blood center personnel for the pooling and leukoreduction of interim platelet units. This occurs in a laboratory setting.

Adverse Events

The Disposable Kits described here contain phthalates (DEHP). Patient groups including pregnant or nursing women and children are considered to be the most at risk of harmful effects due to DEHP exposure. Based on five decades of successful use in medical devices without a reported adverse event related to DEHP, the risk of not doing a needed procedure is far greater than the risk associated with exposure to DEHP.

CLASSIFICATION OF MEDICAL DEVIICE

The Platelet Pooling Set is a CE-marked class IIb medical device complying with the Medical Device Directive (93/42/EEC) as amended.

GMDN Code: 44034 - Blood Donor Set – Double Pack

Platelet pooling set has no direct contact with human body

TECHNICAL SPECIFICATIONS OF MEDICAL DEVICE

Design and appearance

The **Platelet Pooling Set** is an ethylene oxide sterilized single use disposable set. It contains 6 pooling arms and seperate tubing for sterile connection of platelet additive solution or plasma, a platelet leukoreduction filter, a pooled platelet storage bag and a sampling pouch.

Structure and list of accessories

Catalog Number	Description
CONSUMABLE	
41910	Platelet Pooling Set

Opportunities and ways of integration with other MD

The Platelet Pooling Set has been validated for use with the Atrius 3 Component Integrated Processing Set and the Reveos 3 Component Blood Bag Set and is compatible with the platelet units that would be collected from those sets.

The whole blood-derived platelet units can be pooled and leukoreduced using the Platelet Pooling Set. The platelet pooling set also has a connection for adding plasma or a platelet additive solution.

There are no other known compatible medical devices.

MAIN PARAMETERS AND CHARACTERISTICS OF THE MEDICAL DEVICE


Rated technical characteristics

General Specifications	
Bag Capacity	• Platelet bag: 100mL to 400 mL
Disposable Primary Packaging Dimension	▪ Disposable Tubing Set in header bag: 406mm × 254mm (L × W)

Tubing specifications of Platelet Pooling Set

Item	Wall thickness	Outer diameter	Inner diameter
Blue striped tubing	0,60mm	4,15mm	2,95mm
Clear platelet pooling tubing	0,60mm	4,15mm	2,95mm
Tubing to pooled platelets storage bag	0,56mm	4,06mm	2,95mm

*Note: The tubing dimensions listed on this document are nominal values and are intended for use in determining compatibility with various types of laboratory equipment. The nominal values are specified target dimensions; however, because they are based on non-sterile tubing measurements, and due to variations in the manufacturing process, the actual dimensions may be slightly different.

Position	Pre-Sterile Length (mm)		Post Sterile Length	
	Length	Tolerance	Tubing	Between Components
Platelet Pooling Lines				
1	135 mm	+10 / - 5 mm	132 mm	120 mm
Primary Tubing				
2	80 mm	+10 / - 5 mm	80 mm	62 mm
4	80 mm	+10 / - 5 mm	75 mm	57 mm
6	280 mm	+10 / - 5 mm	266 mm	255 mm
7	110 mm	+10 / - 5 mm	104 mm	83 mm
Blue Striped PAS Line				
3	385 mm	+10 / - 5 mm	380 mm	363 mm
Red Striped Restrictive Tubing Line				
5	229 mm	± 6 mm	230 mm	216 mm

Shelf life

The product expires 2 years from date of manufacture. Refer to product labelling for expiry date

METHODS AND CONDITIONS OF STERILIZATION, THE TERM OF MAINTAINING STERILITY

The Platelet Pooling Set is ethylene oxide sterilized.
This product has been processed with ETO to achieve a 10⁻⁶ sterility assurance level in accordance with ISO 11135-1:2007 *Sterilization of health care products – Ethylene Oxide- Part 1: Requirements for development and routine control of a sterilization process for medical devices*.
This product has been validated to achieve less than 6mg EtO, meeting the requirements of ISO 10993-7:2008 *Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals* prior to release and less then 1ppm in the collect bag at time of use.

If a sterile packaging is damaged, it’s necessary to dispose of the medical device in accordance with local regulations as a class A waste.

OPERATION CONDITIONS

Operating conditions of a MD
The pooling and leukoreduction of platelet units into a platelet unit for transfusion using the Platelet Pooling Set occurs in a laboratory setting under ambient room temperature conditions.

Requirements for installation and set-up
Not applicable

Information on technical maintenance
The Platelet Pooling Set is a single-use medical device and it does not need technical maintenance

Information on verification
Not applicable, the Platelet Pooling Set does not incorporate a measuring function

Requirements for maintenance or repair
Not applicable, the Platelet Pooling Set is a single-use medical device and it does not need maintenance or repair

CONDITIONS OF TRANSPORTATION AND STORAGE

Type of storage, temperature, moisture, illuminance etc.
Protect from freezing. Avoid excessive heat.

Table: Storage Recommendations for platelets in the Pooled Platelets Storage Bag

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WARRANTY LIABILITIES

When purchasing a Terumo BCT product, the general terms and conditions are applicable. This includes the following warranty statement:

LIMITED WARRANTY. Warranty information applicable to the equipment or disposables specified on this invoice is set forth in the price list and any contract for such equipment or disposables and shall be provided by Terumo BCT upon request. For a period of ninety (90) days from the date of the invoice, Terumo BCT will repair or replace Terumo BCT spare parts or optional feature components that Terumo BCT finds defective in material or workmanship. Terumo BCT is not responsible for, and shall have no obligation with respect to any failure caused by normal wear and tear, misuse, unauthorized alterations, accident, neglect, use of non-standard accessory attachments, and/or improper customer maintenance. Terumo BCT does not verify the safety or efficacy of, and makes no warranties, express or implied, with respect to any non-Terumo BCT components or non-Terumo BCT spare parts included at Purchaser's request or any components or products used other than as expressly intended by their manufacturer. Failure to perform preventive maintenance procedures will invalidate this warranty. Terumo BCT is not responsible for any events related to non-Terumo BCT parts installed in the equipment. **THIS LIMITED WARRANTY CONTAINS THE PURCHASER'S EXCLUSIVE WARRANTIES AND REMEDIES.**

INDICATION, IF REQUIRED, OF SPECIAL PRECAUTION MEASURES TO BE TAKEN WHEN DESTROYING UNUSED MEDICAL DEVICES

The Platelet Pooling Set should be disposed of in accordance with local regulations. Waste class "A".

PROCEDURE OF DISPOSABLE (UTILIZATION) AND DESTRUCTION OF THE USED MEDICAL DEVICE

The product packaging complies with European Directive 94/62/EC for packaging and packaging waste. Medical device does not influence to environment.

After use the Platelet Pooling Set may be potentially biohazard and should be disposed of in accordance with local regulations.

Class "A" medical waste in accordance with SanPin 2.1.7.2790-10 for unused products.

Disposal of the Class A medical waste.

Collecting of Class A medical waste can be carried out in reusable containers or in disposable bags. The color of the packages does not matter, however, for this class it is not allowed to use yellow and red packages.

Packages with garbage are located inside reusable containers or on special carts. Containers for waste and carts must be marked as "Waste. Class A". As waste of this class is not potentially dangerous, they can be stored on ordinary landfills or dumped into sewerage.

Class "B" medical waste in accordance with SanPin 2.1.7.2790-10 for the used product.

Disposal of the Class B medical waste.

To collect solid and liquid Class B wastes, disposable, moisture-proof, non-permeable containers must be used. Containers should have a tight fitting lid, which precludes the possibility of spontaneous opening of the container.

Collecting of Class B waste into one-off packages is also permitted. All containers for this type of waste must be yellow, or have a yellow marking "Waste. Class B". The rules for decontamination of medical waste of this class are such that, for this type of waste, decontamination is required, after which they can be burned or buried.

CLAIM

When purchasing a Terumo BCT product, the general terms and conditions are applicable. This includes the following liability statement:

Limited liability: To the extent permitted under applicable law, Purchaser's exclusive remedy with respect to any claim, whether in contract, in tort, under statute, under warranty or otherwise shall be limited solely to the refund of the purchase price or replacement of all Products shown to be other than as warranted. Terumo BCT's liability for losses which do not affect the Product (e.g. indirect, incidental, consequential losses, damages caused by delay) shall be excluded. To the extent the foregoing limitation is not permissible under applicable law, the applicable legal limitation of liability under statute or otherwise if there shall be one, shall apply.

Representative office in Russia:

Terumo Russia LLC
BC Northern Tower, 13th Floor,
10 Testovskaya Street, Moscow 123317,
Russia
Tel: +7-495-988-4740
Fax: +7-495-988-4739

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©2016 Terumo BCT, Inc.