Abbott

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Postfach 1303, 65011 Wiesbaden

Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG., hereby certify that the attached documents that keep the number of pages indicated below are the operational documentation. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: Zaman Khan
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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
INTENDED USE

For use in the calibration of the Microalbumin assay on the ARCHITECT cSystems.

CONTENTS / MATERIALS PROVIDED

REF 2K98-03 Microalbumin Calibrators contain human serum albumin. Sodium azide (0.09%) is added as a preservative.

Contents	Concentration ($\mu\text{g/mL}$ or mg/L)	Quantity
CAL 1	5	1 x 2 mL
CAL 2	25	1 x 2 mL
CAL 3	100	1 x 2 mL
CAL 4	300	1 x 2 mL
CAL 5	500	1 x 2 mL

Results for the Microalbumin assay can be reported as $\mu\text{g/mL}$ or mg/L , which are equivalent.

STANDARDIZATION

The values of the Microalbumin Calibrators are determined by turbidimetric immunoassay using Certified Reference Material 470 (CRM 470)¹ as the primary standard.

PRECAUTIONS

- IVD
- For *In Vitro* Diagnostic Use
- Do not use components beyond the expiration date.
- Do not mix components from different kit lot numbers.

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS / MATERIALS PROVIDED** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these calibrators and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.²⁻⁵

The human albumin used in CAL 1-5 has been tested and found to be nonreactive for HBsAg, anti-HIV-1/HIV-2, and anti-HCV.

The following warnings and precautions apply to: CAL 1-5	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

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STORAGE AND STABILITY

- Unopened Microalbumin Calibrators are stable until the expiration date printed on the label when stored at 2 to 8°C. Do not use past the expiration date.
- Opened calibrators are stable for six months.
- Do not freeze calibrators. Improper storage of calibrators can affect assay performance.

PREPARATION OF CALIBRATOR

The calibrators are supplied in liquid form and are ready to use. Invert bottles several times before use, avoiding the formation of foam.

INSTRUCTIONS FOR USE

Calibration is performed by running a water blank and the Microalbumin Calibrators set. Water for the blank is provided by the instrument.

1. Verify that the correct calibrator values have been entered into the calibration file.
2. Mix bottles several times by gentle inversion.
3. Open the bottles, place an appropriate amount of each calibrator in a separate sample cup, and place in the assigned positions.
4. Cap bottles tightly and return to refrigerated storage immediately after use. Always return each cap to its original bottle.
5. Perform calibration as indicated in the ARCHITECT System Operations Manual, Section 6.
6. Verify calibration by running two levels of controls according to the instructions in the ARCHITECT System Operations Manual.

INDICATIONS OF INSTABILITY OR DETERIORATION

Microalbumin Calibrators are clear, colorless, odorless solutions. Instability or deterioration should be suspected if there are visible signs of precipitates, leakage, turbidity, microbial growth, or if calibration does not meet the appropriate Microalbumin reagent package insert claims and/or the ARCHITECT System Operations Manual criteria.

LIMITATIONS OF THE PROCEDURE

Accurate and reproducible results are dependent upon properly functioning instruments and reagents, storage of product as directed, and good laboratory technique.

BIBLIOGRAPHY

1. Ledue TB, Johnson AM. Commutability of serum protein values: persisting bias among manufacturers using values assigned from the certified reference material 470 (CRM 470) in the United States. *Clin Chem Lab Med* 2001;39(11):1129–1133.
2. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
3. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
4. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
5. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Calibrators 1 through 5
	Calibrator 1
	Calibrator 2
	Calibrator 3
	Calibrator 4
	Calibrator 5
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Distributed in the USA by
	Do not freeze
	Identifies products to be used together
	Information needed for United States of America only
	In Vitro Diagnostic Medical Device
	Lot Number
	Manufactured for
	Product of Japan
	List Number

The ARCHITECT cSystem family of instruments consists of c4000, c8000, and c16000 instruments. ARCHITECT, c4000, c8000, c16000, and cSystem are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

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MULTIGENT Microalbumin Calibrators – 15 months from date of manufacturing.

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POTENTIAL USER

MULTIGENT Microalbumin Calibrators is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the ARCHITECT cSystem.

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MATERIALS OF ANIMAL ORIGIN

Materials of animal origin a
at www.abbottdiagnostics
INDICATIONS A
Not a

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No available data.

COLOR AND ODOR

Calibrator 1-5: Colorless and odorless

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1. Ledue TB, Johnson AM. Commutability of serum protein values: persisting bias among manufacturers using values assigned from the certified reference material 470 (CRM 470) in the United States. *Clin Chem Lab Med* 2001;39(11):1129-1133.
2. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
3. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
4. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
5. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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Signed: Za
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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
INTENDED USE

For use in the quality control of the Microalbumin assay on the ARCHITECT cSystems.

CONTENTS / MATERIALS PROVIDED

REF 2K98-11 Microalbumin Controls contain human serum albumin. Sodium azide (0.09%) is added as a preservative.

Contents	Concentration ($\mu\text{g/mL}$ or mg/L)	Quantity
CONTROL 1	20 to 40	2 x 2 mL
CONTROL 2	75 to 105	2 x 2 mL

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-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS / MATERIALS PROVIDED** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these controls and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴

The human albumin used in CONTROL 1-2 has been tested and found to be nonreactive for HBsAg, anti-HIV-1/HIV-2, and anti-HCV.

The following warnings and precautions apply to: CONTROL 1-2	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

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STORAGE AND STABILITY

- Unopened Microalbumin Controls are stable until the expiration date printed on the label when stored at 2 to 8°C. Do not use past the expiration date.
- Opened controls are stable for six months.
- Do not freeze controls. Improper storage of controls can affect assay performance.

PREPARATION OF CONTROLS

The controls are supplied in liquid form and are ready to use. Invert bottles several times before use, avoiding the formation of foam.

INSTRUCTIONS FOR USE

1. Verify that the correct control values have been entered into the control file.
2. Mix bottles several times by gentle inversion.
3. Open the bottles, place an appropriate amount of each control in a separate sample cup, and place in the assigned positions.
4. Cap bottles tightly and return to refrigerated storage immediately after use. Always return each cap to its original bottle.
5. Run controls according to the instructions in the ARCHITECT System Operations Manual, Section 5.
6. Verify control results meet acceptance criteria before reporting patient results.

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Microalbumin Controls are clear, colorless, odorless solutions. Instability or deterioration should be suspected if there are visible signs of precipitates, leakage, turbidity, microbial growth, or if control does not meet the appropriate Microalbumin reagent package insert and/or the ARCHITECT System Operations Manual criteria.

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BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL 1	Control 1
CONTROL 1-2	Control 1-2
CONTROL 2	Control 2
DISTRIBUTED IN THE USA BY	Distributed in the USA by
DO NOT FREEZE	Do not freeze
EC REP	Authorized Representative in the European Community
FOR USE WITH	Identifies products to be used together
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
MANUFACTURED FOR	Manufactured for
PRODUCT OF JAPAN	Product of Japan
REF	List Number

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ARCHITECT, c4000, c8000, c16000, and cSystem are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
 Max-Planck-Ring 2
 65205 Wiesbaden
 Germany
 +49-6122-580

DISTRIBUTED IN THE USA BY

Abbott Laboratories
 Abbott Park, IL 60064 USA

MANUFACTURED FOR

Abbott Laboratories

Customer Service: Contact your local representative
 or find country-specific contact information on
www.abbottdiagnostics.com

Revised February 2017.

©2016, 2017 Abbott Laboratories

ADDENDUM

to the Instruction for Use MULTIGENT Microalbumin Controls for Russia

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MULTIGENT Microalbumin Controls – 15 months from date of manufacturing.

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POTENTIAL USER

MULTIGENT Microalbumin Controls is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the ARCHITECT cSystem.

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Not applicable.

TECHNICAL CHARACTERISTICS

- Control 1-2:
 - o pH-value at 20 °C: 7
 - o Density at 20 °C: 1g/cm³

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Signed: ZK
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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Microalbumin
2K98
307202/R03
B2K9R0

Revised February 2017.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

Microalbumin

INTENDED USE

The Microalbumin assay is used for the quantitative measurement of albumin in human urine on the ARCHITECT cSystems.

SUMMARY AND EXPLANATION OF THE TEST

Microalbuminuria is a condition characterized by increased urinary excretion of albumin in the absence of overt nephropathy, and can be used to predict diabetic nephropathy.^{1,2} Diabetic nephropathy is a major cause of death in individuals with insulin-dependent diabetes; and because it is accompanied by irreversible kidney damage and persistent proteinuria, it is a major factor in the decision to initiate hemodialysis.^{3,4}

Early detection of glomerular damage, when it is minimal and reversible, is extremely important. Monitoring urinary microalbumin is an important component of treatment for both Type I and Type II diabetes mellitus.³ Methods for monitoring microalbuminuria include measurement of protein excretion in 24-hour, timed, or overnight collections, and determination of the albumin:creatinine ratio in an untimed "spot" urine specimen. Twenty-four hour and timed urine collections may be associated with collection errors including improper timing, missed samples, and incomplete bladder emptying. The concentration of protein in a spot urine sample provides an estimate of the protein excretion rate, but is affected by patient hydration. The ratio of protein or albumin to creatinine in a spot urine sample corrects for variations in hydration and avoids the sources of error associated with 24-hour and timed urine collection.⁵

PRINCIPLES OF THE PROCEDURE

The Microalbumin assay is a turbidimetric immunoassay that uses polyclonal antibodies against human albumin. When a specimen is mixed with the reagents, albumin in the specimen combines with the anti-human albumin antibody (goat) in the reagent to yield an insoluble aggregate that causes increased turbidity in the solution. The degree of turbidity is proportional to the concentration of albumin in the specimen, and can be measured optically.

Methodology: Turbidimetric/Immunoassay

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Kit Contents

Microalbumin 2K98.

Supplied as a liquid, ready-to use two-reagent kit.

REF	2K98-21
Σ	500*
R1	2 x 53 mL
R2	2 x 12 mL

*Calculation is based on the minimum reagent fill volume per kit.

Reactive Ingredients	Concentration
R1	
Good's buffer	1.03%
Sodium chloride	1.48%
Sodium hydroxide	< 0.14%
R2	
TRIS buffer	1.17%
Anti-human albumin antibody (goat)	0.17%
Sodium chloride	1.14%
Hydrochloric acid	< 0.9%

Inactive Ingredients: R1 and R2 contain 0.09% sodium azide.

Warnings and Precautions

- IVD
- For *In Vitro* Diagnostic Use
- R2 Contains nonsterile goat polyclonal antibodies.

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.⁶⁻⁹

The following warnings and precautions apply to: R1 and R2	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagents beyond the expiration date.
- Do not pool reagents within a kit or between kits.
- Do not use components from one lot with components from another lot.
- When either the R1 or R2 reagent cartridge becomes empty, replace both cartridges.
- Before use, invert reagent cartridges several times to mix well avoiding the formation of bubbles.
- Remove any air bubbles present in the reagents with a new applicator stick, or allow the reagents to sit at the appropriate storage temperature to allow the bubbles to dissipate. To minimize volume depletion, do not use a transfer pipette to remove bubbles.

CAUTION: Bubbles may interfere with proper detection of reagent level in the cartridge and cause insufficient reagent aspiration which could impact results.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2-8°C	Until expiration date	Do not freeze.
On board	System temperature	28 days (672 hours)	After 28 days (672 hours), the reagent kit must be discarded. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

Reagents may be stored on or off the ARCHITECT cSystem. If reagents are removed from the system, store at 2-8°C (with replacement caps) in their original boxes. When reagent is placed back on the system, run controls and if appropriate criteria are not met, recalibration may be required. For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage or contamination, turbidity, or if calibration or controls do not meet the appropriate package insert and/or ARCHITECT System Operations Manual criteria.

For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Microalbumin assay file must be installed on the ARCHITECT cSystem prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters or for a detailed description of system procedures, refer to the ARCHITECT System Operations Manual, Section 5.

Alternate Result Units

Results for the Microalbumin assay can be reported as µg/mL or mg/L, which are equivalent. In this package insert, results are expressed as µg/mL.

When converting to units other than those listed, refer to the ARCHITECT System Operations Manual, Section 2.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

Verified specimen types to be used with this assay.

Specimen Type	Collection Vessel
Urine (Collect urine as a 24 hour, timed, or random midstream sample (spot collection))	Clean plastic or glass container

Other specimen types and collection tube types/anticoagulants have not been verified with this assay.

The following are acceptable preservatives: 6N hydrochloric acid, acetic acid, chloroform, formalin, toluene, and xylene.

A urine creatinine result is required to calculate a spot collection microalbumin-to-creatinine ratio.

The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Preparation for Analysis

Stored specimens must be inspected for particulates. If particulates are present, mix and centrifuge the specimen to remove particulates prior to testing.

For total sample volume requirements, refer to the ARCHITECT System Operations Manual, Section 5.

Specimen Storage

Analyze fresh specimens if possible.

Avoid repeated freeze/thaw cycles.

Specimen Type	Storage Temperature	Maximum Storage Time
Urine	4°C	2 weeks
	-70°C prior to testing ³	5 months
	4°C	6 days ¹⁰

For timed or 24 hour samples, analyze specimens as soon as possible following collection.

Specimens may be stored for up to 2 weeks at 4°C or up to 5 months at -70°C prior to testing.³

For spot samples, which require the analysis of both albumin and creatinine, analyze fresh specimens if possible; if not, specimens may be stored for up to 6 days at 4°C.¹⁰

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PROCEDURE

Materials Provided

2K98 Microalbumin Reagent Kit

Materials Required but not Provided

- REF 2K98-03 Microalbumin Calibrators
- REF 2K98-11 Microalbumin Controls
- Saline (0.85% to 0.90% NaCl) for specimen dilution

For information on materials required for operation of the instrument, refer to the ARCHITECT System Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the ARCHITECT System Operations Manual, Section 5.

For sample volume requirements, refer to the ASSAY PARAMETERS section of this package insert. Minimum sample volume is calculated by the system and printed on the Orderlist report. Ensure adequate sample volume is present prior to running the test.

Specimen Dilution Procedures

The ARCHITECT cSystems have an automatic dilution feature; for additional information, refer to the ARCHITECT System Operations Manual, Section 2.

Specimens with albumin concentrations values exceeding the 500 µg/mL are flagged and may be diluted by following the Manual Dilution Procedure, or the Automatic Dilution Protocol provided in the assay parameters. If an Automatic Dilution Protocol is not provided, refer to the ARCHITECT System Operations Manual, Section 2 for configuration information, and verify results according to your laboratory's standard operating procedures.

Automated Dilution Protocol

When using the Automated Dilution Protocol, the system performs a 1:4 dilution of the specimen and automatically corrects the concentration by multiplying the result by the appropriate dilution factor.

Manual Dilution Procedure

- 1. Dilute the specimen with saline (0.85% to 0.90% NaCl).
- 2. Enter the dilution factor in the Patient or Control order screen.
The system uses this dilution factor to automatically correct the concentration by multiplying the result by the entered factor. If the operator does not enter the dilution factor, the result must be manually multiplied by the appropriate dilution factor before reporting the result.

Manual Dilution Factor = (Volume of Specimen + Volume of Dilution Reagent) / Volume of Specimen

If a diluted specimen result is flagged indicating it is less than the linear low limit, do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

Calibration

The Microalbumin assay must be calibrated using the Microalbumin Calibrators [REF] 2K98-03 (supplied separately).

CAUTION: The Microalbumin assay uses polyclonal antibodies against human albumin that can react with closely related proteins such as bovine albumin. **Therefore, calibrators containing non-human protein should not be used with this assay.** The Microalbumin Calibrators contain only human albumin. Calibration is stable for approximately 28 days (672 hours), but is required with each change in reagent lot. Verify calibration curve with at least two levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary. For specific calibration instructions and information on calibrator standardization, refer to the Microalbumin Calibrators package insert. For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

Quality control for the Microalbumin assay must be verified using the Microalbumin Controls [REF] 2K98-11 (supplied separately). **CAUTION:** The Microalbumin assay uses polyclonal antibodies against human albumin that can react with closely related proteins such as bovine albumin. This cross-reactivity may vary among reagent lots and cause shifts in control values. **Therefore, control material containing non-human protein should not be used with this assay.** The Microalbumin Controls contain only human albumin and are not affected by cross-reactivity. As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Two levels of controls spanning the medical decision range are to be run once during each day of assay use.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, patient values may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

RESULTS

Representative performance data are given in the EXPECTED VALUES and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert. Results obtained in individual laboratories may be different from the data presented.

Calculation

For additional information on results calculations, refer to the ARCHITECT System Operations Manual, Appendix C.

Interpretation of Results

As with all analyte determinations, the microalbumin value should be used in conjunction with information available from clinical evaluation and other diagnostic procedures.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

LIMITATIONS OF THE PROCEDURE

In some instances, falsely high or low results occur due to non-specific turbidity. Prozone may be observed at albumin concentrations greater than 10,000 µg/mL. If a result is questionable, dilute the sample and repeat the analysis. Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert.

EXPECTED VALUES

According to the American Diabetes Association (ADA) microalbumin levels can be measured from 24 hour, timed, or random urine specimens. Normal albumin levels, microalbuminuria, and macro (clinical) albuminuria levels for each specimen type are listed below.^{11, 12}

	Spot collection ^a		24 hour collection ^b	Timed collection ^c
	Microalbumin:Creatinine ratio			
Units	(µg/mg)	(mg/mmol)	(mg/24 hr)	(µg/min)
	or (mg/g)			
Normal	< 30	< 2.5 (male) < 3.5 (female)	< 30	< 20
Microalbuminuria	30 to 299	2.5 to 29 (male) 3.5 to 29 (female)	30 to 299	20 to 199
Macro (clinical) albuminuria	≥ 300	> 30 (male) > 30 (female)	≥ 300	≥ 200

a. Spot collection = (uAlb µg/mL ÷ Crea-U mg/dL) × 100 mL / dL

or uAlb mg/L ÷ Crea-U mmol/L

b. 24 hour collection = (uAlb µg/mL × Volume mL) ÷ 1,000 µg / mg

c. Timed collection = (uAlb µg/mL × Volume mL) ÷ Time (min)

Due to variability in urinary albumin excretion, at least two of three test results measured within a 6-month period should show elevated levels before a patient is designated as having microalbuminuria.¹³ Exercise within 24 hours, infection, fever, congestive heart failure, marked hyperglycemia, and marked hypertension may elevate urinary albumin excretion over baseline values.¹¹

SPECIFIC PERFORMANCE CHARACTERISTICS

Reportable Range

The reportable range of the Microalbumin assay is 5 to 500 µg/mL (5 to 500 mg/L) based on accuracy by recovery.

Accuracy by Recovery

Based on guidance from the Clinical and Laboratory Standards Institute (CLSI) protocol NCCLS EP6-A,¹⁴ samples were assayed using the Microalbumin reagent. Samples were prepared by dilution of pooled urine stock with saline or by mixing combinations of calibrators. The reportable range was based on mean results within 10% of the theoretical concentration. Data from this study are summarized below.

Theoretical Concentration	Mean Recovered Concentration	Variation	
(µg/mL)	(µg/mL)	(µg/mL)	(%)*
5.00	5.30	0.30	6.00
11.67	12.20	0.53	4.54
18.33	19.05	0.72	3.93
46.60	45.10	-1.50	-3.22
93.19	90.90	-2.29	-2.46
186.39	195.20	8.81	4.73
279.58	279.58	0.00	0.00
372.77	397.60	24.83	6.66
465.97	482.50	16.53	3.55

* % Variation = $\frac{\text{Mean Recovered Concentration} - \text{Theoretical Concentration}}{\text{Theoretical Concentration}} \times 100$

Limit of Detection (LOD)

The LOD of the Microalbumin assay is 1 µg/mL (1 mg/L). LOD was determined as the concentration equivalent to the mean absorbance of 21 replicates of an analyte-free sample + 2 SD, where SD = the pooled within-run standard deviation.

Interference

A study based on CLSI protocol NCCLS EP7-P¹⁵ was performed using urine specimens containing approximately 27 µg/mL microalbumin supplemented with the individual compounds listed below.

Potential interference in the Microalbumin assay is ≤ 10% for the urine preservatives and interferents listed below when used in the amounts shown.

Preservatives and Interferents	Concentration	
	Conventional Units	SI Units
6N Hydrochloric acid	4 mL/dL	40 mL/L
Acetic acid	1,000 µL/dL	10 mL/L
Chloroform	1,000 µL/dL	10 mL/L
Formalin	1,000 µL/dL	10 mL/L
Thymol	50 mg/dL	500 mg/L
Toluene	1,000 µL/dL	10 mL/L
Xylene	1,000 µL/dL	10 mL/L
Ascorbic acid	400 mg/dL	22,712 µmol/L
Bilirubin	25 mg/dL	427.5 µmol/L
Calcium	400 mg/dL	100 mmol/L
Creatinine	400 mg/dL	35,360 µmol/L
Glucose	4,000 mg/dL	222 mmol/L
Hemoglobin	500 mg/dL	5 g/L
Hippuric acid	400 mg/dL	4 g/L
Inorganic phosphorus	400 mg/dL	130 mmol/L
Potassium chloride	1,000 mg/dL	10 g/L
Sodium chloride	2,000 mg/dL	20 g/L
Urea nitrogen	400 mg/dL	142.8 mmol/L
Uric acid	100 mg/dL	5.9 mmol/L
Urobilinogen	20 mg/dL	33.8 µmol/L

Urine samples ranging in pH levels from 3 to 10 demonstrated ≤ 10% interference.

Precision

The Microalbumin assay precision is ≤ 5% CV. Studies were performed using the Microalbumin assay based on guidance from CLSI protocol NCCLS EP5-A.¹⁶ Two levels of microalbumin controls were used, and each level of control was assayed in duplicate twice a day for 21 days. Representative data from these studies are summarized below.

Control		Level 1	Level 2
N		84	84
Mean (µg/mL)		30.15	90.93
Within Run	SD	0.2219	0.6764
	%CV	0.74	0.74
Between Day	SD	0.000	0.000
	%CV	0.00	0.00
Total	SD	1.1429	1.8850
	%CV	3.79	2.07

Method Comparison

Correlation studies were performed based on CLSI protocol NCCLS EP9-A.¹⁷

Results from the Microalbumin assay on the AEROSSET System were compared with the results from a commercially available turbidimetric/immunoturbidimetric methodology.

Results from the Microalbumin assay on an ARCHITECT cSystem were compared with the results on the AEROSSET System.

	AEROSSET vs. Comparative	
	Method	ARCHITECT vs. AEROSSET
N	53	53
Y-Intercept	2.89	1.30
	95% CI: 0.384 to 5.389	95% CI: -0.297 to 2.897
Correlation Coefficient	0.9984	0.9994
	1.023	1.019
Slope	95% CI: 1.007 to 1.039	95% CI: 1.009 to 1.029
Range (µg/mL)	5.3 to 459.4	7.6 to 458.2
Average Bias	8%	9.5%
Standard Error of the Estimate (Sy.x)	6.99	4.41

95% CI = 95% Confidence Interval

BIBLIOGRAPHY

1. Mogensen CE, Christensen CK, Vittinghus E. The stages in diabetic renal disease. With emphasis on the stage of incipient diabetic nephropathy. *Diabetes* 1983;32 Suppl 2:64-78.
2. Viberti GC, Mackintosh D, Bilous RW, et al. Proteinuria in diabetes mellitus: role of spontaneous and experimental variation of glycemia. *Kidney Int* 1982;21:714-720.
3. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders, 1999:798-799.
4. Rosenstock J, Raskin P. Early diabetic nephropathy: assessment and potential therapeutic interventions. *Diabetes Care* 1986;9: 529-545.
5. NKF K/DOQI Guidelines Web Site. K/DOQI clinical practice guidelines for chronic kidney disease: evaluation, classification, and stratification. Part 5. Evaluation of laboratory measurements for clinical assessments of kidney disease. 2002. Available at: http://www.kidney.org/professionals/kdoqi/guidelines_ckd/p5_lab_g5.htm. Accessed October 11, 2007.
6. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
7. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
8. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
9. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

10. Guder WG, Narayanan S, Wisser H, et al. List of analytes—preanalytical variables. Annex In: *Samples: From the Patient to the Laboratory*. Darmstadt, Germany: GIT Verlag; 1996:Annex 12–3, 36–37.
11. Molitch ME, DeFronzo RA, Franz MJ, et al. Nephropathy in Diabetes. *Diabetes Care* 2004;27 Suppl 1:S79–83.
12. McIntosh A, Hutchinson A, Marshall S, et al. Background to renal disease in type 2 diabetes. In: *Clinical Guidelines for Type 2 Diabetes*. Sheffield, UK: University of Sheffield; 2002:27–36.
13. American Diabetes Association. Standards of medical care for patients with diabetes mellitus. *Diabetes Care* 2003;26 Suppl 1: S33–50.
14. Tholen DW, Kroll M, Astles JR, et al. Evaluation of the Linearity of Quantitative Measurement Procedures: *A Statistical Approach; Approved Guideline (EP6-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards; 2003.
15. National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
16. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline*. NCCLS Document EP5-A. Wayne, PA: NCCLS; 1999.
17. National Committee for Clinical Laboratory Standards (NCCLS). *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline*. NCCLS Document (EP9-A). Wayne, PA: NCCLS; 1995.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
DISTRIBUTED IN THE USA BY	Distributed in the USA by
DO NOT FREEZE	Do not freeze
EC REP	Authorized Representative in the European Community
FOR USE WITH	Identifies products to be used together
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
MANUFACTURED FOR	Manufactured for
PRODUCT OF JAPAN	Product of Japan
R1	Reagent 1
R2	Reagent 2
REF	List Number
SN	Serial number

The ARCHITECT cSystem family of instruments consists of c4000, c8000, and c16000 instruments.

AEROSET, ARCHITECT, c4000, c8000, c16000, cSystem, and SmartWash are trademarks of Abbott Laboratories in various jurisdictions.

All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

MANUFACTURED FOR

Abbott Laboratories

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised February 2017.

©2016, 2017 Abbott Laboratories

ASSAY PARAMETERS

Microalbumin—Conventional and SI Units

Configure assay parameters — General			
☒ General	☐ Calibration	☐ SmartWash	☐ Results
Assay: uAlb		Type: Photometric	Version: †
Number: 2839			
Run controls for onboard reagents by: Lot			
☒ Reaction definition	☐ Reagent / Sample	☐ Validity checks	
Reaction mode: End up			
Primary		Secondary	Read times
Wavelength: 340 / 700		Main: 31 – 33	
Last required read: 33			
Absorbance range: — —		Color correction: — —	
Sample blank type: Self		Blank: 14 – 16	

☐ Reaction definition		☒ Reagent / Sample		☐ Validity checks	
Reagent: UALB0		Reagent volume: 180		R1	R2
Diluent: Saline		Water volume: —			
Diluent dispense mode: Type 0		Dispense mode: Type 0		Type 0	Type 0
Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor
STANDARD	6.0	—	—	—	1:1.00
DII 1	25.0	6.0	75	—	1:4.00

☐ Reaction definition		☐ Reagent / Sample		☒ Validity checks	
Reaction check: None					
Maximum absorbance variation: —					

Configure assay parameters — Calibration			
☐ General	☒ Calibration	☐ SmartWash	☐ Results
Assay: uAlb		Calibration method: Spline	
☒ Calibrators	☐ Volumes	☐ Intervals	☐ Validity checks
Calibrator set: uAlbCal			
Blank: Water		Calibrator level: 0††	
Cal 1: uAlbCal1		Concentration: ‡	
Cal 2: uAlbCal2		‡	
Cal 3: uAlbCal3		‡	
Cal 4: uAlbCal4		‡	
Cal 5: uAlbCal5		‡	
Replicates: 3 [Range 1 – 3]			

☐ Calibrators		☒ Volumes		☐ Intervals		☐ Validity checks	
Calibrator: uAlbCal							
Blank: Water		Sample	Diluted sample	Diluent	Water		
Cal 1: uAlbCal1		6.0	—	—	—		
Cal 2: uAlbCal2		6.0	—	—	—		
Cal 3: uAlbCal3		6.0	—	—	—		
Cal 4: uAlbCal4		6.0	—	—	—		
Cal 5: uAlbCal5		6.0	—	—	—		

☐ Calibrators		☐ Volumes		☒ Intervals		☐ Validity checks	
Calibration intervals:							
Full interval: 672		(hours)					
Calibration type:							
Adjust type: None							

☐ Calibrators		☐ Volumes		☐ Intervals		☒ Validity checks	
Blank absorbance range: — —							
Span: Blank		Blank					
Span absorbance range: — —							
Expected cal factor: 0.00							
Expected cal factor tolerance %: 0							

Configure assay parameters — SmartWash				
☐ General	☐ Calibration	☒ SmartWash	☐ Results	☐ Interpretation
Assay: uAlb				
COMPONENT	REAGENT / ASSAY	WASH	Volume	Replicates
R1	AMIK9	Detergent A	345	1
R1	VANCO	Detergent A	345	1
R1	GENT9	Detergent A	345	1
R1	DIG00	10% Detergent B	345	1
R1	TOBRA	10% Detergent B	345	1
R1	DGT0B	10% Detergent B	345	1
R2	AMIK9	Detergent A	345	1
R2	VANCO	Detergent A	345	1
R2	GENT9	Detergent A	345	1
R2	DIG00	10% Detergent B	345	1
R2	TOBRA	10% Detergent B	345	1
R2	DGT0B	10% Detergent B	345	1
R2	CREAC	Detergent A	345	1
R2	FEPOS	Detergent A	345	1
Sample Probe		Detergent A		
Cuvette	Trig**	10% Detergent B	345	

**Not required for ARCHITECT Software version 7.00 and above.

Microalbumin—Conventional Units

Configure assay parameters — Results			
☐ General	☐ Calibration	☐ SmartWash	☒ Results
Assay: uAlb		Assay number: 2839	
Dilution default range:		Result units: µg/mL	
Low-Linearity: 5			
High-Linearity: 500			
Gender and age specific ranges:			
GENDER	AGE (UNITS)	NORMAL	EXTREME

Configure result units	
Assay: uAlb	
Version: †	
Result units: µg/mL	
Decimal places: 0	[Range 0 – 4]
Correlation factor: 1.0000	
Intercept: 0.0000	

Microalbumin—SI Units

Configure assay parameters — Results			
☐ General	☐ Calibration	☐ SmartWash	☒ Results
Assay: uAlb		Assay number: 2839	
Dilution default range:		Result units: mg/L	
Low-Linearity: 5			
High-Linearity: 500			
Gender and age specific ranges:			
GENDER	AGE (UNITS)	NORMAL	EXTREME

Configure result units	
Assay: uAlb	
Version: †	
Result units: mg/L	
Decimal places: 0	[Range 0 – 4]
Correlation factor: 1.0000	
Intercept: 0.0000	

† Due to differences in instrument systems and unit configurations, version numbers may vary.

†† Displays the number of decimal places defined in the decimal places parameter field.

‡ Refer to the concentration specified on calibrator label or package insert. In ARCHITECT Software version 5.00 and above, these values are defined on the Configure calibrator set screen.

ADDENDUM

to the Instruction for Use MULTIGENT Microalbumin for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for devices are 2-8°C. Medical devices must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Medical devices must be carried out in single-use thermal containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in refrigerated trucks using thermal indicators.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. Manufacturer set the maximum shelf life of

MULTIGENT Microalbumin – 24 months from date of manufacturing.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste". Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

EN ISO13485:2016; ISO13485:2003; EN ISO14971:2012; (EC) No 1272/2008; DIN EN 18113-1-5:2011; DIN EN 13612:2002; DIN EN 15193:2009; DIN EN 17511:2003; DIN EN 23640:2015; DIN EN 62366:2008; EN 980:2008; EN ISO 15223-1:2012; ISO 13612-1:2014; EN ISO 14155:2011

POTENTIAL USER

MULTIGENT Microalbumin is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the ARCHITECT cSystem.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRAINDICATIONS FOR USE

Not applicable.

TECHNICAL CHARACTERISTICS

- Reagent 1 (R1): pH-value at 20 °C: 7.4
- Reagent 2 (R2): pH-value at 20 °C: 7.4

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

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1. Mogensen CE, Christensen CK, Vittinghus E. The stages in diabetic renal disease. With emphasis on the stage of incipient diabetic nephropathy. *Diabetes*, 1983;32 Suppl 2:64-78.
2. Viberti GC, Mackintosh D, Bilous RW, et al. Proteinuria in diabetes mellitus: role of spontaneous and experimental variation of glycemia. *Kidney Int* 1982;21:714-720.

3. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders, 1999:798-799.
4. Rosenstock J, Raskin P. Early diabetic nephropathy: assessment and potential therapeutic interventions. *Diabetes Care* 1986;9: 529-545.
5. NKF K/DOQI Guidelines Web Site. K/DOQI clinical practice guidelines for chronic kidney disease: evaluation, classification, and stratification. Part 5. Evaluation of laboratory measurements for clinical assessments of kidney disease. 2002. Available at: http://www.kidney.org/professionals/kdoqi/guidelines_ckd/p5_lab_g5.htm. По состоянию на 11 октября 2007 г.
6. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
7. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
8. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
9. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
10. Guder WG, Narayanan S, Wisser H, et al. List of analytes—preanalytical variables. Annex In: *Samples: From the Patient to the Laboratory*. Darmstadt, Germany: GIT Verlag; 1996:Annex 12-3, 36-37.
11. Molitch ME, DeFronzo RA, Franz MJ, et al. Nephropathy in Diabetes. *Diabetes Care* 2004;27 Suppl 1:S79-83.
12. McIntosh A, Hutchinson A, Marshall S, et al. Background to renal disease in type 2 diabetes. In: *Clinical Guidelines for Type 2 Diabetes*. Sheffield, UK: University of Sheffield; 2002:27-36.
13. American Diabetes Association. Standards of medical care for patients with diabetes mellitus. *Diabetes Care* 2003;26 Suppl 1: S33-50.
14. Tholen DW, Kroll M, Astles JR, et al. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; *Approved Guideline (EP6-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 2003.
15. National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
16. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline*. NCCLS Document EP5-A. Wayne, PA: NCCLS; 1999.
17. National Committee for Clinical Laboratory Standards (NCCLS). *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline*. NCCLS Document (EP9-A). Wayne, PA: NCCLS; 1995.

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