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Alinity i Insulin Reagent Kit

Read Important Changes: Revised April 2014


Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions. The insulin Assay value of a given specimen is determined with assay from different manufacturers, can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the insulin assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining insulin is changed, additional sequential testing should be performed to confirm

Results. This assay shows no cross-reactivity with proinsulin (1-3-33) or IGF analogs. Assayable proinsulin interference is thought to be the insulin antibodies, which (quantity) is unknown, caused by insulin or proinsulin insulin. Supplementary information first, for the diagnosis of diabetes mellitus and second, for differential diagnosis of fasting hypoglycemia, in these applications, the ratio of immunoreactive insulin to blood glucose (IG) may be more valuable than the insulin level alone. Furthermore, a single random blood sample may provide insufficient information due to wide variations in the time responses of insulin

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Подпись _____
Limited liability company
Abbott Laboratories

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Abbott GmbH & Co. KG
DENKA SEIKEN CO., LTD. Tokyo, Japan