Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Postfach 1303, 65011 Wiesbaden

Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

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Signed: [Signature]
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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Warnings and Precautions

- IVD
- For In Vitro Diagnostic Use
- Rx ONLY

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Representative performance data are provided in this section.

Results obtained in individual laboratories may vary.

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

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The Alinity c analyzer and the ARCHITECT c System and AEROSET System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity c analyzer.

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Potential interference in the Direct Bilirubin assay from 62 mg/dL (0.62 g/L) hemoglobin, 125 mg/dL (1.25 g/L) Intralipid, or 0.50 mmol/L Indican (indoxyl sulfate) is $\leq 10\%$ or ± 0.1 mg/dL, whichever is greater, at the medical decision level of the analyte.

This study was performed on the ARCHITECT c System and the AEROSET System.

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	Units	n	Correlation Coefficient		Concentration Range	
			Intercept	Slope		
Alinity c Direct Bilirubin vs ARCHITECT Direct Bilirubin	Serum mg/dL	132	1.00	0.00	1.00	0.2 to 14.7
	μmol/L	132	1.00	0.00	1.00	3.4 to 250.6





BIBLIOGRAPHY

1. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1136.
2. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
3. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
4. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
5. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
6. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1169.
7. Guder WG, Narayanan S, Wissner H, et al. List of Analytes—preanalytical variables. Annex In: *Samples: From the Patient to the Laboratory*. Darmstadt: GIT Verlag; 1996:Annex 8–9.
8. US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*. 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
9. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press; 1997:3-88.
10. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
11. Donnachie EM, Seccombe DW, Urquhart NI, et al. Indocyanine green interference in the Kodak Ektachem determination of total bilirubin. *Clin Chem* 1989;35(5):899-900.
12. Meijer DKF, Weert B, Vermeer GA. Pharmacokinetics of biliary excretion in man. VI. Indocyanine green. *Eur J Clin Pharmacol* 1988;35:295-303.
13. Data on file at Abbott Laboratories.
14. National Committee for Clinical Laboratory Standards (NCCLS). *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline*. NCCLS Document C28-A. Villanova, PA: NCCLS; 1995.
15. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
16. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
17. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
18. National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
19. Taki K, Tsuruta Y, Niwa T. Indoxyl sulfate and atherosclerotic risk factors in hemodialysis patients. *Am J Nephrol* 2007; 27:30–35.
20. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3-105–3-108.
21. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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Key to Symbols

ISO 15223 Symbols	
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
	In Vitro Diagnostic Medical Device
	Lot Number
	List Number
	Serial number
Other Symbols	
	Distributed in the USA by
	Information needed for United States of America only
	Produced for Abbott by
	Product of Canada
	Reagent 1
	Reagent 2
	For use by or on the order of a physician only (applicable to USA classification only).

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Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown
Prince Edward Island
C1E 2B9 Canada

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MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottiagnostics.com or contact your local representative.

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No available data.

COLOUR AND ODOUR

Alinity c Direct Bilirubin Reagent 1: Odourless and Colourless
Alinity c Direct Bilirubin Reagent 2: Odourless and Colourless

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- Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1136.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1169.
- Guder WG, Narayanan S, Wisser H, et al. List of Analytes—preanalytical variables. Annex In: *Samples: From the Patient to the Laboratory*. Darmstadt: GIT Verlag; 1996:Annex 8–9.
- US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*. 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
- Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press; 1997:3-88.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Donnachie EM, Seccombe DW, Urquhart NI, et al. Indocyanine green interference in the Kodak Ektachem determination of total bilirubin. *Clin Chem* 1989;35(5):899-900.
- Meijer DKF, Weert B, Vermeer GA. Pharmacokinetics of biliary excretion in man. VI. Indocyanine green. *Eur J Clin Pharmacol* 1988;35:295-303.
- Data on file at Abbott Laboratories.
- National Committee for Clinical Laboratory Standards (NCCLS). *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline*. NCCLS Document C28-A. Villanova, PA: NCCLS; 1995.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
- Taki K, Tsuruta Y, Niwa T. Indoxyl sulfate and atherosclerotic risk factors in hemodialysis patients. *Am J Nephrol* 2007; 27:30–35.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3-105–3-108.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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Signed: Zaman Khan
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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G71326R04

B7P970

Read Highlighted Changes: Revised March 2018.

REF 07P9720

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NAME

Alinity c Direct Bilirubin Reagent Kit (also referred to as Bili D)

INTENDED USE

The Alinity c Direct Bilirubin assay is used for the quantitative analysis of direct bilirubin in human serum or plasma on the Alinity c analyzer.

SUMMARY AND EXPLANATION OF THE TEST

Red blood cells at the end of their circulating life are broken down in the reticuloendothelial system, mainly the spleen. The resulting heme, once the iron is removed, is then converted to bilirubin. This process accounts for about 80% of the 500 μmol (300 mg) of bilirubin formed daily. Other sources of bilirubin include the breakdown of myoglobin and cytochromes and the catabolism of immature red blood cells in the bone marrow.

Once formed, bilirubin is transported to the liver bound to albumin. This fraction of bilirubin is referred to as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid (mono- and diglucuronides) to form conjugated bilirubin by the enzyme uridyl diphosphate glucuronyl transferase. Conjugated bilirubin or direct bilirubin is excreted via the biliary system into the intestine, where it is metabolized by bacteria to a group of products known collectively as stercobilinogen. Elimination is almost complete and serum levels are normally negligible.

Direct bilirubin is the sum of the conjugated fractions. Direct bilirubin is elevated in conditions causing hepatic obstruction, hepatitis, cirrhosis, several inherited enzyme deficiencies, and inherited defects in canalicular excretion.

PRINCIPLES OF THE PROCEDURE

Bilirubin determination is generally based on the reaction of bilirubin with a diazotized sulfanilic acid, described by Ehrlich.¹ In this method, direct (conjugated fractions) bilirubin couples with a diazonium salt in the presence of sulfamic acid to form the colored compound azobilirubin. The increase in absorbance at 548 nm due to azobilirubin is proportional to the direct bilirubin concentration.

Methodology: Diazo Reaction

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity c Direct Bilirubin Reagent Kit 07P97

Volumes (mL) listed in the table below indicate the volume per cartridge.

REF	07P9720
Tests per cartridge	360
Number of cartridges per kit	4
Tests per kit	1440
R1	68.1 mL
R2	21.0 mL
R1	Active ingredient: Sulfamic acid (9.7 g/L).
R2	Active ingredients: 2, 4-dichloroaniline (< 0.1 g/L), Sodium nitrite (< 0.1 g/L), HCl (33.54 g/L).

Warnings and Precautions

- IVD
- For *In Vitro* Diagnostic Use
- Rx ONLY

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.²⁻⁵

The following warnings and precautions apply to: R1	
DANGER	Contains sulfamic acid.
H314	Causes severe skin burns and eye damage.
Prevention	
P260	Do not breathe mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower.
P310	Immediately call a POISON CENTER or doctor / physician.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to: R2	
	
DANGER	Contains hydrochloric acid.
H314	Causes severe skin burns and eye damage.
H332	Harmful if inhaled.
H290	May be corrosive to metals.
Prevention	
P260	Do not breathe mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P271	Use only outdoors or in a well-ventilated area.
P280	Wear protective gloves / protective clothing / eye protection.
P234	Keep only in original container.
Response	
P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P310	Immediately call a POISON CENTER or doctor / physician.
P390	Absorb spillage to prevent material damage.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position.
Onboard	System Temperature	28 days	
Opened	2 to 8°C	Until expiration date	Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity c Direct Bilirubin assay file must be installed on the Alinity c analyzer prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

Alternate Result Units

Edit assay parameter "Result Units" to select an alternate unit.

Conversion formula:

$$(\text{Concentration in Default result unit}) \times (\text{Conversion factor}) = (\text{Concentration in Alternate result unit})$$

Default Result Unit	Conversion Factor	Alternate Result Unit
mg/dL	17.1	μmol/L

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay.

Other specimen types, collection tube types, and anticoagulants have not been verified with this assay.

Specimen Type	Collection Vessel	Special Conditions
Serum	Serum tubes (with or without gel barrier)	
Plasma	Collection tubes Acceptable anticoagulants are: Lithium heparin (with or without gel barrier) Sodium heparin EDTA	The use of tubes containing sodium fluoride/potassium oxalate is not recommended due to the potential of hemolysis formation with this anticoagulant.

- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.
- Abbott Laboratories has not verified the assay performance characteristics with neonatal specimens.
NOTE: Abbott Laboratories recommends the use of sample interference indices in the semi-quantitative mode to assist in the determination of sample integrity for all specimens. Refer to the Sample Interference Indices (HIL) application sheets.

Specimen Conditions

- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- For accurate results, plasma specimens should be free of platelets and other particulate matter. Ensure centrifugation is adequate to remove platelets.
- Use serum or plasma specimens without visible hemolysis or lipemia. Refer to the SPECIFIC PERFORMANCE CHARACTERISTICS section of this package insert.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross-contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Specimen Storage

Specimens should be protected from bright light as bilirubin is photolabile.⁶ Direct bilirubin is stable in serum and plasma as follows:

Specimen Type	Temperature	Maximum Storage Time
Serum/Plasma	20 to 25°C	2 days ⁷
	2 to 8°C	7 days ^{7, 8}
	-20°C	3 months ⁹
	-80°C	3 months ⁹

Avoid multiple freeze/thaw cycles.

Guder et al. suggest storage of frozen specimens at -20°C for no longer than the time intervals cited above.⁷

Each laboratory may establish a range around -20°C and/or -80°C from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

Stored specimens must be inspected for particulates. If present, mix with a low speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PROCEDURE

Materials Provided

07P97 Alinity c Direct Bilirubin Reagent Kit

Materials Required but not Provided

- Alinity c Direct Bilirubin assay file
- 08P6101 Alinity c Bilirubin Calibrator Kit
- Commercially available controls containing direct bilirubin
- Saline (0.85% to 0.90% NaCl) for specimen dilution

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
 - Sample volume for single test: 5.0 µL.

NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the Alinity ci-series Operations Manual, Section 4.

- Refer to the Alinity c Bilirubin Calibrator Kit package insert and/or commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples with a direct bilirubin value exceeding 15.0 mg/dL (256.5 µmol/L) are flagged with the code "> 15.0 mg/dL" (> 256.5 µmol/L) and may be diluted with either the Automated Dilution Protocol or the Manual Dilution Procedure.

Automated Dilution Protocol

If using an automated dilution protocol, the system performs a dilution of the sample and automatically calculates the concentration by multiplying the result by the dilution factor. For details on configuring automated dilutions, refer to the Alinity ci-series Operations Manual, Section 2.

Manual Dilution Procedure

Dilute the sample with saline (0.85% to 0.90% NaCl).

The operator must enter the dilution factor in the Specimen or Control tab of the Create Order screen. The system will use this dilution factor to automatically calculate the concentration of the sample and report the result.

If the operator does not enter the dilution factor, the result must be manually multiplied by the appropriate dilution factor before reporting the result. If a diluted sample result is less than the lower value of the measuring interval of 0.1 mg/dL (1.7 µmol/L), do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration is stable for approximately 14 days (336 hours), but is required with each change in reagent lot. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Quality Control Procedures

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Two levels of controls (normal and abnormal) are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Commercial controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

Quality Control Guidance

Refer to "Basic QC Practices" by James O Westgard, Ph.D. for guidance on laboratory quality control practices.¹⁰

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

RESULTS

Calculation

The Alinity c Direct Bilirubin assay utilizes the Linear data reduction method to generate a calibration and results.

For information on alternate result units, refer to the INSTRUMENT PROCEDURE, Alternate Result Units section of this package insert.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in mg/dL (µmol/L) which meets the limits of acceptable performance for linearity, imprecision, and bias.

The measuring interval of the Alinity c Direct Bilirubin assay is 0.1 to 15.0 mg/dL (1.7 to 256.5 µmol/L).

LIMITATIONS OF THE PROCEDURE

Some specimens may give a direct bilirubin result slightly greater than the total bilirubin result. During internal testing at Abbott Laboratories, specimens with total bilirubin concentrations of 0.2 mg/dL (3.4 µmol/L) or less occasionally gave a direct bilirubin result that slightly exceeded their respective total bilirubin result. This may be observed when nearly all reacting bilirubin is direct bilirubin.

For patients undergoing evaluations involving the administration of indocyanine green (ICG), it is recommended that samples are drawn after ICG has been eliminated. See the Interference section for additional information.^{11, 12}

Abbott Laboratories has not verified the assay performance characteristics with neonatal specimens.

Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert.

EXPECTED VALUES

This study was performed on the ARCHITECT c System.

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics. Serum¹³

	Range (mg/dL)	Range (µmol/L)
Adult	0.0 to 0.5	0.0 to 8.6

A study was conducted using 136 serum samples from volunteers ranging in age from 25 to 66 years. Data were analyzed as described by NCCLS C28-A.¹⁴ From this study, 95% of all specimens fell within 0.0 to 0.5 mg/dL, with samples ranging from 0.1 to 0.6 mg/dL.

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

The Alinity c analyzer and the ARCHITECT c System and AEROSET System utilize the same reagents and sample/reagent ratios. Unless otherwise specified, all studies were performed on the Alinity c analyzer.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2.¹⁵ Testing was conducted using 1 lot of the Alinity c Direct Bilirubin Reagent Kit, 1 lot of the Alinity c Direct Bilirubin Calibrator Kit, and 1 lot of commercially available controls and 1 instrument. Three controls and 2 human serum panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.

Sample	n	Mean (mg/dL)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
			SD	%CV	SD	%CV
Control Level 1	120	0.4	0.02	3.9	0.02	3.9
Control Level 2	120	1.2	0.02	1.4	0.06	4.8
Control Level 3	120	2.3	0.02	0.8	0.09	4.0
Panel 1	120	3.5	0.04	1.3	0.06	1.7
Panel 2	120	8.6	0.08	0.9	0.11	1.3

^a Includes within-run, between-run, and between-day variability.

Sample	n	Mean (μmol/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
			SD	%CV	SD	%CV
Control Level 1	120	6.8	0.28	4.1	0.28	4.1
Control Level 2	120	21.4	0.29	1.4	1.03	4.8
Control Level 3	120	39.7	0.32	0.8	1.58	4.0
Panel 1	120	60.5	0.76	1.3	1.03	1.7
Panel 2	120	147.6	1.37	0.9	1.88	1.3

^a Includes within-run, between-run, and between-day variability.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.¹⁶ Testing was conducted using 3 lots of the Alinity c Direct Bilirubin Reagent Kit on each of 2 instruments over a minimum of 3 days. The maximum observed Limit of Blank (LoB), Limit of Detection (LoD), and Limit of Quantitation (LoQ) values are summarized below.

	mg/dL	μmol/L
LoB ^a	0.0	0.0
LoD ^b	0.1	1.7
LoQ ^c	0.1	1.7

^a The LoB represents the 95th percentile from $n \geq 60$ replicates of zero-analyte samples.

^b The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on $n \geq 60$ replicates of low-analyte level samples.

^c The LoQ was determined from $n \geq 60$ replicates of low-analyte level samples and is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met.

Linearity

A study was performed based on guidance from CLSI EP06-A.¹⁷ This assay is linear across the measuring interval of 0.1 to 15.0 mg/dL (1.7 to 256.5 μmol/L).

Interference

Potentially Interfering Substances

Potential interference in the Direct Bilirubin assay from 62 mg/dL (0.62 g/L) hemoglobin, 125 mg/dL (1.25 g/L) Intralipid, or 0.50 mmol/L Indican (indoxyl sulfate) is $\leq 10\%$ or ± 0.1 mg/dL, whichever is greater, at the medical decision level of the analyte. This study was performed on the ARCHITECT c System and the AEROSET System.

Potential interference in the Direct Bilirubin assay from 6.3 mg/L (8.1 μmol/L) indocyanine green is $\leq 10\%$ or ± 0.1 mg/dL, whichever is greater, at the medical decision level of the analyte.

This study was performed on the ARCHITECT c System.

These studies were performed based on guidance from NCCLS EP7-P.¹⁸

Interference effects were assessed by Dose Response and Paired Difference methods at the medical decision level of the analyte.

Potentially Interfering Substance	Interferent Level		Direct Bilirubin	
	Default Units	Alternate Units	Target Level (mg/dL)	Difference from Target (mg/dL)
Hemoglobin	31 mg/dL	0.31 g/L	0.4	-0.1
	62 mg/dL	0.62 g/L	0.4	-0.1
	125 mg/dL	1.25 g/L	0.4	-0.2
	250 mg/dL	2.50 g/L	0.4	-0.2
	500 mg/dL	5.00 g/L	0.4	-0.2
Human triglyceride	519 mg/dL	5.86 mmol/L	0.4	-0.1
	1034 mg/dL	11.68 mmol/L	0.4	0.3
Intralipid	125 mg/dL	1.25 g/L	0.4	-0.1
	250 mg/dL	2.50 g/L	0.4	0.1
	500 mg/dL	5.00 g/L	0.4	0.4
Indocyanine Green	6.3 mg/L	8.1 μmol/L	0.3	+0.1
	12.5 mg/L	16.1 μmol/L	0.3	+0.3
	18.8 mg/L	24.2 μmol/L	5.1	+0.4
	25.0 mg/L	32.3 μmol/L	5.1	+0.5

Taki et al. reported indoxyl sulfate concentrations up to 8.62 mg/dL (0.40 mmol/L), with an average of 3.52 mg/dL (0.17 mmol/L), in 224 hemodialysis (HD) patients.¹⁹ Indoxyl sulfate does not cause significant interference with this direct bilirubin method. Testing at Abbott Laboratories demonstrated that addition of 12.57 mg/dL (0.50 mmol/L) 3-indoxyl sulfate potassium salt to specimens increased the direct bilirubin concentration by a maximum of 0.1 mg/dL.

Indocyanine green solutions at the above concentrations were prepared by the individual addition of indocyanine green to two pools of plasma, one with a high concentration of bilirubin and one with a low concentration of bilirubin.

Interferences from medication or endogenous substances may affect results.²⁰

Method Comparison

A study was performed based on guidance from CLSI EP09-A3 using the Passing-Bablok regression method.²¹

	Units	n	Correlation Coefficient		Intercept Slope		Concentration Range
Alinity c Direct Bilirubin vs ARCHITECT Direct Bilirubin	Serum	mg/dL	132	1.00	0.00	1.00	0.2 to 14.7
		μmol/L	132	1.00	0.00	1.00	3.4 to 250.6

BIBLIOGRAPHY

1. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1136.
2. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
3. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
4. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
5. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
6. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1169.
7. Guder WG, Narayanan S, Wisser H, et al. List of Analytes—preanalytical variables. Annex In: *Samples: From the Patient to the Laboratory*. Darmstadt: GIT Verlag; 1996:Annex 8–9.
8. US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*. 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
9. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press; 1997:3–88.
10. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
11. Donnachie EM, Seccombe DW, Urquhart NI, et al. Indocyanine green interference in the Kodak Ektachem determination of total bilirubin. *Clin Chem* 1989;35(5):899–900.
12. Meijer DKF, Weert B, Vermeer GA. Pharmacokinetics of biliary excretion in man. VI. Indocyanine green. *Eur J Clin Pharmacol* 1988;35:295–303.
13. Data on file at Abbott Laboratories.
14. National Committee for Clinical Laboratory Standards (NCCLS). *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline*. NCCLS Document C28-A. Villanova, PA: NCCLS; 1995.
15. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
16. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
17. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
18. National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
19. Taki K, Tsuruta Y, Niwa T. Indoxyl sulfate and atherosclerotic risk factors in hemodialysis patients. *Am J Nephrol* 2007; 27:30–35.
20. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3–105–3–108.
21. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

ISO 15223 Symbols	
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
REF	List Number
SN	Serial number
Other Symbols	
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
PRODUCED FOR ABBOTT BY	Produced for Abbott by
PRODUCT OF CANADA	Product of Canada
R1	Reagent 1
R2	Reagent 2
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).

Alinity, ARCHITECT, and AEROSET are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown
Prince Edward Island
C1E 2B9 Canada

DISTRIBUTED IN THE USA BY
Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised March 2018.

©2016, 2018 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity c Direct Bilirubin Reagent Kit for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for devices are 2-8°C. Medical devices must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Medical devices must be carried out in single-use thermal containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in refrigerated trucks using thermal indicators.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. Manufacturer set the maximum shelf life of

Alinity c Direct Bilirubin Reagent Kit – 14 months from date of manufacturing.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste". Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, 1272/2008/EC (CLP), 1907/2006/EC (REACH), EN ISO 9001, ISO 13485, EN ISO 14971, EN 13641 EN ISO 15223-1, EN ISO 15223-2, EN 980 EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, ISO 80000-1, 1999/45/EC

POTENTIAL USER

The Alinity c Direct Bilirubin Reagent Kit is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the Alinity c analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRAINDICATIONS FOR USE

Not applicable.

TECHNICAL CHARACTERISTICS

- Alinity c Direct Bilirubin Reagent 1: pH-value at 20 °C: 1-1.4
Density at 20 °C: 1.012 g/cm³
- Alinity c Direct Bilirubin Reagent 2: pH-value at 20 °C: ≤0.5
Density at 20 °C: 1.019 g/cm³

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

COLOUR AND ODOUR

Alinity c Direct Bilirubin Reagent 1: Odourless and Colourless

Alinity c Direct Bilirubin Reagent 2: Odourless and Colourless

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- Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:136.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1169.
- Guder WG, Narayanan S, Wisser H, et al. List of Analytes—preanalytical variables. Annex In: *Samples: From the Patient to the Laboratory*. Darmstadt: GIT Verlag; 1996:Annex 8–9.
- US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*. 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
- Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press; 1997:3-88.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Donnachie EM, Seccombe DW, Urquhart NI, et al. Indocyanine green interference in the Kodak Ektachem determination of total bilirubin. *Clin Chem* 1989;35(5):899-900.
- Meijer DKF, Weert B, Vermeer GA. Pharmacokinetics of biliary excretion in man. VI. Indocyanine green. *Eur J Clin Pharmacol* 1988;35:295-303.
- Data on file at Abbott Laboratories.
- National Committee for Clinical Laboratory Standards (NCCLS). *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline*. NCCLS Document C28-A. Villanova, PA: NCCLS; 1995.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
- Taki K, Tsuruta Y, Niwa T. Indoxyl sulfate and atherosclerotic risk factors in hemodialysis patients. *Am J Nephrol* 2007; 27:30–35.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3-105–3-108.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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