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
Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

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I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG., hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: Zaman Khan
Zaman Khan
Associate Director Regulatory
Affairs EMEA and Pakistan
Abbott GmbH & Co. KG

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2018
N. [Signature]
[Circular Stamp]

Key to Symbols

	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Calibrator Low
	Calibrator High
	Control Number
	Concentration
	In Vitro Diagnostic Medical Device
	Lot Number
	Produced for Abbott by
	Product of USA
	List Number
	For use by or on the order of a physician only (applicable to USA classification only).
	Serial Number

Alinity is a trademark of Abbott Laboratories in various jurisdictions.
All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Fisher Diagnostics,
A Div. of Fisher Scientific Company, LLC
A Part of Thermo Fisher Scientific, Inc.
8365 Valley Pike Middletown VA 22645 USA

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Nt

ADDENDUM

to the Instruction for Use Ainity c ICT Serum Calibrator Kit (for Russia)

SHIPPING

Recommended shipping conditions for device is 2-8°C. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

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MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MO are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

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Signed: 
Zaman Khan
Associate Director Regulatory
Affairs EMEA and Pakistan
Abbott GmbH & Co. KG

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S8P700

Revised November 2016.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

Alinity c ICT Urine Calibrator Kit (also referred to as ICT Urine)

INTENDED USE

For use in the calibration of the urine Sodium, Potassium, and Chloride assays.

For additional information, refer to the Alinity c ICT Sample Diluent package insert and the Alinity ci-series Operations Manual.

CONTENTS

The Alinity c ICT Urine Calibrator is prepared in an aqueous solution.

Calibrator	Quantity	CONC		
		Sodium	Potassium	Chloride
URINE CAL L	5 x 2.9 mL	50 mmol/L	9.0 mmol/L	50 mmol/L
URINE CAL H	5 x 2.9 mL	180 mmol/L	90.0 mmol/L	180 mmol/L

MATERIALS REQUIRED BUT NOT PROVIDED

- 04R1001 Alinity ci-series Calibrator/Control Replacement Caps

STANDARDIZATION

Using gravimetric methods, each level of the aqueous-based calibrator is prepared with ACS grade sodium chloride, potassium chloride, and sodium nitrate. The concentration of sodium, potassium, and chloride is determined using a reference method (flame photometry, titration using silver) calibrated against NIST SRM 918 and NIST SRM 919.

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PREPARATION FOR USE

- This product is liquid ready-to-use.
- Prior to each use, mix by gentle inversion.

STORAGE

- Protect from strong light for optimum stability.
- Do not use past expiration date.
- This product may be stored on the instrument.
- If stored onboard, mixing is not required for each additional use.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	15 to 30°C	Until expiration date	
Onboard	System Temperature	6 days (144 hours)	
Opened	2 to 8°C	7 days after opening	Store tightly capped with new replacement cap. Return to refrigerated storage after use.

The analyzer will track In-use Stability, which is the time the calibrator is outside of refrigerated storage while on the analyzer. The analyzer will not allow the use of the calibrator if the In-use Stability has been exceeded. Maximum In-use Stability **CAN** be found in the Assay Parameter Report. For additional information on calibrator In-use Stability, refer to the Alinity ci-series Operations Manual, Section 5.

For additional information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

INSTRUMENT PROCEDURE

- Calibrator lots may be configured using the bar code label on the calibrator carton.
- Calibration is performed by running the Alinity c ICT Urine Calibrator set.
- Verify control results are within acceptable limits before reporting patient results.
- For information on configuring calibrator data, refer to the Alinity ci-series Operations Manual, Section 2.
- For instructions on ordering and loading calibrators on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are visible signs of leakage, if calibration does not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria, or if controls do not meet the appropriate criteria.

LIMITATIONS OF THE PROCEDURE

Accurate and reproducible results are dependent upon properly functioning instruments and reagents, storage of product as directed, and good laboratory technique.

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Abbott

Key to Symbols

	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Control Number
	Concentration
	<i>In Vitro</i> Diagnostic Medical Device
	Lot Number
	Produced for Abbott by
	Product of Canada
	List Number
	For use by or on the order of a physician only (applicable to USA classification only).
	Serial Number
	Urine Calibrator Low
	Urine Calibrator High

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

 Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

 **PRODUCED FOR ABBOTT BY**
Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown,
Prince Edward Island
C1E 2B9 Canada

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com

Revised November 2016.
©2016 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity c ICT Urine Calibrator Kit (for Russia)

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WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 13 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

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ICT Diluent
07P53
G71456R01
B7P530

REF 07P5320

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NAME

Alinity c ICT Sample Diluent (also referred to as ICT Diluent).

INTENDED USE

The Alinity c ICT (Integrated Chip Technology) is used for the quantitation of sodium, potassium, and chloride in human serum, plasma, or urine on the Alinity c analyzer.

SUMMARY AND EXPLANATION OF THE TEST

Sodium is the major cation of extracellular fluid; it plays an essential role in the normal distribution of water and in the maintenance of osmotic pressure in extracellular fluid compartments. Decreased levels of sodium may be caused by an excessive use of diuretics, prolonged vomiting, a decrease in the intake of sodium in the diet, and metabolic acidosis. Increased levels of sodium may be found in Cushing's syndrome, severe dehydration, or in high levels of salt intake without an adequate supply of water.^{1,2}

Potassium is the major intracellular cation. The concentration of potassium in the erythrocytes is approximately 23 times the concentration in plasma. For this reason, only unhemolyzed samples must be used. Decreased levels of extracellular potassium are characterized by weakness in the muscles, irritability, paralysis, accelerated heartbeat, and eventually cardiac arrest, and may be caused by a poor intake of potassium in the diet, by a redistribution of extracellular potassium, and by an increased loss of body fluids rich in potassium. Abnormally elevated levels of extracellular potassium produce mental confusion, general weakness, numbness, flaccid paralysis in the extremities, a slowed heart rate, and eventually collapse of the peripheral vascular system and cardiac arrest. Causes of increased potassium levels may be linked to inappropriate intravenous therapy, dehydration, shock, diabetic ketoacidosis, and severe burns.^{1,2}

Chloride is the major extracellular anion. The majority of ingested chloride is absorbed, and the excess is excreted along with other ions into the urine. Low levels of chloride are observed in the case of prolonged vomiting accompanied by the loss of hydrochloric acid (HCl), in metabolic alkalosis, in critical cases of Addison's disease, and in kidney disease resulting in loss of salt. Elevated levels of chloride are observed in metabolic acidosis associated with prolonged diarrhea and with loss of sodium bicarbonate (NaHCO_3), and in the case of renal tubular diseases in which there is a decreased excretion of hydrogen ion (H^+), which causes in turn a decrease in the reabsorption of bicarbonate ion (HCO_3^-).^{1,2}

PRINCIPLES OF THE PROCEDURE

Ion-selective electrodes (ISE) for sodium, potassium, and chloride utilize membranes selective to each of these ions. An electrical potential (voltage) is developed across the membranes between the reference and measuring electrodes in accordance with the Nernst equation. The voltage is compared to previously determined calibrator voltages and converted into ion concentration.

Methodology: Ion-selective electrode diluted (Indirect)

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity c ICT Sample Diluent 07P53

Volumes (mL) listed in the table below indicate the volume per cartridge.

REF	07P5320
Tests per cartridge	935
Number of cartridges per kit	10
Tests per kit	9350*
R1	68.2 mL
R1	Active Ingredient: buffer.

*One test per sample can generate one to three results (Na^+ , K^+ , and Cl^-).

Warnings and Precautions

- **IVD**
- For In Vitro Diagnostic Use
- **Rx ONLY**

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.³⁻⁶

The following warnings and precautions apply to: ICT Sample Diluent*	
H402	Harmful to aquatic life.
Prevention	
P273	Avoid release to the environment.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

* Not applicable where regulation EU 1272/2008 (CLP) has been implemented.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Abbott

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	15 to 30°C	Until expiration date	Store in upright position.
Onboard	System Temperature	30 days	
Opened	15 to 30°C	Until expiration date	Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 15 to 30°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

Acetic acid odor from the diluent is normal.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity c Sodium, Potassium, and Chloride assay files must be installed on the Alinity c analyzer prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay. Other specimen types, collection tube types, and anticoagulants have not been verified with this assay.

Specimen Type	Collection Vessel	Special Conditions
Serum	Serum tubes (with or without gel barrier)	For Potassium, hemolyzed specimens must not be used.
Plasma	Collection tubes Acceptable anticoagulants are: Lithium heparin (with or without gel barrier) Sodium heparin (full draw)	
Urine (random; 24-hour)	Without preservatives ¹	

- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- For accurate results, plasma specimens should be free of platelets and other particulate matter. Ensure centrifugation is adequate to remove platelets.
- NOTE: Multiple myeloma and lipid samples are known to give low results on diluted ISE systems due to the high level of proteins/lipids present in the sample.^{7,8}
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross-contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter.
- NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Specimen Storage

Serum/Plasma^{9,10}

Temperature	Maximum Storage Time		
	Sodium	Potassium	Chloride
20 to 25°C	2 weeks	1 week	7 days
2 to 8°C	2 weeks	1 week	7 days
-20°C	1 year	1 year	> 1 year

Urine^{9,11}

Temperature	Maximum Storage Time		
	Sodium	Potassium	Chloride
20 to 25°C	45 days	45 days	7 days
2 to 8°C	45 days	2 months	7 days
-20°C	1 year	1 year	7 days

Avoid multiple freeze/thaw cycles.

Guder et al. suggest storage of frozen specimens at -20°C for no longer than the time intervals cited above.⁹

Each laboratory may establish a range around -20°C from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

Stored specimens must be inspected for particulates. If present, mix with a low speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PRC

Material
07P53 /

Material

- Alin
- 09C
- 08F
- 08F
- 08F
- Col
- an

For info
refer to
For info
refer to

Assay

For a
ci-seri

- II
- O
- pr
- T
- ci
- M

- F
- A
- a
- L
- F
- C
- F
- r
- !
- !

Cali

For i
serie

Calil
requ

any
calit

or p
acci

your
rece

This
part

perf

Qu

As
pro

con

PROCEDURE

Materials Provided

07P53 Alinity c ICT Sample Diluent

Materials Required but not Provided

- Alinity c Sodium, Potassium, and Chloride assay files
- 09D2804 Alinity ICT Module
- 08P7640 Alinity c-series ICT Reference Solution
- 08P6901 Alinity c ICT Serum Calibrator Kit
- 08P7001 Alinity c ICT Urine Calibrator Kit
- Commercially available controls containing sodium, potassium, and chloride

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
 - Sample volume for single test: 15 µL (serum/plasma/urine).
NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the Alinity ci-series Operations Manual, Section 4.
- Refer to the Alinity c ICT Serum Calibrator Kit package insert, Alinity c ICT Urine Calibrator Kit package insert and commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration is stable for approximately 1 day (24 hours), but is required with each change in diluent lot. The laboratory may choose any calibration interval up to 24 hours. The use of a particular calibration time interval is dependent on individual laboratory policy or preference. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Quality Control Procedures

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Two levels of controls (normal and abnormal) are to be run every 8 hours and following calibration.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.

- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of ICT Sample Diluent, ICT Reference Solution, or calibrator lot.

Commercial controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

Quality Control Guidance

Refer to "Basic QC Practices" by James O Westgard, Ph.D. for guidance on laboratory quality control practices.¹²

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

RESULTS

Calculation

The Alinity c Sodium, Potassium, and Chloride assays utilize the Potentiometric data reduction method to generate a calibration and results.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in mmol/L which meets the limits of acceptable performance for linearity, imprecision, and bias.

Serum/Plasma

The measuring interval of the Alinity c Sodium assay is 100 to 200 mmol/L.

The measuring interval of the Alinity c Potassium assay is 1.0 to 10.0 mmol/L.

The measuring interval of the Alinity c Chloride assay is 50 to 150 mmol/L.

Urine

The measuring interval of the Alinity c Sodium assay is 20 to 400 mmol/L.

The measuring interval of the Alinity c Potassium assay is 1.0 to 300.0 mmol/L.

The measuring interval of the Alinity c Chloride assay is 20 to 300 mmol/L.

LIMITATIONS OF THE PROCEDURE

Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert.

The Alinity c Sodium assay serum application is susceptible to interference effects from conjugated bilirubin at > 30 mg/dL, hemoglobin at > 500 mg/dL, sodium salicylate at > 35 mg/dL, and benzalkonium chloride at > 1 mg/dL.

The Alinity c Potassium assay serum application is susceptible to interference effects from hemoglobin (hemolysis) at > 100 mg/dL and benzalkonium chloride at > 5 mg/dL.

The Alinity c Chloride assay serum application is susceptible to interference effects from hemoglobin (hemolysis) at > 1000 mg/dL, lithium bromide at > 20 mg/dL, and lithium iodide at > 25.4 mg/dL.

The Alinity c Sodium assay urine application demonstrated greater than 10% interference with sodium carbonate at 1.25 g/dL, sodium fluoride at 400 mg/dL, and sodium oxalate at 60 mg/dL for samples with sodium concentrations ≤ 74 mmol/L.

The Alinity c Chloride assay urine application demonstrated greater than 10% interference from 6N hydrochloric acid at 2.5 mL/dL.

EXPECTED VALUES

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

Reference Range

Serum/Plasma¹³

		Range (mmol/L)
Sodium	Premature, Cord	116 to 140
	Premature, 48 hours	128 to 148
	Newborn, Cord	126 to 166
	Newborn	133 to 146
	Infant	139 to 146
	Child	138 to 145
	Adult	136 to 145
	> 90 years	132 to 146
Potassium	Premature, Cord	5.0 to 10.2
	Premature, 48 hours	3.0 to 6.0
	Newborn, Cord	5.6 to 12.0
	Newborn	3.7 to 5.9
	Infant	4.1 to 5.3
	Child	3.4 to 4.7
	Adult	3.5 to 5.1
	Plasma, Male	3.5 to 4.5
Chloride	Plasma, Female	3.4 to 4.4
	Cord	96 to 104
	Premature	95 to 110
	0 to 30 days	98 to 113
	Adult	98 to 107
	> 90 years	98 to 111

Urine¹³

			Range (mmol/day)
Sodium	6 to 10 years	Male	41 to 115
		Female	20 to 69
	10 to 14 years	Male	63 to 177
		Female	48 to 168
	Adult	Male	40 to 220
		Female	27 to 287
Potassium	6 to 10 years	Male	17 to 54
		Female	8 to 37
	10 to 14 years	Male	22 to 57
		Female	18 to 58
	Adult		25 to 125
Chloride	Infant		2 to 10
	Child < 6 years		15 to 40
	6 to 10 years	Male	36 to 110
		Female	18 to 74
	10 to 14 years	Male	64 to 176
		Female	36 to 173
	Adult		110 to 250
	> 60 years		95 to 195

For Sodium, Potassium, and Chloride, results expressed in mmol/L are equivalent to mEq/L; results expressed in mmol/day are equivalent to mEq/day.

24-Hour Urinary Excretion

Sodium, Potassium, and Chloride

To convert results from mmol/L to mmol/day (24-hour urinary excretion)

$$24\text{-hour excretion} = [(V \times c) + 1000] \text{ mmol/day}$$

Where:

V = 24-hour urine volume (mL)

c = analyte concentration (mmol/L)

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

All studies were performed on the Alinity c analyzer.

Precision

Within-Laboratory Precision

Serum

A study was performed based on guidance from CLSI EP05-A2.

Testing was conducted using 2 lots of the Alinity c ICT Sample Diluent, 2 lots of the Alinity c ICT Serum Calibrator Kit, 2 lots of commercially available controls and 2 instruments. Three controls and 2 human serum panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.¹⁴

Sodium

Sample	Control Lot	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
				SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	1	249	125	0.5	0.4	0.9 (0.8 - 1.0)	0.7 (0.6 - 0.8)
	2	249	125	0.5	0.4	0.9 (0.8 - 0.9)	0.7 (0.6 - 0.8)
Control Level 2	1	252	144	0.5	0.4	0.9 (0.8 - 1.0)	0.6 (0.6 - 0.7)
	2	249	144	0.5	0.4	0.9 (0.8 - 1.0)	0.6 (0.5 - 0.7)
Control Level 3	1	249	161	0.5	0.3	1.0 (0.9 - 1.1)	0.6 (0.6 - 0.7)
	2	249	161	0.6	0.4	1.1 (0.9 - 1.2)	0.7 (0.6 - 0.7)
Panel A	N/A	496	112	0.4	0.4	0.8 (0.6 - 0.9)	0.7 (0.6 - 0.8)
Panel B	N/A	498	190	0.8	0.4	1.4 (1.2 - 1.6)	0.7 (0.6 - 0.9)

Potassium

Sample	Control Lot	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
				SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	1	246	2.8	0.03	0.9	0.04 (0.04 - 0.04)	1.4 (1.3 - 1.5)
	2	246	2.8	0.01	0.3	0.01 (0.00 - 0.02)	0.4 (0.0 - 0.6)
Control Level 2	1	246	4.0	0.02	0.6	0.03 (0.03 - 0.03)	0.8 (0.8 - 0.9)
	2	246	4.0	0.01	0.3	0.02 (0.02 - 0.02)	0.5 (0.5 - 0.5)
Control Level 3	1	246	6.8	0.03	0.4	0.05 (0.04 - 0.05)	0.7 (0.6 - 0.7)
	2	246	6.8	0.03	0.5	0.05 (0.04 - 0.05)	0.7 (0.6 - 0.7)
Panel A	N/A	491	1.6	0.02	1.1	0.03 (0.00 - 0.04)	1.7 (0.0 - 2.7)
Panel B	N/A	492	9.4	0.04	0.4	0.06 (0.05 - 0.07)	0.7 (0.5 - 0.8)

Chloride

Sample
Control
Level 1

Control
Level 2

Control
Level 3

Panel A

Panel B

^a Inclu

^b Maxi

instrum

Urine

A stud

Testin

Diluen

comm

and u

replica

Sodiu

Sampl

Control

Level

Control

Level

Panel

Panel

Pota

Sampl

Control

Level

Control

Level

Panel

Panel

Panel

Chloride

Sample	Control Lot	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
				SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	1	243	85	0.5	0.6	0.7 (0.6 - 0.8)	0.8 (0.7 - 0.9)
	2	243	84	0.4	0.5	0.6 (0.5 - 0.7)	0.7 (0.5 - 0.9)
Control Level 2	1	243	95	0.4	0.4	0.6 (0.6 - 0.7)	0.7 (0.6 - 0.8)
	2	243	94	0.4	0.5	0.6 (0.5 - 0.7)	0.6 (0.5 - 0.7)
Control Level 3	1	243	110	0.5	0.5	0.8 (0.6 - 0.9)	0.7 (0.5 - 0.9)
	2	243	109	0.5	0.5	0.7 (0.5 - 0.9)	0.7 (0.5 - 0.8)
Panel A	N/A	486	55	0.3	0.6	0.6 (0.4 - 0.6)	1.0 (0.8 - 1.2)
Panel B	N/A	485	132	0.6	0.4	1.0 (0.7 - 1.3)	0.8 (0.5 - 0.9)

^a Includes within-run, between-run, and between-day variability.

^b Maximum and minimum SD or %CV for each reagent lot and instrument combination.

Urine

A study was performed based on guidance from CLSI EP05-A2. Testing was conducted using 2 lots of the Afinity c ICT Sample Diluent, 2 lots of the Afinity c ICT Urine Calibrator Kit, 2 lots of commercially available controls and 2 instruments. Two controls and up to 3 human urine panels were assayed in a minimum of 2 replicates at 2 separate times per day for a minimum of 20 days.¹⁴

Sodium

Sample	Control Lot	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
				SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	1	240	92	0.6	0.7	1.0 (0.9 - 1.1)	1.1 (1.0 - 1.2)
	2	240	92	0.5	0.5	0.7 (0.7 - 0.8)	0.8 (0.7 - 0.9)
Control Level 2	1	240	161	0.7	0.4	1.0 (0.9 - 1.1)	0.6 (0.6 - 0.7)
	2	240	161	0.6	0.4	1.5 (1.0 - 1.8)	0.9 (0.7 - 1.1)
Panel A	N/A	480	21	0.5	2.3	0.6 (0.5 - 0.8)	2.9 (2.3 - 3.7)
Panel B	N/A	480	383	1.5	0.4	3.9 (2.5 - 4.7)	1.0 (0.7 - 1.2)

Potassium

Sample	Control Lot	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
				SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	1	240	16.6	0.06	0.3	0.13 (0.13 - 0.14)	0.8 (0.8 - 0.9)
	2	240	16.6	0.05	0.3	0.12 (0.11 - 0.12)	0.7 (0.7 - 0.7)
Control Level 2	1	240	58.1	0.17	0.3	0.34 (0.29 - 0.39)	0.6 (0.5 - 0.7)
	2	240	58.2	0.17	0.3	0.49 (0.29 - 0.62)	0.8 (0.5 - 1.1)
Panel A	N/A	480	1.7	0.02	1.4	0.04 (0.04 - 0.04)	2.4 (2.3 - 2.6)
Panel B	N/A	480	127.7	0.35	0.3	0.68 (0.50 - 0.83)	0.5 (0.4 - 0.6)
Panel C	N/A	479	284.5	0.82	0.3	1.91 (1.65 - 2.12)	0.7 (0.6 - 0.7)

Chloride

Sample	Control Lot	n	Mean (mmol/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
				SD	%CV	SD (Range ^b)	%CV (Range ^b)
Control Level 1	1	240	103	0.5	0.4	0.9 (0.8 - 1.0)	0.9 (0.8 - 0.9)
	2	240	103	0.5	0.5	0.8 (0.8 - 0.8)	0.8 (0.8 - 0.8)
Control Level 2	1	240	193	0.7	0.4	1.2 (1.1 - 1.3)	0.6 (0.6 - 0.7)
	2	240	193	0.7	0.4	1.9 (1.3 - 2.4)	1.0 (0.7 - 1.2)
Panel A	N/A	479	24	0.3	1.1	0.4 (0.1 - 0.5)	1.6 (0.4 - 2.2)
Panel B	N/A	480	273	0.9	0.3	2.2 (1.9 - 2.3)	0.8 (0.7 - 0.8)

^a Includes within-run, between-run, and between-day variability.

^b Maximum and minimum SD or %CV for each reagent lot and instrument combination.

Linearity

A study was performed based on guidance from CLSI EP06-A.¹⁵

Serum

The sodium assay is linear across the measuring interval of 100 to 200 mmol/L.

The potassium assay is linear across the measuring interval of 1.0 to 10.0 mmol/L.

The chloride assay is linear across the measuring interval of 50 to 150 mmol/L.

Urine

The sodium assay is linear across the measuring interval of 20 to 400 mmol/L.

The potassium assay is linear across the measuring interval of 1.0 to 300.0 mmol/L.

The chloride assay is linear across the measuring interval of 20 to 300 mmol/L.

Interference

Potentially Interfering Substances

A study was performed based on guidance from CLSI EP07-A2.¹⁶

Serum

Sodium

For sodium serum a bias of > 2% was considered significant interference.

Potentially Interfering Substance	Interferent Level		Sodium	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Unconjugated Bilirubin	30 mg/dL	513 µmol/L	115	100.4
			144	100.7
	60 mg/dL	1026 µmol/L	115	100.4
			145	100.0
Conjugated Bilirubin	20 mg/dL	287 µmol/L	112	100.9
			141	100.7
	30 mg/dL	356 µmol/L	113	100.9
			141	101.4
Hemoglobin	60 mg/dL	712 µmol/L	112	100.6
			141	102.5
	500 mg/dL	5.0 g/L	128	101.6
			157	101.0
	1000 mg/dL	10.0 g/L	125	101.6
			151	102.6
	2000 mg/dL	20.0 g/L	116	103.9
			141	104.6

Potentially Interfering Substance	Interferent Level		Sodium	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Intralipid	1000 mg/dL	10.0 g/L	111	100.9
			140	100.0
			106	100.9
Ascorbic Acid	6 mg/dL	341 µmol/L	128	100.0
			154	99.4
Acetaminophen	20 mg/dL	1323 µmol/L	128	100.0
			152	100.0
Ibuprofen	50 mg/dL	2424 µmol/L	125	100.0
			150	100.0
Acetylcysteine	167 mg/dL	10.2 mmol/L	127	100.0
			153	100.0
Acetylsalicylic Acid	66 mg/dL	3.7 mmol/L	127	100.0
			153	100.0
Sodium Salicylate	35 mg/dL	2.2 mmol/L	127	102.4
			151	101.3
			128	103.5
Benzalkonium Chloride	1 mg/dL	0.01 g/L (0.001%)	116	100.9
			146	100.0
			116	102.6
	5 mg/dL	0.05 g/L (0.005%)	116	102.6
			145	102.1
			115	106.1
	10 mg/dL	0.1 g/L (0.01%)	115	106.1
			143	106.3

Potassium

For potassium serum a bias of > 10% was considered significant interference.

Potentially Interfering Substance	Interferent Level		Potassium	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Unconjugated Bilirubin	30 mg/dL	513 µmol/L	3.0	100.0
			4.9	100.0
			2.9	103.4
	60 mg/dL	1026 µmol/L	2.9	103.4
			4.8	100.0
			4.8	100.0
Conjugated Bilirubin	20 mg/dL	237 µmol/L	2.9	100.0
			4.8	100.0
			3.0	101.7
	30 mg/dL	356 µmol/L	3.0	101.7
			4.9	98.0
			2.9	100.0
	60 mg/dL	712 µmol/L	2.9	100.0
			4.8	100.0
			4.8	100.0
Hemoglobin	100 mg/dL	1.00 g/L	3.3	109.1
			5.2	105.8
			3.3	112.1
	125 mg/dL	1.25 g/L	3.3	112.1
			5.2	107.7
			3.3	121.2
	250 mg/dL	2.50 g/L	3.3	121.2
			5.1	115.7
			5.1	115.7
Intralipid	1000 mg/dL	10.0 g/L	2.9	100.0
			4.7	102.1
			2.8	100.0
	2000 mg/dL	20.0 g/L	2.8	100.0
			4.5	102.2
			4.5	102.2
Ascorbic Acid	6 mg/dL	341 µmol/L	3.2	100.0
			5.1	100.0
Acetaminophen	20 mg/dL	1323 µmol/L	3.2	100.0
			5.1	100.0
Ibuprofen	50 mg/dL	2424 µmol/L	3.2	100.0
			5.1	100.0
Acetylcysteine	167 mg/dL	10.2 mmol/L	3.2	100.0
			5.1	100.0
Acetylsalicylic Acid	66 mg/dL	3.7 mmol/L	3.2	100.0
			5.1	100.0

Potentially Interfering Substance	Interferent Level		Potassium	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Sodium Salicylate	70 mg/dL	4.4 mmol/L	3.2	100.0
			5.1	100.0
Benzalkonium Chloride	1 mg/dL	0.01 g/L (0.001%)	3.0	103.3
			4.9	102.0
			3.0	110.0
	5 mg/dL	0.05 g/L (0.005%)	3.0	110.0
			4.9	108.2
			2.9	127.6
	10 mg/dL	0.1 g/L (0.01%)	2.9	127.6
			4.8	118.8

Chloride

For chloride serum a bias of > 10% was considered significant interference.

Potentially Interfering Substance	Interferent Level		Chloride	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Unconjugated Bilirubin	30 mg/dL	513 µmol/L	79	100.0
			114	100.0
			79	100.0
	60 mg/dL	1026 µmol/L	79	100.0
			113	100.9
			113	100.9
Conjugated Bilirubin	20 mg/dL	237 µmol/L	79	100.0
			113	100.0
			79	100.0
	30 mg/dL	356 µmol/L	79	100.0
			114	99.1
			79	99.4
	60 mg/dL	712 µmol/L	79	99.4
			113	100.0
			113	100.0
Hemoglobin	1000 mg/dL	10.0 g/L	90	106.7
			116	106.9
			84	113.1
	2000 mg/dL	20.0 g/L	84	113.1
			110	110.0
			110	110.0
Intralipid	1000 mg/dL	10.0 g/L	78	101.3
			112	100.0
	2000 mg/dL	20.0 g/L	75	100.0
			107	100.9
Ascorbic Acid	6 mg/dL	341 µmol/L	85	100.0
			112	100.9
Acetaminophen	20 mg/dL	1323 µmol/L	84	100.6
			112	99.1
Ibuprofen	50 mg/dL	2424 µmol/L	78	100.0
			106	100.0
Acetylcysteine	16.7 mg/dL	1.0 mmol/L	86	100.0
			113	100.9
Acetylsalicylic Acid	66 mg/dL	3.7 mmol/L	84	101.2
			113	100.0
Sodium Salicylate	70 mg/dL	4.4 mmol/L	85	98.8
			113	100.0
Benzalkonium Chloride	1 mg/dL	0.01 g/L (0.001%)	81	100.0
			112	100.0
	5 mg/dL	0.05 g/L (0.005%)	81	100.0
			112	100.0
	10 mg/dL	0.10 g/L (0.010%)	81	100.0
			111	100.3
Lithium Bromide	20 mg/dL	2.3 mmol/L	80	106.9
			111	106.8
	40 mg/dL	4.6 mmol/L	79	117.1
			109	116.1
Lithium Iodide	25.4 mg/dL	1.9 mmol/L	81	106.8
			111	107.2
	50.8 mg/dL	3.8 mmol/L	79	112.7
			109	111.9
Sodium Azide	163 mg/dL	25 mmol/L (0.163%)	80	105.0
			110	103.6
	325 mg/dL	50 mmol/L (0.325%)	79	107.6
			109	105.5

Urine

Sodium

For sodium urine a bias of > 10% was considered significant interference.

Potentially Interfering Substance	Interferent Level		Sodium	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Protein	50 mg/dL	0.5 g/L	33	101.5
			204	100.2
Glucose	1000 mg/dL	55.5 mmol/L	33	100.0
			208	99.8
Ascorbate	200 mg/dL	11.4 mmol/L	33	100.0
			208	100.0
8.5N Acetic Acid	6.25 mL/dL	531 mmol/L	31	100.0
			195	100.5
Boric Acid	250 mg/dL	40.4 mmol/L	33	100.0
			208	99.8
6N Hydrochloric Acid	2.5 mL/dL	150 mmol/L	33	100.0
			205	100.7
6N Nitric Acid	5.0 mL/dL	300 mmol/L	31	100.0
			195	100.8
Sodium Oxalate	60 mg/dL	4.5 mmol/L	33	127.3
			71	112.7
			154	106.5
			206	103.6
Sodium Carbonate	1.25 g/dL	117.9 mmol/L	33	798.5
			206	216.5
Sodium Fluoride	400 mg/dL	95.3 mmol/L	33	381.8
			205	146.2
Acetaminophen	20 mg/dL	1323 µmol/L	39	100.0
			211	100.0
Ibuprofen	50 mg/dL	2424 µmol/L	38	100.0
			207	99.8
Acetylcysteine	167 mg/dL	10.2 mmol/L	39	100.0
			210	100.0

Potassium

For potassium urine a bias of > 10% was considered significant interference.

Potentially Interfering Substance	Interferent Level		Potassium	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Protein	50 mg/dL	0.5 g/L	23.7	100.4
			117.9	100.2
Glucose	1000 mg/dL	55.5 mmol/L	23.8	100.0
			118.5	100.0
Ascorbate	200 mg/dL	11.4 mmol/L	23.8	100.0
			118.7	100.0
8.5N Acetic Acid	6.25 mL/dL	531 mmol/L	22.5	100.9
			112.2	100.1
Boric Acid	250 mg/dL	40.4 mmol/L	23.8	100.0
			118.4	100.0
6N Hydrochloric Acid	2.5 mL/dL	150 mmol/L	23.7	100.8
			118.3	100.6
6N Nitric Acid	5.0 mL/dL	300 mmol/L	22.5	100.9
			112.6	100.2
Sodium Oxalate	60 mg/dL	4.5 mmol/L	23.7	100.0
			118.9	100.0
Sodium Carbonate	1.25 g/dL	117.9 mmol/L	23.8	100.0
			119.2	99.5
Sodium Fluoride	400 mg/dL	95.3 mmol/L	23.9	99.6
			118.8	99.9
Acetaminophen	20 mg/dL	1323 µmol/L	24.4	100.0
			132.0	96.1
Ibuprofen	50 mg/dL	2424 µmol/L	24.3	100.0
			126.7	100.0

Potentially Interfering Substance	Interferent Level		Potassium	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Acetylcysteine	167 mg/dL	10.2 mmol/L	24.3	100.0
			126.6	100.1

Chloride

For chloride urine a bias of > 10% was considered significant interference.

Potentially Interfering Substance	Interferent Level		Chloride	
	Default Units	Alternate Units	Target Level (mmol/L)	Recovery (% of Target)
Protein	50 mg/dL	0.5 g/L	102	101.0
			234	100.2
Glucose	1000 mg/dL	55.5 mmol/L	103	100.0
			237	99.8
Ascorbate	200 mg/dL	11.4 mmol/L	103	101.0
			236	100.0
8.5N Acetic Acid	6.25 mL/dL	531 mmol/L	98	99.0
			222	100.0
Boric Acid	250 mg/dL	40.4 mmol/L	103	100.0
			235	100.0
6N Hydrochloric Acid	2.5 mL/dL	150 mmol/L	103	250.5
			236	165.3
6N Nitric Acid	5.0 mL/dL	300 mmol/L	98	96.9
			223	98.2
Sodium Oxalate	60 mg/dL	4.5 mmol/L	103	101.0
			236	100.0
Sodium Carbonate	1.25 g/dL	117.9 mmol/L	104	101.9
			238	101.3
Sodium Fluoride	400 mg/dL	95.3 mmol/L	103	96.1
			236	100.4
Acetaminophen	20 mg/dL	1323 µmol/L	114	100.0
			258	100.0
Ibuprofen	50 mg/dL	2424 µmol/L	107	100.9
			252	100.0
Acetylcysteine	16.7 mg/dL	1.0 mmol/L	115	100.9
			261	100.8

Interferences from medication or endogenous substances may affect results.¹⁷

Method Comparison

A study was performed based on guidance from CLSI EP09-A3 using the Passing-Bablok regression method.¹⁸

Sodium

	Units	n	Correlation	Concentration			
			Coefficient	Intercept	Slope	Range	
Alinity c	Serum	mmol/L	141	1.00	-0.50	1.00	101 - 197
Sodium vs ARCHITECT Sodium	Urine	mmol/L	101	1.00	0.88	0.99	22 - 386

Potassium

	Units	n	Correlation Coefficient		Concentration Range		
			Intercept	Slope			
Alinity c Potassium vs ARCHITECT Potassium	Serum	mmol/L	122	1.00	0.00	1.00	1.3 - 9.4
	Urine	mmol/L	107	1.00	-1.23	1.05	4.3 - 266.1

Chloride

	Units		n	Correlation Coefficient	Intercept	Slope	Concentration
							Range
Alinity c Chloride vs ARCHITECT Chloride	Serum	mmol/L	120	1.00	0.00	1.00	52 - 148
	Urine	mmol/L	112	1.00	-0.30	1.00	25 - 299

BIBLIOGRAPHY

1. Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. St. Louis, MO: Elsevier Saunders; 2012:807-810, 812-813.
2. Jacobs DS, DeMott WR, Oxley DK, et al. *Laboratory Test Handbook*. Hudson, OH: Lexi-Comp; 2001:144-147, 258-260, 275-279.
3. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
4. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office; December 2009.
5. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
6. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
7. Kaplan LA, Pesce AJ, Kazmierczak SC. *Clinical Chemistry – Theory, Analysis, and Correlation*, 3rd ed. St. Louis, MO: Mosby, 1996:461–463.
8. McPherson RA, Bidorpeh EK, Castellani WJ, et al. *Analysis of Body Fluids in Clinical Chemistry; Approved Guideline (C49-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 2007.
9. Guder WG, da Fonseca-Wollheim F, Hell W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2001:22,40,44,52.
10. US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*, 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
11. Pesce AJ, Kaplan LA, editors. *Methods in Clinical Chemistry*. St. Louis, MO: CV Mosby, 1987:75.
12. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
13. Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. St. Louis, MO: Elsevier Saunders; 2012:2139,2164,2168.
14. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
15. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
16. Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP07-A2. Wayne, PA: CLSI; 2005.
17. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press, 2000:3-176, 3-644, 3-714.
18. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
	In Vitro Diagnostic Medical Device
	Lot Number
	Produced for Abbott by
	Product of USA
	Reagent 1
	List Number
	For use by or on the order of a physician only (applicable to USA classification only).
	Serial number

Alinity and ARCHITECT are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

Fisher Diagnostics,
A Div. of Fisher Scientific Company, LLC
A Part of Thermo Fisher Scientific, Inc.
8365 Valley Pike
Middletown VA 22645 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Created September 2016.
©2016 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity c ICT Sample Diluent (for Russia)

COMPOSITION:

TRIS buffer, magnesium acetate. Preservatives: Silver Probin Nitrate 150

pH: 6.5

SHIPPING

Recommended shipping conditions for device is ambient. Device is stable until the expiration date indicated on the package. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 14 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

EN ISO 9001, ISO 13485

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

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- Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. St. Louis, MO: Elsevier Saunders; 2012:807-810, 812-813.
- Jacobs DS, DeMott WR, Oxley DK, et al. *Laboratory Test Handbook*. Hudson, OH: Lexi-Comp; 2001:144-147, 258-260, 275-279.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Kaplan LA, Pesce AJ, Kazmierczak SC. *Clinical Chemistry – Theory, Analysis, and Correlation*, 3rd ed. St. Louis, MO: Mosby, 1996:461–463.
- McPherson RA, Bldkorpeh EK, Castellani WJ, et al. *Analysis of Body Fluids in Clinical Chemistry; Approved Guideline (C49-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards; 2007.
- Guder WG, da Fonseca-Woltheim F, Heil W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2001:22,40,44,52.
- US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*, 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
- Pesce AJ, Kaplan LA, editors. *Methods in Clinical Chemistry*. St. Louis, MO: CV Mosby, 1987:75.
- Westgard JO. *Basic QC Practices*, 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 5th ed. St. Louis, MO: Elsevier Saunders; 2012:2139,2164,2168.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP07-A2. Wayne, PA: CLSI; 2005.

- Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press, 2000:3-176, 3-644, 3-714.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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