NC0K770IND

NOTARIAL CERTIFICATE

TO ALL TO WHOM these presents shall come

I, Wee Soon Keng, NOTARY PUBLIC duly admitted, authorised to practise in the Republic of Singapore, DO HEREBY CERTIFY

that Chien Chiang for Signet Armorlite Amera Pte Ltd had on this day subscribed his signature to the document titled INSTRUCTIONS FOR USE, hereunto attached.

The said Chien Chiang affirms that he has the power and authority to execute the document on behalf of the company.

IN FAITH AND TESTIMONY whereof I the said notary have subscribed my name and set and affixed my seal of office at Singapore, this 14th day of September 2020.

**NOTARY PUBLIC
SINGAPORE**

By virtue of Rule 8(3)(c) of the Notaries Public Rules, a Notarial Certificate must be authenticated by the Singapore Academy of Law in order to be valid.

Verification code:
01022253

This **Authentication Certificate** only certifies the authenticity of the signature and the capacity of the person who signed the Notarial Certificate.

This **Authentication Certificate** is not valid if the seal of the Singapore Academy of Law is removed or altered in any way whatsoever. This Certificate does not authenticate or confirm the content of the Document attached to the annexed Notarial Certificate.

Authentication

1. Country:	Singapore
2. This public document has been signed by:	Wee Soon Keng
3. Acting in the capacity of:	Notary Public
4. Bears the seal/stamp:	Notary Public

Certified

5. Authentication Cert No.:	AC0K7802ME
6. At:	Singapore Academy of Law
7. The:	16th September 2020
8. By:	Melissa Goh, Deputy Director, SAL
9. SAL Certification Seal:	10. Signature:

Melissa Goh


Singapore
14.09.2020

INSTRUCTIONS FOR USE

Spectacle lenses made of polymer material: clear, tinted, photochromic, afocal, bifocal, transfocal, single-vision, polarizing with antireflex or reinforcing coating or without coating , assorted

Company Signet Armorlite Amera Pte . Ltd. . 201 Kallang Bahru . # 04-00, Essilor Building 1 Singapore 339338 / Signet Armorlite Amera Pte Company. Ltd., 201 Callang Bahru, no. 04-00, Essilor Building 1 , Singapore. 339338 .

Chien Chiang
Managing Director, Amera
Signet Armorlite Amera Pte Ltd

The Kodak trademark and trade dress are used by Signet Armorlite under license from Kodak ©2019

1. PURPOSE AND DESIGN OF THE MEDICAL DEVICE.

Spectacle lenses made of polymer material: clear, colored, photochromic, afocal, bifocal, transfocal, single-vision, polarizing with antireflex or reinforcing coating or without coating (hereinafter referred to as lenses, or product) is an optical product designed for the production of individual corrective glasses intended for the vision correction of various degrees of ametropia.

All types of lenses in a pair are combined with afocal lenses (supplied in a pair of lens + afocal lens) for **amblyopia**.

Scope of application - ophthalmological offices of hospitals, clinics, health institutions.

Indications

Afocal lenses - improve contrast, support accommodation function, disrupt the binocular system, etc.

Single vision lenses - myopia, hyperopia, astigmatism, presbyopia, binocular system disorders.

Bifocal lenses - myopia, hyperopia, astigmatism in presbyopia in the absence of the need for vision at medium distances.

Transfocal (progressive) lenses - myopia, hyperopia, astigmatism in presbyopia.

Contraindications

Absent

Side Effects: The only potential physical hazard is risk

laceration from broken lenses. However, the lenses have excellent impact resistance, but they can break with a strong application of physical force. Broken lenses should be disposed of properly.

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The Kodak trademark and trade dress are used by Signet Armorlite under license from Kodak ©2019

1.2. Assortment

Spectacle lenses made of polymer material: clear, tinted, photochromic, afocal, bifocal, transfocal, single-vision, polarized with antireflex or reinforcing coating or without coating, assorted options:

Spectacle lenses with a refraction from -40 to +30 diopters, diameter 50-80 mm:

- 1 Spectacle lenses made of polymer material, clear single vision
- 2 Spectacle lenses made of polymer material, clear bifocal
- 3 Clear transfocal spectacle lenses made of polymer material
- 4 Spectacle lenses made of polymer material colored single vision
- 5 Spectacle lenses made of polymer material, colored bifocal
- 6 Spectacle lenses made of polymer material, colored transfocal
- 7 Spectacle lenses made of polymer material, photochromic single vision
- 8 Spectacle lenses made of polymer material, photochromic bifocal
- 9 Spectacle lenses made of polymer material photochromic transfocal
- 10 Spectacle lenses made of polymer material polarizing single vision
- 11 Spectacle lenses made of polymer material polarizing bifocal
- 12 Spectacle lenses made of polymer material, polarizing zoom
- 13 Spectacle lenses made of polymer material, polarizing, photochromic, single-vision
- 14 Spectacle lenses made of polymeric material polarizing photochromic transfocal
- 15 Spectacle lenses made of polymer material polarizing photochromic bifocal
- 16 Spectacle lenses made of polymer material, polarized, colored, single-vision
- 17 Spectacle lenses made of polymer material polarized colored transfocal
- 18 Spectacle lenses made of polymer material polarized colored bifocal

Spectacle lenses with a refraction of 0 diopters, diameter 50-80 mm:

- 1 Spectacle lenses made of polymer material, clear afocal
- 2 Spectacle lenses made of polymer material, colored afocal
- 3 Spectacle lenses made of polymer material photochromic afocal
- 4 Spectacle lenses made of polymer material polarizing afocal
- 5 Spectacle lenses made of polymer material polarizing photochromic afocal

The Kodak trademark and trade dress are used by Signet Armorlite under license from Kodak ©2019

6 Spectacle lenses made of polymer material polarized colored afocal

2. TECHNICAL CHARACTERISTICS

Products are supplied non-sterile.

Spectacle lens characteristic name	Spectacle lens type and value of characteristics			
	Afocal	Single vision	Bifocal	Transfocal (progressive)
spherical (S) refraction value, diopters	0	-40 to +30, step 0.25	-40 to +30, step 0.25	-40 to +30, step 0.25
value of cylindrical (C) refraction, diopters	0	0 to 8.0, step 0.25	0 to 8.0, step 0.25	0 to 8.0, step 0.25
position of the cylinder axis / prism base (°)	0	from 0 to 360, step 1.0	from 0 to 360, step 1.0	from 0 to 360, step 1.0
value of prismatic action (prism), diopters	from 0 to 5, Step 0, 25	from 0 to 20, step 0, 25	from 0 to 20, step 0, 25	from 0 to 20, step 0, 25
nominal size, diameter, mm	50-80, step 1mm	50-80, step 1mm	50-80, step 1mm	50-80, step 1mm
additional refraction for near, diopters	-	-	from 0.25 to 4.0, step 0.25	-
prismatic action of the lens (prism), PDPT.	from 0 to 5, Step 0, 25	from 0.25 to 20.0, step 0.25	from 0.25 to 20.0, step 0.25	from 0.25 to 20.0, step 0.25
prismatic action of progressive spectacle lenses, diopters	-	-	-	from 0.25 to 10.0, step 0.25
additional refraction for near progressive spectacle lenses, diopters	-	-	-	from 0.75 to 3.5, step 0.25
information about lens thickness in the center, mm	from 1 to 5, step 0.1	from 1 to 20, step 0.1	from 1 to 20, step 0.1	from 1 to 20, step 0.1

Signet Armorlite

201 Kallang Bahru, #04-00, Singapore 339338

T: +65 6589 9399 F: +65 6292 0245

information about the thickness of the lenses along the edge, mm	from 0.5 to 7, step 0.1	from 0.5 to 18, step 0.1	from 0.5 to 18, step 0.1	from 0.5 to 18, step 0.1
--	-------------------------	--------------------------	--------------------------	--------------------------

The maximum permissible deviations of the refraction of single-vision lenses and zones for the distance of multifocal lenses from the nominal values should correspond to those given in Table 1.

Table 1 In diopters

Surface refraction at the <i>second</i> major meridian	Maximum deviation on the <i>first</i> major meridian	Maximum deviation of the absolute value of the astigmatic difference (cylinder)			
		from 0.00 to 0.75 incl.	from 0.75 to 4.00 incl.	from 4.00 to 6.00 incl.	over 6.00
0.00 to 3.00 incl.	± 0.09	± 0.09	± 0.12	± 0.18	-
3.00 to 6.00 incl.	± 0.12	± 0.12	± 0.12	± 0.18	± 0.25
6.00 to 9.00 incl.	± 0.12	± 0.12	± 0.18	± 0.18	± 0.25
From 9.00 to 12.00 incl.	± 0.18	± 0.12	± 0.18	± 0.25	± 0.25
12.00 to 20.00 incl.	± 0.25	± 0.18	± 0.25	± 0.25	± 0.25
More than 20.00	± 0.37	± 0.25	± 0.25	± 0.37	± 0.37

The maximum permissible deviations of the refraction of the zone for given progressive spectacle lenses *from the nominal values* should correspond to those given in table 2.

Table 2 In diopters

Surface refraction at the <i>second</i> major meridian	Maximum deviation on the <i>first</i> major meridian	Maximum deviation of the absolute value of the astigmatic difference (cylinder)			
		from 0.00	from	from 4.00	over 6.00

The Kodak trademark and trade dress are used by Signet Armorlite under license from Kodak ©2019

		to 0.75 incl.	0.75 to 4.00 incl.	to 6.00 incl.	
0.00 to 6.00 incl.	± 0.12	± 0.12	± 0.18	± 0.18	± 0.25
6.00 to 9.00 incl.	± 0.18	± 0.18	± 0.18	± 0.18	± 0.25
9.00 to 12.00 incl.	± 0.18	± 0.18	± 0.18	± 0.25	± 0.25
12.00 to 20.00 incl.	± 0.25	± 0.18	± 0.25	± 0.25	± 0.25
More than 20.00	± 0.37	± 0.25	± 0.25	± 0.37	± 0.37

The maximum permissible deviations of the cylinder axis position from the nominal values (refer to multifocal, progressive and single vision spectacle lenses with a predetermined orientation, for example, the position of the prism base), should correspond to the values indicated in Table 3.

Table 3

The absolute value of the astigmatic difference (cylinder), diopters	Less than 0.50	0.50 to 0.75 incl.	0.75 to 1.50 incl.	More than 1.50
Axis deviation limit	$\pm 7^\circ$	$\pm 5^\circ$	$\pm 3^\circ$	$\pm 2^\circ$

The maximum permissible deviations of the additional refraction of the zone for the vicinity of multifocal and progressive spectacle lenses from the nominal values should correspond to those given in Table 4.

Table 4 In diopters

The value of the additional refractive zone for near	Less than 4.00	More than 4.00
Maximum permissible deviation	± 0.12	± 0.18

At the base point for given, the maximum permissible deviations of the resulting prismatic action of the prescribed prism and, if present, a thinning prism, should not exceed the values given in Table 5. (For spectacle lenses without a specified prismatic action, the data in Table 5 represent the maximum permissible values of the undesirable prismatic action caused by the deviation of the position of the constructive base point from the calculated one.)

Table 5 In prism diopters

Prismatic action	Lenses		
	Single vision	Multifocal and progressive	
		Horizontally	Vertically
0.00 to 2.00	$\pm (0.25 + 0.1)$	$\pm (0.25 + 0.1)$	$\pm (0.25 + 0.05)$
2.00 to 10.00	$\pm (0.37 + 0.1)$	$\pm (0.37 + 0.1)$	$\pm (0.37 + 0.05)$
More than 10.00	$\pm (0.50 + 0.1)$	$\pm (0.50 + 0.1)$	$\pm (0.50 + 0.05)$
Note - is the highest absolute value of refraction at the main meridians.			

Note - An example of the application of the *zone* tolerances indicated in the table for the distance of a multifocal lens with a prescription : +0.50; -2.50; 20 ° with a prismatic action of no more than 2.00 pdpt.

With regard to spectacle lenses without a given prismatic action, the data in Table 5 represent the maximum permissible values of undesirable prismatic action caused by the deviation of the position of the constructive base point from the calculated one. For this recipe, the refractions on the main meridians are +0.50 diopters and -2.00 diopters. With the highest absolute value of refraction equal to 2.00 diopters , the horizontal deviation of the prismatic action is $\pm (0.25 + 0.1 \times 2.00) = \pm 0.45$ pdpt, the vertical deviation of the prismatic action is $\pm (0.25 + 0.05 \times 2, 00) = \pm 0.35$ pdpt.

Dimensional tolerances:

dimensional tolerances (for diameter) are calculated by the formula:

effective size: from nominal size -1 to usable size +2 mm;

useful size: more useful size -2mm.

Thickness tolerance (at the reference point of the front surface for distance normal to this surface) : The thickness should not differ from the nominal value by more than ± 0.3 mm.

Segment tolerances for multifocal lenses:

Any of the segment dimensions (width, depth and depth of the intermediate zone) must not differ from the nominal value by more than ± 0.5 mm.

When a matched pair of spectacle lenses is sold, any of the segment dimensions

(width, depth and depth of the intermediate zone) must not differ from each other by more than ± 0.7 mm.

The quality of processing of the refractive surfaces of the lens must meet the following requirements:

- The deviation of the lens surface shape (waves), distorting the image of the object under consideration, is not allowed.
- The width of the scratches (b) and their total length (l) must not exceed the values given in table 5:

Table 5

Lens area	b, mm	l, mm
Central, diameter 30mm	Over 0.006 to 0.01 incl.	ten
Regional	Over 0.006 to 0.02-incl.	ten

Scratches up to 0.006mm and dots with a diameter of 0.05mm are allowed if their area in a limited area with a diameter of 5mm does not exceed 0.1mm².

- The number (N) of bubbles, dots and other foreign inclusions, allowed only when the distance between them is more than 5 mm, and their diameter (d) should not exceed the values indicated in table 6:

Table 6

Lens area	d, mm	N, no more
Central, diameter 30mm	0.05 to 0.1 incl.	1
Regional	0.1 to 0.2 incl.	2

- In the annular zone 2 mm wide along the edge of the lens, defects, except for cracks, are not standardized. Cracks are not allowed.

The lenses must withstand the impact of a steel ball with a diameter of 22 mm, to which a force of (100 ± 2) N is applied.

UVB-100%.

The light transmittance of lenses at the design reference point (for uncolored lenses with coatings) must be at least 3%.

The light transmittance of spectacle lenses for driving at night must be at least 75%.

Depending on the color intensity of the lens, they are divided into the following categories 0,1,2,3,4, photochromic lens.

Not suitable for driving
Not suitable for night driving
universal

TRANSMISSION COEFFICIENT CATEGORY	LIMITATIONS ON USE OF THE LENS		

Lenses in category 0	X		
Lenses in category 1		X	
Lenses in category 2		X	
Lenses in category 3 if TV > 8.0%		X	
Lenses in category 4 if TV < 8.0%			X
Photochromic lens	X		

Lenses should be made from the specified materials (under the table, the designation of materials in Russian):

Material	Index	Abbe ne	Specific gravity	UV cut	Chemical components
1.5 (CR39)**	1.502	58	1.31 g/cm ³	355nm	
1.5 Photochromic (CR607)**	1.499	60	1.27 g/cm ³	383nm	diethylene glycol bis
1.5 Photochromic (CR39)	1.502	58	1.31 g/cm ³	393nm	(allyl carbonate)
1.5 Polarized (CR39)	1.502	58	1.31 g/cm ³	380nm	
1.6 (MR8)	1.601	42	1.30 g/cm ³	398nm	polythiourethanes
1.6 Photochromic (MR8)	1.601	42	1.30 g/cm ³	400nm	produced by a reaction
1.6 Polarized (MR8)	1.601	42	1.30 g/cm ³	400nm	of a polythiol with di
1.67 (MR7)	1.665	42	1.36 g/cm ³	395nm	isocyanates
1.67 Photochromic (MR7)	1.665	42	1.36 g/cm ³	400nm	
1.74 (N19)	1.710	31	1.40 g/cm ³	400nm	episulphides
1.74 Photochromic (N19)	1.710	31	1.40 g/cm ³	400nm	
Polycarbonate (PC)	1.584	31	1.20 g/cm ³	380nm	
Polycarbonate Photochromic (PC)	1.584	31	1.20 g/cm ³	390nm	polycarbonate
Polycarbonate Polarized (PC)	1.584	31	1.20 g/cm ³	390nm	

The Kodak trademark and trade dress are used by Signet Armorlite under license from Kodak ©2019

Signet Armorlite

201 Kallang Bahru, #04-00, Singapore 339338
T: +65 6589 9399 F: +65 6292 0245

using dyes: PhysioBrown brands (color class 0 (15% chromaticity), 1 (35% chromaticity), 2 (65% chromaticity), 3 (85% chromaticity), 4 (95% chromaticity), 2 gradient, 3 gradient), PhysioGrey (color class 0 (15% chroma), 1 (35% chroma), 2 (65% chroma), 3 (85% chroma), 4 (95% chroma), 2 gradient, 3 gradient), Yellow (class colors 0 (15% chromaticity), 1 (35% chromaticity), 2 (65% chromaticity), 3 (85% chromaticity), 4 (95% chromaticity) The chromaticity of photochromic lenses is from 5% to 85%.

Outer lens coatings consist of (under the table the designation of materials in Russian):

1.5 UV CLEAR	1,5	59	2.48 g/cm ³	silicon dioxide (SiO ₂), sodium oxide (Na ₂ O) from sodium carbonate (Na ₂ CO ₃), calcium oxide (CaO)
1.5 UNICROWN	1,5	58	2.61 g/cm ³	
1.6/41 White	1,6	42	2.63 g/cm ³	
1.7/35 White	1,7	34	3.21 g/cm ³	
1.7/42 White	1,7	42	3.21 g/cm ³	
1.8/35 White	1,8	34	3.65 g/cm ³	
1.9/31 White	1,9	30	3.99 g/cm ³	
UV DG 37	1,5	60	2.48 g/cm ³	
PHOTOGRAY EXTRA	1,5	57	2.41 g/cm ³	
PHOTOBROWN EXTRA	1,5	46	2.41 g/cm ³	
FUSIBLE PHOTOGRAY 16	1,6	42	2.70 g/cm ³	
XDF DARK GRAY	1,5	56	2.41 g/cm ³	
XDF LIGHT GRAY	1,5	57	2.41 g/cm ³	

All grades of colorants and coatings are from Company Signet Armorlite Amera Pte . Ltd .

3 . OPERATING CONDITIONS AND RULES.


3 .1. Lenses must be suitable for use at air temperatures from -45 to + 40 ° C, with relative humidity up to 80%.

3 .2. Operating rules:

From the lens in an electronic faceting machine (bevelling method), an eyeglass lens is machined into the shape of a spectacle frame.


The method of decorative processing of the edge of the lens by cutting the edge for a given frame.


After that, the frame is straightened in accordance with STB GOST, as well as the patient's anthropological data.

The basis for the manufacture of glasses is a prescription from an ophthalmologist.

An automatic dioptrimeter controls the glasses being manufactured, and prints a check based on the result of the check, which confirms that the glasses are made in accordance with the prescription.

Glasses touch the skin, we wear glasses in smoky rooms, we cook food in them and look into a boiling saucepan, we leave them without a case and dust settles on them. If you look closely, all these traces of wearing can be seen on the nose pads, in the places where the lens touches the frame.

According to this, once in 2 years, you need to change the frame or make ultrasonic cleaning of glasses and lenses.

Ultrasonic cleaning is a high-tech procedure that allows you to clean the lenses and frames of glasses from various types of dirt. This procedure will not do any damage to your lenses.

3.3. Care rules.

Reusable Microfiber Lens Care Cloths

They are composed of ultra-fine fibers. With a small diameter of just 0.06 mm,

microfibers can penetrate the micropores of the lens, removing any liquid or particulate matter faster than any other tissue. Microfiber fibers literally suck dust and dirt into their matrix and hold them securely. It can absorb and retain dirt and moisture up to 7 times its own weight. The main rule of effective use of these wipes is to remember to keep them clean, as a dirty wipe can scratch the lens coating.

Disposable wet wipes

They are packed in sealed bags, which is very convenient for those who save time. These wipes contain a cleaner and a volatile alcohol, and are also composed of a thin non-woven material to prevent scratching of the lens coating.

Cleaning sprays

The cleaning ability of sprays is due to the content of surfactants - special chemical compounds. The molecules of these substances readily come into contact with contaminants, for example, fats and oils that are difficult to wash, and remove them without difficulty. After spray application, the lens should be wiped with a reusable microfiber cloth. It is important to note that lens cleaning sprays do not contain alkaline additives that can damage the surface of plastic, polycarbonate, and highly refractive organic materials.

In addition to cleaning, lenses require daily care. There are a few simple rules to follow:

1. Store lenses only in their original packaging. Never put it in any packaging next to other things, as this can scratch the lenses.
2. Do not put lenses down. This way you will protect them from scratches.
3. Make sure that loose lenses are not exposed to direct sunlight. Heating the lens in direct sunlight can cause a fire.
4. Do not throw your lenses.
5. It is not allowed to lay lenses with their convex side on a hard surface;
6. It is not allowed to apply significant mechanical stress to the lenses;
7. It is not allowed to heat lenses to temperatures above + 60 ° C, including using them in a bath, sauna.

4. MAINTENANCE INFORMATION

The product is warranted for a period of 18 months from the date of purchase, it is guaranteed that the product is free from defects in workmanship and quality of manufacturing materials. Make sure that any malfunction does not occur due to non-compliance with the rules of operation. In the event of a malfunction, while observing the operating requirements in accordance with the instructions for use, the manufacturer replaces the product. For all questions please contact: LLC "USRM": Legal entity, address / fact. address: 127254 Moscow, Ogorodny proezd 20A, building 3, room 20, TIN 7732505764, checkpoint 771501001, OGRN 1057748396023, OKPO 78966933, tel / fax +7 (495) 799-90-91, ru-regulatory@essilor.com

5. ENVIRONMENTAL PROTECTION REQUIREMENTS.

The main type of possible hazardous impact on the environment is the pollution of atmospheric air in populated areas, soil and water as a result of:

- emergency leaks (placers) of production materials;
- unorganized waste disposal on the territory of the manufacturer or outside it;
- their arbitrary dump in places not intended for these purposes .

The materials used in the manufacture of the product should not pose a danger to life, health of people and the environment both during operation and after the end of its life.

6. LABELING AND PACKAGING

The product is supplied non-sterile, 1 and 2 pcs. in the packaging bag.

6.1. The following information must be printed on each package:

- optical effect of the lens: the value of spherical (S) and cylindrical (C) refraction, diopters
- nominal size, mm;
- color;
- coverage identification;
- trade name and location of the manufacturer;
- refractive index;
- **number and date of the registration certificate;**

For multifocal spectacle lenses, the following additional data should be indicated:

- additional refraction for near, diopters;
- type and width of the segment, mm, for example, D ;
- prismatic action of the lens (prism), PDPT;
- the method of measuring additional refraction for near, if it is not measured from the side of the segment;
- the position of the base point for distance in aspherical multifocal lenses.

For progressive spectacle lenses, the following additional data should be indicated:

- additional refraction for near, diopters;
- indication of the place of installation of the lens in the frame (right or left);
- the method of measuring the additional refraction for near, if it is not measured from the progressive side.

6.2. Requirements for marking a progressive spectacle lens:

Permanent (permanent) markings

The following permanent markings must be applied to the spectacle lens:

- a) two alignment points located at a distance of 34 mm from each other, equidistant from the vertical plane passing through the center of the alignment crosshair or the base point of the prism;
- b) the value of additional refraction for near, diopters;
- c) manufacturer's or supplier's logo or trade name or trade mark.

Additional (removable) markings The following removable markings are recommended:

- a) designation of adjustment points;
- b) designation of the base point for distance;
- c) designation of a constructive base point for near;
- d) designation of the installation point (installation crosshair);
- e) designation of the base point of the prism .

The product must be placed in the appropriate slot of **individual** packaging, and fixed in it during transportation (**packaging material is cardboard, for example of type T-23, the material of the retainer (for fixation during transportation) is cardboard, for example of type T-23**: Each package with the product along with operational documentation for delivery , should be packed in a box or wrapped in paper. At the customer's request, operational documentation can be provided in electronic form. The box or bag must be tied with twine or pasted over with paper tape, glue on a paper base or polyethylene tape with an adhesive layer so that the packaging cannot be opened without breaking the integrity.

For transportation, canned products in packages with operational documentation must be packed in boxes made of sheet wood materials. Board boxes should be lined with waterproof material inside.

Transport markings depending on the types of vehicles used. Marking is applied directly to containers or labels: plywood, metal, paper or fiberboard. In this case, the manipulation signs applied to the container must correspond to the values: "Fragile. Caution", "Protect from moisture".

Symbols that may appear on the packaging, including handling marks :

Symbol	Description
	Batch code
	Manufacturer
	date of manufacture

The Kodak trademark and trade dress are used by Signet Armorlite under license from Kodak ©2019

	Not sterile
	Fragile, handle with care
	Keep away from moisture

7. TRANSPORTION AND STORAGE

The product is transported by all types of covered vehicles in accordance with the requirements and rules for the carriage of goods in force for each type of transport. Transportation of cargo by sea transport must be carried out in accordance with the "Safety Rules for the Sea Carriage of General Cargo".

Temperature limits: -50 to + 50 ° C

Relative humidity up to 100 %, with no condensation moisture.

Lenses are stored in warehouses with conditions assuring protection from the atmospheric precipitation and aggressive environment, at a distance of at least 1 m from heating and heating devices for heated warehouses.

Storage temperature from +5 to + 40 ° C

Relative humidity up to 80 %, with no condensation moisture.

Service life is not less than 5 years.

8. MANUFACTURER'S WARRANTIES .

8.1. The manufacturer guarantees lens high quality, in accordance with the associated requirements and technical documentation when conditions of transportation, storage and usage are appropriate.

8.2. The warranty period is 18 months from the date of receipt of the lenses by the customer.

9 . REPAIR INFORMATION . INFORMATION ABOUT CLAIMS.

In the event of a factory defect it must be claimed and replaced.

The complaints to be sent to:

LLC "EGRM": Legal address: 127254 Moscow, Ogorodny proezd 20A, building 3, room 20, TIN 7732505764, IN 771501001, OGRN 1057748396023, OKPO 78966933, tel / fax +7 (495) 799-90-91, ru-regulatory@essilor.com

10 . RECYCLING INFORMATION.

When disposing of the products comply with the requirements of SanPiN 2.1.7.2790-10 , Class A . Also, for disposal, please contact the official representative of the manufacturer:

LLC "EGRM": Legal address: 127254 Moscow, Ogorodny proezd 20A, building 3, room 20, TIN 7732505764 , IN 771501001 , OGRN 1057748396023, OKPO 78966933 , tel / fax +7 (495) 799-90-91, ru-regulatory@essilor.com

11 . THE PRODUCT MEETS STANDARDS :

ISO 8980-1: 2004 "Ophthalmic optics - Uncut finished spectacle lenses - Part 1: Specifications for single-vision and multifocal lenses", MOD 1: 2006;

ISO 8980-2: 2004 "Ophthalmic optics - Uncut finished spectacle lenses - Part 2: Specifications for progressive power lenses", MOD 1: 2006;

ISO 14889: 2003 "Ophthalmic optics - Spectacle lenses - Fundamental requirements for uncut finished lenses", MOD

Signet Armorlite

201 Kallang Bahru, #04-00, Singapore 339338
T: +65 6589 9399 F: +65 6292 0245

EXAMPLE OF A WARRANTY CARD:

Name _____

Model _____

Date of purchase _____ l.m.

Malfunction _____

Warranty center mark

model	malfunction	decision	result

_____ Warranty center stamp

_____ Warranty center stamp

_____ Warranty center stamp

The Kodak trademark and trade dress are used by Signet Armorlite under license from Kodak ©2019

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