USER'S MANUAL

**Mobile elevator Rifton TRAM for vertical adjustment and movement of patients,
with accessories**

Clement Huley

Authorised Signature

Kirk Wareham

**Kirk Wareham, VP and General Manager
Community Products, LLC, USA**

Contents:

- 1. Warnings and important information**
- 2. Prescription**
- 3. Product information**
 - 3.1. Regulatory**
 - 3.2. Design of mobile elevator Rifton TRAM**
 - 3.3. Technical data**
 - 3.4. Accessories**
- 4. Initial assembly**
- 5. Directions for use**
- 6. Operation and transfers**
 - 6.1. Actions to be taken before and during each transfer**
 - 6.2. Transfer to sitting position using thigh straps**
 - 6.3. Using walking saddle**
 - 6.4. Toileting and the TRAM**
- 7. Troubleshooting**
 - 7.1. Checks**
 - 7.2. Current repair**
- 8. Maintenance conditions and repair conditions. Conditions of storage, operation and transportation.**
 - 8.1. Maintenance conditions and repair conditions.**
 - 8.2. Conditions of storage, operation and transportation.**
- 9. Cleaning**
- 10. Utilization**
- 11. Guarantee certificates / Service life**

1. Warnings and important information

WARNINGS!!!

- Thoroughly read and understand the information in this product manual before attempting to use this product. If the procedures and instructions in this manual are not followed, serious injury could occur.
- The use of this product is subject to approval from qualified medical personnel and must be supervised by specialist at all times.
- This product is intended for use by children and adults with cognitive and corporal defects.
- This product is intended for indoor use only. It must be used on flat level surface only. This product must not be used in and around water, near staircases, swimming-pools, etc.
- Straps and supports are provided for the safety of the user and must be carefully adjusted for comfort and security. The use of straps and supports must be supervised by medical personnel at all times.
- Do not leave the patient in the elevator unsupervised.
- Slippery clothing, patient's low muscular tone, patient's specific body structure may cause the TRAM to slide up on the patient, making a safe transfer difficult. In this case immediately stop lifting.
- Do not lift the patient higher than necessary for transfer.
- Do not use the elevator in vertically confined space, always ensure there is free space over patient's head before lifting.
- The initial and destination points of transfer should be as close as possible.
- Do not use the elevator if battery status indicator is red.
- Before lifting keep your hands and feet away from the lift piston.
- Before extending or retracting the base frame keep your hands and feet away at the safe distance.
- Transmission devices (e.g. cell phones) may create disturbances and should not be used during the use of elevator.

IMPORTANT: Please save this product manual for future use.

2. Prescription

Elevator Rifton TRAM is a mobile transfer device. Zero lifting enables a single caregiver to transfer a patient weighing up to 350 lbs. (160 kg) to and from a wheelchair, bed, chair, toilet or gait trainer. Elevator Rifton TRAM allows transferring the patient to standing position, rendering support when standing or moving.

3. Product information

3.1. Regulatory

Mobile elevator Rifton TRAM with accessories for vertical adjustment and movement of patients conforms to the EU Medical Device Directive 93/42/EEC. This device conforms to the requirements of EN 12182.

The Quality management system of the Company Community Products, LLC, USA conforms to ISO 13485:2003.

Mobile elevator Rifton TRAM with accessories for vertical adjustment and movement of patients is produced according to the Quality management system of Company Community Products, LLC, USA; Risk Analysis Corporate standard RS-107; Failure Mode and Effects Analysis. The Risk Analysis program of the Company Community Products, LLC, USA conforms to ISO 14971.

3.2. Design of mobile elevator Rifton TRAM



Figure 1. Dimensions of mobile elevator Rifton Rifton TRAM

3.3. Technical data

Technical data of mobile elevator Rifton TRAM	
Patient's max. height	193 cm
Minimum distance from the floor to armpit	77 cm
Waist circumference	56 to 152 cm
Overall length	116 cm
Overall width	70 cm
Width with legs extended	102 cm
Overall height with support handles	110 to 146 cm
Distance from the floor to the top of base frame (for passage under furniture)	17,5 cm

Turning diameter	127 cm
Overall weight of the device	32 kg
Patient's max. weight	160 kg
Lifting/lowering speed (no load)	4 cm/s
Battery	Detachable, sealed lead-acid gel-type batteries 24 V, 2,9 Ah
Battery charger	Wall-mounted charger, 100-240 VAC, max 650 mA
Motor	24V, 6A, permanent magnet motor
Emergency lowering	Electrical
Wheels	Front: 100mm dual with brake. Rear: 100mm dual
Duty cycle	One minute continuous use followed by 19 minutes idle
Protection degree	IPX6

3.4. Accessories

Frame of mobile elevator Rifton TRAM is a metal structure with electrical lifting device, intended for fixing the body support. The frame of mobile elevator Rifton TRAM consists of base frame, lifting column and handle of base frame. Frame of mobile elevator Rifton TRAM may be optionally fitted with display. Built-in display shows the load on patient's body during transfer.

Control box is built in the elevator frame and allows controlling the elevator: lower or lift its moving part as necessary. The control is brought out to the right elevator handle.

Battery – Detachable, sealed lead-acid gel-type battery 24 V, 2.9 Ah. A fully charged battery allows approximately 70 full duty lifting cycles. Service life of the battery is 5 years.

Additional battery - Detachable, sealed lead-acid gel-type battery 24 V, 2.9 Ah. This battery is required to use the elevator all day long, if the main battery is discharged. Service life of the battery is 5 years.

Battery charger is designed to charge the battery. This is a wall-mounted battery charger, 100-240 VAC, maximum 650 mA.

Additional battery charger is designed to charge the backup (additional) battery in parallel with charging the main battery.

User's manual is provided in hardcopy.

Forearm supports, small and big size are used to support the weight of the patient when using the elevator as auxiliary device for walking training, supporting the body leaning forward, for controlling body and head, for positioning of the patient, for additional support to the patients who have low tone or need arm support. Forearm supports can be adjusted as follows:

- for height,
- horizontally: back and forth, left and right to patient's body,
- turning: left and right, closer and farther from patients body,
- up and down at any angle.

Forearm support, small size: Soft surface area 11x18 cm, adjustment range from elbow to hand – 20 to 32 cm.

Forearm support, big size: Soft surface area 13x22 cm, adjustment range from elbow to hand – 27 to 41 cm.

Body support ensures torso support during transfer and walking. Body support is installed on the elevator frame, includes body support pads and ring clips to fix the straps

Dimensions of body support: overall width 55 cm, support circumference is adjustable – from 56 to 152 cm.

Walking saddle, small and big size is used to transfer the patient from a sitting to a standing position and for support during movement (walking training). The walking saddle can be easily installed using rings and body support parts with colored marking. The basic parameter for determining the size of walking saddle: half of thigh circumference. If the half of patient's thigh circumference is 10 to 20 cm, then walking saddle of small size should be used. If the half of patient's thigh circumference is 20 to 35 cm then walking saddle of big size should be used.

Thigh straps are fixed to body support and used to transfer the patient to a sitting position.

Wheel direction lock ensures straightforward movement of Rifton elevator. The wheel direction lock can be installed on any wheel. To ensure straightforward

movement of Rifton elevator the locks should be installed at least on two wheels. If unlocked, the wheels rotate freely, which is necessary for maneuverability in narrow places.

Wheels ensure mobility of elevator frame. Front wheels: 100 mm, with a brake, double; rear wheels: 100 mm, double. Material – nylon wheels, tires filled with polyurethane foam.

Network cable – to plug the battery charger into wall outlet.

Body support pads – these pads can wrap up the patient's body and prevent excessive compression of the body when tightening the safety belt.

Ring clips are used to fix thigh straps to body support.

Patient handle. During lifting the patient can hold by the patient handle. The patient handle is fixed to the elevator frame.

Safety belt with a buckle allows fixing the body and preventing the patient from sliding during the transfer. The safety belt is fitted with a dual action safety buckle to prevent accidental release.

4. Initial assembly

The caregiver push handle must be attached to the Rifton TRAM frame with four bolts: With the Allen wrench remove the four bolts (Figure 2).

Figure 2


Hold the caregiver push handle against the attachment point as shown on Figure 3, insert the bolts in the holes, and starting with the top bolts, tighten them snugly with the Allen wrench (Figure 3).

Figure 3

Finally, remove the cover on the lower front of the battery box, and insert the electrical jack plug into its receptacle on the actuator control box (Figure 4). The release tab on the plug must be facing towards you. Replace the cover.

Figure 4


5. Directions for use

Caregiver push handle

The caregiver push handle enables the caregiver to maneuver the TRAM. It includes the up/down switch, located on top of the right hand-grip, through which the caregiver can adjust the height of the body support. Push the rocker switch up to raise the body support, down to lower it. (Figure 5).

WARNING! To prevent pinching or crushing, ensure that all hands and feet are clear of the vertical lifting column before activating the up/down switch.


Figure 5

Battery

The battery box houses the microprocessor and rechargeable battery which power and control the TRAM's electric actuator (see Figure 6).

Figure 6

To insert the battery, push the red, elliptical battery eject/emergency stop button (see Figure 7) and lift the lid, slide the battery down into the box and close the lid until the latch engages.

To remove the battery, push the red, elliptical battery eject/emergency stop button and slide battery up and out of the box (see Figure 8).

Figure 7



Figure 8


WARNING! To prevent inadvertent operation of the up/down switch when the TRAM is not in use, press the elliptical, red battery eject button to disconnect the battery.

CAUTION! To prevent long-term battery damage, remove the battery from the battery box and place it in the charger every night, even if the battery level indicator is still green.

Battery level indicator

The battery level light indicates the charge remaining in the battery. A green light means the charge level is adequate for use, yellow indicates that the battery should be charged soon, and red indicates very low charge (see Figure 9). An audible beep indicates that the battery is entirely depleted. A fully charged battery will give approximately 70 complete lift cycles.

Figure 9

To prevent damage to the battery, charge it when the light turns yellow. The indicator light will extinguish 10 minutes after the last activation of the up/down switch. It will turn on again if the up/down switch is used or if a battery is inserted into the battery box.

Notice: If the battery level becomes critically low, indicated by the battery level light turning red, the microprocessor will disable the up/down switch. However, the emergency lowering button will continue to function.

CAUTION! Charge or change the battery immediately if the battery level light turns red; do not attempt further lifts or transfers with a red light.

CAUTION! Prior to every use of the TRAM, activate the up/down button to turn on the battery level light, and verify that the light is green. If it is not, change or charge the battery.

Emergency stop button

To stop the TRAM motor in the event of an emergency, press the elliptical red emergency stop button, located on the front of the battery box. To reset it, close the lid of the battery box (see Figure 7).

Emergency lowering button

The emergency lowering button (A), located on the front of the battery box, provides a means of lowering the body support if the up/down switch fails. The emergency lowering button will continue to function even if the battery level has dropped below the threshold for operation of the up/down switch. It overrides all other inputs.

CAUTION! All caregivers should become familiar with the location and operation of the emergency stop button and the emergency lowering button before using the TRAM.

Battery charger and charging the batteries

Installation

1. Attach the battery charger to the wall near an outlet, using two screws (A) (see Figure 10). Two screws are provided with the charger, however, they may not be suitable for every situation.
2. Place battery in battery charger (Figure 11).
3. Plug the battery charger cord into wall outlet.

Figure 10



Figure 11

Charging

Charge batteries 24 hours before first use.

When the battery needs charging, an indicator light on the hand pendant will go on.

Charge batteries every night, even if indicator light is green.

To charge the battery, remove the battery from the control box, and place it onto the wall-mounted charger.

The battery charger and indicator light will shut off automatically when charging is complete.

Charging normally takes approximately six hours.

During long periods of inactivity or storage, batteries will lose charge.

Allowing batteries to deep cycle (become nearly or completely dead) will destroy them. To prevent this, store batteries in the battery charger.

An additional charging device is required to charge the additional (backup) battery.

WARNING: Do not charge batteries in a wet area.

Old batteries should be disposed of properly at an appropriate recycling facility.

Display

The display has two purposes: it can be used to measure and track a patient's weight, and to measure the amount of weight a patient is bearing during ambulation. (Figure 12).

Figure 12

To use the display, turn on the display and use the lb/kg button to display either pounds or kilograms. Before approaching the patient who will be weighed, attach all the components you will need for this particular patient, (leg straps, arm supports etc.) and then press the ZERO button to zero the display. Proceed to lift the patient. Once a patient is completely supported by the TRAM, the patient's weight will be displayed.

Notice: To obtain an accurate weight, the patient must be lifted clear of all weight-bearing surfaces, and the feet or legs must not be in contact with the base frame of the TRAM.

There are two ways to measure weight-bearing during ambulation:

Method 1: Turn on and zero the display as previously described. Once the patient is lifted and in the walking position, the weight displayed is the amount of weight that the TRAM is bearing.

Method 2: Using the thigh supports and the seated transfer procedure, lift the patient with the TRAM. With the patient's full weight displayed on the display, press the zero button. Next, re-position the patient for ambulation. The weight displayed will now be the amount of weight the patient is bearing (ignore the minus sign on the display).

To replace the display batteries, push outward on the small tab at the rear of the display housing and lift the battery cover (a small, flat head screwdriver inserted into the slot at the back of the lid may be used to gently pry it up). Insert fresh AA Alkaline batteries; the cells must be oriented as shown by the inscription in each battery receptacle (Figure 13).

Figure 13


Notice:

The display uses four AA Alkaline batteries. The battery life is approximately 100 hours of continuous use, or approximately 500 weight recordings. The display will automatically turn off if it does not sense a change in weight for more than 10 minutes.


Direction lock of a wheel

The direction lock enables the caregiver to prevent one or more of the TRAM's wheels from swiveling. Locking a front wheel makes the TRAM easier for medical personnel to maneuver during transfers. Locking a rear wheel can stabilize the TRAM when used for ambulation.

To install a direction lock on a wheel, insert the prongs of the direction lock into the slot found at the ends of each base frame with the red plunger pull handle up (Figure 14). The prongs are fully inserted in the slot when it clicks.

Figure 14

To engage the direction lock, align the wheel with the base frame, and push the plunger all the way down between the wheel's tires (see Figure 15).


Figure 15

Pull the red plunger up to disengage the direction lock.

Base frame

The handle of base frame adjusts the width of the base frame from 70 cm to 102 cm). It expands so that patients can be lifted from wide chairs and wheelchairs, and retracts for maneuverability, and to enable the TRAM to pass through narrow doorways. Swing the handle of base frame to the right to expand the frame, and to the left to close it (see Figure 16).

Figure 16


WARNING: To avoid injury, ensure that all hands and feet are clear of the handle of base frame and the base legs before expanding or retracting the base.

Notice: Disengage direction locks prior to moving the handle of base frame.

Body support


The TRAM's body support includes the body support pads, patient hand grips, ring clips, and the safety belt and buckle. (Figure 17).

Figure 17

The back belt has a dual action safety buckle to prevent accidental release. Press the two release tabs simultaneously to release (Figure 18).

Figure 18

The following can also be attached to the body support, depending on how the TRAM is to be used:

- Thigh straps (Figure 19)
- Walking saddle (Figure 20)
- Forearm supports (Figure 21)

Figure 19




Figure 20

Forearm supports

Forearm supports give extra lifting support for patients who have low tone and provide positioning support during ambulation. Forearm supports adjust in many directions and angles to accommodate different positioning requirements (Figure 21).

Figure 21


Installation of forearm supports

Notice: If forearm supports were purchased initially with the TRAM they will be attached and ready for use.


If you want to remove the forearm supports, unscrew the wing knob and the screw fixing them (Figure 22).

Figure 22


To prevent loss, screw the wing knob back on the threaded stud (Figure 23).

Figure 23

With no mounting brackets on the frame of mobile elevator Rifton TRAM, install the oval shaped hole cover by inserting the threaded stud through center hole and tightening the wing knob inside the body support channel onto the stud (Figure 24).

Figure 24.

If you purchased forearm supports separately you will need to install the mounting brackets onto the body support.

Remove the plastic hole cover from the body support channel by unscrewing wing knob inside the channel. To prevent loss, screw the wing knob back on the threaded stud.

With the large, lipped end up, install the mounting bracket on the body support channel and securely tighten the attachment knob on the threaded stud, which will now be protruding inside the body support channel.

Adjusting the forearm supports on the frame of mobile elevator Rifton TRAM

Height adjustment:

Press button (A) and slide post to desired position. Pressing A butting allows sliding the pad on the slide post backward and forward along the elevator (Figure 25).

Loosen knob (B) to:

- Slide arm pad toward or away from the user.
- Rotate up or down.
- Rotate in or out.

- Move the arm pad backward or forward.

Loosening knob B the user can adjust the position of the pads according to individual needs of each patient, without changing the axis.


Loosen knob C to adjust the patient handle (Figure 25):

1. Slide patient handle forward or back for different forearm lengths.
2. Rotate the patient handle from side to side.

Arm strap (Figure 25-D) and wrist strap (Figure 25-E) secure the patient's arm in the forearm support. The wrist strap prevents the patient's arm from coming out of the forearm support.

The wrist strap is located before the pad on which the patient's arm rests when using the forearm supports, the arm strap is located in the middle of the pad.

Figure 25.

6. Operation and transfers

6.1. Actions to be taken before and during each transfer

Before each transfer

- Thoroughly inspect the TRAM body support for tears, rips and worn areas, giving special attention to the straps and buckles. Remove the product from service if any condition develops that might make use unsafe.
- Move the origin and the destination of the transfer as close to each other as possible.
- Ensure that the body support is positioned comfortably, correctly and securely, and that the safety buckle is properly engaged.
- Lock the wheels.

During each transfer

- Lift the patient only high enough to perform the transfer.
- Check that the straps are snug and the body support does not slide up or cause discomfort.

6.2. Using thigh straps for a seated transfer

Clip one end of each thigh strap onto the yellow clips at the front of the TRAM's body support (Figure 26).

Figure 26


With the patient sitting up and arms raised slightly, position the TRAM's body support pads around the patient's rib cage a few inches below the armpits (Figure 27)

Figure 27.

Secure snugly with the buckle (Figure 28)

Figure 28

Attach the thigh straps (gray side up) to one of the colored clips along each side of the body support. For patients with lower tone, choose a clip further back for greater support. Pull the straps snugly, making sure both thigh pads are positioned under the patient's thigh and adjusted equally (Figure 29)

Figure 29

Raise patient. Watch closely to make sure the patient remains comfortable and does not sag. Move patient to transfer location. The TRAM base legs should be retracted for travel. (Figure 30).

Figure 30.

Gently lower patient onto seat. Avoid lowering too far so body support doesn't push down on hips. Unclip the back ring on the thigh straps and pull straps out from under patient. Then release back buckle (Figure 31).

Figure 31

6.3. Using walking saddle

Position walking saddle under patient, either by lifting patient, or by folding saddle, shifting the patient to one side and positioning saddle under patient one half at a time. (Figure 32)

Figure 32.

Attach walking saddle rings to the colored clips on the body support (Figure 33).

Figure 33.


For some patients positioning will be improved by crossing the front and rear straps as shown (Figure 34)

Figure 34

As you raise the patient, pull the TRAM gently toward you to mimic the natural sit-to-stand arc. Watch closely to make sure the patient is comfortable and well supported throughout the lift (Figure 35)

Figure 35

Adjust body support height so patient can walk comfortably; body support should not be positioned too high under the armpits (shorten walking saddle straps as necessary). Check to make sure the saddle is positioned under the patient as shown (Figure 36). If desired, use forearm supports for additional support or positioning.

Figure 36


6.4. Toileting and the TRAM

Mobile elevator Rifton TRAM may be used for toileting by the patients able to carry the weight load.


Position the patient in the TRAM over the toilet seat, leaving enough room to adjust clothing (Figure 37)

Figure 37


Loosen and remove the supporting straps one at a time so the patient's feet are touching the floor and bearing some weight (Figure 38)

Figure 38


Adjust clothing for toileting (Figure 39)

Figure 39

Using the up/down switch, lower the patient onto the seat. Lock the brakes on the TRAM's front wheels (Figure 40)

Figure 40

7. Troubleshooting

7.1. Checks

Daily check

- Visual damage or wear of elevator Rifton TRAM and accessories.
- Cracks, breaks, missing or loose details and/or improper operation.
- Correct functioning of emergency stop.

Remove the product from service if any condition develops that might make use unsafe.

The product may only be repaired using the original parts.

7.2. Current repair

Lift does not go up or down:

1. Make sure that the red emergency stop button on the control box is not pushed.
Rotate the button clockwise to start again.
2. Make sure the electrical cables are secure.
3. Make sure the battery box lid is closed.
4. Make sure the battery is charged and inserted correctly.
5. Contact Community Products, LLC, USA or its authorised representative RehabMedical, LLC, Russia.

Battery does not charge:

1. Check that the battery contact plate is not damaged or broken.
2. Check that wall charger is plugged securely into the wall socket.

Display does not turn on:

1. Remove and reinsert batteries, making sure that each cell is correctly oriented in polarity.
2. Batteries may be dead and need replacing.

3. Contact Community Products, LLC, USA or its authorised representative RehabMedical, LLC, Russia.

Spare parts replacement

If necessary, the manufacturer will supply any spare parts and repair manuals for the product. The customer is responsible to make sure that the repairs or modifications are correctly and safely completed. To place the order for spare parts please identify the serial number (Item №) of your product. Contact us at the numbers listed in paragraph 11 of this manual. Use spare parts manufactured only by Community Products, LLC, USA.

8. Maintenance conditions and repair conditions. Conditions of storage, operation and transportation

8.1. Maintenance conditions and repair conditions

- Temperature from +14 °C to +40 °C
- Relative humidity up to 90%
- Atmospheric pressure from 680 to 790 Mm hg.

8.2. Conditions of storage, operation and transportation

Operating conditions for medical device

- Temperature from +14°C to +40 °C
- Humidity up to 90%
- Atmospheric pressure from 680 to 790 Mm hg.

Storage conditions

Elevator Rifton TRAM with accessories is stored in protected from direct sunshine dry place at air temperature from - 15°C to +45 °C, humidity to 90%, atmospheric pressure from 680 to 790 Mm hg.

Conditions for transportation

Elevator Rifton TRAM with accessories is transported in covered vehicles, at temperature from -15 °C to +45 °C, at humidity to 90%.

9. Cleaning

Elevator Rifton TRAM should be cleaned with a moist rag using soap solution, mild detergent or 10% bleach solution. Mild disinfectants may be used to clean the forearm supports, body support, walking saddles, thigh straps and body support pads.

10. Utilization

Elevator Rifton TRAM and accessories thereto do not contain any substances dangerous to human life, health and environment. All materials used for packaging Elevator Rifton TRAM and accessories thereto are fully reusable. Packaging materials should be disposed of in an environment-friendly manner. Elevator Rifton TRAM and accessories thereto should be disposed of in accordance with the applicable law and local regulations and instructions. The run-down battery packs should be disposed of in accordance with clause 273 40 Code of Federal Regulations (CFR).

11. Guarantee certificates / Service life

Warranty: 1 year after delivery.

Subject to correct use the service life of the product is 5 years. After 5 years the service life of the product may be extended once the battery is replaced.

Should any questions, failures or operating problems arise, feel free to contact Community Products, LLC, 2032 Route 213, Rifton, NY 12471, USA, + 1 845-658-7700 / + 1 845-658-3144 or its authorised representative RehabMedical, LLC, Russia, 20a/1, Novoalekseevskaya st, Moscow, 129626, Russia, phone / fax: +7 (495) 686-40-47. Community Products, LLC, USA and its authorised representative RehabMedical, LLC, Russia cannot accept responsibility for any modifications or alterations made to mobile elevators Rifton after they leave manufacturer's premises. Community Products, LLC, USA and its authorised representative RehabMedical, LLC, Russia do not accept responsibility for use of elevator Rifton jointly with products of other manufacturers.

USER'S MANUAL

**Mobile elevator Rifton SoloLift for vertical adjustment and movement of
patients, with accessories**

Authorised Signature

**Kirk Wareham, VP and General Manager
Community Products, LLC, USA**

Contents:

- 1. Warnings and important information**
- 2. Prescription**
- 3. Product information**
 - 3.1. Regulatory**
 - 3.2. Design of mobile elevator Rifton SoloLift**
 - 3.3. Technical data**
 - 3.4. Accessories**
- 4. Initial assembly**
- 5. Directions for use**
- 6. Operation and transfers**
 - 6.1. Actions to be taken before and during each transfer**
 - 6.2. Using SoloVest**
 - 6.3. Transfer to sitting position using thigh straps**
 - 6.4. Transfer to standing position without thigh straps**
 - 6.5. Transfer from floor to wheelchair and back.**
- 7. Troubleshooting**
 - 7.1. Checks**
 - 7.2. Current repair**
- 8. Maintenance conditions and repair conditions. Conditions of storage, operation and transportation.**
 - 8.1. Maintenance conditions and repair conditions.**
 - 8.2. Conditions of storage, operation and transportation.**
- 9. Cleaning**
- 10. Utilization**
- 11. Guarantee certificates / Service life**

1. Warnings and important information

WARNINGS!!!

- Thoroughly read and understand the information in this product manual before attempting to use this product. If the procedures and instructions in this manual are not followed, serious injury could occur.
- The use of this product is subject to approval from qualified medical personnel and must be supervised by specialist at all times.
- This product is intended for use by children and adults with cognitive and corporal defects.
- This product is intended for indoor use only. It must be used on flat level surface only. This product must not be used in and around water, near staircases, swimming-pools, etc.
- Straps and supports are provided for the safety of the user and must be carefully adjusted for comfort and security. The use of straps and supports must be supervised by medical personnel at all times.
- Do not leave the patient in the elevator unsupervised.
- Slippery clothing, patient's low muscular tone, patient's specific body structure may cause the SoloVest to slide up on the patient, making a safe transfer difficult. In this case immediately stop lifting.
- Use only the Rifton SoloVest with the Rifton SoloLift. Elevator Rifton SoloLift and SoloVest are incompatible with any other elevator or vest.
- Do not lift the patient higher than necessary for transfer.
- Do not use the elevator in vertically confined space, always ensure there is free space over patient's head before lifting.
- Move the origin and the destination of the transfer as close to each other as possible.
- Do not use the elevator if battery status indicator is red.
- Before lifting keep your hands and feet away from the Lifting column.

- Before extending or retracting the base frame keep your hands and feet away at the safe distance.
- Transmission devices (e.g. cell phones) may create disturbances and should not be used during the use of elevator.

IMPORTANT: Please save this product manual for future use.

2. Prescription

Elevator Rifton SoloLift is a mobile transfer device. Zero lifting enables a single caregiver to transfer a patient weighing up to 350 lbs. (160 kg) to and from a wheelchair, bed, chair, toilet or gait trainer. Elevator Rifton SoloLift allows transferring the patient to standing position, rendering support when standing or moving.

Using elevator Rifton SoloLift together with SoloVest to be put above the waist creates opportunities for easier transfer of patients.

3. Product information

3.1. Regulatory

Mobile elevator Rifton SoloLift with accessories for vertical adjustment and movement of patients conforms to the EU Medical Device Directive 93/42/EEC. This device conforms to the requirements of EN 12182.

The Quality management system of the Company Community Products, LLC, USA conforms to ISO 13485:2003.

Mobile elevator Rifton SoloLift with accessories for vertical adjustment and movement of patients is produced according to the Quality management system of Company Community Products, LLC, USA; Risk Analysis Corporate standard RS-107; Failure Mode and Effects Analysis. The Risk Analysis program of the Company Community Products, LLC, USA conforms to ISO 14971.

3.2. Design of mobile elevator Rifton SoloLift



Figure 1. Dimensions of mobile elevator Rifton SoloLift

3.3. Technical data of mobile elevator Rifton SoloLift

Patient's max. height	193 cm
Minimum distance from the floor to armpit	48 cm
Waist circumference	56 to 152 cm
Overall length	129,2 cm
Overall width	78,7 cm
Width with legs extended	118,1 cm
Overall height with support handles	96,7 to 171,8 cm
Distance from the floor to the top of base frame (for passage under furniture)	14,6 cm
Turning diameter	142 cm
Overall weight of the device	64 kg
Patient's max. weight	160 kg
Lifting/lowering speed (no load)	5 cm/s
Battery	Detachable, sealed lead-acid gel-type batteries 24 V, 2,9 Ah
Battery charger	Wall-mounted charger, 100-240 V, max 650 mA
Motor	24V, 10A, permanent magnet motor
Emergency lowering	Mechanical or electrical
Wheels	Front: 100mm dual. Rear: 100mm dual with brake
Duty cycle	Two minutes continuous use followed by 18 minutes idle
Protection degree	IPX6

3.4. Accessories

Frame of mobile elevator Rifton SoloLift is a metal structure with electrical lifting device, intended for fixing the vest. Frame of mobile elevator Rifton SoloLift consists of base frame, Lifting column, lift arm, swing arm and rotation lock. Frame of mobile elevator Rifton SoloLift may be optionally fitted with display. Built-in display shows the load on patient's body during transfer. Material –ASTM powder-coated steel, A500 Grade B, latex, lead and phthalates free.

Control box is built in the elevator frame and allows controlling the elevator: lower or lift its moving part as necessary. The control is brought out to hand pendant.

Battery – Detachable, sealed lead-acid gel-type battery 24 V, 2.9 Ah. A fully charged battery allows approximately 70 full duty lifting cycles. Service life of the battery is 5 years.

Additional battery - Detachable, sealed lead-acid gel-type battery 24 V, 2.9 Ah. This battery is required to use the elevator all day long, if the main battery is discharged. Service life of the battery is 5 years.

Battery charger is designed to charge the battery. This is a wall-mounted battery charger, 100-240 V, maximum 650 mA.

Additional battery charger is designed to charge the backup (additional) battery in parallel with charging the main battery.

User's manual is provided in hardcopy.

Thigh straps are fixed to SoloVest and designed to transfer the patient to the sitting position.

Wheel direction lock ensures straightforward movement of Rifton elevator. The wheel direction lock can be installed on any wheel. To ensure straightforward movement of Rifton elevator the locks should be installed at least on two wheels. If unlocked, the wheels rotate freely, which is necessary for maneuverability in narrow places.

Three sizes of the SoloVest are available: small, medium and big. SoloVest is intended for supporting the patient's body during transfer.

Parameters	SoloVest sizes		
	Small size	Medium size	Big size
Figures			
Color	green	blue	dark blue
Vest height (front / lateral), cm	24 / 19	30 / 23	30 / 25
Weight, kg	1,6	2	2,5
Girth, cm	56 - 86	71 - 102	91 - 152

The SoloVest adapter strap makes it possible to use the SoloVest with an overhead track system that uses a spreader bar. Overhead track system is a construction used in medical facilities and some social security facilities. SoloVest may be used with overhead track systems when the patients start to walk by themselves.

Wheels ensure mobility of elevator frame. Front wheels: 100 mm, double; rear wheels: 100 mm, double, with a brake. Material – nylon wheels, tires filled with polyurethane foam.

Network cable – to plug the battery charger into wall outlet.

Rings are used to fix the thigh straps to the patient handle. The rings are fixed to blue clips on patient handle.

Patient handle. During lifting in the elevator the patient can hold onto the handle. The patient handle is fixed to the elevator frame and is fitted with blue clips which are used to fix thigh straps.

Hand pendant is used to lift and lower the patient with the elevator.

Foot pedal is used to change the width of frame and fixed to the base frame.

4. Initial assembly

Mobile elevator Rifton SoloLift is supplied assembled, ready for use.

5. Directions for use


Rotation lock: Push the blue rotation lock lever up to unlock and rotate patient. Pull lever down to lock after rotation is complete.

Figure 2.

Emergency mechanical lowering: If electrical operation fails, lower the patient by pulling the red emergency-lowering regulator straight up. (This will only work when the SoloLift is loaded.)

Figure 3.

Hand pendant: Operate the lift by using the up and down buttons on the hand pendant


Figure 4

Electrical emergency lowering: If the hand pendant fails to operate, lower the patient by using a pointed object (such as a ball-point pen) to push the hole marked "EMERGENCY" on the control box.


Emergency stop: To stop the SoloLift in an emergency, push the red emergency stop button on the control box. Reset by turning the button clockwise.

Figure 5

Extending the legs: Step on the upper foot pedal to extend or retract the base frame.

Figure 6

Battery charger and charging the batteries

Installation

1. Attach the charger to the wall near an outlet, using two screws (A) (see Figure 7). Two screws are provided with the charger, however, they may not be suitable for every situation.
2. Place battery in battery charger (see Figure 8).
3. Plug the charger cord into wall outlet.

Figure 7



Figure 8

Charging

Charge batteries 24 hours before first use.

When the battery needs charging, an indicator light on the hand pendant will go on.

Charge batteries every night, even if indicator light is green.

To charge the battery, remove the battery from the control box, and place it onto the wall-mounted charger.

The battery charger and indicator light will shut off automatically when charging is complete.

Charging normally takes approximately six hours.

During long periods of inactivity or storage, batteries will lose charge.

Allowing batteries to deep cycle (become nearly or completely dead) will destroy them. To prevent this, store batteries in battery charger.

An additional battery charger is required to charge the additional (backup) battery.

WARNING: Do not charge batteries in a wet area.

Old batteries should be disposed of properly at an appropriate recycling facility.

6. Operation and transfers

6.1. Actions to be taken before and during each transfer

Before each transfer

- Thoroughly check SoloVest for presence of breaks, cuts or worn out parts, pay specific attention to straps and buckles. Remove the product from service if any defects have been identified.
- Ensure the SoloVest is positioned correctly and adjusted tightly.
- Lock wheels on patient equipment during transfer.

During each transfer


- At the beginning of each transfer, while the patient is still seated and the straps are snug, check that the SoloVest cannot slide up or cause discomfort.
- Lift the patient only high enough to perform the transfer.

WARNING! Supervision by medical personnel is required at all times.

6.2. Using SoloVest

Position the SoloVest well below the patient's arms and the body support pads of the vest should overleap in front. Fix the buckles and tight the straps.

Figure 9. Correct size of SoloVest. Front view.

CAUTION: Make sure that the SoloVest size is chosen correctly. The SoloVest is too big for the patient if the body support pads overlap too far and the straps cannot be tightened.

WARNING: Do not transfer the patient, if the vest slips upwards during transfer.

6.3. Transfer to sitting position using thigh straps

Position the SoloVest around the patient and pull the straps tight.

Roll the SoloLift up to the patient with the swing arm in the raised position (Figure 10). Extend base frame if necessary. Once the SoloLift is in position, rotate the swing arm down behind the patient.

Figure 10.


Attach the SoloVest clips through the holes on the SoloLift swing arm (Figure 11)

Figure 11

Attach the thigh straps by hooking the silver rings over the blue hooks on the patient handle (Figure 12). Adjust the length of the thigh strap as necessary.

Figure 12.

Once the patient is secure and comfortable, activate the lift by using the hand pendant, verifying that the SoloVest is secure and is not sliding up (Figure 13). If the SoloVest slides up, stop lifting immediately. Lift the patient only high enough to perform the transfer.

Figure 13.

6.4. Transfer to standing position without thigh straps.

Position the SoloVest around the patient and pull the straps until tight.

Roll the SoloLift up to the patient with the swing arm in the raised position (Figure 10). Extend base frame if necessary. Once the SoloLift is in position, rotate the swing arm down behind the patient.

Attach the SoloVest clips through the holes on the SoloLift swing arm (Figure 11).

Leave the thigh straps unattached to lift patient into a standing position.

Once the patient is secure and comfortable, activate the lift by using the hand pendant, verifying that the SoloVest is secure and is not sliding up (Figure 14). If the SoloVest slides up, stop lifting immediately. Lift the patient only high enough to perform the transfer.

Figure 14.


Push the blue rotation lock lever up to position the patient into a standing position. Push lever back down to lock after rotation is complete. (Figure 15). Push lever back down to lock after rotation is complete. (Figure 16).

Figure 15.

Figure 16.

Secure patient in the intended object of transfer before removing the SoloVest. (Figure 17).

Figure 17.

6.5. Transfer from floor to wheelchair and back.

Sit the patient up.

Position the SoloVest around the patient and pull the straps until tight.

Roll the SoloLift up to the patient with the swing arm in the raised position (Figure 10). Extend base frame if necessary. Once the SoloLift is in position, rotate the swing arm down behind the patient.

Attach the SoloVest clips through the holes on the SoloLift swing arm (Figure 11). Leave the thigh straps unattached to lift patient into a standing position.

Once the patient is secure and comfortable, activate the lift by using the hand pendant, verifying that the SoloVest is secure and is not sliding up (Figure 13). If the SoloVest slides up, stop lifting immediately. Lift the patient only high enough to perform the transfer.

As the patient rises, move the lift forward slowly so that the patient's feet remain in the same position during the entire lift to standing (Figure 18).

Figure 18.

Once standing, turn patient 180°, bring wheelchair in behind patient and lower patient into wheelchair (Figure 19).

Figure 19.

To transfer the patient from wheelchair to floor, repeat the above steps in the reverse order.

7. Troubleshooting

7.1. Checks

Daily check

- Visual damage or wear of elevator Rifton SoloLift and accessories.
- Cracks, breaks, missing or loose details and/or improper operation.
- Correct functioning of emergency stop.

Remove the product from service if any condition develops that might make use unsafe.

The product may only be repaired using the original parts.

7.2. Current repair

Lift does not go up or down:

1. Make sure that the red emergency stop button on the control box is not pushed.
Rotate the button clockwise to start again.
2. Make sure the electrical cables are secure.
3. Make sure the battery box lid is closed.
4. Make sure the battery is charged and inserted correctly.
5. Contact Community Products, LLC, USA or its authorised representative RehabMedical, LLC, Russia.

Battery does not charge:

1. Check that the battery contact plate is not damaged or broken.
2. Check that wall charger is plugged securely into the wall socket.

Display does not turn on:

1. Remove and reinsert batteries, making sure that each cell is correctly oriented in polarity.
2. Batteries may be dead and need replacing.
3. Contact Community Products, LLC, USA or its authorised representative RehabMedical, LLC, Russia.

Spare parts replacement

If necessary, the manufacturer will supply any spare parts and repair manuals for the product. The customer is responsible to make sure that the repairs or modifications are correctly and safely completed. To place the order for spare parts please identify the serial number (Item №) of your product. Contact us at the numbers listed in paragraph 11 of this manual. Use spare parts manufactured only by Community Products, LLC, USA.

8. Maintenance conditions and repair conditions. Conditions of storage, operation and transportation

8.1. Maintenance conditions and repair conditions

- Temperature from $+14^{\circ}\text{C}$ to $+40^{\circ}\text{C}$
- Relative humidity up to 90%
- Atmospheric pressure from 680 to 790 Mm hg.

8.2. Conditions of storage, operation and transportation

Operating conditions for medical device

- Temperature from $+14^{\circ}\text{C}$ to $+40^{\circ}\text{C}$
- Humidity up to 90%
- Atmospheric pressure from 680 to 790 Mm hg.

Storage conditions

Elevator Rifton SoloLift with accessories is stored in protected from direct sunshine dry place at air temperature from -15°C to $+45^{\circ}\text{C}$, humidity to 90%, atmospheric pressure from 680 to 790 Mm hg.

Conditions for transportation

Elevator Rifton SoloLift with accessories is transported in covered vehicles, at temperature from -15°C to $+45^{\circ}\text{C}$, at humidity to 90%.

9. Cleaning

Elevator Rifton SoloLift should be cleaned with a moist rag using soap solution, mild detergent or 10% bleach solution. Mild disinfectants may be used to clean the straps and body support pads. SoloVest may not be washed in washing machine or immersed in water.

10. Utilization

Elevator Rifton SoloLift and accessories thereto do not contain any substances dangerous to human life, health and environment. All materials used for packaging Elevator Rifton SoloLift and accessories thereto are fully reusable. Packaging materials should be disposed of in an environment-friendly manner. Elevator Rifton SoloLift and accessories thereto should be disposed of in accordance with the applicable law and local regulations and instructions. The run-down battery packs should be disposed of in accordance with clause 273 40 Code of Federal Regulations (CFR).

11. Guarantee certificates / Service life

Warranty: 1 year after delivery.

Subject to correct use the service life of the product is 5 years. After 5 years the service life of the product may be extended once the battery is replaced.

Should any questions, failures or operating problems arise, feel free to contact Community Products, LLC, 2032 Route 213, Rifton, NY 12471, USA, + 1 845-658-7700 / + 1 845-658-3144 or its authorised representative RehabMedical, LLC, 20a/1, Novoalekseevskaya st, Moscow, 129626, Russia, phone / fax: +7 (495) 686-40-47.

Community Products, LLC, USA and its authorised representative RehabMedical, LLC, Russia cannot accept responsibility for any modifications or alterations made to mobile elevators Rifton after they leave manufacturer's premises.

Community Products, LLC, USA and its authorised representative RehabMedical, LLC, Russia do not accept responsibility for use of elevator Rifton jointly with products of other manufacturers.

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