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
NOTARY PUBLIC

On this 25th day of August 2015 before me, the undersigned Notary Public, personally appeared Denise Oliveira, who proved to me through satisfactory evidence of identification (RI State License) to be the person whose name is signed on this document.

Keyi Sino

VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System

Low Profile Bioresorbable Coated Permanent Mesh with Positioning System
Designed for Laparoscopic Ventral Hernia Repair

Includes: Introducer Tool* Inflation Assembly

Instructions for Use

Product Description

VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System (the device) is a low profile, bioresorbable, coated, permanent mesh, with a pre-attached removable positioning system, designed for the reconstruction of soft tissue deficiencies during laparoscopic ventral hernia repair.

VENTRALIGHT™ ST Mesh is a dual-component (absorbable and nonabsorbable) sterile prosthesis designed for the reconstruction of soft tissue deficiencies. The low profile mesh facilitates laparoscopic deployment and the pre-sized shapes offer maximum ready-to-use benefits. VENTRALIGHT™ ST Mesh is co-knitted using polypropylene (PP) and polyglycolic acid (PGA) fibers to result in a two-sided mesh with a PP surface and a PGA surface. The mesh is coated on the PGA surface with a bioresorbable, chemically modified sodium hyaluronate (HA). The carboxymethylcellulose (CMC) and polyethylene glycol (PEG) based hydrogel. The fascial side of the mesh allows for a prompt fibroblastic response through the interstices of the mesh, allowing for complete tissue ingrowth, similar to polypropylene mesh alone. The visceral side of the mesh is a bioresorbable coating that separates the mesh from underlying tissue and visceral organ surfaces to minimize tissue attachment to the mesh. Shortly after hydration, the biopolymer coating becomes a hydrated gel that is resorbed from the site in less than 30 days.

The ECHO PS™ Positioning System is made up of a thermoplastic coated nylon balloon with an attached inflation tube which is pre-attached to the VENTRALIGHT™ ST Mesh with small connectors. Once inserted, the balloon is inflated with the included inflation assembly to help facilitate laparoscopic deployment, including unrolling, positioning, and placement of the mesh. Once initial mesh fixation is complete, the ECHO PS™ Positioning System is completely removed from the body. (See Section 6, Directions for Using VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System.)

* The sizes of VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System that are packaged with an Introducer Tool are indicated in the chart below. The Introducer Tool consists of a metal tines/T-cap assembly which is used to roll VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System in order to facilitate laparoscopic introduction (see Section 4, Directions for Using the Introducer Tool). Products packaged without an introducer tool are designed to be introduced with a grasper (see Section 5, Directions for Use with Graspers).

Product Code (REF)	Mesh Size	Description	Introducer Tool Included	Recommended Minimum Trocar/Incision Size*
5955450	4.5" x 11.4 cm	Circle	No	10 mm
5955460	4" x 6" / 10.2 cm x 15.2 cm	Ellipse	No	10 mm
5955600	6" / 15.2 cm	Circle	Yes	12 mm
5955680	6" x 8" / 15.2 cm x 20.3 cm	Ellipse	Yes	12 mm
5955610	6" x 10" / 15.2 cm x 25.4 cm	Oval	Yes	12 mm
5955790	7" x 9" / 17.8 cm x 22.9 cm	Ellipse	Yes	12 mm
5955800	8" / 20.3 cm	Circle	Yes	12 mm
5955810	8" x 10" / 20.3 cm x 25.4 cm	Ellipse	Yes	15 mm
5955113	10" x 13" / 25.4 cm x 33.0 cm	Ellipse	Yes	12 mm trocar incision site
5955124	12" x 14" / 30.5 cm x 35.6 cm	Rectangle	No	12 mm trocar incision site
5955000	Reorder Code for ECHO PS™ Positioning System Inflation Assembly (EICS)			

*If a proximal cap is available on the trocar, removing the proximal cap can help facilitate deployment. Deployment capability may vary depending on rolled patch size and graspers/trocars used.

Indications

VENTRALIGHT™ ST Mesh is indicated for use in the reconstruction of soft tissue deficiencies, such as for the repair of hernias.

The ECHO PS™ Positioning System is intended to be used to facilitate the delivery of soft tissue prostheses during laparoscopic hernia repair.

Contraindications

- Do not use the VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System in infants or children whereby future growth will be compromised by use of such material.
- Do not use VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System for the reconstruction of cardiovascular defects.
- Literature reports there is a possibility for adhesion formation when the polypropylene is placed in direct contact with the bowel or viscera.

Warnings

- This device is supplied sterile. Inspect the packaging to be sure it is intact and undamaged prior to use.
- This device has been designed for single use only. Reuse, resterilization, reprocessing and/or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device and may lead to device failure which may result in injury to the patient. Reuse, reprocessing, resterilization, or repackaging may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient or end user.

- If unused prosthesis has been in contact with instruments or supplies used on a patient or contaminated with body fluids, discard with care to prevent risk of transmission of viral infections.
- Ensure proper orientation; the coated side of VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System should be oriented against the bowel or sensitive organs. Do not place the polypropylene side against the bowel. There is a possibility for adhesion formation when the polypropylene side is placed in direct contact with the bowel or viscera (reference Surface Orientation section below).
- The use of any permanent mesh or patch in a contaminated or infected wound could lead to infection, fistula formation and/or extrusion of the prosthesis.
- If an infection develops, treat the infection aggressively. Consideration should be given regarding the need to remove the prosthesis. An unresolved infection may require removal of the prosthesis.
- To prevent recurrences when repairing hernias, it is recommended that the prosthesis be large enough to extend at least 3 to 5 cm beyond the margins of the defect.
- Do not apply sharp, heat emitting, or ultrasonic tools (such as scissors, needles, tackers, diathermic tools, etc.) to the ECHO PS™ Positioning System.
- The ECHO PS™ Positioning System should not be used with any other hernia prosthesis aside from those with which it comes pre-attached/packaged.
- VENTRALIGHT™ ST Mesh is the only permanent implant component of the device. The inflation adapter and syringe are to be kept external to the patient and discarded after use. The ECHO PS™ Positioning System (including the balloon, all connectors, and inflation tube) is to be removed from the patient and appropriately discarded as it is not part of the permanent implant.
- Discard Introducer Tool and all components of the ECHO PS™ Positioning System (including the inflation adapter and syringe) after use. This product may be a potential biohazard. Handle and dispose in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Precautions

- Please read all instructions prior to use.
- Only physicians qualified in the appropriate surgical techniques should use this device.
- The safety and effectiveness of VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System has not been evaluated in clinical studies for the presence of malignancies in the abdominopelvic cavity.
- Visualization should be maintained throughout the course of the entire procedure. Additionally, laparoscopic removal of the ECHO PS™ Positioning System must be performed under sufficient visualization of the entire device and surrounding anatomy to ensure proper removal.
- Do not trim the mesh. This will affect the interface between the mesh and positioning system.

INSTRUCTIONS FOR USE

Note on Transfascial Sutures: The ECHO PS™ Positioning System eliminates the need for transfascial orientation sutures. If transfascial sutures are to be used for fixation, place sutures after removing the ECHO PS™ Positioning System from the body and after completing all mechanical fixation.

1) Preparation

It is recommended that VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System be completely immersed in sterile saline for no more than 1–3 seconds immediately prior to placement in order to maximize the flexibility of the prosthesis. The safety and effectiveness of VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System in combination with solutions other than saline have not been tested.

2) Surface Orientation

It is extremely important that this product be oriented correctly to function as intended. The bioresorbable coated side of the VENTRALIGHT™ ST Mesh (which contains the ECHO PS™ Positioning System) must always be positioned against those surfaces where minimal tissue attachment is desired, i.e., against bowel or other visceral structures. The bioresorbable coated surface is designed to temporarily separate tissue surfaces and minimize tissue attachment to the mesh. The uncoated polypropylene mesh side should face the surface where tissue ingrowth is desired, and never be placed against the bowel or other visceral structures.


3) Trocar Size Recommendations

VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System should be hydrated for no more than 1–3 seconds just prior to laparoscopic placement. The VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System must be rolled immediately after hydration. Using the Introducer Tool (See Section 4, Directions for Using the Introducer Tool) or grasper (See Section 5 Directions for Use with Graspers) roll the VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System parallel to the BAW® logo with the polypropylene on the outside and bioresorbable coated side with the ECHO PS™ Positioning System on the inside. Insert into the abdomen through the recommended minimum trocar or trocar incision site. Do not force the device through the trocar. If the VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System will not easily deploy down the trocar, remove the trocar and insert through the next largest available size trocar or through the trocar incision site and reinsert trocar. See table in Product Description Section for recommended minimum trocar size.


4) Directions for Using the Introducer Tool

Note: If using REF 5955450, REF 5955460, or REF 5955124 see Section 5, Directions for Use with Graspers.

- Remove the VENTRALIGHT™ ST Mesh with ECHO PS™ Positioning System (the device) and included Introducer Tool from the sterile pouch. Set inflation assembly pouch aside.
- Hydrate the device for no more than 1–3 seconds just prior to rolling (Fig. 1).
- Place the device between the metal tines parallel to the BAW® logo approximately 2.5 cm to 5 cm from the edge of the device. Ensure that the device is centered on the tines and that the inflation tube is lying flat and facing the proximal end of the tool in parallel to the rolling tines (Fig. 2).
- Place T-cap on the end of the tines (Fig. 3).
- Grasp the center of the device/tines to provide counter pressure against the device/tines. With one hand roll the device, polypropylene side out and bioresorbable coated side with ECHO PS™ Positioning System on the inside, by turning the tines handle until the device is completely wrapped around the tines. Ensure the inflation tube is not wrapped around the mesh (Fig. 4).
- Remove the T-cap. Ensure the device is positioned such that at least ¼ cm of the tines extends beyond the mesh edge (Fig. 5).
- Deliver the device through the trocar under sufficient visualization of the device and the surrounding anatomy.

- As the device is being deployed through the trocar, rotate the tines handle in the direction that the device was rolled. This will keep the device tight around the tines, thus facilitating deployment (Fig. 6).
- Rotate the tines handle approximately $\frac{1}{2}$ turn in the opposite direction the device was rolled and partially slide the tines handle out of the device. Do not completely remove the tines from the device until the device has passed through the trocar entirely (Fig. 7, Fig. 8).
- Note: If the tines are removed from the device before it is completely deployed through the trocar, use the laparoscope to push the mesh through the trocar, or use a grasper from an opposing trocar location to pull the mesh into the abdomen.
- Under visualization, continue to advance the device and tines through the trocar. Repeat the previous steps and rotate the tines handle in the direction the device was rolled to completely deploy the patch through the trocar (Fig. 9).
- After the device has cleared the trocar, remove the tines from the trocar and discard appropriately.

5) Directions for Use with Graspers

Note: For use with REF 5955450, REF 5955460, and REF 5955124 only.

- Remove the VENTRAUGHT™ ST Mesh with Echo PS™ Positioning System (the device) from the sterile pouch. Set inflation assembly pouch aside.
- Hydrate the device for no more than 1–3 seconds just prior to rolling (Fig. 1).
- Roll the device tightly on itself starting at the edge parallel to the BARD® logo and continue across with the polypropylene side on the outside and bioresorbable coated side with the Echo PS™ Positioning System on the inside (Fig. 10).
- While holding the rolled device, grasp the leading edge with the grasper as shown (Fig. 11a). For optimal results, it is important to ensure that the balloon is protected and does not come in direct contact with the trocar or trocar incision site. Ensure that the inflation tube is lying flat and facing the proximal end of the grasper. Place through the trocar or trocar incision site (Fig. 11b).


Note: If the grasper slips from the device before the device is completely deployed through the trocar incision site, use a grasper from an opposing trocar location to pull the mesh into the abdomen.


6) Directions for Using VENTRAUGHT™ ST Mesh with Echo PS™ Positioning System

To activate and inflate the device:


- Once the VENTRAUGHT™ ST Mesh with Echo PS™ Positioning System (the device) is inserted into the abdomen, use a grasper to locate the blue retrieval loop on the inflation tube, ensuring that the blue inflation tube is not wrapped around the mesh and is clearly visible (Fig. 12).
- Pass a suture passer device through healthy skin at the center of the hernia defect (avoid going directly through the umbilicus). Grasp the blue retrieval loop and pull the retrieval loop and inflation tube out of the abdominal cavity (Fig. 13).

- Place an atraumatic clamp or hemostat on the inflation tube at the level of the skin to temporarily hold the device in place (Fig. 14a). Cut the inflation tube on the dashed line between the retrieval loop and the yellow anchor with surgical scissors (Fig. 14b) to open the tube for inflation. Discard the retrieval loop (Fig. 14c).
- Remove the inflation assembly from the sterile pouch and screw the inflation adapter tightly to the syringe (Fig. 15). Connect the inflation tube and inflation assembly as follows (Fig. 16):
 - Ensure the clear cap of the adapter is pushed downward to open the inflation tube channel.
 - Insert the inflation tube as far as possible into the cap opening.
 - Pull the clear cap of the adapter upward to lock in place.


- To inflate the device, release the clamp or hemostat and pull the inflation tube upward to lift the mesh off the viscera (Fig. 17). Always grab and pull the tube directly. Do not lift up using the inflation adapter/assembly. Inflate the device by pumping the syringe until the balloon is fully inflated (Fig. 18). The syringe may be filled while attached to the inflation assembly. One to three pumps will be required to fully inflate the device depending on the size of the mesh. A slight high pitch sound may occur; this is normal and indicates the inflation assembly is working properly. Note: The mesh should now be unrolled, opened and ready to be properly positioned.


To properly position and fixate the mesh:

- Raise the inflation tube to properly adjust the mesh to the desired position then clamp the inflation tube to hold the device in place (Fig. 19). Once the device is inflated, if desired, remove the syringe by unscrewing from the adapter and setting aside (Fig. 20).
- If necessary, use a grasper to orient the VENTRAUGHT™ ST Mesh in relation to the defect (Fig. 21). On oval and ellipse sizes, the BARD® logo and dark shaded arrows on the device represent the long axis of the mesh. The two marked areas at both ends of the long axis indicate the midline of the VENTRAUGHT™ ST Mesh.
- Once the device has been properly positioned and before the Echo PS™ Positioning System is deflated, ensure that no tissue is entrapped between the device and abdominal wall. Fixate around the entire patch perimeter of VENTRAUGHT™ ST Mesh with fasteners placed 1–2 cm apart (Fig. 22). Ensure that no fasteners are placed through the Echo PS™ Positioning System (Fig. 27).


To deflate and remove the Echo PS™ Positioning System:

- To deflate the Echo PS™ Positioning System, release the clamp on the inflation tube, cut the tube as close to the skin as possible, and then discard the excess tube (Fig. 23).
- Begin removal of the Echo PS™ Positioning System by grasping one of the two removal points marked by the dark arrows adjacent to the BARD® logo (Fig. 24). Begin pulling the positioning system off the mesh in one smooth motion.
- Continue removing the Echo PS™ Positioning System by pulling it up to the tip of the trocar (Fig. 25). Remove both the Echo PS™ Positioning System and trocar simultaneously (Fig. 26).

Verify that the Echo PS™ Positioning System including the balloon, all mesh connectors, and the inflation tube (as illustrated below) is fully intact after removal. The number of mesh connectors per balloon can be found in the adjacent table. Once the positioning system is verified to be fully intact, discard the Echo PS™ Positioning System appropriately.


Product Code (REF)	Mesh Size	Number of Mesh Connectors
5955450	4.5" / 11.4 cm	2
5955460	4" x 6" / 10.2 cm x 15.2 cm	2
5955600	6" / 15.2 cm	4
5955680	6" x 8" / 15.2 cm x 20.3 cm	4
5955610	6" x 10" / 15.2 cm x 25.4 cm	4
5955790	7" x 9" / 17.8 cm x 22.9 cm	4
5955800	8" / 20.3 cm	4
5955610	8" x 10" / 20.3 cm x 25.4 cm	4
5955113	10" x 13" / 25.4 cm x 33.0 cm	8
5955124	12" x 14" / 30.5 cm x 35.6 cm	8

- Reinsert the trocar, place any additional fixation as desired, and complete the procedure. Care should be taken to ensure that the patch is adequately fixated to the abdominal wall. If necessary, additional fasteners and/or sutures should be placed.

Fixation

BARD® fixation devices or nonabsorbable monofilament sutures are recommended to properly secure VENTRAUGHT™ ST Mesh. If other fixation devices are used, they must be indicated for use in hernia repair. If transfascial sutures are to be used for fixation, place sutures after removing the Echo PS™ Positioning System from the body and after completing all mechanical fixation. Care should be taken to ensure that the patch is adequately fixated to the abdominal wall. If necessary, additional fasteners and/or sutures should be placed.

Fig. 27
Reference Step #8 above.

Adverse Reactions

Possible complications include seroma, adhesions, hematomas, inflammation, extrusion, fistula formation, infection, allergic reaction, and recurrence of the hernia or soft tissue defect.

Traceability

Traceability labels which identify the type, size, and lot number of the device are attached to every package. A label should be affixed to the patient's permanent medical record to clearly identify the device which was implanted.

If you experience a product failure, please contact Davol Inc. at 1-800-556-6275 for instructions on returning the product.

Storage

VENTRAUGHT™ ST Mesh with Echo PS™ Positioning System should be stored in a clean, dry area at ambient room temperature until ready for use.

Bard, Davol, Echo PS, and Ventrailight are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 2010 - 2013 C. R. Bard, Inc. All Rights Reserved.

	ENGLISH	FRANÇAIS	DEUTSCH	ITALIANO	ESPAÑOL
	Contents	Contenu	Inhalt	Contenuto	Contenido
	Do not use if package is damaged	Ne pas utiliser si l'emballage est endommagé	Nicht verwenden, wenn die Verpackung beschädigt ist	Non utilizzare se la confezione è danneggiata	No usar si el envase está dañado
	Actual Size	Taille réelle	Tatsächliche Größe	Dimensioni reali	Tamaño real
	Partially Absorbable	Partiellement résorbable	Teilweise resorbierbar	Parzialmente assorbibile	Parcialmente reabsorbible
	Circle	Cercle	Kreisförmig	Circolo	Círculo
	Ellipse	Ellipse	Ellipse	Ellissi	Elipse
	Oval	Ovale	Oval	Ovale	Óvalo
	Rectangle	Rectangle	Rechteck	Rettangolo	Rectángulo
	Patents Pending	Brevets en instance	Patente angemeldet	Brevetti in corso di registrazione	Patentes pendientes
	Package 1 of 2	Colis 1 sur 2	Packung 1 von 2	Confezione 1 di 2	Envase 1 de 2
	Package 2 of 2	Colis 2 sur 2	Packung 2 von 2	Confezione 2 di 2	Envase 2 de 2
	U. S. Federal law restricts this device to sale by or on the order of a physician.	Conformément à la loi fédérale américaine, ce dispositif ne peut être vendu que par un médecin ou sur prescription médicale.	Nach dem US-Bundesgesetz darf die Abgabe dieses Produkts nur durch einen Arzt oder auf ärztliche Anordnung erfolgen.	Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo ai soli medici o su prescrizione degli stessi.	La legislación federal estadounidense restringe este dispositivo a su venta por un médico o bajo prescripción médica.
	For use with	À utiliser avec	Zur Verwendung mit	Per uso con	Para utilizar con

	NEDERLANDS	PORTUGUÊS	ΕΛΛΗΝΙΚΑ	DANSK	SVENSKA
	Inhoud	Conteúdo	Περιεχόμενα	Indhold	Innehåll
	Het product niet gebruiken als de verpakking beschadigd is	Não utilizar se a embalagem estiver danificada	Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία έχει υποστεί ζημιά	Må ikke bruges, hvis emballagen er beskadiget	Använd inte om förpackningen är skadad
	Ware grootte	Tamanho real	Πραγματικό μέγεθος	Faktisk størrelse	Verklig storlek
	Gedeeltelijk resorbeerbaar	Parcialmente absorvível	Μερικώς απορροφήσιμο	Delvist resorberbar	Delvis absorberbar
	Cirkel	Círculo	Κυκλικό	Cirkel	Cirkel
	Ellips	Elipse	Ελλειψοειδές σχήμα	Ellipse	Ellips
	Ovaal	Oval	Οβάλ	Oval	Oval
	Rechthoek	Retângulo	Ορθογώνιο σχήμα	Rektangel	Rektangel
	Octrooien aangevraagd	Patentes pendentes	Αναμένεται η έγκριση διπλωμάτων ευρεσιτεχνίας	Patentanmeldt	Patentsökt
	Verpakking 1 van 2	Embalagem 1 de 2	Συσκευασία 1 από 2	Pakke 1 af 2	Förpackning 1 av 2
	Verpakking 2 van 2	Embalagem 2 de 2	Συσκευασία 2 από 2	Pakke 2 af 2	Förpackning 2 av 2
	Volgens de Amerikaanse wetgeving mag dit product uitsluitend worden verkocht door of op voorschrift van een arts.	A legislação federal dos EUA restringe a venda deste dispositivo a um médico ou por encomenda de um médico.	Η ομοσπονδιακή νομοθεσία των Η.Π.Α. περιορίζει την πώληση της συσκευής αυτής μόνο σε ιατρούς ή κατόπιν εντολής ιατρού.	I henhold til amerikansk lovgivning (USA) må dette udstyr kun sælges eller ordineres af en læge.	Enligt amerikansk federal lag får detta instrument endast säljas av läkare eller på läkares ordination.
	Voor gebruik samen met	Para utilização com	Για χρήση με	Til brug sammen med	Används tillsammans med

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 **EC REP**
Bard Limited
Crawley, UK
RH11, 9BP

CE
0086


Medical Services & Support
Clinical Information Line
1-800-562-0027

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
Denise Oliveira
Regulatory Affairs Manager
Davol Inc., Subsidiary of C.R. Bard, Inc.

NOTARY PUBLIC

On this 22nd day of December 2014 before me, the undersigned Notary Public, personally appeared Denise Oliveira, who proved to me through satisfactory evidence of identification (Rhode Island State License) to be the person whose name is signed on this document.

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