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Instructions for use for the medical device

“Apex locator device AirPex”

Name of the medical device: "Apex locator device AirPex"

Manufacturer of the medical device: Changzhou Sifary Medical Technology Co., Ltd (China).

Address: No. 99, Qingyang Road, Xuejia County, Xinbei District, Changzhou City, 213000 Jiangsu, China.

Phone number: +86-0519-85962691

Email: ivy@sifary.com

Medical device purpose and indication:

Medical device purpose: Location of the root canal apex during the endodontic treatment.

Medical device indication: endodontic treatment (root canal treatment).

Indication: dentistry. This medical device should be used only in specialized medical institutions by qualified medical personnel.

Possible side effects and risks:

Possible side effects: not found

Risks: The main cumulative residual risk: non-compliance with the operation instructions described in the operating manual. The main way to manage residual risks in order to reduce them to a permissible level: It is necessary to follow the instructions described in the Operation Manual. Before each use the endodontic device must be switched on and checked and tested to guarantee its uninterrupted operation.

Parts Identification:

Charge base

Apex locator

File hook (2 pcs)


Touch probe

File clip (2 pcs)

Extension cord

Adapter

Tester.

APX Locator(1pcs) 	Charge Base (1pcs) 	Adapter (1pcs) 
File Hook (2pcs) 	Touch Probe (1pcs) 	File Clip (2pcs) 
Extension Cord (2pcs) 	Tester (1pcs) 	

Symbols used in the User Manual and marking

	If the instructions are not followed properly, operation may lead to hazards for the product or the user/patient.
	Additional information, explanation of operation and performance.
	Serial number
	Catalogue number
	Manufacturer
	Date of manufacture
	Class II equipment
	Type BF applied part
	CE marking
	Direct current
	Dispose of in accordance with the WEEE directive
	Keep dry
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified
	Authorized Representative in the European Community
	Temperature limitation
	Humidity limitation
	Manufacturer's LOGO
	Follow instructions for use
	Shell protection

7. Contraindications

Do not use this unit in conjunction with an electric scalpel or on patients who have a pacemaker.

Blocked canals cannot be accurately measured.

Safety and effectiveness have not been established in pregnant women and children.

WARNING

Read the following warnings before use:

The device must not be placed in humid surroundings or anywhere where it can come into contact with any type of liquids.

Do not expose the device to direct or indirect heat sources. The device must be operated and stored in a safe environment.

The device requires special precautions with regard to electromagnetic compatibility (EMC) and must be installed and operated in strict compliance with the EMC information. In particular, do not use the device in the vicinity of fluorescent lamps, radio transmitters, remote controls and do not use this system near the active HF Surgical equipment in the hospital. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the AirPex, including cables supplied by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Do not use, operate or store at high temperatures. Comply with the specified operating and storage conditions.

Gloves and a rubber dam are compulsory during treatment.

If irregularities occur in the device during treatment, switch it off. Contact the agency.

No modification of this device is allowed. Never open or repair the device yourself, otherwise, void the warranty.

Installing the AirPex

Connecting file clip, lip hook, extension cord and clip

Connect file clip, extension cord and lip hook to APEX locator as shown in the picture. Also, connect both extension cords according to actual situation. Then place APEX locator in clip for fixation.

Use of touch probe

The equipped touch probe can replace the use of the clip.

When using the touch probe, connect the lip hook, touch probe, extension cord to APEX Locator as shown in the below picture.

NOTE

Please use the original file clip and lip hook that manufactured by Sifary. Because the size of the unoriginal file clip and lip hook are different, it may damage the APEX locator or cause deviation of measurement accuracy.

Please check the connection of the device before use to make sure the device functions well.


Connecting charge base

Plug the USB of adapter into the charge base, plug the other end into a power outlet.

Put the APEX locator in the groove in the middle of the charge base. The power LED on charge base will light up. And the interface of APEX locator will light up to show that is charging.

Power LED


NOTE

Put the APEX locator on the charge base in the right place, otherwise the APEX locator will not be charged.

NOTE

Only the original adapter can be used. Do not use the device while charging. The APEX Locator power connector can only be used to connect the original adapter cord for charging purpose.

Use Interface

Set key ⑤
 Power key ⑥
 Reference range
 Apical area display
 Reference point
 Measuring bar
 Measured value
 Volume display
 Power display

Turn Power On/Off

Press ⑥ to turn on. Long press ⑥ more than 2 seconds, or no operation for 3 minutes to turn off.

Volume control

During standby state, shot press ⑥ to cycle the volume through the minor to the maximum.

Setting the reference point

During standby state, press ⑤ to set the reference point between 0~1. Seven points can be selected circularly. The cursor flashing position indicates the selected reference point.

Power display

Display the remaining power through the number of grid.

Reference range

The flash of the measuring bar is the current measured value, represents the estimated distance from the apical foramen in millimeters.

Display reversing



During standby state, press ⑥ and ⑤ together to reverse the display.


Operation and functional characteristics

1 Charge


	Display the present remaining amount of the battery. Less than 15% remains, please charge.
	NOTE If the power is less than 15%, the device must be recharged within 30 days, otherwise the battery will be damaged.
	Charging indication appears on the screen, and flashes slowly, when battery is fully charged or in a state near full charge, the flash will stop. It takes about 4-5 hours for full charge, depending on residual battery power and battery state. It can be recharged 300-500 times, depending on the operating conditions of the device. WARNING Do not change the battery. Only trained technicians or distributors can change the battery. The electronic parts will be damaged if use a wrong battery or installed with a wrong way.

2 Function checking of APEX locator

	- After turning on, insert the tester into the APEX locator. - Clamp the groove of tester with file clip or touch the groove of tester with touch probe. - The measuring bar flashes at the point when the displayed measured value is 02, 03 or 04 - Recommend to test the APEX locator with tester once a week. NOTE If the measurements are not expected, check whether the tester is connected properly. If the connection is normal but the screen still does not show the expected value, please stop using the device and contact the local dealer for processing.
	- Before each use, make the file clip touch the lip hook, or use the touch probe to touch the lip hook to confirm the device condition (short circuit). - Make sure that the tester is not installed on the APEX locator. Then connect the file clip, lip hook and extension cord according to chapter 4.1. Finally, make the lip hook touch the metal part of file clip. The measured value on the screen should be shown as -2. WARNING If the measurement shown is not -2, check if the connection is normal. If the connection is normal but the screen still does not show the expected value, please stop using the device and contact your local dealer.

Operation and not suitable condition


- Press the back cover of the file clip to make the hook of the file clip stick out. And hook the metal handle of the root canal file. Release the pressure and use the elasticity of the file clip to complete the connection between the file clip and the dental file.


NOTE

- When connecting the root canal file, make sure that the file clip and the root canal file handle are basically perpendicular, otherwise the chuck of the file clip will be easily damaged.
- This equipment does not include the root canal files. Please select suitable root canal files according to the clinical needs. The metal part of the root canal file should be well conductive.

- When the file clip can't enter the patient's mouth, the file clip can be replaced by the extension cord with the touch probe. Press the touch probe on the metal handle of the root canal file to complete the connection between the touch probe and the root canal file.

- Hook the lip hook to the patient's lip. Make sure that the lip hook contacts with the lip fully. Then slowly insert the root canal file into the prepared root canal.
- If the patient is fitted with a metal crown or other conductive device, the root canal file and the metal part of file clip should not be in contact with it, to avoid causing wrong measurement results.
- APEX locator should be fixed in patient's collar with the clip.

NOTE

- To prevent measurement errors caused by conduction between the gums or adjacent root canals, dry the pulp chamber floor with a cotton pellet or other means before testing.
- Use the root canal file with correct number and taper. Make the file fully contact with the canal wall, which facilitates accurate measurements.

Fig.1



Fig.2

- As the file progresses in the root canal, subsequent measuring bar lights up gradually.
- When displayed as shown in figure 1, it indicates that the distance from the apical foramen is about 2mm at this time. Meanwhile, the APEX locator makes a "di-di" alarm sound at a long interval.
- When displayed as shown in figure 2, it means that the distance from the apical foramen is close, about 1.5 mm, and the time interval of "di-di~" alarm sound made by the APEX locator becomes shorter.
- When displayed as shown in figure 3, it means that it reaches the apical foramen. The measured value is 00. Meanwhile, the Apex Locator emits a long beep sound with no interval.
- When displayed as shown in figure 4, it means that the root canal file has penetrated the apical foramen, and the main unit makes a very urgent "di~~" alarm sound.

NOTE

- The position of apical foramen (measured value: 00) measured by this device is the Major/Anatomic apical foramen. In clinical practice, in order to prevent surgical failure caused by piercing the root apical orifice, 0.5-1.0 mm is usually taken from the measured value, which is the Minor/Physiological apical foramen prepared for root canal.
- The value of the reference distance is only an estimated value, not a clinical basis.

Fig.3



Fig.4

- The measured value does not represent the distance. It simply indicates the file progression towards the apex.
- ⚠ **WARNING**
- During measuring, insert the file slowly to prevent penetrating the apical foramen.
- The APEX locator is used to detect the apex of root canal. In clinical use, it must be combined with X-ray and other means to determine the working length of the root canal.
- The device should be used by dentists with knowledge of dental root canal length and skill in operation.

Unsuitable situation of root canals for Electric Measurement

cannot obtain precise measurements if the root canal conditions as below

Root canal with a large apical foramen

The root canal cannot be accurately measured because of the lesion or incomplete development of the apical foramen. The results may show that the length measured is shorter than the actual one.


Root canal blood overflow from the opening

If blood spills from the root opening and contacts the gums, it will cause leakage of electricity, which cannot be accurately measured. Wait for the bleeding to stop completely. Clean the root canal and the opening, completely empty the root canal blood, and then measure it.

The root canal uses a chemical solution to flow out from the opening


If a chemical solution flows out of the root canal, it is impossible to get an accurate measurement.

It is important to remove the overflow from the opening.

Broken crown

If the crown is broken, a segment of the gingival tissue enters the lumen, and the contact between the gingival tissue and the root file causes electrical leakage, which cannot be accurately measured. In this case, the appropriate material should be used to isolate the gingival tissue.

The crack tooth Leakage through branch of the root canal


Broken teeth can cause electrical leakage and cannot be accurately measured.

Branch tubes can also cause leakage.






Retreatment canal which was filled with gutta-percha

The gutta-percha must be completely removed to eliminate its insulation, then pass a small file all the way through the apical foramen and then put a little saline in the canal, but do not let it overflow the canal opening.

Crown or metal prosthesis that touches gingival tissue

Accurate measurement cannot be obtained if the file touches a metal prosthesis that is touching gingival tissue. In this case, widen the opening at the top of the crown so that the file will not touch the metal prosthesis before taking a measurement.

	Cutting debris on tooth Pulp inside canal Remove all cutting debris on the tooth. Remove all the pulp inside the canal. Otherwise an accurate measurement cannot be obtained.
	Caries touching the gums In this case, electrical leakage through the caries infected area to the gums are impossible to obtain an accurate measurement.
	Blocked canal The meter will not run if the canal is blocked. Opening the canal all the way to the apical construction to measure it.
	Extremely dry canal If the canal is extremely dry, the meter may not work until it is quite close to the apex. In this case, try to moisten the canal with oxydol or saline.
Difference measuring result between Apex locator reading and Radiography Sometimes the reading of the apex locator reading does not correspond to the X-ray image. this does mean inaccurate of apex locator or X-ray, depending on the angle of the X-ray beam, the root tip may not be displayed correctly. The position of the root tip seems to differ from its true position.	
	The X-ray photo shows that the actual apex of the root canal is not the same as the anatomic end. In fact, the apical foramen is located at the coronal end. in this case, X-ray may indicate that the file needle has not reached the apical foramen, even if it has actually reached the apical foramen.

Cleaning, Disinfection and Sterilization

Preword

For hygiene and sanitary safety purpose, the file clip, lip hook, touch probe and extension cord must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well use the subsequent uses. Comply with your national guidelines, standards and requirements for cleaning, disinfection and sterilization.

Reprocessing procedures have only limited implications to this dental instrument. The limitation of the numbers of reprocessing procedures is therefore determined by the function / wear of the device. From the processing side there is no maximum number of allowable reprocessing. The device should no longer be reused in case of signs of material degradation. In case of damage, the device should be reprocessed before sending back to the manufacturer for repair.

General recommendations

- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments, where applicable after sterility.
- For your own safety, please wear personal protective equipment (gloves, safety glasses, etc.).
- Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, and FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer.
- The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.
- Thoroughly clean and wash the components before autoclaving.
- Do not use bleach or chloride disinfectant materials.

Autoclavable Components

File clip	Lip hook	Touch probe	Extension cord
			

WARNING

Only the components above can be autoclaved.

Before first use and after each use, sterilize the above components

Preparation at the Point of Use: Disconnect the components (Lip hook, file clip, touch probe and extension cord) from the main unit. Remove gross contaminations from the components with code water (40°C) immediately after use. Don't use a fixating detergent or hot water (>40°C) as this can cause the formation of residuals which may influence the result of the reprocessing process. Store the instruments in a humid surrounding.

WARNING

Do not submerge the components or wipe them with any of the following functional water (acidic electrolyzed water, strong alkaline solution, or ozone water), medical agents (glutaral, etc.), or any other types of water or commercial cleaning liquids. Such liquids may result in metal corrosion and fixation of the residual medical agents to the components.

Transportation: Safe storage and transportation to the reprocessing area to avoid any damage and contamination to the environment.

Preparation for Decontamination: The devices must be reprocessed in a disassembled state.

WARNING

Do not fail to take out the file before cleaning the file clip.

Observe suitable personal protective measures.

Cleaning: Do a manual pre-cleaning, until the components are visually clean. Submerge the components in a cleaning solution. Clean the surfaces with a soft bristle brush.

Cleaning: Regarding cleaning/disinfection, rinsing and drying, it is to distinguish between manual and automated reprocessing methods. Preference is to be given to automated reprocessing methods, especially due to the better standardizing potential and industrial safety.

Automated Cleaning:

Use a washer-disinfector meeting the requirements of the ISO 15883 series.

Carefully put the instrument into the washer-disinfector on a tray and set the parameters as follows start the program:

- 4 min pre-washing with cold water (<40°C)
- emptying
- 5 min washing with a mild alkaline cleaner at 55°C
- emptying
- 3 min neutralising with warm water (40°C)
- emptying
- 5 min intermediate rinsing with warm water (40°C)
- emptying

The automated cleaning processes have been validated by using 0.5% neodisher MediClean forte (Weigert).

Note Acc. to EN ISO 17664 no manual reprocessing methods are required for these devices. If a manual reprocessing method has to be used, please validate it prior to use.

WARNING

Use only approved washer-disinfectors according to EN ISO 15883, maintain and calibrate it regularly.

Follow instructions and observe concentrations given by the manufacturer (see general recommendations).

Disinfection: Automated Thermal Disinfection in washer/disinfector under consideration of national requirements in regards to A0 value (see EN ISO 15883).

Disinfection cycle of 5 min disinfection at 93°C has been validated for the device to achieve an A0 value 3000.

After manual cleaning, the instrument should be automated disinfected or sterilized immediately. A manual disinfection is not recommended.

Drying:

Automated Drying:

Drying of outside of instrument through drying cycle of washer/disinfector. If needed, additional manual drying can be performed through lint free towel. Insufflate cavities of instruments by using sterile compressed air.

Functional Testing, Maintenance: Visual inspection for cleanliness of the components and reassembling. Functional testing according to instructions of use. If necessary, perform reprocessing process again until the component is visibly clean.

Before packaging and autoclaving, make sure that the device has been maintained acc. to the manufacturer's instruction.

Packaging: Pack the instruments in an appropriate packaging material for sterilization.

WARNING

Check the validity period of pouch given by the manufacturer to determine the shelf life.

Use pouches which resist to a temperature up to 141°C and in accordance with EN ISO 11607

Sterilization: Sterilization of instruments by applying a fractionated pre-vacuum steam sterilization process (according to EN 285/EN 13060/EN ISO 17665) under consideration of the respective country requirements.

Minimum requirements: 5 min at 134 °C

Maximum sterilization temperature: 137°C

Drying time: 8 min

Flash sterilization is not allowed on lumen instruments!

WARNING

Use only approved autoclave devices according to EN 13060 or EN 285.

Use a validated sterilization procedure according to EN ISO 17665.

Respect the maintenance procedure of the autoclave device given by the manufacturer.

Use only this recommended sterilization procedure.

Control the efficiency (packaging integrity, no humidity, color change of sterilization indicators, physicochemical integrators, digital records of cycles parameters).

The sterilization procedure must comply with EN ISO 17665.

Wait for cooling before touching.

Storage: Storage of sterilized instruments in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.

WARNING

Integrity cannot be guaranteed if packaging is open, damaged or wet.
Check the packaging before using (packaging integrity, no humidity and validity period).

WARNING

The instructions provided above have been validated by the manufacturer of the medical device as being suitable for preparing a medical device for use. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

Disinfection components

APEX locator

Charge Base

Adapter

Water

Disinfect all the surfaces with a cloth lightly moistened with Ethanol for Disinfection (Ethanol 70 to 80 vol%) at least 2min, repeat for 5 times.

NOTE

Do not use anything except Ethanol for Disinfection (Ethanol 70 to 80 vol%).
Do not use too much ethanol as it's going into machine and damage the components inside.

Troubleshooting

When trouble is found, check the following points before contacting your distributor. If none of these are applicable or the trouble is not remedied even after action has been taken, the product may have failed. Contact your distributor.

Problem	Cause	Solution
The power is not turned on.	The battery is flat.	Charge the battery.
	Press the power switch too short time.	Long press the power switch.
The charge indicator flash on the piece screen.	Put the APEX locator on the charge base in the wrong location.	Check the location.
	Charging is completed.	Checking the instructions of the battery.
	The charge base is broken.	Contact your distributor.
No sound.	Beep volume is set to 0.	Set beep volume to 1, 2 or 3.

Operating time at the full charge and charging time of AirPex

Continuous operating time at the full charge: not less than 9 hours. Charging time of the battery: 4 hours $\pm$ 0.5 hours.

Label dimensions:

Device marking: label dimensions: 35x25mm $\pm$ 0.5mm

Box marking: Label dimensions: 90x70mm $\pm$ 0.5mm

Weight/size characteristics of the product and its components:

AirPex

Component	Weight	Size mm ($\pm 5\%$)
PEX Locator	14.5 $\pm$ 5g	28 x 48.5 x 15.3
Tip hook	2.5 $\pm$ 1g	64 x 24.6, d2
Touch probe	2.2 $\pm$ 0.5g	82, d2,5
File clip	7.6 $\pm$ 2g	416, d12.8
Master	2.2 $\pm$ 0.5	34 x 8 x 10, d2
Charge base	38.5 $\pm$ 5	60 x 60 x 14,4
Extension Cord	6.3 $\pm$ 1.5	386, d2.5
Adapter	82 $\pm$ 8g	78,5 x 74,5 x 29, cable 1200

Information on the range and accuracy of measurement of the distance to the apex, as well as the voltage, current strength and the pulse form created by the device during the apex localization.

Accuracy of the root canal measurement: ± 0.5 mm.

Current strength: is not an outcoming characteristic for the device

Voltage: 70 $\pm$ 10 mV

Pulse form: rectangular

Technical Data

Manufacturer	Changzhou Sifary Medical Technology Co.,Ltd
Model	AirPex
Dimensions	13cm x 11cm x 8cm±1cm (Outer box)
Weight	0.35kg±10% (in packaging)
Power supply	Li-Polymer Battery: 3.7V, 110mAh
Battery charge	Cordless from the charger
Charger supply	5V  1A From the adapter: Input (100 ÷ 240) V, 50/60 Hz, (0,5 ÷ 0,2) A Output 5 V dc, 1A
Level of protection against electrical shock	Device with an internal power supply and the working part BF type. Class II for battery charge
Working regime	Continuous
Shell protection	IPX0
Software	Version: AP1.1.1.010, Date: 2019.9.10, Class B.
Operating conditions	Use: in enclosed spaces Ambient temperature: 5°C ~ 40 °C Relative humidity: <80% Operating altitude < 3000m above sea level Atmospheric pressure: 70kPa~106kPa
Transport and storage conditions	Ambient temperature: -20 °C ~ +55 °C Relative humidity: 20% ~ 80 % Atmospheric pressure: 70kPa~106kPa

Guidance and manufacturer's declaration – electromagnetic emissions		
The AirPex is intended for use in the electromagnetic environment specified below. The customer or user of the AirPex should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
Emissions SPR 11	Group 1	The AirPex uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Emissions SPR 11	Class B	
Harmonic emissions C61000-3	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	The AirPex is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity			
The AirPex is intended for use in the electromagnetic environment specified below. The customer or user of the AirPex should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transients bursts IEC 61000-4-	±2kV 100kHz repetition frequency	±2kV 100kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage IEC 61000-4-	Line to line: ±0.5kV, ±1kV Line to earth: ±0.5kV, ±1kV, ±2kV	Line to line: ±0.5kV, ±1kV Line to earth: ±0.5kV, ±1kV, ±2kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips IEC 61000-4-	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0° 0% UT; 250/300 cycle	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0° 0% UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of devices require continued operation during power mains interruptions, it is recommended that devices be powered from an uninterruptible power supply or a battery
Power frequency magnetic field IEC 61000-4-	30 A/m 50Hz or 60Hz	30 A/m 50Hz or 60Hz	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT: rated voltage(s); E.g. 25/30 cycles means 25 cycles at 50Hz or 30 cycles at 60Hz

Guidance and manufacturer's declaration – electromagnetic immunity

The AirPex is intended for use in the electromagnetic environment specified below. The customer or user of the AirPex should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Induced disturbances produced by RF fields IEC 61000-4-6	3 V 0.15 MHz – 80 MHz, 6 V in ISM bands between 0.15 MHz and 80 MHz, 80 % AM at 1 kHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the AirPex, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended minimum separation distances See the RF wireless communication equipment table in "Recommended minimum separation distances"
Induced RF fields IEC 61000-4-3	3 V/m, 80 MHz – 2.7 GHz, 80 % AM at 1 kHz	3V/m	
Immunity fields in RF wireless communication equipment IEC 61000-4-3	See the RF wireless communication equipment table in "Recommended minimum separation distances"	Complies	

Recommended minimum separation distances

Today, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. AirPex has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipments and the AirPex as recommended below.

Test frequency (Hz)	Band (MHz)	Service	Modulation	Max power (W)	Distance (m)	Immunity test level (V/m)
30-50	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
40-60	430-470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
70-80	704-787	LTE Band 13 17	Pulse modulation 217Hz	0.2	0.3	9
80-100	800-960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18Hz	2	0.3	28
120-140	1700-1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation 217Hz	2	0.3	28
240-260	2400-2570	Bluetooth	Pulse modulation	2	0.3	28

		WLAN, 802.11 b/g/n, RFID, 2450, LTE, Band 7	217Hz			
240 500 785	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9

WARNING

Use of accessories and cables other than those specified or provided by the manufacturer of **AirPex** could result in increased electromagnetic emissions or decreased electromagnetic immunity of **AirPex** and result in improper operation.

Cable information:

Cable	Max. cable length shielded/unshielded	
Adapter cable	1.2m	Unshielded
Measuring cable	0.8m	Unshielded

Use of **AirPex** adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, **AirPex** and the other equipment should be observed to verify that they are operating normally.

Statement

Service Life

The service life is 3 years since it is sold (operation life started)

Warranty

The manufacturer warranty is 12 months since the device is sold (operation life started). Within 12 months since it is sold (operation life started) the manufacturer or its legal representative commits to fix the device for free or replace it with a new one if the quality of the device results poor.

Disposal

The package should be recycled. Metal parts of the device are disposed as scrap metal. Synthetic materials, electrical components, and printed circuit boards are disposed as electrical scrap. The lithium batteries are disposed as special refuse taking into consideration the local environmental protection laws and regulation. It is necessary to count with the local rules, laws and regulations when organizing the disposal of the medical device.

Product labels sample

Sample of device label (label size: 35x25mm)

APEX LOCATOR

Model: AirPex REF

Input: 5V = 1A Battery: DC 3.7V/110mAh

Changzhou Sifary Medical Technology Co., Ltd.
No.99, Qingyang Road, Xuejia County, Xinbei District,
Changzhou City, 213000 Jiangsu, China
Lins Service & Consulting GmbH
Tel: +49 175 4870819
Add: Obere Seegasse 34/2, 69124,
Heidelberg, Germany

Sample of outside box label (label size: 90x70mm)

APEX LOCATOR

Model: AirPex

REF

SN

Box Size: 132mmX113mmX78mm

Input: 5V = 1A

Battery: DC 3.7V/110mAh

Changzhou Sifary Medical Technology Co., Ltd.
No.99, Qingyang Road, Xuejia County, Xinbei District,
Changzhou City, 213000 Jiangsu, China
Lins Service & Consulting GmbH
Tel: +49 175 4870819
Add: Obere Seegasse 34/2, 69124, Heidelberg, Germany

EC REP

Applied standards:

- ISO 1640:2009
- ISO 13485:2016
- ISO 14971:2012
- ISO 60601-1:2006+A1:2013
- ISO 60601-1-2: 2015
- ISO 62304:2006+A1:2015
- ISO 10993-1:2018
- ISO 10993-5:2009
- ISO 10993-10:2010
- ISO 60601-1-6:2010 + A1: 2015
- ISO 62366-1:2015
- ISO 17665-1:2006
- ISO 1041:2008+A1:2013
- ISO 17664:2017
- ISO 15223-1:2016
- ISO 7405:2018

EN 80601-2-60:2015

EN ISO 780:2015

IEC 60601-2-60:2015

IEC 60601-2-60:2015

IEC 60601-2-60:2015

1. Information for identification of the medical device in terms of combination:

The "Apex locator device AirPex" is used with endodontic instruments (files) for cleaning and shaping root canals.

2. Last update IFU information:

This IFU has been last updated on May 01, 2021.

3. Contact and legal information:

Manufacturer: Changzhou Sifary Medical Technology Co., Ltd (Китай)

Address: No. 99, Qingyang Road, Xuejia County, Xinbei District, Changzhou City, 213000 Jiangsu, China.

Phone number: +86-0519-85962691

Email: ivy@sifary.com

Website: www.eighteenth.com

CE REP

Legal representative in EU: Llins Service & Consulting GmbH

Phone number: +49 175 4870819

Address: Obere Seegasse 34/2, 69124, Heidelberg, Germany

Email: Llins.Service@gmail.com

Legal representative in Russian Federation:

Dentalink LTD

Address: 123592 Moscow, Kulakova st, 20 building 1A, Unit XVII, room 64, floor 5.

Phone number: (495) 766-16-71

Email: info@dentalink.ru

Website: www.dentalink.ru

Signature

Wen Tao

General Manager




Attachment 1.


Maintenance Manual: AirPex

Warning

This product must only be repaired by manufacturer approved partner, Otherwise the warranty will be invalid.

Disassembly


Sketches	Note
	Remove the black rubber pads on the bottom
	Remove four screws from the bottom of the base with a cross screwdriver, then remove the upper housing
	Remove two screws (marked with red arrow) from the circuit board with cross screwdriver. Then separate the circuit board and charge base lower housing
	Melt the welding spot (marked by red arrow) of the charging coil and separate the charging coil and circuit board


Installation

Diagrammatic Sketch


Note

Weld the charging coil to the circuit board, then fix the charging coil with heat-resist tape

Assemble the circuit board and charge base lower housing in the figure, and fix with two ST2x4 screws

Assemble the upper housing, fasten with four ST2.2x6 screws, then attach the rubber pad at the bottom of the base

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(Translation)

(2021) SCCZZ Zi, No. 6655

Applicant: Changzhou Sifary Medical Technology Co., Ltd., Domicile: No. 99, Qingyang Road, Xuejia County, Xinbei District, Changzhou City.

Legal representative: Tao Huan, male, born on April 14, 1990, ID No. 320483199004146532.

Issue under notarization: Signature

This is to certify that Tao Huan, the legal representative of Changzhou Sifary Medical Technology Co., Ltd., personally appeared in my office on May 17, 2021, and signed in my presence the Instructions for use for the medical device "Apex locator device AirPex" attached hereto.

Notary: Lu Chunhong (Signature)
Changzhou Notary Public Office (Seal)
Changzhou City, Jiangsu Province
The People's Republic of China
May 24, 2021

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(Translation)

(2021) SCCZZ Zi, No. 6656

Applicant: Changzhou Sifary Medical Technology Co., Ltd., Domicile: No. 99, Qingyang Road, Xuejia County, Xinbei District, Changzhou City.

Legal representative: Tao Huan, male, born on April 14, 1990, ID No. 320483199004146532.

Issue under notarization: The translated copy in conformity with the original copy

This is to certify that the English translated copy of NOTARIAL CERTIFICATE (2021) SCCZZ Zi, No. 6655 attached hereto is in conformity with the original Chinese copy.

Notary: Lu Chunhong (Signature)

Changzhou Notary Public Office (Seal)

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
1. Name of the medical device:

“Apex locator device E-PEX”

2. Manufacturer of the medical device:

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Phone number: +86-0519-85962691

Email: ivy@sifary.com

3. Medical device purpose and indication:

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Field: dentistry. This medical device should be used only in specialized medical institutions by qualified medical personnel.



4. Possible side effects and risks:

Possible side effects: not found


Risks: The main cumulative residual risk: non-compliance with the operation instructions described in the operating manual. The main way to manage residual risks in order to reduce them to a permissible level: It is necessary to follow the instructions described in the Operation Manual. Before each use the endodontic device must be switched on, checked and tested to guarantee its uninterrupted operation.


5. Parts Identification:

1. Apex locator (main unit)
2. Measuring wire
3. File clip (2 pcs)
4. Lip hook (2 pcs)
5. Tester
6. Adapter

Apex Locator (1pcs) 	File clip (2pcs) 	Tester(1pcs) 
Adapter (1pcs) 	Measuring Wire (1pcs) 	Lip Hook (2pcs) 

6. Symbols used in the User Manual and marking

 WARNING	If the instructions are not followed properly, operation may lead to hazards for the product or the user/patient.
 NOTE	Additional information, explanation of operation and performance.
 SN	Serial number
 REF	Catalogue number
	Manufacturer
	Date of manufacture
	Class II equipment
	Type BF applied part
	CE marking
	Direct current
	Dispose of in accordance with the WEEE directive
	Keep dry
	Consult instructions for use
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified
	Authorized Representative in the European Community
	Temperature limitation
	Humidity limitation
	Manufacturer's LOGO
IPX0	Shell protection

7. Contraindications

Do not use this unit in conjunction with an electric scalpel or on patients who have a pacemaker. Blocked canals cannot be accurately measured. Safety and effectiveness

have not been established in pregnant women and children.

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
- 1、 The device must not be placed in humid surroundings or anywhere where it can come into contact with any type of liquids.
- 2、 Do not expose the device to direct or indirect heat sources. The device must be operated and stored in a safe environment.
- 3 、 The device requires special precautions with regard to electromagnetic compatibility (EMC) and must be installed and operated in strict compliance with the EMC information. In particular, do not use the device in the vicinity of fluorescent lamps, radio transmitters, remote controls, portable or mobile RF communication devices and do not charge, operate or store at high temperatures. Comply with the specified operating and storage conditions.
- 4、 Gloves and a rubber dam are compulsory during treatment.
- 5、 If irregularities occur in the device during treatment, switch it off. Contact the agency.
- 6、 Never open or repair the device yourself, otherwise, void the warranty.


8. Installing the E-PEX

Install the E-PEX

Insert the measuring wire into the socket as shown in left picture, make sure connect properly.

Connect the file clip, measuring wire and lip hook as shown in the picture.

NOTE

When installing the measuring wire, please pay attention to the orientation of the slots in the attachment part and do not apply too much force while adapting it.

Incorrect connection will result in inaccurate measurement, even the device cannot be used.


Connection Operation

Make sure E-CONNECT in standby.

Open rubber cover, plug data transfer cable into E-CONNECT.

Turn on the E-PEX, and insert the other end of data transfer cable into E-PEX.


After connect the cable, the screen of the E-CONNECT will display "CONNECTED !" indicating that the connection is properly.



CONNECTED !

E-PEX can only connect to E-CONNECT manufactured by Sifary.


After connecting E-CONNECT and E-PEX, do the below steps to make sure the device is working normally.

1. Insert the file into the contra angle.
2. Make the file touch the lip hook (short circuit)
3. Press the main switch of E-CONNECT. All the indicator bars in the screen will light up. That means the system is working normally.

After confirming the system can working normally, user can hang the lip hook into the patient's mouth, and start the treatment.


E-PEX Charging


When the power indicator flashes, please stop using the device and charge it immediately. We suggest the user to charge the device when there is only one bar left.

Connect the Apex Locator main unit with the power adapter.

When the power indicator is as shown below, it indicates that the device is in charging.

WARNING

Keep the device away from the heat source and make sure that there is no combustible surrounding.

When battery is low charge the device fully. Charging frequently in low power state for short time will reduce the battery life.

Do not use other power adapter to charge the device, otherwise it will damage the device.

Do not charge the device while using it.
Do not use other battery for the device, otherwise it will damage the device.

9. Operation and Functional characteristics

9.1 Functions setting


Function Checking

1. Press the Power switch to turn the device on. The display will show measuring interface.


(The device will automatically shut down if it is not used for 10 minutes.)

2. Check that the measuring wire, file clip, lip hook and APEX LOCATOR main unit are properly connected. Touch the metal part of the file clip with the lip hook (short circuit).


3. Observe the E-PEX display. All the meter indicator bars on the display will light up, and a rapid beep sound will be generated at the same time. The "APEX" sign will be flashed, which means that the E-PEX is working normally.


Volume control

The E-PEX's volume of the key and alarm sounds can be adjusted. Press the volume keys to cycle the volume through the minor to the maximum.

Setting the Reference point

Press SET switch to set the reference point (between 0~1).

Press SET to adjust reference point

The point will be automatically saved.

9.2 Display screen

Instruction

1. When the file reaches the front region of the apical foramen, the screen displays the white indicator bars (As shown in picture 1).


Pic. 1

2. When the file reaches the position near by the apical foramen, the screen displays the green indicator bars (As shown in picture 2).

Pic. 2

3. When the red indicator bars light up, it means that the file has exceeded the apical foramen. A rapid beep sound will be generated at the same time (As shown in picture 3).

Pic.3

WARNING

Avoid using apex locator for working length determination in the following conditions:

1. Open apex cases.
2. Draining canals.
3. Poor isolation from oral environment (avoid seepage of oral fluids into access cavity).
4. Root fractures / perforation.
5. Gutta percha filled canals;

Please use the original accessories, otherwise the device may measure inaccurately or not even function.

NOTE

The green part "00" display means major apical foramen (not the minor apical foramen). Hence it is recommended to reduce the working length by 0.5-1 mm.

The device's screen does not show the actual length of the root canal, the number reducing only means a trend that file is progressing apically.


The gingival crevicular fluid / saliva / gingival polyp will interfere with device functioning. Hence it is recommended to isolate the tooth.

The accessories which contact with patient (file clip and lip hook) can be reused and should be sterilized by high temperature before first use and after each use.

Display the root canal on


E-CONNECT

1. The white band on handpiece screen displays the progression of the file into the root canals.
2. The closer the file tip reaches the apical foramen, the more rapid the beep sound makes.
3. After connection, it will activate the advanced setting in chapter 9.5.

Combination Function

Set "ON" to choose the combination function.

The position of the reference point is automatically set with the E-PEX, and the cursor is displayed on the E-CONNECT screen.

When the file reaches the reference point, E-CONNECT will start Apical Reverse, Apical Slow Down and Apical Torque Reduction function (If the function is activated).

WARNING

Do not use a non - specified data transfer cable, otherwise it will damage the device.

Do not hit device and splash liquids.

NOTE

Make sure to connect the two devices with right position.

After connecting the two devices with the cable, gently push and pull the interface to ensure that the connection is stable, otherwise the data transmission may not be accurate.

In certain cases, for example when the canal is blocked, the measurement may be unable.

The device will not be able to perform a precise measurement for every time, especially in cases of abnormal or unusual morphology of the root canal. The user needs to coordinate with x-ray to check the results of the measurement.

If the meter does not move when you insert the file, it is possible that the device is not working normally, therefore, stop using.

10. Cleaning, Disinfection and Sterilization

Foreword

For hygiene and sanitary safety purpose, the components (file clip, lip hook) must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well use the subsequent uses. Comply with your national guidelines, standards and requirements for cleaning, disinfection and sterilization.

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General recommendations

- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments, where applicable after sterility.
- For your own safety, please wear personal protective equipment (gloves, safety glasses, etc.).
- Use only a disinfecting solution which is approved for its efficacy (VAH/DGHEM-listing, CE marking, and FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer.
- The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.
- Thoroughly clean and wash the components before autoclaving.
- Do not use bleach or chloride disinfectant materials.

Autoclavable Components	
	
File clip	Lip hook

WARNING

- Only the components above can be autoclaved.
- Before first use and after each use, sterilize the above components heat sources.

Autoclave Procedure:



Reprocessing Instructions	
Preparation at the Point of Use:	Disconnect the components (Lip hook and file clip) from the main unit. Remove gross contaminations from the components with code water (<40°C) immediately after use. Don't use a fixating detergent or hot water (>40°C) as this can cause the fixation of residuals which may influence the result of the reprocessing process. Store the instruments in a humid surrounding.  WARNING Do not submerge the components or wipe them with any of the following functional water (acidic electrolyzed water, strong alkaline solution, or ozone water), medical agents (glutaral, etc.), or any other special types of water or commercial cleaning liquids. Such liquids may result in metal corrosion and adhesion of the residual medical agents to the components.
Transportation:	Safe storage and transportation to the reprocessing area to avoid any damage and contamination to the environment.
Preparation for Decontamination:	The devices must be reprocessed in a disassembled state.  WARNING - Do not fail to take out the file before cleaning the file clip. - Observe suitable personal protective measures.
Pre-Cleaning:	Do a manual pre-cleaning, until the components are visually clean. Submerge the components in a cleaning solution and flush the lumens with a water jet pistol with cold tap water for at least 10 seconds. Clean the surfaces with a soft bristol brush.
Cleaning:	Regarding cleaning/disinfection, rinsing and drying, it is to distinguish between manual and automated reprocessing methods. Preference is to be given to automated reprocessing methods, especially due to the better standardizing potential and industrial safety. Automated Cleaning: Use a washer-disinfector meeting the requirements of the ISO 15883 series. Carefully put the instrument into the washer-disinfector on a tray and set the parameters as follows and start the program: 4 min pre-washing with cold water (<40°C) - emptying 5 min washing with a mild alkaline cleaner at 55°C

	• emptying • 3 min neutralising with warm water (>40°C) • emptying • 5 min intermediate rinsing with warm water (>40°C) • emptying <i>The automated cleaning processes have been validated by using 0.5% neodisher MediClean forte (Dr. Weigert).</i> Note Acc. to EN ISO 17664 no manual reprocessing methods are required for these devices. If a manual reprocessing method has to be used, please validate it prior to use.  WARNING • Use only approved washer-disinfectors according to EN ISO 15883, maintain and calibrate it regularly. • Follow instructions and observe concentrations given by the manufacturer (see general recommendations).
Disinfection:	Automated Thermal Disinfection in washer/disinfector under consideration of national requirements in regards to A0 value (see EN ISO 15883). A disinfection cycle of 5 min disinfection at 93°C has been validated for the device to achieve an A0 value of 3000. After manual cleaning, the instrument should be automated disinfected or sterilized immediately. A manual disinfection is not recommended.
Drying:	Automated Drying: Drying of outside of instrument through drying cycle of washer/disinfector. If needed, additional manual drying can be performed through lint free towel. Insufflate cavities of instruments by using sterile compressed air.
Functional Testing, Maintenance:	Visual inspection for cleanliness of the components and reassembling. Functional testing according to instructions of use. If necessary, perform reprocessing process again until the component is visibly clean. Before packaging and autoclaving, make sure that the device has been maintained acc. to the manufacturer's instruction.
Packaging:	Pack the instruments in an appropriate packaging material for sterilization.  WARNING • Check the validity period of pouch given by the manufacturer to determine the shelf life. • Use pouches which resist to a temperature up to 141°C and in accordance with EN ISO 11607.

Sterilization	Sterilization of instruments by applying a fractionated pre-vacuum steam sterilization process (according to EN 285/EN 13060/EN ISO 17665) under consideration of the respective country requirements. Minimum requirements: 3 min at 134 °C (in EU: 5 min at 134 °C) Maximum sterilization temperature: 137°C Flash sterilization is not allowed on lumen instruments!  WARNING ● Use only approved autoclave devices according to EN 13060 or EN 285. ● Use a validated sterilization procedure according to EN ISO 17665. ● Respect the maintenance procedure of the autoclave device given by the manufacturer. ● Use only this recommended sterilization procedure. ● Control the efficiency (packaging integrity, no humidity, color change of sterilization indicators, physicochemical integrators, digital records of cycles parameters). ● The sterilization procedure must comply with EN ISO 17665. ● Wait for cooling before touching.
Storage:	Storage of sterilized instruments in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.  WARNING Sterility cannot be guaranteed if packaging is open, damaged or wet. Check the packaging before using (packaging integrity, no humidity and validity period).
Reprocessing validation study information:	The above-mentioned reprocessing process (cleaning, disinfection, sterilization) has been successfully validated. Refer to test reports: - Changzhou Sifary _Cleaning Disinfection Validation Report - Changzhou Sifary _Sterilization Validation Report_File clip
 WARNING The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for use. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility,	

achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

Disinfection

Wipe with Ethanol for Disinfection Ethanol 70 to 80 vol%		Wipe all the surfaces with a cloth lightly moistened with Ethanol for Disinfection (Ethanol 70 to 80vol%) at least 2 min, repeat for 5 times.
		
Charger	Measuring wire	 NOTE Do not use anything except Ethanol for Disinfection (Ethanol 70 to 80 vol%). Do not use too much ethanol as it's going into machine and damage the components inside.

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Problem	Cause	Solution
The power is not turned on.	The battery is flat.	Charge the battery.
	Press the power switch too short time.	Long press the power switch.
No charge indicator flash on handpiece screen.	Put the APEX locator on the charge base in the wrong location.	Check the location.
	Charging is completed.	Checking the instructions of the battery.
	The charge base is broken.	Contact your distributor.
No sound.	Beep volume is set to 0.	Set beep volume to 1, 2 or 3.

12. Operating time at the full charge and charging time of E-PEX

Continuous operating time at the full charge: not less than 10 hours. Charging time of the battery: 4 hours $\pm$ 0.5 hour.

13. Label dimensions:

Device marking: label dimensions: 50x25mm $\pm$ 0.5mm

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14. Weight and size characteristics of the product and its components:

E-PEX

Component	Weight	Size mm ($\pm$ 5%)
Apex Locator (main unit)	142 $\pm$ 10g	67.9 x 97 x 46.6 mm
Measuring wire	21 $\pm$ 5g	1500, d3.7
File clip	7 $\pm$ 2g	250, d9
Lip hook	2.5 $\pm$ 1g	64 x 24.6, d2
Tester	2 $\pm$ 0.5	29 x 8.4, d10
Adapter	82 $\pm$ 8g	78.5 x 74.5 x 29, cable 1200

15. Information on the range and accuracy of measurement of the distance to the apex, as well as the voltage, current strength and the pulse form created by the device during the apex localization.

Accuracy of the root canal measurement: $\pm$ 0.5mm.

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Pulse form: rectangular

16. Technical Data

Manufacturer	Changzhou Sifary medical technology Co., Ltd
Model	E-PEX
Dimensions	20cm x 17.5cm x 7cm $\pm$ 1cm (Outer box)
Weight	600g $\pm$ 10% (in packaging)
Display	3.5' color LCD
Power supply	Lithium ion battery: 3.7V, 1500mAh
Battery Charger	5V $\pm$ 1A From the adapter: Input (100 $\div$ 240) V, 50/60 Hz, (0,5 $\div$ 0,2) A Output 5 V dc, 1A
Type of protection against electrical shock	Device with an internal power supply and the working part BF type. Class II for battery charge
Working regime	Continuous
Shell protection	IPX0
Software	Version: EP 1.1.1.1, Date: 2018.5.12, Class B
Operation conditions	Use: in enclosed spaces Ambient temperature: 5°C ~ 40 °C Relative humidity: <80% Operating altitude < 3000m above sea level Atmospheric pressure: 70kPa ~ 106kPa
Transport and storage conditions	Ambient temperature: -20 °C ~ +55 °C Relative humidity: 20% ~ 80 % Atmospheric pressure: 70kPa ~ 106kPa

17. EMC Tables

Guidance and manufacturer's declaration – electromagnetic emissions		
The E-PEX is intended for use in the electromagnetic environment specified below. The customer or the user of the E-PEX should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The E-PEX uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	The E-PEX is suitable for use in all

CISPR 11		establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The E-PEX is intended for use in the electromagnetic environment specified below. The customer or the user of the E-PEX should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transients/bursts IEC 61000-4-4	±2kV 100kHz repetition frequency	±2kV 100kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Line to line: ±0.5kV, ±1kV Line to earth: ±0.5kV, ±1kV, ±2kV	Line to line: ±0.5kV, ±1kV Line to earth: ±0.5kV, ±1kV, ±2kV	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips IEC 61000-4-11	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of devices require continued operation during power mains interruptions, it is recommended that devices be powered from an uninterruptible power supply or a battery
	0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0°	0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0°	
Voltage interruptions IEC 61000-4-11	0% UT; 250/300 cycle	0% UT; 250/300 cycle	
Rated Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz or 60Hz	30 A/m 50Hz or 60Hz	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: UT: rated voltage(s); E.g. 25/30 cycles means 25 cycles at 50Hz or 30 cycles at 60Hz			

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The **E-PEX** is intended for use in the electromagnetic environment specified below. The customer or the user of the **E-PEX** should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted dis-turbances induced by RF fields IEC 61000-4-6	3 V 0.15 MHz – 80 MHz, 6 V in ISM bands between 0.15 MHz and 80 MHz, 80 % AM at 1 kHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the E-PEX , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF EM fields IEC 61000-4-3	3 V/m, 80 MHz – 2,7 GHz, 80 % AM at 1 kHz	3V/m	Recommended minimum separation distances See the RF wireless communication equipment table in "Recommended minimum separation distances"

Proximity fields from RF wireless communication equipment IEC 61000-4-3	See the RF wireless communication equipment table in "Recommended minimum separation distances"	Complies	
---	---	----------	--

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Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. The **E-PEX** has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipments and the **E-PEX** as recommended below.

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900;	Pulse modulation 217Hz	2	0.3	28
1845						
1970						

		DECT; LTE Band 1, 3, 4, 25; UMTS				
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9
5500						
5785						

WARNING

1. Use of accessories and cables other than those specified or provided by the manufacturer of **E-PEX** could result in increased electromagnetic emissions or decreased electromagnetic immunity of **E-PEX** and result in improper operation.

Cable information:

Cable Name	Cable Length (m)	Shielded or not	Remark
Adapter Cable	1.2	No	/

2. Use of **E-PEX** adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, **E-PEX** and the other equipment should be observed to verify that they are operating normally.

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Disposal

The package should be recycled. Metal parts of the device are disposed as scrap metal. Synthetic materials, electrical components, and printed circuit boards are disposed as electrical scrap. The lithium batteries are disposed as special refuse taking into consideration the local environmental protection laws and regulation. It is necessary to count with the local rules, laws and regulations when organizing the refusal of the medical device.

19. Product labels sample

1) Sample of device label (label size: 50x25mm)

2) Sample of outside box label (label size: 90x70mm)

20. Applied standards:

EN 1640:2009
EN ISO 13485:2016
EN ISO 14971:2012
EN 60601-1:2006+A1:2013
EN 60601-1-2: 2015
EN 62304:2006+A1:2015
ISO 10993-1:2018
ISO 10993-5:2009
ISO 10993-10:2010
EN 60601-1-6:2010 + A1: 2015
EN 62366-1:2015
EN ISO 17665-1:2006

EN 1041:2008+A1:2013
EN ISO 17664:2017
EN ISO 15223-1:2016
EN ISO 7405:2018
EN 80601-2-60:2015
EN ISO 780:2015
MEDDEV 2.7/1 rev.4
MDD93/42/EEC+2007/47/EC
MEDDEV 2.12-1 Rev. 8

21. Information for identification of the medical device in terms of combination:

The "Apex locator device E-PEX" is used only in combination with the medical device "Endomotor for working with NiTi instruments E-CONNECT" manufactured by Changzhou Sifary Medical Technology Co., Ltd (China)

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Email: info@dentalink.ru

Website: www.dentalink.ru



Signature
Huan Tao
General Manager



Attachment 1.

Maintenance Manual: E-PEX




WARNING This product must only be repaired by manufacturer approved partner,
Otherwise the warranty will be invalid.



Disassembly



Step	Diagrammatic Sketch	Description
1		Tear off the back label
2		Remove two screws under the label counterclockwise with the cross screwdriver
3		Take the left black foam pad and unplug the battery connector with tweezers to remove the battery

4		Remove the four screws on the back side counterclockwise with the cross screwdriver, then turn the whole back, and please note that the front cover is off. You need to turn it around with the front cover altogether
5		Remove the two screws from the circuit board as shown in the figure, Lift the circuit board from the charging port side.

Installation

Step	Diagrammatic Sketch	Description
1		Install the circuit board and screw it
2		Install front cover and screw it
3		Install battery ,Please pay attention to the battery joint boss needs to be aligned with the mounting joint groove, then insert black foam pad under the bottom of battery and put on the battery cover

4		After install the battery cover and screw it up
5		Put the original label on the battery cover

4		After install the battery cover and screw it up
5		Put the original label on the battery cover

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This is to certify that Tao Huan, the legal representative of Changzhou Sifary Medical Technology Co., Ltd., personally appeared in my office on May 17, 2021, and signed in my presence the Instructions for use for the medical device "Apex locator device E-PEX" attached hereto.

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This is to certify that the English translated copy of NOTARIAL CERTIFICATE (2021) SCCZZ Zi, No. 6653 attached hereto is in conformity with the original Chinese copy.

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