Medtronic
Further Together


Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA

ATTESTATION


STATE OF MINNESOTA

HENNEPIN COUNTY

As a responsible member of Medtronic Inc., I hereby attest that I believe the Recharger 37751 Charging System User Manual attached herein to be a true and accurate copy.


Ashley Luloff
Regulatory Affairs Specialist

Official Notary Signature

United States of America
State of Minnesota
Hennepin County

Subscribed and sworn to before me
This 22 day of March 2017

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432, USA
Tel: +1 763 514 4000

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CE0123
2005

M924449A048 Rev A 2015-04-15

Printing instructions: Refer to the "Patient Manual" category in doc# A00002 for Neuro Core European Printing Instructions.

Medtronic® is a registered trademark of Medtronic, Inc.

[USA] FCC Information

The following is communications regulation information on the Model 37751 Recharger.

FCC ID: LF537751

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

IMPORTANT: Changes or modifications to this product not authorized by Medtronic, Inc., could void the FCC Certification and negate your authority to operate this product.

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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M924449A048 Rev A 2015-04-15

Label symbols

Explanation of symbols on products and packaging. Refer to the appropriate product to see symbols that apply.

Consult instructions for use

Manufacturer

Temperature limitation

Keep dry

IP21

Ingress protection rating IP21, per 60601-1-11

IP22

Ingress protection rating IP22, per 60601-1-11

Serial number

PIN number

Magnetic Resonance (MR) Unsafe

Non-ionizing electromagnetic radiation

IEC 60601-1/EN60601-1, Type BF Equipment

Label symbols

Chinese Standard (SJ/T11364-2006)
Logo: Electronic Information Products
Pollution Control Symbol. (The date in this
logo means the environmental protection
use period of the product.)

System meets the applicable Canadian
(CAN/CSA-C22.2 No. 60601-1) and US
(UL 60601-1:2003) electrical safety
standard requirements.

The model 37761 meets the applicable
Canadian (CAN/CSA-C22.2 No. 60601-1)
and US (UL 60601-1) electrical safety
standard requirements.

Conformité Européenne (European
Conformity). This symbol means that the
device fully complies with AIMD Directive
90/385/EEC (NB 0123) and R&TTE
Directive 1999/5/EC.

Authorized Representative in the
European Community

For USA audiences only


Do not dispose of this product in the
unsorted municipal waste stream.
Dispose of this product according to local
regulations. See [http://
recycling.medtronic.com](http://recycling.medtronic.com) for instructions
on proper disposal of this product.

Label symbols

Filename Date Time
UC200xxxxxx EN
4.625" x 6.0" inches (117 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.2: 03-27-2015

Rechargeable battery

Class II device

Label symbols

37751 2015-04-15 English

M924449A048 Rev A 2015-04-15

Filename Date Time
UC200xxxxxx EN
4.625" x 6.0" inches (117 mm x 152 mm)

Medtronic Confidential
PPManual.xml - PatientProgrammerTemplate.fm
Template version 6.2: 03-27-2015

Label symbols

10 English 37751 2015-04-15

M924449A048 Rev A 2015-04-15

Table of contents

Label Symbols 7

Glossary 14

1 Introduction 18

- Purpose of the device 18
- Warnings 18
- Precautions 19
- How to use this manual 23
- Components of your charging system 25
- Recharger keys 28
- Recharger screens 29

2 Using your charging system 34

- Battery and charging information 34
 - Recharger and neurostimulator 34
 - Recharger 35
 - Neurostimulator 35
- Neurostimulator battery status 36
- Charging your neurostimulator battery 40

Table of contents

Accessories 47

Using the recharger pouch 47

Using the charging system carrying
case 48

Identification label 48

Additional features (available for some
neurostimulator models) 48

Turning the neurostimulator on and off
while charging 48

3 Troubleshooting 52

Recharger screens 52

Warning screens 52

Information screens 57

Possible problems and solutions 65

User assistance 74

4 Maintenance 76

Recharger maintenance tasks 76

Checking the recharger battery 76

Charging your recharger battery 78

Cleaning and Care 82

Safety and technical checks 83

Battery and device disposal 83

Declaration of conformity 84
Specifications 84

Index 87

Glossary

Caution - A statement describing actions that could result in damage to or improper functioning of a device.

Clinician - A healthcare professional such as a doctor or nurse.

Contraindication - A condition or circumstance when a person should not have a neurostimulation system.

Discharged battery - The rechargeable battery is depleted and should be charged as soon as possible. When the battery is in a discharged state, therapy is not available.

Indication - The purpose of the neurostimulation system and the medical condition for which it may be implanted.

Neurostimulator - The power source of a neurostimulation system. It contains the battery and electronics that control the stimulation.

Overdischarge - The rechargeable battery cannot be recharged when it is overdischarged. Once a battery is in a discharged state, if it is not recharged it continues to lose charge.

Patient Programmer - A hand-held device that allows you to check your neurostimulator battery daily. It is also used to adjust some of the stimulation settings.

Precaution - See Caution.

Recharger - The component of the neurostimulation system that is used to charge your neurostimulator battery.

Template - The cut-out area in the plastic part of the charging system belt.

Warning - A statement describing an action or situation that could harm the patient.

Warning screen - A screen displayed on the recharger that alerts you to a problem with the charging system.

Filename Date Time
UC200xxxxxx EN
4.625" x 6.0" inches (117 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.2: 03-27-2015

Glossary

16 English 37751 2015-04-15

M924449A048 Rev A 2015-04-15

Filename Date Time
UC200xxxxxx EN
4.625" x 6.0" inches (117 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.2: 03-27-2015

1 Introduction

M924449A048 Rev A 2015-04-15

Purpose of the device

The Medtronic Model 37751 Charging System is designed to charge Medtronic rechargeable neurostimulators. The indications, contraindications, risks and benefits, warnings, precautions, and individualization of treatment information for your specific neurostimulation system are located in the patient programmer user manual supplied with your patient programmer.

Warnings

Wound contact—DO NOT use the recharger on an unhealed wound. The recharging system is not sterile, and contact with the wound may cause an infection.

Power supply—Do not use the power supply near water or during a lightning storm. Electric shock, serious injury, or death may result.

Precautions

Low battery charge level—Charge the neurostimulator when you see a **Low battery** screen displayed on the patient programmer or recharger. The battery charge level is at 0 and your therapy will stop soon unless you recharge the battery.

If you see the following screen (Figure 1.1), the neurostimulator battery is in a discharged state and your therapy has stopped.


Figure 1.1 Neurostimulator battery is in a discharged state.

You should charge the battery as soon as possible. If the battery remains in a discharged state, it continues to lose charge and it will become overdischarged.

If the neurostimulator battery is allowed to overdischarge, you may not be able to

recharge or communicate with the neurostimulator; however, your doctor may be able to restore the battery function.

Allowing the neurostimulator battery to overdischarge will permanently affect the neurostimulator in one of the following ways:

- Battery function is restored; however, charging sessions may be more frequent because battery capacity has been reduced.
- Battery function is not restored and the neurostimulator must be surgically replaced. Battery function is not restored because:
 - the neurostimulator battery is permanently damaged.
 - the neurostimulator battery has been overdischarged and restored twice before. The third time the battery is overdischarged, the neurostimulator will reach end of service. Surgery is required to replace the neurostimulator.

Battery charge level—Check battery status and charge the battery regularly. Battery capacity will continue to slowly drain, even when the neurostimulator is off. If the battery becomes completely discharged, therapy will be lost.

Storage temperature: recharger and antenna—Do not expose the recharger or antenna to extremely hot [eg, hot cars, saunas, or direct sunlight; 40 °C (104 °F)] or cold [eg, walk-in freezers; –20 °C (–4 °F)] environments. The recharger and antenna are sensitive to temperature extremes and can become damaged.

Component compatibility (recharger)—Use only the Medtronic Neurological charging system components to charge your recharger. Use of components other than Medtronic charging system components to charge the recharger can damage your charging system components.

Patient device use—When operating a patient control device (eg, external neurostimulator, patient programmer,

controller, recharging system), use special care near flammable or explosive atmospheres. An interaction between the flammable or explosive atmospheres and the battery in the device could occur. The consequences of using a battery-powered device near flammable or explosive atmospheres are unknown.

Equipment modification—Do not modify this equipment. Modification of this equipment can result in damage to the device, causing the device to malfunction or become unusable.

Interference from EMI—When using the recharger to communicate with your neurostimulator, move away from equipment that may generate electromagnetic interference (EMI) or turn off the likely source of EMI. EMI may disrupt communication between the recharger and neurostimulator. Examples of EMI sources are computer monitors, cellular telephones, and motorized wheelchairs.

Recharger use—Check for skin irritation or redness near the neurostimulator during recharging. Do not sit or lie on the antenna, or apply excessive pressure to the antenna. Take periodic breaks during prolonged recharging. Although no direct cause and effect has been established, some patients have reported heating sensation, discomfort, blistering not caused by heating, skin irritation, or redness near the implanted neurostimulator during or after recharging.

For additional information, refer to "Possible problems and solutions" on page 65.

How to use this manual

This manual is designed to help you understand your charging system. Ask your clinician to explain anything that is unclear.

- Chapter 1 "Introduction", describes the recharger keys and screens and how the system works.

- Chapter 2 "Using your charging system", describes the recharger and how to perform specific tasks.
- Chapter 3 "Troubleshooting", describes recharger warning and information screens, how to solve possible charging system problems, and who to contact if your recharger is lost or broken.
- Chapter 4 "Maintenance", describes how to care for your recharger and contains system specifications and battery information.
- A glossary is provided at the beginning of this manual.

Components of your charging system


Figure 1.2 Components of your charging system.

Recharger: The recharger is a hand-held device used to charge the battery in your neurostimulator. The display on the recharger provides information on your charging system.

Recharger battery: The recharger is used to charge your implanted neurostimulator battery. The recharger has its own rechargeable battery that is similar to the

rechargeable batteries commonly found in electronic devices. The recharger battery typically lasts for approximately 5 years, depending upon your therapy needs and how often you use your recharger.

Rechargeable batteries lose some of their power every day, even when not used. Charging sessions using the AC power supply and power cord restore power to the recharger battery. You can charge the recharger battery anytime, including during a recharging session for your neurostimulator.

AC power supply and power cord (Model 37761): This device is used to charge the recharger using AC power. If you need to charge your neurostimulator but your recharger battery is low, connect your recharger to the AC power supply to conduct a charging session. The power cord plugs into an electrical outlet and connects the AC power supply to the recharger.

Note: Use only the power supply provided by Medtronic with your charging system.

Antenna: The antenna is connected to the recharger and is used to establish communication with your neurostimulator. You can adjust a dial on the antenna to improve the charging efficiency between the recharger and the neurostimulator.

Belt: If your neurostimulator is implanted in your abdominal area, then you will be given a belt to allow you to keep the antenna directly over your neurostimulator during a charging session.

Holster: If your neurostimulator is implanted in your upper chest area, then you will be given a holster to allow you to keep the antenna directly over your neurostimulator during a charging session.

Pouch: The recharger fits into a pouch that can be attached to a belt around your waist.

System carrying case: The carrying case is used to store the recharger, the charging system components, and your patient programmer.

Recharger keys

Three keys on the recharger allow you to control the recharger.


Figure 1.3 Recharger keys.

- **Audio control** : Turns your recharger on. When your recharger is on, pressing the **Audio control** key turns audible status tones on and off.
- **Start charge/test** : Begins charging the neurostimulator.
- **Stop charge** : Ends or interrupts charging the neurostimulator.

Recharger screens

The recharger screen displays information about your charging session (see Figure 1.4).


Figure 1.4 Recharger screens.

Status row — Displays icon information for the recharger and neurostimulator.

During a neurostimulator or recharger charging session, you will usually see the following icons on the **Status row**:

- Neurostimulator is on  or off . If the stimulation status is unknown, this area of the screen will be blank.
- Audio (sound) is on  or off .

- Battery charge status of the recharger
 (icons show from “needs recharging” to “full,” respectively).
- The recharger is plugged into the AC power supply . If the recharger is not plugged into a power supply, this area of the screen will be blank.

Charge progress row — Indicates the progress of the charging session.

During a neurostimulator charging session, you will see the following icon on the **Charge progress row**: a flashing **Battery** icon , indicates the neurostimulator battery is charging.

During a recharger charging session, you will see the following icons on the **Charge progress row**:

- A flashing **Battery** icon  indicates the recharger is charging.
- The **Connection** icon , indicates that the power supply is connected to an operational electrical outlet.

Charge efficiency row — Indicates the efficiency of the charging signal between the neurostimulator and recharger. During a neurostimulator charging session, the more solid black boxes ■■■■■■ that appear on the **Charge efficiency row**, the more efficient the charging signal and the shorter the charging session will be. If you see 6 or more empty boxes ■■■■■■, adjust the antenna to improve the signal strength between the neurostimulator and the recharger.

Filename Date Time
UC200xxxxxx EN
4.625" x 6.0" inches (117 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.2: 03-27-2015


Introduction 1

32 English 37751 2015-04-15

M924449A048 Rev A 2015-04-15

Filename Date Time
UC200xxxxxx EN
4.625" x 6.0" inches (117 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.2: 03-27-2015

2 Using your charging system

M924449A048 Rev A 2015-04-15

This chapter includes information on:

- the recharger battery and the charging procedure.
- the neurostimulator battery status icons.
- charging your neurostimulator.
- the charging system accessories.
- the charging system optional features.

Battery and charging information

Recharger and neurostimulator

- Charge the recharger and neurostimulator even if you are not using them for an extended period of time.
- You can charge the neurostimulator and recharger at the same time. Follow the instructions for both procedures ("Charging your neurostimulator battery" on page 40 and "Charging your recharger battery" on page 78) at the same time.

Using your charging system 2

- Charge the batteries in your recharger and neurostimulator before they are low.

Recharger

- The recharger battery might have run low since it was manufactured and shipped. You might need to charge it before you use it the first time.
- As your recharger ages, you will need to charge the battery more often.

Neurostimulator

- Your ability to locate your neurostimulator and the depth of the neurostimulator determine the efficiency of the charging session. With practice, you will find the best location on your skin for placing your recharger and antenna.
- If your neurostimulator is on while you are charging the battery, your charging session might last longer.
- Depending on your therapy needs, charging your neurostimulator could take

from one hour to more than 12 hours per week.

- You can partially charge your neurostimulator if you do not have enough time for a complete charging session.
- You can charge your neurostimulator at any time; you do not have to wait for a low battery message.
- Depending on where your neurostimulator is implanted, you will be given a belt that is placed around your waist or a holster that is placed over your shoulder to hold the antenna in place during your charging session.



Neurostimulator battery status

Use your patient programmer to check your neurostimulator battery charge level every day. Refer to the manual supplied with the patient programmer.

The recharger will also indicate if your neurostimulator battery needs to be recharged. If either message in Table 2.1 is

displayed on the recharger, recharge the neurostimulator battery immediately.

Table 2.1 When to charge the neurostimulator

When you see this	It means
	The neurostimulator battery charge level is low, and soon your therapy will not be available. Charge the neurostimulator battery. Press the gray Audio key on the recharger to clear the message.
	The neurostimulator battery charge level is depleted; the battery is in a discharged state. Therapy has stopped. Charge the neurostimulator battery now. Press the green Start charge/test key to clear the screen.

△ **Caution:** Charge the neurostimulator when you see a **Low battery** screen displayed on the patient programmer or recharger. The battery charge level is at 0

Using your charging system 2

and your therapy will stop soon unless you recharge the battery.

If you see the following screen (Figure 2.1), the neurostimulator battery is in a discharged state and your therapy has stopped.


Figure 2.1 Neurostimulator battery is in a discharged state.

You should charge the battery as soon as possible. If the battery remains in a discharged state, it continues to lose charge and it will become overdischarged.

If the neurostimulator battery is allowed to overdischarge, you may not be able to recharge or communicate with the neurostimulator; however, your doctor may be able to restore the battery function.

Allowing the neurostimulator battery to overdischarge will permanently affect the neurostimulator in one of the following ways:

- Battery function is restored; however, charging sessions may be more frequent because battery capacity has been reduced.
- Battery function is not restored and the neurostimulator must be surgically replaced. Battery function is not restored because:
 - the neurostimulator battery is permanently damaged.
 - the neurostimulator battery has been overdischarged and restored twice before. The third time the battery is overdischarged, the neurostimulator will reach end of service. Surgery is required to replace the neurostimulator.

Using **your** charging system 2

Charging your neurostimulator battery

Check with your clinician prior to using your recharger for the first time. When your neurostimulator needs to be charged, follow the steps below.

Notes:

- You might need to charge the recharger battery before you charge your neurostimulator battery. Refer to "Recharger maintenance tasks" on page 76.
- Based on the position of your implanted neurostimulator, you might need to use more than one template (Figure 2.2) to ensure the antenna remains positioned correctly over your neurostimulator.
- Based on the location of your implanted neurostimulator, you will use either a belt around your waist or a holster over your shoulder to hold the antenna in place.

- Do not place the antenna over bulky clothing or bandages.
 - The antenna may be placed directly on your skin. However, if the antenna irritates your skin, you might want to wear a thin undergarment to protect your skin.
1. Press your skin where your neurostimulator is located to determine the best place for the antenna.
 2. Put on your belt or holster, with the template or opening centered over your neurostimulator (Figure 2.2) and secure it in place. Reposition the template or opening, if necessary.


Figure 2.2 Place the belt or holster over your neurostimulator.

3. Place the antenna in the template or opening, and press it into the belt or holster until it snaps into place (Figure 2.3).


Figure 2.3 Attach the antenna.

4. Press the green **Start charge/test**  key on the recharger.
 - a. The **Recharger wait** screen will appear on the recharger display (Figure 2.4). The **Recharger wait** screen indicates that the neurostimulator and recharger are trying to communicate.


Figure 2.4 Recharger wait screen.

- b. When charging begins, the **Neurostimulator charging session** screen will appear. This screen indicates the neurostimulator battery charge level and efficiency of the charging session (Figure 2.5).


Figure 2.5 Neurostimulator charging session screen.

- c. The **Battery charge level** updates throughout the charging session (Figure 2.6). Charging sessions can take from a few minutes to more than

12 hours, depending on the charging efficiency, how often you charge, and your therapy needs.

Figure 2.6 Battery charge levels.

Notes:

- If you need to interrupt your charging session, press the **Stop charge**  on the recharger.
 - You might hear periodic clicking noises during recharging, and you might see the screen refresh.
5. Check the progress of the charging session by reviewing the icons on your recharger screen. If you see 6 or more empty boxes ■■■■■■, reposition the antenna or use the dial on the antenna to try to make the session more efficient.

Using your charging system 2

- a. Turn the dial on the antenna one quarter turn in either direction (Figure 2.7).


Figure 2.7 Adjusting the antenna dial.

- b. Press the green **Start charge/test**  key. The recharger will display an updated efficiency status.
- c. Repeat step a and b above until the recharger screen indicates a more efficient charging session.

Note: Charging efficiency is affected by factors such as implant depth and location, and efficiency will vary between patients. You and your clinician can work

Using **your** charging system 2

together to optimize your recharging efficiency.

6. Periodically check the progress of the charging session by reviewing the icons on your recharger screen.
7. When your neurostimulator has a charge sufficient for use, you will see the **Neurostimulator charge sufficient** screen (Figure 2.8)


Figure 2.8 Neurostimulator charge sufficient screen.

8. You can press the **Stop charge** key to stop your charging session, if desired. You will hear multiple tones if the audio is on. If you wait for the system to fully charge your neurostimulator, the **Neurostimulator charge complete** screen will appear on the recharger

display (Figure 2.9), and you will hear multiple tones, if the audio is on.


Figure 2.9 Neurostimulator charge complete screen.

9. Carefully release the antenna from the template or opening and remove the antenna from the belt or holster.
10. Remove the belt or holster and store it, the antenna, and the recharger in the charging system carrying case.

Accessories

Using the recharger pouch

The recharger comes with a small pouch that you can wear on a belt or holster to hold the recharger during a charging session.

Using your charging system 2

Using the charging system carrying case

The recharger comes with a carrying case to store the recharger and accessories. Using the case keeps the components clean and free from damage.

Identification label

Place an identification label on the back of your recharger in case the recharger is lost.

Additional features (available for some neurostimulator models)

Turning the neurostimulator on and off while charging

If you have an implantable neurostimulator for pain therapy, you can use the recharger to turn your neurostimulator on and off during a charging session. Use the two keys on the side of the recharger to turn your neurostimulator on and off.


Figure 2.10 Neurostimulator ON and OFF keys.

To turn the neurostimulator on while charging:

1. Press the **Neurostimulator ON**  key.
2. If audio is enabled, listen for the confirmation tone.
3. Verify that the **Neurostimulator ON** icon  is displayed on the recharger screen.

To turn the neurostimulator off while charging:

1. Press the **Neurostimulator OFF**  key.


Using your **charging** system 2

2. If audio is enabled, listen for the confirmation tone.
3. Verify that the **Neurostimulator OFF** icon  is displayed on the recharger screen.

Note: If you want to adjust stimulation using your patient programmer before your neurostimulator is completely charged, press the **Stop charge**  key on the recharger. You might need to remove the antenna so your patient programmer can communicate with your neurostimulator. Adjust stimulation with your patient programmer. Resume the charging session by pressing the **Start charge/test**  key.

Filename Date Time
UC200xxxxxx EN
4.625" x 6.0" inches (117 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.2: 03-27-2015

3 Troubleshooting

M924449A048 Rev A 2015-04-15

This section will help you solve problems with your charging system. It also provides information on when to call your clinician or Medtronic.

Note: If you cannot solve a problem, or if a problem is not described here, contact your clinician.

Recharger screens

The recharger displays **Warning** ⚠ and **Information** ⓘ screens to provide you with information about your system, alert you to a problem with your system, or to guide you during recharger use.

Warning screens

Warning screens ⚠ indicate a problem with the recharger, the antenna, or the neurostimulator. If the audio is on, multiple tones alert you to the screen. Table 3.1 describes **Warning** screens and provides instructions (**see blue text**) on how to resolve the problem.

△ **Caution:** Charge the neurostimulator when you see a **Low battery** screen displayed on the patient programmer or recharger. The battery charge level is at 0 and your therapy will stop soon unless you recharge the battery.

If you see the following screen (Figure 3.1), the neurostimulator battery is in a discharged state and your therapy has stopped.


Figure 3.1 Neurostimulator battery is in a discharged state.

You should charge the battery as soon as possible. If the battery remains in a discharged state, it continues to lose charge and it will become overdischarged.

If the neurostimulator battery is allowed to overdischarge, you may not be able to




recharge or communicate with the neurostimulator; however, your doctor may be able to restore the battery function.

Allowing the neurostimulator battery to overdischarge will permanently affect the neurostimulator in one of the following ways:

- Battery function is restored; however, charging sessions may be more frequent because battery capacity has been reduced.
- Battery function is not restored and the neurostimulator must be surgically replaced. Battery function is not restored because:
 - the neurostimulator battery is permanently damaged.
 - the neurostimulator battery has been overdischarged and restored twice before. The third time the battery is overdischarged, the neurostimulator will reach end of




service. Surgery is required to replace the neurostimulator.

Table 3.1 Recharger Warning screens

Warning screen	Cause and action
Charge your neurostimulator battery 	The neurostimulator battery charge level is depleted; the battery is in a discharged state. Therapy has stopped. Press the green Start charge/test key to clear the screen. Charge the neurostimulator battery now.
Charge your recharger battery 	Your recharger battery charge is insufficient to charge your neurostimulator. Connect the recharger to the AC power supply to clear the screen. Charge your recharger now.
Incompatible neurostimulator 	Your neurostimulator is incompatible with the recharger, and the recharger cannot be used to charge your neurostimulator. Call your clinician.

Troubleshooting 3

Table 3.1 Recharger Warning screens
(continued)

Warning screen	Cause and action
Call doctor 	Your neurostimulator has reached end of service. Stimulation has stopped. Call your doctor immediately.
	
Call clinician 	Your system is not working correctly. Stimulation might have stopped. Call your clinician and provide the specific error code shown on the recharger display.
	

Troubleshooting 3

Information screens

Information screens ⓘ provide information and instructions on how to proceed. If the audio is on, multiple tones alert you to the message. Press the **Audio** key on the recharger to clear the information screen. Table 3.2 describes information screens.

Table 3.2 Recharger Information screens

Information screen	Description and action
Reposition antenna	The communication between your recharger and your neurostimulator was interrupted or could not be established. Reposition your antenna over your neurostimulator, and press the green Start charge/test key.


Table 3.2 Recharger Information screens
(continued)


Information screen	Description and action
Check antenna 	The antenna connection to the recharger is loose or your antenna is not connected. Check that the antenna and recharger are making full contact. Connect or tighten the connection and retry. The antenna or recharger might require service. Call Medtronic. Refer to the addresses and phone numbers on the inside back cover of this manual.

Table 3.2 Recharger Information screens
(continued)



Information screen	Description and action
Move to a cooler environment; ensure your system is not too hot	Your antenna is too hot to work properly. Allow the antenna to cool. Store the charging system components in a cool, dry place.
	Your skin is too warm for a successful charging session.
	Move to a cool location, remove any excess clothing, allow your skin temperature to cool.
	The antenna may require service.
	Contact Medtronic. Refer to the addresses and phone numbers on the inside back cover of this manual.

Table 3.2 Recharger Information screens
(continued)


Information screen	Description and action
Note: This screen might appear during a recharging session, even if the antenna surface does not feel warm to you. Recharging will stop when this screen is displayed because the antenna is too warm for the system to work properly. Recharging will begin again once the surface temperature of the recharger antenna is within the operating range of the system.	
Physician recharger information 	An incorrect combination of keys was pressed and an error has occurred. The recharger cannot be used to charge your neurostimulator. Contact your physician.

Table 3.2 Recharger Information screens
(continued)



Information screen	Description and action
Diagnostics screens   	Diagnostics screens are used for troubleshooting only. These screens should only be accessed under the supervision of a Medtronic representative or your clinician. If you see a Diagnostics screen, use the following procedure to exit. Press the Stop charge key on the front of the recharger to clear this message. ⚠ Warning: Ensure that the key you press is the Stop charge key. If you accidentally press the Neurostimulator OFF key, you will turn off your stimulation. Use the patient programmer to check that your neurostimulator is on or off as desired. If you have any problems, contact Medtronic. Refer to the addresses and phone numbers inside the back cover of this manual.

Table 3.2 Recharger Information screens
(continued)










Information screen	Description and action
Charge your neurostimulator 	The neurostimulator battery charge level is low and stimulation will stop soon. Charge your neurostimulator battery using the recharger.
	
Charge your recharger 	The recharger battery charge level is low. Charge your recharger battery using the AC power supply.

Table 3.2 Recharger Information screens
(continued)

Information screen	Description and action
Finish neurostimulator charging session 	The neurostimulator battery is only partially charged. You did not completely charge the neurostimulator battery. Finish charging your neurostimulator battery using the recharger, if desired.
Finish recharger charging session 	The recharger battery is only partially charged. You did not completely charge the recharger battery. Finish charging your recharger battery using the AC power supply, if desired.
Neurostimulator charging complete 	The neurostimulator battery contains sufficient charge for use. Press the Audio key to clear the message.


Troubleshooting 3

Table 3.2 Recharger Information screens
(continued)

Information screen	Description and action
Recharger charging complete 	The recharger battery is fully charged. Press the Audio key to clear the message.
Call clinician 	There is a problem with a component of your system. Call your clinician and provide the specific error code shown on the recharger display.
Call clinician 	Your neurostimulator needs to be replaced, because stimulation will stop soon. Call your clinician and provide the message shown on the recharger display.

Troubleshooting 3

Table 3.2 Recharger Information screens
(continued)

Information screen	Description and action
POR 	Your neurostimulator has undergone a power-on-reset. Call your clinician and provide the message shown on the recharger display.
Reposition antenna  	Your neurostimulator can be turned on and off with the recharger. You tried to turn your neurostimulator on or off, but the recharger did not communicate successfully with the neurostimulator. Reposition your antenna over your neurostimulator. Try again to turn your neurostimulator on or off.

Possible problems and solutions

The following table will help you solve problems or identify when to call your

clinician. Problems are described in the left column (**bold black text**). The right column lists possible causes of the problem (plain text) and how to correct the problem (**bold blue text**).

Note: If a problem is not resolved after several attempts or if a problem is not described in the following table, contact your clinician.

Table 3.3 Troubleshooting

Problems	Causes and actions
Your recharger screen does not turn on when you try to charge your neurostimulator	The recharger battery is low. Attach the AC power supply to the recharger. The recharger is broken. Contact your doctor.
You spill fluid on your recharger, antenna, or AC power supply	If the AC power supply is plugged into an electrical outlet, interrupt the flow of electricity to the outlet. Unplug the power supply. Dry the component with a towel.

Troubleshooting 3

Table 3.3 Troubleshooting (continued)

Problems	Causes and actions
Your AC power supply will not turn on (the light on the supply does not turn green when you plug it in)	The power cord might not be plugged into the AC power supply or electrical outlet. Connect the power cord to the electrical outlet. The connection to the power cord might be loose. Reseat the electrical cord on the AC power supply. The power cord connected to the AC power supply might be incompatible with the electrical outlet. Obtain an adapter for the type of electrical outlet you are using.

Table 3.3 Troubleshooting (continued)

Problems	Causes and actions
Your recharger or antenna falls into water	1. If the AC power supply is plugged into an electrical outlet, interrupt the flow of electricity to the outlet. Unplug the power supply. 2. Remove the recharger or antenna from the water. 3. If the recharger or antenna is dirty, wipe it with a towel dampened with clean water. 4. Shake the recharger or antenna. 5. If you see or hear water inside, wipe it dry with a clean towel. Repeat steps 4 and 5. 6. When the recharger or antenna is dry, try to use it.

Troubleshooting 3

Table 3.3 Troubleshooting (continued)

Problems	Causes and actions
Your AC power supply falls into water	1. If the AC power supply is plugged into an electrical outlet, interrupt the flow of electricity to the outlet. Unplug the power supply. 2. Remove the AC power supply from water. 3. If it is dirty, wipe it with a towel dampened with clean water. 4. Leave the AC power supply on a dry towel for a few hours. 5. Shake the AC power supply. If you see or hear water inside, repeat steps 4 and 5 until the power supply is thoroughly dry. 6. Plug the power cord into a grounded electrical outlet. A green light on the top of the power supply will come on if it is working properly.

70 English 37751 2015-04-15

M924449A048 Rev A 2015-04-15

Table 3.3 Troubleshooting (continued)

Problems	Causes and actions
Your recharger, antenna, or AC power supply falls off a cabinet or table	The recharger, antenna, and AC power supply are designed to withstand falls of short distances. Try the recharger, antenna, or AC power supply; they should work.
Your recharger is beeping	You cannot resolve the problem or you see a screen other than those described in this section. Contact your clinician to resolve the problem.
Your recharger screen is blank	The Recharger wait screen did not appear when you pressed the Start charge/test key. Connect your recharger to the AC power supply and follow the directions for charging your recharger on page 78. Note: If the Recharger charging screen does not appear after you connect the recharger to the AC power supply, contact your clinician.

Troubleshooting 3

72 English 37751 2015-04-15

M924449A048 Rev A 2015-04-15

Table 3.4 Possible problems

Problems	Causes and actions
You experience heating sensation, discomfort, blistering not caused by heating, skin irritation, or redness near the implanted neurostimulator during or after recharging.	Your surgical wound is not fully healed. Stop recharging. If skin irritation persists, contact your clinician. The belt or holster is too tight and creating too much pressure on your skin from the antenna. Loosen the belt or holster so that it is comfortable and only tight enough to secure the antenna in place. You are sitting or lying on the antenna. Move to a position where you are not sitting or lying on the antenna. You experience skin irritation caused by direct contact with the antenna. Wear a thin undergarment or place a thin piece of fabric between your skin and the antenna. You may also want to try using more than one template to move the antenna further away from the skin.

Table 3.4 Possible problems (continued)

Problems	Causes and actions
You experience heating sensation, discomfort, blistering not caused by heating, skin irritation, or redness near the implanted neurostimulator during or after recharging. (continued)	You have been recharging for a long time without a break and are feeling discomfort. If skin redness occurs, take a break from recharging long enough for your skin to recover. During longer recharging sessions, take periodic breaks and remove the recharger from your skin. If any of these symptoms are severe, stop recharging and contact your clinician.
You experience uncomfortable warmth.	There is not enough air flow across the back of the antenna to keep it cool. Ensure that the back of the antenna is exposed to open air during recharging.

Troubleshooting 3

37751 2015-04-15 English 73

User assistance

Your recharger has been designed and tested to provide trouble-free service. If repair or service is needed, contact your clinician or a Medtronic sales office. Refer to the list of Medtronic contacts at the end of this manual.

The serial number is located on the back of the recharger. This number identifies each recharger. If you contact Medtronic about your recharger, refer to the serial number.

If your recharger stops working - First try the steps in Table 3.3. Otherwise, contact Medtronic.

If you lose your recharger - Contact your clinician to order a new recharger.

! USA To register the recharger for service covered by the warranty, complete and mail the warranty registration.

Filename Date Time
UC200xxxxxx EN
4.625" x 6.0" inches (117 mm x 152 mm)

Medtronic Confidential
PPManual.xsl - PatientProgrammerTemplate.fm
Template version 6.2: 03-27-2015

4 Maintenance

M924449A048 Rev A 2015-04-15

This section describes how to care for and dispose of your recharger and accessories.

Recharger maintenance tasks

Maintaining the battery in your recharger ensures that you will be able to charge your neurostimulator battery when needed. For information on charging and the batteries within the recharger, see "Battery and charging information" on page 34.

Checking the recharger battery

Press the gray **Audio** key on the recharger to turn on your recharger and check the recharger battery charge level. The **Recharger battery charge level** screen displays the current battery charge level (Figure 4.1).



Figure 4.1 Recharger battery charge level screen and icons.

Charge your recharger with the AC power supply when you see either of the screens in Table 4.1.

Table 4.1 When to charge the recharger

When you see this	It means
	The recharger battery charge level is low and will lose power soon. Charge the recharger battery. Press the Audio key to clear the message, and charge the recharger battery.

Table 4.1 When to charge the recharger
(continued)

When you see this	It means
	The recharger battery lost or is about to lose power. Charge the recharger battery immediately.

Charging your recharger battery

Follow the steps below to charge the recharger battery using the AC power supply.

1. Plug the AC power supply into an AC outlet (Figure 4.2). A green light indicates it has power.

 **Warning:** Do not use the AC power supply near water or during a lightning storm. Electric shock, serious injury, or death may result.


Figure 4.2 Plug in the AC power supply.

2. Connect the recharger to the AC power supply (Figure 4.3). The recharger will automatically start to charge when connected to the AC power supply, and the **Recharger charging session** screen will be displayed (Figure 4.4).

Notes:

- The charging session could take several hours, depending on the charge level of the batteries.

- To interrupt a recharger charging session, disconnect the recharger from the AC power supply.

Align triangles on both the cable and the recharger


Figure 4.3 Connect the recharger to the AC power supply.


Figure 4.4 Recharger charging session screen.

3. When the charging session is complete, the **Recharger charging complete** screen (Figure 4.5.) will appear.


Figure 4.5 Recharger charging complete screen.

4. Unplug the AC power supply from the electrical outlet.
5. Disconnect the AC power supply from the recharger.

6. Return the recharger, power cord, and AC power supply to the system carrying case.

Cleaning and Care

Follow these guidelines to ensure that the recharger system and accessories function properly:

- Keep the recharger, antenna, and AC power supply out of the reach of children and pets.
- Use the recharger, antenna, and AC power supply only as instructed by your clinician and this manual.
- Handle the recharger, antenna, and AC power supply with care. Do not drop, strike or step on the components.
- Do not dismantle or tamper with the recharger, antenna, or AC power supply.
- Clean the outside of the recharger, antenna, and AC power supply with a damp cloth when necessary. Mild

household cleaners will not damage the components.

- The recharger, antenna, and AC power supply are not waterproof. Do not allow moisture to get inside the components.
- The belt or holster, with the antenna removed, can be washed with a mild laundry detergent in a washing machine and dried in a dryer on a low heat setting.

Safety and technical checks

Periodic safety and technical checks or periodic maintenance of the recharger are not required.

The recharger contains no user-serviceable parts. If repair or service is needed, contact your clinician or a Medtronic sales office. Refer to the list of Medtronic contacts at the end of this manual.

Battery and device disposal

Dispose of worn out devices according to local regulations. Do not incinerate. If you no

longer need your recharger and would like to donate it, contact your clinician.

Declaration of conformity

Medtronic declares that this product is in conformity with the essential requirements of Directive 1999/5/EC on Radio and Telecommunications Terminal Equipment and Directive 90/385/EEC on Active Implantable Medical Devices.

For additional information, contact Medtronic. Refer to the list of Medtronic contacts at the end of this manual.

Specifications

The following table lists the specifications for your recharger.

Table 4.2 Recharger specifications

Parameter	Specification
Power source	2 non-removable lithium-ion batteries (rechargeable)

Maintenance 4

Table 4.2 Recharger specifications (continued)

Parameter	Specification
Operating temperature	10 °C to 35 °C (50 °F to 95 °F)
Temperature limitation	–20 °C to 40 °C (–4 °F to 104 °F)
Ingress protection	IP22 rating for solid objects greater than or equal to 12.5 mm, and for vertically dripping water when the device is tilted 15 degrees, per 60601-1-11.
Recharger size (approximate)	13 cm x 8 cm x 3.18 cm (5 in x 3 in x 1.25 in)
Recharger weight, including battery (approximate)	227 g (8 oz)
Battery life	5 years (approximate life of device) ^a
Inputs	0.85 A, 8.4 V DC (charge) 0.5 A, 10 V DC (base)

^a The longevity of your recharger battery is dependent upon how often you use your recharger.

Table 4.3 AC power supply and power cord specifications

Parameter	Specification
Operating temperature	10 °C to 35 °C (50 °F to 95 °F)
Temperature limitation	–20 °C to 40 °C (–4 °F to 104 °F)
Ingress protection	IP21 rating for solid objects greater than or equal to 12.5 mm, and for vertically dripping water, per 60601-1-11.
Power supply voltage	100–240 V (auto-ranging)
Rated frequency	50 Hz, 60 Hz
Input current	0.35 A
Outputs	0.85 A, 8.4 V DC (charge) 0.5 A, 10 V DC (base)
Mode of operation	Continuous

Index

Accessories 47
AC power supply
 description 26
 troubleshooting 68, 70, 71
Antenna
 description 27
 troubleshooting 69, 71
Antenna dial
 about 27
 adjusting 44
Audio control key 28
Audio tones 47, 52
Battery (neurostimulator)
 see Neurostimulator battery
Battery (recharger)
 see Recharger battery
Belt 27, 36
Care 82
Carrying case 27, 48
Charge efficiency row 31

Charge progress row 30
Charging efficiency 31, 44-46
Charging neurostimulator and recharger at
same time 34
Charging session
neurostimulator battery 43-47
recharger battery 79-81
Charging system
cleaning and care 82
parts 25-27
Charging the neurostimulator battery
procedure 40-47
turning neurostimulator on and off while
charging 48-50
Charging the recharger battery 78-81
Checking recharger battery 76-78
Cleaning 82
Dial (antenna)
about 27
adjusting 44
Disposing of devices 84
Holster 36

Index

88 English 37751 2015-04-15

M924449A048 Rev A 2015-04-15

Identification label 48
Information screens 57-65
Interrupting
 neurostimulator battery charging
 sessions 44
 recharger battery charging session 80
Keys, recharger 28
Low battery level screen
 neurostimulator 52
Maintenance 76-82
Neurostimulator battery
 battery level messages 37, 55
 charging 40-47
 checking status 36
 interrupting charging session 44
 see Neurostimulator
Neurostimulator battery (rechargeable)
 overdischarge 52
Neurostimulator
 charging information 34-35
 turning on and off while charging 48-50
 see Neurostimulator battery

Neurostimulator OFF key 49
Neurostimulator ON key 49
Overdischarge 52
Patient programmer 36
Pouch 27, 47
Power cord 26
Power supply
 see AC power supply
Recharger battery
 about 25
 charging 78-81
 checking charge level 76-78
 interrupting charging session 80
 maintaining 76-82
 See also Recharger 34
Recharger keys 28
Recharger
 charging information 34
 description 25
 possible problems 72
 troubleshooting 67-71
 see Recharger battery

Recharger screens

about 29-31

information screens 57-65

Neurostimulator charging session
screen 29

Recharger charging session screen 29

warning screens 52-56

Recharger system

see Charging system

Repair 83

Screens

see Recharger screens

Serial number 74

Specifications 84

Start charge/Test key 28

Status row 29

Stop charge key 28, 44

Template 41

Test (Start charge/Test) key 28

Tones, audio 47, 52

Troubleshooting 52-71
 assistance locating neurostimulator for
 charging 61
Turning neurostimulator on and off while
 charging 48-50
Warning screens 52-56
Warranty registration 74

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1. PATIENT PROGRAMMER MyStim, model 97740

Components	Material	Material contacts human tissue
Bottom Housing	PC/ABS	Yes (skin)
Battery Door	PC/ABS	Yes (skin)
Screen	PC	Yes (skin)
Button	Silicone	Yes (skin)
Battery Door Label	Clear Polycarbonate	Yes (skin)
Carrying case	Material	NEOPRENE
	Zippers	Black Plastic
	Velcro	Acrylic Copolymer

2. CHARGING SYSTEM, model 97754

Components	Material	Material contacts human tissue
Exterior Casing	PC/ABS	Yes (skin)
Belt	polyester	Yes (skin)
The power cord	Polyvinylchloride	No
Bag	Material	Nylon
	Zippers	Black Plastic
	Velco	Acrylic Copolymer
	Metal elements	Stainless Steel
The carrying case system	Material	NEOPRENE
	Snapper	Stainless Steel
	Velcro	Acrylic Copolymer
Spacer	PC/ABS	Yes (skin)


3. Detachable antenna, model 37092

Components	Material	Material contacts human tissue
Exterior Casing	PC/ABS	Yes (skin)

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
Item	Specification
Power source	2 AAA alkaline batteries (nonrechargeable, LR03) 1.5V voltage
Operating temperature	9°C to 43°C (49°F to 110°F)
Storage temperature	-34°C to 57°C (-30°F to 135°F)
Size (approximate)	9.4cm × 5.6 cm × 2.8 cm (3.7 in ± 0.1in × 2.2 in ± 0.1in × 1.1 in ± 0.1in)
Weight, including batteries (approximate)	111 g (3.9 oz) ± 0.5g
Battery life	2 months (average) for alkaline batteries
Mode of operation	Continuous
Input characteristics of the control panel channel, which provides reception of signals from the neuro stimulator	All inputs and outputs are conveyed via the IrDA connection, the Patient programmer user interface, or the telemetry interface
Characteristics of control panel signals used for programming of neurostimulants	The IrDA physical layer shall be consistent with the IrDA standard for signaling rates up to and including 38.4kb/s in low-power to low-power transmission mode
Parameters of acoustic sound signals are not specified when changing settings: force (loudness), acoustic frequency, duration	Force (loudness) = 70 dB Acoustic frequency / = 2700 Hertz Duration: 3 short beeps in 1 second during telemetry attempts 1 short beep when adjust STIM up or down.

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
2. CHARGING SYSTEM, model 97754

-Charger, model 37751 antenna

Parameter	Specification
Power source	2 non-removable lithium-ion batteries (rechargeable) - Voltage: 8.4 V at full charge - Capacity: Nominal Capacity = 170mAh - The charging time? Starting at 6.0 Volts, charge time to 8.4 V should be less than 4.5 hours
Operating temperature	10°C to 35°C (50°F to 95°F)
Storage temperature	-20°C to 40°C (-4°F to 104°F)

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Operating/storage relative humidity	30% to 95%
Operating/storage atmospheric pressure	70 kPa to 106 kPa (20.7 in Hg to 31.3 in Hg)
Recharger size (approximate)	13 cm (+/- 0,5 cm) x 8 cm (+/- 0,5 cm) x 3.18 cm (+/- 0,15 cm) (5 in ± 0.1in x 3 in ± 0.1in x 1.25 in ± 0.01in)
Rechargerweight, including battery (approximate)	227 g (8 oz) ± 1g
Battery life	5 years (approximate life of device)
Battery life	The longevity of your recharger battery is dependent upon how often you use your recharger.
Inputs	0.85 A, 8.4 V DC (charge) 0.5 A, 10 V DC (base)
Outputs: -Carrier frequency -Flux density of the antenna power radiated	175 khz 30 dBuA/m OOK burst

- Belt with template

Weight:	100 g +/- 1 g
Strap:	
Length	60 inches (152, 30 cm) +/- .5in
Width	2 inches (5.08 cm) +/- .5in
Mount:	
Height	.425 of an inch (1.08 cm) +/- .01in
Width	3.513 inches (8,92 cm) +/- .5in
Length	5.852 inch (14, 86cm) +/- .5in

- Spacer

Size (approximate)	3.66 inches +/- .01 in width 5.29 inches +/- .01 in lenght
Weight, (approximate)	25 g +/- 2 g

- Power supply, model 37761

Operating temperature	10°C to 35°C (50°F to 95°F)
Temperature limitation	-20°C to 40°C (-4°F to 104°F)
Power supply voltage	100-240 V (auto-ranging)
Rates frequency	50 Hz, 60 Hz
Input current	0.35 A
Outputs	0.85 A, 8.4 V DC (ccharge) 0.5 A, 10 V DC (base)
Mode of operation	Continuous
Diameter of the power cord	10mm +/- 1mm

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

Software:

Software Version:	not later 2.0C
Safety class (МЭК 62304):	Class C

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replacement surgery, help minimize potential kinking of the extension or lead, and minimize interference with telemetry.

- To use the nonsterile clinician programmer in a sterile field, place a sterile barrier between the patient and the programming head to prevent infection. Do not sterilize any part of the clinician programmer. Sterilization may damage the programmer.

Before explanting and replacing an existing neurostimulator, print a report from the existing neurostimulator that shows implanted and any abandoned component information. When adding, changing, or removing neurostimulators, leads, extensions, and accessories, always program up-to-date component model numbers, implant locations, and any abandoned component information to re-establish MRI scan eligibility.

If this information is not updated or is entered incorrectly, MRI conditionally safe (MRICS) scan eligibility data will be inaccurate, and the patient is at risk for one of the following:

- The patient is allowed to have an MRI scan inappropriate for the implanted components, which could cause tissue heating, resulting in tissue damage or serious patient injury.

- The patient is unnecessarily restricted from having an MRI scan.

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All Rights Reserved