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Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG. hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed:

Date: 18-Oct-2018

Zaman Khan
Associate Director Regulatory
Affairs EMEA and Pakistan
Abbott GmbH & Co. KG

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
Alanine Aminotransferase Reagent Kit

en
ALT
07P98
G71193R04
B7P980

Read Highlighted Changes: Revised February 2018.

REF 07P9820

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NAME

Alinity c Alanine Aminotransferase Reagent Kit (also referred to as ALT)

INTENDED USE

The Alinity c Alanine Aminotransferase assay is used for the quantitation of alanine aminotransferase in human serum or plasma on the Alinity c analyzer.

SUMMARY AND EXPLANATION OF THE TEST

Alanine Aminotransferase (ALT), also referred to as glutamate pyruvate transaminase (GPT), is an enzyme involved in amino acid metabolism. It is found in many tissues, but the highest levels are found in liver and kidney tissues. Tissue destruction leads to the release of the intracellular enzyme into the circulating blood. Markedly elevated serum ALT levels may be found in a variety of diseases which involve the liver, such as hepatitis, mononucleosis, and cirrhosis. These very high levels of ALT are not usually observed in other disease processes, e.g., myocardial infarction; thus, ALT is regarded as a reasonably specific indicator of liver disease.

PRINCIPLES OF THE PROCEDURE

ALT present in the sample catalyzes the transfer of the amino group from L-alanine to α -ketoglutarate, forming pyruvate and L-glutamate. Pyruvate in the presence of NADH and lactate dehydrogenase (LD) is reduced to L-lactate. In this reaction NADH is oxidized to NAD. The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due to the oxidation of NADH to NAD.

Methodology: Enzymatic: NADH (without P-5-P)

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity c Alanine Aminotransferase Reagent Kit 07P98

Volumes (mL) listed in the table below indicate the volume per cartridge.

REF	07P9820
Tests per cartridge	360
Number of cartridges per kit	10
Tests per kit	3600
R1	68.1 mL
R2	21.0 mL
1 Active Ingredients: β -NADH (0.16 mg/mL), Lactate dehydrogenase (2.57 U/mL), L-Alanine (392 mmol/L).	
2 Active Ingredients: α -Ketoglutaric acid (77 mmol/L), L-Alanine (1000 mmol/L).	

Warnings and Precautions

IVD

For In Vitro Diagnostic Use

Rx ONLY

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴

The following warnings and precautions apply to R2:	
Contains tris hydroxymethyl aminomethane.	
H316	Causes mild skin irritation.
Response	
P332+P313	If skin irritation occurs: Get medical advice / attention.

* Not applicable where regulation EU 1272/2008 (CLP) or OSHA Hazard Communication 29CFR 1910.1200 (HCS) 2012 have been implemented.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2-8°C	Until expiration date	Store in upright position.
Onboard	System Temperature	27 days	
Opened	2-8°C	Until expiration date	Store in upright position. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance.

Abbott

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity c Alanine Aminotransferase assay file must be installed on the Alinity c analyzer prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay. Other specimen types, collection tube types, and anticoagulants have not been verified with this assay.

Specimen Types	Collection Tubes	Special Conditions
Serum	Serum tubes (with or without gel barriers)	
Plasma	Collection tubes Acceptable anticoagulants are: Lithium heparin (with or without gel barrier) Sodium heparin EDTA	Do not use ammonium heparin. ⁵

Hemolysis in serum or plasma can increase test results.

CAUTION: Erythrocytes contain approximately 3 to 5 times more ALT than does serum.⁶

- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- For accurate results, plasma specimens should be free of platelets and other particulate matter. Ensure centrifugation is adequate to remove platelets.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross-contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Specimen Storage

Numerous publications have defined storage conditions for ALT.⁶⁻¹⁵ Examples are shown below.

Specimen Type	Temperature	Maximum Storage Time	Special Instructions
Serum/ Plasma	30°C	3 days ¹⁰	Remove serum or plasma from the clot, red blood cells, or separator gel.
	2 to 8°C	7 days ¹⁰	
	-40°C	60 days ¹⁵	

It is recommended that specimens be assayed on the day of collection.^{16, 17}

When samples were stored at -20°C for 8 days, an 11% reduction in ALT activity was observed; a 20% reduction in ALT activity was observed when specimens were stored at -20°C for 1 month.¹⁸

Avoid multiple freeze/thaw cycles.

Stored specimens must be inspected for particulates. If present, mix with a low speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PROCEDURE

Materials Provided

07P98 Alinity c Alanine Aminotransferase Reagent Kit

Materials Required but not Provided

- Alinity c Alanine Aminotransferase assay file
- Commercially available controls containing alanine aminotransferase
- Saline (0.85% to 0.90% NaCl) for specimen dilution

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

Ass

For
ci-se

- I
- J
- K
- L
- M
- N
- O
- P
- Q
- R
- S
- T
- U
- V
- W
- X
- Y
- Z

- Re
- in:
- Fo
- Op
- For
- ma
- Ma
- req

Sampl

Sample
are flag,
either th
Procedu

Automat

The syst
calculate
factor.

Manual I

Dilute the
The opera
Control ta
dilution fa
sample an
If the opera
manually n
the result,
the measur
an appropri
For detail
series Ope

Calibratio

For instruct
series Oper
Calibration
is required
at least 2 le
control requ
acceptable
This assay r
parts or sub
performed.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
 - Sample volume for single test: 5.3 µL.

NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the Alinity ci-series Operations Manual, Section 4.
- Refer to the commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples with an alanine aminotransferase value exceeding 3899 U/L are flagged with the code "> 3899 U/L" and may be diluted with either the Automated Dilution Protocol or the Manual Dilution Procedure.

Automated Dilution Protocol

The system performs a 1:5 dilution of the sample and automatically calculates the concentration by multiplying the result by the dilution factor.

Manual Dilution Procedure

Dilute the sample with saline (0.85% to 0.90% NaCl).

The operator must enter the dilution factor in the Specimen or Control tab of the Create Order screen. The system will use this dilution factor to automatically calculate the concentration of the sample and report the result.

If the operator does not enter the dilution factor, the result must be manually multiplied by the appropriate dilution factor before reporting the result. If a diluted sample result is less than the lower value of the measuring interval of 5 U/L, do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration is stable for approximately 27 days (648 hours), but is required with each change in reagent lot. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Quality Control Procedures

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Two levels of controls (normal and abnormal) are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Commercial controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

Quality Control Guidance

Refer to "Basic QC Practices" by James O Westgard, Ph.D. for guidance on laboratory quality control practices.¹⁹

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

RESULTS

Calculation

The Alinity c Alanine Aminotransferase assay utilizes the Factor data reduction method to generate a calibration curve and results. The calibration factor for the Alinity c Alanine Aminotransferase assay is 8141.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

Measuring Interval

Measuring Interval is defined as the range of values in U/L which meets the limits of acceptable performance for linearity, imprecision, and bias.

The measuring interval of the Alinity c Alanine Aminotransferase assay is 5 to 3899 U/L.

LIMITATIONS OF THE PROCEDURE

Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert.

EXPECTED VALUES

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

Reference Range

Serum/Plasma^{20, 21}

	Range (U/L)
Adult	0 to 55

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

The Alinity c analyzer, and the ARCHITECT c System and AEROSET System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity c analyzer.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2. Testing was conducted using 1 lot of Alinity c Alanine Aminotransferase Reagent Kit, water calibrator, and 1 lot of commercially available controls, and 1 instrument. Three controls were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.²²

Sample	n	Mean (U/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
			SD	%CV	SD	%CV
Control Level 1	120	30	0.5	1.8	0.9	2.8
Control Level 2	120	108	0.7	0.6	1.0	1.0
Control Level 3	120	230	0.7	0.3	1.4	0.6

^a Includes within-run, between-run, and between-day variability.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2 and was replicated on 3 reagent lots and 2 instruments over a minimum of 3 days on each instrument/reagent lot combination. The maximum observed Limit of Blank (LoB), Limit of Detection (LoD) and Limit of Quantitation (LoQ) values are reported in the table.²³

	U/L
LoB ^a	1
LoD ^b	2
LoQ ^c	5

^a The LoB represents the 95th percentile from $n \geq 60$ replicates of zero-analyte samples.

^b The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on $n \geq 60$ replicates of low-analyte level samples.

^c The LoQ was determined from $n \geq 60$ replicates of low-analyte level samples and is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met.

Linearity

A study was performed based on guidance from CLSI EP06-A.²⁴

This assay is linear across the measuring interval of 5 to 3899 U/L.

Interference

This study was performed on the AEROSET System.

Potentially Interfering Endogenous Substances

Interference studies were conducted using NCCLS EP7-P.²⁵

Interference effects were assessed by Dose Response and Paired Difference methods, at the medical decision level of the analyte.

Potentially Interfering Substance	Interferent Level		Target U/L	Recovery (% of Target)
	Default Units	Alternate Units		
Bilirubin	30 mg/dL	513 μ mol/L	53.1	95.3
	60 mg/dL	1026 μ mol/L	53.1	88.1
Hemoglobin	750 mg/dL	7.5 g/L	47.4	107.9
	1000 mg/dL	10.0 g/L	47.4	111.0
Intralipid	550 mg/dL	5.5 g/L	50.6	97.2
	625 mg/dL	6.25 g/L	50.6	96.8

The following drugs were tested on the ARCHITECT system for interference at the concentrations indicated using an acceptance criteria of $\pm 10\%$ from the target value.

Potentially Interfering Substance	Interferent Level		Target U/L	Recovery (% of Target)
	Default Units	Alternate Units		
Sulfapyridine	60 mg/L	241.0 μ mol/L	43.1	95.1
Sulfasalazine	20 mg/L	50.3 μ mol/L	43.1	92.8
Tamoxifen	20 mg/L	103.1 μ mol/L	58.2	105.0

Interferences from medications or endogenous substances may affect results.²⁶

Method Comparison

A study was performed based on guidance from CLSI EP9-A3 using the Passing-Bablok regression method.²⁷

	Units	n	Correlation Coefficient			Concentration Range	
			Intercept	Slope			
Alinity c Alanine Aminotransferase vs ARCHITECT Alanine Aminotransferase	Serum U/L	130	1.00	1.11	0.97	6-3727	

BIBLIOGRAPHY

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 2nd ed. Philadelphia, PA: WB Saunders; 1994:795-797.
- Henry RJ, Cannon DC, Winkelman JW. *Clinical Chemistry Principles and Techniques*, 2nd ed. Hagerstown, MD: Harper and Row; 1974:888.
- Ruby SG, Reiber NE, Lonser RE. Preanalytical variation in alanine aminotransferase. *Clin Chem* 1988;34(4):744-745.
- Heins M, Heil W, Withold W. Storage of serum and whole blood samples? Effects of time and temperature on 22 serum analytes. *Eur J Clin Chem Clin Biochem* 1995;33:231-238.
- Dale JC, Pruett SK. Phlebotomy—a minimalist approach. *Mayo Clin Proc* 1993;68(3):249-255.
- Young D. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press; 1997:3-12.
- Elfath D, Cooney J, McDaniel R, et al. Effect of frozen storage of serum on the level of 22 chemistry analytes [Poster 0102 presented at AACC 43rd National Meeting in Washington, DC: July 30, 1991]. *Clin Chem* 1991;37(6):931.
- Faulkner AM, Lukes-Hall AM, White GW. Evaluation of the Grenier plasma separator blood tube. *Ann Clin Biochem* 1990;27:386-387.
- Wilding P, Silva JA, Wilde CE. Transport of specimens for clinical chemistry analysis. *Ann Clin Biochem* 1977;14:301-306.
- Schwartz MK. Interferences in diagnostic biochemical procedures. *Adv Clin Chem* 1973;16:10.
- Mosley JW, Goodwin RF. Stability of serum glutamic pyruvic transaminase activity on storage. *Am J Clin Pathol* 1965;44:591-595.
- Williams K, Williams A, Kline L, et al. Stability of serum alanine aminotransferase activity. *Transfusion* 1987;27(5):431-433.
- Cuccherini B, Nussbaum S, Seeff L, et al. Stability of aspartate aminotransferase and alanine aminotransferase activities. *J Lab Clin Med* 1983;102(3):370-376.
- Donnelly JG, Soldin SJ, Nealon DA, et al. Stability of twenty-five analytes in human serum at 22°C, 4°C, and -20°C. *Pediatr Pathol Lab Med* 1995;15:869-874.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Kaplan LA, Pesce AJ, editors. *Clinical Chemistry Theory, Analysis, and Correlation*, 2nd ed. St Louis, MO: CV Mosby; 1989:895-898.
- Sherman KE, Dodd RY, et al. Alanine aminotransferase levels among volunteer blood donors: geographic variation and risk factors. *J Infect Dis* 1982;145(3):383-386.

22. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
23. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
24. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
25. National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
26. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3-16–3-22.
27. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

■ Key to Symbols

ISO 15223 Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
	In Vitro Diagnostic Medical Device
	Lot Number
	List Number
	Serial number

Other Symbols

	Distributed in the USA by
	Information needed for United States of America only
	Produced for Abbott by
	Product of Canada
	Reagent 1
	Reagent 2
	For use by or on the order of a physician only (applicable to USA classification only).

Alinity, ARCHITECT, and AEROSSET are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown
Prince Edward Island
C1E 2B9 Canada

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60054 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised February 2018.

©2016, 2018 Abbott Laboratories

ADDENDUM

to the Instruction for Use "Alinity c Alanine Aminotransferase" (for Russia)

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for device is 2-8°C. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

Medical device must be stored at a temperature of 2°C to 8°C in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

SHELF LIFE

18 months at the recommended storage condition (2 to 8°C).

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. Manufacturer set the warranty shelf life of device which corresponds to period between the date of manufacture and expiration date that are indicated on the label and/or in accompanying documentation. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required

REPAIRS

Not required

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-11 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste". Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

ISO 9001, ISO 13485, ISO 13485 (including CMDCAS)

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

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1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
5. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 2nd ed. Philadelphia, PA: WB Saunders; 1994:795–797.
6. Henry RJ, Cannon DC, Winkelman JW. *Clinical Chemistry Principles and Technics*, 2nd ed. Hagerstown, MD: Harper and Row; 1974:888.
7. Ruby SG, Reiber NE, Lonser RE. Preanalytical variation in alanine aminotransferase. *Clin Chem* 1988;34(4):744–745.
8. Heins M, Heil W, Withold W. Storage of serum and whole blood samples? Effects of time and temperature on 22 serum analytes. *Eur J Clin Chem Clin Biochem* 1995;33:231–238.
9. Dale JC, Pruett SK. Phlebotomy—a minimalist approach. *Mayo Clin Proc* 1993;68(3):249–255.
10. Young D. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press; 1997:3–12.
11. Elfath D, Cooney J, McDaniel R, et al. Effect of frozen storage of serum on the level of 22 chemistry analytes [Poster 0102 presented at AACC 43rd National Meeting in Washington, DC: July 30, 1991]. *Clin Chem* 1991;37(6):931.
12. Faulkner AM, Lukes-Hall AM, White GW. Evaluation of the Grenier plasma separator blood tube. *Ann Clin Biochem* 1990;27:386–387.
13. Wilding P, Silva JA, Wilde GE. Transport of specimens for clinical chemistry analysis. *Ann Clin Biochem* 1977;14:301–306.
14. Schwartz MK. Interferences in diagnostic biochemical procedures. *Adv Clin Chem* 1973;16:10.
15. Mosley JW, Goodwin RF. Stability of serum glutamic pyruvic transaminase activity on storage. *Am J Clin Pathol* 1965;44:591–595.
16. Williams K, Williams A, Kline L, et al. Stability of serum alanine aminotransferase activity. *Transfusion* 1987;27(5):431–433.
17. Cuccherini B, Nussbaum S, Seeff L, et al. Stability of aspartate aminotransferase and alanine aminotransferase activities. *J Lab Clin Med* 1983;102(3):370–376.
18. Donnelly JG, Soldin SJ, Nealon DA, et al. Stability of twenty-five analytes in human serum at 22°C, 4°C, and -20°C. *Pediatr Pathol Lab Med* 1995;15:869–874.
19. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
20. Kaplan LA, Pesce AJ, editors. *Clinical Chemistry Theory, Analysis, and Correlation*, 2nd ed. St Louis, MO: CV Mosby; 1989:895–898.
21. Sherman KE, Dodd RY, et al. Alanine aminotransferase levels among volunteer blood donors: geographic variation and risk factors. *J Infect Dis* 1982;145(3):383–386.
22. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
23. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
24. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
25. National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Proposed Guideline*. NCCLS Document EP7-P. Villanova, PA: NCCLS; 1986.
26. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3–16–3–22.
27. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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