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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

01271929

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Translation from English into Russian

From ZHEJIANG BIOREAR BIOTECHNOLOGY CO., LTD. ("ZHEJIANG BIOREAR BIOTECHNOLOGY CO., LTD."), I certify that the document provided is accurate, correct and valid.

Approved: October 28, 2022

Signature/Name:  (Bian Li)

Position: Vice President

INSTRUCTIONS FOR USE OF THE MEDICAL DEVICE

Deep well plate 96 wells, colorless, polypropylene, 2.2 ml volume, V -bottom, non-sterile, compatible with automated nucleic acid extraction systems , in the composition

Version 28102022 dated October 28, 2022

1. Name and composition of the medical device:

Deep well plate 96 wells, colorless, polypropylene, 2.2 ml volume, V-shaped bottom, non-sterile, compatible with automated nucleic acid extraction systems, in the composition

- deep well plate - 10 packs of 5 pieces per box;
- instructions for use - 1 pc.

2. Indications and contraindications for the use of a medical device.

Indications for use

" Deep well plate 96 wells, colorless, polypropylene, 2.2 ml, V -bottom, non-sterile, compatible with automated nucleic acid extraction systems, in the composition" is to be used once when carrying out automatic extraction of nucleic acids from any type of human biological material, regardless from the analyzed analyte .

Contraindications for use, possible side effects

There are no contraindications for use, possible side effects within the established purpose.

3. Information about potential consumers of the medical product.

Potential consumers are the personnel of clinical diagnostic laboratories (specialists in the field of clinical laboratory diagnostics). For professional use only.

4. Specifications. Description of the main functional elements.

Production material	polypropylene
Color	colorless _
Number of holes, pcs.	96
Hole bottom type	V -shaped
Nominal well volume, ml	2.2
non- sterile	

The product is biologically inert. Free from DNA, RNases , DNases , PCR inhibitors.

The dimensions of the plate are shown in Figure 1.


Figure 1 - Dimensions of the plate, mm

5. Functional purpose, principles of operation and application procedure

Intended use

" Deep well plate 96 wells, colorless, polypropylene, 2.2 ml, V-shaped bottom, non-sterile, compatible with automated nucleic acid extraction systems, in the composition" is designed to contain liquids (reagents and samples of biological material) of the required volume at isolation of nucleic acids using compatible automated systems for purification of nucleic acids in the process of simultaneous multichannel processing of biomaterial (including magnetic separation). The product is intended for single use by a professional user when performing diagnostic procedures in vitro .

Functional use

Deep-well plate 96 wells, colorless, polypropylene, 2.2 ml, V -shaped bottom, non-sterile, compatible with automated nucleic acid isolation systems, included, used in clinical diagnostic laboratories as a consumable (auxiliary) for diagnostics in vitro .

Operating principle

The success of the subsequent analysis depends on the qualitative and reproducible purification of nucleic acids, proteins and cells. Nucleic acid extraction systems provide quality results with minimal manual preparation time, helping to automate much of your workflow.

Deep well plates are resistant to interactions with biomolecules and are ideal for nucleic acid isolation, including for working with magnetic particles. When using plates and tip combs specially designed for nucleic acid extraction systems, the

quality of protein and nucleic acid extraction and the yield of the target product are significantly improved (see Figure 2).

The deep-well plate contains 96 wells and allows you to place the required amount of biomaterial in it. The tablet material is chemically resistant to alcohols and organic solvents.

When working with magnetic particles, plates are used in conjunction with a comb of tips for magnets with deep wells. In such nucleic acid extraction systems, magnetic rods with disposable tips collect, transfer, and mix the magnetic particles.


Figure 2 - Principle of operation

Application procedure.

deep well plates is given in the operating documentation of compatible systems for nucleic acid extraction.

Before installation in the system holders, inspect the product for any damage.

6. Operational / functional characteristics and methods of control.

6.1 Compatibility of plates with the corresponding positions (nests) of the nucleic acid isolation system :

The plates should be freely installed in the appropriate position (socket) of the nucleic acid extraction system and located in it (it) without a visible slope and displacement.

Compliance is assessed visually in 8-10 repetitions when running on a compatible nucleic acid isolation system.

6.2 Mutual geometric compatibility of deep well plates and tip combs:

The tip combs should be free to dip into the appropriate wells of the plate. The comb of tips placed in the plate should be located in it without a visible slope and displacement. The comb tips should be free to remove from the wells of the plate.

Compliance is assessed visually by aligning plates and tip combs by testing 8-10 pieces of each unit.

Deep well plates are geometrically compatible with tip combs having the dimensions shown in Figure 3:


Figure 3 - Dimensional characteristics of compatible tip combs

6.3 Preservation of the integrity of the walls of the deep well plate during the extraction of nucleic acids:

At the end of the release cycle, there should be no traces of liquid visible to the eye on the outer surface of the tablet. After 3-3.5 minutes after placing the tablet on a sheet of filter paper, there should be no traces of liquid.

Compliance is assessed visually in 3 repetitions.

7. Marking .

On the labeling of packages (primary and transport), symbols are used, the decoding of which is presented in table 1.

Table 1

Symbol name	Symbol
Batch code	LOT

Symbol name	Symbol
Use up to	
Manufacturer (name and address of manufacturer)	
Date of manufacture	
Medical device for diagnostics in vitro	
Reuse prohibition	
Do not use if packaging is damaged	
Protect from moisture / Keep dry	
Avoid exposure to sunlight / Keep away from heat sources	
Temperature limit	

8. Expiry date.

The manufacturer guarantees that the products meet the requirements of the technical documentation. Warranty period - 36 months from the date of manufacture.

9. Transportation and storage.

9.1 The product is stored in a transport package in closed, clean and dry warehouses protected from moisture and direct sunlight, at a distance of at least 1 m from heating devices on racks at least 5 cm from the floor in conditions that exclude the impact of aggressive environments (acidic, alkaline, etc.), as well as flammable and combustible liquids at a temperature of 5 to 30 ° C in a dry place.

9.2 The product is transported by all means of transport in covered vehicles in accordance with the rules for the carriage of goods in force for this type of transport.

9.3 During storage and transportation, it is not recommended to stack transport packages on top of each other.

9.4 It is not allowed to transport products without transport packaging.

9.5 Environmental conditions during transportation:

- air temperature from minus 40°C to 60°C,
- relative air humidity $\leq 95 \%$;
- atmospheric pressure 84.0-106.7 kPa.

10. Disposal.

Used products are classified as epidemiologically hazardous waste. Unused products are epidemiologically safe waste. Disposal must be carried out in accordance with the rules in force in the country of the consumer, observing all established restrictions at the time of disposal.

In the Russian Federation, the product is disposed of in accordance with the requirements of SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures. Waste should be disposed of through organizations licensed for this type of activity.

11. Instructions for use.

The use of deep well plates should be carried out in accordance with the instructions for use of the device itself and compatible systems for nucleic acid extraction, for example, KingFisher magnetic particle processor for purification of nucleic acids, cells and proteins (FSZ 2009/05562), Microlab STARlet automatic dosing station with accessories (RZN 2018/6981), Freedom EVO automatic modular analyzer with accessories (FSZ 2008/03047).

Permissible operating conditions:

Air temperature from 5 to 30°C.

The order of transportation, storage and use, set out in the instructions for use, must be observed.

12. Requirements for the maintenance and repair of a medical device .

Deep well plate 96 wells, colorless, polypropylene, volume 2.2 ml, V -shaped bottom, non-sterile, compatible with automated nucleic acid extraction systems, as a part, is a single-use product, not subject to maintenance.

13. Applied standards.

Standard	Name of the standard
ISO 13485	Medical products. Quality management systems. Requirements for Regulatory Purposes
ISO 14971	Medical products. Application of risk management to medical devices
ISO 15223-1	Medical products. Symbols used with information provided by the manufacturer. Part 1. General requirements
EN IS O 23640	Medical products for in vitro diagnostics. Evaluation of the stability of reagents for in vitro diagnostics
EC Directive 98/79/EC	Directive 98/79/ EC of the European Parliament and of the Council of 27 October 1998 on medical devices for in -diagnosis vitro

14. Medical device manufacturer:

ZHEJIANG BIOREAR BIOTECHNOLOGY CO., LTD.
("ZHEJIANG BIOREAR BIOTECHNOLOGY CO., LTD.")

Address of the manufacturer of the medical device:

1st _ Floor , Building 2, No.3535 , Linggongtang Road , Nanhu district , Jiaxing , Zhejiang , China / 1 Floor, Building 2, No .3535, Linggongtan Road, Nanhu District, Jiaxing, Zhejiang, China

Phone: + 86 057382618891

Email: sales@biorear.com

15. Authorized representative of the manufacturer:

Cosmofarm Limited Liability Company (Cosmofarm LLC) ,

Address of the authorized representative of the manufacturer of the medical device:

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