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(ETHICON, LLC)

Manager Global Regulatory Affairs

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

Instruction for Use
PROLENE Polypropylene Mesh in variants

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INSTRUCTION FOR USE

PROLENE Polypropylene Mesh in variants:

PROLENE Polypropylene Mesh: 6 x 11 cm;
 PROLENE Polypropylene Mesh: 30 x 30 cm;
 PROLENE Polypropylene Mesh: 15 x 15 cm;
 PROLENE Polypropylene Mesh: 15 x 10 cm;
 PROLENE Soft Polypropylene Mesh: 2,5 x 10,2 cm;
 PROLENE Soft Polypropylene Mesh: 6,4 x 11,4 cm;
 PROLENE Soft Polypropylene Mesh: 7,6 x 15,2 cm;
 PROLENE Soft Polypropylene Mesh: 15,2 x 15,2 cm;
 PROLENE Soft Polypropylene Mesh: 25,4 x 25,4 cm;
 PROLENE Soft Polypropylene Mesh: 30,5 x 35,6 cm;
 PROLENE Soft Polypropylene Mesh: 50 x 50 cm;

Hereinafter may be referred to as mesh, PROLENE Mesh, PROLENE Polypropylene Mesh, product, device, MD.

INTRODUCTION

This package insert is designed to provide instructions for use of PROLENE Mesh.

It is not a comprehensive reference to surgical technique for repair of abdominal wall hernias.

PROLENE Mesh is intended for use only by physicians who are trained in the surgical procedures and techniques required for hernia repairs and the implantation of synthetic meshes. The selection of mesh for any given patient is a function of numerous factors including, but not limited to, the patient's past medical and surgical history, current medical condition (i.e., comorbidities), surgical technique, and size and location of the defect hernia. The physician is advised to consult the medical literature regarding techniques, complications, and adverse reactions before selecting a mesh.

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DESCRIPTION

PROLENE Mesh is comprised of polypropylene filaments. The polypropylene polymer is identical to the material used in PROLENE Suture.

PROLENE Soft Mesh is comprised of nonabsorbable polypropylene filaments. The undyed and dyed (phthalocyanine blue, Color Index No. 74160) polypropylene polymers are identical to the material used in PROLENE Suture.

INDICATIONS

PROLENE Mesh is indicated for the repair of abdominal wall hernias and abdominal wall deficiencies that require the addition of a reinforcing material to obtain the desired surgical result.

CONTRAINDICATIONS

PROLENE Mesh should not be used intraperitoneally.

WARNINGS

- Failure to properly follow instructions may result in improper functioning of the device and could lead to injury. Please read all information carefully.
- If this device is used in patients with the potential for growth or tissue expansion (such as infants, children, or women who may become pregnant), the surgeon should be aware that the device will not stretch significantly as the patient grows.
- This device is indicated for abdominal wall hernia repair and not for gynecologic procedures. Gynecologic procedures should be performed using devices indicated for gynecologic repairs.
- It is recommended that the device not be used in a contaminated field, because contamination of the device may lead to infection that may require removal of the device.
- As with any implant, an acute and permanent foreign body response will occur. In some patients, this response can result in one or more of the adverse reactions listed below.
- The device is a permanent implant that is designed to integrate into the tissue. In cases in which the device needs to be removed, in part or in whole, significant dissection may be required.
- Insufficient overlap on any side of the defect may increase the risk of

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postoperative complications, including recurrence. Consult Application/Instructions for Use section.

- Insufficient or improper fixation may increase the risk of postoperative complications, including recurrence. Consult Application/Instructions for Use section.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of bloodborne pathogens to patients and anyone coming in contact with the device.
- Inspect the mesh carefully before implantation. Do not use the device if it is damaged.

PRECAUTIONS

- The safety and effectiveness of pretreating PROLENE Surgical Mesh with solutions (e.g., saline, medications) prior to implantation have not been studied.

ADVERSE REACTIONS

Potential adverse reactions are those typically associated with surgically implantable materials, including infection, inflammation, seroma formation, acute or chronic pain, foreign body sensation, hematoma, nerve damage, soft tissue injury, adhesion formation, fistula formation, extrusion/erosion, excessive contraction or shrinkage of the tissue surrounding the mesh, and mesh failure/hernia recurrence.

One or more revision surgeries may be necessary to treat the above-mentioned adverse reactions. Revision surgery may not resolve the adverse reactions and may pose a risk of additional adverse reactions.

APPLICATION / INSTRUCTIONS FOR USE

Sizing and Placement

PROLENE Mesh must always be placed extraperitoneally. PROLENE Mesh can be trimmed according to the surgeon's discretion, while providing the necessary overlap to reduce the likelihood of recurrence. The use of thermal

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cutting devices is not recommended, because it has not been tested.

Ventral/Incisional and Inguinal Hernia Repair

In ventral or incisional hernia repair, the mesh should extend at least 3 to 5 cm beyond the margins of the defect, unless at the surgeon's discretion, additional overlap onto healthy tissue is needed. When used in inguinal hernia repair, the mesh should provide sufficient overlap of the fascial defect on all sides. The colored stripes (PROLENE Soft Mesh) on the mesh can be used for orientation and alignment purposes. Place the mesh so that it lies flat to the tissue.

Fixation

The method of securing the implant (e.g., nonabsorbable/absorbable sutures or tackers) to provide for adequate tissue fixation and to reduce the risk of recurrence should be determined at the surgeon's discretion, based on the individual patient's needs and the nature of the repair.

Spacing and distribution between fixation points and technique should be determined at the surgeon's discretion to provide adequate mesh fixation, to reduce the risk of mesh migration, and to optimize mesh-to-tissue contact to foster tissue ingrowth.

Preclinical data and reported experience suggest having the fixation points at least 1 cm from the edge of the mesh.

PERFORMANCE / ACTIONS

An animal study showed that implantation of PROLENE Mesh elicits an inflammatory reaction that stimulates the deposition of a thin fibrous layer of tissue, which can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The material is not absorbed.

STERILITY

PROLENE Mesh is sterilized using ethylene oxide gas. Do not resterilize. Do not use if the package is opened or damaged.

STORAGE

No special conditions required. Do not use after expiry date.

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Transportation conditions:

Temperature: Do not exceed -60 or + 60 °C

Relative humidity: up to 100% (non-condensing at 25 °C)

Operation conditions:

The medical device is intended for use at the standard conditions of hospital operating room environment

The conditions in the rooms where the medical device is used must meet the following requirements:

Temperature from -10°C to 50°C (75% non-condensing humidity equivalent to 30°C) when utilized immediately after opening sterile packaging.

HOW SUPPLIED

PROLENE Mesh is available in single-use, sterile packets in a variety of sizes. Each box contains instructions for use

Product Code	Packaging Piece/s per box
PML1	1
PMM3	3
PMN3	3
PMS3	3
SPMXS	6
SPMS	6
SPMΠ	6
SPMII	6
SPMLI	3
SPMXXL	3
SPM3XL	3

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TRACEABILITY

A traceability label that identifies the type, size, expiry date and batch code of the device is included in every package. This label should be affixed or electronically added to the patient's permanent medical record to clearly identify the device that was implanted.

Interpretation of symbols used on the package.

Symbol	Meaning
	Manufacturer
	CE-mark and Identification number of the Notified Body. The product meets the essential requirements of Medical Device Directive 93/42/EEC.
	Caution
	Batch code
	Catalogue number
	Do not reuse
	Do not re-sterilize
	Do not use if package is damaged
	Authorized Representative in the European Community
	Quantity
	Use-by date
	Sterilized using ethylene oxide
	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed healthcare practitioner

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INTENDED USE

This mesh is intended for the repair of abdominal wall hernias and abdominal wall deficiencies that require the addition of a reinforcing material to obtain the desired surgical result.

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60 months

MD DISPOSAL RULES AND CONDITIONS

Products coming into contact with blood and/or biological fluids are class B medical devices under medical wastes classification and may be recycled and / or disposed in accordance with applicable national regulations.

Unused products (not having contact with blood and/or biological fluids), including those with expired shelf life, belong to class A medical devices under medical waste classification and may be recycled and/or disposed in accordance with applicable national regulations

MRI STATEMENT

The materials are non-metallic and does not contain ferromagnetic or conductive materials that would interact with any magnetic energy/field

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created by MRI and thus the device can be considered MR safe.

INFORMATION ABOUT IMPLANT REMOVAL AND EXTRACTION

The type of contact is permanent, so removal and extraction is not foreseen.

WARRANTY

Manufacturer guarantees to the end user absence of defects and conformity to product specifications for the entire shelf life of the product under the conditions of adherence to storage, transportation and package integrity requirements. There are no warranties that extend beyond the description on the face hereof. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use

INFORMATION ABOUT MANUFACTURER

Legal Manufacturer:

Ethicon LLC, Highway 183 Km. 8.3, San Lorenzo, Puerto Rico 00754, USA

Physical Manufacturer:

Ethicon Inc., 655 Ethicon Circle, Cornelia, Georgia, 30531, USA
Ethicon LLC, Highway 183 Km. 8.3, San Lorenzo, Puerto Rico 00754, USA

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER

For additional information related to the use of the medical product, please contact an authorized representative of the manufacturer, "Johnson & Johnson" LLC, 121614, Russia, Moscow, st. Krylatskaya, building 17/2,
phone: 8 (495) 580 77 77,
fax: 8 (495) 580 78 78.

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