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


Note: The reusable ARTH 263300 42 tubing set
is available as an optional extra
(see 'Recommended accessories').

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Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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Device description

The ARTHROPUMP® is a microprocessor-controlled reel pump for controlled extension of the joint in arthroscopic examinations and operations. The speed-controlled reel pump permits a pre-selectable filling rate of up to 1000 ml/min. An electronic pressure sensor ensures constant intraarticular pressure and additional safety is guaranteed by an integrated overpressure valve. Operation and monitoring of the unit have been simplified as far as possible. Due to a new design of pump head insertion of the hoses is considerably simplified; correct positioning of the hose system is also monitored electronically. Parallel displays of the rated value and actual value permit optimal monitoring of the current operating status. An electronic safety shutdown system stops flow if there are continuing departures from set value and warns the user via audible signals.

An electronic autocheck system which checks the various system components when the unit is put into operation and reports any error statuses contributes to the high safety standard of the ARTHROPUMP®.

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User qualification

The ARTHROPUMP® may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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**Indications for Use/
Contraindications for Use****Indications for Use**

The ARTHROPUMP® is used for the infusion of sterile irrigation solution into the knee, shoulder ankle, elbow, wrist or temporal mandibular joint (TMJ). The delivery of pressurized irrigation solution through the arthroscope aids in the visualization and distention of the joint, rinsing or removal of deposits, blood clots or excised tissue from the surgical field for accurate diagnosis and therapeutic treatment during arthroscopic procedures.

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
Contraindications

The ARTHROPUMP® should not be used when arthroscopy is contraindicated.

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Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations. It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants. The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using it (see the section "Testing the unit for proper operation").

During treatment using the ARTHROPUMP® the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

Any treatment may only be performed if there is visual observation of the action of the unit.

Increase in the backward-rotation of the reel pump (Arthropump as of serial number 688)

If the actual pressure exceeds the set pressure during operation, the reel pump rotates backward 20 turns instead of 3/4 turns.

Warning: The pressure measuring chamber/tubing set may be contaminated by the increased backward-rotation. For this reason, the tubing set must be replaced before every application.

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Installation and operating instructions

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Unpacking the equipment

Carefully unpack the ARTHROPUMP® and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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Basic equipment

- 1 ARTHROPUMP® 283300 20-1
- 3 mtp Disposable tubing kit
- 1  connecting cord 200901 70
- 1 Power cord 400 A
- 1 Instruction manual
- 1 Wallet for instruction manual
- 1 Operator's quick reference guide (self-adhesive)
- 6 Label strips

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
Installation

Note: The ARTHROPUMP® including the accessories connected may be used only in medical with electrical installations conforming to the applicable national, state, and local electrical regulations.

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Warning: The unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Fitting a label strip to the front panel

Undo the left-hand front side cover with a screwdriver.

Note: The appropriate controls/indicators on the unit are designated on the label strip.

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
Remove any unmarked label strip. Select a label strip ② (currently 6 languages available) and slide it between the plastic and the metal front.

Fold the strip along the broken line to an angle of 90° and approx. 2 mm further on (toward the end of the strip) fold it back to an angle of 90° so that the side cover can be refitted without any obstructions.

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
The unit is equipped with a connector ⑤ for attaching a ground line. The unit's ground line should be installed by a qualified electrician.

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Cleaning and care of the unit

Warning: Always pull out power plug before cleaning!

Caution: Avoid allowing fluids to enter the unit.

The exterior surfaces of the unit can be disinfected with a disposable cloth moistened with disinfectant (preferably alcohol-free).

Note: For cleaning, disinfecting and sterilizing the reusable tubing set ARTH 26 3300 42, observe the instruction manual delivered with the tubing set.

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Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year. For detailed instructions please refer to the version of the Service Manual currently valid, catalog number SV 3377 (English version).

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important Information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Standard compliance (for 283300 20-1)

According to IEC 60601-1, UL 2601.1,
CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type BF 

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the appendix (p. 29-41).

Directive compliance (for 283300 20-1)

According to Medical Device Directive (MDD):

This medical product belongs to Class II b

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

Note: The code number after the CE mark indicates the responsible notified body.

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

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We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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