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Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG., hereby certify that the attached documents that keep the number of pages indicated below are the operational documentation. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: ZQ
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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
Highlighted Changes: Revised December 2017.

NAME

Alinity i Syphilis TP Controls (also referred to as Syphilis Ctrls)

INTENDED USE

The Alinity i Syphilis TP Controls are for the estimation of test precision and the detection of systematic analytical deviations of the Alinity i analyzer when used for the qualitative detection of antibody to *Treponema pallidum* (TP) in human serum and plasma.

For additional information, refer to the Alinity i Syphilis TP reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

CONTROL - prepared with recalcified human plasma. **CONTROL +** prepared in recalcified human plasma (inactivated); reactive for anti-TP. Preservatives: sodium azide and other antimicrobial agents.

The controls are at the following ranges:

Control	Quantity	Color	RANGE S/CO
CONTROL -	1 x 8.0 mL	Natural	≤ 0.40
CONTROL +	1 x 8.0 mL	Blue ^a	1.25 - 3.75

^a Dye: Acid Blue No. 9

NOTE: The insert ranges for the controls are not lot specific and represent the total range of values which may be generated throughout the life of the product. It is recommended that each laboratory establish its own means and acceptable ranges which should fall within the package insert ranges. Sources of variation that can be expected include:

- Calibration
- Control lot
- Reagent lot
- Calibrator lot
- Instrument

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

Safety Precautions

- **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that this product and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human-sourced material used in the Negative Control is nonreactive for anti-TP, anti-HCV, HBsAg, HIV-1 RNA or HIV-1 Ag, and anti-HIV-1/HIV-2.
- The human-sourced material used in the Positive Control is reactive for anti-TP and nonreactive for anti-HCV, HBsAg, HIV-1 RNA or HIV-1 Ag, and anti-HIV-1/HIV-2.

The following warnings and precautions apply to: **CONTROL -** and

CONTROL +

Contains sodium azide.

EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- This product may be used immediately after removal from 2 to 8°C storage.
- Prior to each use, mix by gentle inversion (5 to 10 times).

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	Until expiration date	Store tightly capped. Store in an upright position. Return to refrigerated storage after use.

INSTRUMENT PROCEDURE

- To obtain the recommended volume requirements for the controls, hold the bottle vertically, and dispense 5 drops of the negative control and 5 drops of the positive control into each sample cup in the assigned position.
- For instructions on ordering and loading controls on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, or if controls do not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

ISO 15223 Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONTROL -	Negative Control
CONTROL +	Positive Control
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
REF	List Number

Other Symbols

CN	Control Number
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information needed for United States of America only
PRODUCT OF GERMANY	Product of Germany
RANGE	Range
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).

Alinity is a trademark of Abbott Laboratories in various jurisdictions.
All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com

Revised December 2017.

©2016, 2017 Abbott Laboratories

ADUM

Instruction for Use Alinity i Syphilis TP Controls (for Russia)

PNIG

Recommended shipping conditions for device is 2-8 in single-use thermo containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in autorefrigerators using thermal indicators. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 12 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

90/79/EC, 1272/2008/EC (CLP) as amended, 1907/2006/EC (REACH) as amended, Regulation (EC) 1069/2009, Regulation (EC) 1432/2011, CLASS 2013 (as amended), GHS (as amended), WHSR Part 7.1 (as amended), EN ISO 9001, ISO 13485, EN ISO 14971, EN 13641, EN 13612, EN ISO 18223-1, EN 980, EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, EN ISO 23940, EN ISO 20776-1, EN 13675, EN ISO 17511

POTENTIAL USER

The Alinity i Syphilis TP Controls is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals and laboratories, whose personnel have been trained in the use of the Alinity i analyser.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SD6. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRA-INDICATIONS FOR USE

Not applicable.

STORAGE AND USE CONDITION

Storage of Medical device at a temperature of 2 °C to 8 °C should be carried out in cold rooms or in refrigerators providing a regulated temperature regime with a daily temperature record.

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- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

DISTRIBUTED IN THE USA BY
Abbott Laboratories
Abbott Park, IL 60064 USA

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