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Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Michael Fritz, Manager Regulatory/Quality Europe, Middle East, Africa (EMEA) and Pakistan, International Regulatory Affairs & Division Labeling, Abbott GmbH & Co. KG, hereby certify that attached documents are true and accurate copies of the manufacturer's instructions of usage with addendums and also their translation into Russian language in accordance with number of pages indicated below.

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#	Title	Number of pages
1.	ARCHITECT Methotrexate Reagent Kit (package insert) including ADDENDUM to the package insert of ARCHITECT Methotrexate Reagent Kit	9
2.	ARCHITECT Methotrexate Calibrators (package insert) including ADDENDUM to the package insert of ARCHITECT Methotrexate Calibrators	3
3.	ARCHITECT Methotrexate Controls (package insert) including ADDENDUM to the package insert of ARCHITECT Methotrexate Controls	3

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Michael Fritz,
Manager Regulatory/Quality Europe, Middle East,
Africa (EMEA) and Pakistan, International Regulatory
Affairs & Division Labeling

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden-Delkenheim
Germany

2016-02-11

Abbott

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden

RCHITECT**SYSTEM****en**

Methotrexate

REF 2P49**B2P490****611-031 9/14**

Methotrexate

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Key to Symbols

Consult instructions for use.

Manufacturer

Store at 2-8°C

Sufficient for

Use by/Expiration Date

ASSAY SPECIFIC DILUENT

Assay Specific Diluent

CONJUGATE

Conjugate

CONTAINS: CYANATE

Contains sodium thiocyanate.
Contact with acids liberates
very toxic gas.

CONTROL NO.

Control Number

GTIN

Global Trade Item Number

INFORMATION FOR USA ONLY

Information needed for
United States of America only
In Vitro Diagnostic Medical
Device

IVD

Lot Number

LOT**MICROPARTICLES**

Microparticles

MULTI-ASSAY MANUAL DILUENT

Multi-Assay Manual Diluent

PRE-TRIGGER SOLUTION

Pre-Trigger Solution

PRODUCED FOR ABBOTT BY

Produced for Abbott by

PRODUCT OF USA

Product of USA

REACTION VESSELS

Reaction Vessels

REAGENT LOT

Reagent Lot

REF

List Number

REPLACEMENT CAPS

Replacement Caps

Rx ONLY

For use by or on the order of
a physician only (applicable to
USA classification only)

SAMPLE CUPS

Sample Cups

SEPTUM

Septum

SN

Serial Number

TRIGGER SOLUTION

Trigger Solution

WARNING: SENSITIZER

Warning: May cause an
allergic reaction.

WASH BUFFER

Wash Buffer

**Abbott**

NAME

ARCHITECT Methotrexate

INTENDED USE

The ARCHITECT Methotrexate assay is a chemiluminescent microparticle immunoassay (CMIA) for the quantitative determination of methotrexate in human serum and plasma on the ARCHITECT i System with STAT protocol capability. The measurements obtained are used in monitoring levels of methotrexate to ensure appropriate therapy.

SUMMARY AND EXPLANATION OF TEST

Methotrexate (MTX, amethopterin) is a folate antimetabolite. It is an analog of aminopterin, which is also derived from folic acid. The molecular structure of methotrexate differs from folic acid in that it has a hydroxyl group in place of the 4-amino group on the pteridine ring and there is no methyl group at the N position.

Methotrexate is an antineoplastic drug used solely or in combination with other antineoplastic drugs for the treatment of leukemia and other diseases.^{1,2} Relatively low doses of methotrexate have been used in the treatment of nonmalignant diseases such as severe psoriasis, asthma, rheumatoid arthritis, sarcoidosis, and transplantation therapy.³⁻⁸ Intermediate to high doses of methotrexate with Leucovorin rescue have been used with favorable results in the treatment of osteogenic sarcoma, leukemia, non-Hodgkin's lymphoma, lung and breast cancer.^{9,13}

Methotrexate levels are monitored to avoid excessive toxic effects of the drug and to determine when to intervene with counter-acting 'rescue' therapy. Adverse reactions include myelosuppression, stomatitis, nausea, vomiting, convulsions, liver and renal abnormalities, anemia, leukopenia, thrombocytopenia, osteoporosis, neurotoxicity, and leukoencephalopathy.² There are no current indications for monitoring low dose methotrexate therapies.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT Methotrexate assay is a one-step immunoassay for the quantitative determination of methotrexate in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

Sample, anti-methotrexate coated paramagnetic microparticles, assay specific diluent, and methotrexate acridinium-labeled conjugate are combined to create a reaction mixture. The anti-methotrexate coated microparticles bind to the methotrexate present in the sample and to the acridinium-labeled conjugate. After washing, Pre-Trigger and Trigger Solutions are added to the reaction mixture. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is an inverse relationship between the amount of methotrexate in the sample and the RLUs detected by the ARCHITECT i System optics.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Reagent Kit, 100 Tests

ARCHITECT Methotrexate Reagent Kit (2P49)

- **MICROPARTICLES** 1 Bottle (6.6 mL) Anti-methotrexate (mouse monoclonal) coated microparticles in MES buffer with non-protein stabilizers. Minimum concentration: 0.073% solids. Preservative: ProClin 300.
- **CONJUGATE** 1 Bottle (15.0 mL) Acridinium-labeled methotrexate conjugate in citrate buffer with non-protein stabilizers. Minimum concentration: 22 ng/mL. Preservative: ProClin 300.
- **ASSAY SPECIFIC DILUENT** 1 Bottle (5.9 mL) Assay Specific Diluent containing MES buffered saline diluent. Preservative: ProClin 300.

Assay Diluent

ARCHITECT Multi-Assay Manual Diluent (7D82-50)

- **MULTI-ASSAY MANUAL DILUENT** 1 Bottle (100 mL) ARCHITECT Multi-Assay Manual Diluent containing phosphate buffered saline solution. Preservative: antimicrobial agent.

Other Reagents

ARCHITECT Pre-Trigger Solution

- **PRE-TRIGGER SOLUTION** Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

ARCHITECT Trigger Solution

- **TRIGGER SOLUTION** Trigger Solution containing 0.35N sodium hydroxide.

ARCHITECT Wash Buffer

- **WASH BUFFER** Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

WARNINGS AND PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

Safety Precautions

- **CAUTION:** This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹⁴ Biosafety Level 2¹⁵ or other appropriate biosafety practices^{16,17} should be used for materials that contain or are suspected of containing infectious agents.

- The following warnings and precautions apply to these components:

- Microparticles
- Conjugate

WARNING:	Contains methylisothiazolones
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

- The following warnings and precautions apply to the Assay Specific Diluent:

WARNING:	Contains methylisothiazolones and sodium thiocyanate.
H317	May cause an allergic skin reaction.
EU032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.
- For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagent kits beyond the expiration date.
- **Do not pool reagents within a kit or between reagent kits.**
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**

To avoid contamination, wear gloves and use a septum on an uncapped reagent bottle. Once a septum has been used, do not invert the bottle as this will compromise assay results. Over time, residual liquid in the bottle may be typically dried salts.

For a detailed discussion of safety precautions, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

2-8°C

The ARCHITECT i System reagent kit must be stored at 2-8°C. When stored and used within the expiration date.

The ARCHITECT i System reagent kit must be stored at 2-8°C. When stored and used within the expiration date.

Reagents may be removed from the system, it and boxes does not storage

For info Oper

Indicate

When a

deteric

result

be n

Syst

IN

To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.

- Once a septum has been placed on an open reagent bottle, do not invert the bottle as this will result in reagent leakage and may compromise assay results.
- Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

-  The ARCHITECT Methotrexate Reagent Kit must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.
- When stored and handled as directed, the reagents are stable until the expiration date.
- The ARCHITECT Methotrexate Reagent Kit may be stored on board the ARCHITECT i System for a maximum of 30 days. After 30 days, the reagent kit must be discarded. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.
- Reagents may be stored on or off the ARCHITECT i System. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.**
- For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

- The ARCHITECT Methotrexate assay file must be installed on the ARCHITECT i System with STAT protocol capability from an ARCHITECT i System Assay CD-ROM prior to performing the assay.
- When the Methotrexate (MTX) assay file is installed, calculated assays MTX20ug, MTX20uM, MTX400ug, MTX400uM, MTX8000ug, and MTX8000uM are automatically installed in order to perform calculations on diluted samples. These calculated assays are used in place of entering manual dilution factors.
- **Three of the calculated assays must be uninstalled prior to running the Methotrexate assay.** Each laboratory should determine which calculated assays to uninstall based on their individual unit of measure requirements.
 - IF laboratory results are reported in $\mu\text{mol/L}$, THEN uninstall calculated assays 3030 (MTX20ug), 3031 (MTX400ug), and 3032 (MTX8000ug) per ARCHITECT System Operations Manual, Section 2, as these will not be used.
 - IF laboratory results are reported in $\mu\text{g/mL}$, THEN configure the instrument to report the alternate units of $\mu\text{g/mL}$ according to the ARCHITECT System Operations Manual, specifically the "Change the result units setting" of Section 2, AND uninstall calculated assays 3027 (MTX20uM), 3028 (MTX400uM), and 3029 (MTX8000uM) per ARCHITECT System Operations Manual, Section 2, as these will not be used.
- For detailed information on assay file installation and deletion, and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.
- For specific details on ordering a calculated assay, refer to the ARCHITECT System Operations Manual, Section 5.
- For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.
- For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

Note: Do not use collection tubes with separator gel.

- Human serum or plasma samples tested for methotrexate should be protected from light.
- The specimen collection tubes listed below were verified for use with the ARCHITECT Methotrexate assay. Other specimen collection tubes have not been tested with this assay.
 - Human serum collected in tubes without separator gel (red top tubes only)
 - Human plasma collected in tubes without separator gel
 - Dipotassium (K_2) EDTA
 - Sodium Heparin
 - Lithium Heparin
- Performance has not been established for the use of cadaveric specimens or body fluids other than human serum or plasma.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Human serum or plasma samples tested for methotrexate should be protected from light.
- Do not use specimens with the following conditions:
 - heat-inactivated
 - pooled
 - grossly hemolyzed
 - obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use disposable pipettes or pipette tips.
 - Replace disposable pipettes or pipette tips between samples.
 - Replace pipette tips between serial dilution steps for samples and controls.
- Replace gloves if contamination is suspected.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for serum and plasma collection tubes. Gravity separation is not sufficient for specimen preparation.
- Mix thawed specimens thoroughly by low-speed vortexing or by inverting 10 times. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous. If samples are not mixed thoroughly, inconsistent results may be obtained.
- To ensure consistency in results, specimens must be transferred to a centrifuge tube and centrifuged at a Relative Centrifugal Force (RCF) $\geq 10,000 \times g$ for 10 minutes before testing if
 - they contain fibrin, red blood cells, or other particulate matter,
 - they require repeat testing, or
 - they were frozen and thawed.

Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

- Human serum or plasma samples tested for methotrexate should be protected from light.
- Specimens may be stored on or off the clot or red blood cells
 - for up to 1 day at room temperature
 - for up to 2 days at 2-8°C
- If testing will be delayed more than 2 days, remove serum or plasma from the clot or red blood cells and store at -10°C or colder.
- Avoid more than 3 freeze/thaw cycles.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- It is recommended that specimens be removed from the clot or red blood cells.
- Ship at 2-8°C (wet ice) or frozen (dry ice). Protect from light.
- Do not exceed the storage time limitations listed above.

PROCEDURE

Materials Provided

- 2P49 ARCHITECT Methotrexate Reagent Kit

Materials Required but not Provided

- ARCHITECT i System
- ARCHITECT Methotrexate Assay File, may be obtained from:
 - ARCHITECT i System e-Assay CD-ROM found on www.abbottdiagnostics.com
 - ARCHITECT i System Assay CD-ROM
- 2P49-01 ARCHITECT Methotrexate Calibrators (A-F)
- 2P49-10 ARCHITECT Methotrexate Controls (L, M, H, and X) or other commercially available controls
- 2P49-15 ARCHITECT Methotrexate Extended Range Controls (Y and Z)
- 7D82-50 ARCHITECT **MULTI-ASSAY MANUAL DILUENT**
- ARCHITECT **PRE-TRIGGER SOLUTION**
- ARCHITECT **TRIGGER SOLUTION**
- ARCHITECT **WASH BUFFER**
- ARCHITECT **REACTION VESSELS**
- ARCHITECT **SAMPLE CUPS**
- ARCHITECT **SEPTUM**
- ARCHITECT **REPLACEMENT CAPS**
- Microcentrifuge tubes or equivalent
- Pipettes or pipette tips (optional) to deliver the specified volumes.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
 - **Invert the microparticle bottle 30 times.**
- Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
- **If the microparticles do not resuspend, DO NOT USE. Contact your Abbott representative.**
- Once the microparticles have been resuspended, place a septum on the bottle. For instructions about placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT i System with STAT protocol capability.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles.
- Order calibration, if necessary.
 - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
 - For information on ordering patient specimens and controls, and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test. Maximum number of replicates sampled from the same sample cup: 10.
 - Priority:
 - Sample volume for first test: 60 µL
 - Sample volume for each additional methotrexate test from the same sample cup: 10 µL
 - ≤ 3 hours on-board:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional methotrexate test from the same sample cup: 10 µL
 - > 3 hours on-board: Replace with a fresh sample (patient specimens, controls, and calibrators).
 - If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT Methotrexate Calibrators and Controls.
 - Mix calibrators and controls by gentle inversion before use.
 - Hold bottles **vertically**, and dispense recommended volumes into each respective sample cup.
 - Recommended volumes:
 - for each Calibrator (A-F): 5 drops
 - for each Control (Low, Medium, and High): 5 drops

- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Specimen Dilution Procedures

- Specimens with a methotrexate concentration greater than 1.500 µmol/L (0.682 µg/mL) will be flagged as "> 1.500 µmol/L" ("> 0.682 µg/mL") and must be diluted with the Manual Dilution Procedure.

Note: Specimens with a very high Methotrexate concentration may be flagged with the following Error Code: "1232/ Result cannot be calculated, final RLU read is outside the specification of the highest calibrator." Dilute the patient sample using the Manual Dilution Procedure below and rerun using the appropriate calculated assays.

- **Caution:** In order to prevent cross contamination, **do not remove dropper tips** from the ARCHITECT Methotrexate control bottles. Use caution when handling control bottles and caps.
- Replace pipette tips between samples, controls, and serial dilution steps.

Manual Dilution Procedure

- The suggested dilution for the ARCHITECT Methotrexate assay is 1:20. Up to 3 serial dilutions are made to the sample (1:20, 1:400, and 1:8000).
- **If a sample requires dilution, run one replicate of the Medium control and the diluted X, Y, and Z controls along with the diluted samples.**
- All samples in bold font below will be tested. Some tubes prepared are for intermediate steps in the serial dilution procedure and will not be tested.

Note: For all manual dilutions:

Diluent = ARCHITECT Multi-Assay Manual Diluent
DF = Dilution Factor

Serial Dilution of Specimen > 1.500 µmol/L (> 0.682 µg/mL)

1. Label 3 microcentrifuge tubes to be used for the dilution of the specimen:
 - Tube #1 DF=20
 - Tube #2 DF=400
 - Tube #3 DF=8000
2. Add 950 µL Diluent to each of the 3 tubes.
3. Add 50 µL of the patient specimen to Tube #1.
 - a. Vortex the sample for 10-30 seconds to mix.
 - b. Tube #1 DF=20 is ready for testing.
4. Remove 50 µL from Tube #1 (Step 3b) and add it to Tube #2.
 - a. Vortex the sample for 10-30 seconds to mix.
 - b. Tube #2 DF=400 is ready for testing.
5. Remove 50 µL from Tube #2 (Step 4b) and add it to Tube #3.
 - a. Vortex the sample for 10-30 seconds to mix.
 - b. Tube #3 DF=8000 is ready for testing.

Serial Dilution of Controls X, Y, and Z

6. Label 6 microcentrifuge tubes to be used for the dilution of the controls:
 - Tube X DF=20
 - Tube Y1
 - Tube Y2 DF=400
 - Tube Z1
 - Tube Z2
 - Tube Z3 DF=8000

Note: Only the tubes labeled with "DF" will be tested.

7. Add 950 µL Diluent to each of the 6 tubes.
8. Dispense 5 drops of Control X into a separate tube.
 - a. Remove 50 µL of Control X from this tube and add it to Tube X DF=20.
 - b. Vortex the sample for 10-30 seconds to mix.
 - c. Tube X DF=20 is ready for testing.
9. Dispense 5 drops of Control Y into a separate tube.
 - a. Remove 50 µL of Control Y from this tube and add it to Tube Y1.
 - b. Vortex the sample for 10-30 seconds to mix.
 - c. Tube Y1 will not be tested.
10. Remove 50 µL from Tube Y1 (Step 9c) and add it to Tube Y2 DF=400.
 - a. Vortex the sample for 10-30 seconds to mix.
 - b. Tube Y2 DF=400 is ready for testing.

Dispense 5 drops of Control Z into a separate tube.

- Remove 50 µL of Control Z from this tube and add it to Tube Z1.
- Vortex the sample for 10-30 seconds to mix.
- Tube Z1 will not be tested.

- Remove 50 µL from Tube Z1 (Step 11c) and add it to Tube Z2.

- Vortex the sample for 10-30 seconds to mix.
- Tube Z2 will not be tested.

- Remove 50 µL from Tube Z2 (Step 12b) and add it to Tube Z3 DF=8000.

- Vortex the sample for 10-30 seconds to mix.
- Tube Z3 DF=8000 is ready for testing.

- Test the Medium Control, Control X (1:20 dilution), Control Y (1:400 dilution), and Control Z (1:8000 dilution). Refer to the appropriate table below.

- Test diluted patient samples using the calculated assays files per the appropriate table below.

- FOR laboratory results reported in µmol/L:

To Order	Order Assay File
• Medium Control (to evaluate run validity)	MTX
• Control X (1:20 dilution)	MTX20uM*
• Patient Sample (1:20 dilution)	
• Control Y (1:400 dilution)	MTX400uM*
• Patient Sample (1:400 dilution)	
• Control Z (1:8000 dilution)	MTX8000uM*
• Patient Sample (1:8000 dilution)	

*Note: Do not enter the manual dilution factor in the order screen. The assay file calculates the result based on the appropriate dilution factor.

- FOR laboratory results reported in µg/mL:

To Order	Order Assay File
• Medium Control (to evaluate run validity)	MTX
• Control X (1:20 dilution)	MTX20ug*
• Patient Sample (1:20 dilution)	
• Control Y (1:400 dilution)	MTX400ug*
• Patient Sample (1:400 dilution)	
• Control Z (1:8000 dilution)	MTX8000ug*
• Patient Sample (1:8000 dilution)	

*Note: Do not enter the manual dilution factor in the order screen. The assay file calculates the result based on the appropriate dilution factor.

- IF results are obtained for 2 dilutions of the same patient sample, THEN report the result with the lowest dilution factor per the appropriate table below.

- FOR laboratory results reported in µmol/L:

If results are obtained from	Then report result from
MTX20uM and MTX400uM	MTX20uM
MTX400uM and MTX8000uM	MTX400uM

- FOR laboratory results reported in µg/mL:

If results are obtained from	Then report result from
MTX20ug and MTX400ug	MTX20ug
MTX400ug and MTX8000ug	MTX400ug

Note: If your laboratory uses the "Result Details Report", then the calculated result is provided next to the carrier/position (C / P) and is circled in the Result Details Report example below. Do NOT report the result next to the Constituent Assay MTX.

C / P	Result	Units	Range	Flags	Code
A765/3	27.200	µmol/L			

Module	Constituent assay	Result	Flags
1	MTX	0.068 µmol/L	

Calibration

- Test Calibrators A-F in duplicate. The calibrators should be priority loaded.
- A single sample of each ARCHITECT Methotrexate Low, Medium, and High Control must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Calibration Range: 0.000 to 1.500 µmol/L (0.000 to 0.682 µg/mL)
- Once an ARCHITECT Methotrexate calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

QUALITY CONTROL PROCEDURES

The recommended control requirement for the ARCHITECT Methotrexate assay is that a single sample of each Low, Medium, and High Control be tested once every 24 hours each day of use. It is recommended that Controls X, Y, and Z be diluted and tested with diluted samples.

For detailed dilution procedure instructions and quality control, refer to the **PROCEDURE, Specimen Dilution Procedures** section of this package insert.

If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures. Additional controls may be tested in conformance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

Each laboratory should establish control ranges to monitor the acceptable performance of the assay. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be indicated.

If a diluted control is out of its specified range, the associated specimen dilution test results are invalid. Specimens and controls must be prepared again. Refer to the **PROCEDURE, Specimen Dilution Procedures** section of this package insert.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

The ARCHITECT Methotrexate assay belongs to method group 6.

RESULTS

Calculation

The ARCHITECT Methotrexate assay uses a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration curve.

Alternate Result Units

- The default unit for the ARCHITECT Methotrexate assay is µmol/L.
- When the alternate result unit, µg/mL, is selected, the following conversion factor is used by the system:
- Conversion Formula:
(Concentration in µmol/L) x (0.45444) = µg/mL

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in µmol/L which meets the limits of acceptable performance for both imprecision and bias for an undiluted sample. For the studies described in this package insert, the range is 0.040 µmol/L (Limit of Quantitation - LoQ) to 1.500 µmol/L (0.018 to 0.682 µg/mL).

LIMITATIONS OF THE PROCEDURE

- Do not use collection tubes with separator gel.
- Samples to be tested for methotrexate should be protected from light.
- Results should be used in conjunction with other data; e.g., symptoms, results of other tests, and clinical impressions.
- If the methotrexate results are inconsistent with clinical evidence, additional testing is recommended.

- Specimens from patients who have received glucarpidase (carboxypeptidase G₂) as a high-dose methotrexate rescue therapy should not be tested with the ARCHITECT Methotrexate assay for at least 48 hours following the last dose of glucarpidase.¹⁸ These specimens have increased serum levels of 4-[[2,4-diamino-6-(pteridinyl)methyl]-methylamino]-benzoic acid (DAMPA) that result from metabolism of methotrexate by glucarpidase.¹⁹ DAMPA crossreacts with the methotrexate antibody used in this assay. Oncologists on the clinical team should notify the laboratory when glucarpidase is administered.²⁰
- Aminopterin is derived from folic acid and crossreacts with the methotrexate antibody used in this assay.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA).²¹ Such specimens may show anomalous values when tested with assay kits such as ARCHITECT Methotrexate that employ mouse monoclonal antibodies. Additional information may be required for diagnosis.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays.²² Patients routinely exposed to animals or to animal serum products can be prone to this interference, and anomalous results may be observed. Additional information may be required for diagnosis.

EXPECTED VALUES

No precise relationship between methotrexate serum levels and antineoplastic efficacy has been established, although levels below approximately 0.02 $\mu\text{mol/L}$ were seen as necessary for resumption of DNA synthesis.²³ The correlation between serum methotrexate drug concentration and duration of tumor cell exposure has been shown to predict methotrexate toxicity.

Following high-dose methotrexate treatment with Leucovorin rescue, patients with a 24-hour serum methotrexate concentration $>10 \mu\text{mol/L}$, a 48-hour level $>1.0 \mu\text{mol/L}$, and a 72-hour level $>0.1 \mu\text{mol/L}$ are at an increased risk of toxicity.²⁴

Toxicity is typically present in the form of myelosuppression, stomatitis, nausea, vomiting, convulsions, and liver and renal abnormalities.² Anemia, leukopenia, thrombocytopenia, osteoporosis, skin and mucosal involvement with a fatal outcome have also been reported. Neurotoxicity and leukoencephalopathy are also reported as toxic effects involving methotrexate.²

SPECIFIC PERFORMANCE CHARACTERISTICS

Data in the SPECIFIC PERFORMANCE CHARACTERISTICS section were generated using the ARCHITECT i 2000_{SR} System.

Assay results obtained in individual laboratories may vary from data presented.

Precision

The ARCHITECT Methotrexate assay is designed to have within-laboratory (total) imprecision of $\leq 7.5\%$ CV for samples with methotrexate concentrations between 0.040 $\mu\text{mol/L}$ and $\leq 12.500 \mu\text{mol/L}$ and $\leq 10\%$ CV for samples with methotrexate concentrations $> 12.500 \mu\text{mol/L}$.

A study was performed based on guidance from the National Committee for Clinical Laboratory Standards (NCCLS) document EP5-A2.²⁵ The ARCHITECT Methotrexate Controls, Extended Range Controls, and 5 levels of human serum panels (Panels 1, 3, and 5-7) and 2 levels of human plasma panels (Panels 2 and 4) were assayed. All samples were tested in replicates of 2 at 2 separate times per day for 20 days. Additional replicates for Panels 6 and 7 were performed to assess precision at the measuring interval boundaries. The data are summarized in the following table.

Sample	N	Mean ($\mu\text{mol/L}$)	Within-Laboratory (Total)			
			SD	%CV	SD	%CV
Low Control	80	0.072	0.0023	3.2	0.0029	4.0
Medium Control	80	0.436	0.0134	3.1	0.0148	3.4
High Control	80	0.949	0.0401	4.2	0.0402	4.2
Control X	80	10.011	0.4060	4.1	0.4642	4.6
Control Y	80	55.645	1.9365	3.5	2.4197	4.4
Control Z	80	554.900	20.0998	3.6	39.0768	7.0
Panel 1	80	0.065	0.0026	4.1	0.0029	4.6
Panel 2	80	0.101	0.0044	4.4	0.0053	5.3
Panel 3	80	1.103	0.0475	4.3	0.0547	5.0
Panel 4	80	7.966	0.2683	3.4	0.3084	3.9
Panel 5	80	1679.100	60.8687	3.6	83.4626	5.0
Panel 6	160	0.047	0.0020	4.4	0.0025	5.4
Panel 7	160	1.347	0.0554	4.1	0.0590	4.4

Recovery

The ARCHITECT Methotrexate assay is designed to have a recovery of $100 \pm 10\%$ for samples with methotrexate concentrations from 0.040 $\mu\text{mol/L}$ to 1.500 $\mu\text{mol/L}$.

A study was performed with 5 serum and 5 plasma samples that spiked with known amounts of methotrexate to create samples across the measuring interval. The concentration of methotrexate was determined using the ARCHITECT Methotrexate assay and the resulting percent recovery was calculated. The mean percent recovery was 97% for serum samples and 98% for plasma samples.

Linearity

The ARCHITECT Methotrexate assay is designed to have a deviation from linearity within $\pm 10\%$ for samples with methotrexate concentrations from 0.040 $\mu\text{mol/L}$ to 1.500 $\mu\text{mol/L}$.

A study was performed based on guidance from the NCCLS document EP6-A.²⁶ Serum and plasma dilution sets with methotrexate concentrations ranging from 0.010 $\mu\text{mol/L}$ to 2.001 $\mu\text{mol/L}$ were prepared gravimetrically and evaluated. The observed deviation from linearity was within $\pm 10\%$ for samples from 0.040 $\mu\text{mol/L}$ to 1.500 $\mu\text{mol/L}$.

An additional study was performed using serum and plasma samples with methotrexate concentrations ranging from 1.2 $\mu\text{mol/L}$ to 2500 $\mu\text{mol/L}$. The samples were prepared gravimetrically and diluted using serial 1:20 dilutions into the measuring interval with ARCHITECT Multi-Assay Manual Diluent. Regression of assayed methotrexate concentrations was linear throughout the range.

Sensitivity

Limit of Quantitation (LoQ)

The ARCHITECT Methotrexate assay is designed to have a Limit of Quantitation (LoQ) of $\leq 0.040 \mu\text{mol/L}$.

Based on guidance from the Clinical and Laboratory Standards Institute (CLSI) document EP17-A2,²⁷ a study was performed with 4 zero-level samples and 8 low-level methotrexate samples (targeted to approximately 0.010, 0.020, 0.030, 0.040, 0.050, 0.080, 0.100, and 0.200 $\mu\text{mol/L}$). These samples were tested in 6 separate runs over 4 days using 2 reagent lots and 2 instruments. The LoQ is defined as the lowest amount of analyte in a sample that can be quantitatively determined with a Percent Total Allowable Error $< 25\%$. The LoQ for the ARCHITECT Methotrexate assay was 0.020 $\mu\text{mol/L}$.

Limit of Blank and Limit of Detection

In the same study, the Limit of Blank (LoB) and Limit of Detection (LoD) were determined. The LoB was 0.005 $\mu\text{mol/L}$ and the LoD was 0.009 $\mu\text{mol/L}$.

Interference

Potentially Interfering Drugs

A study was performed based on guidance from the CLSI document EP7-A2.²⁸ Potentially interfering drugs were evaluated to determine whether methotrexate concentrations were affected when using the ARCHITECT Methotrexate assay. The drugs listed below were spiked into samples with methotrexate concentrations at the medical decision points of approximately 0.050 $\mu\text{mol/L}$ and 1.000 $\mu\text{mol/L}$. The samples were assayed, and the methotrexate concentrations of the spiked samples were compared to reference samples. The ARCHITECT Methotrexate Assay did not crossreact with the major metabolite 7-Hydroxymethotrexate.

The drugs were tested at 1000 $\mu\text{mol/L}$ except where noted. The data are summarized in the following table.

The compounds Aminopterin and DAMPA were also tested in serum at 1000 $\mu\text{mol/L}$ in the absence of methotrexate and the cross-reactivity was 61% and 46%, respectively.

have a reactivity
interfering from
that were
across
determined
percent
for serum
from
from
ument
ations
ically
% for
with
'L
'D
'I

Potentially Interfering Substances	% Cross-Reactivity ^a	
	Methotrexate Concentration	
	0.050 µmol/L	1.000 µmol/L
Amphotericin	0	0
Minopterin ^b (5 µmol/L)	43	72
Cyclophosphamide (1500 µmol/L)	0	0
Cytosine	0	0
DAMPA ^b (5 µmol/L)	46	83
Dihydrofolic Acid	0	0
5-Fluorouracil (3000 µmol/L)	0	0
Folic Acid	0	0
Folinic Acid (Leucovorin)	0	0
7-Hydroxymethotrexate	0	0
6-Mercaptopurine	0	0
Methopterin	0	0
5-Methyltetrahydrofolic Acid	0	0
DL-8-Methyl-5,6,7,8-Tetrahydropterine	0	0
Prednisolone	0	0
Pyrimethamine	0	0
Sulfamethoxazole (1600 µmol/L)	0	0
Tetrahydrofolic Acid	0	0
Triamterene	0	0
Trimethoprim	0	0
Vinblastine	0	0
Vincristine	0	0

^a % Cross-reactivity = $\frac{\text{Mean Test Concentration} - \text{Mean Control Concentration}}{\text{Therapeutic Concentration}} \times 100$

^b The highest concentration that could be quantitated was 5 µmol/L.

Potentially Interfering Substances

A study was performed based on guidance from the CLSI document EP7-A2.²⁸ Potentially interfering substances were evaluated to determine whether methotrexate concentrations were affected when using the ARCHITECT Methotrexate assay. The substances listed below were spiked into samples with methotrexate concentrations at the medical decision points of approximately 0.050 µmol/L and 1.000 µmol/L. The samples were assayed and the methotrexate concentrations of the spiked samples were compared to reference samples. The data are summarized in the following table.

Potentially Interfering Substances	Interferent Concentration	% Difference ^a	
		Methotrexate Concentration 0.050 µmol/L	1.000 µmol/L
Bilirubin (Unconjugated)	60 mg/dL	4	4
Bilirubin (Conjugated)	60 mg/dL	2	3
Cholesterol	700 mg/dL	4	3
Hemoglobin	1000 mg/dL	4	-3
Human Anti-Mouse Antibodies	1500 ng/mL	6	3
Monoclonal IgM	0.386 g/dL	7	-5
Rheumatoid Factor	1500 IU/mL	-4	8
Total Protein (albumin)	12 g/dL	-7	-6
Total Protein (gamma globulin)	12 g/dL	-5	-6
Triglycerides	3000 mg/dL	4	-3
Uric Acid	30 mg/dL	-4	2

^a % Difference = $\frac{\text{Mean Test Concentration} - \text{Mean Control Concentration}}{\text{Mean Control Concentration}} \times 100$

Method Comparison

The ARCHITECT Methotrexate assay is designed to have a slope of 1.00 ± 0.10 and a correlation coefficient (r) of ≥ 0.95 for specimens ranging from 0.040 µmol/L to 1.000 µmol/L methotrexate when compared to the TDx method and ranging from 0.040 µmol/L to 1.500 µmol/L methotrexate when compared to the Liquid Chromatography-Tandem Mass Spectrometry (LC/MS/MS) method.

Studies were performed using serum specimens based on guidance from the CLSI document EP9-A3²⁹ using the Passing-Bablok regression method to compare the ARCHITECT Methotrexate assay to the TDx Methotrexate method and to the LC/MS/MS method. The data are summarized in the following tables including data from a study with specimens above the measuring interval requiring dilution.

ARCHITECT Methotrexate vs. TDx Methotrexate (N = 86)						
Concentration Range (µmol/L)		Correlation Coefficient				
ARCHITECT	TDx	(r)	Intercept	95% CI ^a	Slope	95% CI ^a
0.040-0.993	0.05-0.99	0.9982	0.005	(0.001,0.009)	0.948	(0.924,0.970)

^a CI = Confidence Interval

ARCHITECT Methotrexate vs. LC/MS/MS (N = 101)						
Concentration Range (µmol/L)		Correlation Coefficient				
ARCHITECT	LC/MS/MS	(r)	Intercept	95% CI ^a	Slope	95% CI ^a
0.040-1.438	0.028-1.588	0.9980	0.016	(0.011,0.022)	0.923	(0.908,0.944)

^a CI = Confidence Interval

ARCHITECT Methotrexate vs. TDx Methotrexate (N = 119)						
Concentration Range (µmol/L)		Correlation Coefficient				
ARCHITECT	TDx	(r)	Intercept	95% CI ^a	Slope	95% CI ^a
0.040-888.000	0.05-855.50	0.9983	-0.004	(-0.007,-0.001)	1.018	(0.995,1.037)

^a CI = Confidence Interval

BIBLIOGRAPHY

1. Cole PD, Kamen BA, Bertino JR. Folate Antagonists. In: Hong WK, Bast RC Jr, Hait WN, *et al.*, editors. *Holland-Frei Cancer Medicine*. 8th ed. Shelton, CT: People's Medical Publishing House; 2010:611-620.
2. Methotrexate Injection USP [package insert]. New York, NY: Pfizer Labs; 2012.
3. Tung JP, Maibach HI. The practical use of methotrexate in psoriasis. *Drugs* 1990;40(5):697-712.
4. Guss S, Portnoy J. Methotrexate treatment of severe asthma in children. *Pediatrics* 1992;89(4):635-639.
5. Brooks PJ, Sprull WJ, Parish RC, *et al.* Pharmacokinetics of methotrexate administered by intramuscular and subcutaneous injections in patients with rheumatoid arthritis. *Arthritis Rheum* 1990;33(1):91-94.
6. Baughman RP, Lower EE. The effect of corticosteroid or methotrexate therapy on lung lymphocytes and macrophages in sarcoidosis. *Am Rev Respir Dis* 1990;142(6, pt 1):1268-1271.
7. Hosenpud JD, Hershberger RE, Ratkovec RR, *et al.* Methotrexate for the treatment of patients with multiple episodes of acute cardiac allograft rejection. *J Heart Lung Transplant* 1992;11(4, pt 1):739-745.
8. Ringden O, Klaesson S, Sundberg B, *et al.* Decreased incidence of graft-versus-host disease and improved survival with methotrexate combined with cyclosporin compared with monotherapy in recipients of bone marrow from donors other than HLA identical siblings. *Bone Marrow Transplant* 1992;9(1):19-25.
9. Sæter G, Alvegård TA, Elomaa I, *et al.* Treatment of osteosarcoma of the extremities with the T-10 protocol, with emphasis on the effects of preoperative chemotherapy with single-agent high-dose methotrexate: A Scandinavian Sarcoma Group Study. *J Clin Oncol* 1991;9(10):1766-1775.
10. Abramowitz M, Ochs J, Pui C-H, *et al.* High-dose methotrexate improves clinical outcome in children with acute lymphoblastic leukemia: St. Jude Total Therapy Study X. *Med Pediatr Oncol* 1988;16(5):297-303.
11. Hann IM, Eden OB, Barnes J, *et al.* 'MACHO' chemotherapy for stage IV B cell lymphoma and B cell acute lymphoblastic leukaemia of childhood. *Br J Haematol* 1990;78(3):359-364.
12. Wheeler CA, Shulman LN, Ervin T, *et al.* Cisplatin, continuous-infusion 5-fluorouracil, and intermediate-dose methotrexate in the treatment of unresectable non-small cell carcinoma of the lung. *Cancer* 1991;67(4):892-895.
13. Powles TJ, Jones AL, Judson IR, *et al.* A randomised trial comparing combination chemotherapy using mitomycin C, mitozantrone and methotrexate (3M) with vincristine, anthracycline and cyclophosphamide (VAC) in advanced breast cancer. *Br J Cancer* 1991;64(2):406-410.
14. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
15. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
16. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.

17. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.
18. VORAXAZE (glucarpidase) [package insert]. Brentwood, TN: BTG International Inc.; 2012.
19. Buchen S, Ngampolo D, Melton RG, *et al*. Carboxypeptidase G₂ rescue in patients with methotrexate intoxication and renal failure. *Br J Cancer* 2005;92(3):480-487.
20. Al-Turkmani MR, Law T, Naria A, *et al*. Difficulty measuring methotrexate in a patient with high-dose methotrexate-induced nephrotoxicity. *Clin Chem* 2010;56(12):1792-1796.
21. Schroff RW, Foon KA, Beatty SM, *et al*. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
22. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
23. Young RC, Chabner BA. An in vivo method for monitoring differential effects of chemotherapy on target tissues in animals and man: Correlation with plasma pharmacokinetics. *J Clin Invest* 1973;52(6) abstract 340: 92a.
24. Nirenberg A, Mosende C, Mehta BM, *et al*. High-dose methotrexate with citrovorum factor rescue: Predictive value of serum methotrexate concentrations and corrective measures to avert toxicity. *Cancer Treat Rep* 1977;61(5):779-783.
25. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
26. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. NCCLS Document EP6-A. Wayne, PA: NCCLS; 2003.
27. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
28. Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP7-A2. Wayne, PA: CLSI; 2005.
29. Clinical and Laboratory Standards Institute (CLSI). *Method Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP9-A3. Wayne, PA: CLSI; 2013.

ARCHITECT, TDx and Chemiflex are trademarks of Abbott Laboratories in various jurisdictions.

All other trademarks are property of their respective owners.

ABBOTT
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Fujirebio Diagnostics Inc.
Malvern, PA 19355, USA

Distributed by Abbott Laboratories
Abbott Park, IL 60064 USA

September 2014
©2014 Abbott Laboratories

ADDENDUM for Russia

to the package insert of ARCHITECT Methotrexate Reagent Kit

SHIPPING CONDITIONS

ARCHITECT Methotrexate Reagent Kit should be shipped in accordance with the storage requirements mentioned in the package insert: vertical position of reagent bottles, shipping temperature under 2-8°C.

WARRANTY OBLIGATIONS

The manufacturer guarantees the quality of the valid results of the test for the period of time mentioned on the package label as expiry date. This date correspond to 8 months from the date of Reagent Kit manufacturing on condition that all requirements of this package insert are followed including the storage of ARCHITECT Methotrexate Reagent Kit under 2-8°C.

ARCHITECT SYSTEM

Methotrexate Calibrators

INTENDED USE

The ARCHITECT Methotrexate Calibrators are for the calibration of the ARCHITECT *i* System when used for the quantitative determination of methotrexate in human serum and plasma.

Refer to the ARCHITECT Methotrexate reagent package insert and ARCHITECT System Operations Manual for additional information.

CONTENTS

6 Bottles (4.0 mL each) of ARCHITECT Methotrexate Calibrators. Calibrators A through F are prepared with recalcified human plasma. Calibrators B through F contain USP Grade Methotrexate. Preservatives: sodium azide and ProClin 950.

The calibrators are at the following concentrations:

Calibrator	Concentration CONC	
	μmol/L	μg/mL
CAL A	0.000	0.000
CAL B	0.040	0.018
CAL C	0.100	0.045
CAL D	0.250	0.114
CAL E	0.700	0.318
CAL F	1.500	0.682

STANDARDIZATION

Internal reference standards are manufactured for ARCHITECT Methotrexate using Methotrexate Reference Standard (USP). Methotrexate calibrators are manufactured gravimetrically and tested against these internal reference standards.

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

-  **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.
- The calibrators contain recalcified human plasma that was tested and found to be nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HIV-1/HIV-2, and anti-HCV.

- The following warnings and precautions apply to these components:
 - Calibrators A through F

WARNING: Contains methylisothiazolones and sodium azide.
H317 May cause an allergic skin reaction.
EUH032 Contact with acids liberates very toxic gas.

Prevention

P261 Avoid breathing mist / vapors / spray.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves / protective clothing / eye protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of water.
P333+P313 If skin irritation or rash occurs: Get medical advice / attention.
P362+P364 Take off contaminated clothing and wash it before reuse.

Disposal

P501 Dispose of contents / container in accordance with local regulations.

- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

STORAGE

- ARCHITECT Methotrexate Calibrators are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- ARCHITECT Methotrexate Calibrators must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.

QUALITY CONTROL PROCEDURES

A single sample of each Low, Medium, and High Control must be tested to evaluate the assay calibration. For information on ordering controls, refer to the ARCHITECT System Operations Manual, Section 5.

- Ensure that assay control values are within the ranges specified in the control package insert.

Once an ARCHITECT Methotrexate calibration is accepted and stored all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used
- Controls are out of range

PREPARATION FOR USE

- The ARCHITECT Methotrexate Calibrators may be used immediately after removal from 2-8°C storage.
- Mix the calibrators by gentle inversion before use.
- To perform a calibration, test the calibrators in duplicate. T calibrators should be priority loaded.
- To obtain the recommended volume requirements for the calibrator hold the bottles **vertically** and dispense a minimum of 5 drops of e calibrator into each respective sample cup.
- After each use, tightly close the caps and return the calibrator 2-8°C storage.
- For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.

en

thotrexate
P49-01
490
7/14

LIOGRAPHY

US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.

ents:
ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. All trademarks are property of their respective owners.

 **ABBOTT**
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Fujirebio Diagnostics Inc.,
Malvern, PA 19355, USA


DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

July 2014
©2014 Abbott Laboratories

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Store at 2-8°C
	Use by/Expiration date
	Calibrators A, B, C, D, E, and F
	Concentration
	Contains sodium azide. Contact with acids liberates very toxic gas.
	Distributed in the USA by
	Global Trade Item Number
	Information needed for United States of America only
	In Vitro Diagnostic Medical Device
	Lot Number
	Produced for Abbott by
	Product of USA
	List Number
	For use by or on the order of a physician only (applicable to USA classification only)
	Warning: May cause an allergic reaction.

ADDENDUM for Russia

to the package insert of ARCHITECT Methotrexate Calibrators

SHIPPING CONDITIONS

ARCHITECT Methotrexate Calibrators should be shipped in accordance with the storage requirements mentioned in the package insert: vertical position of reagent bottles, shipping temperature under 2-8°C.

WARRANTY OBLIGATIONS

The manufacturer guarantees the quality of the valid results of the test for the period of time mentioned on the package label as expiry date. This date correspond to 8 months from the date of Calibrators manufacturing on condition that all requirements of this package insert are followed including the storage of ARCHITECT Methotrexate Calibrators under 2-8°C.

ARC

lethr

NTENDE
The ARC
the ester
deviation
determi
Method
perfor
Refer
ARC
CC
4 f
M
p

Methotrexate Controls

en

Methotrexate

REF 2P49-10

C2P490

611-016 9/14

INTENDED USE

The ARCHITECT Methotrexate Low, Medium, and High Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT *i* System when used for the quantitative determination of methotrexate in human serum and plasma. ARCHITECT Methotrexate Control X is used to provide assurance of accuracy when performing manual sample dilutions.

Refer to the ARCHITECT Methotrexate reagent package insert and ARCHITECT System Operations Manual for additional information.

CONTENTS

4 Bottles (8.0 mL each) of ARCHITECT Methotrexate Controls. Controls L, M, H, and X are prepared with USP Grade Methotrexate added to recalcified human plasma. Preservatives: sodium azide and ProClin 950. The controls are at the following concentrations:

Control	Target Concentration / Range			
	µmol/L		µg/mL	
	CONC	RANGE	CONC	RANGE
CONTROL L	0.070	0.048 - 0.092	0.032	0.022 - 0.042
CONTROL M	0.450	0.304 - 0.596	0.204	0.138 - 0.271
CONTROL H	1.000	0.675 - 1.325	0.454	0.307 - 0.602
CONTROL X	10.000	6.750 - 13.250	4.544	3.067 - 6.021

Each laboratory should establish its own concentration ranges for new control lots at each control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days. Sources of variation that can be expected should be included in this study in order to be representative of future system performance. These may include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules
- Data points collected at different times of the day
- Multiple technicians performing manual dilution procedure

These results should be applied to your laboratory's quality control practices.

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

CAUTION: This product contains human sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.

- The controls contain recalcified human plasma that was tested and found to be nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HIV-1/HIV-2, and anti-HCV.

- The following warnings and precautions apply to these components:

- Controls L, M, H, and X

WARNING: Contains methylisothiazolones and sodium azide.

H317 May cause an allergic skin reaction.

EUH032 Contact with acids liberates very toxic gas.

Prevention

P261 Avoid breathing mist / vapors / spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves / protective clothing / eye protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of water.

P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

P362+P364 Take off contaminated clothing and wash it before reuse.

Disposal

P501 Dispose of contents / container in accordance with local regulations.


- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

Handling Precautions

- **Caution:** In order to prevent cross contamination, do not remove dropper tips from the ARCHITECT Methotrexate Control bottles. Use caution when handling control bottles and caps.

STORAGE

- ARCHITECT Methotrexate Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- ARCHITECT Methotrexate Controls must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.

QUALITY CONTROL PROCEDURES

The recommended control requirement for the ARCHITECT Methotrexate assay is that a single sample of each Low, Medium, and High Control be tested once every 24 hours each day of use. It is recommended that Control X be diluted and tested with diluted samples. For detailed dilution procedure instructions and quality control, refer to the ARCHITECT Methotrexate reagent package insert, Specimen Dilution Procedures section.

If your laboratory quality control procedures require more frequent use of controls to verify test results, follow those procedures. Additional controls may be tested in conformance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy. The control values must be within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated specimen test results are invalid and must be retested. Recalibration may be indicated. If a diluted control is out of its specified range, the associated specimen dilution test results are invalid and must be retested.

PREPARATION FOR USE

- The ARCHITECT Methotrexate Controls may be used immediately after removal from 2-8°C storage.
- Mix the controls by gentle inversion before use.
- To obtain the recommended volume requirements for the Low, Medium, and High Controls, hold the bottles **vertically** and dispense a minimum of 5 drops of the controls into each respective sample cup.
- The ARCHITECT Methotrexate Control X must be diluted per the ARCHITECT Methotrexate reagent package insert, Specimen Dilution Procedures section, and then run as a control with diluted samples.
 - Replace pipette tips between serial dilution steps.
- Replace gloves if contamination is suspected.
- After each use, tightly close the caps and return the controls to 2-8°C storage.
- For information on ordering controls, refer to the ARCHITECT System Operations Manual, Section 5.
- To troubleshoot control values that fall outside the control range, refer to the ARCHITECT System Operations Manual, Section 10.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. All trademarks are property of their respective owners.

ABBOTT
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Fujirebio Diagnostics Inc.,
Malvern, PA 19355, USA

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

September 2014
©2014 Abbott Laboratories

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Store at 2-8°C
	Use by/Expiration date
	Concentration
	Contains sodium azide. Contact with acids liberates very toxic gas.
	Controls L, M, H, and X
	Distributed in the USA by
	Global Trade Item Number
	Information needed for United States of America only
	In Vitro Diagnostic Medical Device
	Lot Number
	Produced for Abbott by
	Product of USA
	Range
	List Number
	For use by or on the order of a physician only (applicable to USA classification only)
	Warning: May cause an allergic reaction.

WARRANTY for R
the package ins
IPPING COND
ARCHITECT Ma
in the storage
position of reas
WARRANTY C
The manufac
the period c
date corre
on contri
including

ADDENDUM for Russia

the package insert of ARCHITECT Methotrexate Controls

SHIPPING CONDITIONS

ARCHITECT Methotrexate Controls should be shipped in accordance with the storage requirements mentioned in the package insert: vertical position of reagent bottles, shipping temperature under 2-8°C.

WARRANTY OBLIGATIONS

The manufacturer guarantees the quality of the valid results of the test for the period of time mentioned on the package label as expiry date. This date correspond to 8 months from the date of Controls manufacturing on condition that all requirements of this package insert are followed including the storage of ARCHITECT Methotrexate Controls under 2-8°C.

ARCHITECT

SYSTEM

Methotrexate Extended Range Controls

INTENDED USE

The ARCHITECT Methotrexate Extended Range Controls are used to provide assurance of accuracy when performing manual sample dilutions. Refer to the ARCHITECT Methotrexate reagent package insert and ARCHITECT System Operations Manual for additional information.

CONTENTS

2 Bottles (8.0 mL each) of ARCHITECT Methotrexate Extended Range Controls. Controls Y and Z are prepared with USP Grade Methotrexate added to recalcified human plasma. Preservatives: sodium azide and ProClin 950.

The controls are at the following concentrations:

Control	Target Concentration / Range			
	$\mu\text{mol/L}$		$\mu\text{g/mL}$	
	CONC	RANGE	CONC	RANGE
CONTROL Y	50.000	30.000 - 70.000	22.722	13.633 - 31.811
CONTROL Z	500.000	320.000 - 700.000	227.220	145.421 - 318.108

Each laboratory should establish its own concentration ranges for new control lots at each control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days. Sources of variation that can be expected should be included in this study in order to be representative of future system performance. These may include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules
- Data points collected at different times of the day
- Multiple technicians performing manual dilution procedure

These results should be applied to your laboratory's quality control practices.

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use
- **Rx ONLY**

-  **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.

- The controls contain recalcified human plasma that was tested and found to be nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HIV-1/HIV-2, and anti-HCV.

Methotrexate

REF 2P49-15

C2P4Q0

611-046 9/14

- The following warnings and precautions apply to these components:
 - Extended Range Controls Y and Z

WARNING: Contains methylisothiazolones and sodium azide.

H317 May cause an allergic skin reaction.

EUH032 Contact with acids liberates very toxic gas.

Prevention

P261 Avoid breathing mist / vapors / spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves / protective clothing / eye protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of water.

P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

P362+P364 Take off contaminated clothing and wash it before reuse.

Disposal

P501 Dispose of contents / container in accordance with local regulations.

- Safety Data Sheets are available at www.abbottddiagnostics.com or contact your local representative.

Handling Precautions

- **Caution:** In order to prevent cross contamination, do not remove dropper tips from the ARCHITECT Methotrexate Extended Range Control bottles. Use caution when handling control bottles and caps.

STORAGE

- ARCHITECT Methotrexate Extended Range Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.
- ARCHITECT Methotrexate Extended Range Controls must be stored at 2-8°C in an upright position and may be used immediately after removal from 2-8°C storage.

2-8°C

QUALITY CONTROL PROCEDURES

The recommended control requirement for the ARCHITECT Methotrexate assay is that a single sample of each Low, Medium, and High Control be tested once every 24 hours each day of use. It is recommended that Control Y and Control Z be diluted and tested with diluted samples. For detailed dilution procedure instructions and quality control, refer to the ARCHITECT Methotrexate reagent package insert, Specimen Dilution Procedures section.

If your laboratory quality control procedures require more frequent use of controls to verify test results, follow those procedures. Additional controls may be tested in conformance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

The control values must be within the acceptable ranges specified in the control package insert. If a diluted control is out of its specified range, the associated specimen dilution test results are invalid and must be retested.

en
Methotrexate
2P40-15
4Q0
1/14
nts:

PREPARATION FOR USE

The ARCHITECT Methotrexate Extended Range Controls may be used immediately after removal from 2-8°C storage.

Mix the controls by gentle inversion before use.

The ARCHITECT Methotrexate Control Y and Control Z must be diluted per the ARCHITECT Methotrexate reagent package insert, Specimen Dilution Procedures section, and then run as controls with diluted samples.

- Replace pipette tips between serial dilution steps.
- Replace gloves if contamination is suspected.
- After each use, tightly close the caps and return the controls to 2-8°C storage.
- For information on ordering controls, refer to the ARCHITECT System Operations Manual, Section 5.
- To troubleshoot control values that fall outside the control range, refer to the ARCHITECT System Operations Manual, Section 10.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. All trademarks are property of their respective owners.

 **ABBOTT**
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

PRODUCED FOR ABBOTT BY

Fujirebio Diagnostics Inc.,
Malvern, PA 19355, USA

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

September 2014
©2014 Abbott Laboratories

Key to Symbols

Caution

Consult instructions for use

Manufacturer

Store at 2-8°C

Use by/Expiration date

CONC

Concentration

CONTAINS: AZIDE

Contains sodium azide. Contact with acids liberates very toxic gas.

CONTROL Y **CONTROL Z**

Controls Y and Z

DISTRIBUTED IN THE USA BY

Distributed in the USA by

GTIN

Global Trade Item Number

INFORMATION FOR USA ONLY

Information needed for United States of America only

IVD

In Vitro Diagnostic Medical Device

LOT

Lot Number

PRODUCED FOR ABBOTT BY

Produced for Abbott by

PRODUCT OF USA

Product of USA

RANGE

Range

REF

List Number

Rx ONLY

For use by or on the order of a physician only (applicable to USA classification only)

WARNING: SENSITIZER

Warning: May cause an allergic reaction.

ADDENDUM for Russia

to the package insert of ARCHITECT[®] Methotrexate Extended Range Controls

SHIPPING CONDITIONS

ARCHITECT[®] Methotrexate Extended Range Controls should be shipped in accordance with the storage requirements mentioned in the package insert: vertical position of reagent bottles, shipping temperature under 2-8°C.

WARRANTY OBLIGATIONS

The manufacturer guarantees the quality of the valid results of the test for the period of time mentioned on the package label as expiry date. This date correspond to 8 months from the date of Extended Range Controls manufacturing on condition that all requirements of this package insert are followed including the storage of ARCHITECT[®] Methotrexate Extended Range Controls under 2-8°C.

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1. Cole PD, Kamen BA, Bertino JR. Folate Antagonists. In: Hong WK, Bast RC Jr, Hait WN, *et al.*, editors. *Holland-Frei Cancer Medicine*. 8th ed. Shelton, CT: People's Medical Publishing House; 2010:611-620.
2. Methotrexate Injection USP [package insert]. New York, NY: Pfizer Labs; 2012.
3. Tung JP, Maibach HI. The practical use of methotrexate in psoriasis. *Drugs* 1990;40(5):697-712.
4. Guss S, Portnoy J. Methotrexate treatment of severe asthma in children. *Pediatrics* 1992;89(4):635-639.
5. Brooks PJ, Sprull WJ, Parish RC, *et al.* Pharmacokinetics of methotrexate administered by intramuscular and subcutaneous injections in patients with rheumatoid arthritis. *Arthritis Rheum* 1990;33(1):91-94.
6. Baughman RP, Lower EE. The effect of corticosteroid or methotrexate therapy on lung lymphocytes and macrophages in sarcoidosis. *Am Rev Respir Dis* 1990;142(6, pt 1):1268-1271.
7. Hosenpud JD, Hershberger RE, Ratkovec RR, *et al.* Methotrexate for the treatment of patients with multiple episodes of acute cardiac allograft rejection. *J Heart Lung Transplant* 1992;11(4, pt 1):739-745.
8. Ringdén O, Klaesson S, Sundberg B, *et al.* Decreased incidence of graft-versus-host disease and improved survival with methotrexate combined with cyclosporin compared with monotherapy in recipients of bone marrow from donors other than HLA identical siblings. *Bone Marrow Transplant* 1992;9(1):19-25.
9. Sæter G, Alvegård TA, Elomaa I, *et al.* Treatment of osteosarcoma of the extremities with the T-10 protocol, with emphasis on the effects of preoperative chemotherapy with single-agent high-dose methotrexate: A Scandinavian Sarcoma Group Study. *J Clin Oncol* 1991;9(10):1766-1775.
10. Abramowitch M, Ochs J, Pui C-H, *et al.* High-dose methotrexate improves clinical outcome in children with acute lymphoblastic leukemia: St. Jude Total Therapy Study X. *Med Pediatr Oncol* 1988;16(5):297-303.
11. Hann IM, Eden OB, Barnes J, *et al.* 'MACHO' chemotherapy for stage IV B cell lymphoma and B cell acute lymphoblastic leukaemia of childhood. *Br J Haematol* 1990;76(3):359-364.
12. Wheeler CA, Shulman LN, Ervin T, *et al.* Cisplatin, continuous-infusion 5-fluorouracil, and intermediate-dose methotrexate in the treatment of unresectable non-small cell carcinoma of the lung. *Cancer* 1991;67(4):892-895.
13. Powles TJ, Jones AL, Judson IR, *et al.* A randomised trial comparing combination chemotherapy using mitomycin C, mitozantrone and methotrexate (3M) with vincristine, anthracycline and cyclophosphamide (VAC) in advanced breast cancer. *Br J Cancer* 1991;64(2):406-410.
14. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
15. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
16. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
17. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.
18. VORAXAZE (glucarpidase) [package insert]. Brentwood, TN: BTG International Inc.; 2012.
19. Buchen S, Ngampolo D, Melton RG, *et al.* Carboxypeptidase G₂ rescue in patients with methotrexate intoxication and renal failure. *Br J Cancer* 2005;92(3):480-487.
20. Al-Turkmani MR, Law T, Naria A, *et al.* Difficulty measuring methotrexate in a patient with high-dose methotrexate-induced nephrotoxicity. *Clin Chem* 2010;56(12):1792-1796.
21. Schroff RW, Foon KA, Beatty SM, *et al.* Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
22. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
23. Young RC, Chabner BA. An in vivo method for monitoring differential effects of chemotherapy on target tissues in animals and man: Correlation with plasma pharmacokinetics. *J Clin Invest* 1973;52(6) abstract 340: 92a.
24. Nirenberg A, Mosende C, Mehta BM, *et al.* High-dose methotrexate with citrovorum factor rescue: Predictive value of serum methotrexate concentrations and corrective measures to avert toxicity. *Cancer Treat Rep* 1977;61(5):779-783.
25. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
26. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. NCCLS Document EP6-A. Wayne, PA: NCCLS; 2003.
27. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
28. Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP7-A2. Wayne, PA: CLSI; 2005.
29. Clinical and Laboratory Standards Institute (CLSI). *Method Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP9-A3. Wayne, PA: CLSI; 2013.

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- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.

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