Medtronic

DLP®

Femoral Arterial and Venous Cannula/Placement Kits
Canule artérielle et veineuse fémorale/Kits de positionnement
Femoralarterien- und -venenkanule/Positionierungskits
Canula venosa y arterial femoral/Equipos de colocación
Femorale arteriele en veneuze canule/plaatsingskits
Cannula arteriosa e venosa femorale/Kit di posizionamento
Femoral arterie- og venekanyler / innføringssett
Reisivaltimo- ja -laskimokanyyli / asetuspakkaukset
Sats för placering av femoral artar- och venkanyler
Μηριαία αρτηριακή και φλεβική κανούλα/Kit τοποθέτησης
Femoral arterie- og venekanyler/anlæggelsessæt
Udowa kaniula tętnicza i żylna/zestawy do wprowadzania
Canula arterial e venosa femoral/Kits de colocação
Femoral Arteriyel ve Venöz Kanul/Yerleştirme Kitleri
Femorální arteriální a venózní kanyla / sada pro umístění
Femoralis artérias es vénas kanul/benhelyezőkészlet
Бедренная артериальная и венозная канюля/Комплекты для установки

As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of the
original document. The original is filed at Medtronic Inc.

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432, USA
Tel: +1 763 514 4000

Handwritten: Aug 3, 2017

Handwritten: Reg Affairs Spec.

Handwritten: Gayatri Ghadge

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Käyttöohjeet
Bruksanvisning • Οδηγίες χρήσης • Bruksanvisning • Instrukcja użytkowania
Instruções de utilização • Kullanım Talimatları • Návod k použití • Használati útmutató
Инструкция по эксплуатации

USA Caution: Federal law (USA) restricts this device to sale by or on the
order of a physician.

2011/SEP/14 at 9:48 a.m. Doc number: M934624A001 Rev. A [NO17]
Printspec c - 5.5 x 8.5 inches

Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf dem Verpackungsetikett / Explicación de los símbolos que aparecen en el etiquetado del envase / Verklaring van de symbolen op de verpakkingslabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på emballasjen / Pakkauksen symbolien selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στην ετικέτα της συσκευασίας / Forklaring af symboler på emballagens mærkater / Objasnienie symboli znajdujących się na etykietach opakowania / Explicação dos símbolos nas etiquetas da embalagem / Ambalaj etiketindeki sembollerin açıklamaları / Vysvětlení symbolů uvedených na štítcích balení / A csomagoláson látható szimbólumok jelentése / Объяснение символов на этикетках упаковки

Refer to the outer package label to see which symbols apply to this product / Se référer à l'étiquette sur l'emballage extérieur pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt gelten, entnehmen Sie bitte dem Etikett auf der Verpackung / Consulte la etiqueta del envase exterior para comprobar que símbolos son aplicables a este producto / Controleer het label op de buitenverpakking om te zien welke symbolen op dit product van toepassing zijn / Vedere l'etichetta sulla confezione esterna per sapere quali simboli si applicano a questo prodotto / Se på etiketten utenpå pakningen for å fastslå hvilke symboler som gjelder for dette produktet / Katso ulkopakkauksesta, mitkä symbolit koskevat tätä tuotetta / Se etiketten på ytterförpackningen för de symboler som gäller denna produkt / Ανοίξετε στην ετικέτα της εξωτερικής συσκευασίας για να δείτε ποια σύμβολα αφορούν αυτό το προϊόν / Se mærkaten på den ydre emballage for de symboler, der gælder for dette produkt / Należy zapoznać się z etykietą na zewnętrznym opakowaniu, aby określić, które symbole dotyczą niniejszego produktu / Consultar as etiquetas exteriores da embalagem para ver que símbolos se aplicam a este produto / Hangi szimbólok bu ürün için geçerli olduğunu görmek için ambalaj etiketine başvurun / Podle štítků na vnější straně balení zjistíte, které symboly se týkají tohoto produktu / A termékre vonatkozó szimbólumok a külső csomagoláson találhatók / См. на упаковке, какие символы применимы к данному продукту

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DLP®

Femoral Arterial and Venous Cannula/Placement Kits

Product Description

The cannula and dilator components are polyurethane and have tapered proximal tips to aid in diameter transition from the dilator. The dilator components utilize tip indicator bands for minimum tip alignment of the component parts. The cannula body is marked for depth insertion. The components are designed to fit over each other to progressively dilate the vessel during insertion. The dilator allows a 0.1-cm (0.038-in) guidewire to pass through it to aid in insertion. A soft, tapered transition fitting terminates in a 0.95-cm (3/8-in) connection site or barbed connector.

Sterile, nonpyrogenic, single use, radiopaque.

Indications for Use

These cannulae are intended for rapid femoral venous and arterial access for short-term cardiopulmonary bypass. This product is intended for use up to six hours or less.

Contraindications

Do not use if the patient has severe peripheral atherosclerosis or severe arterial dissection. Use of these cannulae in smaller sized patients (ie, children and small women) is contraindicated due to ischemic risk.

This device is not intended for use except as indicated above.

Warnings

Bending, kinking or clamping the white polyurethane cannula tubing may weaken the cannula and cause collapse during its use. Only use clamp at recommended clamp site (Figure 1).

Patient and procedure selection is a medical responsibility and the outcome is dependent on many variables, including patient pathology, surgical and perfusion procedures.

These cannulae are for short-term, temporary use only.

Ischemia may occur distal to the cannula access site. Measures (such as cannulating the distal portion of the femoral artery) should be taken to avoid ischemic injury.

Insure that the appropriate cannula has entered the appropriate vessel using blood pressures or blood samples as indicators. When connecting the cannula to the venous or arterial line a tubing clamp may be required to clamp the cannula at the clamp site.

Removal should always be done under direct visualization to repair the opening in the vessel.

Precautions

Note: Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Surgeons must evaluate the appropriateness of the procedure based on their own medical training and experience, and the type of surgical procedure.

Please tighten luer caps before using.

These cannulae are designed and intended for single use only. DO NOT REUSE.

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Instructions for Use

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the cannula onto the sterile field utilizing an aseptic technique.
2. The method of access is left to the discretion of the physician. Insertion under direct vision via cut down or standard Seldinger technique utilizing a series of dilators is recommended.
3. Prior to insertion, the cannula components must be assembled and properly aligned. Moisten each component with sterile saline. For models 96328 and 96428, insert the dilator into the sheath. Insert the sheath into the cannula. Minimum tip alignment of each component occurs when the black indicator band on the proximal end of the dilator is positioned at the beveled end of the sheath, and the black indicator band on the proximal end of the sheath is positioned at the end of the 0.95-cm (3/8-in) connector site of the cannula (Figure 2). Minimum tip alignment is necessary to facilitate insertion.
4. Prior to insertion, position the venous cannula over the patient and determine approximate depth of insertion by noting the relative distance from the right atrium to the femoral vein insertion site. Depth markings will aid in determining appropriate insertion depth.

5. Using the cut down technique or Seldinger technique, insert the guidewire into the femoral vein.
Caution: Insure that the vessel is of adequate size to accommodate the venous cannula. Insert the assembled cannula into the femoral vessel, beginning with the dilator. Advance the dilator sufficiently to allow the cannula to follow the natural pathway up to the desired location.
Note: Insertion technique will vary depending on patient anatomy. Use care during insertion to follow the natural anatomical pathway of the vessel.
6. The markings on the cannula and the cannula dilator indicate 5-cm intervals. Twenty centimeters is indicated by 2 closely placed bands, 30 cm by 3 bands, etc. The bold band near the cannula side holes indicates a 5-cm length between the band and the first side hole. When this band is visible upon removal, the physician is alerted that the side holes will soon be removed from the patient. When the bold band on the dilator is visible just outside of the cannula connector, the cannula has been inserted to the same depth as the dilator.
7. When using Models 96328 and 96428, progressively advance the dilator until the dilator reaches the right atrium. While holding the dilator, advance the dilator sheath up to the right atrium. Insure cannula tip position in the right atrium by referring to the predetermined insertion depth (see step 4). Verify positioning with tactile feel or fluoroscopy. Refer to cannula depth markings to aid in tip positioning.
8. When proper positioning of cannula is verified, slowly withdraw the dilator components and guidewire from the cannula. When the cannula is fully primed, clamp the cannula at the designated clamping site (Figure 3). Place a 0.95-cm x 0.95-cm (3/8-in x 3/8-in) or 0.95-cm x 1.27-cm (3/8-in x 1/2-in) connector into the connection site (Model 96328 only). Connect the cannula to the appropriate line from the extracorporeal circuit.
Warning: Clamping or kinking the white polyurethane cannula tubing may weaken the cannula and cause collapse during its use. CLAMP ONLY AT THE RECOMMENDED CLAMP SITE (Figure 3).
9. Remove clamp prior to initiation of the pump run.
10. Upon completion of cardiopulmonary bypass, remove the cannula from the vessel under direct visualization.
11. Upon removal of the cannula, repair the opening in the vessel using appropriate surgical techniques.
12. Dispose of products properly according to hospital policy with the other blood contaminated products.

The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE FEMORAL ARTERIAL AND VENOUS CANNULA/PLACEMENT KITS, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE FEMORAL ARTERIAL AND VENOUS CANNULA/PLACEMENT KITS, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of

competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Medtronic Inc.
710 Medtronic Pkwy
Minneapolis, MN 55432
USA
1:54

Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, **having** its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby swear (or affirm) that the attached documentation:

- 1- Appendix to IFU - 97 pages;
Are true and unaltered copies of the original document.

Huda Yusuf
Sr. Regulatory Affairs Specialist

Name
Title

September 22, 2017
Date


Diane Marie Larson Sept 22, 2017

Bio-Medicus™

Adult Cannula and Introducer
 Canule et introducteur pour adulte
 Kanüle und Einführhilfe für Erwachsene
 Cânula e introductor para adultos
 Canule voor volwassenen en introducer
 Cannula e Introduttore per adulti
 Kanyle og innfører for voksne
 Aikuisten kanyyli ja ohjain
 Kanyl och introducer för vuxna
 Κάνουλα ενηλίκων και εισαγωγέας
 Kanyle og indfører til voksne
 Kaniula i introduktor dla dorosłych
 Cânula e introdutor para adulto
 Yetişkin Kanül ve İntrodüseri
 Kanyla pro dospěle a zaváděc
 Kanul és bevezető felhőttek számára
 Канюля для взрослых и интродьюсер

As a registered Notary Public for the State of Minnesota,
 I hereby certify that this is a true and accurate copy of the
 original document. The original is filed at Medtronic Inc.

Haw Aug 3, 2017

Medtronic Inc.
 710 Medtronic Parkway
 Minneapolis, MN 55432, USA
 Tel: +1 763 514 4000

Reg. Affairs Spec
 Gayatri Ghadge

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
 Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Kayttoohjeet
 Bruksanvisning • Οδηγίες χρήσης • Brugsanvisning • Instrukcja użytkowania
 Instruções de utilização • Kullanım Talimatları • Návod k použití • Használati útmutató
 Инструкция по эксплуатации

Caution: Federal law (USA) restricts this device to sale by or on the
 order of a physician

master:20140730, alldocMod:20140911, alldocMod2:20140911, alldocStyle:20140911, alldocStyle_graph:20140911,
 master_pageset:20140911, languageLookup:20140514, lable:20140620, characterJoin:20110211 3of9:20110211, 3of9_svg:
 20110211

Trademarks may be registered and are the property of their respective owners.
Les marques commerciales mentionnées peuvent être déposées et appartiennent à leurs propriétaires respectifs.
Marken sind eventuell eingetragen und sind das Eigentum ihrer jeweiligen Inhaber.
Las marcas comerciales pueden estar registradas y pertenecen a sus respectivos propietarios.
Handelsmerken kunnen zijn geregistreerd en zijn het eigendom van de desbetreffende eigenaars.
I marchi di fabbrica possono essere registrati e sono di proprietà dei rispettivi titolari.
Varumärken kan vara registrerade och tillhör respektive ägare.
Tavaramärgid võivad olla registreeritud ja ne kuuluvad vastavalt omajärgi.
Tavumärken kan vara registrerade och tillhör respektive ägare.
Також торгові знаки можуть бути зареєстровані та належати відповідним суб'єктам.
Varumärker kan vara registrerade og tillhører deres respektive ejere.
Znaki towarowe mogą być zastrzeżone i są własnością odpowiednich podmiotów.
As marcas comerciais podem ser registadas e são propriedade dos respectivos detentores.
Ticarî markalar teschi olabir ve işi sahiplerinin mülkiyetindedir.
Ochránené známky (mohou byť registrované) sú majetkom príslušných vlastníkov.
A védjegyek a megfélelő védjegytulajdonosok tulajdonát képezhetik.
Товарные знаки могут быть зарегистрированными и являются собственностью их владельцев.

Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf der Verpackung / Explicación de los símbolos que aparecen en la documentación del envase / Verklaring van de symbolen op de verpakkingslabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på pakningen / Pakkausetikettin symbolien selitykset / Förklaring av symboler på förpackningsmärkningen / Επεξήγηση των συμβόλων στη σήμανση της συσκευασίας / Forklaring til symbolerne på emballagens mærkater / Objaśnienie symboli znajdujących się na etykietach opakowania / Explicação dos símbolos na documentação da embalagem / Ambalaj etiketi üzerindeki sembollerin açıklaması / Vysvětlení symbolů uvedených na štítcích balení / A csomagoláson látható szimbólumok jelentése / Объяснение символов на этикетках упаковки

Refer to the outer package label to see **which** symbols apply to this product / Se referer à l'étiquette sur l'emballage extérieur pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt zutreffen, entnehmen Sie bitte dem Etikett auf der Verpackung. / Consulte la etiqueta del envase exterior para comprobar qué símbolos son aplicables a este producto / Controleer het label op de buitenverpakking om te zien welke symbolen op dit product van toepassing zijn / Fare riferimento all'etichetta sulla confezione esterna per sapere quali simboli si applicano a questo prodotto. / Se på etiketten utenpå pakningen for å fastslå hvilke symboler som gjelder for dette produktet / Katso ulkopakkauksen etiketistä, mitkä symbolit koskevat tätä tuotetta / Se etiketten på ytterförpackningen för vilka symboler som gäller denna produkt / Ανατρέξτε στην ετικέτα της εξωτερικής συσκευασίας για να δείτε ποια σύμβολα ισχύουν για το παρόν προϊόν / Se mærkaten på den ydre emballage for de symboler der gælder for dette produkt / Należy zapoznać się z etykietą na zewnętrznym opakowaniu, aby określić, które symbole dotyczą niniejszego produktu / Consultar a etiqueta exterior da embalagem para ver que símbolos se aplicam a este produto / Bu üründe hangi sembollerin geçerli olduğunu görmek için dış ambalaj etiketine başvurun / Symboly, které se vztahují k tomuto výrobku, naleznete na štítku na vnější straně obalu. / A termékre vonatkozó szimbólumok a külső csomagoláson találhatók. / См на внешней этикетке упаковки, какие символы применимы к данному продукту.

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1. Product Description

Adult cannula and introducer.

Sterile, nonpyrogenic, disposable, single use only.

2. Indications for Use

These devices are to be used by a trained physician only. Cannulae are used to cannulate vessels, perfuse vessels or organs, and/or connect with accessory extracorporeal equipment. The cannula introducer is intended to facilitate proper insertion and placement of the appropriately-sized cannula within the vessel for cardiopulmonary bypass. The Bio-Medicus™ cannula (18-cm [7.09-in] tip length models) may be used in either the femoral position as an arterial delivery cannula or in the jugular position as a venous return cannula. This product is intended for use up to 6 hours.

3. Contraindications

Alone, the cannula and introducer are not medical treatment devices. The cannula introducer is only to be used with the appropriately-sized Bio-Medicus cannula. These devices are not intended for use, except as indicated above. Do not use if the patient has severe peripheral atherosclerosis or severe arterial dissection.

4. Warnings and Precautions

- Read all warnings, precautions, and Instructions for Use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- Verify Use By date.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness, or death.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken in any way.
- Do not use alcohol or alcohol-based fluids for lubrication.
- Thoroughly inspect the cannula prior to use to verify that the lumen and side holes (if applicable) are patent and that the cannula has not been damaged or kinked.
- Use caution when assembling or manipulating the introducer within the cannula, or moving the suture ring.
- Avoid false lumen cannulation.
- Never advance, position, or place the cannula without the introducer as this may kink and damage the cannula.
- Use care to protect the side hole areas from kinking as they are passed into the vessel.
- Manipulate cannula in such a manner as to prevent kinks or any restrictions that may alter flow.
- Manipulation of the heart with cannula in place may cause flow restrictions.
- Avoid kinking when attaching tubing and when anchoring or untethering the cannula body.
- Caution should be used when attaching the cannula to other devices in the extracorporeal circuit so as not to damage the device being connected or the cannula.
- Clamping too close to connector may damage connector end.
- Do not clamp wire-reinforced sections. Clamping on reinforced sections may result in permanent wall distortion and/or lumen collapse.
- Since the cannula wall is opaque due to the presence of the spring, be sure to completely remove the air from the cannula during the priming procedure.
- Care must be taken to avoid overinsertion of the cannula and possible impingement of the cannula tip against the vessel/ventricular/atrial walls.
- The use of cold solutions and cardioplegic agents, administered externally or via the coronary circulation, may increase the stiffness of urethane plastic and contribute to ventricular perforation under conditions of cannula tip pressure during manipulation of the heart.
- Do not attempt to cross-clamp the cannula with the introducer in place.
- Carefully monitor for signs of occlusion on both the inlet and outlet cannulae.
- User assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- As with all medical devices, this device is to be used by a trained physician only.
- Initiate adequate anticoagulation therapy prior to patient cannulation.
- Only physicians trained and experienced in using percutaneous cannulation techniques (such as the Seldinger technique) and cardiopulmonary bypass should use this device.

- Ischemia may occur distal to the cannula placement site.
- Cannula tip placement should be confirmed by using fluoroscopy.
- Extreme care must be exercised to ensure that the arterial cannula and venous cannula are not confused while being connected to the perfusion circuit.
- In some cases, suture rings may not be detectable on x-ray due to the variances in anatomical densities and x-ray techniques. Modifying the x-ray angle or distance may aid in detection.
- When utilizing vacuum for augmented drainage, do not exceed -100 mmHg.
- Maximum pressure is 300 mmHg.
- Do not expose to MRI while cannula is inserted in patient.
- Due to the presence of phthalates in the product, the clinician must weigh the medical benefits of product use against the drawbacks of phthalate exposure for male children and pregnant or nursing women.

5. Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, embolic events and dislocation, infections, blood loss, and thrombus formation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

6. Complications

Possible complications include those normally associated with cardiopulmonary bypass, anticoagulation, and large bore venous and arterial catheterization.

Note: The benefits of cannulation for extracorporeal circulation must be weighed against the risk of systemic anticoagulation and subsequent propensity for hemorrhage.

7. Sterilization

All Bio-Medicus cannulae are sterilized (EtO) prior to shipment. The Lot Number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Resterilization is not recommended under any circumstances as this device is intended for single use only.

8. Instructions for Use

Note: Sterile cannula should be placed into the vascular system according to standard surgical or percutaneous techniques. Various sizes of dilators and cannulae are available in order to select the best cannula for each specific procedure. When it is determined that extracorporeal bypass is no longer required, the cannulae are carefully removed from the vessels and their placement sites are repaired according to standard surgical practice.

Note: The introducer has been designed to accommodate a 0.038-in (0.97-mm) diameter guidewire (not supplied with the device).

1. Remove the adult cannula and introducer from the shipping carton and inspect the protective, sterile package for damage.
2. Using aseptic technique, open the protective sterile package containing the cannula and introducer.
Note: Use a Bio-Medicus™ insertion kit when performing Seldinger or percutaneous cannulation technique.
3. Insert the introducer through the hemostasis cap and ensure the introducer hub is fully inserted into the hemostasis cap.
4. Assemble the cannula/introducer. Insert introducer into the lumen of the cannula and advance it until the hemostasis cap is pushed firmly onto the barb connector.
Note: Tighten the luer fitting cap on the cannula connector (if present) to ensure that the luer fitting is completely sealed.
5. Insert the cannula using a guidewire placement Seldinger technique. Determine the length of cannula to be introduced and note the insertion markings on the cannula.
Note: Some cannulae are labeled with length markings in 10 cm (3.94 in) increments for depth reference during cannula placement.
Note: Confirm the anatomical location of the guidewire before inserting the cannula. Anatomical location of guidewire is critical before cannula insertion.
Caution: Ensure that the vessel is of adequate size for cannula size selected to allow distal perfusion of the limb when the cannula is in place. Improper cannula size may be difficult to advance. The vessel must be large enough to ensure perfusion and venous return.
Caution: Do not kink guidewire during insertion, vessel dilation, and manipulation or it may be difficult to remove.
6. Position the adjustable radiopaque suture ring to the desired depth.

- 7 Pass the cannula/introducer assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.
- 8 Insert tip of the cannula/introducer into the vessel approximately 3 cm (1.18 in) to 5 cm (1.97 in). Grasp the guidewire and hold stationary **while** advancing the cannula to the desired depth. Blood leakage from the proximal drainage holes may occur during catheter insertion until these holes are fully advanced into the vessel.
Note: Advancing the cannula and introducer with a rotating motion will make cannula placement easier.
Note: When the bold band closest to the introducer hub is in alignment with the edge of the hemostasis cap, the introducer tip and cannula tip are at the same depth (Refer to Figure 1). For venous cannula the alignment of the cannula and introducer ends is indicated by the tip of the triangle-shaped symbol.
- 9 When advancing the venous cannula into the superior vena cava (SVC), verify positioning with tactile feel, indirect visualization method, or echocardiographic imaging.
Note: For venous models, graduated marks on the hub end of the introducer allow users to approximate the introducer length beyond the tip of the cannula.
- 10 For proper intra-cardiac advancement, once the side holes are intra-vascular, withdraw the introducer to the base of the triangle-shaped symbol representing the tapered tip of the introducer.
Warning: Do not withdraw to the point of the triangle shaped symbol or there is risk of SVC, innominate vein perforation or vascular cardiac trauma during intra-vascular advancement. Never force the cannula. It should advance easily.
- 11 When desired cannula placement has been achieved, remove the guidewire and retract the introducer clear of the clamp site.
Caution: To maintain hemostasis do not completely remove the introducer or hemostasis cap until the clamp is applied.
Warning: Ensure clamp is properly positioned on cannula tubing between tubing connector and wire-reinforced sections. Clamping too close to connector may damage connector end. Clamping on wire-reinforced sections may result in permanent wall distortion and/or lumen collapse.
Warning: When the bold band closest to the vascular end of the introducer is visible at the edge of the hemostasis cap, the cannula can be safely cross-clamped. Do not attempt to cross-clamp the introducer.
Note: Perfusion tubing clamps must be used. Many clamps are too weak to ensure complete clamping.
- 12 Remove the hemostasis cap, exposing the 3/8-in (0.95-cm) barbed connector.
Note: Venous caps are color-coded blue. Since 18-cm tip length models may be used as either an arterial delivery cannula or as a venous return cannula, those hemostasis caps are clear, not color-coded.
Note: If needed, an adapter is included with 23 Fr and larger venous models for connection to 1/2-in (1.27-cm) tubing.
- 13 Slowly fill the remainder of the cannula tube with sterile priming fluid, displacing any air out through the tube connector end before connecting fluid-filled tubing.
- 14 Attach the extracorporeal circuit, taking care to connect the venous cannula to the venous return line (inlet to pump) and the arterial cannula to the arterial return line (outlet from membrane oxygenator), ensuring air-free connection.
Note: Venous cannula can be identified by an arrow, in the clamp area next to the Bio-Medicus cannula logo, indicating flow out of the connector. See Figure 1.
- 15 The tubing may be further secured with cable ties at the barbed connector. Secure the device to the patient following standard medical practice.
Note: In order to avoid accidental displacement during extracorporeal circulation, suture may be wrapped around the suture ring to secure the cannula and avoid kinking the body of the device.
- 16 The retention clip may be used to secure cannula and/or tubing to the patient or drapes. The clip may be initially placed using adhesive tape. To ensure that the cannula is secure during the procedure, clamp or suture the retention clip through the holes.
Warning: Do not rely on the adhesive alone to secure the clip.
- 17 When it is determined that extracorporeal bypass is no longer required, clamp the non wire-reinforced tube section, carefully remove the cannulae from the vessels and repair their placement sites according to standard surgical practice.
Warning: This product has been qualified for use up to 6 hours. Suitability for long term use has not been determined.
- 18 After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

9. The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE ADULT CANNULA AND INTRODUCER, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

10. The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE ADULT CANNULA AND INTRODUCER, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Bio-Medicus™

Pediatric Cannula and Introducers

1. Product Description

Pediatric cannula and introducers. Sterile, nonpyrogenic, disposable single use only.

2. Indications for Use

These devices are to be used by a trained physician only.

Cannulae are used to cannulate vessels, perfuse vessels or organs, and/or connect with **ACCESSORY** extracorporeal equipment. The cannula introducer is intended to facilitate proper insertion and placement of the appropriately-sized cannula within the vessel for cardiopulmonary bypass. These products are intended for use up to 6 hours.

3. Contraindications

Alone, the cannula and introducer are not medical treatment devices. The introducer is to be used **only with the** appropriately-sized Bio-Medicus™ cannula. These devices are not intended for use except as indicated above.

4. Warnings and Precautions

- Read all warnings, precautions, and Instructions for Use carefully prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- Verify Use By date.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken.
- Do not use alcohol or alcohol-based fluids for lubrication.
- Thoroughly inspect the cannula prior to use to verify that the lumen and side holes (if applicable) are patent and that the cannula has not been damaged or kinked.
- Use caution when assembling or manipulating the introducer within the cannula, or moving the suture ring.
- Avoid false lumen cannulation.
- Never advance, position, or place the cannula without the introducer as this may kink and damage the cannula.
- Use care to protect the side hole areas from kinking as they are passed into the vessel.
- Manipulate cannula in such a manner as to prevent kinks or any restrictions that may alter flow.
- Manipulation of the heart with cannula in place may cause flow restrictions.
- Avoid kinking when attaching tubing and when anchoring or untethering the cannula body.
- Caution should be used when attaching the cannula to other devices in the extracorporeal circuit so as not to damage the cannula or the device being connected.
- Clamping too close to connector may damage connector end.
- Do not clamp wire-reinforced sections. Clamping on reinforced sections may result in permanent wall distortion and/or lumen collapse.
- Since the cannula wall is opaque due to the presence of the spring, be sure to completely remove air from the cannula during the priming procedure and before clamping.
- Care must be taken to avoid overinsertion of the cannula and possible impingement of the cannula tip against the vessel/ventricular/atrial walls.
- The use of cold solutions and cardioplegic agents, administered externally or via the coronary circulation may increase the stiffness of urethane plastic and contribute to ventricular perforation under conditions of cannula tip pressure during manipulation of the heart.
- Do not attempt to cross-clamp the cannula with the introducer in place.
- Carefully monitor for signs of occlusion on both the inlet and outlet cannulae.
- User assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- As with all medical devices, this device is to be used by a trained physician only.
- Initiate adequate anticoagulation therapy prior to patient cannulation.
- Only physicians trained and experienced in using percutaneous cannulation techniques (such as the Seldinger technique) and cardiopulmonary bypass should use this device.
- Ischemia may occur distal to the cannula placement site.
- Cannula tip placement should be confirmed by using fluoroscopy.

Instructions for Use English 7

master.20140730, alldocMod.20140911, alldocMod2.20140911, alldocStyle.20140911, **alldocStyle_graph.20150220**
master_pageSet.20140911, languageLookup.20140514, table.20140620, characterJoin.20110211, 3of9.20110211, 3of9_svg.
20110211

- Extreme care must be exercised to ensure that the arterial cannula and venous cannula are not confused while being connected to the perfusion circuit.
- In some cases, suture rings may not be detectable on x-ray due to the variances in anatomical densities and x-ray techniques. Modifying the x-ray angle or distance may aid in detection.
- When utilizing vacuum for augmented drainage, do not exceed -100 mmHg.
- Maximum pressure is 300 mmHg.
- Do not expose to MRI while cannula is inserted in patient.

5. Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to: embolic events and dislocation, infections, blood loss and thrombus formation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

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7. Sterilization

All Bio-Medicus cannulae are sterilized (EIO) prior to shipment. The Lot Number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Resterilization is not recommended under any circumstances as this device is intended for single use only.

8. Instructions for Use

Note: Sterile cannula should be placed into the vascular system according to standard surgical or percutaneous technique. Various sizes of dilators and cannulae are available in order to select the best cannula for each specific procedure. When it is determined that extracorporeal bypass is no longer required, the cannulae are carefully removed from the vessels and their placement sites are repaired according to standard surgical practice.

Note: The introducer has been designed to accommodate a 0.025-in (0.635-mm) diameter guidewire (not supplied with the device).

1. Remove the cannula and introducers from the shipping carton and inspect the protective, sterile package for damage.
2. Using aseptic technique, open the protective sterile package containing the cannula and the introducers.
3. Select either the standard solid core introducer for traditional placement method or the hollow core guidewire-capable introducer for Seldinger technique placement method.
Note: Use a pediatric-sized Bio-Medicus™ insertion kit when performing Seldinger or percutaneous cannulation.
4. Insert the introducer through the hemostasis cap and ensure the introducer hub is fully inserted into the hemostasis cap. Insert the introducer into the lumen of the cannula through the connector (Figure 1).
5. Firmly place hemostasis cap onto the connector of cannula.
6. Insert the cannula using a guidewire placement Seldinger technique. Determine the length of cannula to be introduced and note the insertion markings on the cannula.

Note: Some cannulae are labeled with length markings in 10-cm (3.94-in) increments for depth reference during cannula placement.

Note: Confirm the anatomical location of the guidewire before inserting the cannula. Anatomical location of guidewire is critical before cannula insertion.

Caution: Ensure that the vessel is of adequate size for cannula size selected to allow distal perfusion of the limb when the cannula is in place. Improper cannula size may be difficult to advance. The vessel must be large enough to ensure perfusion and venous return.

Caution: Do not kink guidewire during insertion, vessel dilation and manipulation or it may be difficult to remove.

7. Position the adjustable radiopaque suture ring to the desired depth.
8. Pass the cannula/introducer assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.
9. Insert tip of the cannula/introducer into the vessel. Grasp the guidewire and hold stationary while advancing the cannula to the desired depth.

10. Blood leakage from the proximal drainage holes of the venous cannulae models may occur during catheter insertion until these holes are intravascular.
- Note:** Advancing the cannula and introducer with a rotating motion will make cannula placement easier.
- Note:** When the bold band closest to the introducer hub is in alignment with the edge of the hemostasis cap, the cannula tip and the introducer tip are at the same depth (Refer to Figure 1).
11. When the desired cannula placement has been achieved, remove the introducer (and guidewire if applicable).
- Warning:** To maintain hemostasis, do not completely remove the introducer or hemostasis cap until the clamp is applied.
- Warning:** Ensure clamp is properly positioned on cannula tubing between tubing connector and wire-reinforced sections. Clamping too close to the connector may damage the connector end. Clamping on wire-reinforced sections may result in permanent wall distortion and/or lumen collapse.
- Warning:** When the bold band closest to the vascular end of the introducer is visible at the edge of the hemostasis cap, the cannula can be safely cross-clamped. Do not attempt to cross-clamp the introducer.
- Note:** Perfusion tubing clamps must be used. Many clamps are too weak to ensure complete clamping.
12. Remove the hemostasis cap exposing the 1/4-in (0.635-cm) barbed connector.
- Note:** Arterial hemostasis caps are color-coded red. Venous caps are color-coded blue.
13. Slowly fill the remainder of the cannula tube with sterile priming fluid, displacing any air out through the tube connector end before connecting the fluid-filled tubing.
14. Attach the extracorporeal circuit, taking care to connect the venous cannula to the venous return line (inlet to pump) and the arterial cannula to the arterial return line (outlet from membrane oxygenator), ensuring air-free connection. The tubing may be further secured with cable ties at the barbed connector.
- Note:** Venous cannula can be identified by an arrow, in the clamp area next to the Bio-Medicus cannula logo, indicating flow out of the connector. See Figure 1.
15. Secure the device to the patient following standard medical practice. In order to avoid accidental displacement during extracorporeal circulation, suture may be wrapped around the suture ring to secure the cannula and avoid kinking the body of the device.
16. The retention clip may be used to secure cannula and/or tubing to the patient or drapes. The clip may be initially placed using the adhesive tape. To ensure that the cannula is secure during the procedure, clamp or suture the retention clip through the holes.
- Warning:** Do not rely on the adhesive alone to secure the clip.
17. When it is determined that extracorporeal bypass is no longer required, clamp the non wire-reinforced tube section, carefully remove the cannulae from the vessels and repair their placement sites according to standard surgical practice.
- Warning:** This product has been qualified for use up to 6 hours. Suitability for long-term use has not been determined.
18. After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

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Disclaimer of Warranty

ALTHOUGH THE PEDIATRIC CANNULA AND INTRODUCERS, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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Disclaimer of Warranty

ALTHOUGH THE PEDIATRIC CANNULA AND INTRODUCERS, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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9. Cannula femoral pediatric arterial Bio-Medicus, effective length 100 mm, external diameter 2.7 mm, composed of:

- Cannula femoral pediatric arterial Bio-Medicus, effective length 100 mm, outer diameter 2.7 mm,
- The introducer with a solid conical tip,
- Introducer,
- Hemostatic cap,
- Clamp for fixing the cannula,
- Instruction manual.

10. Cannula femoral pediatric arterial Bio-Medicus, effective length 105 mm, outer diameter 3.3 mm, composed of:

- Cannula femoral pediatric arterial Bio-Medicus, effective length 105 mm, outer diameter 3.3 mm,
- The introducer with a solid conical tip,
- Introducer,
- Hemostatic cap,

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- Clamp for fixing the cannula,
- Instruction manual.

11. Cannula femoral pediatric arterial Bio-Medicus, effective length 110 mm, outer diameter 4.0 mm, composed of:

- Cannula femoral pediatric arterial Bio-Medicus, effective length 110 mm, outer diameter 4.0 mm,
- The introducer with a solid conical tip,
- Introducer,
- Hemostatic cap,
- Clamp for fixing the cannula,
- Instruction manual.

12. Cannula femoral pediatric arterial Bio-Medicus, effective length 115 mm, external diameter 4.7 mm, composed of:

- Cannula femoral pediatric arterial Bio-Medicus, effective length 115 mm, outer diameter 4.7 mm,
- The introducer with a solid conical tip,
- Introducer,
- Hemostatic cap,
- Clamp for fixing the cannula,
- Instruction manual.

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	connection, sizes 4.7 mm, 5.7 mm, 7.0 mm	environment of the body, blood Collar, Clamp site, Vented connector, Luer Cap, Solvent: Device does not contact with the patient (Medical staff works in gloves).	manufactured by Lubrizol; Depth markings: Ink, Permanent, 0.1mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd.; Ink, Permanent, 0.5 mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd; Collar: PCV 3006/1-85 CLEAR 0003, manufactured by AlphaGary Corp.; Colorant 8538-E1 red, manufactured by AMERICHEM, Inc; Clamp site: PVC 3006/1-85 CLEAR 0003, manufactured by AlphaGary Corp. Vented connector: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. Luer Cap: Polypropylene homopolymer Pro-fax PD626, manufactured by Basell USA Inc.; Polypropylene Color Concentrate CH-01010 SP, manufactured by FERRO Corp. -Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
	Femoral arterial cannula DLP, sizes 5.7 mm, 7.0 mm	Cannula tube, Depth markings: Short-term contact with the internal environment of the body, blood Collar, Clamp site, Connector, Solvent: Device does not contact with the patient (Medical staff works in gloves).	Cannula tube: Liquid silicone 360 Medical Fluid 100 CST, manufactured by DOW CORNING; Polyether aliphatic Tecoflex EG-100A-B40, manufactured by Lubrizol; Depth markings: Ink, Permanent, 0.1mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd.; Ink, Permanent, 0.5 mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd; Collar: PCV 3006/1-85 CLEAR 0003, manufactured by AlphaGary Corp.; Colorant 8538-E1 red, manufactured by AMERICHEM, Inc ; Clamp site: PVC 3006/1-85 CLEAR 0003, manufactured by AlphaGary Corp. Connector: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. -Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich

	Femoral venous cannula DLP, sizes 5.7 mm, 7.0 mm	Cannula tube, Depth markings: Short-term contact with the internal environment of the body, blood Collar, Clamp site, Connector, Solvent: Device does not contact with the patient (Medical staff works in gloves).	Cannula tube: Liquid silicone 360 Medical Fluid 100 CST, manufactured by DOW CORNING; Polyether aliphatic Tecoflex EG-100A-B40, manufactured by Lubrizol; Depth markings: Ink, Permanent, 0.1mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd.; Ink, Permanent, 0.5 mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd; Collar: PCV 3006/1-85 CLEAR 0003, manufactured by AlphaGary Corp.; Colorant 8533-E1 dark blue, manufactured by AMERICHEM, Inc; Clamp site: PVC 3006/1-85 CLEAR 0003, manufactured by AlphaGary Corp. Non-vented connector: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. -Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
	Polyurethane dilator	Polyurethane dilator tube, Depth markings: Short-term contact with the internal environment of the body, blood	Polyurethane dilator tube: Liquid silicone, 360 Medical Fluid, 100 CST, manufactured by DOW CORNING; Polyether aliphatic, Tecoflex EG-93A-B40, manufactured by Lubrizol Depth markings: Ink, Permanent, 0.1mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd. ; Ink, Permanent, 0.5 mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd



	Femoral venous cannula DLP, size 9.3 mm	Cannula tube, Depth markings: Short-term contact with the internal environment of the body, blood Collar, Clamp site, Transition fitting, Solvent: Device does not contact with the patient (Medical staff works in gloves).	Cannula tube: Liquid silicone 360 Medical Fluid 100 CST, manufactured by DOW CORNING; Polyether aliphatic Tecoflex EG-100A-B40, manufactured by Lubrizol; Depth markings: Ink, Permanent, 0.1mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd.; Ink, Permanent, 0.5 mm FinePoint System Pen «PFS» (2-Methoxy ethanol - 23,7%, Methanol – 63,7%, Keton Resin – 7,3%, C.I. Solvent Black 27 – 5,3 %), manufactured by ITOYA of America, Ltd; Collar: PCV 3006/1-85 CLEAR 0003, manufactured by AlphaGary Corp.; Colorant 8533-E1 dark blue, manufactured by AMERICHEM, Inc.; Clamp site: Polypropylene Pro-fax 6823, manufactured by Basell Canada Inc. - Transition fitting: PVC 3006-77 CLEAR 0003, manufactured by AlphaGary Corp. -Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
	Sheath for Femoral venous cannula DLP, size 9.3 mm	Short-term contact with the internal environment of the body, blood	Liquid silicone, 360 Medical Fluid, 100 CST, manufactured by DOW CORNING; Polyether aliphatic, Tecoflex EG-100A-B40, manufactured by Lubrizol


	Guidewire diameter 0.1 cm for Femoral venous cannula DLP, size 9.3 mm	Guidewire body: Short-term contact with the internal environment of the body, blood Dispenser, Track clip: Device does not contact with the patient (Medical staff works in gloves).	Guidewire body: Stainless Steel Type 304 (iron – 70%, chromium - 18,5%, nickel - 9%, manganese - 1%, cobalt - 0,2%, others – 1,3%); Copolymer Polypropylene Titanpro SM546, manufactured by Titan Petchem -Dispenser: Polyethylene Marlex, manufactured by Phillips Chemical Company; - Track clip: PVC 3006-77 CLEAR 0003, manufactured by AlphaGary Corp.; Colorant Brilliant Ultra Blue VF156, manufactured by Nazdar Ltd.; Colorant Permanent Blue VF159, manufactured by Nazdar Ltd.
	Femoral arterial cannula Bio-Medicus, sizes 5.0 mm, 5.7 mm, 6.3 mm, 7.0 mm, 7.7 mm, 8.3 mm	Cannula tube, Cannula body, Depth markings, Adjustable radiopaque suture ring: Short-term contact with the internal environment of the body, blood Connector, Luer Cap, Solvent: Device does not contact with the patient (Medical staff works in gloves).	-Cannula tube: Polyurethane - Pellethane 5863-91A, manufactured by Lubrizol Corp.; - Cannula body: Stainless Steel Type 304 (iron – 70%, chromium - 18,5%, nickel - 9%, manganese - 1%, cobalt - 0,2%, others – 1,3%); - Depth markings: L-219 Ink (Methyl Ethyl Ketone – 80%, Isopropyl Alcohol – 10%, Cellulose Nitrate – 10%) - Connector: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. - Adjustable radiopaque suture ring: Silicone - MED-4860 Part A, manufactured by NuSil Technology LLC; Silicone - MED-4860 Part B, manufactured by NuSil Technology LLC Radiopaque marking in adjustable radiopaque suture ring: Barium Sulfate - 100%; - Luer Cap: Polypropylene Pro-fax 6524, manufactured by Basell Canada Inc. -Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
	Introducer for Femoral arterial	Introducer tube, Depth markings: Short-term contact with the internal environment of the body, blood	- Introducer tube: PVC 9066M-RAD14, manufactured by Colorite; - Depth markings: L-219 Ink (Methyl Ethyl Ketone – 80%, Isopropyl

cannula Bio-Medicus, sizes 5.0 mm, 5.7 mm, 6.3 mm, 7.0 mm, 7.7 mm, 8.3 mm	Hub, Solvent: Device does not contact with the patient (Medical staff works in gloves).	Alcohol – 10%, Cellulose Nitrate – 10%) - Hub: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. - Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
Hemostatic Cap for Femoral arterial cannula Bio-Medicus, sizes 5.0 mm, 5.7 mm, 6.3 mm, 7.0 mm, 7.7 mm, 8.3 mm	Device does not contact with the patient (Medical staff works in gloves).	Silicone LIM 6050-D2, manufactured by Momentive Performance Materials LLC
Cannula retention clip for Femoral arterial cannula Bio-Medicus, sizes 5.0 mm, 5.7 mm, 6.3 mm, 7.0 mm, 7.7 mm, 8.3 mm	Device does not contact with the patient, it is the primary packaging (designed to fix the product during storage and transportation). Medical staff works in gloves.	Polyethylene Marlex, manufactured by Phillips Chemical Company
Femoral venous cannula Bio-Medicus, sizes 5.0 mm, 5.7 mm, 6.3 mm, 7.0 mm, 7.7 mm, 8.3 mm, 9.0 mm, 9.7 mm	Cannula tube, Canula body, Depth markings, Adjustable radiopaque suture ring: Short-term contact with the internal environment of the body, blood Connector, Solvent: Device does not contact with the patient (Medical staff works in gloves).	- Cannula tube: Polyurethane - Pellethane 5863-91A, manufactured by Lubrizol Corp.; - Cannula body: Stainless Steel Type 304 (iron – 70%, chromium - 18,5%, nickel - 9%, manganese - 1%, cobalt - 0,2%, others – 1,3%) ; - Depth markings: L-219 Ink (Methyl Ethyl Ketone – 80%, Isopropyl Alcohol – 10%, Cellulose Nitrate – 10%) - Connector: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. - Adjustable radiopaque suture ring: Silicone - MED-4860 Part A, manufactured by NuSil Technology LLC; Silicone - MED-4860 Part B, manufactured by NuSil Technology LLC; Radiopaque marking in adjustable radiopaque suture ring: Barium Sulfate - 100% ; - Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl

			Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
	Introducer for Femoral venous cannula Bio-Medicus, sizes 5.0 mm, 5.7 mm, 6.3 mm, 7.0 mm, 7.7 mm, 8.3 mm, 9.0 mm, 9.7 mm	Introducer tube, Depth markings: Short-term contact with the internal environment of the body, blood Hub, Solvent: Device does not contact with the patient (Medical staff works in gloves).	- Introducer tube: PVC 9066M-RAD14, manufactured by Colorite; - Depth markings: L-219 Ink (Methyl Ethyl Ketone – 80%, Isopropyl Alcohol – 10%, Cellulose Nitrate – 10%) - Hub: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. - Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
	Hemostatic Cap for Femoral venous cannula Bio-Medicus, sizes 5.0 mm, 5.7 mm, 6.3 mm, 7.0 mm, 7.7 mm, 8.3 mm, 9.0 mm, 9.7 mm	Device does not contact with the patient (Medical staff works in gloves).	Silicone LIM 6050-D2, manufactured by Momenive Performance Materials LLC; Colorant 8533-E1 dark blue, manufactured by AMERICHEM, Inc
	Cannula retention clip for Femoral venous cannula Bio-Medicus, sizes 5.0 mm, 5.7 mm, 6.3 mm, 7.0 mm, 7.7 mm, 8.3 mm, 9.0 mm, 9.7 mm	The product does not contact with the patient, (designed to fix the product during storage and transportation). Medical staff works in gloves.	Polyethylene Marlex, manufactured by Phillips Chemical Company

Pediatric femoral arterial cannula Bio-Medicus, sizes 2.7 mm, 3.3 mm, 4.0 mm, 4.7 mm	Cannula tube, Cannula body, Depth markings, Adjustable radiopaque suture ring: Short-term contact with the internal environment of the body, blood Connector, Solvent: Device does not contact with the patient (Medical staff works in gloves).	- Cannula tube: Polyurethane - Pellethane 5863-91A, manufactured by Lubrizol Corp.); - Cannula body: Stainless Steel Type 304 (iron – 70%, chromium - 18,5%, nickel - 9%, manganese - 1%, cobalt - 0,2%, others – 1,3%) ; - Depth markings: L-219 Ink (Methyl Ethyl Ketone – 80%, Isopropyl Alcohol – 10%, Cellulose Nitrate – 10%) - Connector: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. - Adjustable radiopaque suture ring: Silicone - MED-4860 Part A, manufactured by NuSil Technology LLC; Silicone - MED-4860 Part B, manufactured by NuSil Technology LLC Radiopaque marking in adjustable radiopaque suture ring: Barium Sulfate - 100%; -Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
Introducer with solid tapered tip for Pediatric femoral arterial cannula Bio-Medicus, sizes 2.7 mm, 3.3 mm, 4.0 mm, 4.7 mm; Pediatric femoral venous cannula Bio-Medicus, sizes 2.7 mm, 3.3 mm, 4.0 mm, 4.7 mm	Introducer tube, Depth markings: Short-term contact with the internal environment of the body, blood Hub, Solvent: Device does not contact with the patient (Medical staff works in gloves).	- Introducer tube: PVC 9066M-RAD14, manufactured by Colorite; - Depth markings: L-219 Ink (Methyl Ethyl Ketone – 80%, Isopropyl Alcohol – 10%, Cellulose Nitrate – 10%) - Hub: PVC 2212C2B/H-118 Clear 0003, manufactured by AlphaGary Corp. -Solvent: PVC Homopolymer Suspension Resin, manufactured by Oxy Vinyl Tetrahydrofuran, manufactured by Sigma-Aldrich Cyclohexanone, manufactured by Sigma-Aldrich
Hemostatic Cap for Pediatric femoral arterial cannula Bio-Medicus, sizes 2.7 mm, 3.3 mm, 4.0 mm, 4.7 mm; Pediatric femoral	Device does not contact with the patient (Medical staff works in gloves).	For Pediatric femoral arterial cannula Bio-Medicus: Silicone LIM 6050-D2, manufactured by Momentive Performance Materials LLC; Colorant 8538-E1 red, manufactured by AMERICHEM, Inc For Pediatric femoral venous cannula Bio-Medicus: Silicone LIM 6050-D2, manufactured by Momentive Performance Materials LLC ; Colorant 8533-E1 dark blue, manufactured by

	venous cannula Bio-Medicus, sizes 2.7 mm, 3.3 mm, 4.0 mm, 4.7 mm		AMERICHEM, Inc
	Cannula retention clip for Pediatric femoral arterial cannula Bio-Medicus, sizes 2.7 mm, 3.3 mm, 4.0 mm, 4.7 mm; Pediatric femoral venous cannula Bio-Medicus, sizes 2.7 mm, 3.3 mm, 4.0 mm, 4.7 mm	The product does not contact with the patient, (designed to fix the product during storage and transportation). Medical staff works in gloves.	Polyethylene Marlex, manufactured by Phillips Chemical Company
10. Данные о маркировке медицинского изделия и его упаковке/ Labelling and Packaging of medical device	Для канюль бедренных артериальных DLP (первичная упаковка – слева, вторичная упаковка - справа) Femoral arterial cannula DLP (Primary package left, secondary package right):  		

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	Lot Number
	Quantity
	Open Here
	Pull Tab to Open
	Caution, Consult Accompanying Documents
	Consult Instructions for Use
	For US Audiences Only
	Keep Dry
	Upper Limit of Temperature
	Fragile, Handle with Care
	Catalog Number
	Manufacturer
	Date of Manufacture
	Contains Phthalate DEHP
	Do Not Use if Package is Damaged
	Authorized Representative in the European Community
	Manufactured In

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Information in Russian sticker The final packaging for the Russian market is accompanied with additional sticker. Devices supplied in packaging marked with the following possible information: • Name of medical device in Russian language; • Name of the manufacturer; • Country of origin; • Code UPN (unique product number); • Code CFN (customer facing number); • Date and number of Registration Certificate; • Sign of conformity in Russia; • Warning and informing symbols and inscriptions (if applicable). Sample of the sticker	<table><tr><td rowspan="5"> Non-toxic DEHP EC REP</td><td colspan="2">Name of medical device Cannulas of femoral arterial and venous DLP and Bio-Medicus, variants of execution</td></tr><tr><td>Manufacturer: Medtronic Inc.</td><td>Code UPN</td></tr><tr><td>Country of origin: USA</td><td>Code CFN</td></tr><tr><td colspan="2">Manufacturer's authorized representative: LLC «Medtronic»</td></tr><tr><td colspan="2">Regestretion Certificate № ____ от ____</td></tr></table>	 Non-toxic DEHP EC REP	Name of medical device Cannulas of femoral arterial and venous DLP and Bio-Medicus, variants of execution		Manufacturer: Medtronic Inc.	Code UPN	Country of origin: USA	Code CFN	Manufacturer's authorized representative: LLC «Medtronic»		Regestretion Certificate № ____ от ____	
 Non-toxic DEHP EC REP	Name of medical device Cannulas of femoral arterial and venous DLP and Bio-Medicus, variants of execution											
	Manufacturer: Medtronic Inc.		Code UPN									
	Country of origin: USA		Code CFN									
	Manufacturer's authorized representative: LLC «Medtronic»											
	Regestretion Certificate № ____ от ____											

	Date of	See IFU
--	---------	---------

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© 2003, 2011 Medtronic