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VP OPERATIONS DPS DACH & AMERICAS OPERATIONS

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Senior Quality Operations manager

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Johnson & Johnson CAMPUS NEUCHÂTEL

LEGALIZATION

The undersigned, Bernard Lamboley, notary public of the Canton of Neuchâtel, Switzerland, with office at La Chaux-de-Fonds, certifies that the signatures set on this Approval are the true and genuine signatures of **Mr. Thomas Alwin Karl SCHÄDLER**, from Germany, in Feldbrunnen-St. Niklaus, Switzerland, acting as VP Operations DPS DACH & Americas and of **Mr. Nicolas Pierre HAINARD**, from Val-de-Travers, in Cudrefin, Switzerland, acting as Senior Quality Operations Manager who are competent to act on behalf of Medos International SARL, in Le Locle, Switzerland. -----
La Chaux-de-Fonds, this twenty-seventh day of April two thousand twenty (April, 27, 2020) -
RGV 62 No 131

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
Instruction for use

The ENVOY DA Guiding Catheter with introducer (for use in the Russian Federation)

Medos International SÀRL
Chemin-Blanc 38, 2400
Le Locle, Switzerland
41 32 933 83 00

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA:
Johnson & Johnson LLC,
121614, Moscow, Krylatskaya str., 17/2.
Tel. +7 (495) 580-7777, Fax: +7 (495) 580-7878

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IMPORTANT INFORMATION

Please Read Before Use

The ENVOY DA Guiding Catheter with introducer**Medical Device name, models or variations of MD**

The ENVOY DA Guiding Catheter with introducer

Product Code	Description
67125805D	ENVOY DA, 6F, 105cm, MPD
67125895D	ENVOY DA, 6F, 95 cm, MPD
67126005D	ENVOY DA, 6F, 105cm, Straight
67126095D	ENVOY DA, 6F, 95 cm, Straight
67125605D	ENVOY DA, 6F, 105cm, MPC
67126005DB	ENVOY DA XB, 6F, 105cm, Straight
67125805DB	ENVOY DA XB, 6F, 105cm, MPD
67125605DB	ENVOY DA XB, 6F, 105cm, MPC

Description

The **ENVOY DA** Guiding Catheter is a variable stiffness, braided catheter with a large non-tapered lumen that facilitates the intravascular passage of interventional devices. The distal segment is flexible for navigation into distal vasculature. The catheter has an outer hydrophilic coating that reduces friction during manipulation in the vessel. The lubricious PTFE lined inner lumen is designed to facilitate movement of the guide wires and other devices. A luer fitting located on the end of the catheter hub can be used to attach accessories. The distal section of the catheter is radiopaque to aid visualization under fluoroscopy, and the distal tip is clearly distinguished by a radiopaque marker. The catheter is available with preshaped tips to facilitate positioning. A peel away introducer is included to facilitate insertion into the sheath.

IndicationsThe **ENVOY DA** Guiding Catheter is intended for use in the peripheral,

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Indications

The ENVOY DA Guiding Catheter is intended for use in the peripheral, coronary, and neuro vasculature for the intravascular introduction of interventional/diagnostic devices.

Contraindications

None known.

Warnings

Discard the guiding catheter after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused. Do not use with Ethiodol™ or Lipiodol™ contrast media, or other such contrast media which incorporates the components of these agents.

Precautions

Inspect the sterile package carefully. Do not use if:

- The package or seal appears damaged,
- Contents appear damaged, or
- The expiry date has passed.
- Store in a cool, dark, dry place.
- Exposure to temperatures above 54°C (130°F) may damage the catheter.
- Do not expose to organic solvents.
- Inspect the guiding catheter before use to verify that its size, shape, and condition are suitable for the specific procedure. If strong resistance is met during manipulation, discontinue the procedure and determine the cause for the resistance before proceeding. If the cause of the resistance cannot be determined, withdraw the catheter.
- Extreme care must be taken to avoid damage to the vasculature through which the catheter passes. The guiding catheter may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.

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- Torqueing the catheter excessively while kinked may cause damage which could result in possible separation along the catheter shaft. Should the guiding catheter shaft become severely kinked, withdraw the entire system (guiding catheter, guidewire, and catheter sheath introducer).

- Advancement, manipulation, and withdrawal of the guiding catheter should always be performed under fluoroscopic guidance.

- Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

Adverse Events

Possible complications include:

- acute occlusion
- air embolism
- death
- distal embolization
- emboli
- false aneurysm formation
- hematoma or hemorrhage at puncture site
- infection
- intracranial hemorrhage
- ischemia
- neurological deficits including stroke
- vessel spasm, thrombosis, dissection, or perforation

Recommended Procedure

1. Gently remove the ENVOY DA Guiding Catheter from its packaging.
2. Inspect the catheter upon removal from packaging to verify that it is undamaged.

WARNING: Do not use a catheter that has been damaged in any way. If damage is detected, replace with another ENVOY DA Guiding Catheter that is not damaged.

3. Prior to use, flush the lumen of the Guiding Catheter with a heparinized saline solution. Immerse 10 cm of the distal end of the catheter in the heparinized saline to wet the hydrophilic coating.
4. Use the supplied peel away introducer to aid in the insertion of the catheter

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through the catheter sheath introducer.

5. Introduce the guiding catheter into the vasculature through the catheter sheath introducer, and over a .038" (0.97 mm) or .035" (0.89 mm) guidewire using a percutaneous entry technique of choice.

6. Set up a continuous heparinized saline flush through the sidearm of a rotating hemostatic valve attached to the Luer lock fitting of the guiding catheter hub.

Note: It is recommended that continuous heparinized saline flush be maintained between the ENVOY DA Guiding Catheter and any intraluminal device passed coaxially through it.

7. Under fluoroscopic guidance, advance the guiding catheter over the guidewire until the desired position is attained.

8. Remove the guidewire prior to the introduction of other intravascular devices or infusion of contrast agents.

How Supplied

The ENVOY DA Guiding Catheter is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

The manufacturer is not responsible for any product subjected to re-sterilization, and also does not exchange any product that has been opened but has not been used. As long as the inner unit is not opened or damaged, the product is sterile and nonpyrogenic.

Warranty

The manufacturer does not assume any additional warranty obligations, except those that are legally approved.

The manufacturer warrants that this medical device is free from defects

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in both materials and workmanship. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.

Disposal

Used products, as well as packaging, must be disposed of in accordance with the standards provided for the disposal of potentially polluting medical waste.

In the Russian Federation, guiding catheters are to be disposed of at the end of their service life in accordance with the requirements of SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for the management of medical waste."

The hazard class of medical waste is based on the characteristics of the morphological composition: class B (epidemiologically hazardous waste). In the case of expiration date or damage to the packaging, the products are also subject to disposal.

Conditions of transportation and storage

The medical device is transported by all means of transport in closed vehicles in accordance with the rules for the carriage of goods in force on this type of transport.

Store in a cool, dark, dry place. Do not expose to temperatures in excess of 131 ° F (55 ° C).

Transportation conditions:

Designation	Range
Temperature	8°C-30°C
Relative humidity	No more than 80%

Storage conditions:




Designation	Range
Temperature	8°C-30°C
Relative humidity	No more than 80%

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Interpretation of symbols used on the package

REF	Catalogue Number
 YYYY-MM-DD	Manufacturing date (YYYY- MM-DD)
 YYYY-MM-DD	Use by (YYYY-MM-DD)
STERILE EO	Sterilized using ethylene oxide
	Do not reuse
	Caution! Consult instructions for use
RADIOPAQUE	Radiopaque
QTY	Quantity
Made in	Made in
Rx Only	Federal Law (United States of America) restricts this Device to sale by, or on the order of a physician.
DO NOT AUTOCLAVE	Do not autoclave

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6. LIST OF STANDARDS

Standard/ Guidance/ Directive	Description
Quality System	
BS EN ISO 13485: 2012	Quality Management Systems – Medical devices – System requirements for regulatory purposes
BS EN ISO 14971: 2012	Medical Devices – Application of risk management to medical devices
European Council Directive 93/42/EEC of 14 th of June 1993 Amended 2007, Annex II	Medical Device Directive
Packaging	
BS EN ISO 11607- 1: 2009	Packaging for Terminally Sterilized Medical Devices – Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems
Sterilization	
BS EN ISO 11135-1: 2007	Sterilization of Healthcare products Ethylene Oxide: Part 1 Requirements for development, validation and routine control of a sterilization process for medical devices
BS EN ISO 10993-7: 2008	Biological Evaluation of medical devices: Ethylene oxide sterilization residuals – Part 7
Product Safety/Performance	
BS EN ISO 10555-1: 2009*	Sterile, single use intravascular catheters Part 1: General requirements; Sterile
ISO 594-1: 1986 (E)	Conical fitting with a 6% (Luer) taper for syringes, needles and certain other medical equipment: Part 1 – General Requirements

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ISO 594-2: 1998 (E)	Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment: Part 2 – Local fittings
AAMI / ANSI HE75: 2009	Human Factors Engineering – Design of Medical Devices
BS EN 62366: 2008	Medical Devices. Application of usability engineering to Medical Devices
Biocompatibility	
BS EN ISO 10993-1: 2009	Biological evaluation of medical devices: Evaluation & Testing – Part 1
BS EN ISO 10993-5: 2009	Biological evaluation of medical devices Part 5: Tests for in vitro cytotoxicity
BS EN ISO 10993-4: 2009	Biological evaluation of medical devices Part 4: Selection for tests for interactions with blood
ISO 10993-10: 2010	Biological evaluation of medical devices Part 10: Tests for irritation and skin sensitization
ISO 10993-11: 2009	Biological evaluation of medical devices Part 11: Test for systemic toxicity
ISO 10993-3: 2009	Biological evaluation of medical devices Part 3: Test for genotoxicity, carcinogenicity, and reproductive toxicity
Labeling	
BS EN ISO 15223-1: 2012	Medical Device – Symbols to be used with medical device labels, labeling and information to be supplied
BS EN 1041: 2008	Information supplied by the manufacturer with medical devices
Clinical	
MEDDEV 2.7.1 December 2009	Evaluation of Clinical Data: A Guide for Manufacturers and Notified Bodies
MEDDEV 2.12-2, January 2012	Guidelines on Post Market Clinical Follow-Up

«...gap analysis was completed and it was determined that testing that was conducted meets or exceeds the requirements defined in BS EN ISO 10555-1: 2013 Sterile, single use intravascular catheters Part 1: General requirements; Sterile.

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A gap analysis was completed and it was determined that testing that was conducted meets or exceeds the requirements defined in BS EN ISO 10555-1: 2013 Sterile, single use intravascular catheters Part 1: General requirements; Sterile.

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Description of sterilization method

The ENVOY DA Guiding Catheters are disposable and are supplied in a sterile package. Packaged products are sterilized with ethylene oxide (EO) in accordance with DIN EN ISO 11135 BS EN ISO 11135-1: 2007. The minimum guaranteed sterility level (SAL) is 10^{-6}

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Shelf life of the ENVOY DA Guiding Catheters – 3 years, after sterilization.

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The guaranteed shelf life of the product is 3 years from the date of manufacture (sterilization).

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